

Mobilität wird einfach!

Im Hofstätt 11
D- 72359 Dotternhausen
Phone +49 (0)7427 9480-0
Fax +49 (0)7427 9480-7025
E-Mail info@proactiv-gmbh.de
www.proactiv-gmbh.com



Quotation date: _____

Order date: _____

according to GTC PRO ACTIV GmbH

Medical suppliers/company address:

Cust. no.: _____
Name: _____
Street: _____
Postcode, city: _____
Phone: _____
Fax: _____
E-Mail: _____
Contact: _____

Invoicing address, other (see below)

Cust. no.: _____
Name: _____
Street: _____
Postcode, city: _____

Shipping address, other (see below)

Cust. no.: _____
Name: _____
Street: _____
Postcode, city: _____

Customer name/commission:

Name: _____
Street: _____
Postcode, city: _____
Phone: _____
Fax: _____
E-Mail: _____
DoB: _____
Paying organisation: _____
Indication: _____



FREAK youth bike

Maximum user weight limit 100 kg

- girder chassis with reinforcement brace
- drive unit, can be decoupled via the eccentric quick-lock system
- drive unit natural anodised with 20" drive wheel, optional 26"
- A-shaped handles with integrated switching fittings
- Shimano V-brake rim brake
- infinitely adjustable backrest
- wheel camber 6°, rear wheel size 20"
- safety flag

9113000013 € 3,390.00

Services:

Delivery, instruction and handing over at the customer by PRO ACTIV

9990000003 € 425.00
plus shipping costs

Body size

Body height: _____ cm

Weight: _____ kg

Frame colour – One-colour design

Information on current colour designs and example images under:
www.proactiv-gmbh.de/farben.php

White RAL 9016

- gloss 9000100170-023
 matt 9000100170-123 **€ 33.30**

Brilliant silver

- gloss 9000100170-011

Anthracite RAL 7016

- gloss 9000100170-001
 matt 9000100170-101 **€ 33.30**

Black RAL 9005

- gloss 9000100170-003
 matt 9000100170-019 **€ 33.30**

Deep black, silver effect

- gloss 9000100170-007

Pure orange RAL 2004

- gloss 9000100170-022

Fire Red

- gloss 9000100170-018
 matt 9000100170-118 **€ 33.30**

Pink magenta RAL 4010

- gloss 9000100170-037
 matt 9000100170-137 **€ 33.30**

Sky blue

- gloss 9000100170-017

Pastel blue RAL 5024

- gloss 9000100170-036
 matt 9000100170-136 **€ 33.30**

Marine blue

- gloss 9000100170-045
 matt 9000100170-145 **€ 33.30**

Apple green

- gloss 9000100170-014
 matt 9000100170-114 **€ 33.30**

Lemon green

- gloss 9000100170-005
 matt 9000100170-105 **€ 33.30**

Violet

- gloss 9000100170-044

Yellow RAL 1023

- gloss 9000100170-046
 matt 9000100170-146 **€ 33.30**

Red RAL 3001

- gloss 9000100170-047
 matt 9000100170-147 **€ 33.30**

Hub drive system

Shimano "Nexus"

7-speed hub gear

- 9113000014

Gear ratio: Sprocket 32 teeth
Cassette 20 teeth (with bottom bracket support L – 19 teeth)
with rear pedal brake
Shimano V-brake brake lever, left side

Shimano "Nexus"

8-speed hub gear

- 9113000015 **€ 138.50**

"Grip Shift" twist grip shifter
Gear ratio: Sprocket 32 teeth,
Cassette 20 teeth (with bottom bracket support L – 19 teeth)
with rear pedal brake
Shimano V-brake brake lever, left side

Shimano "Alfine"

11-speed hub gear

- 9113200026 **€ 639.60**

Handle "Rapid Fire"
Gear ratio: Sprocket 32 teeth,
Cassette 20 teeth
(with bottom bracket support L – 19 teeth)



Rohloff "Speedhub 500/14"

14-speed hub gear

- 9113000016 **€ 1,549.60**

"Grip Shift" twist grip shifter
Gear ratio: Sprocket 32 teeth,
Cassette 16 teeth
with hydraulic disk brake (Magura MT2)
brake lever, right side
and Shimano V-brake brake handle, left side

Gear change accessories hub gear

Mountain drive 2.5: 1

Hill gear ratio in combination with hub gears

- 9100200003 **€ 499.80**

Planetary gear set for doubling of the number of gears.
All gears (the hub gear and the gear reduction for hills) can be changed when stationary!
Only in combination with grip size 40 or 45 cm!



Dérailleur

Shimano Deore 3x9 Rapide Fire 27-speed dérailleur

9113000017 € 609.40

Gear ratio: sprocket 22-32-44 teeth, cassette 11-13-15-17-20-23-26-30-32 teeth with hydraulic disk brake (Magura MT2) brake lever, right side and Shimano V-brake, brake lever left side

Shimano Deore XT 1x11 Rapid Fire 11- speed chain shift dérailleur

9113200034 € 690.00

mechanical dérailleur operation from the handle. 11-way cassette 11-42 teeth

Shimano Deore XT 1x11 Di2 Button operation 11- speed chain shift dérailleur

9113200035 € 1,725.00

electronic dérailleur operation from the handle. 11-way cassette 11-42 teeth

Shimano Deore XT 1x11 Di2 Rapid Fire 11- speed chain shift dérailleur

9113200032 € 1,688.00

electronic dérailleur operation from the handle. 11-way cassette 11-42 teeth

Shimano Deore XT 1x11 Di2 chin operation 11- speed chain shift dérailleur

9113200033 € 1,795.00

electronic dérailleur with chin operation (Execution Tetra). 11-way cassette 11-42 teeth



Brake system

Powerful HD brake system "Magura"

9100300101 € 78.60

Recommended from a body weight of 80 kg




Chassis

Chassis geometry

9113100000-035

Length 102 cm
Seat width approx. 35 cm



Dimensions


Pedal bearing shaft

Grip width 40 cm

9100500100-040

Grip width 45 cm

9100500100-045




Fixed bottom bracket support

XS (X/Y 113/130)

9100501100-004

S (X/Y 67/153)

9100501100-003



Adjustable bottom bracket support

9100500001 € 182.90

Adjustable angle via clamped hinge on the fork bridge and in their length.

For the pre-setting, please specify the bottom bracket support size that comes as close to your desired position.



Size: _____

If not specified, the adjustable bottom bracket support size S will be delivered.

Pedal crank

Length 135 mm


9100600100-135

Length 155 mm

9100600100-155

Length 175 mm

9100600100-175



Customer-specific crank length according to dimension _____ mm (minimum length 125 mm)

9100600101 € 155.80

Leg rest

Leg rest separated

9100100008

Adaptation on the fork for short legs

Only possible up to body size 165 cm!



Gear accessories for chain shift

Mountain drive 2.5: 1 gear reduction for hills in combination with chain shift

9100200006 € 371.30

Planetary gear set
Replaces the 3-fold sprocket and allows the gear reduction for hills to be switched when stationary!

Only in combination with grip size 40 or 45 cm!



"Grip Shift" twist grip shifter 3-way, mounting left side

9110200017 € 43.30

SRAM Centera for 3-way dérailleur
In combination with Shimano Deore chain shift!



"Grip Shift" twist grip shifter 9-way, mounting right side

9110200016 € 43.30

SRAM Centera for 9-fold dérailleur
In combination with Shimano Deore chain shift!



Foot stopper for leg rest

9100100012 € 71.90



Thigh anchorage

9100100013 € 127.90



Back equipment

Backrest with filler and 1.5 cm foam padding

9100800000



Body contour back padding

9100800010 € 123.80

with back padding 3 cm
Cover: black/breathable



Seat upholstery

Open belt system seat upholstery

9100801010



Body contour seat upholstery

9100801025
breathable, cushioning effect



Seat cushion

"Standard" seat cushion

9113000190 € 88.30

For 102 cm chassis
Seat width: 35 cm
Height: 4 cm
Cover: black, breathable



Seat cushion is existing

Brand/Type _____

Cushion height _____

Wheels

20x1" back wheels (406)

9100203050
Black rim

20x1" front wheel (406)

9100203052
Black rim

26x1" front wheel with double-wall rim

9100203051
(ETRTO 559)
Silver rim

Tyres

20x1.1" Schwalbe "Durano"

9113105001
(ETRTO: 28-406)



20x2" Schwalbe "Big Apple"

9110105125
(ETRTO: 50-406)



20x1,75" Schwalbe „Marathon Plus“

9113105120 € 58.80
(ETRTO: 47-406)



26x2,0" Schwalbe Marathon Plus Tour

9100105140 € 58.80
(ETRTO: 50-559)



Accessories

Plug-in protection plate

9100502100 € 45.80



Leg guard

9113900000 € 141.50

Guard for MTB tyres

Collision guard clamped on axle tube

9100800100 € 171.80
with 40 mm mounting clip



Neck rest

9100800150 € 206.00
height-adjustable, angle-adjustable



Hip safety belt

9100800200 € 93.50



**SIGMA BC 1200
Bike Computer**

9110901000 € 55.30
wireless, with bracket



Your additional notes:

**SIGMA BC 1606L DTS
Bike Computer**

9110900018 € 135.20
with pedal frequency set
wireless



**Bracket ø 22 mm x 42 mm for fixing
of accessories
(speedometer, lighting or similar)**

9110903271 Item price €
18.30

Colour: natural anodised

_____ Items



**Lighting set SIGMA Lightster LED
with bracket**

9100900019 € 102.50

with battery charger and
set of batteries (4 x NICD)
Front and rear lights



Options for users with reduced function

**Tetra hand support, pivotable,
adjustable grip distance**

9110601101 € 140.30



Tetra shell handles

9110601100 € 135.90



Tetra backhand support

9110601102 € 158.60

pivotable, adjustable distance to the bike
handle



I herewith confirm the details listed above

Location, Date

Signature (End customer)

Location, Date

Signature (Dealer)