

Mobilität wird einfach!

Im Hofstätt 11
D- 72359 Dotternhausen
Phone +49 (0)7427 9480-0
Fax +49 (0)7427 9480-7025
E-Mail info@proactiv-gmbh.de
www.proactiv-gmbh.com



Quotation date: _____

Order date: _____

according to GTC PRO ACTIV GmbH

Medical suppliers/company address:

Cust. no.: _____

Name: _____

Street: _____

Postcode, city: _____

Phone: _____

Fax: _____

E-Mail: _____

Contact: _____

Invoicing address, other (see below)

Cust. no.: _____

Name: _____

Street: _____

Postcode, city: _____

Shipping address, other (see below)

Cust. no.: _____

Name: _____

Street: _____

Postcode, city: _____

Customer name/commission:

Name: _____

Street: _____

Postcode, city: _____

Phone: _____

Fax: _____

E-Mail: _____

DoB: _____

Paying organisation: _____

Indication: _____



NJ1 el. Compact bike RS

Maximum user weight limit 100 kg

5 different drive versions for different purposes,
equipped with:

- girder chassis with integrated rechargeable battery holder
- drive unit, can be decoupled via the eccentric quick-lock system
- A-shaped handles with integrated switching and brake fittings
- infinitely adjustable backrest
- 6° wheel camber
- hydraulic disc brake system
- secondary V-brake (rim brake) with locking mechanism
- safety flag

Services:

1. Delivery, instruction and handing over at the customer by PRO ACTIV

9990000003 € 425.00

plus shipping costs

Body size

Body height: _____ cm

Weight: _____ kg

Drive unit colour

Silver,
anodised drive unit

9110100170-000

Black Shadow,
deep black drive unit
RAL 9005, smooth gloss

9110100170-001 € 103.80

- fork tube R and L
- fork bridge, top and bottom
- bottom bracket support

Neodrives

Neodrives
"The powerful hub drive for max.
performance with TFT display" and
27-speed dérailleur

Torque 12/40 Nm, power 250 Watt



Wheel size: 26"

9100100010-026 € 6,540.00

Neodrives battery

The ranges specified depend on the range
profile.

Battery neodrives G1 36V 10.75 Ah /
387Wh, range approx. 70 km

9100900105



Battery neodrives G2 36V 14.25 Ah /
513Wh, range approx. 70-90 km

9100900106 € 119.00



Double rechargeable battery,
neodrives G1

2 pcs 36V 10.75Ah / 387Wh,
range approx. 100 km

9100900107 € 1,017.00

Quick-change fitting for easy change between
batteries using a magnetic plug

Double rechargeable battery,
neodrives G2

2 pcs 36V 14.25Ah / 513Wh,
range approx. 130 km

9100900108 € 1,232.00

Quick-change fitting for easy change between
batteries using a magnetic plug

Drive, BionX

BionX PL 250HT IGH2
"The comfortable drive with
integrated 2-speed automatic free-
running shifting hub", drive
assistance "0 km/h"

Torque 9/40 Nm, power 250 Watt

Wheel size: 26"

9100100020-026 € 6,520.00

BionX PL 250HT IGH3
"The comfortable drive with
integrated 3-speed
free-running shifting hub and twist
grip shifter"

Torque 9/40 Nm, power 250 Watt

Wheel size: 26"

9100100030-026 € 6,400.00

BionX PL 250 HT
"The powerful drive for permanent
support" with 27-speed dérailleur

Torque 9/40 Nm, power 250 Watt

Wheel size: 26"

9100100040-026 € 6,198.00

BionX D series
"The powerful drive for demanding
users" with 27-speed dérailleur

Torque 25/50 Nm, power 250 Watt

Wheel size: 26"

9100100050-026 € 6,820.00

Only in combination with powerful HD brake
system and second independent disk brake
system.

Not in combination with V-brake rim brakes.

BionX, rechargeable batteries and accessories

The ranges specified depend on the range
profile.

BionX, rechargeable battery 48V 8.8
Ah 423 Wh, range: approx. 55 km

9100900100

Battery BionX 48V 11.6 Ah 556Wh,
range approx. 70 km

9100900102 € 208.00

BionX, double rechargeable battery
48V 8.8 Ah 423 Wh including
exchange plug connection,
range: approx. 110 km

9100900103 € 1,133.00

Changer plug connection for easy change
between the batteries.

BionX, double rechargeable battery
48V 11.6 Ah
556Wh including change plug connection,
 range. approx. 140 km

9100900104 € 1,320.00

Changer plug connection for easy change between the batteries.

Thumb operating switch for the starting assistant mounted on the handle

9116200000 € 127.00

Assembly:

- left side
 right side

Only in combination with BionX drive, not in combination with Neodrives Z10!



Gear change accessories

Mountain drive 2.5: 1 gear reduction for hills in combination with BionX PL 250HT IGH2 and IGH3 to expand to double number of gears

9116200003 € 499.80

Planetary gear set
 All gears (3 gears of the the hub gear and the gear reduction for hills) can be changed when stationary!
 In combination with grip width 40 or 45 cm!



Mountain drive 2.5: 1 gear reduction for hills for BionX PL 250HT and neodrive – 18 gear –

9116200004 € 371.30

Planetary gear set
 Replaces the 3-fold sprocket and allows the gear reduction for hills to be switch when stationary!
 In combination with grip width 40 or 45 cm!



Dérailleur

Shimano Deore 3x9 Rapide Fire 27-speed dérailleur

9100200018

Gear ratio:
 sprocket 22-32-44 teeth,
 cassette 11-13-15-17-20-23-26-30-32 teeth
 Not in combination with BionX PL 250HT IGH2 and BionX PL 250HT IGH3!



Shimano Deore XT 1x11 Rapid Fire 11- speed chain shift dérailleur

9110200034 € 151.00

mechanical dérailleur operation from the handle. 11-way cassette 11-42 teeth



Shimano Deore XT 1x11 DI2 Button operation 11- speed chain shift dérailleur

9110200035 € 1,175.00

electronic dérailleur operation from the handle. 11-way cassette 11-42 teeth



Shimano Deore XT 1x11 DI2 Rapid Fire 11- speed chain shift dérailleur

9110200032 € 1,138.00

electronic dérailleur operation from the handle. 11-way cassette 11-42 teeth



Shimano Deore XT 1x11 DI2 chin operation 11- speed chain shift dérailleur

9110200033 € 1,245.00

electronic dérailleur with chin operation (Execution Tetra).
 11-way cassette 11-42 teeth



Brake system

Powerful HD brake system "Magura"

9100300101 € 78.60

Recommended from a body weight of 80 kg
 Not required with hydraulic back pedal disk brake selected



Second independent disk brake system "Magura"

9100300100 € 183.90

Shimano V-Brake rim brake omitted!
 Only in combination with the powerful HD brake system "Magura"!

Hydraulic back-peddalling disc brake with crank release function for moving backwards

9110300003 € 618.80


Compatible with chain shift and free-running shifting hub
 In combination with:
 - pedal bearing shaft grip widths 45 to 55 cm
 - pedal crank crank lengths 155 or 175 mm up to 120 kg body weight
 Not in combination with Mountaindrive 2.5:1 gear reduction for hills!
 Not in combination with the bottom bracket support XS!



Chassis

Chassis geometry

- Length 104 cm
Seat width approx. 39 cm
9100100000-001
- Length 116 cm
Seat width approx. 39 cm
9100100000-002
- Length 116 cm **€ 161.20**
Seat width approx. 43 cm
9100100000-003
- Heavy Duty version **€ 335.50**
up to 120 kg body weight
9100800013
Recommended for off-road use



Leg rest

Leg rest separated

- 9100100008
Adaptation on the fork for short legs
Only possible up to body size 165 cm!




Leg rest wrapping

- 9100100009
Adaptation on the fork



Leg rest adapted to the chassis

- 9100100010 **€ 400.40**
Legs not influenced by the steering motion, increases the vehicle width!
Not in combination with the leg shield.



Dimensions

Pedal bearing shaft

- Grip width 40 cm
9100500100-040
- Grip width 45 cm
9100500100-045
- Grip width 50 cm
9100500100-050




Foot stopper for leg rest

- 9100100012 **€ 71.90**



Fixed bottom bracket support

- XS (X/Y 113/130)
9100501100-004
- S (X/Y 67/153)
9100501100-003
- M (X/Y 90/250)
9100501100-001



Thigh anchorage

- 9100100013 **€ 127.90**



Adjustable bottom bracket support

- 9100500001 **€ 182.90**

Adjustable angle via clamped hinge on the fork bridge and in their length.

For the pre-setting, please specify the bottom bracket support size that comes as close to your desired position.

Size: _____

If not specified, the adjustable bottom bracket support size S will be delivered.



Back equipment

Backrest with filling plate

- 9100800000
Foam padding 1.5 cm



Body contour back padding

- 9100800010 **€ 123.80**
with back padding 3 cm
Cover: black/breathable



Pedal crank

- Length 155 mm
9100600100-155
- Length 175 mm
9100600100-175
- Length 195 mm
9100600100-195



Seat upholstery

Open belt system seat upholstery

- 9100801010



Body contour seat upholstery

- 9100801025



Seat cushion

shape dimensions in cm	seat cushion "Standard" 9100000190	seat cushion with upper leg guidance 9100800103
Chassis: 104 SB: 39 ST: 30	-001 <input type="checkbox"/>	-001 <input type="checkbox"/>
Chassis: 116 SB: 39 ST: 40	-002 <input type="checkbox"/>	-002 <input type="checkbox"/>
Chassis: 116 SB: 43 ST: 40	-003 <input type="checkbox"/>	-003 <input type="checkbox"/>
price	€ 88,30	€ 161,10

Cover: black, breathable, open fabric surface
height: 4cm

seat cushion with upper leg guidance only in combination with hip safety belt.

Tyres

Schwalbe Marathon

- 26x1,5" (559-35)
9100105100

Touring tyres



Schwalbe Marathon Plus Tour

- 26x2,0" (559-50) € 58.80
9100105140



Schwalbe Smart Sam

- 26x2,1" (559-54) € 58.80
9100105150



Bike accessories

Plug-in protection plate

- 9100502100 € 45.80



Fast-Charger Neodrives

- 9116900006 € 66.40

Leg guard

- 9100100231 € 141.50

Screwed onto the fork tube
Colour: black matt
Recommended for off-road use or
MTB tyres



Collision guard clamped on the axle tube

- 9100800100 € 171.80



Neck rest

- 9100800150 € 206.00

height-adjustable
angle-adjustable



Hip safety belt

- 9100800200 € 93.50




Bracket ø 22 mm x 42 mm for fixing of accessories (speedometer, lighting or similar)

- 9110903271 Item price € 18.30

Colour: natural anodised

_____ Items



Lighting set SIGMA Lightster LED with bracket

- 9100900019 € 102.50

with battery charger and set of batteries
(4 pcs. NICD) front and rear lights



Options for users with reduced function

Tetra hand support, pivotable, adjustable grip distance

- 9110601101 € 140.30



Tetra shell handles

- 9110601100 € 135.90



Tetra backhand support

- 9110601102 € 158.60

pivotable, adjustable distance to the bike handle



Chassis colour

Information on current colour designs and example images under:
www.proactiv-gmbh.de/farben.php

White RAL 9016

- gloss 9000100170-023
 matt 9000100170-123 € 33.30

Brilliant silver

- gloss 9000100170-011

Anthracite RAL 7016

- gloss 9000100170-001
 matt 9000100170-101 € 33.30

Black RAL 9005

- gloss 9000100170-003
 matt 9000100170-019 € 33.30

Deep black, silver effect

- gloss 9000100170-007

Pure orange RAL 2004

- gloss 9000100170-022

Fire Red

- gloss 9000100170-018
 matt 9000100170-118 € 33.30

Pink magenta RAL 4010

- gloss 9000100170-037
- matt 9000100170-137 **€ 33.30**

Sky blue

- gloss 9000100170-017

Pastel blue RAL 5024

- gloss 9000100170-036
- matt 9000100170-136 **€ 33.30**

Marine blue

- gloss 9000100170-045
- matt 9000100170-145 **€ 33.30**

Apple green

- gloss 9000100170-014
- matt 9000100170-114 **€ 33.30**

Lemon green

- gloss 9000100170-005
- matt 9000100170-105 **€ 33.30**

Violet

- gloss 9000100170-044

Transparent, powder coated

- gloss 9000100170-029

Yellow RAL 1023

- gloss 9000100170-046
- matt 9000100170-146 **€ 33.30**

Red RAL 3001

- gloss 9000100170-047
- matt 9000100170-147 **€ 33.30**

Your additional notes:

I herewith confirm the details listed above

Location, Date

Signature (End customer)

Location, Date

Signature (Dealer)