

Mobilität wird einfach!

Im Hofstätt 11
D- 72359 Dotternhausen
Phone +49 (0)7427 9480-0
Fax +49 (0)7427 9480-7025
E-Mail info@proactiv-gmbh.de
www.proactiv-gmbh.com



Quotation date: _____

Order date: _____

according to GTC PRO ACTIV GmbH

Medical suppliers/company address:

Cust. no.: _____
Name: _____
Street: _____
Postcode, city: _____
Phone: _____
Fax: _____
E-Mail: _____
Contact: _____

Invoice address different (see below)

Cust. no.: _____
Name: _____
Street: _____
Postcode, city: _____

Shipping address different (see below)

Cust. no.: _____
Name: _____
Street: _____
Postcode, city: _____

Customer name/commission:

Name: _____
Street: _____
Postcode, city: _____
Phone: _____
Fax: _____
E-Mail: _____
DoB: _____
Paying organisation: _____
Indication: _____



Raptor "4all"

Maximum user weight limit 100 kg

- torsion-resistant aluminium tube chassis
- drive unit "black" with adjustable pedal position
- Shimano Deore 3x9 Rapid Fire chain gear
- A-shaped handles with integrated switching and brake fittings
- Magura hydraulic disk brake system, brake lever, right side
- Shimano V-Brake rim brake, brake handle, left side
- infinitely adjustable backrest
- 1° wheel camber
- safety flag

9115810000 € 4,120.00

Services:

1. Delivery, instruction and handing over at the customer by PRO ACTIV

9990000003 € 425.00
plus shipping costs

Body size

Body height: _____ cm

Weight: _____ kg

Chassis colour

Information on current colour designs and example images under:
www.proactiv-gmbh.de/farben.php

White RAL 9016

- gloss 9000100170-023
- matt 9000100170-123 **€ 33.30**

Brilliant silver

- gloss 9000100170-011

Anthracite RAL 7016

- gloss 9000100170-001
- matt 9000100170-101 **€ 33.30**

Black RAL 9005

- gloss 9000100170-003
- matt 9000100170-019 **€ 33.30**

Deep black, silver effect

- gloss 9000100170-007

Pure orange RAL 2004

- gloss 9000100170-022

Fire Red

- gloss 9000100170-018
- matt 9000100170-118 **€ 33.30**

Pink magenta RAL 4010

- gloss 9000100170-037
- matt 9000100170-137 **€ 33.30**

Sky blue

- gloss 9000100170-017

Pastel blue RAL 5024

- gloss 9000100170-036
- matt 9000100170-136 **€ 33.30**

Marine blue

- gloss 9000100170-045
- matt 9000100170-145 **€ 33.30**

Apple green

- gloss 9000100170-014
- matt 9000100170-114 **€ 33.30**

Lemon green

- gloss 9000100170-005
- matt 9000100170-105 **€ 33.30**

Violet

- gloss 9000100170-044

Yellow RAL 1023

- gloss 9000100170-046
- matt 9000100170-146 **€ 33.30**

Red RAL 3001

- gloss 9000100170-047
- matt 9000100170-147 **€ 33.30**

Powder coating, transparent

- gloss 9000100170-029

Dérailleur

Shimano Deore 3x9 Rapid Fire 27-speed dérailleur

- 9115200002
- Gear ratio: sprocket 22-32-44 teeth,
cassette 11-13-15-17-20-23-26-30-32 teeth

Shimano Deore XT 1x11 Rapid Fire 11- speed chain shift dérailleur

- 9110200034 **€ 151.00**
- mechanical derailleur operation from the
handle. 11-way cassette 11-42 teeth



Shimano Deore XT 1x11 Di2 Rapid Fire 11- speed chain shift dérailleur

- 9110200032 **€ 1,138.00**
- electronic derailleur operation from the
handle. 11-way cassette 11-42 teeth

Hub drive system

Rohloff "Speedhub 500/14" 14-speed hub gear

- 9115200004 **€ 1,038.00**
- "Grip Shift" twist grip shifter
Gear ratio: sprocket 32 teeth, cassette 16
teeth

Gear accessories for chain shift

Competition gear ratio


- 9115200010 **€ 193.90**

Sprocket 52/42/30
 - 5-arm chain wheel holder
 - chain wheels Shimano Ultegra
 incl. custom aluminium chain guard

Chassis

Chassis geometry

- 9115810001
Seat width 38 cm length
122 cm
- 9115810002
Seat width 42 cm length
122 cm



Dimensions

Pedal bearing shaft

- Grip width 40 cm
9115500100-040
- Grip width 45 cm
9115500100-045
- Grip width 50 cm
9115500100-050



Adjustable bottom bracket support

911550001

Adjustable angle via clamped hinge on the fork bridge and in their length.

For the pre-setting, please specify the bottom bracket support size that comes as close to your desired position.

Size: _____

If not specified, the adjustable bottom bracket support size S will be delivered.


Pedal crank

Length 175 mm
9115600100-175

Length 195 mm
9115600100-195



Drive unit welded permanently



Dimension: X=_____/Y=_____

Note:

Position of the pedal bracket is welded fixed and cannot be changed later.

We therefore recommend defining the dimensions for the drive unit with adjustable pedal support in advance!

9115400000 € 717.60

Brake system

Powerful HD brake system "Magura"

9100300101 € 78.60

Recommended from a body weight of 80 kg



2nd independent V-brake rim brake system

9115300000 € 130.00

Individual checking required as risk of collision with the chain depending on the chain guide/pedal position. Not in combination with hydraulic disk brake system.

Leg rest

Leg rest wrapping

9115100009
Adaptation on the fork



Foot stopper for leg rest

9115100012 € 71.90



Back equipment

Backrest adjustable with filling plate

9115800001
Infinitely adjustable backrest with 3 cm removable back cushion, black, breathable



Body contour back padding

9115800010 € 113.40

with back cushion, 3 cm thick
Cover: black, breathable

Seat upholstery

Open belt system seat upholstery

9115801010



Seat cushion

"Standard" seat cushion

9115000190 € 88.30

For 116 cm chassis
Seat width: 39 cm
Depth: 39 cm
Height: 4 cm



Cover: black/breathable

Running wheels

Running wheels with double-wall rim

26" (559)
9115104260

Tyres

Schwalbe "Durano"

26x1.1" (559-28)
9115105000

Racing-Slic

Schwalbe "Marathon"

26x1,5" (559-40)
9100105100

Touring tyres



Bike accessories

Your additional notes:

Plug-in protection plate

9115502100 € 45.80



Collision guard clamped on the axle tube

9115800100 € 171.80

with 40 mm mounting clip



Neck rest

9115800150 € 206.00

height-adjustable
angle-adjustable



Hip safety belt

9115800200 € 93.50



SIGMA BC 1200 Bike Computer

9115901000 € 55.30

wireless
with bracket



SIGMA BC 1606L DTS Bike Computer

9115900018 € 135.20

with pedal frequency set
wireless



Lighting set SIGMA Lightster LED with bracket

9100900019 € 102.50

with battery charger and
set of batteries (4 x NICD)
Front and rear lights



Options for users with reduced function

Tetra hand support, pivotable, adjustable grip distance

9110601101 € 140.30



Tetra shell handles

9115601100 € 135.90



Tetra backhand support

9110601102 € 158.60

pivotable, adjustable distance to the bike handle



I herewith confirm the details listed above

Location, Date

Signature (End customer)

Location, Date

Signature (Dealer)