

Mobilität wird einfach!

Im Hofstätt 11
D- 72359 Dotternhausen
Phone +49 (0)7427 9480-0
Fax +49 (0)7427 9480-7025
E-Mail info@proactiv-gmbh.de
www.proactiv-gmbh.com



Raptor "4you"

Maximum user weight limit 100 kg

- torsion-resistant aluminium tube chassis
- drive unit "black" with adjustable pedal position
- Shimano Deore 3x9 Rapid Fire chain gear
- A-shaped handles with integrated switching and brake fittings
- Magura hydraulic disk brake system, brake lever, right side
- Shimano V-Brake rim brake, brake handle, left side
- infinitely adjustable backrest
- 1° wheel camber
- safety flag

9115004000 € 4,660.00

Services:

1. Delivery, instruction and handing over at the customer by PRO ACTIV

- 9990000003 € 400.00
plus shipping costs

Quotation date: _____

Order date: _____

according to GTC PRO ACTIV GmbH

Medical suppliers/company address:

Cust. no.: _____
Name: _____
Street: _____
Postcode, city: _____
Phone: _____
Fax: _____
E-Mail: _____
Contact: _____

Invoicing address, other (see below)

Cust. no.: _____
Name: _____
Street: _____
Postcode, city: _____

Shipping address, other (see below)

Cust. no.: _____
Name: _____
Street: _____
Postcode, city: _____

Customer name/commission:

Name: _____
Street: _____
Postcode, city: _____
Phone: _____
Fax: _____
E-Mail: _____
DoB: _____
Paying organisation: _____
Indication: _____

Body size

Body height: _____ cm

Weight: _____ kg

Chassis colour

Information on current colour designs and example images under:
www.proactiv-gmbh.de/farben.php

White RAL 9016

- gloss 9000100170-023
- matt 9000100170-123 **€ 33.30**

Brilliant silver

- gloss 9000100170-011

Anthracite RAL 7016

- gloss 9000100170-001
- matt 9000100170-101 **€ 33.30**

Black RAL 9005

- gloss 9000100170-003
- matt 9000100170-019 **€ 33.30**

Deep black, silver effect

- gloss 9000100170-007

Pure orange RAL 2004

- gloss 9000100170-022

Fire Red

- gloss 9000100170-018
- matt 9000100170-118 **€ 33.30**

Pink magenta RAL 4010

- gloss 9000100170-037
- matt 9000100170-137 **€ 33.30**

Sky blue

- gloss 9000100170-017

Pastel blue RAL 5024

- gloss 9000100170-036
- matt 9000100170-136 **€ 33.30**

Marine blue

- gloss 9000100170-045
- matt 9000100170-145 **€ 33.30**

Apple green

- gloss 9000100170-014
- matt 9000100170-114 **€ 33.30**

Lemon green

- gloss 9000100170-005
- matt 9000100170-105 **€ 33.30**

Violet

- gloss 9000100170-044

Yellow RAL 1023

- gloss 9000100170-046
- matt 9000100170-146 **€ 33.30**

Red RAL 3001

- gloss 9000100170-047
- matt 9000100170-147 **€ 33.30**

Powder coating, transparent

- gloss 9000100170-029

Dérailleur

Shimano Deore 3x9 Rapid Fire 27-speed dérailleur

- 9115200002
- Gear ratio: sprocket 22-32-44 teeth,
cassette 11-13-15-17-20-23-26-30-32 teeth

Shimano Deore XT 1x11 Rapid Fire 11- speed chain shift dérailleur

- 9110200034 **€ 151.00**
- mechanical dérailleur operation from the
handle. 11-way cassette 11-42 teeth

Shimano Deore XT 1x11 Di2 Rapid Fire 11- speed chain shift dérailleur

- 9110200032 **€ 1,138.00**
- electronic dérailleur operation from the
handle. 11-way cassette 11-42 teeth



Hub drive system

Rohloff "Speedhub 500/14" 14-speed hub gear

- 9115200004 **€ 1,038.00**

"Grip Shift" twist grip shifter
Gear ratio: sprocket 32 teeth, cassette 16
teeth

Gear accessories for chain shift

Competition gear ratio


- 9115200010 **€ 193.90**

Sprocket 52/42/30
- 5-arm chain wheel holder
- chain wheels Shimano Ultegra
incl. custom aluminium chain guard

Chassis

Chassis geometry individual

- Dimensional
sketch, see
appendix



Dimensions

Pedal bearing shaft

- Grip width 40 cm
9115500100-040
- Grip width 45 cm
9115500100-045
- Grip width 50 cm
9115500100-050



Adjustable bottom bracket support

911550001

The angle and length can be adjusted via the clamped rotary point on the fork bridge.

For the pre-setting, please specify the bottom bracket support size that comes as close to your desired position.

Size: _____

If no details are provided, the adjustable bottom bracket support shall be supplied in S


Pedal crank

Length 175 mm
9115600100-175

Length 195 mm
9115600100-195



Drive unit permanently welded



Dimension: X=_____/Y=_____

Note:

The position of the pedal bracket is welded fixed and cannot be changed later. We therefore recommend defining the dimensions for the drive unit with adjustable pedal support in advance!

9115400000 € 717.60

Brake system

Powerful HD brake system "Magura"

9100300101 € 78.60

Recommended from a body weight of 80 kg



2nd independent V-brake rim brake system

9115300000 € 130.00

Individual checking required as risk of collision with the chain depending on the chain guide/pedal position. Not in combination with hydraulic disk brake system.

Leg rest

Leg rest wrapping

9115100009
Adaptation on the fork



Foot stopper for leg rest

9115100012 € 71.90



Back equipment

without back rest

9115800013
Prepared for customer's installation of a complete back/seat system.



Backrest adjustable with filling plate

9115800001
Infinitely adjustable backrest with 3 cm removable back cushion, black, breathable



Body contour back padding

9115800010 € 113.40

with back cushion, 3 cm thick
Cover: black, breathable

Seat upholstery

Open belt system seat upholstery

9115801010



Seat cushion

"Standard" seat cushion

9115000190 € 88.30

For 116 cm chassis
Seat width: 39 cm
Depth: 39 cm
Height: 4 cm



Cover: black, breathable

Running wheels

Running wheels with double-wall rim

26" (559)
9115104260

Tyres

Schwalbe "Durano"

26x1.1" (559-28)
9115105000

Racing-Slic



Schwalbe "Marathon"

26" (559-35)
9100105100
Touring tyres



Tetra backhand support

9110601102 € 158.60
pivotal, adjustable distance to the bike handle



Bike accessories

Plug-in protection plate

9115502100 € 45.80



Collision guard clamped on axle tube

9115800100 € 171.80
with 40 mm mounting clip



Neck rest

9115800150 € 206.00
height-adjustable
angle-adjustable



Hip safety belt

9115800200 € 93.50



SIGMA BC 1200 Bike Computer

9115901000 € 55.30
wireless
with bracket



SIGMA BC 1606L DTS Bike Computer

9115900018 € 135.20
with pedal frequency set
wireless



Lighting set SIGMA Lightster LED with bracket

9100900019 € 102.50
with battery charger and set of batteries (4 x NICD)
Front and rear lights



Options for users with reduced function

Tetra hand support, pivotal, adjustable grip distance

9110601101 € 140.30



Tetra shell handles

9115601100 € 135.90



Your additional notes:

I herewith confirm the details listed above

Location, Date

Signature (End customer)

Location, Date

Signature (Dealer)

