

Im Hofstätt 11  
D- 72359 Dotternhausen  
Tel +49 (0)7427 9480-0  
Fax +49 (0)7427 9480-7025  
e-mail info@proactiv-gmbh.de  
www.proactiv-gmbh.de



**Quotation date:**

**Order date:**

according to GTC PRO ACTIV GmbH

**Medical suppliers / Company address:**

Cust.no.:

Name: \_\_\_\_\_

Street: \_\_\_\_\_

Postcode, city: \_\_\_\_\_

Phone: \_\_\_\_\_

Fax: \_\_\_\_\_

E-Mail: \_\_\_\_\_

Contact: \_\_\_\_\_

**Invoicing address, other** (see below)

Cust.no.:

Name: \_\_\_\_\_

Street: \_\_\_\_\_

Postcode, city: \_\_\_\_\_

**Shipping address, other** (see below)

Cust.no.:

Name: \_\_\_\_\_

Street: \_\_\_\_\_

Postcode, city: \_\_\_\_\_

**Customer name / Commission:**

Name: \_\_\_\_\_

Street: \_\_\_\_\_

Postcode, city: \_\_\_\_\_

Phone: \_\_\_\_\_

Fax: \_\_\_\_\_

E-Mail: \_\_\_\_\_

Date of birth.: \_\_\_\_\_

Paying organisation: \_\_\_\_\_

Indication: \_\_\_\_\_



**Wheel-e**

up to 120 kg trailer load and 10 kg load capacity  
2 different basic models for different purposes,  
equipped with:

- docking plate for quick coupling to the "Quick-Snap Adapter"
- Drive signal transmitter and display screen with support up to 6km/h
- hydraulic disc brake system, disc Ø180mm
- parking stand
- balance weight (only single battery)
- Handle bar sports edition, mounted on the angle adjustable stem
- lighting set Sigma with rechargeable battery charger, front and rear lights
- bell

**Services:**

**1. Consultation, definition of dimensions and testing at the customer by:**

- 9990000000 **PRO ACTIV**
- 9990000004 **Dealer**

**2. Assembly of the frame adapter at PRO ACTIV**

- 9990000001  
on newly ordered PRO ACTIV wheelchairs (assembly included)  
Prod. no.: \_\_\_\_\_
- 9990000002 **€ 125,00**  
on existing PRO ACTIV wheelchairs  
Prod. no.: \_\_\_\_\_
- 9990000005  
on third party product by the dealer (compatibility check by the dealer)  
Product: \_\_\_\_\_

**3. Delivery, instruction and handing over at the customer by PRO ACTIV**

- 9990000003 **€ 400,00**  
Plus shipping costs

**Body size**

Body height: \_\_\_\_\_ cm

upper body height: \_\_\_\_\_ cm

Weight: \_\_\_\_\_ kg

Information are required for processing!

## Colour

**Silver,  
anodised drive unit**

9110100170-000

**Black Shadow,  
deep black drive unit  
RAL 9005, smooth gloss**

9110100170-001 **€ 99,80**

- Fork tube R and L  
- Fork bridge, top and bottom  
- Bottom bracket support

## Handlebar

### Steering column

angle adjustable at the triple clamp

9130600001  
length of the steering column  
220mm

9130600002  
length of the steering column  
280mm

9130600003  
length of the steering column  
340mm



## 2nd independent disc brake system "Magura"

9130300100 **€ 176,80**

Shimano V-brake rim brake omitted

### Parking brake

9130300102 **€ 64,00**

mounted on the steering column.  
V-Brake omitted



## BionX D-Series

"The powerful drive for demanding users"

6 km/h: torque 25/50 Nm, Power 250 Watt  
Wheel size: 20", double brake system  
"Magura MT4"

9130100020-020 **€ 4.405,00**  
20" (ETRTO 406)

### rim colour silver

9130200105

drive wheel silver matt coated

### rim colour orange RAL 2004

9130200100

drive wheel orange gloss coated

### Schwalbe "Knobby"

Profiled tyres

9130105126  
20x2" (406-54)



### Schwalbe "Big Ben"

Profiled tyres

9130105127 **€ 22,90**  
20x2,15" (406-55)



## Powerful double Brake system "Magura MT5" with 2x 4-piston-brake caliper and brake disc Ø203mm

9130300105 **€ 156,00**

Recommended from a body weight of 80 kg

## BionX

### BionX PL 250 HT

"The comfortable drive for permanent support" with 27-speed derailleur

6km/h: torque 9/40 Nm, Power 250 Watt  
Wheel size: 20", rim colour: black  
Additional brake V-Brake, double brake system "Magura MT4"

9130100010-020 **€ 3.500,00**  
20" (ETRTO 406)



### Schwalbe "Marathon"

Touring tyres

9130105120  
20x1,5" (406-40)



### Schwalbe "Knobby"

Profiled tyres

9130105126  
20x2" (406-54)



### Schwalbe "Marathon Plus"

unpuncturable, with anti-puncture belt  
Smart Guard

9130105128 **€ 22,90**  
20x1,75" (406-47)



## Powerful HD brake system "Magura MT5" with 4-piston-brake caliper

9130300101 **€ 82,50**

Recommended from a body weight of 80 kg



## BionX, accessories and rechargeable batteries

### Console / Display for 15 km/h

9130900100 **€ 108,00**

removable



### Console / Display for 24 km/h

9130900105 **€ 108,00**

removable



The ranges specified depend on the range profile.

### Rechargeable battery BionX

**48V 8,8Ah 423Wh,**  
range approx. 30-50km

9130900010  
Assembly to the right on the fork



### Rechargeable battery BionX

**48V 11,6 Ah 556Wh,**  
range approx. 40-60km

9130900011      **€ 200,00**  
Assembly to the right on the fork



### Double rechargeable battery BionX

**48V 8,8Ah 423Wh including  
exchange plug connection,**  
range approx. 50-70km

9130900012      **€ 1.090,00**  
Assembly on both sides of the fork, exchange  
plug connection mounted on the fork to allow  
the simple exchange between the  
rechargeable batteries.



### Double rechargeable battery BionX

**48V 11,6 Ah  
555Wh including  
exchange plug connection,**  
range approx.. 60-80km

9130900013      **€ 1.270,00**  
Assembly on both sides of the fork, exchange  
plug connection mounted on the fork to allow  
the simple exchange between the  
rechargeable batteries.



## Accessoires

### SKS fender

9130500003      **€ 79,00**  
removable



### fender MX

9130500002      **€ 129,00**  
fixed assembled  
Colour: orange



### fender MX

9130500004      **€ 129,00**  
fixed assembled  
Colour: black



### Manoeuvring rollers for stands

9110900000      **€ 92,50**  
Manoeuvring rollers mounted to the parking  
stand via dual clamping bracket. (figure  
similar)



### Aluminium luggage rack "cargo systems"

9110903290      **€ 77,10**  
Colour: black  
Maximum load weight of 10 kg on the  
loading platform,  
max. 15 kg in combination with luggage  
bags with low loading height.



### Bottle holder with drinking bottle

9130900001      **€ 18,00**  
Mounted on the rechargeable battery rail,  
right side.  
Not in combination with double  
rechargeable battery.



### rearview mirror

9130900002      **€ 74,20**  
left side

### rearview mirror

9130900003      **€ 74,20**  
right side

### LED headlight

9130900004      **€ 97,00**  
Supply of the LED front light via the  
rechargeable battery of the drive system.  
Back light as plug-on light supplied by a  
battery  
Lighting set omitted

## Options for users with reduced function

### Quick release clip lever extension for fixed frame adapter

9110700004      **€ 78,80**  
to make the opening and closing of the  
quick release clips easier

### Eccentric connection lever extension for Quick Snap connection

9110700005      **€ 55,90**  
to make the closing or opening of the Quick  
Snap adapter easier. Tetra version.  
Only in combination with fixed frame  
adapter.



### Centring column adaptation assistance

9110700002      **€ 41,70**  
in order to make inserting the central  
adapter easier

## Frame adapters

### SPEEDY fixed frame adapter with "Quick-Snap Adaption"

9110700141 € 392,20

with quick release clips and centring column mount on the frame  
Colour: black



### TRAVELER folding frame adapter with "Quick-Snap Adaption"

9110701300 € 419,20

with quick release clips and centring column mount on the folding mechanism  
Colour: black



### BUDDY fixed frame adapter with "Quick-Snap Adaption"

9110701400 € 392,20

with quick release clips and centring column mount on the frame  
Colour: black



### Wheelchair front adapter with universal clamping and "Quick-Snap Adaption"

9110701500 € 446,80

Adapter mount on the wheelchair front end

mounting situation for PA Products:  
Straight front end from SB (seat width)  
31cm to 54cm  
Front end V-shape from SB  
(seat width)30cm to 50cm



#### PLEASE NOTE:

If the tube diameter is not specified, as standard we shall deliver an adapter for a tube diameter of 23-27 mm!

Please make the following selection:

9110701525

Tube diameter 23-27mm

9110701530

Tube diameter 28-32mm

9110701535

Tube diameter 33-37mm


## Wheelchair-dependent accessories

### Extended Bike plate

For Speedy A1/A2  
9110100100-002 € 280,50

For Buddy  
9110100100-003 € 280,50


For Traveler  
9110100100-001 € 280,50



4° camber integrated, independently adjustable track

### Wheelbase extension removable via "grid mechanism with pulling operation", with continuous axle tube

For Speedy 4all /  
4teen / 4you /  
Traveler 4you  
9110102001 € 428,00



4° / 6° camber integrated, independently adjustable track

### Wheelbase extension removable via "Quick-Pins", with continuous axle tube

Für Speedy F2  
9015101000 € 428,00

Für Speedy F4  
9110101000 € 428,00

4° / 6° camber integrated, independently adjustable track.

Attachment possibility for F2/F4 must be checked by PRO ACTIV!



I herewith confirm the details listed above

Location, Date Signature (End customer)

Location, Date Signature (Dealer)