

31 Appendix: Medical product passport / record of training

Product specifications:

Serial number: _____ Key number: _____

Customer data:

Surname, forename: _____

Street: _____

Postcode, city: _____

Phone: _____

Paying organisation: _____

Training carried out by:

Medical supplies dealer

PROACTIV field representative

Stamp / Date / Dealer's signature

Record of training

I was / we were instructed in accordance with the associated hand-over certificate about the operation of the product listed and informed about possible operator errors. I was / we were also advised about situations where the assistance of another person is required. The operating instructions were handed to me / us.

Instructor

Name, date, signature _____

1. Person being trained

Name, date, signature _____

2. Person being trained

Name, date, signature _____

3. Person being trained

Name, date, signature _____

For minors, or persons who are not responsible for their actions, legal guardians / supervisors / responsible persons are to be trained in the use, this is confirmed by their signature. The data is recorded in the feedback system of PROACTIV Reha-Technik GmbH, as the manufacturer of the above named product. It will be managed in accordance with Section 16 BDSG (Federal Data Protection Law).

32 Appendix: Hand-over certificate

32.1 Required compliance criteria to authorise use

| Topics | Completed / fulfilled | Remarks |
|--|-----------------------|---------|
| The product is suitable for the customer based on their own judgement and the customer information received regarding the disability-related restrictions. | | |
| The use intended by the customer is fully consistent with the intended use as described in the operating instructions (see the Product description / intended use chapter). | | |
| The product's equipment is suitable to allow the customer safe use with maximum reduction of risks (see check list on the following page). | | |
| The customer was informed about the current / applicable regulations in accordance with the road traffic regulations. | | |
| The customer's driving ability was checked during a test drive in difficult driving situations and found to be appropriate (see the check list on the following page). | | |
| The user, according to their own statements, or those of the legal representative or guardian and the assessment of the person providing the training, is able to meet the requirements of public traffic in full and to act accordingly. This ability to act, which is the basis for reducing the risk for the user and other road users to an acceptable level, is also completely achievable taking current illnesses / disabilities into full account. | | |
| The operating instructions, and explicitly all of the warning and safety instructions contained therein, were discussed during the training in detail and understood by the user. The user was then handed these operating instructions. | | |

32.2 Check list for training the user

| Topics | Completed / fulfilled |
|--|-----------------------|
| Advised of the applicable legal regulations when driving on public roads. | |
| All mechanical and electrical functional controls were explained and their function demonstrated. | |
| Adaptation of the adapter to the wheelchair and removing the adapter from the wheelchair was demonstrated and then performed by the user themselves and / or their assistant. | |
| Adaptation and uncoupling the product to / from the wheelchair was demonstrated and then performed by the user themselves and / or their assistant. | |
| Adaptation of the wheelbase extension on the wheelchair – if present – was demonstrated and then performed by the user themselves and / or their assistant. | |
| Repositioning the wheelchair drive wheels from the standard position in the wheelbase extension sockets – if fitted – was demonstrated and then performed by the user themselves and / or their assistant. | |
| How the parking stand works was demonstrated and then performed by the user themselves and / or their assistant. | |
| Operation and basic settings on the display / operating console were demonstrated and then performed by the user themselves and / or their assistant. | |
| How the starting assistant / pushing air – if fitted – was demonstrated and then performed by the user themselves and / or their assistant. | |
| Operation of the drive system and the reaction of the drive system to the various settings was demonstrated and then performed by the user themselves and / or their assistant. | |
| Removal and insertion of the operating console and the rechargeable batteries as well as operation of the exchange plug connection as switching solution – if fitted – was demonstrated and then performed by the user themselves and / or their assistant. | |
| Handling and charging the batteries as well as the charger functions were demonstrated and then performed by the user themselves and / or their assistant. The instructions about charging the batteries during a prolonged period of non-use / storage of the product are important here. | |
| Use of the parking brake and the service brakes was demonstrated and then performed by the user themselves and / or their assistant. | |
| The operation and function of the gears was demonstrated and then performed by the user themselves and / or their assistant. | |
| The operation of the lights – if fitted – was demonstrated and then performed by the user themselves and / or their assistant. | |
| Test drive: Forwards and, if required, backwards travel through 4 cones spaced at 1.5 m or 2 m | |
| Test drive: Driving on the level, uphill and downhill in the direction of travel | |
| Test drive: Emergency stop from maximum speed | |
| Information for care, cleaning and maintenance of the product have been provided and understood by the user and / or assistant. | |
| Information on the wheel with regard to inflation pressure and tread depth have been provided and understood by the user and / or assistant. | |
| Information on regular checks of the brakes have been provided and understood by the user and / or assistant. | |
| Information on checking the gears including cables and lines and the maintenance of the chain have been provided and understood by the user and / or assistant. | |
| The contents of the operating instructions from PRO ACTIV and the other component manufacturers were completely worked through based on the product training and were understood by the user and / or assistant. | |

The use of the product is only permitted when all topics listed in "Required compliance criteria for those permitted to use" have been met by the user and all the points have been ticked off in the "Check list for training the user".