

Serial number:

**Subsequent inspection:** After another year after the last inspection or already beforehand following particular loading

	OK / carried out	not OK	resolved
Check all screws / fastening elements are firmly seated	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Clean and oil / grease all pivot points and bearings	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Visual inspection of the frame and attachments for crack formation, deformation, etc.	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Functional and safety check of the brakes and, where necessary, replacement of the brake pads, brake cables	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Check the spoke tension of the product wheel and, if required, correct the tension / re-centring	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Functional and safety check of the product wheel and, where necessary, replacement of the tyre on the product	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Functional and safety check of the steering, adaptation to the product and the adapted wheelchair / wheelchairs	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Test drive / functional test	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

resolved = the fault was corrected

**Comments:**

Stamp:

\_\_\_\_\_

Date / Signature