

Im Hofstätt 11  
D-72359 Dotternhausen  
Phone +49 (0)7427 9480-0  
Fax +49 (0)7427 9480-7025  
e-mail info@proactiv-gmbh.de  
www.proactiv-gmbh.com



**Quotation date:** \_\_\_\_\_  
**Consultation was carried out by**  
**Staff at the medical supplies shop:** \_\_\_\_\_  
**PRO ACTIV field representative:** \_\_\_\_\_

**Order date:** \_\_\_\_\_  
According to GTC PRO ACTIV Reha-Technik GmbH

### Medical suppliers/company address:

Cust. no.: \_\_\_\_\_  
Name: \_\_\_\_\_  
Street: \_\_\_\_\_  
Postcode / city: \_\_\_\_\_  
Phone: \_\_\_\_\_  
Fax: \_\_\_\_\_  
E-Mail: \_\_\_\_\_  
Contact: \_\_\_\_\_

**Invoicing address, other**(see below)

Cust. no.: \_\_\_\_\_  
Name: \_\_\_\_\_  
Street: \_\_\_\_\_  
Postcode / city: \_\_\_\_\_

**Shipping address, other** (see below)

Cust. no.: \_\_\_\_\_  
Name: \_\_\_\_\_  
Street: \_\_\_\_\_  
Postcode / city: \_\_\_\_\_

### Customer name/commission:

Surname,  
forename: \_\_\_\_\_  
Street: \_\_\_\_\_  
Postcode / city: \_\_\_\_\_  
Phone: \_\_\_\_\_  
D.O.B: \_\_\_\_\_  
Paying  
organisation: \_\_\_\_\_  
Indication: \_\_\_\_\_



## NJ1 el. compact bike

Maximum user weight limit 100 kg

- girder chassis with integrated rechargeable battery holder
- drive unit, can be decoupled via the eccentric quick-lock system
- A-shaped handles with integrated switching and brake fittings
- infinitely adjustable backrest
- 6° wheel camber
- hydraulic disc brake system
- Shimano V-brake (rim brake) with locking mechanism
- safety pennon with holder
- leg guard
- black drive unit

### Neodrives

**"The powerful hub drive for maximum performance  
with TFT display and 11 speed dérailleur"**

Torque 12/40 Nm, power 250 Watt

Wheel size: 26"

A9100100011-026 € 8,610.00

### Body size

Body height: \_\_\_\_\_ cm

Weight: \_\_\_\_\_ kg

## Chassis colour

Information about current colour designs and example images at: [www.proactiv-gmbh.de/farben.php](http://www.proactiv-gmbh.de/farben.php)

### White RAL 9016

- gloss A9000100170-023  
 matt A9000100170-123 € 40.10

### Brilliant silver

- gloss A9000100170-011

### Anthracite RAL 7016

- gloss A9000100170-001  
 matt A9000100170-101 € 40.10

### Black RAL 9005

- gloss A9000100170-003  
 matt A9000100170-019 € 40.10

### Deep black, silver effect

- gloss A9000100170-007

### Pure orange RAL 2004

- gloss A9000100170-022

### Fire Red

- gloss A9000100170-018  
 matt A9000100170-118 € 40.10

### Pink magenta RAL 4010

- gloss A9000100170-037  
 matt A9000100170-137 € 40.10

### Sky blue

- gloss A9000100170-017

### Pastel blue RAL 5024

- gloss A9000100170-036  
 matt A9000100170-136 € 40.10

### Marine blue

- gloss A9000100170-045  
 matt A9000100170-145 € 40.10

### Electric green

- gloss A9000100170-015  
 matt A9000100170-115 € 40.10

### Lemon green

- gloss A9000100170-005  
 matt A9000100170-105 € 40.10

### Violet

- gloss A9000100170-044

### Powder coating transparent

- gloss A9000100170-029

### Yellow RAL 1023

- gloss A9000100170-046  
 matt A9000100170-146 € 40.10

### Red RAL 3001

- gloss A9000100170-047  
 matt A9000100170-147 € 40.10

### Lagoon blue

- gloss A9000100170-052  
 matt A9000100170-152 € 40.10

### Athletic grey

- gloss A9000100170-050  
 matt A9000100170-150 € 40.10

### Muddy grey

- gloss A9000100170-051  
 matt A9000100170-151 € 40.10

## Neodrives battery

The ranges specified depend on the range profile, user weight and speed.

### Neodrives DT rechargeable battery

**36V 16Ah 576Wh**  
range approx. 75 km

- A9100900109



### Neodrives DT double rechargeable battery

**36V 16Ah 576Wh**  
range approx. 150 km

- A9100900110 € 1,344.40

Quick-change fitting for easy change between batteries using a magnetic plug

## Dérailleur

### Shimano Deore 1x11 Rapid Fire 11-speed dérailleur

- A9110200032  
mechanical dérailleur gear with operation on the handle  
11-speed cassette 11-42 teeth  
top chain ring with chain guards

#### Chain ring selection:

- 32 teeth  
 38 teeth



### Shimano Deore XT 1x11 Rapid Fire 11-speed dérailleur

- A9110200034 € 182.90

Mechanical dérailleur gear with operation on the handle. 11-speed cassette 11-42 teeth

#### Chain ring selection:

- 32 teeth  
 38 teeth

### Shimano Deore XT 1x11 Di2 with button operation 11-speed dérailleur

- A9110200035 € 1,414.50

electronic dérailleur with operation on the handle. 11-speed cassette 11-42 teeth

#### Chain ring selection:

- 32 teeth  
 38 teeth



**Shimano Deore XT 1x11 DI2**  
**With chin operation**  
**11-speed dérailleur**

A9110200033      € 1,552.50

Electronic dérailleur gear with chin operation (tetra equipment).  
11-speed cassette 11-42 teeth

Chain ring selection:

32 teeth

38 teeth



## Gear change accessories

**Mountain drive 2.5: 1 gear reduction for hills for neodrive**

A9116200003      € 690.00

To expand to double number of gears  
Planetary gear set permits switching to the gear reduction for hills when stationary.  
In combination with grip width 40, 45, 48 or 53 cm




## Brake system

**Powerful brake system with brake disc Ø 203 mm**

A9110300105      € 79.10

Recommended from a body weight of 80 kg.  
Not required with hydraulic back pedal disk brake selected



**Double disc brake system "Magura" brake disc Ø 203 mm**

A9100300100      € 222.10

V-brake rim brake omitted  
Only in combination with powerful brake system.  
Recommended from a body weight of 100 kg

**Mechanical backpedal disk brake with crank release function**

A9110300004      € 702.80

Compatible with dérailleur and free-running shifting hub.  
In combination with:  
- pedal bearing shaft grip widths 45 to 55 cm  
- pedal crank crank lengths 155 or 175 mm up to 120 kg body weight  
Not in combination with Mountaindrive 2,5:1 gear reduction for hills!

## Dimensions

### Crank type and grip width









**Standard crank**

- Grip width 40 cm  
A9110500100-040
- Grip width 45 cm  
A9110500100-045
- Grip width 50 cm  
A9110500100-050

**Tetra pedal crank (flared)**

With flared crank arms for max. space for upper body

- Grip width 48 cm (176 mm)  
A9110500100-048
- Grip width 53 cm (226 mm)  
A9110500100-053
- Grip width 58 cm (276 mm)  
A9110500100-058

<p><b>Standard pedal crank</b></p> <ul style="list-style-type: none"> <li><input checked="" type="checkbox"/> Length 155 mm A9110600100-155</li> <li><input checked="" type="checkbox"/> Length 175 mm A9110600100-175</li> <li><input checked="" type="checkbox"/> Length 195 mm A9110600100-195</li> </ul> 	<p><b>Tetra pedal crank</b></p> <ul style="list-style-type: none"> <li><input type="checkbox"/> Length 155 mm A9100600105-155      <b>€ 129.50</b></li> <li><input type="checkbox"/> Length 175 mm A9100600105-175      <b>€ 129.50</b></li> <li><input type="checkbox"/> Individual length 195 mm to 140 mm      <b>€ 246.30</b> A9100600100</li> </ul> 	
<p><b>Fixed bottom bracket support</b></p> <ul style="list-style-type: none"> <li><input checked="" type="checkbox"/> S (X/Y 67/153) A9100501100-003</li> <li><input checked="" type="checkbox"/> M (X/Y 90/250) A9100501100-001</li> </ul> 	<p><b>Adjustable bottom bracket support</b></p> <p>the angle and length can be adjusted via the clamped rotary point on the fork bridge.</p> <p>For the pre-setting, please specify the bottom bracket support size that comes as close to your desired position.</p> <ul style="list-style-type: none"> <li><input type="checkbox"/> S (180 mm)      <b>€ 220.90</b> A9100500001-002</li> <li><input type="checkbox"/> M (240 mm)      <b>€ 220.90</b> A9100500001-003</li> </ul> <p>If no details are provided, the adjustable bottom bracket support shall be supplied in S.</p> 	
<b>Leg rest</b>		
<p><b>Leg rest separated</b></p> <ul style="list-style-type: none"> <li><input checked="" type="checkbox"/> A9100100008 Adaptation on the fork for short legs Only possible up to body size 165 cm</li> </ul> 	<p><b>Leg rest wrapping</b></p> <ul style="list-style-type: none"> <li><input checked="" type="checkbox"/> A9100100009 Adaptation on the fork</li> </ul> 	<p><b>Foot stopper for leg rest</b></p> <ul style="list-style-type: none"> <li><input type="checkbox"/> A9100100012      <b>€ 86.80</b></li> </ul> 
<p><b>Thigh anchorage</b></p> <ul style="list-style-type: none"> <li><input type="checkbox"/> A9100100013      <b>€ 204.70</b></li> </ul> 		

## Back equipment

### Backrest with filling plate and back cushion

A9100800002  
Cover: Black, breathable, open surface

### Body contour backrest upholstery

A9100800010 € 85.10  
with back cushion 3 cm  
Cover: Black, breathable, open surface




## Chassis

### Chassis geometry

Length 116 cm  
Seat width approx.  
39 cm  
A9100100000-005

Length 116 cm  
Seat width approx.  
43 cm  
A9100100000-006 € 194.70

Heavy Duty version  
up to 120 kg body  
weight  
A9100800013 € 405.10  
Recommended for off-road  
use



## Seat upholstery

### Open belt system seat upholstery

Chassis: 116 cm  
SB 39 cm, ST 40 cm

A9100801010-002

Chassis: 116 cm  
SB 43 cm, ST 40 cm

A9100801010-003




### Body Contour upholstery

Chassis: 116 cm  
SB 39 cm, ST 40 cm

A9100801025-002

Chassis: 116 cm  
SB 43 cm, ST 40 cm

A9100801025-003



## Seat cushion

### Standard seat cushion


Cover: Black, breathable, open surface  
Seat cushion height: 4 cm

Chassis: 116 cm  
SB 39 cm, ST 40 cm

A9100000190-002 € 106.60

Chassis: 116 cm  
SB 43 cm, ST 40 cm

A9100000190-003 € 106.60



### Seat cushion with thigh guide

Cover: Black, breathable, open surface  
Seat cushion height: 4 cm

Only in combination with thigh anchorage

Chassis: 116 cm  
SB 39 cm, ST 40 cm

A9100800103-002 € 194.60

Chassis: 116 cm  
SB 43 cm, ST 40 cm

A9100800103-003 € 194.60



## Tyres

### 26x1.5" Schwalbe "Marathon"

A9100105100  
(ETRTO 559-40)  
Touring tyres



### 26x2" Schwalbe "Marathon Plus Tour"

A9100105140 € 67.60  
(ETRTO 559-50)  
Touring tyres



### 26x2.1" Schwalbe "Smart Sam"

A9100105150 € 71.00  
(ETRTO 559-54)



## Bike accessories

### Plug-in protection plate

A9100502100 € 55.30



### Neodrives rapid charger

A9116900006 € 80.20  
Replaces the standard charger



### Collision guard clamped on axle tube

A9100800100 € 207.50



### Neck rest

A9100800150 € 248.40  
height-adjustable  
angle-adjustable



### Hip safety belt

A9100800200 € 112.90



### Bracket ø 22 mm x 42 mm for fixing of accessories (speedometer, lighting or similar)

A9110903271 Item price € 22.10

Colour: Natural anodised

\_\_\_\_\_ Items



### SIGMA lighting set

A9120900019 € 117.90

Consisting of front light and rear lights, chargeable via Micro USB socket, incl. Micro USB charging cable. Lighting set compliant with road traffic regulations



## Options for users with reduced function

### Palm rest on the handle

A9110601200 € 97.80

Padded, bolted, angle-adjustable. Not in combination with other handles.



### Tetra hand support, pivotable, adjustable grip distance

A9110601101 € 161.30




### Tetra hand support, pivotable, adjustable grip width

A9110601102 € 182.40

Very stiff outer support as round tube ø 20 mm with foam rubber cover



<p><b>Tetra shell handles</b></p> <p><input type="checkbox"/> A9110601100      <b>€ 164.10</b></p> 	<p><b>Operating satellite mounted on the handle</b></p> <p>Assembly: Right <input type="checkbox"/> Left <input type="checkbox"/></p> <p>with existing support line (laid along a brake line)</p> <p><input type="checkbox"/> A9116200000      <b>€ 96.30</b></p> <p>without an existing support line (an additional dummy line is required)</p> <p><input type="checkbox"/> A9116200001      <b>€ 136.50</b></p>	
--	---	--

**Notes and remarks:**

**Location, Date**

**Signature End customer/dealer**

The signature field is intended for the consultant and insured person. This confirms that the intended equipment is intended for the insured person.

- 7 - 2023-02-16 / The prices are plus statutory value-added tax and shipping costs     = Included in the basic price     = Position required  
Order sheet NJ1 el. compact bike Rev. 02 UK