



**Quotation date:** \_\_\_\_\_

**Consultation was carried out by**  
**Staff at the medical supplies shop:** \_\_\_\_\_

**Enquiry no.:** \_\_\_\_\_

**Order date:** \_\_\_\_\_

according to GTC PRO ACTIV Reha-Technik GmbH

**Order no.:** \_\_\_\_\_

### Medical suppliers / company address:

**Cust. no.:** \_\_\_\_\_

**Name:** \_\_\_\_\_

**Street:** \_\_\_\_\_

**Postcode / city:** \_\_\_\_\_

**Phone:** \_\_\_\_\_

**Fax:** \_\_\_\_\_

**E-Mail:** \_\_\_\_\_

**Contact:** \_\_\_\_\_

**Different invoicing address (see below)**

**Cust. no.:** \_\_\_\_\_

**Name:** \_\_\_\_\_

**Street:** \_\_\_\_\_

**Postcode / city:** \_\_\_\_\_

**Different shipping address (see below)**

**Cust. no.:** \_\_\_\_\_

**Name:** \_\_\_\_\_

**Street:** \_\_\_\_\_

**Postcode / city:** \_\_\_\_\_

### Customer name / commission:

**ID / Comm. No.:** \_\_\_\_\_

**Surname, forename:** \_\_\_\_\_

**Street:** \_\_\_\_\_

**Postcode / city:** \_\_\_\_\_

**Phone:** \_\_\_\_\_

**E-Mail:** \_\_\_\_\_

**D.o.B.:** \_\_\_\_\_

**Paying organisation:** \_\_\_\_\_

**Indication:** \_\_\_\_\_



## HUSK-E adaptive bike

**powerful Shimano version**

**up to 120 kg load capacity**

- For users with reduced strength and greater impairments and for everyone who needs a powerful drive for a strong helping hand that can be accessed at any time
- HS version for best possible response behaviour
- Suitable for uneven terrain, hills, climbs and challenging routes
- derailleur and hub gear systems possible
- Disk brakes can be implemented with a brake lever or a back-pedal disk brake

### basic equipment

- Shimano middle motor| 85Nm
- Docking plate for quick attachment to the „Quick-Snap adapter“
- adjustable and foldable bottom bracket support to minimize the packing size
- Mechanical parking disc brake
- Parking stand with longer operating lever
- LED lighting set for Shimano system
- Mudguard
- 20" drive wheel
- Small display steps
- Support level switch with ON/OFF button and light button

A9117100031-020

€ 7,390.00

### Body size

Body height: \_\_\_\_\_ cm

Upper body height: \_\_\_\_\_ cm

Weight: \_\_\_\_\_ kg

This information is necessary for checking!

## Heavy Duty version

### Special clearance user weight

A9117100100 € 507.00

for user weights from 120 kg to 150 kg

## Display

### Small display steps

included in the basic price

Features:

- time
- support mode
- battery charge level
- speed
- Distance and travel time
- total kilometers
- Estimated range and range overview
- pedal cadence
- gear display for Di2 gear systems

Only in combination with the green-backgrounded circuits



### Large display steps

A9117200110 € 61.30

replaces the smaller colour display, which is installed as standard.

Features:

- good readable display, large sharp writing
- many display contents on the same screen

Only in combination with the green-backgrounded circuits



### Large display E6100 recommended for Tetra configuration


A9117200100 € 130.00

replaces the smaller colour display EM800, which is installed as standard.

Features:

- on and off button, system start of the display
- light button (on/off at the display)
- mode button
- good readable display, large sharp writing
- many display contents on the same screen

Only in combination with the orange-backgrounded circuits for Tetra configuration



## Transmission

### Chainring on bottom bracket with 28 teeth

A9117200050  
for a higher top speed

### Chainring on bottom bracket with 22 teeth

A9117200051  
for a very light response

## Derailleur gears

### Shimano Deore 1x11 Rapid Fire 11-speed derailleur gear

A9117200026  
mechanical derailleur gear system with control on the right handle.  
11-speed cassette 11-42 teeth and a gear range of 381 % with hydraulic disc brake Ø 203 mm



SHIMANO  
**DEORE**

### Shimano CUES 1x10 Di2 with push-button operation 10-speed Linkglide derailleur gear

A9117200035 € 292.00  
electronic derailleur gear system with control on the handlebar via a button and automatic function.

Operation: right side  left side   
If not specified, the operation will be mounted on the „right side“.

10-speed cassette 11-43 teeth and a gear range of 390 % with hydraulic disc brake Ø 203 mm




SHIMANO  
**CUES**  
Di2

### Nexus 5-speed gear hub Di2 with push-button operation

A9117200036 € 496.00  
electronic hub gear with automatic shift function and manual shift function via control button.

Operation: right side  left side   
If not specified, the operation will be mounted on the „right side“.

with hydraulic disc brake Ø 180 mm.  
with additional mechanical parking disc brake.




## Gear hub for tetra configuration

### Shimano Deore XT 1x11 DI2 with chin control

A9117200023 € 1,224.30


electronic derailleur with chin control (tetra equipment).  
11-speed cassette 11-42 teeth  
and a gear ratio range of 381 %.



### Nexus 5-speed hub gear DI2 with chin control

A9117200033 € 984.40

electronic hub gear with automatic shifting function and manual shifting function via chin control  
button mounted centrally on the bottom bracket



## Brake system

### Double hydraulic disc brake system „Magura“

A9117300100 € 168.10

Two Magura 2-piston callipers  
Parking brake via locking bar on the handle  
Omitting the mechanical parking disk brake.

### Mechanical back-pedal disc brake with crank release function

A9117300110 € 738.40

compatible with derailleur gears and freewheel hub.

Incl. second mechanical disc brake designed as a parking brake.

The hydraulic disc brake is omitted.

In combination with:

- grip width 45 cm, 50 cm, 53 cm or 58 cm
- crank length 155 mm or 175 mm
- up to 120 kg body weight



### Driving brake - supplement for the mechanical parking brake

A9117300120 € 116.10

The parking disc brake is supplemented with a brake handle on the left-hand handle.

This means that the mechanical disc brake can be used as a driving brake and parking brake.

## Batteries and accessories

The specified ranges depend on the distance profile, user weight and speed.

### Battery Shimano 36V 14Ah 504Wh

Range approx. 45 km - 75 km

A9117900150  
mounted on the left of the fork

### Battery Shimano 36V 17.5Ah 630Wh

Range approx. 55 km - 85 km

A9117900155 € 186.60  
mounted on the left of the fork



### Double Battery Shimano 14Ah 504Wh

Range approx. 90 km - 150 km

A9117900160 € 1,467.80

mounted on both sides of the fork,  
interchangeable connector mounted on the fork  
for easy switching between the batteries.  
Bracket swivels to the right for removing the  
battery.

Only in combination with frame adapter with  
right-hand operation.

### Double Battery Shimano 17.5Ah 630Wh

Range approx. 110 km - 170 km

A9117900165 € 1,872.70

mounted on both sides of the fork,  
interchangeable connector mounted on the fork  
for easy switching between the batteries.  
Bracket swivels to the right for removing the  
battery.

Only in combination with frame adapter with  
right-hand operation.



## Colour mudguard

If not specified, the mudguard and chain guard will be delivered in anthracite RAL 7016 matt.

### Mudguard, chain guard and motor cover anthracite RAL 7016 matt

A911750004

### Mudguard, chain guard and motor cover pure orange RAL 2004 gloss

A911750006



### Mudguard, chain guard and motor cover lemon green gloss

A911750007



### Mudguard, chain guard and motor cover customized colour selection

A9117500010 € 153.50

from colour palette one-colour design  
powder colours wheelchair

A9000100170-\_\_\_\_

Colour: \_\_\_\_\_

## Tyres

### 20x2.15" Schwalbe „Marathon 365“

A9117105103  
(ETRTO 406-55)  
profile tyre  
puncture protection level 5



### 20x2.15" Schwalbe „Big Ben Plus“

A9117105102 € 28.70  
(ETRTO 406-55)  
profile tyre  
puncture protection level 6



## Dimensions

### Crank type and grip width

#### Standard crank

- grip width 40 cm  
A9117500100-040
- grip width 45 cm  
A9117500100-045
- grip width 50 cm  
A9117500100-050

#### Tetra crank (flared)

with extended crank arms for maximum upper body space

- grip width 48 cm  
A9117500100-048
- grip width 53 cm  
A9117500100-053
- grip width 58 cm  
A9117500100-058

#### Bottom bracket support

maximum length from the pivot point in the fork to the bottom bracket shaft

- length 520 mm (short)  
using shortened holding tubes € 155.30  
A9117500001  
recommendation: up to approx. 165 cm body height
- length 570 mm (standard)  
A9117500008  
recommendation: from approx. 165 cm to 185 cm body height
- length 620 mm (long)  
using a longer bottom bracket distance € 204.50  
A9117500003  
recommendation: from approx. 185 cm body height

<p><b>Standard pedal crank</b></p> <ul style="list-style-type: none"> <li><input type="checkbox"/> length 155 mm A9110600100-155</li> <li><input type="checkbox"/> length 175 mm A9110600100-175</li> <li><input type="checkbox"/> length 195 mm A9110600100-195</li> </ul> 	<p><b>Tetra pedal crank</b></p> <ul style="list-style-type: none"> <li><input type="checkbox"/> length 155 mm A9100600105-155      <b>€ 136.10</b></li> <li><input type="checkbox"/> length 175 mm A9100600105-175      <b>€ 136.10</b></li> <li><input type="checkbox"/> Individual length 140 mm to 195 mm A9100600100      <b>€ 258.70</b></li> </ul> 	
--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------	------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------	--

## Handles

<p><b>Palm rest on the handle</b></p> <ul style="list-style-type: none"> <li><input type="checkbox"/> A9110601200      <b>€ 102.80</b> upholstered, bolted, angle-adjustable.</li> </ul> 	<p><b>Tetra hand support, pivotable, adjustable grip distance</b></p> <ul style="list-style-type: none"> <li><input type="checkbox"/> A9110601101      <b>€ 169.40</b></li> </ul> 	<p><b>Tetra backrest swivel and adjustable grip width</b></p> <ul style="list-style-type: none"> <li><input type="checkbox"/> A9110601102      <b>€ 191.60</b> very stiff external support as round tube Ø 20 mm with foam rubber coating</li> </ul> 
<p><b>Tetra ergonomic handle, freestanding</b></p> <ul style="list-style-type: none"> <li><input type="checkbox"/> A9100600001      <b>€ 164.60</b> crank handles arranged vertically can be combined with all other handle options</li> </ul> 	<p><b>Tetra shell handles</b></p> <ul style="list-style-type: none"> <li><input type="checkbox"/> A9110601100      <b>€ 172.40</b></li> </ul> 	

## Accessories Bike

<p><b>Maneuvering castors for parking stand</b></p> <ul style="list-style-type: none"> <li><input type="checkbox"/> A9110900017      <b>€ 122.10</b> Maneuvering castors mounted on the parking stand via inlet</li> </ul> 	<p><b>Luggage bags including mounts for the fork tubes</b></p> <ul style="list-style-type: none"> <li><input type="checkbox"/> A9117900220      <b>€ 315.40</b> two panniers, each with a volume of 5.8 liters Colour of mount: silver Colour of bags: dark sand Max. load per side: 5 kg</li> </ul> 	<p><b>Operating satellite on the handle</b></p> <p>Mounting: Right side <input type="checkbox"/> Left side <input type="checkbox"/></p> <p>with existing support line (laid along a brake line)</p> <ul style="list-style-type: none"> <li><input type="checkbox"/> A9117200040      <b>€ 89.50</b></li> </ul> <p>without existing support line (an additional blind line is required)</p> <ul style="list-style-type: none"> <li><input type="checkbox"/> A9117200041      <b>€ 126.90</b></li> </ul>
--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------	-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------	--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

## Options for users with reduced functions

### Lever extension excenter lock for Quick Snap connection

A9110700006 € 73.90

simplifies the opening and closing of the Quick Snap adapter  
Tetra version, curved



### Adaptation help for centering pin

A9110700002 € 55.00  
to simplify the insertion of the central adapter.



### Chest strap for upper body fixation

- size S (105 cm) A9000800016 € 241.80
- size M (125 cm) A9000800013 € 241.80
- size L (145 cm) A9000800014 € 241.80
- size XL (165 cm) A9000800017 € 241.80

consisting of elastic webbing material in combination with neoprene with Velcro and operating tab.  
Colour: black



## Frame adapter

### SPEEDY rigid frame adapter right-hand operation with „Quick-Snap-Adaption“

A9110700142 € 537.20

with quick-release clamp and centering pin holder  
Colour: black

Please select accordingly:

**Centring pin for Ø 25 mm**  
for A1 and A2

A9110700022

**Centring pin for Ø 40 mm**  
for 4all, 4all Ergo, 4you Ergo, F2, F3 and 4teen

A9110700023



### TRAVELER folding frame adapter right-hand operation with „Quick-Snap-Adaption“

A9110701302 € 537.20

with quick-release clamp and centering pin holder on the folding mechanism  
Colour: black



### BUDDY rigid frame adapter right-hand operation with „Quick-Snap-Adaption“

A9110701401 € 537.20

with quick-release clamp and centering pin mount on the frame  
Colour: black



**Wheelchair front adapter  
right-hand operation  
with universal clamp  
and „Quick-Snap-Adaption“**

A9110701502 € 572.30

Adapter attachment to the wheelchair stem.  
Operation on the right-hand side.

Mounting situation for PA products:  
Straight front end from SB 31 cm to 54 cm  
V-shaped front end from SB 30 cm to 50 cm

**Please select the frame clamp below:**

clamping range 95 mm  
A9110701525  
Tube diameter 23-27 mm

clamping range 95 mm  
A9110701530  
Tube diameter 28-32 mm

clamping range 95 mm  
A9110701535  
Tube diameter 33-37 mm

clamping range 70 mm € 76.50  
A9110900072  
Tube diameter 23-27 mm

clamping range 70 mm € 76.50  
A9110900078  
Tube diameter 28-32 mm

clamping range 70 mm € 76.50  
A9110900079  
Tube diameter 33-37 mm

clamping range 75 mm € 268.10  
A9199000080  
Tube diameter Kuschall 2.0

**PLEASE NOTE:**

If the exact tube diameter is not specified, an adapter for a tube diameter of 23 - 27 mm and a clamping range of 95 mm will be delivered as standard.

**Please select the clamping sleeve below:**

clamping sleeve [80 mm]  
included as standard  
for close clearance between footrest and drive wheel.

clamping sleeve [110 mm] € 46.70  
A9199000117  
for medium distance between footrest and drive wheel.

clamping sleeve [170 mm] € 56.60  
A9199000154  
for a wide distance between footrest and drive wheel.

**PLEASE NOTE:**

If not specified, the clamping sleeve with 80 mm length will be delivered.

**Clamp lever with 3-way star handle for actuating the frame clamping**

A9199701507 € 106.10

mounted on the right and left for easy terminal operation with limited manual function  
Only in combination with front adapter



**Wheelchair-dependent accessories**

**Extended bike plate**

for SPEEDY A1 / A2 € 370.00  
A9110100100-002

for BUDDY € 370.00  
A9110100100-003

for TRAVELER € 370.00  
A9110100100-001  
4° camber integrated  
Track independently adjustable



**Wheelbase extension via „grid mechanism with pulling operation“**

for SPEEDY 4all / 4teen / 4you / TRAVELER 4you € 564.70  
A9110102001

removable, with continuous axle tube, track width independently adjustable

**Please select the wheel camber:**

- 1° drive wheel camber
- 4° drive wheel camber
- 6° drive wheel camber



**Wheelbase extension angle-adjustable and removable**

for SPEEDY F2 / F3 € 746.80  
A9015100005

with continuous axle tube, height adjustment via gear joint adjustable in 6° angular increments wheel position displacement approx. 14 cm adjustable track

Attachment possibility must be checked by PRO ACTIV.

Only in combination with base plate. Up to 120 kg load capacity.

**Please select the wheel camber below:**

- 1° drive wheel camber
- 4° drive wheel camber
- 6° drive wheel camber



## Services

### Consultation, clarification of the installation situation and measurement through PRO ACTIV

A9990000000 € 354.50

### Handover and instruction of the handbike / traction device by PRO ACTIV

A9990000011 € 354.50

### Attachment of the handbike / traction device to the wheelchair by PRO ACTIV

#### at specialist dealer

#### on a PRO ACTIV wheelchair

A9990000012 € 230.00

product no.: \_\_\_\_\_

#### on a third-party product

A9990000013 € 366.90

manufacturer: \_\_\_\_\_

model: \_\_\_\_\_

frame tube Ø \_\_\_\_ mm

### Attachment of the handbike / traction device to the wheelchair through specialist dealer

including mounting of the frame adapter

#### on a PRO ACTIV wheelchair

A9990000008

product no.: \_\_\_\_\_

#### on a third-party product

A9990000005

manufacturer: \_\_\_\_\_

model: \_\_\_\_\_

frame tube Ø \_\_\_\_ mm

### Notes and remarks:

### Location, Date

The signature field is intended for the consultant and insured person. This confirms that the intended equipment is intended for the insured person.  
- 8 - 2026-03-03 / The prices are plus statutory value-added tax and shipping costs

### Signature End customer/dealer

= included in the basic price  = Position selected