

Im Hofstätt 11
 D-72359 Dotternhausen
 Tel +49 (0)7427 9480-0
 Fax +49 (0)7427 9480-7025
 e-mail info@proactiv-gmbh.de
 www.proactiv-gmbh.de



BUDDY

Active easy running rigid frame wheelchair

Maximum user weight limit 100 kg

A9041130001 € 4,080.00

Quotation date: _____
Consultation was carried out by
Staff at the medical supplies shop: _____
Enquiry no.: _____

Order date: _____
 according to GTC PRO ACTIV Reha-Technik GmbH

Order no.: _____

Medical suppliers / company address:

Cust. no.: _____
Name: _____
Street: _____
Postcode / city: _____
Phone: _____
Fax: _____
E-Mail: _____
Contact: _____

Shipping address: same as company address
 Other (see below)

Name: _____
Street: _____
Postcode / city: _____


Customer name / commission:

ID / Comm. No.: _____
Surname, forename: _____
Street: _____
Postcode / city: _____
Phone: _____
Fax: _____
E-Mail: _____
D.o.B.: _____
Paying organisation: _____

Invoicing address: same as company address
 Other (see below)

Cust. no.: _____
Name: _____
Street: _____
Postcode / city: _____
Phone: _____
Fax: _____
E-Mail: _____
Contact: _____

Body size:




Body height: _____ cm
Weight: _____ kg
Lower leg length: _____ cm
 (anatomic dimension: hollow of the knee to the foot sole)

Indication: _____

This information is necessary for checking. If the information is missing, a check of the relevant items cannot be made!

Basic size



Seat width **SB**

<input type="checkbox"/> 20	<input type="checkbox"/> 21	<input type="checkbox"/> 22	<input type="checkbox"/> 23	<input type="checkbox"/> 24	<input type="checkbox"/> 25	<input type="checkbox"/> 26	<input type="checkbox"/> 27	<input type="checkbox"/> 28	<input type="checkbox"/> 29
<input type="checkbox"/> 30	<input type="checkbox"/> 31	<input type="checkbox"/> 32	<input type="checkbox"/> 33	<input type="checkbox"/> 34	<input type="checkbox"/> 35	<input type="checkbox"/> 36	<input type="checkbox"/> 37	<input type="checkbox"/> 38	

Size M (Frame length 36 cm (to the bend) and start of abduction at ST 27 cm)

<input type="checkbox"/> 22	<input type="checkbox"/> 23	<input type="checkbox"/> 24	<input type="checkbox"/> 25	<input type="checkbox"/> 26	<input type="checkbox"/> 27	<input type="checkbox"/> 28	<input type="checkbox"/> 29	<input type="checkbox"/> 30	<input type="checkbox"/> 31
-----------------------------	-----------------------------	-----------------------------	-----------------------------	-----------------------------	-----------------------------	-----------------------------	-----------------------------	-----------------------------	-----------------------------

Seat depth **ST**

Size L (Frame length 39 cm (to the bend) and start of abduction at ST 30 cm)

<input type="checkbox"/> 26	<input type="checkbox"/> 27	<input type="checkbox"/> 28	<input type="checkbox"/> 29	<input type="checkbox"/> 30	<input type="checkbox"/> 31	<input type="checkbox"/> 32	<input type="checkbox"/> 33	<input type="checkbox"/> 34	<input type="checkbox"/> 35
-----------------------------	-----------------------------	-----------------------------	-----------------------------	-----------------------------	-----------------------------	-----------------------------	-----------------------------	-----------------------------	-----------------------------

Size XL (Frame length 43 cm (to the bend) and start of abduction at ST 34 cm)

<input type="checkbox"/> 31	<input type="checkbox"/> 32	<input type="checkbox"/> 33	<input type="checkbox"/> 34	<input type="checkbox"/> 35	<input type="checkbox"/> 36	<input type="checkbox"/> 37	<input type="checkbox"/> 38	<input type="checkbox"/> 39	<input type="checkbox"/> 40
-----------------------------	-----------------------------	-----------------------------	-----------------------------	-----------------------------	-----------------------------	-----------------------------	-----------------------------	-----------------------------	-----------------------------

Back height **RH**

<input type="checkbox"/> 20	<input type="checkbox"/> 21	<input type="checkbox"/> 22	<input type="checkbox"/> 23	<input type="checkbox"/> 24	<input type="checkbox"/> 25	<input type="checkbox"/> 26	<input type="checkbox"/> 27	<input type="checkbox"/> 28	<input type="checkbox"/> 29
<input type="checkbox"/> 30	<input type="checkbox"/> 31	<input type="checkbox"/> 32	<input type="checkbox"/> 33	<input type="checkbox"/> 34	<input type="checkbox"/> 35	<input type="checkbox"/> 36	<input type="checkbox"/> 37	<input type="checkbox"/> 38	<input type="checkbox"/> 39
<input type="checkbox"/> 40	<input type="checkbox"/> 41	<input type="checkbox"/> 42	<input type="checkbox"/> 43	<input type="checkbox"/> 44	<input type="checkbox"/> 45	<input type="checkbox"/> 46	<input type="checkbox"/> 47	<input type="checkbox"/> 48	

Lower leg length **USL**

<input type="checkbox"/> 20	<input type="checkbox"/> 21	<input type="checkbox"/> 22	<input type="checkbox"/> 23	<input type="checkbox"/> 24	<input type="checkbox"/> 25	<input type="checkbox"/> 26	<input type="checkbox"/> 27	<input type="checkbox"/> 28	<input type="checkbox"/> 29
<input type="checkbox"/> 30	<input type="checkbox"/> 31	<input type="checkbox"/> 32	<input type="checkbox"/> 33	<input type="checkbox"/> 34	<input type="checkbox"/> 35	<input type="checkbox"/> 36	<input type="checkbox"/> 37	<input type="checkbox"/> 38	<input type="checkbox"/> 39
<input type="checkbox"/> 40	<input type="checkbox"/> 41	<input type="checkbox"/> 42	<input type="checkbox"/> 43	<input type="checkbox"/> 44	<input type="checkbox"/> 45	<input type="checkbox"/> 46	<input type="checkbox"/> 47		

Seat height back **SHh**

with 20" drive wheels						
<input type="checkbox"/> 31	<input type="checkbox"/> 32	<input type="checkbox"/> 33	<input type="checkbox"/> 34	<input type="checkbox"/> 35	<input type="checkbox"/> 36	<input type="checkbox"/> 37
with 22" drive wheels						
<input type="checkbox"/> 33.5	<input type="checkbox"/> 34.5	<input type="checkbox"/> 35.5	<input type="checkbox"/> 36.5	<input type="checkbox"/> 37.5	<input type="checkbox"/> 38.5	<input type="checkbox"/> 39.5
with 24" drive wheels						
<input type="checkbox"/> 36	<input type="checkbox"/> 37	<input type="checkbox"/> 38	<input type="checkbox"/> 39	<input type="checkbox"/> 40	<input type="checkbox"/> 41	<input type="checkbox"/> 42

Seat height front **SHv**
(with 5" caster wheel)

<input type="checkbox"/> 40	<input type="checkbox"/> 41	<input type="checkbox"/> 42	<input type="checkbox"/> 43	<input type="checkbox"/> 44	<input type="checkbox"/> 45	<input type="checkbox"/> 46	<input type="checkbox"/> 47
-----------------------------	-----------------------------	-----------------------------	-----------------------------	-----------------------------	-----------------------------	-----------------------------	-----------------------------

Tipping point **KP**
(inner edge of rear tube to
centre of drive wheel axle)

<input type="checkbox"/> <u>5</u>	<input type="checkbox"/> 6	<input type="checkbox"/> 7	<input type="checkbox"/> 8	<input type="checkbox"/> 9	<input type="checkbox"/> 10
-----------------------------------	----------------------------	----------------------------	----------------------------	----------------------------	-----------------------------

Please tick the required dimensions (all dimensions are in cm)

The underlined values are the standard settings. If not specified, the standard value will be delivered.

The SHh and SHv values might be a bit different from the specified values,
depending on the wheel size and camber.

We recommend a seat height difference (difference between SHv and SHh) between 5 cm and 12 cm.

Frame colour - One-colour design

Information on current colour designs and example images can be found at: <https://www.proactiv-gmbh.com/colours/>

Frame colour: Selection from colour palette

Frame colour will be specified in the order

A9000100170-000

Raw, blasted, not coated

A9000100170-020

Note: In combination with moisture e.g. hand perspiration, can affect dark aluminum discoloration.

Powder colours

white RAL 9016

gloss A9000100170-023

matt A9000100170-123 € 43.40

brilliant silver

gloss A9000100170-011

anthracite RAL 7016

gloss A9000100170-001

matt A9000100170-101 € 43.40

black RAL 9005

gloss A9000100170-003

matt A9000100170-019 € 43.40

pure orange RAL 2004

gloss A9000100170-022

matt A9000100170-122 € 43.40

pink magenta RAL 4010

gloss A9000100170-037

matt A9000100170-137 € 43.40

pastel blue RAL 5024

gloss A9000100170-036

matt A9000100170-136 € 43.40

lemon green

gloss A9000100170-005

matt A9000100170-105 € 43.40

yellow RAL 1023

gloss A9000100170-046

matt A9000100170-146 € 43.40

red RAL 3001

gloss A9000100170-047

matt A9000100170-147 € 43.40

lagoon blue

gloss A9000100170-052

matt A9000100170-152 € 43.40

athletic grey RAL 7046

gloss A9000100170-050

matt A9000100170-150 € 43.40

muddy grey RAL 7033

gloss A9000100170-051

matt A9000100170-151 € 43.40

peach RAL 3022

gloss A9000100170-038

matt A9000100170-138 € 43.40

Powder colours „dormant“

deep black effect silver

gloss A9000100170-007

sky blue

gloss A9000100170-017

navy blue

gloss A9000100170-045

matt A9000100170-145 € 43.40

electric green

gloss A9000100170-015

matt A9000100170-115 € 43.40

violet

gloss A9000100170-044

fire red

gloss A9000100170-018

matt A9000100170-118 € 43.40

Frame colour - One-colour design special colours „liquid coating“

€ 349.60

storm (ice blue)

gloss A9000100170-500

matt A9000100170-501

prairie (rose gold)

gloss A9000100170-504

matt A9000100170-505

berry

gloss A9000100170-526

matt A9000100170-527

blue ocean

gloss A9000100170-522

matt A9000100170-523

emerald

gloss A9000100170-514

matt A9000100170-515

chili

gloss A9000100170-516

matt A9000100170-517

Frame colour - Two-colour design

Note: In the two-colour design, both colours are the same gloss level as the design colour.

Two-colour design with one base colour (on the rear frame) and one design colour (on the front frame) € 268.50

Please select a colour below:

Base colour (on the rear frame)

white RAL 9016 <input type="checkbox"/> A9000100170-200	brilliant silver <input type="checkbox"/> A9000100170-201	anthracite RAL 7016 <input type="checkbox"/> A9000100170-202
black RAL 9005 <input type="checkbox"/> A9000100170-203	pastel blue RAL 5024 <input type="checkbox"/> A9000100170-206	lemon green <input type="checkbox"/> A9000100170-207
yellow RAL 1023 <input type="checkbox"/> A9000100170-208	lagoon blue <input type="checkbox"/> A9000100170-212	athletic grey RAL 7046 <input type="checkbox"/> A9000100170-210
muddy grey RAL 7033 <input type="checkbox"/> A9000100170-211		

Design colour (on the front frame)

black RAL 9005 <input type="checkbox"/> gloss A9000100170-259 <input type="checkbox"/> matt A9000100170-260	lemon green <input type="checkbox"/> gloss A9000100170-250 <input type="checkbox"/> matt A9000100170-251	chili <input type="checkbox"/> gloss A9000100170-253 <input type="checkbox"/> matt A9000100170-254
pure orange RAL 2004 <input type="checkbox"/> gloss A9000100170-255 <input type="checkbox"/> matt A9000100170-256	emerald <input type="checkbox"/> gloss A9000100170-274 <input type="checkbox"/> matt A9000100170-275	storm (ice blue) <input type="checkbox"/> gloss A9000100170-264 <input type="checkbox"/> matt A9000100170-265
prairie (rose gold) <input type="checkbox"/> gloss A9000100170-266 <input type="checkbox"/> matt A9000100170-267	berry <input type="checkbox"/> gloss A9000100170-284 <input type="checkbox"/> matt A9000100170-285	blue ocean <input type="checkbox"/> gloss A9000100170-280 <input type="checkbox"/> matt A9000100170-281

Colour - Attachment parts

Note: If not specified, the caster fork and bearing block will be mounted in black.

Drive wheel: aluminum rim, hub and spokes

Black A9000100171-003

Caster fork and bearing block

Silver A9000100173-011

Black A9000100173-003

Colour design of attachment parts according to customer request:

Drive wheel: aluminum rim and design inlet¹ according to colour selection / hub and spokes black

frame colour A9000100171-999 € 249.60

Selection from one-colour design (different from the frame colour) A9000100171-998 € 288.90

Colour: _____

A9000100170-

Caster wheel: 5" aluminum rim² (6-spoke design)

frame colour A9000100174-999 € 75.20

Selection from one-colour design (different from the frame colour) A9000100174-998 € 87.10

Colour: _____

A9000100170-

Special colour request (Frame colour - One-colour design)

Frame colour for single purchasing

Selection from RAL colour palette

A9000100170-100 € 414.40

gloss level

gloss

matt

colour designation

RAL number

Note: Attachments **not** available in design colour options

¹ Only in combination with Liteline hub

² Only in combination with caster wheel 5" aluminum rim 6-spoke-design good shock absorption A900050121*, as well as with suspension fork 5-spoke-design A9000500160*

Frame

Front end abducted shape

A9041130030-K*
Knee angle 105°
Bearing block firmly welded



Push handles

Equipment without push handles

A9000401000



Safety push handles back-positioned

Length: 45 cm
A9041401075-001
max. height = SHh + approx. 53 cm

Length: 55 cm
A9041401075-002
max. height = SHh + approx. 63 cm
continuous height-adjustable
push handle holder back-positioned



Aluminum push handles, back-positioned with side handle offset 145 mm

A9041401080-K* € 89.50

continuous height-adjustable
Safety push handle holder mounted on the back cross bar
Length: 55 cm
Up to SB 38 cm or body weight 60 kg



Aluminium central push handle Foldable handle area

A9041401086-K € 102.60

height-adjustable and foldable in the handle area via toothed joint, holder welded on the middle on the back cross bar.
Up to SB 38 cm or body weight 50 kg
Not in combination with headrest mount
Please note: Not suitable for carrying the wheelchair — use the back crossbar for this purpose.



Central push handle for brake handles Foldable grip area

A9041401070-K € 141.90

closed, with continuous height adjustment, with foam rubber handle cover
Up to SB 38 cm or body weight 50 kg
Not in combination with headrest mount
Only in combination with parking lever brake for accompanying person or drum brake.



Aluminium push handle closed, foldable via toothed joint

A9041401069-K € 389.90

ontinuous height-adjustable
push handles shifted to the back
Up to SB 38 cm or body weight 60 kg



Backrest

Fixed backrest

A9041420020
screwed to aluminium side parts
Back angle 90° to the floor
Backrest cannot fold down
Side parts in frame colour



Folding backrest

A9000400023 € 306.30

7-way angle adjustable,
with integrated plastic-coated aluminum side parts



Back equipment

If not specified, the „Adaptable back with back padding“ will be delivered.

Adaptable back with back padding

A9000802010-K001

Cover: Black, breathable, open surface
Equipment with back cross bar



Adaptable back with back padding „light“

A9000802320-K001

Cover: Black, breathable, open surface
Cushion thickness: 9 mm
Equipment with back cross bar

A9000802320-K002

Cover: Black, breathable, open surface
Cushion thickness: 9 mm
Equipment without back cross bar



Adaptable back with back padding „Comfort“

A9000815101-K € 131.40

Cover: Black, breathable, open surface
Front padding: 20 mm
Rear padding: 5 mm
Equipment with back cross bar

Ergonomic back shell aluminum from RH 24 cm

A9000802210-K € 441.60

with lumbar support.
Colour: Anthracite matt
(If no colour is specified under „Select colour“)
If required, please select the back cushion separately.



Back cushion for ergonomic back shell

A9000802220-K € 221.30


cushion thickness: 2 cm

Cover:

Black, breathable, open surface

Black, breathable, closed surface

If not specified, the cover material „Black, breathable, open surface“ will be delivered.



ERGO Back backrest bar for ergonomic back shell from RH 24 cm and SB 30 cm

A9000802019-K € 242.60

with recessed handle area
Colour: Anthracite matt
(If no colour is specified under „Select colour“)
Only compatible with central push handle
Item no. A9041401082
If required, please select the back shell and the back cushion separately.



Equipment without back upholstery and back padding or back system

A9000802011 € - 53.00

Equipment with back cross bar



Back crossbar flared

by 20 mm A9000400835-020 € 109.30

by 30 mm A9000400835-030 € 109.30

by 50 mm A9000400835-050 € 140.80

by 80 mm A9000400835-080 € 155.30

Only in combination with back-positioned push handles




Seat upholstery

If not specified, the „Body Contour“ seat upholstery will be delivered.

Body Contour

A9000801025-K



Seat plate aluminum

A9000801212-K € 292.70

screwed onto the frame tubes
Colour: anthracite matt
(If no colour is specified under „Select colour“)



Equipment without seat upholstery

A9000801050

Colour Back equipment and Seat plate

Information on current colour designs and availability at www.proactiv-gmbh.com

If not specified, the seat and backrest system will be delivered in the standard colour anthracite matt.

anthracite RAL 7016 gloss

- A9000100179-995
- Ergonomic back shell
- ERGO Back backrest
- Seat plate aluminum

in frame colour

- A9000100179-999 **€ 92.30**
- Ergonomic back shell
- ERGO Back backrest
- Seat plate aluminum

For frame coating in two-colour design



Please specify if the desired option should be coated in the base colour or the design colour.

- A9000100179-999 **€ 92.30**
- Ergonomic back shell
- base colour
- design colour
- ERGO Back backrest
- base colour
- design colour
- Seat plate aluminum
- base colour
- design colour

If a different colour is required, please select it from the wheelchair colour palette:

- A9000100179-998 **€ 149.40**
- Please enter the last 3 digits of the desired colour:
A9000100170-
- Ergonomic back shell
- ERGO Back backrest
- Seat plate aluminum

Seat cushion



Seat cushion in rectangular contour Inlet: Standard

- A9000805180-K* **€ 118.30**

Seat cushion with radius contour Inlet: Standard

- A9000800300-K* **€ 163.10**

Seat cushion height:

- 3 cm 5 cm
- 4 cm 6 cm

Seat cushion in rectangular contour Inlet: Anti-Decubitus

- up to 30x30 cm A9000805204-K* **€ 301.70**
- >30x30 up to 50x50 cm A9000805200-K* **€ 437.70**

Seat cushion with radius contour Inlet: Anti-Decubitus

- up to 30x30 cm A9000800311-K* **€ 356.10**
- >30x30 up to 50x50 cm A9000800310-K* **€ 482.90**

Seat cushion height:

- 4 cm 6 cm

Cover:

- Black, breathable, open surface
- Black, breathable, closed surface

If not specified, the cover material „Black, breathable, open surface“ will be delivered.

- Seat cushion is provided**

Clothing guard

Aluminum clothing guard

fixed to the side panel
Only in combination with fixed backrest



in frame colour (base colour for two-colour design)

- A9000420101-K **€ 311.10**

Selection from colour palette „One-colour design“

- A9000420103-K **€ 454.40**

Colour: _____

gloss level: gloss matt

Aluminum clothing guard removable via lock function adjustable following wheel arch Only in combination with folding backrest					
in anthracite matt <input type="checkbox"/> A9000420215-K € 243.20		in frame colour (base colour for two-colour design) <input type="checkbox"/> A9000420216-K € 311.10		in design colour (only possible with frame colour in two-colour design) <input type="checkbox"/> A9000420217-K € 345.10	
Carbon clothing guard removable via lock function adjustable following wheel arch Only in combination with folding backrest 					
<input type="checkbox"/> A9000420013-K € 547.30					
Knee lever brake system, brake lever integrated in the clothing guard see „Brake systems“.					
Armrests					
Armrest, height adjustable, integrated in the clothing guard Including aluminum clothing guard with arm padding 25 cm long Removable via lock function Only in combination with folding backrest					
in anthracite matt <input type="checkbox"/> A9000420014-K € 569.60		in frame colour (base colour for two-colour design) <input type="checkbox"/> A9000420015-K € 621.90		in design colour (only possible with frame colour in two-colour design) <input type="checkbox"/> A9000420016-K € 648.10	
Armrest, height adjustable, integrated in the clothing guard <input type="checkbox"/> A9000420017-K € 862.00 Including carbon clothing guard with arm padding 25 cm long Removable via lock function Only in combination with folding backrest					
Axle plate Standard equipment with 12 mm sockets / axles					
Axle plate <input checked="" type="checkbox"/> A9000220000 quick-release axle Ø 12 mm 		Extended bike plate <input type="checkbox"/> A9000220030 € 315.20 adjustable track 			
Wheel camber *Total width with tight-mounted handrim.					
1° drive wheel camber <input checked="" type="checkbox"/> A9030202010 Total width = approx. SB + 20 cm*		2,5° drive wheel camber <input checked="" type="checkbox"/> A9030202020 Total width = approx. SB + 23 cm*		4° drive wheel camber <input checked="" type="checkbox"/> A9030202030 Total width = approx. SB + 26 cm*	
6° drive wheel camber <input checked="" type="checkbox"/> A9030202040 Total width = approx. SB + 31 cm*		8° drive wheel camber <input type="checkbox"/> A9041202080 € 77.00 Total width = approx. SB + 35 cm*		10° drive wheel camber <input type="checkbox"/> A9041202100 € 136.90 Total width = approx. SB + 39 cm*	

E-Adaptations

Preparation for the „e-motion“ / „e-motion DuoDrive“ mounting

A9000220200 € 502.90

Equipment: e-motion torque arm, 1/2" axle and socket length 50.8 mm (2")
PRO ACTIV and Alber recommend mounting two anti-tipping supports. Please select under „Adaptations“ or order Alber tipping supports directly from Alber.

Preparation for the „e-fix“ mounting

A9000220210 € 502.90

Equipment: e-motion torque arm, 1/2" axle and socket length 50.8 mm (2")
PRO ACTIV and Alber recommend mounting two anti-tipping supports. Please select under „Adaptations“ or order Alber tipping supports directly from Alber, as well as recording BUDDY.

Brake system

If not specified, the „Brake system push to lock with standard brake handle“ will be delivered.

Recommendation: For KP 7 cm or lower, we recommend the „Brake system push to lock with brake handle 20° outwards“.

Knee lever brake with standard brake lever

A9000701010-001



Knee lever brake with long brake lever

A9000701010-002



Knee lever brake system Brake lever angled 20°

standard brake handle A9000700002

long brake handle A9000701010-004



Brake lever with M8 adapter and plastic ball

A9000706001 € 170.50

only in combination with knee lever brake system, standard or long brake handle



Brake system push to lock Brake handle to flip-down

A9000706008 € 158.60

Length: 9 cm



Knee lever brake with one-hand operation

Operation right side A9000707020-001 € 392.30

Operation left side A9000707020-002 € 392.30
via Bowden cable



Knee lever brake system, brake lever integrated in the clothing guard

Brake lever attached to the lower frame tube.

Actuation via cable pull.

Clothing guard cannot be removed.



in schwarz glanz with fixed backrest
in anthracite matt with foldable backrest

A9041720015-K € 643.00

in frame colour (base colour for two-colour design)

A9041720016-K* € 710.30

in design colour (only possible with frame colour in two-colour design)

A9041720017-K* € 744.30

Knee lever brake with lock

Standard brake handle A9000700214 € 158.60

Long brake handle A9000700215 € 158.60

Locking block for brake lever
Prevents braking during pushing by accompanying persons
Not in combination with knee lever brake with one-hand operation, brake lever angled 20° as well as parking knee lever brake for accompanying person



Parking knee lever brake for accompanying person

A9000701020 € 464.10

Push handle operation via Bowden cable, not in combination with 1° wheel camber



Drum brake

A9000720011-K* € 474.00

Operating from the push handle
Total wheelchair width (see selection „Wheel camber“) + 3 cm

Please note:

- Only in combination with drive wheel aluminum 36-spoke rim wheel
- not in combination with one-hand drive, e-fix or e-motion
- not in combination with 1° wheel camber



Footrest

If not specified, the „continuous footrest” will be delivered.

Continuous footrest

A9041601010-K*
with tubular resting bracket
height adjustable
possible from USL 26 cm to USL 47 cm
If you want to equip the footrest with a „rubber coated footplate” please select this under „Footrest accessories”.



High-distanced footrest

A9041606010-K
height adjustable
Possible from USL 20 cm to 25 cm.
Rubber coated footplate necessary, please select this under „Footrest accessories”.



Footrest with variable positioning

A9047600001-K € 281.10
adjustable in height, depth and angle
including footplate cover with anti-slip film
up to USL 37 cm



Footrest folds up to the rear with spring locking mechanism

A9047602005-K* € 347.20
with tubular resting bracket
up to USL 32 cm
including steering wheel track width extension
20 mm per side
If you want to equip the footrest with a „rubber coated footplate” please select this under „Footrest accessories”.



Footrest with variable positioning and folding up to the rear

A9047600002-K € 373.20
adjustable in height, depth and angle
including footplate cover with anti-slip film
up to USL 37 cm




Footrest folding up to one side

A9041601300-K* € 173.90
from USL 25 cm:
with tubular resting bracket
height adjustable +/- 3 cm
cannot be removed
from USL 29 cm:
can also be swivelled outwards
If you want to equip the footrest with a „rubber coated footplate” please select this under „Footrest accessories”.



Divided footrest

A9041603200-K € 336.80
with rubber coated footplate,
height adjustable
can be pivoted outwards and folded up to the
side
cannot be removed
possible from USL 29 cm and SB 26 cm



Footrest accessories

Rubber coated footplate with 14 cm depth

A9000610000-K € 125.90
bolted, angle-adjustable



Rubber coated footplate with 18 cm depth

A9000610005-K € 140.10
bolted, angle-adjustable



Footplate side stop

A9000610100 € 59.70
adjustable
Only in combination with rubber coated
footplate



Deflector rollers mounted on the footplate

A9000600006 € 117.40
Only in combination with rubber coated
footplate



Adaptations

Anti-tipping support

- right side
A9041301160-K001 € 199.00
- left side
A9041301160-K002 € 199.00
- with tension spring
can be folded up
removable via Quick-pin



Outdoor wheel attachment with steering fork

- A9041900100 € 845.70
- equipped with straight-running lock,
can be removed from the wheelchair via quick
release lever.
With swivelling pneumatic tyre 12.1/2 x 2.1/4"
caster wheel that can be removed from the
support tube via thru axle



3" caster wheel

3" skate wheel (80 mm)

- A9000501110 € 70.80
- puncture-proof



4" caster wheel

4" caster wheel „Standard“ solid rubber black

- A9000501011-004
(82A), puncture-proof



4" skate wheel (100 mm)

- A9000501120 € 70.80
- puncture-proof



4" safety light-up wheel (100 mm)

- A9000501220 € 94.80
- puncture-proof



4" CorMORun caster wheel (100 mm)

- A9000501042 € 173.00
- with replaceable rubber 6 mm,
puncture-proof
Rim colour: Black
Incl. hair guard spacer (5.3 mm)



4" Frog Legs caster wheel (100 mm)

- A9000501230 € 206.30
- with replaceable rubber 6 mm,
puncture-proof
Colour: black, silver rim
Incl. hair guard spacer (5.3 mm)



5" caster wheel

5" caster wheel „Standard“ PU black

- A9000501036-005
(82A), puncture-proof



5" caster wheel „Standard“ PU grey

- A9000501031-005
(82A), puncture-proof



5" safety light-up wheel (125 mm)

- A9000501221 € 94.80
- puncture-proof



<p>5" aluminum rim 3-spoke design „good shock absorption“</p> <p><input type="checkbox"/> silver anodised A9000501210 € 179.80</p> <p><input type="checkbox"/> black A9000501211 € 179.80</p> <p>with replaceable rubber 6 mm medium hard, puncture-proof Incl. hair guard spacer (3 mm)</p> 	<p>5" CorMORun caster wheel (125 mm)</p> <p><input type="checkbox"/> A9000501041 € 181.30</p> <p>with replaceable rubber 6 mm, puncture-proof Rim colour: Black Incl. hair guard spacer (5.3 mm)</p> 	<p>5" Frog Legs caster wheel (125 mm)</p> <p><input type="checkbox"/> A9000501232 € 206.30</p> <p>with replaceable rubber 6 mm, puncture-proof Colour: black, silver rim Incl. hair guard spacer (5.3 mm)</p> 
--	---	--

6" caster wheel

<p>6" caster wheel „Standard“ PU grey</p> <p><input checked="" type="checkbox"/> A9000501010-006 (82A), puncture-proof</p> 	<p>6" caster wheel „pneumatic tyre“</p> <p><input type="checkbox"/> A9000501050 € 57.20</p> <p>Front seat height min. 42 cm</p> 	
--	--	--

Caster fork

For caster fork colours, see „Colour of add-on parts“
If not specified, the „Standard caster fork with screw axle“ will be delivered.

<p>Caster fork Standard</p> <p><input checked="" type="checkbox"/> A9000500100-K001 with screwed axle</p> <p><input type="checkbox"/> A9000500100-K002 € 167.40 with quick-release axle</p> 	<p>Caster fork Heavy Duty for high loads or with the use of additional electric drives</p> <p><input type="checkbox"/> A9000510111-K001 € 54.00 with screwed axle</p> <p><input type="checkbox"/> A9000510111-K002 € 218.90 with quick-release axle</p> 	<p>Suspension fork inkl. 5" aluminum rim „good shock absorption“</p> <p><input type="checkbox"/> A9000500160-K002 € 616.90 with quick-release axle</p> 
--	--	---

<p>Caster fork spare wheel „Standard“</p> <p><input type="checkbox"/> A9000500300-K € 296.70</p> <p>Caster fork Standard with quick-release axles</p> 	<p>Caster fork spare wheel „Heavy Duty“</p> <p><input type="checkbox"/> A9000500301-K € 337.90</p> <p>Caster fork Heavy Duty with quick-release axles</p> 	<p>The spare caster fork does not have a caster wheel. In case a caster wheel is requested, please put a cross on a second selection in the category „Caster wheels“.</p>
---	--	---

Track distance

Caster wheel track extension



A9000500400*
10 mm per side

€ 59.10

A9000500401*
15 mm per side

€ 59.10

A9000500410* € € 59.10
20 mm per side
included in „Footrest folds up to the rear with
spring locking mechanism“
item no. A9047602005-K*

Drive wheels

If not specified, the „drive wheel 36 spoke aluminum rim“ will be delivered.

Drive wheel 18 spoke aluminum rim

20" (ETRTO 451-19)
A9000203052

22" (ETRTO 489-19)
A9000203062

24" (ETRTO 540-19)
A9000203122

with Liteline hub, radially spoked

not in combination with:

- Puncture-proof tires
- One-hand drive
- Drum brake
- Transport restraint connection
- Handbike

up to 70 kg body weight:

Gewicht	Sturz				
	1°	4°	6°	8°	10°
bis 50 Kg	✓	✓	✓	✓	✓
51 - 60 Kg	✓	✓	✓	✓	
61 - 70 Kg	✓	✓			



Drive wheel 36 spoke aluminum rim

22" (ETRTO 489-19)
A9000203060-010

24" (ETRTO 540)
A9000203000-010
with Liteline hub (left image),
radially spoked

20" (ETRTO 451-19)
A9000203050

22" (ETRTO 451-19)
A9000203060-020

24" (ETRTO 540)
A9000203000-020
with aluminum front hub (right image),
radially spoked
in combination with one-hand drive



Drive wheel double rim 36-spokes

24" (ETRTO 540)
A9000203100 € € 157.80

with aluminum front hub,
radially spoked
for use with high loads



Spinergy wheel „Spox Sport“ with 24 PBO spokes; rim, hub and spokes in black

24" (ETRTO 540)
A9000203550 € € 784.10



Spinergy „Everyday Light Extreme“ with 12 PBO spokes; rim and hub in black, spokes in white

24" (ETRTO 540)
A9000203554 € € 957.20

Recommendation: For heavy loads
(e.g. higher body weight), we recommend the
Spinergy „Spox Sport“ wheel with 24 PBO
spokes.



Spinergy „Everyday Light Extreme“ with 12 PBO spokes; rim, hub and spokes in black


24" (ETRTO 540)
A9000203552 € € 957.20

Recommendation: For heavy loads
(e.g. higher body weight), we recommend the
Spinergy „Spox Sport“ wheel with 24 PBO
spokes.



Handrims

x = dimension handrim attachment



Handrim attachment
general:

- x = 11 mm
- x = 19 mm

Handrim attachment
Curve + Gekko handrim:

- x = 10 mm
- x = 18 mm

If not specified, the handrim is mounted tightly (11 mm / 10 mm).

Stainless steel Ø 20 mm

- 20" handrim for 20" wheel
A9000205010-001
- 22" handrim for 22" wheel
A9000205010-002
- 24" handrim for 24" wheel
A9000205010-004



Aluminum Ø 19 mm silver anodised

- 22" handrim for 22" wheel
A9000205019-002
- 24" handrim for 24" wheel
A9000205019-004



Aluminum Ø 16 mm powder coated in frame colour

(basic colour with two-colour design)

- = Gloss surface
- = Matt surface

Note: Matt coating is not recommended due to the low grip.

- 20" handrim for 20" wheel
A9000205021-001 **€ 68.30**
- 22" handrim for 22" wheel
A9000205021-002 **€ 68.30**



Aluminum Ø 19 mm powder coated in frame colour

(basic colour with two-colour design)

- = Gloss surface
- = Matt surface

Note: Matt coating is not recommended due to the low grip.

- 20" handrim for 20" wheel **€ 68.30**
A9000205020-001
- 22" handrim for 22" wheel **€ 68.30**
A9000205020-002
- 24" handrim for 24" wheel **€ 68.30**
A9000205020-004



Aluminum Ø 16 mm polished and hard-anodised

- 20" handrim for 20" wheel **€ 164.90**
A9000205031-001
- 22" handrim for 22" wheel **€ 164.90**
A9000205031-002



Aluminum Ø 19 mm polished and hard-anodised

- 20" handrim for 20" wheel **€ 164.90**
A9000205030-001
- 22" handrim for 22" wheel **€ 164.90**
A9000205030-002
- 24" handrim for 24" wheel **€ 164.90**
A9000205030-004



Titanium Ø 19 mm

- 24" handrim for 24" wheel **€ 335.60**
A9000205040-004



Silicone handrim cover

- for handrim 20" **€ 232.80**
A9000205100-001
 - for handrim 22" **€ 232.80**
A9000205100-002
 - for handrim 24" **€ 232.80**
A9000205100-004
- Surface without structure, compatible with handrim Ø 19 mm or Ø 20 mm.



Ultra-Grip „very grippy“

- 20" handrim for 20" wheel **€ 402.40**
A9000205070-001
- 22" handrim for 22" wheel **€ 402.40**
A9000205070-002
- 24" handrim for 24" wheel **€ 402.40**
A9000205070-004



<p>Ergo-Grip</p> <p><input type="checkbox"/> 20" handrim for 20" wheel A9000205080-001 € 402.40</p> <p><input type="checkbox"/> 22" handrim for 22" wheel A9000205080-002 € 402.40</p> <p><input type="checkbox"/> 24" handrim for 24" wheel A9000205080-004 € 402.40</p> <p>Attachment fitting to the rim</p> <p>Note: Due to the complex geometry, concentricity deviations of +/- 2-3 mm are possible.</p> 	<p>Maxgrepp with cellular and natural rubber coating very grippy, soft</p> <p><input type="checkbox"/> 20" handrim for 20" wheel A9000200788 € 402.40</p> <p><input type="checkbox"/> 22" handrim for 22" wheel A9000200789 € 402.40</p> <p><input type="checkbox"/> 24" handrim for 24" wheel A9000200787 € 402.40</p> 	<p>Curve black anodised, ergonomic handrim profile</p> <p><input type="checkbox"/> 24" handrim for 24" wheel Profile height 24 mm A9000205130-004 € 402.40</p> <p><input type="checkbox"/> 24" handrim for 24" wheel Profile height 31 mm Size L A9000205140-004 € 402.40</p> 
<p>Curve silver anodised ergonomic handrim profile</p> <p><input type="checkbox"/> 24" handrim for 24" wheel Profile height 31 mm Size L A9000205160-004 € 402.40</p> 	<p>Curve powder-coated, ergonomic handrim profile</p> <p><input type="checkbox"/> = frame colour (gloss surface)</p> <p><input type="checkbox"/> = black (gloss surface)</p> <p><input type="checkbox"/> 24" handrim for 24" wheel Profile height 24 mm A9000205135-004 € 467.60</p> <p><input type="checkbox"/> 24" handrim for 24" wheel Profile height 31 mm Size L A9000205145-004 € 467.60</p> 	<p>Curve Grip, with Tetra-Grip coating, ergonomic handrim profile</p> <p><input type="checkbox"/> 24" handrim for 24" wheel Profile height 24 mm A9000205155-004 € 520.70</p> 
<p>Quadro with Tetra-Grip coating</p> <p><input type="checkbox"/> 24" handrim for 24" wheel A9000205170-004 € 520.70</p> 	<p>Gekko black anodised</p> <p><input type="checkbox"/> 24" handrim for 24" wheel A9000205175-004 € 520.70</p> <p>with a robust and replaceable silicone lip, fully integrated into the profile</p> <p>Profile height 31 mm</p> 	<p>Gekko silver anodised</p> <p><input type="checkbox"/> 24" handrim for 24" wheel A9000205185-004 € 520.70</p> <p>with a robust and replaceable silicone lip, fully integrated into the profile</p> <p>Profile height 31 mm</p> 
<p>Gekko S black anodised</p> <p><input type="checkbox"/> 22" handrim for 22" wheel A9000205195-002 € 520.70</p> <p><input type="checkbox"/> 24" handrim for 24" wheel and 25" wheel A9000205195-004 € 520.70</p> <p>with a robust and replaceable silicone lip, fully integrated into the profile</p> <p>Profile height 24 mm</p> 	<p>Equipment without handrims</p> <p><input checked="" type="checkbox"/> A9000205990</p>	

Spoke guard

Finger- and spoke guard „transparent“

Not in combination with Spinergy drive wheels



for 20" wheel
A9000206300-001

€ 207.10

for 22" wheel
A9000206300-002

€ 207.10

for 24" wheel
A9000206300-004

€ 207.10

Tyres

If not specified, the „Schwalbe RightRun in black“ tyres will be delivered.

Schwalbe „RightRun“ High-pressure tyres

22x1" (489-25)
A9000204001-002

24x1" (540-25)
A9000204001-004

Colour: Black
7.5 bar



Schwalbe „RightRun“ High-pressure tyres

20x1" (451-25)
A9000204000-001

22x1" (489-25)
A9000204000-002

24x1" (540-25)
A9000204000-004

Colour: Grey, black wall
7.5 bar



Schwalbe „DownTown“ Profiled tyres

24x1" (540-25)
A9000204010-004
7.5 bar



Schwalbe „Marathon Plus Evolution“ with „smart guard“ puncture protection and reflector strips

22x1" (489-25)
A9000204030-002

€ 49.20

24x1" (540-25)
A9000204031-004

€ 49.20

Colour: Black
10 bar



Schwalbe „Marathon Plus Evolution“ with „smart guard“ puncture protection

24x1" (540-25)
A9000204030-014

€ 49.20

Colour: Black
10 bar



Schwalbe „MarathonPlus“ Schwalbe „Marathon Plus“ with „smart guard“ puncture protection

24x1" (540-25)
A9000204030-004

€ 49.20

Colour: Grey, black wall
10 bar



Schwalbe „One“

24x1" (540-25)
A9000204123-004

€ 93.00

Colour: Black
9 bar



Puncture-proof „solid rubber tyres“

22x1" (489-25)
A9000204200-002

€ 122.70

Colour: Grey

24x1" (540-25)
A9000204201-004

€ 122.70

Colour: Black
Not in combination with 18 spoke or Spinergy
drive wheels.



Accessories

24" MTB spare wheel set

with double rim (ETRTO 507), aluminum front hub, quick-release axle and stainless steel handrim
Mountain bike tyres, selectively:

A9000203410 € 703.40

24x1.75" Schwalbe "Marathon Plus"
Brake setting compatible with 24x1" wheel



A9000203413 € 703.40

24x2.1" Schwalbe "SMART SAM Sport"
Brake setting not compatible with 24x1" wheel




NOTE: We charge a surcharge for the use of drum brakes

A9000200004 € 232.20

Through axle with assistance, Large push button

A9000200096 € 105.80

Interlock dimension: 122 / 110 mm
Axle diameter: 12 mm
Not in combination with one-hand drive



Hand air compressor with pressure gauge

A9000901100 € 106.30

up to 14 bar



Crutch holder with flexible nylon cup

A9000902310* € 137.40



Length-adjustable safety belt with buckle

A9000902100 € 80.60



Adjustable safety belt with hook and fastener

A9000902110 € 64.50



SIGMA lighting set

A9000903270 € 107.90

Consisting of front light and rear lights, chargeable via Micro USB socket, incl. Micro USB charging cable
Lighting set compliant with road traffic regulations



Transport bag for drive wheels

for 24" wheel
A9000903530-004 € 77.80
Ø 76 cm x 2 cm



Therapy table

width 350 mm
(up to SB 37 cm) € 480.20
A9000902800-001

width 400 mm
(SB 38 - 41 cm) € 480.20
A9000902800-002

width 420 mm
(from SB 42 cm) € 480.20
A9000902800-003

removable table plate from Lexan, transparent, with front edging, with 2 positioning stops
Ø 25 mm (for adaptation to frame tube diameter Ø 25 mm)



Separate order sheets

Category „Rehabilitation technology accessories”

e.g. Seat & backrest systems, one-hand drive or transport restraint systems can be found at:
www.proactiv-gmbh.com/downloads/

Additional products

such as hand bikes or traction & attachment devices can be found at:
www.proactiv-gmbh.com/downloads/

Additional sheets

such as short abduction beginning etc. can be found at:
www.proactiv-gmbh.com/downloads/

Services

Consultation, clarification of the installation situation and measurement through PRO ACTIV

A9990000000 € 354.50

Handover and instruction on the wheelchair by PRO ACTIV

A9990000006 € 354.50

Notes and remarks:

Location, Date

The signature field is intended for the consultant and insured person. This confirms that the intended equipment is intended for the insured person.

- 18 - 2026-03-03 / The prices are plus statutory value-added tax and shipping costs

Signature End customer/dealer

= included in the basic price = Position selected