

Im Hofstätt 11
D-72359 Dotternhausen
Phone +49 (0)7427 9480-0
Fax +49 (0)7427 9480-7025
e-mail info@proactiv-gmbh.de
www.proactiv-gmbh.com



Quotation date: _____

Consultation was carried out by
Staff at the medical supplies shop: _____

PRO ACTIV field representative: _____

Order date: _____

According to GTC PRO ACTIV Reha-Technik GmbH

Medical suppliers/company address:

Cust. no.: _____

Name: _____

Street: _____

Postcode / city: _____

Phone: _____

Fax: _____

E-Mail: _____

Contact: _____

Invoicing address, other(see below)

Cust. no.: _____

Name: _____

Street: _____

Postcode / city: _____

Shipping address, other (see below)

Cust. no.: _____

Name: _____

Street: _____

Postcode / city: _____

Customer name/commission:

Surname,
forename: _____

Street: _____

Postcode / city: _____

Phone: _____

E-Mail: _____

D.O.B: _____

Paying
organisation: _____

Indication: _____



WHEEL-E

Up to 150 kg trailer load and 10 kg load capacity,
equipped with:

- Docking plate for quick coupling to the "Quick-Snap adapter"
- Indicator display for driving speed, support level and charge status display
- Driving speed up to 6 km/h with reverse gear, adjustable via drive signal transmitter on the right
- With hydraulic disc brake system and additional mechanical disc brake with parking function
- Parking stand with manoeuvring rollers
- Handlebars, attached to the angle-adjustable front end
- Sigma LED battery lighting set with front and rear light
- Bell
- Drive wheel 20" (ETRTO 406), black rim
- Black drive unit

**"The powerful, high-performance
drive for demanding conditions"**

Climbing power up to 15%

A9130100031-020

€ 6,400.00

Please note: When purchasing a WHEEL-E in combination with a wheelchair, the discounted VAT rate applies for the complete unit.

Body size

Body height: _____ cm

Upper body
height: _____ cm

Weight: _____ kg

Information essential to process the order.

Brake system

Parking disk brake, locking lever mounted to the right on the handlebar

A9130300106 € 153.40

Not in conjunction with double disk brake system

Double disc brake system "Magura" brake disk Ø 203 mm

A9130300105 € 193.20

Recommended from a towed weight of 80 kg and steep slopes, brake-lock function via locking bar



Rechargeable batteries and accessories

Battery 48 V 10.5 Ah 504 Wh

A9130900016

Range: approx. 25-45 km

Mounting the battery on the right, balancing weight mounted on the left



Double battery 48 V 10.5 Ah 504 Wh

A9130900017 € 1,483.50

Range: approx. 50-90 km

Positioned on the steering column with changeover plug connection via magnetic connector
incl. 2. chargers



Activation max. speed 15 km/h

A9130900101 € 138.00

Use only on private property, outside the scope of application of road traffic regulations

Activation max. speed 24 km/h

A9130900106 € 138.00

Use only on private property, outside the scope of application of road traffic regulations

LEGAL WARNING for the WHEEL-E towing device.

Equipping the WHEEL-E with a console for a maximum speed of 6 km/h.

The factory settings of the WHEEL-E PRO ACTIV are designed for an electric drive system that supports a maximum speed of 6 km/h.

Consoles for a top speed of 15 km/h or 24 km/h are optionally available. PRO ACTIV points out that the specific laws of your country must be observed for the use of the 15 km/h and 24 km/h consoles.

When ordering, please confirm with your signature that you have taken note of this fact.

Tyres

20x2.25" Schwalbe "Jumpin'Jack"

(Low availability – only as long as stocks last. Alternative option: Schwalbe "Smart Sam")

A9130105131
(ETRTO 406-57)



20x2.35" Schwalbe "Smart Sam"

A9117105101
(ETRTO 406-60)
Profiled tyres
Puncture protection level 3.5



20x2.40" Schwalbe "Super Moto-X"

A9130105132 € 27.40
(ETRTO 406-62)
Puncture protection level 6



20x4.0" BFT tyres

- A9130105130
(ETRTO 406-100)
Profile tyres, Fat Tire

€ 41.90




Handlebar

Steering column

- A9130600001
Steering column length
220 mm
- A9130600002
Steering column length
280 mm
- A9130600003
Steering column length
340 mm

Adjustable angle via clamped hinge on the fork bridge



Accessories

SKS mudguard

- A9130500003
removable



MX mudguard Colour: Orange

- A9130500002
Permanently mounted

€ 78.20



MX mudguard Colour: black

- A9130500004
Permanently mounted

€ 78.20



Aluminium luggage rack "cargo systems"

- A9110903290 € 179.70
Colour: black
Maximum load weight of 10 kg on the loading platform,
max. 15 kg in combination with luggage bags with low loading height.









Bottle holder with drinking bottle

- A9130900001 € 22.70
mounted on the rechargeable battery rail left, side side
Not in combination with double rechargeable battery



LED headlight

- A9130900004 € 116.00
Supply of the LED front light via the rechargeable battery of the drive system.
Back light as plug-on light supplied by a battery. Lighting set omitted.

Rear view mirror left <input type="checkbox"/> A9130900002 € 88.80	Rear view mirror right <input type="checkbox"/> A9130900003 € 88.80	
Options for users with reduced function		
Rotary throttle Triride <input type="checkbox"/> A9130200019 € 44.40 Instead of the thumb throttle lever, the rotary throttle is installed. Mounted on the right-hand side of the handlebar. 	Gas assistance for WHEEL-E <input type="checkbox"/> A9130200017 € 26.70 Can be clamped without tools to the WHEEL-E rotary grip 	Tetra gas actuation <input type="checkbox"/> A9130400100 € 327.80 right side 
Tetra brake actuation <input type="checkbox"/> A9130400101 € 357.70 Brake actuation side left, with brake actuation via rotary function with hand rest 	Chest strap for upper body fixation <input type="checkbox"/> Size S (105 cm) A9000800016 € 220.80 <input type="checkbox"/> Size M (125 cm) A9000800013 € 220.80 <input type="checkbox"/> Size L (145 cm) A9000800014 € 220.80 <input type="checkbox"/> Size XL (165 cm) A9000800017 € 220.80 Consisting of elastic webbing material in combination with neoprene with Velcro and operating tab. Colour: black 	Quick release clip lever extension for fixed frame adapter <input type="checkbox"/> A9110700004 € 98.90 to make the opening and closing of the quick release clips easier Only in combination with rigid frame adapter
Securing the lever extension with locking bolt <input type="checkbox"/> A9199000084 € 82.90 Can only be used for the lever extension on the quick release clip.	Eccentric connection lever extension for Quick Snap connection simplifies opening or closing the Quick Snap Adapter. Tetra version <input type="checkbox"/> A9110700005 € 70.30 Straight version Only in combination with rigid frame adapter <input type="checkbox"/> A9110700006 € 70.30 Curved version Only in conjunction with front adapter 	Centring column adaptation assistance <input type="checkbox"/> A9110700002 € 52.40 in order to make inserting the central adapter easier 
Parking stand with operating lever for high operating comfort <input type="checkbox"/> A9110900004 € 55.90 with longer actuating lever and reduced spring force		

Frame adapters

SPEEDY fixed frame adapter with "Quick Snap Adaptation"

with quick release clips and centring column mount

Colour: black

A9110700141 € 492.50

Heavy Duty version

A9110700143 € 623.90

up to a max. body weight of 200 kg



TRAVELER folding frame adapter with "Quick Snap Adaptation"

A9110701300 € 526.70

with quick release clips and centring column mount on the folding mechanism

Colour: black



BUDDY fixed frame adapter with "Quick Snap Adaptation"

A9110701400 € 492.50

with quick release clips and centring column mount on the frame

Colour: black



Wheelchair front adapter with universal clamping and "Quick-Snap Adaptation"

A9110701500 € 561.10

Adapter mount on the wheelchair front end

Please make the following selection:

A9110701525

Tube diameter 23-27 mm

Installation situation with PA products:

Frame stem G-shape from SB 31 cm to 54 cm

Frame stem V-shape from SB 30 cm to 50 cm

PLEASE NOTE:

If the tube diameter is not specified, as standard we shall deliver an adapter for a tube diameter of 23 - 27 mm.

Wheelchair front adapter Heavy Duty > 120 kg - 200 kg with universal clamping and "Quick-Snap Adaptation"

A9199701510 € 813.80

Adapter with reinforced double tube steering head connection for adaptations with wheelchairs for user weight > 120 kg – 200 kg

Adapter mount on the wheelchair front end



A9110701530

Tube diameter 28-32 mm

A9110701535

Tube diameter 33-37 mm

Clamp lever with 3-way star handle for actuating the frame clamping

A9199701507 € 101.00

mounted on the right and left for easy terminal operation with limited manual function
Only in conjunction with front adapter



Wheelchair-dependent accessories

Extended bike plate

- For SPEEDY A1/A2
A9110100100-002 € 352.20
- For BUDDY
A9110100100-003 € 352.20
- For TRAVELER
A9110100100-001 € 352.20
- 4° camber integrated
Track can be adjusted independently



Wheelbase extension via "grid mechanism with pulling operation"

- For SPEEDY 4all /
4teen / 4you /
TRAVELER 4you
A9110102001 € 537.50
- removable, with continuous axle tube
4°/ 6° camber integrated
Track can be adjusted independently



Wheelbase extension, angle-adjustable and removable

- For SPEEDY F2/F4
A9015100005 € 677.30
- with continuous axle tube,
height adjustment via gear joint
adjustable in 6° angular increments
wheel position displacement approx. 14 cm
4°/6° camber integrated, adjustable track
Attachment possibility must be checked by
PRO ACTIV.
Only in combination with base plate.
Up to 120 kg load capacity.



Front brake

- A9125500120 € 145.20
- as hill start assist with anti-roll back device,
compatible for mounting in combination with
the wheelchair knee lever brake.
Required for WHEEL-E with 15 km/h
individual operating permit.
- For combination with wheelbase extension,
item A9000700000 / A9000700173 Retaining
bars Adjustment range
295 mm is also required



Services

Handover and instruction for the handbike/traction device by the dealer incl. mounting the frame adapter

on a PRO ACTIV wheelchair

A9990000008
Prod. no.: _____

on a third-party product

A9990000005
Manufacturer: _____
Model: _____
Frame tube ø ____ mm

Notes and remarks:

Location, Date

The signature field is intended for the consultant and insured person. This confirms that the intended equipment is intended for the insured person.

Signature End customer/dealer