

Mobilität wird einfach!

Im Hofstätt 11
D- 72359 Dotternhausen
Phone +49 (0)7427 9480-0
Fax +49 (0)7427 9480-7025
E-Mail info@proactiv-gmbh.de
www.proactiv-gmbh.com



NJ1 Adaptive bike RS

with *dérailleurs* or free-wheel hub
and disk brake

up to 120 kg trailer load and 10 kg load capacity

- docking plate for quick coupling to the "Quick-Snap" adapter
- A-shaped handles with integrated switching and brake fittings
- hydraulic disc brake system
- secondary V-brake (rim brake) with locking mechanism
- parking stand
- chain protection tube
- lighting set Sigma with rechargeable battery charger, front and rear lights

9110100000 € 3,046.00

Quotation date: _____

Order date: _____

according to GTC PRO ACTIV GmbH

Medical suppliers/company address:

Cust. no.: _____
Name: _____
Street: _____
Postcode, city: _____
Phone: _____
Fax: _____
E-Mail: _____
Contact: _____

Invoicing address, other (see below)

Cust. no.: _____
Name: _____
Street: _____
Postcode, city: _____

Shipping address, other (see below)

Cust. no.: _____
Name: _____
Street: _____
Postcode, city: _____

Customer name/commission:

Name: _____
Street: _____
Postcode, city: _____
Phone: _____
Fax: _____
E-Mail: _____
DoB: _____
Paying organisation: _____
Indication: _____

Services:

1. Assembly of the frame adapter for PRO ACTIV

9990000001
on newly ordered PRO ACTIV wheelchair (assembly included)

9990000002 € 135.00
on existing PRO ACTIV wheelchair
Prod. no.: _____

2. Delivery, instruction and handing over at the customer by PRO ACTIV

9990000003 € 425.00
plus shipping costs

4. Assembly of the frame adapter by a specialist dealer

9990000008
on existing PRO ACTIV wheelchair
Prod. no.: _____
(assembly and delivery by specialist dealer)

9990000005
onto a third-party product
(compatibility testing, assembly and delivery by specialist dealer)

Body size

Body height: _____ cm

Upper body height: _____ cm

Weight: _____ kg

Information essential to process the order!

Colour

**Silver,
anodised drive unit**

9110100170-000

**Black Shadow,
deep black drive unit
RAL 9005, smooth gloss**

9110100170-001 **€ 103.80**

- fork tube R and L
- fork bridge, top and bottom
- bottom bracket support

Dérailleur

**Shimano Deore 3x9 Rapid Fire
27-speed chain shift dérailleur**

9110200020

Gear ratio:
sprocket 22-32-44 teeth,
cassette 11-13-15-17-20-23-26-30-32 teeth



**Shimano Deore XT 1x11 Rapid Fire
11- speed chain shift derailleur**

9110200034 **€ 151.00**

mechanical derailleur operation from the handle. 11-way cassette 11-42 teeth
Not in combination with Schwalbe "Knobby" 20x2" or Schwalbe "Big Ben" 20x2,15"!



**Shimano Deore XT 1x11 DI2
Button operation
11- speed chain shift derailleur**

9110200035 **€ 1,175.00**

electronic derailleur operation from the handle. 11-way cassette 11-42 teeth
Not in combination with Bereifung Schwalbe "Knobby" 20x2" oder Schwalbe "Big Ben" 20x2,15"!



**Shimano Deore XT 1x11 DI2
Rapid Fire
11- speed chain shift dérailleur**

9110200032 **€ 1,138.00**

electronic derailleur operation from the handle. 11-way cassette 11-42 teeth
Not in combination with Bereifung Schwalbe "Knobby" 20x2" oder Schwalbe "Big Ben" 20x2,15"!



**Shimano Deore XT 1x11 DI2
Kinnbedienung
11- speed chain shift derailleur**

9110200033 **€ 1,245.00**

electronic derailleur with chin operation (Execution Tetra).
11-way cassette 11-42 teeth
Not in combination with Bereifung Schwalbe "Knobby" 20x2" oder Schwalbe "Big Ben" 20x2,15"!



Hub drive system

**Shimano "Alfine"
11-speed hub gear**

9110200026 **€ 88.60**

Handle "Rapid Fire"
Gear ratio: Sprocket 32 teeth,
cassette 20 teeth
(with bottom bracket support L – 19 teeth)



**Shimano "Alfine DI2"
11-speed gear shift**

9100200029 **€ 1,135.00**

Electric gears, gear change using control buttons, with freewheel without back pedal brake

Gear ratio: Sprocket 32 teeth,
Chain sprocket 20 teeth.

Hydraulic back pedal brake possible as an additional option



**Rohloff "Speedhub 500/14"
14-speed hub gear**

9110200027 **€ 1,038.00**

"Grip Shift" twist grip shifter
Gear ratio: Sprocket 32 teeth,
Cassette 16 teeth



Gear change accessories

**Mountain drive 2.5: 1 gear reduction for
hills in combination with hub gear
Expansion to double number of gears**

9110200003 **€ 499.80**

Planetary gear set
All gears (the hub gear and the gear reduction for hills) can be changed when stationary!
In combination with grip width 40 or 45 cm!



**Mountain drive 2.5: 1 gear reduction for
hills in combination with chain shift
Expansion to double number of gears**

9110200006 **€ 371.30**

Planetary gear set
Replaced the 3-way sprocket and permits the switching of the gear reduction for hills when stationary!
In combination with grip width 40 or 45 cm!



Gear accessories for chain shift

**"Grip Shift" twist grip shifter
3-way, mounting left side**

9110200017 **€ 43.30**

SRAM Centera for 3-way dérailleur
In combination with Shimano Deore or XT chain shift!



**"Grip Shift" twist grip shifter
9-way, mounting right side**

9110200016 **€ 43.30**

SRAM Centera for 9-fold dérailleur
In combination with Shimano Deore or XT chain shift!



**Hand lever for direct switching of the
3-way dérailleur**

9110200030 **€ 76.70**

Direct hand operation of the dérailleur on the 3-way sprocket



Switching accessories for free-wheeling shifting hub

The following options are included with the dérailleurs in the delivery! Combined with the free-wheeling shifting hub, the number of gears is increased by a factor of 3!

Deore 3-way chain wheels with thumb switch, left side

9110000014 **€ 385.70**

For Alfine 11-speed and Rohloff 14-speed Hub gears; with front dérailleur to expand to 33 or 42 gears
Chain wheel sizes: Z 22; Z32; Z44
Dérailleur and chain wheel design: Shimano Deore



Deore 3-way chain wheels with dérailleur direct operation

9110000016 **€ 283.50**

For Alfine 11-speed and Rohloff 14-speed Hub gears; with front dérailleur to expand to 33 or 42 gears
Chain wheel sizes: Z 22; Z32; Z44
Direct dérailleur operation,
Dérailleur and chain wheel design: Shimano Deore



Brake system

Powerful HD brake system "Magura MT4"

9110300101 **€ 78.60**

recommended from a body weight of 80 kg



Second independent disk brake system "Magura MT4"

9110300100 **€ 183.90**

Shimano V-brake rim brake omitted
Only in combination with the powerful HD brake system "Magura"!

Parking brake can be operated via the gear lever

9110300002 **€ 66.60**

mounted on the bottom bracket support, 2nd brake handle on the crank handle, left side is omitted



Hydraulic back-peddalling disc brake with crank release function for moving backwards

9110300003 **€ 618.80**

Compatible with dérailleurs and freewheel gear hub in combination with:
- pedal bearing shaft grip widths 45 or 50 cm
- pedal crank crank lengths 155 or 175 mm up to 120 kg body weight
not in combination with
Mountaindrive 2.5:1 gear reduction for hills!



Dimensions

Pedal bearing shaft

Grip width 40 cm
9110500100-040

Grip width 45 cm
9110500100-045

Grip width 50 cm
9110500100-050

Grip width 55 cm **€ 57.20**
9110500100-055
Handle distance of more than 25 mm



Fixed bottom bracket support

M (X/Y 90/250 mm)
9110501100-001

M-L (X/Y 130/310 mm)
9110501100-004

L (X/Y 158/362 mm)
9110501100-002

XL (X/Y 158/410 mm) **€ 52.60**
9110501100-006

XXL (X/Y 170/465 mm) **€ 52.60**
9110501100-007



Adjustable bottom bracket support

9110500001 **€ 182.90**


Adjustable angle via clamped hinge on the fork bridge and in their length.

For the pre-setting, please specify the bottom bracket support size that comes as close to your desired position.

Size: _____

If no details are provided, the adjustable bottom bracket support shall be supplied in XL.

Deliverable up to size XL



Pedal crank

Length 155 mm
9110600100-155

Length 175 mm
9110600100-175

Length 195 mm
9110600100-195



Customer-specific crank length according to dimension

_____ mm (minimum length 125 mm)

9110600100-000 **€ 155.80**

Drive wheel

Drive wheel with double-wall rim

20" (406)
9110203050

24" (507)
9110104259

26" (559)
9110104260

Tyres

Schwalbe "Marathon"

- 20x1.5" (406-40)
9110105120
- 26x1.5" (559-35)
9110105100



Touring tyres

Continental "Ride Tour"

- 24x1.75" (507-47)
9110105127
- Touring tyres



Schwalbe "Big Apple"

- 20x2" (406-50)
9110105125

Balloon tyres



Schwalbe "Knobby"

- 20x2" (406-54)
9110105126

Profiled tyres



Schwalbe "Marathon Plus"

unplattbar, mit Pannenschutzeinlageband
Smart Guard

- 9130105128 € 23.80
20x1,75" (406-47)
- 24x1,75" (507-47) € 23.80
9110105135



Bike accessories

Plug-in protection plate

- 9110502100 € 45.80



SIGMA BC 1200 Bike Computer

- 9110901000 € 55.30
- wireless



SIGMA BC 1606L DTS Bike Computer

- 9110900018 € 133.20
- wireless
with wireless pedal frequency set



Bracket ø 22 mm x 42 mm for fixing of accessories (speedometer, lighting or similar)

- 9110903271 Item price € 18.30

Colour: natural anodised

_____ Items



Aluminium luggage rack "cargo systems"

- 9110903290 € 80.20

Colour: black
Max. load weight 10 kg on the loading platform, max. 15 kg in combination with saddle bags with low loading height.



Manoeuvring rollers for parking stands

- 9110900017 € 96.20

Manoeuvring rollers mounted to the parking stand via dual clamping bracket. (figure similar)



Manoeuvring bar mounted on the chain guard

- 9110500000 € 133.10
- with foam rubber cover



Options for users with reduced function

Tetra hand support, pivotable, adjustable grip distance

- 9110601101 € 140.30



Tetra shell handles

- 9110601100 € 135.90



Tetra backhand support

- 9110601102 € 158.60
- pivotable, adjustable distance to the bike handle



Quick release clip lever extension for fixed frame adapter

- 9110700004 € 82.00

in order to make the opening and closing of the quick release clips easier



Eccentric connection lever extension for Quick Snap connection

- 9110700005 € 58.10

in order to make the closing or opening of the Quick Snap adapter easier. Tetra version



Ratchet with lever extension for eccentric connection for quick snap adaption

- 9110700007 € 99.80

Makes it easier to close respectively open the quick snap adapter. Execution Tetra



Centring column adaptation assistance

- 9110700002 € 43.40

in order to make inserting the central adapter easier.



Frame adapters

SPEEDY fixed frame adapter with "Quick Snap Adaption"

9110700141 € 407.90

with quick release clips and centring column mount on the frame
Colour: black



TRAVELER folding frame adapter with "Quick Snap Adaption"

9110701300 € 436.00

with quick release clips and centring column mount on the folding mechanism
Colour: black



BUDDY fixed frame adapter with "Quick Snap Adaption"

9110701400 € 407.90

with quick release clips and centring column mount on the frame
Colour: black

Wheelchair front adapter with universal clamping and "Quick Snap Adaption"

9110701500 € 464.70

Adapter mount on the wheelchair front end

Installation situation with PA products:
Frame stem G-From from seat width 31 cm to 54 cm
Frame stem V-From from seat width 30 cm to 50 cm



PLEASE NOTE:

if the frame tube diameter is not specified, as standard we shall deliver an adapter for a frame tube diameter of 23-27 mm!

Please make the following selection:

9110701525

Frame tube diameter 23-27 mm

9110701530

Frame tube diameter 28-32 mm

9110701535

Frame tube diameter 33-37 mm

Wheelbase extension removable via "grid mechanism with pulling operation", with continuous axle tube

For Speedy 4all / 4teen / 4you / Traveler 4you 9110102001 € 445.10

4°/6° camber integrated, independently adjustable track


Wheelbase extension removable via "Quick-Pins", with continuous axle tube

For Speedy F2 9015101000 € 445.10

For Speedy F4 9110101000 € 445.10

4°/6° camber integrated, independently adjustable track

Attachment possibility for F2/F4 must be checked by PRO ACTIV!



chest strap for torso fixation

size S (105cm) 9000800016 € 221.50

size M (125cm) 9000800013 € 221.50

size L (145cm) 9000800014 € 221.50

size XL (165cm) 9000800017 € 221.50

Made out of elastic belt material with neoprene and velcro, with lug for easy opening.
Colour: black



your additional notes:


Wheelchair-dependent accessories

Extended bike plate

For Speedy A1/A2 9110100100-002 € 291.70

For Buddy 9110100100-003 € 291.70

For Traveler 9110100100-001 € 291.70



4° camber integrated, independently adjustable track

I herewith confirm the details listed above

Location, Date

Signature (End customer)

Location, Date

Signature (Dealer)