



☐ **Quotation date:** _____

Consultation was carried out by
Staff at the medical supplies shop: _____

Enquiry no.: _____

☐ **Order date:** _____

according to GTC PRO ACTIV Reha-Technik GmbH

Order no.: _____

Medical suppliers / company address:

Cust. no.: _____

Name: _____

Street: _____

Postcode / city: _____

Phone: _____

Fax: _____

E-Mail: _____

Contact: _____

☐ **Different invoicing address (see below)**

Cust. no.: _____

Name: _____

Street: _____

Postcode / city: _____

☐ **Different shipping address (see below)**

Cust. no.: _____

Name: _____

Street: _____

Postcode / city: _____

Customer name / commission:

ID / Comm. No.: _____

Surname, forename: _____

Street: _____

Postcode / city: _____

Phone: _____

E-Mail: _____

D.o.B.: _____

Paying organisation: _____

Indication: _____



HUSK-E adaptive bike

powerful Shimano version

up to 140 kg load capacity*

*load capacity = user weight + wheelchair + items transported on or in the wheelchair

- For users with reduced strength and greater impairments and for everyone who needs a powerful drive for a strong helping hand that can be accessed at any time
- HS version for best possible response behaviour
- Suitable for uneven terrain, hills, climbs and challenging routes
- derailleur and hub gear systems possible
- Disk brakes can be implemented with a brake lever or a back-pedal disk brake

basic equipment

- Shimano middle motor| 85Nm
- Docking plate for quick attachment to the „Quick-Snap adapter”
- adjustable and foldable bottom bracket support to minimize the packing size
- Mechanical parking disc brake
- Parking stand with longer operating lever
- LED lighting set for Shimano system
- Mudguard
- 20" drive wheel
- Small display steps
- Support level switch with ON/OFF button and light button

A9117100031-020

€ 7,390.00

Body size

Body height: _____ cm

Upper body height: _____ cm

User weight: _____ kg

This information is necessary for checking!

Heavy Duty version

Extended trailer load capacity*
from 140 kg to 170 kg

☐ A9117100100 € 507.00

Only possible in combination with an HD adapter.

Display

Small display steps

☐ included in the basic price

Features:

- time
- support mode
- battery charge level
- speed
- Distance and travel time
- total kilometers
- Estimated range and range overview
- pedal cadence
- gear display for Di2 gear systems

Only in combination with the green-backgrounded circuits



Large display steps

☐ A9117200110 € 61.30

replaces the smaller colour display, which is installed as standard.

Features:

- good readable display, large sharp writing
- many display contents on the same screen

Only in combination with the green-backgrounded circuits



Large display E6100

recommended for Tetra configuration

☐ A9117200100 € 130.00

replaces the smaller colour display EM800, which is installed as standard.

Features:

- on and off button, system start of the display
- light button (on/off at the display)
- mode button
- good readable display, large sharp writing
- many display contents on the same screen

Only in combination with the orange-backgrounded circuits for Tetra configuration



Transmission

Chainring on bottom bracket with 28 teeth

☐ A9117200050
for a higher top speed


Chainring on bottom bracket with 22 teeth

☐ A9117200051
for a very light response

Deraillieur gears

Shimano Deore 1x11 Rapid Fire 11-speed deraillieur gear

☐ A9117200026
mechanical deraillieur gear system with control on the right handle.
11-speed cassette 11-42 teeth
and a gear range of 381 %
with hydraulic disc brake Ø 203 mm



SHIMANO
DEORE

Shimano CUES 1x10 Di2 with push-button operation 10-speed Linkglide deraillieur gear

☐ A9117200035 € 292.00

electronic deraillieur gear system with control on the handlebar via a button and automatic function.

Operation: right side ☐ left side ☐

If not specified, the operation will be mounted on the „right side“.

10-speed cassette 11-43 teeth
and a gear range of 390 %
with hydraulic disc brake Ø 203 mm



SHIMANO
CUES
Di2

Nexus 5-speed gear hub Di2 with push-button operation


☐ A9117200036 € 496.00

electronic hub gear with automatic shift function and manual shift function via control button.

Operation: right side ☐ left side ☐

If not specified, the operation will be mounted on the „right side“.

with hydraulic disc brake Ø 180 mm.
with additional mechanical parking disc brake.




Gear hub for tetra configuration

Shimano Deore XT 1x11 Di2 with chin control

☐ A9117200023 € 1,224.30


electronic derailleur with chin control (tetra equipment).
11-speed cassette 11-42 teeth
and a gear ratio range of 381 %.



Nexus 5-speed hub gear Di2 with chin control

☐ A9117200033 € 984.40

electronic hub gear with automatic shifting function and manual shifting function via chin control
button mounted centrally on the bottom bracket



Brake system

Double hydraulic disc brake system „Magura“

☐ A9117300100 € 168.10

Two Magura 2-piston callipers
Parking brake via locking bar on the handle
Omitting the mechanical parking disk brake.

Mechanical back-pedal disc brake with crank release function

☐ A9117300110 € 738.40

compatible with derailleur gears and freewheel hub.
Incl. second mechanical disc brake designed as a parking brake.
The hydraulic disc brake is omitted.
In combination with:
- grip width 45 cm, 50 cm, 53 cm or 58 cm
- crank length 155 mm or 175 mm
up to 120 kg body weight



Driving brake - supplement for the mechanical parking brake

☐ A9117300120 € 116.10

The parking disc brake is supplemented with a brake handle on the left-hand handle.
This means that the mechanical disc brake can be used as a driving brake and parking brake.

Batteries and accessories

The specified ranges depend on the distance profile, user weight and speed.

Battery Shimano 36V 14Ah 504Wh

Range approx. 45 km - 75 km

☐ A9117900150
mounted on the left of the fork

Battery Shimano 36V 17.5Ah 630Wh

Range approx. 55 km - 85 km

☐ A9117900155 € 186.60
mounted on the left of the fork



Double Battery Shimano 14Ah 504Wh

Range approx. 90 km - 150 km

☐ A9117900160 € 1,467.80

mounted on both sides of the fork,
interchangeable connector mounted on the fork
for easy switching between the batteries.
Bracket swivels to the right for removing the battery.
Only in combination with frame adapter with right-hand operation.

Double Battery Shimano 17.5Ah 630Wh

Range approx. 110 km - 170 km

☐ A9117900165 € 1,872.70

mounted on both sides of the fork,
interchangeable connector mounted on the fork
for easy switching between the batteries.
Bracket swivels to the right for removing the battery.
Only in combination with frame adapter with right-hand operation.



Colour mudguard

If not specified, the mudguard and chain guard will be delivered in anthracite RAL 7016 matt.

Mudguard, chain guard and motor cover anthracite RAL 7016 matt

☐ A9117500004

Mudguard, chain guard and motor cover pure orange RAL 2004 gloss

☐ A9117500006



Mudguard, chain guard and motor cover lemon green gloss

☐ A9117500007



Mudguard, chain guard and motor cover customized colour selection

☐ A9117500010 € 153.50

from colour palette one-colour design
powder colours wheelchair

A9000100170-____

Colour: _____

Tyres

20x2.15" Schwalbe „Marathon 365“

☐ A9117105103
(ETRTO 406-55)
profile tyre
puncture protection level 5



20x2.15" Schwalbe „Big Ben Plus“

☐ A9117105102 € 28.70
(ETRTO 406-55)
profile tyre
puncture protection level 6



Dimensions

Crank type and grip width

Standard crank

- ☐ grip width 40 cm
A9117500100-040
- ☐ grip width 45 cm
A9117500100-045
- ☐ grip width 50 cm
A9117500100-050

Tetra crank (flared)

with extended crank arms for maximum upper body space

- ☐ grip width 48 cm
A9117500100-048
- ☐ grip width 53 cm
A9117500100-053
- ☐ grip width 58 cm
A9117500100-058

Bottom bracket support

maximum length from the pivot point in the fork to the bottom bracket shaft

- ☐ length 520 mm (short)
using shortened holding tubes € 155.30
A9117500001
recommendation: up to approx. 165 cm body height
- ☐ length 570 mm (standard)
A9117500008
recommendation: from approx. 165 cm to 185 cm body height
- ☐ length 620 mm (long)
using a longer bottom bracket distance € 204.50
A9117500003
recommendation: from approx. 185 cm body height

Standard pedal crank ☐ length 155 mm A9110600100-155 ☐ length 175 mm A9110600100-175 ☐ length 195 mm A9110600100-195 	Tetra pedal crank ☐ length 155 mm A9100600105-155 € 136.10 ☐ length 175 mm A9100600105-175 € 136.10 ☐ Individual length 140 mm to 195 mm A9100600100 € 258.70 	
--	---	--

Handles

Palm rest on the handle ☐ A9110601200 € 102.80 upholstered, bolted, angle-adjustable. 	Tetra hand support, pivotable, adjustable grip distance ☐ A9110601101 € 169.40 	Tetra backrest swivel and adjustable grip width ☐ A9110601102 € 191.60 very stiff external support as round tube Ø 20 mm with foam rubber coating 
Tetra ergonomic handle, freestanding ☐ A9100600001 € 164.60 crank handles arranged vertically can be combined with all other handle options 	Tetra shell handles ☐ A9110601100 € 172.40 	

Accessories Bike

Maneuvering castors for parking stand ☐ A9110900017 € 122.10 Maneuvering castors mounted on the parking stand via inlet 	Luggage bags including mounts for the fork tubes ☐ A9117900220 € 315.40 two panniers, each with a volume of 5.8 liters Colour of mount: silver Colour of bags: dark sand Max. load per side: 5 kg 	Operating satellite on the handle Mounting: Right side ☐ Left side ☐ with existing support line (laid along a brake line) ☐ A9117200040 € 89.50 without existing support line (an additional blind line is required) ☐ A9117200041 € 126.90
--	--	---

Options for users with reduced functions

Lever extension excenter lock for Quick Snap connection

☐ A9110700006 € 73.90

simplifies the opening and closing of the Quick Snap adapter
Tetra version, curved



Adaptation help for centering pin

☐ A9110700002 € 55.00
to simplify the insertion of the central adapter.



Chest strap for upper body fixation

☐ size S (105 cm) A9000800016 € 241.80

☐ size M (125 cm) A9000800013 € 241.80

☐ size L (145 cm) A9000800014 € 241.80

☐ size XL (165 cm) A9000800017 € 241.80

consisting of elastic webbing material in combination with neoprene with Velcro and operating tab.
Colour: black



Frame adapter

SPEEDY rigid frame adapter right-hand operation with „Quick-Snap-Adaption“

☐ A9110700142 € 537.20

with quick-release clamp and centering pin holder
Colour: black

Heavy Duty version

☐ A9110700144 € 816.70

for a load capacity* of 140 kg – 220 kg

Please select accordingly:

Centring pin for Ø 25 mm
for A1 and A2

☐ A9110700022

Centring pin for Ø 40 mm
for 4all, 4all Ergo, 4you Ergo, F2, F3 and 4teen

☐ A9110700023



TRAVELER folding frame adapter right-hand operation with „Quick-Snap-Adaption“

☐ A9110701302 € 537.20

with quick-release clamp and centering pin holder on the folding mechanism
Colour: black

Heavy Duty version

☐ A9110700146 € 816.70

for a load capacity* of 140 kg – 200 kg



BUDDY rigid frame adapter right-hand operation with „Quick-Snap-Adaption“

☐ A9110701401 € 537.20

with quick-release clamp and centering pin mount on the frame
Colour: black



Wheelchair front adapter right-hand operation with universal clamp and „Quick-Snap-Adaption“ ☐ A9110701502 € 572.30 Adapter attachment to the wheelchair stem. Operation on the right-hand side. Mounting situation for PA products: Straight front end from SB 31 cm to 54 cm V-shaped front end from SB 30 cm to 50 cm Please select the frame clamp below: ☐ clamping range 95 mm A9110701525 Tube diameter 23-27 mm	Wheelchair front adapter HD right-hand operation with universal clamp and „Quick-Snap-Adaption“ ☐ A9199701511 € 871,60 Adapter with reinforced double-tube steering head connection for wheelchair adaptations with a load capacity* of 140 kg – 220 kg Adapter mounting on the wheelchair frame extension	
☐ clamping range 70 mm € 76.50 A9110900072 Tube diameter 23-27 mm ☐ clamping range 75 mm € 268.10 A9199000080 Tube diameter KÜSCHALL 2.0 PLEASE NOTE: If the exact tube diameter is not specified, an adapter for a tube diameter of 23 - 27 mm and a clamping range of 95 mm will be delivered as standard. Please select the clamping sleeve below: ☐ clamping sleeve [80 mm] included as standard for close clearance between footrest and drive wheel. PLEASE NOTE: If not specified, the clamping sleeve with 80 mm length will be delivered.	☐ clamping range 95 mm A9110701530 Tube diameter 28-32 mm ☐ clamping range 70 mm € 76.50 A9110900078 Tube diameter 28-32 mm ☐ clamping sleeve [110 mm] € 46.70 A9199000117 for medium distance between footrest and drive wheel.	☐ clamping range 95 mm A9110701535 Tube diameter 33-37 mm ☐ clamping range 70 mm € 76.50 A9110900079 Tube diameter 33-37 mm ☐ clamping sleeve [170 mm] € 56.60 A9199000154 for a wide distance between footrest and drive wheel.
Clamp lever with 3-way star handle for actuating the frame clamping ☐ A9199701507 € 106.10 mounted on the right and left for easy terminal operation with limited manual function Only in combination with front adapter 		

Wheelchair-dependent accessories

Extended bike plate

- ☐ for SPEEDY A1 / A2
A9110100100-002 € 370.00
- ☐ for BUDDY
A9110100100-003 € 370.00
- ☐ for TRAVELER
A9110100100-001 € 370.00
- 4° camber integrated
Track independently adjustable



Wheelbase extension via „grid mechanism with pulling operation”

- ☐ for SPEEDY 4all /
4teen / 4you /
TRAVELER 4you
A9110102001 € 564.70

removable, with continuous axle tube, track
width independently adjustable
Up to 150 kg user weight.

Please select the wheel camber:

- ☐ 1° drive wheel camber
☐ 4° drive wheel camber
☐ 6° drive wheel camber

Note: The camber of the wheelbase extension
must not be less than that in the normal driving
position. A maximum camber of 6° is permitted.



Wheelbase extension angle-adjustable and removable

- ☐ for SPEEDY F2 / F3
A9015100005 € 746.80

with continuous axle tube,
height adjustment via gear joint
adjustable in 6° angular increments
wheel position displacement approx. 14 cm
adjustable track

Attachment possibility must be checked by
PRO ACTIV.

Only in combination with base plate.
Up to 120 kg user weight.

Please select the wheel camber below:

- ☐ 1° drive wheel camber
☐ 4° drive wheel camber
☐ 6° drive wheel camber

Note: The camber of the wheelbase extension
must not be less than that in the normal driving
position. A maximum camber of 6° is permitted.



Services

Consultation, clarification of the installation situation and measurement through PRO ACTIV

- ☐ A9990000000 € 354.50

Handover and instruction of the handbike / traction device by PRO ACTIV

- ☐ A9990000011 € 354.50

Attachment of the handbike / traction device to the wheelchair by PRO ACTIV at specialist dealer

on a PRO ACTIV wheelchair

- ☐ A9990000012 € 230.00

product no.: _____

on a third-party product

- ☐ A9990000013 € 366.90

manufacturer: _____

model: _____

frame tube Ø ____ mm

Attachment of the handbike / traction device to the wheelchair through specialist dealer

including mounting of the frame adapter

on a PRO ACTIV wheelchair

- ☐ A9990000008

product no.: _____

on a third-party product

- ☐ A9990000005

manufacturer: _____

model: _____

frame tube Ø ____ mm

Notes and remarks: