

# FOXX



## ***BUILT FOR THE WILD. MADE FOR YOU.***

### **FOXX off-road compact bike - adaptive mountain biking rethought.**

Whether steep climbs, technical trails or flowing passages - the FOXX is your reliable companion. With innovative high-pivot kinematics, sophisticated geometry and a powerful drive system, it will take you anywhere - with maximum control, comfort and power.

### **Big on performance. Compact in size.**

Despite its off-road power, the FOXX remains compact and easy to transport with a length of 2 metres, a wheelbase of 1.20 metres and a weight of just 53 kg. The pack size of the FOXX can be adjusted in a few simple steps for easy transport.

### **Maximum performance with high-pivot kinematics.**

The 4-bar linkage system ensures perfect grip, efficiency and comfort. Thanks to the intelligent gear ratio and sophisticated anti-squat behaviour, you can crank powerfully and efficiently, while the sensitive suspension absorbs bumps with ease. Ideal for any terrain - whether uphill, downhill or cross-country.



Front wheel | Rear wheel  
24" | 27,5"



Full suspension  
100 mm | 140 mm



Motors up to  
500 W & 130 Nm



Range  
up to 150 km

# WIDE RANGE OF OPTIONS

## Precision and traction on any surface

The FOXX's innovative independent wishbone suspension guarantees you maximum traction and control. The individually adjustable suspension and the stabiliser between the wheels ensure stable cornering and safety.

## Customised suspension

Choose from high-quality suspension options, perfectly matched to your riding style and area of use. All dampers have three riding modes that can be adjusted depending on the terrain and desired response behaviour.



## Comfort meets control

The customisable back system with adjustable back angle, ergonomic seat cushion and individually adjustable legrests ensure a perfect riding position and maximum riding stability.

## Dual battery & quick change

Two selectable powerful drive systems guarantee pure power. Can be combined with a dual battery system for a range of up to 150 km.

## TECHNICAL DATA

Suspension:	Full suspension with 100 mm travel front and 140 mm travel rear Standard suspension: Front axle: FOX Float DPS Performance   Rear axle: Manitou Mara Inline Offroad suspension: Front axle: FOX Float DPS Factory   Rear axle: Manitou Mara Pro
Mid-motor:	Bafang M560   130 Nm   500 W   3.3 kg Shimano   85 Nm   250 W   2.7 kg
Battery:	Bafang 48 V   14.5 Ah   696 Wh   Range approx. 50 - 75 km   3.77 kg Shimano LiMn 36 V   17.5 Ah   630 Wh   Range approx. 45 - 70 km   3.5 kg
Wheel sizes:	Front wheel size: 24"   Rear wheel size: 27,5"
Ground clearance:	19 cm
Seat width:	37 cm
Weight:	from 53 kg (with battery)
Pack size:	from 110 x 80 x 56 cm

MADE IN GERMANY

# PROACTIV®

PRO ACTIV Reha-Technik GmbH

Im Hofstätt 11  
72359 Dotternhausen

E-mail: [info@proactiv-gmbh.de](mailto:info@proactiv-gmbh.de) | Tel. +49 (0) 7427 9480-0

[www.proactiv-gmbh.com](http://www.proactiv-gmbh.com)

Rev. 05/2025

All information is subject to design changes, technical adjustments, printing errors and colour deviations.



FOLLOW US ON SOCIAL MEDIA

