



BUDDY

**Adaptive rigid frame wheelchair
for children and teenagers
with growing-with-you function***

*Equipment-related replacement parts required

Maximum user weight limit 100 kg

A9041130001 € 4,080.00

☐ **Quotation date:** _____

Consultation was carried out by
Staff at the medical supplies shop: _____

Enquiry no.: _____

☐ **Order date:** _____

according to GTC PRO ACTIV Reha-Technik GmbH

Order no.: _____

Medical suppliers / company address:

Cust. no.: _____
Name: _____
Street: _____
Postcode / city: _____
Phone: _____
Fax: _____
E-Mail: _____
Contact: _____

Invoicing address: ☐ same as company address

☐ Other (see below)

Cust. no.: _____
Name: _____
Street: _____
Postcode / city: _____
Phone: _____
Fax: _____
E-Mail: _____
Contact: _____

Shipping address: ☐ same as company address


☐ Other (see below)

Name: _____
Street: _____
Postcode / city: _____

Customer name / commission:

ID / Comm. No.: _____
Surname, forename: _____
Street: _____
Postcode / city: _____
Phone: _____
Fax: _____
E-Mail: _____
D.o.B.: _____
Paying organisation: _____

Body size:



Body height: _____ cm


Weight: _____ kg

Lower leg length: _____ cm
(anatomic dimension: hollow of the knee to the foot sole)

Indication: _____

This information is necessary for checking. If the information is missing, a check of the relevant items cannot be made!

Basic size



Seat width **SB**

☐ 20	☐ 21	☐ 22	☐ 23	☐ 24	☐ 25	☐ 26	☐ 27	☐ 28	☐ 29
☐ 30	☐ 31	☐ 32	☐ 33	☐ 34	☐ 35	☐ 36	☐ 37	☐ 38	

Size M (Frame length 36 cm (to the bend) and start of abduction at ST 27 cm)

☐ 22	☐ 23	☐ 24	☐ 25	☐ 26	☐ 27	☐ 28	☐ 29	☐ 30	☐ 31
-----------------------------	-----------------------------	-----------------------------	-----------------------------	-----------------------------	-----------------------------	-----------------------------	-----------------------------	-----------------------------	-----------------------------

Seat depth **ST**

Size L (Frame length 39 cm (to the bend) and start of abduction at ST 30 cm)

☐ 26	☐ 27	☐ 28	☐ 29	☐ 30	☐ 31	☐ 32	☐ 33	☐ 34	☐ 35
-----------------------------	-----------------------------	-----------------------------	-----------------------------	-----------------------------	-----------------------------	-----------------------------	-----------------------------	-----------------------------	-----------------------------

Size XL (Frame length 43 cm (to the bend) and start of abduction at ST 34 cm)

☐ 31	☐ 32	☐ 33	☐ 34	☐ 35	☐ 36	☐ 37	☐ 38	☐ 39	☐ 40
-----------------------------	-----------------------------	-----------------------------	-----------------------------	-----------------------------	-----------------------------	-----------------------------	-----------------------------	-----------------------------	-----------------------------

Back height **RH**

☐ 20	☐ 21	☐ 22	☐ 23	☐ 24	☐ 25	☐ 26	☐ 27	☐ 28	☐ 29
☐ 30	☐ 31	☐ 32	☐ 33	☐ 34	☐ 35	☐ 36	☐ 37	☐ 38	☐ 39
☐ 40	☐ 41	☐ 42	☐ 43	☐ 44	☐ 45	☐ 46	☐ 47	☐ 48	

Lower leg length **USL**

☐ 20	☐ 21	☐ 22	☐ 23	☐ 24	☐ 25	☐ 26	☐ 27	☐ 28	☐ 29
☐ 30	☐ 31	☐ 32	☐ 33	☐ 34	☐ 35	☐ 36	☐ 37	☐ 38	☐ 39
☐ 40	☐ 41	☐ 42	☐ 43	☐ 44	☐ 45	☐ 46	☐ 47		

Seat height back **SHh**

with 20" drive wheels									
☐ 31	☐ 32	☐ 33	☐ 34	☐ 35	☐ 36	☐ 37			
with 22" drive wheels									
☐ 33	☐ 34	☐ 35	☐ 36	☐ 37	☐ 38	☐ 39	☐ 40		
with 24" drive wheels									
☐ 36	☐ 37	☐ 38	☐ 39	☐ 40	☐ 41	☐ 42			

Seat height front **SHv**
(with 5" caster wheel)

☐ 40	☐ 41	☐ 42	☐ 43	☐ 44	☐ 45	☐ 46	☐ 47
-----------------------------	-----------------------------	-----------------------------	-----------------------------	-----------------------------	-----------------------------	-----------------------------	-----------------------------

Tipping point **KP**
(inner edge of rear tube to
centre of drive wheel axle)

☐ <u>5</u>	☐ 6	☐ 7	☐ 8	☐ 9	☐ 10
-----------------------------------	----------------------------	----------------------------	----------------------------	----------------------------	-----------------------------

Please tick the required dimensions (all dimensions are in cm)

The underlined values are the standard settings. If not specified, the standard value will be delivered.

The SHh and SHv values might be a bit different from the specified values,
depending on the wheel size and camber.

We recommend a seat height difference (difference between SHv and SHh) between 5 cm and 12 cm.

Frame colour - One-colour design		
Information on current colour designs and example images can be found at: https://www.proactiv-gmbh.com/colours/		
Frame colour: Selection from colour palette		
Frame colour will be specified in the order ☒ A9000100170-000	Raw, blasted, not coated ☒ A9000100170-020 Note: In combination with moisture e.g. hand perspiration, can affect dark aluminum discoloration.	
Powder colours		
white RAL 9016 ☒ gloss A9000100170-023 ☐ matt A9000100170-123 € 43.40	brillant silver ☒ gloss A9000100170-011	anthracite RAL 7016 ☒ gloss A9000100170-001 ☐ matt A9000100170-101 € 43.40
black RAL 9005 ☒ gloss A9000100170-003 ☐ matt A9000100170-019 € 43.40	pure orange RAL 2004 ☒ gloss A9000100170-022 ☐ matt A9000100170-122 € 43.40	pink magenta RAL 4010 ☒ gloss A9000100170-037 ☐ matt A9000100170-137 € 43.40
pastel blue RAL 5024 ☒ gloss A9000100170-036 ☐ matt A9000100170-136 € 43.40	lemon green ☒ gloss A9000100170-005 ☐ matt A9000100170-105 € 43.40	yellow RAL 1023 ☒ gloss A9000100170-046 ☐ matt A9000100170-146 € 43.40
red RAL 3001 ☒ gloss A9000100170-047 ☐ matt A9000100170-147 € 43.40	lagoon blue ☒ gloss A9000100170-052 ☐ matt A9000100170-152 € 43.40	athletic grey RAL 7046 ☒ gloss A9000100170-050 ☐ matt A9000100170-150 € 43.40
muddy grey RAL 7033 ☒ gloss A9000100170-051 ☐ matt A9000100170-151 € 43.40	peach RAL 3022 ☒ gloss A9000100170-038 ☐ matt A9000100170-138 € 43.40	
Powder colours „dormant“		
deep black effect silver ☒ gloss A9000100170-007	sky blue ☒ gloss A9000100170-017	navy blue ☒ gloss A9000100170-045 ☐ matt A9000100170-145 € 43.40
electric green ☒ gloss A9000100170-015 ☐ matt A9000100170-115 € 43.40	violet ☒ gloss A9000100170-044	fire red ☒ gloss A9000100170-018 ☐ matt A9000100170-118 € 43.40
Frame colour - One-colour design special colours „liquid coating“		
€ 349.60		
storm (ice blue) ☐ gloss A9000100170-500 ☐ matt A9000100170-501	prairie (rose gold) ☐ gloss A9000100170-504 ☐ matt A9000100170-505	berry ☐ gloss A9000100170-526 ☐ matt A9000100170-527
blue ocean ☐ gloss A9000100170-522 ☐ matt A9000100170-523	emerald ☐ gloss A9000100170-514 ☐ matt A9000100170-515	chili ☐ gloss A9000100170-516 ☐ matt A9000100170-517

Frame colour - Two-colour design

Note: In the two-colour design, both colours are the same gloss level as the design colour.

☐ **Two-colour design with one base colour (on the rear frame) and one design colour (on the front frame)**

€ 268.50

Please select a colour below:

Base colour (on the rear frame)

white RAL 9016

☐ A9000100170-200

brilliant silver

☐ A9000100170-201

anthracite RAL 7016

☐ A9000100170-202

black RAL 9005

☐ A9000100170-203

pastel blue RAL 5024

☐ A9000100170-206

lemon green

☐ A9000100170-207

yellow RAL 1023

☐ A9000100170-208

lagoon blue

☐ A9000100170-212

athletic grey RAL 7046

☐ A9000100170-210

muddy grey RAL 7033

☐ A9000100170-211

Design colour (on the front frame)

black RAL 9005

☐ gloss A9000100170-259

☐ matt A9000100170-260

lemon green

☐ gloss A9000100170-250

☐ matt A9000100170-251

chili

☐ gloss A9000100170-253

☐ matt A9000100170-254

pure orange RAL 2004

☐ gloss A9000100170-255

☐ matt A9000100170-256

emerald

☐ gloss A9000100170-274

☐ matt A9000100170-275

storm (ice blue)

☐ gloss A9000100170-264

☐ matt A9000100170-265

prairie (rose gold)

☐ gloss A9000100170-266

☐ matt A9000100170-267

berry

☐ gloss A9000100170-284

☐ matt A9000100170-285

blue ocean

☐ gloss A9000100170-280

☐ matt A9000100170-281

Colour - Attachment parts

Note: If not specified, the caster fork and bearing block will be mounted in black.

Drive wheel: aluminum rim, hub and spokes

☐ Black A9000100171-003

Caster fork and bearing block

☐ Silver A9000100173-011

☐ Black A9000100173-003

Colour design of attachment parts according to customer request:

Drive wheel: aluminum rim and design inlet¹ according to colour selection / hub and spokes black

☐ frame colour A9000100171-999 € 249.60

☐ Selection from one-colour design (different from the frame colour) A9000100171-998 € 288.90

Colour: _____

A9000100170- ☐ ☐ ☐

Caster wheel: 5" aluminum rim² (6-spoke design)

☐ frame colour A9000100174-999 € 75.20

☐ Selection from one-colour design (different from the frame colour) A9000100174-998 € 87.10

Colour: _____

A9000100170- ☐ ☐ ☐

Special colour request (Frame colour - One-colour design)

Frame colour for single purchasing Selection from RAL colour palette

☐ A9000100170-100 € 414.40

gloss level

☐ gloss

☐ matt

colour designation

RAL number

Note: Attachments **not** available in design colour options

¹ Only in combination with Liteline hub

² Only in combination with caster wheel 5" aluminum rim 6-spoke-design good shock absorption A9000500121*, as well as with suspension fork 5-spoke-design A9000500160*

Frame

Front end abducted shape

☒ A9041130030-K*
Knee angle 105°
Bearing block permanently welded



Push handles

Equipment without push handles

☐ A9000401000



Safety push handles back-positioned

☐ Length: 45 cm
A9041401075-001
max. height = SHh + approx. 53 cm

☐ Length: 55 cm
A9041401075-002
max. height = SHh + approx. 63 cm
continuous height adjustment
Safety push handle holder mounted on the back cross bar



Aluminum push handles Back positioned with side handle offset 145 mm

☐ A9041401080-K* € 89.50
height adjustment
Safety push handle holder mounted on the back cross bar
Length: 55 cm
Up to seat width SB 38 cm or body weight 50 kg



Aluminium central push handle foldable handle area

☐ A9041401086-K € 102.60
height-adjustable and foldable in the handle area via toothed joint, bracket welded to the center of the back crossbar.
Up to seat width of 38 cm or body weight 50 kg
Not in combination with headrest mount.
Caution: Not suitable for carrying the wheelchair - hold the back crossbar.




Central push handle for brake handles foldable grip area

☐ A9041401070-K € 141.90
closed, with continuous height adjustment, with foam rubber handle cover
Up to seat width SB 38 cm or body weight 50 kg
Not in combination with headrest mount
Only in combination with parking lever brake for accompanying person or drum brake.



Aluminium push handle closed, foldable via toothed joint

☐ A9041401069-K € 389.90
with continuous height adjustment, push handle holder back-positioned
Up to seat width SB 38 cm or body weight 60 kg



Backrest

Fixed backrest

☐ A9041420020
screwed to aluminum side parts
Back angle 90° to the floor
Backrest cannot fold down
Side parts in frame colour



Folding backrest

☐ A9000400023 € 306.30
7-way angle adjustable, with integrated plastic-coated aluminum side parts



Back equipment

If not specified, the „Adaptable back with back padding“ will be delivered.

Adaptable back with back padding

☒ A9000802010-K001
Cover: Black, breathable, open surface
Equipment with back crossbar



Adaptable back with back padding „light“

☒ A9000802320-K001
Cover: Black, breathable, open surface
Cushion thickness: 9 mm
Equipment with back crossbar

☒ A9000802320-K002
Cover: Black, breathable, open surface
Cushion thickness: 9 mm
Equipment without back crossbar



Adaptable back with back padding „Comfort“

☐ A9000815101-K **€ 131.40**
Cover: Black, breathable, open surface
Front padding: 20 mm
Rear padding: 5 mm
Equipment with back crossbar

Ergonomic back shell aluminum from back height RH 24 cm


☐ A9000802210-K **€ 441.60**
with lumbar support.
Colour: Anthracite matt
(If no colour is specified under „Select colour“)
If required, please select the back cushion separately.



Back cushion for ergonomic back shell

☐ A9000802220-K **€ 221.30**
cushion thickness: 2 cm

Cover:
☐ Black, breathable, open surface
☐ Black, breathable, closed surface
If not specified, the cover material „Black, breathable, open surface“ will be delivered.



ERGO Back backrest bar for ergonomic back shell from back height RH 24 cm and seat width SB 30 cm

☐ A9000802019-K **€ 242.60**
with recessed handle area
Colour: Anthracite matt
(If no colour is specified under „Select colour“)
Only compatible with central push handle
Item no. A9041401082
If required, please select the back shell and the back cushion separately.



Equipment without back upholstery and back padding or back system

☐ A9000802011 **€ - 53.00**
Equipment with back crossbar



Back crossbar flared

☐ by 20 mm A9000400835-020 **€ 109.30**
☐ by 30 mm A9000400835-030 **€ 109.30**
☐ by 50 mm A9000400835-050 **€ 140.80**
☐ by 80 mm A9000400835-080 **€ 155.30**
Only in combination with back-positioned push handles




Seat upholstery

If not specified, the „Body Contour“ seat upholstery will be delivered.

Body Contour

☒ A9000801025-K



Seat plate aluminum

☐ A9000801212-K **€ 292.70**
screwed onto the frame tubes
Colour: anthracite matt
(If no colour is specified under „Select colour“)



Equipment without seat upholstery

☒ A9000801050

Colour Back equipment and Seat plate

Information on current colour designs and availability at www.proactiv-gmbh.com

If not specified, the seat and backrest system will be delivered in the standard colour anthracite matt.

anthracite RAL 7016 gloss

- ☐ A9000100179-995
☐ Ergonomic back shell
☐ ERGO Back backrest
☐ Seat plate aluminum

in frame colour

- ☐ A9000100179-999 € 92.30
☐ Ergonomic back shell
☐ ERGO Back backrest
☐ Seat plate aluminum

For frame coating in two-colour design


Please specify if the desired option should be coated in the base colour or the design colour.

- ☐ A9000100179-999 € 92.30
☐ Ergonomic back shell
☐ base colour
☐ design colour
☐ ERGO Back backrest
☐ base colour
☐ design colour
☐ Seat plate aluminum
☐ base colour
☐ design colour

If a different colour is required, please select it from the wheelchair colour palette:


- ☐ A9000100179-998 € 149.40
 Please enter the last 3 digits of the desired colour:
 A9000100170-
☐ Ergonomic back shell
☐ ERGO Back backrest
☐ Seat plate aluminum

Seat cushion



Seat cushion in rectangular contour Inlet: Standard

- ☐ A9000805180-K* € 118.30



Seat cushion with radius contour Inlet: Standard

- ☐ A9000800300-K* € 163.10

Seat cushion height:

- ☐ 3 cm ☐ 5 cm
☐ 4 cm ☐ 6 cm

Seat cushion in rectangular contour Inlet: Anti-Decubitus

- ☐ up to 30x30 cm A9000805204-K* € 301.70
☐ >30x30 up to 50x50 cm A9000805200-K* € 437.70

Seat cushion with radius contour Inlet: Anti-Decubitus

- ☐ up to 30x30 cm A9000800311-K* € 356.10
☐ >30x30 up to 50x50 cm A9000800310-K* € 482.90

Seat cushion height:

- ☐ 4 cm ☐ 6 cm

Cover:

- ☐ Black, breathable, open surface
☐ Black, breathable, closed surface

If not specified, the cover material „Black, breathable, open surface“ will be delivered.

- ☐ Seat cushion is provided

Clothing guard

Aluminum clothing guard

fixed to the side panel
 Only in combination with fixed backrest



in frame colour (base colour for two-colour design)


- ☐ A9000420101-K € 311.10

Selection from colour palette „One-colour design“

- ☐ A9000420103-K € 454.40

Colour: _____

gloss level: ☐ gloss ☐ matt

Aluminum clothing guard removable via lock function adjustable following wheel arch Only in combination with folding backrest			
in anthracite matt ☐ A9000420215-K € 243.20	in frame colour (base colour for two-colour design) ☐ A9000420216-K € 311.10	in design colour (only possible with frame colour in two-colour design) ☐ A9000420217-K € 345.10	
Carbon clothing guard removable via lock function ☐ A9000420013-K € 547.30 adjustable following wheel arch Only in combination with folding backrest 			
Knee lever brake system, brake lever integrated in the clothing guard see „Brake systems”.			
Armrests			
Armrest, height adjustable, integrated in the clothing guard Including aluminum clothing guard with arm padding 25 cm long Removable via lock function Only in combination with folding backrest			
in anthracite matt ☐ A9000420014-K € 569.60	in frame colour (base colour for two-colour design) ☐ A9000420015-K € 621.90	in design colour (only possible with frame colour in two-colour design) ☐ A9000420016-K € 648.10	
Armrest, height adjustable, integrated in the clothing guard ☐ A9000420017-K € 862.00 Including carbon clothing guard with arm padding 25 cm long Removable via lock function Only in combination with folding backrest			
Axle plate Standard equipment with 12 mm sockets / axles			
Axle plate ☒ A9000220000 quick-release axle Ø 12 mm 	Extended bike plate ☐ A9000220030 € 315.20 adjustable track 		
Wheel camber *Total width with tight-mounted handrim.			
1° drive wheel camber ☒ A9030202010 Total width = approx. SB + 20 cm*	2,5° drive wheel camber ☒ A9030202020 Total width = approx. SB + 23 cm*	4° drive wheel camber ☒ A9030202030 Total width = approx. SB + 26 cm*	
6° drive wheel camber ☒ A9030202040 Total width = approx. SB + 31 cm*	8° drive wheel camber ☐ A9041202080 € 77.00 Total width = approx. SB + 35 cm*	10° drive wheel camber ☐ A9041202100 € 136.90 Total width = approx. SB + 39 cm*	

E-Adaptations

Preparation for the „e-motion“ / „e-motion DuoDrive“ mounting

☐ A9000220200 € 502.90

Equipment: e-motion torque arm, 1/2" axle and socket length 50.8 mm (2")
PRO ACTIV and Alber recommend mounting two anti-tipping supports. Please select under „Adaptations“ or order Alber anti-tipping supports directly from Alber. Not in combination with 10° drive wheel camber.

Preparation for the „e-fix“ mounting

☐ A9000220210 € 502.90

Equipment: e-fix torque arm, 1/2" axle and socket length 50.8 mm (2")
PRO ACTIV and Alber recommend mounting two anti-tipping supports. Please select under „Adaptations“ or order Alber anti-tipping supports directly from Alber, as well as recording BUDDY. Only in combination with 1° or 2.5° drive wheel camber.

Brake system

If not specified, the „Brake system push to lock with standard brake handle“ will be delivered.

Recommendation: For KP 7 cm or lower, we recommend the „Brake system push to lock with brake handle 20° outwards“.

Knee lever brake with standard brake lever

☐ A9000701010-001



Knee lever brake with long brake lever

☐ A9000701010-002



Knee lever brake system Brake lever angled out 20°

☐ standard brake lever A9000700002

☐ long brake lever A9000701010-004



Brake lever with M8 adapter and plastic ball

☐ A9000706001 € 170.50

only in combination with knee lever brake system, standard or long brake handle



Brake system push to lock Brake handle to flip-down

☐ A9000706008 € 158.60
length: 9 cm



Knee lever brake with one-hand operation

☐ Operation right side A9000707020-001 € 392.30

☐ Operation left side A9000707020-002 € 392.30
via Bowden cable



Knee lever brake system, brake lever integrated in the clothing guard

Brake unit attached to the lower frame tube.
Actuation via cable pull.
Clothing guard cannot be removed.



glossy black for fixed backrest
matt anthracite for foldable backrest

☐ A9041720015-K € 643.00

in frame colour (base colour for two-colour design)

☐ A9041720016-K* € 710.30

in design colour (only possible with frame colour in two-colour design)

☐ A9041720017-K* € 744.30

Knee lever brake with lock

☐ Standard brake lever A9000700214 € 158.60

☐ Long brake lever A9000700215 € 158.60

Locking block for brake lever
Prevents braking during pushing by accompanying persons
Not in combination with knee lever brake with one-hand operation, brake lever angled out 20° as well as parking knee lever brake for accompanying person



Parking knee lever brake for accompanying person

☐ A9000701020 € 464.10

Push handle operation via Bowden cable, not in combination with 1° wheel camber



Drum brake

☐ A9000720011-K* € 474.00

Operating from the push handle
Total width of wheelchair (see selection „Wheel camber“) + 3 cm
Please note:
- Only in combination with drive wheel aluminum 36-spoke rim wheel
- not in combination with one-hand drive, e-fix or e-motion
- not in combination with 1° wheel camber



Footrest

If not specified, the „continuous footrest“ will be delivered.

Continuous footrest

☐ A9041601010-K*
with tubular resting bracket
height adjustable
possible from lower leg length USL 26 cm to
USL 47 cm
If you want to equip the footrest with a „rubber
coated footplate“ please select this under
„Footrest accessories“.



Footrest for short lower legs

☐ A9041606010-K
height adjustable
From lower leg length USL 20 cm to 25 cm.
Rubber coated footplate necessary, please
select this under „Footrest accessories“.




Footrest with variable positioning

☐ A9047600001-K **€ 281.10**
adjustable in height, depth and angle
including footplate cover with anti-slip film
up to lower leg length USL 37 cm



Footrest folds up to the rear with spring locking mechanism

☐ A9047602005-K* **€ 347.20**
with tubular resting bracket
up to lower leg length USL 32 cm
including caster wheel track extension 20 mm
per side
If you want to equip the footrest with a „rubber
coated footplate“ please select this under
„Footrest accessories“.



Footrest with variable positioning and folding up to the rear

☐ A9047600002-K **€ 373.20**
adjustable in height, depth and angle
including footplate cover with anti-slip film
up to lower leg length USL 37 cm



Footrest folding up to one side

☐ A9041601300-K* **€ 173.90**
from lower leg length USL 25 cm:
with tubular resting bracket
height adjustable +/- 3 cm
cannot be removed
from lower leg length USL 29 cm:
can also be swivelled outwards
If you want to equip the footrest with a „rubber
coated footplate“ please select this under
„Footrest accessories“.



Divided footrest

☐ A9041603200-K **€ 336.80**
with rubber coated footplate,
height adjustable
can be swivelled outwards and folded up to the
side
cannot be removed
possible from lower leg length USL 29 cm and
seat width SB 26 cm



Footrest accessories

Rubber coated footplate, 14 cm deep

☐ A9000610000-K **€ 125.90**
bolted, angle-adjustable



Rubber coated footplate, 18 cm deep

☐ A9000610005-K **€ 140.10**
bolted, angle-adjustable



Footplate side stop

☐ A9000610100 **€ 59.70**
adjustable
Only in combination with rubber coated
footplate



Deflector rollers mounted on the footplate

☐ A9000600006 **€ 117.40**
Only in combination with rubber coated
footplate



Adaptations

Anti-tipping support

☐ right side
A9041301160-K001 € 199.00

☐ left side
A9041301160-K002 € 199.00

with tension spring
can be folded up
removable via Quick-pin



Outdoor wheel attachment with steering fork

☐ A9041900100 € 845.70

equipped with straight-running lock,
can be removed from the wheelchair via quick
release lever.

With swivelling pneumatic tyre 12.1/2 x 2.1/4"
caster wheel that can be removed from the
support tube via thru axle



3" caster wheel

3" skate wheel (80 mm)

☐ A9000501110 € 70.80

puncture-proof



4" caster wheel

4" caster wheel „Standard“ solid rubber black

☒ A9000501011-004
(82A), puncture-proof



4" skate wheel (100 mm)

☐ A9000501120 € 70.80

puncture-proof



4" safety light-up wheel (100 mm)

☐ A9000501220 € 94.80

puncture-proof



4" CorMORun caster wheel (100 mm)

☐ A9000501042 € 173.00

puncture-proof
Rim colour: Black
Incl. hair guard spacer (5.3 mm)



4" Frog Legs caster wheel (100 mm)

☐ A9000501230 € 206.30

with replaceable rubber 6 mm,
puncture-proof
Colour: black, silver rim
Incl. hair guard spacer (5.3 mm)



5" caster wheel

5" caster wheel „Standard“ PU black

☒ A9000501036-005
(82A), puncture-proof



5" caster wheel „Standard“ PU grey

☒ A9000501031-005
(82A), puncture-proof



5" safety light-up wheel (125 mm)

☐ A9000501221 € 94.80


puncture-proof



5" aluminum rim 3-spoke design „good shock absorption“ ☐ silver anodised A9000501210 € 179.80 ☐ black A9000501211 € 179.80 with replaceable rubber 6 mm medium hard, puncture-proof Incl. hair guard spacer (3 mm) 	5" CorMORun caster wheel (125 mm) ☐ A9000501041 € 181.30 puncture-proof Rim colour: Black Incl. hair guard spacer (5.3 mm)  	5" Frog Legs caster wheel (125 mm) ☐ A9000501232 € 206.30 with replaceable rubber 6 mm, puncture-proof Colour: black, silver rim Incl. hair guard spacer (5.3 mm)  
6" caster wheel		
6" caster wheel „Standard“ PU grey ☒ A9000501010-006 (82A), puncture-proof 	6" caster wheel „pneumatic tyre“ ☐ A9000501050 € 57.20 Front seat height min. 42 cm 	
Caster fork For caster fork colours, see „Colour of add-on parts“ If not specified, the „Standard caster fork with screw axle“ will be delivered.		
Caster fork Standard ☒ A9000500100-K001 with screwed axle ☐ A9000500100-K002 € 167.40 with quick-release axle 	Caster fork Heavy Duty for high loads or with the use of additional electric drives ☐ A9000510111-K001 € 54.00 with screwed axle ☐ A9000510111-K002 € 218.90 with quick-release axle 	Suspension fork incl. 5" aluminum rim „good shock absorption“ ☐ A9000500160-K002 € 616.90 with quick-release axle 
Caster fork spare wheel „Standard“ ☐ A9000500300-K € 296.70 with quick-release axles 	Caster fork spare wheel „Heavy Duty“ ☐ A9000500301-K € 337.90 with quick-release axles 	The spare caster fork does not have a caster wheel. In case a caster wheel is requested, please put a cross on a second selection in the category „Caster wheels“.

Track distance

Caster wheel track extension



☐ A9000500400*
10 mm per side

€ 59.10

☐ A9000500401*
15 mm per side

€ 59.10

☐ A9000500410* € € 59.10
20 mm per side
included in „Footrest folds up to the rear with
spring locking mechanism“
item no. A9047602005-K*

Drive wheels

If not specified, the „drive wheel 36 spoke aluminum rim“ will be delivered.

Drive wheel 18 spoke aluminum rim

☐ 20" (ETRTO 451-19)
A9000203052

☐ 22" (ETRTO 489-19)
A9000203062

☐ 24" (ETRTO 540-19)
A9000203122

with Liteline hub, radially spoked

not in combination with:

- Puncture-proof tires
- One-hand drive
- Drum brake
- Transport restraint connection
- Handbike

up to 70 kg body weight:

Sturz	1°	4°	6°	8°	10°
Gewicht					
bis 50 Kg	✓	✓	✓	✓	✓
51 - 60 Kg	✓	✓	✓	✓	✓
61 - 70 Kg	✓	✓			



Drive wheel 36 spoke aluminum rim

☐ 22" (ETRTO 489-19)
A9000203060-010

☐ 24" (ETRTO 540)
A9000203000-010
with Liteline hub (left image),
radially spoked

☐ 20" (ETRTO 451-19)
A9000203050

☐ 22" (ETRTO 451-19)
A9000203060-020

☐ 24" (ETRTO 540)
A9000203000-020
with aluminum front hub (right image),
radially spoked
in combination with one-hand drive



Drive wheel double rim 36-spokes

☐ 24" (ETRTO 540)
A9000203100 € 157.80

with aluminum front hub,
radially spoked
for use with high loads



Spinergy wheel „Spox Sport“ with 24 PBO spokes; rim, hub and spokes in black

☐ 24" (ETRTO 540)
A9000203550 € 784.10



Spinergy „Everyday Light Extreme“ with 12 PBO spokes; rim and hub in black, spokes in white

☐ 24" (ETRTO 540)
A9000203554 € 957.20

Recommendation: For heavy loads
(e.g. higher body weight), we recommend the
Spinergy „Spox Sport“ wheel with 24 PBO
spokes.



Spinergy „Everyday Light Extreme“ with 12 PBO spokes; rim, hub and spokes in black


☐ 24" (ETRTO 540)
A9000203552 € 957.20

Recommendation: For heavy loads
(e.g. higher body weight), we recommend the
Spinergy „Spox Sport“ wheel with 24 PBO
spokes.



Handrims

x = dimension handrim attachment



Handrim attachment
general:

☐ x = 11 mm

☐ x = 19 mm

Handrim attachment
Curve + Gekko handrim:

☐ x = 10 mm

☐ x = 18 mm

If not specified, the handrim is mounted tightly (11 mm / 10 mm).

Stainless steel Ø 20 mm

☒ 20" handrim for
20" wheel
A9000205010-001

☒ 22" handrim for
22" wheel
A9000205010-002

☒ 24" handrim for
24" wheel
A9000205010-004



Aluminum Ø 19 mm silver anodised

☒ 22" handrim for
22" wheel
A9000205019-002

☒ 24" handrim for
24" wheel
A9000205019-004



Aluminum Ø 16 mm powder coated in frame colour

(basic colour with two-colour design)

☐ = Gloss surface

☐ = Matt surface

Note: Matt coating is not recommended due to the low grip.

☐ 20" handrim for
20" wheel
A9000205021-001 € 68.30

☐ 22" handrim for
22" wheel
A9000205021-002 € 68.30



Aluminum Ø 19 mm powder coated in frame colour

(basic colour with two-colour design)

☐ = Gloss surface

☐ = Matt surface

Note: Matt coating is not recommended due to the low grip.

☐ 20" handrim for
20" wheel
A9000205020-001 € 68.30

☐ 22" handrim for
22" wheel
A9000205020-002 € 68.30

☐ 24" handrim for
24" wheel
A9000205020-004 € 68.30



Aluminum Ø 16 mm polished and hard-anodised

☐ 20" handrim for
20" wheel
A9000205031-001 € 164.90

☐ 22" handrim for
22" wheel
A9000205031-002 € 164.90




Aluminum Ø 19 mm polished and hard-anodised

☐ 20" handrim for
20" wheel
A9000205030-001 € 164.90

☐ 22" handrim for
22" wheel
A9000205030-002 € 164.90

☐ 24" handrim for
24" wheel
A9000205030-004 € 164.90



Titanium Ø 19 mm

☐ 24" handrim for
24" wheel
A9000205040-004 € 335.60



Silicone handrim cover

☐ for handrim 20"
A9000205100-001 € 232.80

☐ for handrim 22"
A9000205100-002 € 232.80

☐ for handrim 24"
A9000205100-004 € 232.80

Surface without structure,
compatible with handrim Ø 19 mm or Ø 20 mm.



Ultra-Grip „very grippy“

☐ 20" handrim for
20" wheel
A9000205070-001 € 402.40

☐ 22" handrim for
22" wheel
A9000205070-002 € 402.40

☐ 24" handrim for
24" wheel
A9000205070-004 € 402.40



Ergo-Grip ☐ 20" handrim for 20" wheel A9000205080-001 € 402.40 ☐ 22" handrim for 22" wheel A9000205080-002 € 402.40 ☐ 24" handrim for 24" wheel A9000205080-004 € 402.40 Attachment fitting to the rim Note: Due to the complex geometry, concentricity deviations of +/- 2-3 mm are possible. 	Maxgrepp with cellular and natural rubber coating very grippy, soft ☐ 20" handrim for 20" wheel A9000200788 € 402.40 ☐ 22" handrim for 22" wheel A9000200789 € 402.40 ☐ 24" handrim for 24" wheel A9000200787 € 402.40 	Curve black anodised, ergonomic handrim profile ☐ 24" handrim for 24" wheel Profile height 24 mm A9000205130-004 € 402.40 ☐ 24" handrim for 24" wheel Profile height 31 mm Size L A9000205140-004 € 402.40 
Curve silver anodised ergonomic handrim profile ☐ 24" handrim for 24" wheel Profile height 31 mm Size L A9000205160-004 € 402.40 	Curve powder-coated, ergonomic handrim profile ☐ = frame colour (gloss surface) ☐ = black (gloss surface) ☐ 24" handrim for 24" wheel Profile height 24 mm A9000205135-004 € 467.60 ☐ 24" handrim for 24" wheel Profile height 31 mm Size L A9000205145-004 € 467.60 	Curve Grip, with Tetra-Grip coating, ergonomic handrim profile ☐ 24" handrim for 24" wheel Profile height 24 mm A9000205155-004 € 520.70 
Quadro with Tetra-Grip coating ☐ 24" handrim for 24" wheel A9000205170-004 € 520.70 	Gekko black anodised ☐ 24" handrim for 24" wheel A9000205175-004 € 520.70 with a robust and replaceable silicone lip, fully integrated into the profile Profile height 31 mm 	Gekko silver anodised ☐ 24" handrim for 24" wheel A9000205185-004 € 520.70 with a robust and replaceable silicone lip, fully integrated into the profile Profile height 31 mm 
Gekko S black anodised ☐ 22" handrim for 22" wheel A9000205195-002 € 520.70 ☐ 24" handrim for 24" wheel and 25" wheel A9000205195-004 € 520.70 with a robust and replaceable silicone lip, fully integrated into the profile Profile height 24 mm 	Equipment without handrims ☒ A9000205990	

Spoke guard

Finger- and spoke guard „transparent“

Not in combination with Spinergy drive wheels



☐ for 20" wheel
A9000206300-001

€ 207.10

☐ for 22" wheel
A9000206300-002

€ 207.10

☐ for 24" wheel
A9000206300-004

€ 207.10

Tyres

If not specified, the „Schwalbe RightRun in black“ tyres will be delivered.

Schwalbe „RightRun“ High-pressure tyres

☐ 22x1" (489-25)
A9000204001-002

☐ 24x1" (540-25)
A9000204001-004

Colour: Black
7.5 bar



Schwalbe „RightRun“ High-pressure tyres

☐ 20x1" (451-25)
A9000204000-001

☐ 22x1" (489-25)
A9000204000-002

☐ 24x1" (540-25)
A9000204000-004

Colour: Grey, black wall
7.5 bar



Schwalbe „DownTown“ Profiled tyres

☐ 24x1" (540-25)
A9000204010-004
7.5 bar



Schwalbe „Marathon Plus Evolution“ with „smart guard“ puncture protection and reflector strips

☐ 22x1" (489-25)
A9000204030-002 € 49.20

☐ 24x1" (540-25)
A9000204031-004 € 49.20

Colour: Black
10 bar



Schwalbe „Marathon Plus Evolution“ with „smart guard“ puncture protection

☐ 24x1" (540-25)
A9000204030-014 € 49.20

Colour: Black
10 bar



Schwalbe „MarathonPlus“ with „smart guard“ puncture protection

☐ 24x1" (540-25)
A9000204030-004 € 49.20

Colour: Grey, black wall
10 bar



Schwalbe „One“

☐ 24x1" (540-25)
A9000204123-004 € 93.00

Colour: Black
9 bar



Puncture-proof „solid rubber tyres“

☐ 22x1" (489-25)
A9000204200-002 € 122.70

Colour: Grey

☐ 24x1" (540-25)
A9000204201-004 € 122.70

Colour: Black
Not in combination with 18 spoke or Spinergy
drive wheels.



Accessories

24" MTB spare wheel set

with double rim (ETRTO 507), aluminum front hub, quick-release axle and stainless steel handrim
Mountain bike tyres, selectively:

☐ A9000203410 € 703.40

24x1.75" Schwalbe "Marathon Plus"
Brake setting compatible with 24x1" wheel



☐ A9000203413 € 703.40

24x2.1" Schwalbe "SMART SAM Performance"
Brake setting not compatible with 24x1" wheel




NOTE: We charge a surcharge for the use of drum brakes

☐ A9000200004 € 232.20

Through axle with assistance, Large push button

☐ A9000200096 € 105.80

Interlock dimension: 122 / 110 mm
Axle diameter: 12 mm
Not in combination with one-hand drive



Hand air compressor with pressure gauge

☐ A9000901100 € 106.30

up to 14 bar



Crutch holder with flexible nylon cup

☐ A9000902310* € 137.40



Length-adjustable safety belt with buckle

☐ A9000902100 € 80.60



Adjustable safety belt with hook and fastener

☐ A9000902110 € 64.50



SIGMA lighting set

☐ A9000903270 € 107.90

Consisting of front light and rear lights,
chargeable via Micro USB socket,
incl. Micro USB charging cable
Lighting set compliant with road traffic
regulations



Transport bag for drive wheels

☐ for 24" wheel
A9000903530-004 € 77.80
Ø 76 cm x 2 cm



Therapy table

☐ width 350 mm
(up to SB 37 cm) € 480.20
A9000902800-001

☐ width 400 mm
(SB 38 - 41 cm) € 480.20
A9000902800-002

☐ width 420 mm
(from SB 42 cm) € 480.20
A9000902800-003

removable table plate from Lexan, transparent,
with front edging, with 2 positioning stops
Ø 25 mm (for adaptation to frame tube
diameter Ø 25 mm)



Transport restraint system

Transport restraint connection on rear frame integrated with belt eye

Frame weld-on part for mounting the lap belt and for holding the rear retractor.
If required, please select lap belt with lock (A9000900001 / A9000900007) separately.

☐ A9041910001

€ 365.40



Bearing block intermediate plate

with belt eye for holding the front retractor.

Please note: The caster wheel track becomes 12 mm wider on each side.

☐ Colour: Black
A9000900008-001

€ 181.80

☐ Colour: Silver
A9000900008-002

€ 181.80



Lap belt 1350 mm with lock, connection to transport restraint system

☐ A9000900001 € 169.10

adjustable belt length. With belt links on both sides for mounting the diagonal shoulder belt. Belt colour: Red

Lap belt 1670 mm with lock, connection to transport restraint system


☐ A9000900007 € 169.10

adjustable belt length. With belt links on both sides for mounting the diagonal shoulder belt. For use on wheelchairs with a seat shell. Belt colour: Red

Headrest with segment upholstery

☐ A9000900004 € 415.60

3-part and support tube Ø 25 mm



Headrest holder Ø 25 mm

☐ A9000900005-K001 € 279.60

to accommodate the headrest with segment padding on the back cross bar. Reinforced version with double cross-tube.



Headrest holder with square mount for the support headrest

☐ Tube dimensions of support tube 15 x 15 mm
A9000400850-001 € 289.80

☐ Tube dimensions of support tube 12,7 x 12,7 mm
A9000400850-002 € 289.80

Mounted on the back crossbar. Reinforced version with double cross tube. From SB 23 cm possible. Not in combination with central push handle.

Separate order sheets

Category „Rehabilitation technology accessories“

e.g. Seat & backrest systems, one-hand drive or transport restraint systems can be found at:
www.proactiv-gmbh.com/downloads/

Additional products

such as hand bikes or traction & attachment devices can be found at:
www.proactiv-gmbh.com/downloads/

Additional sheets

such as short abduction beginning etc. can be found at:
www.proactiv-gmbh.com/downloads/

Services

Consultation, clarification of the installation situation and measurement through PRO ACTIV

☐ A9990000000 € 354.50

Handover and instruction on the wheelchair by PRO ACTIV

☐ A9990000006 € 354.50

Notes and remarks:

Location, Date

Signature End customer/dealer

The signature field is intended for the consultant and insured person. This confirms that the intended equipment is intended for the insured person.

- 19 - 2026-07-27 / The prices are plus statutory value-added tax and shipping costs ☐ = included in the basic price ☒ = Position selected
Order sheet BUDDY Rev. 02 UK