



**Quotation date:** \_\_\_\_\_  
**Consultation was carried out by**  
**Staff at the medical supplies shop:** \_\_\_\_\_

**Enquiry no.:** \_\_\_\_\_

**Order date:** \_\_\_\_\_  
according to GTC PRO ACTIV Reha-Technik GmbH

**Order no.:** \_\_\_\_\_

### Medical suppliers / company address:

**Cust. no.:** \_\_\_\_\_  
**Name:** \_\_\_\_\_  
**Street:** \_\_\_\_\_  
**Postcode / city:** \_\_\_\_\_  
**Phone:** \_\_\_\_\_  
**Fax:** \_\_\_\_\_  
**E-Mail:** \_\_\_\_\_  
**Contact:** \_\_\_\_\_

**Different invoicing address (see below)**

**Cust. no.:** \_\_\_\_\_  
**Name:** \_\_\_\_\_  
**Street:** \_\_\_\_\_  
**Postcode / city:** \_\_\_\_\_

**Different shipping address (see below)**

**Cust. no.:** \_\_\_\_\_  
**Name:** \_\_\_\_\_  
**Street:** \_\_\_\_\_  
**Postcode / city:** \_\_\_\_\_

### Customer name / commission:

**ID / Comm. No.:** \_\_\_\_\_  
**Surname, forename:** \_\_\_\_\_  
**Street:** \_\_\_\_\_  
**Postcode / city:** \_\_\_\_\_  
**Phone:** \_\_\_\_\_  
**E-Mail:** \_\_\_\_\_  
**D.o.B.:** \_\_\_\_\_  
**Paying organisation:** \_\_\_\_\_  
**Indication:** \_\_\_\_\_

## SPIKE adaptive bike

up to 75 kg trailer load and 5 kg payload

- docking plate for quick coupling to the „Quick-Snap“ adapter
- drive unit with steering fork
- A-shaped handles with integrated switching and brake fittings
- chain guard in black on both sides
- without gear shift
- bottom bracket support can be adjusted in height and angle
- parking stand
- Sigma lighting set with rechargeable battery charger, front and rear lights

A911400000 € 2,678.00

### basic equipment:

- Freewheel hub with singlespeed kit, 15 teeth
- Hydraulic disc brake with Ø 180 mm disc
- Mech. Shimano disc brake, brake handle left

### Body size

Body height: \_\_\_\_\_ cm

Upper body height: \_\_\_\_\_ cm

Weight: \_\_\_\_\_ kg

This information is necessary for checking.

## Fork colour

Information on current colour designs and example images can be found at: <https://www.proactiv-gmbh.com/colours/>

### white RAL 9016

- gloss A9000100170-023
- matt A9000100170-123      **€ 43.40**

### brilliant silver

- gloss A9000100170-011

### anthracite RAL 7016

- gloss A9000100170-001
- matt A9000100170-101      **€ 43.40**

### black RAL 9005

- gloss A9000100170-003
- matt A9000100170-019      **€ 43.40**

### deep black effect silver

- gloss A9000100170-007

### pure orange RAL 2004

- gloss A9000100170-022

### fire red

- gloss A9000100170-018
- matt A9000100170-118      **€ 43.40**

### pink magenta RAL 4010

- gloss A9000100170-037
- matt A9000100170-137      **€ 43.40**

### sky blue

- gloss A9000100170-017

### pastel blue RAL 5024

- gloss A9000100170-036
- matt A9000100170-136      **€ 43.40**

### navy blue

- gloss A9000100170-045
- matt A9000100170-145      **€ 43.40**

### electric green

- gloss A9000100170-015
- matt A9000100170-115      **€ 43.40**

### lemon green

- gloss A9000100170-005
- matt A9000100170-105      **€ 43.40**

### violet

- gloss A9000100170-044

### yellow RAL 1023

- gloss A9000100170-046
- matt A9000100170-146      **€ 43.40**

### red RAL 3001

- gloss A9000100170-047
- matt A9000100170-147      **€ 43.40**

### lagoon blue

- gloss A9000100170-052
- matt A9000100170-152      **€ 43.40**

### athletic grey RAL 7046

- gloss A9000100170-050
- matt A9000100170-150      **€ 43.40**

### muddy grey RAL 7033

- gloss A9000100170-051
- matt A9000100170-151      **€ 43.40**

### peach RAL 3022

- gloss A9000100170-038
- matt A9000100170-138      **€ 43.40**

## Equipment

### Shimano Deore 1x11 Rapid Fire 11-speed derailleur gears

- A9114200032      **€ 403.70**
- mechanical derailleur gear system with control on the handle.  
11-speed cassette 11-42 teeth

### Shimano „Nexus“ 8-speed hub gear „Grip Shift“ twist grip shifter

- A9114000008      **€ 718.70**
- with V-brake rim brake, brake lever on the left



### Mountaindrive 2.5 : 1 Hill reduction in combination with hub gear


- A9100200003      **€ 797.40**
- planetary gearbox for doubling the gears.  
All gears (the hub gear and the mountain gear) can be shifted while stationary!  
In combination with grip width 40 cm or 45 cm



## Dimensions

### Bottom bracket shaft with standard crank

- grip width 36 cm  
A9114500100-036
- grip width 40 cm  
A9110500100-040
- grip width 45 cm  
A9110500100-045



### Standard pedal crank

- length 135 mm  
A9110600100-135
- length 155 mm  
A9110600100-155
- length 175 mm  
A9110600100-175

## Tyres

### 20x2" Schwalbe „Big Apple“

A9114105125  
(ETRTO 406-50)  
balloon tire



### Schwalbe „Marathon Plus“

20x1.75" (406-47)  
A9130105128  
unpuncturable, with anti-puncture belt  
Smart Guard




### Puncture protection band for pneumatic tyres

A9000206100-K € 55.80  
Only in combination with Schwalbe „Big Apple“.

## Wheelchair-dependent accessories

### Extended bike plate

- for SPEEDY A1 / A2  
A9110100100-002 € 370.00
- for BUDDY  
A9110100100-003 € 370.00
- for TRAVELER  
A9110100100-001 € 370.00  
4° camber integrated  
Track independently adjustable



### Wheelbase extension via „grid mechanism with pulling operation“

- for SPEEDY 4all /  
4teen / 4you /  
TRAVELER 4you € 564.70  
A9110102001  
removable, with continuous axle tube, track  
width independently adjustable
- Please select the wheel camber:
- 1° drive wheel camber
- 4° drive wheel camber
- 6° drive wheel camber



## Frame adapter

### SPEEDY rigid frame adapter left-hand operation with „Quick-Snap-Adaption“

- A9110700141 € 537.20  
with quick-release clamp and  
centering pin holder  
Colour: black

Please select accordingly:

**Centring pin for Ø 25 mm**  
for A1 and A2

- A9110700022

**Centring pin for Ø 40 mm**  
for 4all, 4all Ergo, 4you Ergo, F2, F3 and 4teen

- A9110700023



### TRAVELER folding frame adapter with „Quick-Snap-Adaption“

- A9110701300 € 537.20  
with quick-release clamp and  
centering pin holder on the folding mechanism  
Colour: black



### BUDDY rigid frame adapter with „Quick-Snap-Adaption“

- A9110701400 € 537.20  
with quick-release clamp and  
centering pin mount on the frame  
Colour: black



**Wheelchair front adapter with universal clamp and „Quick-Snap-Adaption“**

A9110701500 € 572.30

Adapter attachment to the wheelchair stem. Operation on the right-hand side.

Mounting situation for PA products:  
Straight front end from SB 31 cm to 54 cm  
V-shaped front end from SB 30 cm to 50 cm

**Please select the frame clamp below:**

clamping range 95 mm  
A9110701525  
**Tube diameter 23-27 mm**

clamping range 95 mm  
A9110701530  
**Tube diameter 28-32 mm**

clamping range 95 mm  
A9110701535  
**Tube diameter 33-37 mm**

**PLEASE NOTE:**

If the exact tube diameter is not specified, an adapter for a tube diameter of 23 - 27 mm will be delivered as standard.

**Please select the clamping sleeve below:**

clamping sleeve [80 mm] included as standard  
for close clearance between footrest and drive wheel.

clamping sleeve [110 mm] € 46.70  
A9199000117  
for medium distance between footrest and drive wheel.

clamping sleeve [170 mm] € 56.60  
A9199000154  
for a wide distance between footrest and drive wheel.

**PLEASE NOTE:**

If not specified, the clamping sleeve with 80 mm length will be delivered.

**Combination adapter PA-Triride, steerable Quick-Snap version**

A9199000083 € 659.70

with Quick-Snap Adaptation for adapting PRO ACTIV handbikes and pulling devices as well as Triride products that are equipped with the PRO ACTIV adapter attachment package (A9140700000).

Frame clamps are NOT included in the delivery. Please select these separately.



**Please select the frame clamp below:**

Ø 25 mm E8140BS-C-106 € 205.20

Ø 30 mm E8140BS-C-107 € 205.20

Ø 32 mm E8140BS-C-108 € 205.20

**Clamp lever with 3-way star handle for actuating the frame clamping**

A9199701507 € 106.10

mounted on the right and left for easy terminal operation with limited manual function  
Only in combination with front adapter



**Accessories**

**Plug-in mudguard**

made of plastic  
including necessary mounting components

A9122502100 € 109.40

**In combination with luggage rack**

A9122502101 € 58.10



**Aluminum mudguard**

including necessary mounting components

A9122500005 € 153.50

**In combination with luggage rack**

A9122500006 € 112.40

**SIGMA BC 8.0 ATS bike computer**

A9110901000 € 70.10  
wireless, with mount



**Chain guard disc in frame colour**

A9114500005 € 85.80

**Aluminum luggage rack „cargo systems“**

A9120900000 € 210.00

Colour: black  
max. load weight on the loading platform 10 kg,  
max. 15 kg in combination with luggage bags  
at lower loading height.



**Maneuvering castors for parking stand**

A9110900017 € 122.10

Maneuvering castors mounted on the parking stand via inlet



## Options for users with reduced functions

### Tetra hand support pivotable and adjustable grip distance

A9110601101 € 169.40



### Tetra shell handles

A9110601100 € 172.40



### Adaptation help for centering pin

A9110700002 € 55.00  
to simplify the insertion of the central adapter.



### Lever extension excenter lock for Quick Snap connection

Makes it easier to open and close the Quick Snap adapter. Tetra version

A9110700005 € 73.90

straight design  
only in combination with rigid frame adapter

A9110700006 € 73.90

curved design  
only in combination with front adapter



### Chest strap for upper body fixation

size S (105 cm) A9000800016 € 241.80

size M (125 cm) A9000800013 € 241.80

size L (145 cm) A9000800014 € 241.80

size XL (165 cm) A9000800017 € 241.80

consisting of elastic webbing material in combination with neoprene with Velcro and operating tab.  
Colour: black



## Services

### Consultation, clarification of the installation situation and measurement through PRO ACTIV

A9990000000 € 354.50

### Handover and instruction of the handbike / traction device by PRO ACTIV

A9990000011 € 354.50

### Attachment of the handbike / traction device to the wheelchair by PRO ACTIV at specialist dealer

#### on a PRO ACTIV wheelchair

A9990000012 € 230.00

product no.: \_\_\_\_\_

#### on a third-party product

A9990000013 € 366.90

manufacturer: \_\_\_\_\_

model: \_\_\_\_\_

frame tube Ø \_\_\_\_ mm

### Attachment of the handbike / traction device to the wheelchair through specialist dealer

including mounting of the frame adapter

#### on a PRO ACTIV wheelchair

A9990000008

product no.: \_\_\_\_\_

#### on a third-party product

A9990000005

manufacturer: \_\_\_\_\_

model: \_\_\_\_\_

frame tube Ø \_\_\_\_ mm

**Notes and remarks:**

---

**Location, Date**

**Signature End customer/dealer**

The signature field is intended for the consultant and insured person. This confirms that the intended equipment is intended for the insured person.  
- 6 - 2026-03-03 / The prices are plus statutory value-added tax and shipping costs     = included in the basic price     = Position selected  
Order sheet SPIKE adaptive bike Rev. 07 UK