

Im Hofstätt 11
D-72359 Dotternhausen
Tel +49 (0)7427 9480-0
Fax +49 (0)7427 9480-7025
e-mail info@proactiv-gmbh.de
www.proactiv-gmbh.de



LIFT solid with mechanical seat height adjustment

The desired work height is set by means of the user-friendly, manual and continuous adjustment of the seat height, with the help of gas pressure springs. With just a little muscle power, the seat lifts up to 30 cm.

Maximum user weight limit 120 kg

A9060100000* € 7,860.00

Quotation date: _____
Consultation was carried out by
Staff at the medical supplies shop: _____
Enquiry no.: _____

Order date: _____
according to GTC PRO ACTIV Reha-Technik GmbH
Order no.: _____

Medical suppliers / company address:

Cust. no.: _____
Name: _____
Street: _____
Postcode / city: _____
Phone: _____
Fax: _____
E-Mail: _____
Contact: _____

Shipping address: same as company address
 Other (see below)

Name: _____
Street: _____
Postcode / city: _____

Customer name / commission:

ID / Comm. No.: _____
Surname, forename: _____
Street: _____
Postcode / city: _____
Phone: _____
Fax: _____
E-Mail: _____
D.o.B.: _____
Paying organisation: _____


Invoicing address: same as company address
 Other (see below)

Cust. no.: _____
Name: _____
Street: _____
Postcode / city: _____
Phone: _____
Fax: _____
E-Mail: _____
Contact: _____

Basic size


Body size:

Body height: _____ cm
 Weight: _____ kg
 Lower leg length: _____ cm
 (anatomic dimension: hollow of the knee to the foot sole)
 Indication: _____



UL This information is necessary for checking. If the information is missing, a check of the relevant items cannot be made!

Wheelchair dimensions:



Seat width SB	<input type="checkbox"/> 35 <input type="checkbox"/> 36 <input type="checkbox"/> 37 <input type="checkbox"/> 38 <input type="checkbox"/> 39 <input type="checkbox"/> 40											
	<input type="checkbox"/> 41 <input type="checkbox"/> 42 <input type="checkbox"/> 43 <input type="checkbox"/> 44 <input type="checkbox"/> 45 <input type="checkbox"/> 46											
	The following dimensions are optionally available subject to surcharge:						<input type="checkbox"/> 30 <input type="checkbox"/> 31 <input type="checkbox"/> 32 <input type="checkbox"/> 33 <input type="checkbox"/> 34 <input type="checkbox"/> 47 <input type="checkbox"/> 48 <input type="checkbox"/> 49 <input type="checkbox"/> 50					
Seat depth ST	Frame size 1		Frame size 2			Frame size 3						
	<input type="checkbox"/> 35 <input type="checkbox"/> 36		<input type="checkbox"/> 39 <input type="checkbox"/> 40 <input type="checkbox"/> 41			<input type="checkbox"/> 44 <input type="checkbox"/> 45 <input type="checkbox"/> 46						
	<input type="checkbox"/> 37 <input type="checkbox"/> 38		<input type="checkbox"/> 42 <input type="checkbox"/> 43			<input type="checkbox"/> 47 <input type="checkbox"/> 48						
						The following dimensions are optionally available subject to surcharge:						
						<input type="checkbox"/> 25 <input type="checkbox"/> 26 <input type="checkbox"/> 27 <input type="checkbox"/> 28 <input type="checkbox"/> 29 <input type="checkbox"/> 30 <input type="checkbox"/> 31 <input type="checkbox"/> 32 <input type="checkbox"/> 33 <input type="checkbox"/> 34						
Back height RH	<input type="checkbox"/> 20 <input type="checkbox"/> 21 <input type="checkbox"/> 22 <input type="checkbox"/> 23 <input type="checkbox"/> 24 <input type="checkbox"/> 25 <input type="checkbox"/> 26 <input type="checkbox"/> 27 <input type="checkbox"/> 28 <input type="checkbox"/> 29											
	<input type="checkbox"/> 30 <input type="checkbox"/> 31 <input type="checkbox"/> 32 <input type="checkbox"/> 33 <input type="checkbox"/> 34 <input type="checkbox"/> 35 <input type="checkbox"/> 36 <input type="checkbox"/> 37 <input type="checkbox"/> 38 <input type="checkbox"/> 39											
	<input type="checkbox"/> 40 <input type="checkbox"/> 41 <input type="checkbox"/> 42 <input type="checkbox"/> 43 <input type="checkbox"/> 44 <input type="checkbox"/> 45 <input type="checkbox"/> 46 <input type="checkbox"/> 47 <input type="checkbox"/> 48											
Lower leg length USL	<input type="checkbox"/> 30 <input type="checkbox"/> 31 <input type="checkbox"/> 32 <input type="checkbox"/> 33 <input type="checkbox"/> 34 <input type="checkbox"/> 35 <input type="checkbox"/> 36 <input type="checkbox"/> 37 <input type="checkbox"/> 38 <input type="checkbox"/> 39 <input type="checkbox"/> 40											
	<input type="checkbox"/> 41 <input type="checkbox"/> 42 <input type="checkbox"/> 43 <input type="checkbox"/> 44 <input type="checkbox"/> 45 <input type="checkbox"/> 46 <input type="checkbox"/> 47 <input type="checkbox"/> 48 <input type="checkbox"/> 49 <input type="checkbox"/> 50											
	The following dimensions are optionally available subject to surcharge:					<input type="checkbox"/> 20 <input type="checkbox"/> 21 <input type="checkbox"/> 22 <input type="checkbox"/> 23 <input type="checkbox"/> 24 <input type="checkbox"/> 25 <input type="checkbox"/> 26 <input type="checkbox"/> 27 <input type="checkbox"/> 28 <input type="checkbox"/> 29						
Seat height back SHh	<input type="checkbox"/> 42 <input type="checkbox"/> 45											
Seat height front SHv	with 4" caster wheel:				<input type="checkbox"/> 44 <input type="checkbox"/> 46 <input type="checkbox"/> 48 <input type="checkbox"/> 50							
	with 5" caster wheel:				<input type="checkbox"/> 46 <input type="checkbox"/> 48 <input type="checkbox"/> 50 <input type="checkbox"/> 52							
	with 6" caster wheel:				<input type="checkbox"/> 48 <input type="checkbox"/> 50 <input type="checkbox"/> 52 <input type="checkbox"/> 54							

Please tick the required dimensions (all dimensions are in cm)

We recommend a seat height difference (difference between SHv and SHh) between 5 cm and 12 cm.

The SHh and SHv values might be a bit different from the specified values,
depending on the wheel size and camber.

Individual dimensions available on request.

Frame colour - One-colour design

Information on current colour designs and example images can be found at: <https://www.proactiv-gmbh.com/colours/>

Frame colour: Selection from colour palette

Frame colour will be specified in the order

A9000100170-000

Raw, blasted, not coated

A9000100170-020

Note: In combination with moisture e.g. hand perspiration, can affect dark aluminum discoloration.

Powder colours

white RAL 9016

gloss A9000100170-023

matt A9000100170-123 € 43.40

brilliant silver

gloss A9000100170-011

anthracite RAL 7016

gloss A9000100170-001

matt A9000100170-101 € 43.40

black RAL 9005

gloss A9000100170-003

matt A9000100170-019 € 43.40

pure orange RAL 2004

gloss A9000100170-022

matt A9000100170-122 € 43.40

pink magenta RAL 4010

gloss A9000100170-037

matt A9000100170-137 € 43.40

pastel blue RAL 5024

gloss A9000100170-036

matt A9000100170-136 € 43.40

lemon green

gloss A9000100170-005

matt A9000100170-105 € 43.40

yellow RAL 1023

gloss A9000100170-046

matt A9000100170-146 € 43.40

red RAL 3001

gloss A9000100170-047

matt A9000100170-147 € 43.40

lagoon blue

gloss A9000100170-052

matt A9000100170-152 € 43.40

athletic grey RAL 7046

gloss A9000100170-050

matt A9000100170-150 € 43.40

muddy grey RAL 7033

gloss A9000100170-051

matt A9000100170-151 € 43.40

peach RAL 3022

gloss A9000100170-038

matt A9000100170-138 € 43.40

Powder colours „dormant“

deep black effect silver

gloss A9000100170-007

sky blue

gloss A9000100170-017

navy blue

gloss A9000100170-045

matt A9000100170-145 € 43.40

electric green

gloss A9000100170-015

matt A9000100170-115 € 43.40

violet

gloss A9000100170-044

fire red

gloss A9000100170-018

matt A9000100170-118 € 43.40

Frame colour - One-colour design special colours „liquid coating“

€ 349.60

storm (ice blue)

gloss A9000100170-500

matt A9000100170-501

prairie (rose gold)

gloss A9000100170-504

matt A9000100170-505

berry

gloss A9000100170-526

matt A9000100170-527

blue ocean

gloss A9000100170-522

matt A9000100170-523

emerald

gloss A9000100170-514

matt A9000100170-515

chili

gloss A9000100170-516

matt A9000100170-517

Frame colour - Two-colour design

Note: In the two-colour design, both colours are the same gloss level as the design colour.

Two-colour design with one base colour (on the rear frame) and one design colour (on the front frame) € 268.50

Please select a colour below:

Base colour (on the rear frame)

white RAL 9016 <input type="checkbox"/> A9000100170-200	brilliant silver <input type="checkbox"/> A9000100170-201	anthracite RAL 7016 <input type="checkbox"/> A9000100170-202
black RAL 9005 <input type="checkbox"/> A9000100170-203	pastel blue RAL 5024 <input type="checkbox"/> A9000100170-206	lemon green <input type="checkbox"/> A9000100170-207
yellow RAL 1023 <input type="checkbox"/> A9000100170-208	lagoon blue <input type="checkbox"/> A9000100170-212	athletic grey RAL 7046 <input type="checkbox"/> A9000100170-210
muddy grey RAL 7033 <input type="checkbox"/> A9000100170-211		

Design colour (on the front frame)

black RAL 9005 <input type="checkbox"/> gloss A9000100170-259 <input type="checkbox"/> matt A9000100170-260	lemon green <input type="checkbox"/> gloss A9000100170-250 <input type="checkbox"/> matt A9000100170-251	chili <input type="checkbox"/> gloss A9000100170-253 <input type="checkbox"/> matt A9000100170-254
pure orange RAL 2004 <input type="checkbox"/> gloss A9000100170-255 <input type="checkbox"/> matt A9000100170-256	emerald <input type="checkbox"/> gloss A9000100170-274 <input type="checkbox"/> matt A9000100170-275	storm (ice blue) <input type="checkbox"/> gloss A9000100170-264 <input type="checkbox"/> matt A9000100170-265
prairie (rose gold) <input type="checkbox"/> gloss A9000100170-266 <input type="checkbox"/> matt A9000100170-267	berry <input type="checkbox"/> gloss A9000100170-284 <input type="checkbox"/> matt A9000100170-285	blue ocean <input type="checkbox"/> gloss A9000100170-280 <input type="checkbox"/> matt A9000100170-281

Colour - Attachment parts

Note: The attachments are always mounted in black.

Drive wheel: aluminum rim, hub and spokes

Black A9000100171-003

Caster fork and bearing block

Black A9000100173-003

Colour design of attachment parts according to customer request:

Drive wheel: aluminum rim and design inlet¹ according to colour selection / hub and spokes black

frame colour A9000100171-999 € 249.60

Selection from one-colour design (different from the frame colour) A9000100171-998 € 288.90

Colour: _____

A9000100170-

Caster wheel: 5" aluminum rim² (6-spoke design)

frame colour A9000100174-999 € 75.20

Selection from one-colour design (different from the frame colour) A9000100174-998 € 87.10

Colour: _____

A9000100170-

Special colour request (Frame colour - One-colour design)

Frame colour for single purchasing Selection from RAL colour palette

A9000100170-100 € 414.40

gloss level

gloss

matt

colour designation

RAL number

Note: Attachments **not** available in design colour options


¹ Only in combination with Liteline hub

² Only in combination with caster wheel 5" aluminum rim 6-spoke-design good shock absorption A900050121*, as well as with suspension fork 5-spoke-design A9000500160*

Frame


Straight front end

A9060100010-K*
Knee angle 105°



Front end V-shape

A9060100020-K*
Knee angle 105°
Note: Footrest area with V-shape
 $x = SB - (0.25 \times USL + 6 \text{ cm})$



Individual production

On request and after checking

Seat width and Seat depth

A9060100100-K € 1,249.50

Seat width 30 cm - 34 cm
and / or

Seat depth 25 cm - 34 cm

If necessary, with forward displacement of the backrest

Seat width 47 cm - 50 cm

A9060100200-K € 1,230.10

Lower leg leght 20 cm - 29 cm

A9060100300-K € 237.60

only for continuous footrest

Back tubes & push handles

Equipment without push handles

A9000401000



Back tube with integrated handles

A9000401051-K € 130.40
straight.



Aluminium push handles fixed in back tube

A9000401050* € 130.40
from RH 23 cm



Horizontal push handles screwed into the back tube

A9000401076 € 130.40
can be unscrewed



Safety push handles with continuous height adjustment

A9000401060* € 216.40
mount in the back tube
from RH 30 cm with adaptable back
from RH 35 cm with lordosis support
from RH 43 cm with ergonomic backshell



Safety push handles back-positioned length 45 cm

A9000401075* € 372.90
continuous height adjustment,
Safety push handle holder mounted on the back cross bar
not in combination with flip-down armrest



Additional option

Back tube with lordosis support

A9000400005-K € 114.40

- 1 cm
- 2 cm
- 3 cm



Backrest

Folding backrest

A9060400011
7-way angle adjustable,
with integrated plastic-coated aluminum side parts



Back equipment

If not specified, the „Adaptable back with back padding“ will be delivered.

Adaptable back with back padding

A9000802010-K001
Cover: Black, breathable, open surface
Equipment with back cross bar



Adaptable back with back padding „light“

A9000802320-K001
Cover: Black, breathable, open surface
Cushion thickness: 9 mm
Equipment with back cross bar

A9000802320-K002
Cover: Black, breathable, open surface
Cushion thickness: 9 mm
Equipment without back cross bar



Adaptable back with back padding „Comfort“

A9000815101-K **€ 131.40**
Cover: Black, breathable, open surface
Front padding: 20 mm
Rear padding: 5 mm
Equipment with back cross bar

Back rise + 3 - 5 cm

A9000802100 **€ 159.10**
for the nominal back height dimension
Only in combination with height-adjustable
push handles



Equipment without back upholstery and back padding or back system


A9000802011 **€ - 53.00**
Equipment with back cross bar

Seat upholstery

If not specified, the „Body Contour“ seat upholstery will be delivered.


Body Contour

A9000801025-K



Open belt system

A9000801010-K **€ 79.90**





Equipment without seat upholstery

A9000801050



Seat cushion



Seat cushion in rectangular contour Inlet: Standard

A9000805180-K* **€ 118.30**

Seat cushion with radius contour Inlet: Standard

A9000800300-K* **€ 163.10**

Seat cushion height:

3 cm 5 cm
 4 cm 6 cm

Seat cushion in rectangular contour Inlet: Anti-Decubitus

up to 30x30 cm
A9000805204-K* **€ 301.70**

>30x30 up to 50x50 cm
A9000805200-K* **€ 437.70**

Seat cushion with radius contour Inlet: Anti-Decubitus

up to 30x30 cm
A9000800311-K* **€ 356.10**

>30x30 up to 50x50 cm
A9000800310-K* **€ 482.90**

Seat cushion height:

4 cm 6 cm

Cover:

Black, breathable, open surface
 Black, breathable, closed surface

If not specified, the cover material „Black, breathable, open surface“ will be delivered.

Seat cushion is provided

Lift mechanism

Mechanical seat height adjustment by 30 cm

A9060300100-K*
via gas springs up to max. 2x1200 N
up to approx. 100 kg

Mechanical seat height adjustment by 30 cm

A9060300200-K* € 985.00
via gas springs up to max. 4x900 N
from approx. 100 kg
incl. lever extension

Armrest

Armrest with aluminum handle bar

A9060400100-K*

Armrest with arm padding 25 cm long

A9060400200-K* € 165.80

Armrest with arm padding 38 cm long

A9060400300-K* € 197.00

Release

Lifting function release mechanism

attached to the armrest
A9060403040
 right
 left

Note: The release mechanism is only possible on the right **OR** on the left.
If not specified, the „release on the right“ will be delivered.

Lever extension for release mechanism with pivoting handle bar

A9060403440-K € 175.20

Release mechanism, additionally equipped with a pressure plate

A9060403450 € 319.20

for force support of the finger function
not in combination with arm padding

Axle plate

Standard equipment with 12 mm sockets / axles

Drive wheel plate


A9060201100
Rear seat height approx. 45 cm
(for 24" wheel)

Axle plate lowered

A9060201120-K € 260.90
Rear seat height approx. 42 cm
(for 24" wheel)

Extended drive wheel plate

A9060201110-K € 292.00
Wheelbase extension by 7.5 cm



Wheel camber

*Total width with tight-mounted handrim.

1° drive wheel camber

A9060202000
Total width = approx. SB + 20 cm*

E-Adaptations

Preparation for the „e-motion" mounting

A9060201200-K € 502.90

Equipment: e-motion torque arm,
1/2" axle and socket length 50.8 mm (2")
Not allowed in combination with „e-motion DuoDrive“.

Preparation for the „e-fix" mounting

E35/E36
(current models) € 502.90
A9060201230-K

E25/E26
(predecessor models for reuse) € 502.90
A9060201220-K

Equipment: 1/2" axle,
socket length 50.8 mm (2")
Torque support e-fix and battery holder.
The e-fix attachment must be carried out by Fa.
Alber, as an additional safety device is required
in accordance with Alber installation guidelines.

Brake system

If not specified, the „Knee lever brake system mounted at the rear” will be delivered.

Knee lever brake system

A9060701010



Knee lever brake system mounted low down

A9060701015 € 156.10

Brake holding rail mounted on the axle plate, straight brake lever version



Dual knee lever brake system

A9060701020 € 373.20

One brake mounted at the rear and another brake mounted at the front on the wheel plate with long brake lever in straight design

Footrest

If not specified, the „continuous footrest” will be delivered.

Continuous footrest

A9016601010-K*
with tubular resting bracket
height adjustable +/- 3 cm
If you want to equip the footrest with a „rubber coated footplate” please select this under „Footrest accessories”.



Footrest folds up to the rear with spring locking mechanism

A9016602005-K* € 287.40

with tubular resting bracket from USL 40 cm only with V-shape and frame size 3
If you want to equip the footrest with a „rubber coated footplate” please select this under „Footrest accessories”.



Footrest folds up to the rear

A9016602010-K* € 307.90

with tubular resting bracket and 20 mm operating lever from frame size 2 from USL 42 cm V-shape from SB 38 cm
If you want to equip the footrest with a „rubber coated footplate” please select this under „Footrest accessories”.



Divided footrest

A9016603200-K € 336.80

with rubber coated footplate, height adjustable +/- 3 cm, angle adjustment, can be pivoted outwards and folded up to the side
Including frame reinforcement V-shape from SB 40 cm



„Swing Away” footrest

A9016604100-K € 404.70

with rubber coated footplate, separated in the middle, can be folded to the side, removable, can be folded up, height-adjustable, angle-adjustable
Including frame reinforcement from USL 34 cm



Removable front frame, foldable with divided footrest

A9065630000-K € 845.70

height-adjustable, with rubber plate
Only with straight front end
From 100 kg user weight only on request and after checking



Footrest accessories

Rubber coated footplate

A9000610000-K € 125.90

bolted, angle-adjustable
Depth: 14 cm



Footplate side stop

A9000610100 € 59.70

adjustable
only in combination with rubber coated footplate



Adaptations

Anti-tipping support, removable

A9060300100
adjustable with Quickpin



3" caster wheel

3" skate wheel (80 mm)

A9000501110 € 70.80
puncture-proof



4" caster wheel

4" caster wheel „Standard“ solid rubber black

A9000501011-004
(82A), puncture-proof



4" skate wheel (100 mm)

A9000501120 € 70.80
puncture-proof



4" safety light-up wheel (100 mm)

A9000501220 € 94.80
puncture-proof



4" CorMORun caster wheel (100 mm)

A9000501042 € 173.00
with replaceable rubber 6 mm,
puncture-proof
Rim colour: Black
Incl. hair guard spacer (5.3 mm)



4" Frog Legs caster wheel (100 mm)

A9000501230 € 206.30
with replaceable rubber 6 mm,
puncture-proof
Colour: black, silver rim
Incl. hair guard spacer (5.3 mm)



5" caster wheel

5" caster wheel „Standard“ PU black

A9000501036-005
(82A), puncture-proof



5" caster wheel „Standard“ PU grey

A9000501031-005
(82A), puncture-proof



5" safety light-up wheel (125 mm)

A9000501221 € 94.80
puncture-proof



<p>5" aluminum rim 3-spoke design „good shock absorption“</p> <p><input type="checkbox"/> silver anodised A9000501210 € 179.80</p> <p><input type="checkbox"/> black A9000501211 € 179.80</p> <p>with replaceable rubber 6 mm medium hard, puncture-proof Incl. hair guard spacer (3 mm)</p> 	<p>5" CorMORun caster wheel (125 mm)</p> <p><input type="checkbox"/> A9000501041 € 181.30</p> <p>with replaceable rubber 6 mm, puncture-proof Rim colour: Black Incl. hair guard spacer (5.3 mm)</p> 	<p>5" Frog Legs caster wheel (125 mm)</p> <p><input type="checkbox"/> A9000501232 € 206.30</p> <p>with replaceable rubber 6 mm, puncture-proof Colour: black, silver rim Incl. hair guard spacer (5.3 mm)</p> 
--	---	--

6" castor wheels available on request after checking, only with V-shape front end.

Caster fork

If not specified, the „Standard caster fork with screw axle“ will be delivered.

Caster fork Standard

- A9000500100-K001 with screwed axle
- A9000500100-K002 € 167.40 with quick-release axle



Caster fork Heavy Duty

- for high loads or for use with additional electric drives
- A9000510111-K001 € 54.00 with screwed axle
- A9000510111-K002 € 218.90 with quick-release axle



Caster fork spare wheel „Standard“

- A9000500300-K € 296.70
- Caster fork Standard with quick-release axles



Caster fork spare wheel „Heavy Duty“

- A9000500301-K € 337.90
- Caster fork Heavy Duty with quick-release axles



The spare caster fork does not have a caster wheel.
In case a caster wheel is requested, please put a cross on a second selection in the category „Caster wheels“.

Track distance

Caster wheel track extension

- A9060500400-K € 73.60
- 13 mm per side



Drive wheels

If not specified, the „drive wheel aluminum rim 36 spokes“ will be delivered.

Drive wheel 36 spoke aluminum rim

24" (ETRTO 540)
A9000203000-010
with Liteline hub (left image),
radially spoked
up to 120 kg body weight

24" (ETRTO 540)
A9000203000-020
with aluminum front hub (right image),
radially spoked
over 120 kg body weight



Drive wheel double rim 36-spokes

24" (ETRTO 540)
A9000203100 € 157.80

25" (ETRTO 559)
A9000203200 € 157.80

26" (ETRTO 590)
A9000203300 € 157.80
with aluminum front hub,
radially spoked
For use with high loads and body weight
>120 kg



Spinergy wheel „Spox Sport“ with 24 PBO spokes; rim, hub and spokes in black

24" (ETRTO 540)
A9000203550 € 784.10


25" (ETRTO 559)
A9000201559 € 784.10

26" (ETRTO 590)
A9000203555 € 784.10



Handrims

x = dimension handrim attachment



Handrim attachment
general:

- x = 11 mm
 x = 19 mm

Handrim attachment
Curve + Gekko handrim:

- x = 10 mm
 x = 18 mm

If not specified, the handrim is mounted wide (19 mm / 18 mm).

Stainless steel Ø 20 mm

24" handrim for
24" wheel and 25" wheel
A9000205010-004

26" handrim for
26" wheel
A9000205010-006



Aluminum Ø 19 mm silver anodised

24" handrim for
24" wheel and 25" wheel
A9000205019-004



Aluminum Ø 19 mm powder coated in frame colour

(basic colour with two-colour design)

- = Gloss surface
 = Matt surface

Note: Matt coating is not recommended due to the low grip.

24" handrim for
24" wheel and 25" wheel € 68.30
A9000205020-004

26" handrim for
26" wheel € 68.30
A9000205020-006



Aluminum Ø 19 mm polished and hard-anodised

24" handrim for
24" wheel and 25" wheel € 164.90
A9000205030-004

26" handrim for
26" wheel € 164.90
A9000205030-006



Titanium Ø 19 mm

24" handrim for
24" wheel and 25" wheel € 335.60
A9000205040-004

26" handrim for
26" wheel € 335.60
A9000205040-006



Silicone handrim cover

for handrim 24"
A9000205100-004 € 232.80

for handrim 26"
A9000205100-006 € 232.80

Surface without structure,
compatible with handrim Ø 19 mm or
Ø 20 mm.



<p>Ultra-Grip „very grippy“</p> <p><input type="checkbox"/> 24" handrim for 24" wheel and 25" wheel € 402.40 A9000205070-004</p> <p><input type="checkbox"/> 26" handrim for 26" wheel € 402.40 A9000205070-006</p> 	<p>Ergo-Grip</p> <p><input type="checkbox"/> 24" handrim for 24" wheel and 25" wheel € 402.40 A9000205080-004</p> <p><input type="checkbox"/> 26" handrim for 26" wheel € 402.40 A9000205080-006</p> <p>Attachment fitting to the rim</p> <p>Note: Due to the complex geometry, concentricity deviations of +/- 2-3 mm are possible.</p> 	<p>Curve black anodised, ergonomic handrim profile</p> <p><input type="checkbox"/> 24" handrim for 24" wheel and 25" wheel € 402.40 Profile height 24 mm A9000205130-004</p> <p><input type="checkbox"/> 24" handrim for 24" wheel and 25" wheel € 402.40 Profile height 31 mm Size L A9000205140-004</p> <p><input type="checkbox"/> 26" handrim for 26" wheel € 402.40 Profile height 31 mm Size L A9000205140-006</p> 
<p>Curve silver anodised ergonomic handrim profile</p> <p><input type="checkbox"/> 24" handrim for 24" wheel and 25" wheel € 402.40 Profile height 31 mm Size L A9000205160-004</p> 	<p>Curve powder-coated, ergonomic handrim profile</p> <p><input type="checkbox"/> = frame colour (gloss surface) <input type="checkbox"/> = black (gloss surface)</p> <p><input type="checkbox"/> 24" handrim for 24" wheel and 25" wheel € 467.60 Profile height 24 mm A9000205135-004</p> <p><input type="checkbox"/> 24" handrim for 24" wheel and 25" wheel € 467.60 Profile height 31 mm Size L A9000205145-004</p> <p><input type="checkbox"/> 26" handrim for 26" wheel € 467.60 Profile height 31 mm Size L A9000205145-006</p> 	<p>Curve Grip, with Tetra-Grip coating, ergonomic handrim profile</p> <p><input type="checkbox"/> 24" handrim for 24" wheel and 25" wheel € 520.70 Profile height 31 mm Size L A9000205150-004</p> 
<p>Quadro with Tetra-Grip coating</p> <p><input type="checkbox"/> 24" handrim for 24" wheel and 25" wheel € 520.70 A9000205170-004</p> <p><input type="checkbox"/> 25" handrim for 25" wheel € 520.70 A9000205170-005</p> <p><input type="checkbox"/> 26" handrim for 26" wheel € 520.70 A9000205170-006</p> 	<p>Gekko black anodised</p> <p><input type="checkbox"/> 24" handrim for 24" wheel and 25" wheel € 520.70 A9000205175-004</p> <p><input type="checkbox"/> 25" handrim for 25" wheel € 520.70 A9000205175-005</p> <p><input type="checkbox"/> 26" handrim for 26" wheel € 520.70 A9000205175-006</p> <p>with a robust and replaceable silicone lip, fully integrated into the profile Profile height 31 mm</p> 	<p>Gekko silver anodised</p> <p><input type="checkbox"/> 24" handrim for 24" wheel and 25" wheel € 520.70 A9000205185-004</p> <p>with a robust and replaceable silicone lip, fully integrated into the profile Profile height 31 mm</p> 

Gekko S black anodised

- 24" handrim for 24" wheel and 25" wheel **€ 520.70**
A9000205195-004

with a robust and replaceable silicone lip, fully integrated into the profile
Profile height 24 mm



Equipment without handrims

- A9000205990

Spoke guard

Finger- and spoke guard „transparent“

- for 24" wheel **€ 207.10**
A9000206300-004

- for 26" wheel **€ 207.10**
A9000206300-006

Not in combination with Spinergy drive wheels.



Tyres

If not specified, the „Schwalbe RightRun in black“ tyres will be delivered.

Schwalbe „RightRun“ High-pressure tyres

- 24x1" (540-25)
A9000204001-004

Colour: Black
7.5 bar



Schwalbe „RightRun“ High-pressure tyres

- 24x1" (540-25)
A9000204000-004

- 25x1" (559-25)
A9000204000-005

- 26x1" (590-25)
A9000204000-006

Colour: Grey, black wall
7.5 bar



Schwalbe „DownTown“ Profiled tyres

- 24x1" (540-25)
A9000204010-004

7.5 bar



Schwalbe „Marathon Plus Evolution“ with „smart guard“ puncture protection

- 24x1" (540-25) **€ 49.20**
A9000204030-014

- 25x1" (559-25) **€ 49.20**
A9000204030-015

- 26x1" (590-25) **€ 49.20**
A9000204030-016

Colour: Black
10 bar



Schwalbe „Marathon Plus“ with „smart guard“ puncture protection

- 24x1" (540-25) **€ 49.20**
A9000204030-004

Colour: Grey, black wall
10 bar



Schwalbe „One“


- 24x1" (540-25) **€ 93.00**
A9000204123-004

Colour: Black
9 bar



<p>Puncture-proof „solid rubber tyres”</p> <p><input type="checkbox"/> 24x1" (540-25) € 122.70 A9000204201-004 Colour: Black Not in combination with Spinergy drive wheels.</p> 	
--	--

Accessories

<p>Through axle with assistance, Large push button</p> <p><input type="checkbox"/> A9000200096 € 105.80 Interlock dimension: 122 / 110 mm Axle diameter: 12 mm</p> 	<p>Hand air compressor with pressure gauge</p> <p><input type="checkbox"/> A9000901100 € 106.30 up to 14 bar</p> 	<p>Crutch holder with flexible nylon cup</p> <p><input type="checkbox"/> A9000902310* € 137.40</p> 
<p>Length-adjustable safety belt with buckle</p> <p><input type="checkbox"/> A9000902100 € 80.60</p> 	<p>Adjustable safety belt with hook and fastener</p> <p><input type="checkbox"/> A9000902110 € 64.50</p> 	<p>SIGMA lighting set</p> <p><input type="checkbox"/> A9000903270 € 107.90 Consisting of front light and rear lights, chargeable via Micro USB socket, incl. Micro USB charging cable Lighting set compliant with road traffic regulations</p> 
<p>Transport bag for drive wheels</p> <p><input type="checkbox"/> for 24" wheel A9000903530-004 € 77.80 Ø 76 cm x 2 cm</p> <p><input type="checkbox"/> for 26" wheel A9000903530-006 € 77.80 Ø 83 cm x 2 cm</p> 	<p>Stand-up aid handles via locking pin</p> <p><input type="checkbox"/> A9060403200 € 470.30 removable, height-adjustable mounted on the seat tube in combination with the „knee lever brake system mounted low down” the height adjustability is limited.</p> 	

Separate order sheets

<p>Category „Rehabilitation technology accessories” e.g. Seat & backrest systems can be found at: www.proactiv-gmbh.com/downloads/</p>	
---	--

Services

<p>Consultation, clarification of the installation situation and measurement through PRO ACTIV</p> <p><input type="checkbox"/> A9990000000 € 354.50</p>	<p>Handover and instruction on the wheelchair by PRO ACTIV</p> <p><input type="checkbox"/> A9990000006 € 354.50</p>	
---	---	--

Notes and remarks:

Location, Date

Signature End customer/dealer

The signature field is intended for the consultant and insured person. This confirms that the intended equipment is intended for the insured person.
- 15 - 2026-03-03 / The prices are plus statutory value-added tax and shipping costs = included in the basic price = Position selected
Order sheet LIFT solid - mechanical Rev. 08 UK