

Customized production - model „SPEEDY A2“

Frame design

Frame design Heavy Duty A1-4 up to 150 kg body weight

A9010100091 € 478.10

resp. from SB 51 cm to SB 58 cm:

- Seat tube with additional reinforcement inner tube
- Additional lower frame cross stiffener tube
- Additional frame diagonal stiffener tubes to the lower rear cross tube
- Back tube reinforcement insert
- Stainless steel backrest locking pin
- Side part thickness 3 mm aluminium
- Reinforced drive wheel plate 2 slots, 12 mm thick

Note:

Please select the „Heavy Duty“ steering wheel fork under Steering wheel fork!

Frame design Heavy Duty A1-4 from 151 kg body weight

A9010100092 € 943.40

according to individual specifications,
SB ____ cm

- Seat tube with additional reinforcement inner tube
- Additional frame cross stiffening tubes
- Additional frame diagonal stiffening tubes to the lower rear cross tube and seat tube
- Reinforced drive wheel plate, 2 slots, 12 mm thick
- Back tube reinforcement insert
- Stainless steel backrest locking bolts
- Stainless steel push handles
- Side part thickness 3 mm aluminium
- Footrest or foot bar, continuous, reinforced design
- Knee lever brake bolt, stainless steel

Note:

Please select the „Heavy Duty“ steering wheel fork under Steering wheel fork!

Frame

Customized frame geometry abducted shape

A9010100305 € 1,281.00

Frame extension per side abducted by
specification

2,5 cm 3,5 cm

Customized frame length (must be checked individually)

A9010100070 € 416.00

Frame length to frame arch ____ cm
(for ST ____ cm)

Customized frame height (must be checked individually)

A9010100110 € 382.60

Frame height: ____ cm

Knee angle

Individual knee angle

A9000100050-003 € 534.50

90° 100° 110° 115°

Services

Consulting and customization by PRO ACTIV according to individual customer requirements

A9990000017 € 612.00

Additional requirements