



FREAK compact bike

approved up to 100 kg payload

- Underframe chassis with reinforcement strut
- chassis length 102 cm, SB 35 cm
- Drive unit can be uncoupled via excenter lock
- black drive unit with 20" drive wheel
- A-shaped handles with integrated gearshift controls
- Shimano V rim brake
- Backrest with continuous angle adjustment
- wheel camber 6°, rear wheel size 20"
- safety flag
- Seat upholstery Body Contour

A9113000013

€ 4,640.00

Drive wheel on Freak with 26"

☐ A9113000014 € 253.60

V-brake replaced by mechanical disc brake.
Not compatible with Nexus 8-speed hub.

☐ **Quotation date:** _____

Consultation was carried out by
Staff at the medical supplies shop: _____

Enquiry no.: _____

☐ **Order date:** _____

according to GTC PRO ACTIV Reha-Technik GmbH

Order no.: _____

Medical suppliers / company address:

Cust. no.: _____

Name: _____

Street: _____

Postcode / city: _____

Phone: _____

Fax: _____

E-Mail: _____

Contact: _____

☐ **Different invoicing address (see below)**

Cust. no.: _____

Name: _____

Street: _____

Postcode / city: _____

☐ **Different shipping address (see below)**

Cust. no.: _____

Name: _____

Street: _____

Postcode / city: _____

Customer name / commission:

ID / Comm. No.: _____

Surname, forename: _____

Street: _____

Postcode / city: _____

Phone: _____

E-Mail: _____

D.o.B.: _____

Paying organisation: _____

Indication: _____

Body size

Body height: _____ cm

Weight: _____ kg

Frame colour - One-colour design

Information on current colour designs and example images can be found at: <https://www.proactiv-gmbh.com/colours/>

white RAL 9016

☐ gloss A9000100170-023

☐ matt A9000100170-123 € 43.40

brilliant silver

☐ gloss A9000100170-011

anthracite RAL 7016

☐ gloss A9000100170-001

☐ matt A9000100170-101 € 43.40

black RAL 9005

☐ gloss A9000100170-003

☐ matt A9000100170-019 € 43.40

deep black effect silver

☐ gloss A9000100170-007

pure orange RAL 2004

☐ gloss A9000100170-022

fire red

☐ gloss A9000100170-018

☐ matt A9000100170-118 € 43.40

pink magenta RAL 4010

☐ gloss A9000100170-037

☐ matt A9000100170-137 € 43.40

sky blue

☐ gloss A9000100170-017

pastel blue RAL 5024

☐ gloss A9000100170-036

☐ matt A9000100170-136 € 43.40

navy blue

☐ gloss A9000100170-045

☐ matt A9000100170-145 € 43.40

electric green

☐ gloss A9000100170-015

☐ matt A9000100170-115 € 43.40

lemon green

☐ gloss A9000100170-005

☐ matt A9000100170-105 € 43.40

violet

☐ gloss A9000100170-044

yellow RAL 1023

☐ gloss A9000100170-046

☐ matt A9000100170-146 € 43.40

red RAL 3001

☐ gloss A9000100170-047

☐ matt A9000100170-147 € 43.40

lagoon blue

☐ gloss A9000100170-052

☐ matt A9000100170-152 € 43.40

athletic grey RAL 7046

☐ gloss A9000100170-050

☐ matt A9000100170-150 € 43.40

muddy grey RAL 7033

☐ gloss A9000100170-051

☐ matt A9000100170-151 € 43.40

peach RAL 3022

☐ gloss A9000100170-038

☐ matt A9000100170-138 € 43.40

Hub gear

Shimano „Nexus“ 8-speed hub gear

☐ A91132000015

„Grip Shift“ twist shift
Gear ratio: 32 teeth sprocket,
cassette 20 teeth
with back-pedal brake
Shimano V-Brake brake handle left

Shimano „Alfine“ 11-speed hub gear

☐ A91132000026

€ 833.70

„Rapid Fire“ gear lever
Gear ratio: 32 teeth sprocket,
cassette 20 teeth
with hydraulic disc brake (Magura MT2) brake
handle right and Shimano V-Brake brake
handle left



Derailleur gears

Shimano Deore 1x11 Rapid Fire 11-speed derailleur

☐ A91132000034

€ 517.50

mechanical derailleur gear with control on the
handle
11-speed cassette 11-42 teeth
with hydraulic disc brake (Magura MT4) brake
lever on the right
and Shimano V-brake brake lever on the left
Chainring at the top with chain guard discs

Shimano Deore XT 1x11 DI2 with push button operation 11-speed derailleur

☐ A91132000035

€ 1,800.20

electronic derailleur gear with operation on the
handle.
11-speed cassette 11-42 teeth
with hydraulic disc brake (Magura MT4) brake
lever on the right
and Shimano V-brake brake lever on the left
Chainring at the top with chain guard discs



Selection of the chainring:

☐ 32 teeth

☐ 38 teeth

Brake system

Powerful braking system with brake disc Ø 203 mm

☐ A9110300105 € 83.10

recommended from 80 kg body weight.



Dimensions

Crank type and grip width

Bottom bracket shaft with standard crank

☐ grip width 40 cm
A9110500100-040

☐ grip width 45 cm
A9110500100-045

Bottom bracket shaft with tetra crank (flared)

with extended crank arms for maximum upper body space

☐ grip width 48 cm (176 mm)
A9110500100-048

Standard pedal crank

☐ length 155 mm
A9110600100-155

☐ length 175 mm
A9110600100-175



Tetra pedal crank

☐ length 155 mm € 136.10
A9100600105-155


☐ length 175 mm € 136.10
A9100600105-175



Fixed bottom bracket support

(Limited availability – only as long as stocks last. Alternative: Adjustable bottom bracket support)

☐ S (X/Y 67/153 mm)
A9100501100-003



Adjustable bottom bracket support

the angle and length can be adjusted via the clamped rotary point on the fork bridge and length.

☐ S (180 mm) € 232.10
A9100500001-002

☐ M (240 mm) € 232.10
A9100500001-003

If not specified, the adjustable bottom bracket support will be delivered in M.



Legrest

Divided legrest

☐ A9100100008
adaptation on the fork
for short legs
Only possible up to body size 165 cm



Foot stop for legrest

☐ A9100100012 € 91.20



Back equipment

Backrest with filling plate and back cushion

(Limited availability – only as long as stocks last.)

☐ A9113800002
Cover: Black, breathable, open surface

Backrest upholstery Body Contour


☐ A9100800010 € 92.40
with back cushion 3 cm
Cover: Black, breathable, open surface



Seat cushion

Seat cushion „standard“

☐ A9113000190 € 115.00
for chassis 102 cm
SB: 35 cm
height: 4 cm
Cover: Black, breathable, open surface



Tyres

20x1.1" Schwalbe „Durano“

☐ A9113105001
(ETRTO 406-28)
road bike tires



20x1.75" Schwalbe „Marathon Plus“

☐ A9113105120 € 71.10
(ETRTO 406-47)
puncture-proof
with Smart Guard puncture protection inserts



Accessories

Plug-in mudguard

☐ A9100502100 € 58.10



Rear guard clamped on the axle tube

☐ A9100800100 € 218.10
with 40 mm mounting clamp



Neck rest

☐ A9100800150 € 261.00
height-adjustable, angle-adjustable



Pelvic fixation belt

☐ A9100800200 € 121.70



SIGMA BC 8.0 ATS bike computer

☐ A9110901000 € 70.10
wireless, with mount



SIGMA BC 14.0 WL CAD bike computer

☐ A9110900018 € 171.50
with pedal frequency set, wireless



Bracket Ø 22 mm x 42 mm for mounting accessories (tacho, lighting, etc.) ☐ A9110903271 pcs. á € 23.20 Colour: Natural anodised _____ piece 	SIGMA lighting set ☐ A9120900019 € 123.90 consisting of front light and rear lights, chargeable via Micro USB socket, incl. Micro USB charging cable. Lighting set compliant with road traffic regulations 	
---	---	--

Options for users with reduced functions

Palm rest on the handle ☐ A9110601200 € 102.80 upholstered, bolted, angle-adjustable. 	Tetra hand support, pivotable, adjustable grip distance ☐ A9110601101 € 169.40 	Tetra backrest swivel and adjustable grip width ☐ A9110601102 € 191.60 very stiff external support as round tube Ø 20 mm with foam rubber coating Not in combination with Tetra hand support pivotable 
Tetra shell handles ☐ A9110601100 € 172.40 		

Services

Consultation, clarification of the installation situation and measurement through PRO ACTIV ☐ A9990000000 € 354.50	Handover and instruction of the handbike / traction device by PRO ACTIV ☐ A9990000011 € 354.50	
--	--	--

Notes and remarks:

Location, Date

Signature End customer/dealer

The signature field is intended for the consultant and insured person. This confirms that the intended equipment is intended for the insured person.

- 5 - 2026-07-15 / The prices are plus statutory value-added tax and shipping costs

☐ = included in the basic price ☒ = Position selected

Order sheet FREAK compact bike Rev. 08 UK