



## One-hand drive with double handrim

Approved up to 120 kg load capacity

**Quotation date:** \_\_\_\_\_  
**Consultation was carried out by**  
**Staff at the medical supplies shop:** \_\_\_\_\_  
**Enquiry no.:** \_\_\_\_\_

**Order date:** \_\_\_\_\_  
according to GTC PRO ACTIV Reha-Technik GmbH  
**Order no.:** \_\_\_\_\_

### Medical suppliers / company address:

Cust. no.: \_\_\_\_\_  
Name: \_\_\_\_\_  
Street: \_\_\_\_\_  
Postcode / city: \_\_\_\_\_  
Phone: \_\_\_\_\_  
Fax: \_\_\_\_\_  
E-Mail: \_\_\_\_\_  
Contact: \_\_\_\_\_

**Shipping address:**  same as company address  
 Other (see below)

Name: \_\_\_\_\_  
Street: \_\_\_\_\_  
Postcode / city: \_\_\_\_\_

### Customer name / commission:


ID / Comm. No.: \_\_\_\_\_  
Surname, forename: \_\_\_\_\_  
Street: \_\_\_\_\_  
Postcode / city: \_\_\_\_\_  
Phone: \_\_\_\_\_  
Fax: \_\_\_\_\_  
E-Mail: \_\_\_\_\_  
D.o.B.: \_\_\_\_\_  
Paying organisation: \_\_\_\_\_  
ID / Comm. No.: \_\_\_\_\_

**Invoicing address:**  same as company address  
 Other (see below)

Cust. no.: \_\_\_\_\_  
Name: \_\_\_\_\_  
Street: \_\_\_\_\_  
Postcode / city: \_\_\_\_\_  
Phone: \_\_\_\_\_  
Fax: \_\_\_\_\_  
E-Mail: \_\_\_\_\_  
Contact: \_\_\_\_\_

## Version, attachment for new wheelchair orders

Please note: For a wheelchair with one-hand drive, the attachments for the wheelchair are only available in **black**.



wheel camber: 1°, 4°, 6°, 8°, 10°, 12° for TRAVELER, SPEEDY A1 / A2  
 wheel camber: 2°, 4°, 6°, 8°, 10° for BUDDY  
 wheel camber: 2° for SPEEDY 4all / 4you / 4all Ergo / 4you Ergo / 4teen

Please select the desired operation side:

Operation:  right side  left side

### TRAVELER

incl. removable drive shaft

A9030208010\* € 1,504.90

### SPEEDY 4all / 4you / 4all Ergo / 4you Ergo / 4teen

A9016208010\* € 1,325.80

### SPEEDY A1 / A2

A9010208010\* € 1,325.80

### BUDDY

A9000220010\* € 1,325.80

## Handrim A

If not specified, the handrim is mounted wide (19 mm).

Attachment:  tight x = 11 mm  wide x = 19 mm

### Handrim stainless steel

for 20" wheel, Ø 452 mm

A9000205010-101

for 22" wheel, Ø 488 mm

A9000205010-102

for 24" and 25" wheel Ø 539 mm

A9000205010-104

for 26" wheel, Ø 590 mm

A9000205010-106

### Handrim aluminum

**powder coated in frame colour**

(basic colour with two-colour design)

= Gloss surface

= Matt surface

Note: Matt coating is not recommended due to the low grip.

for 20" wheel, Ø 452 mm

A9000205020-101 € 34.20

for 22" wheel, Ø 488 mm

A9000205020-102 € 34.20

for 24" and 25" wheel, Ø 538 mm

A9000205020-104 € 34.20

for 26" wheel, Ø 590 mm

A9000205020-106 € 34.20

### Handrim aluminum

**polished and hard-anodised**

for 20" wheel, Ø 452 mm

A9000205030-101 € 82.40

for 22" wheel, Ø 488 mm

A9000205030-102 € 82.40

for 24" and 25" wheel, Ø 538 mm

A9000205030-104 € 82.40

for 26" wheel, Ø 590 mm

A9000205030-106 € 82.40

### Handrim Ultra-Grip

for 20" wheel, Ø 452 mm

A9000205070-101 € 201.30

for 22" wheel, Ø 486 mm

A9000205070-102 € 201.30

for 24" and 25" wheel, Ø 548 mm

A9000205070-104 € 201.30

for 26" wheel, Ø 587 mm

A9000205070-106 € 201.30

## Handrim B with 3-spoke mounting

Please select a handrim B that is at least 2" smaller than handrim A.

If not specified, the handrim is mounted wide (38 mm).

Attachment:  tight y = 28 mm  wide y = 38 mm

### Handrim stainless steel

- 18", Ø 398 mm  
 A9000208010-100
- 20", Ø 452 mm  
 A9000208010-101
- 22", Ø 488 mm  
 A9000208010-102
- 24", Ø 539 mm  
 A9000208010-104

### Handrim aluminum powder coated in frame colour (basic colour with two-colour design)

- = Gloss surface
- = Matt surface
- Note: Matt coating is not recommended due to the low grip.
- 18", Ø 397 mm  
 A9000208020-100 € 34.20
- 20", Ø 452 mm  
 A9000208020-101 € 34.20
- 22", Ø 488 mm  
 A9000208020-102 € 34.20
- 24", Ø 538 mm  
 A9000208020-104 € 34.20

### Handrim aluminum polished and hard-anodised

- 18", Ø 397 mm  
 A9000208030-100 € 82.40
- 20", Ø 452 mm  
 A9000208030-101 € 82.40
- 22", Ø 488 mm  
 A9000208030-102 € 82.40
- 24", Ø 538 mm  
 A9000208030-104 € 82.40

### Handrim Ultra-Grip

- 18", Ø 397 mm  
 A9000208070-100 € 201.30
- 20", Ø 452 mm  
 A9000208070-101 € 201.30
- 22", Ø 486 mm  
 A9000208070-102 € 201.30
- 24", Ø 548 mm  
 A9000208070-104 € 201.30

## Handrim C

If not specified, the handrim is mounted wide (19 mm).

Attachment:  tight x = 11 mm  wide x = 19 mm

### Handrim stainless steel

- for 20" wheel, Ø 452 mm  
 A9000205010-101
- for 22" wheel, Ø 488 mm  
 A9000205010-102
- for 24" and 25" wheel, Ø 539 mm  
 A9000205010-104
- for 26" wheel, Ø 590 mm  
 A9000205010-106

### Handrim aluminum powder coated in frame colour (basic colour with two-colour design)

- = Gloss surface
- = Matt surface
- Note: Matt coating is not recommended due to the low grip.
- for 20" wheel, Ø 452 mm  
 A9000205020-101 € 34.20
- for 22" wheel, Ø 488 mm  
 A9000205020-102 € 34.20
- for 24" and 25" wheel, Ø 538 mm  
 A9000205020-104 € 34.20
- for 26" wheel, Ø 590 mm  
 A9000205020-106 € 34.20

### Handrim aluminum polished and hard-anodised

- for 20" wheel, Ø 452 mm  
 A9000205030-101 € 82.40
- for 22" wheel, Ø 488 mm  
 A9000205030-102 € 82.40
- for 24" and 25" wheel, Ø 538 mm  
 A9000205030-104 € 82.40
- for 26" wheel, Ø 590 mm  
 A9000205030-106 € 82.40

<p><b>Handrim Ultra-Grip</b></p> <p>for 20" wheel, Ø 452 mm  <input type="checkbox"/> A9000205070-101      € 201.30</p> <p>for 22" wheel, Ø 486 mm  <input type="checkbox"/> A9000205070-102      € 201.30</p> <p>for 24" and 25" wheel, Ø 548 mm  <input type="checkbox"/> A9000205070-104      € 201.30</p> <p>for 26" wheel, Ø 587 mm  <input type="checkbox"/> A9000205070-106      € 201.30</p>	<p><b>Equipment without handrim</b></p> <p><input checked="" type="checkbox"/> A9000205990</p>	
<b>Other</b>		
<p><b>Silicone handrim cover Surface without structure</b></p> <p>for handrim 20", Ø 448 mm  <input type="checkbox"/> A9000205100-001      € 116.30/pcs.</p> <p>for handrim 22", Ø 488 mm  <input type="checkbox"/> A9000205100-002      € 116.30/pcs.</p> <p>for handrim 24", Ø 538 mm  <input type="checkbox"/> A9000205100-004      € 116.30/pcs.</p> <p>for handrim 26", Ø 585 mm  <input type="checkbox"/> A9000205100-006      € 116.30/pcs.</p> <p>Please select which handrim you would like to attach the handrim cover to:</p> <p><input type="checkbox"/> handrim A</p> <p><input type="checkbox"/> handrim B</p> <p><input type="checkbox"/> handrim C</p>	<p><b>Knee lever brake with one-hand operation via Bowden cable</b></p> <p><input type="checkbox"/> A9000707020*      € 392.30</p> <p><input type="checkbox"/> Operation right side</p> <p><input type="checkbox"/> Operation left side</p>	<p><b>Parking knee lever brake for accompanying person</b> Push handle operation via Bowden cable</p> <p><b>for TRAVELER</b>  <input type="checkbox"/> A9000701020      € 464.10</p> <p><b>for SPEEDY A1 / A2</b>  <input type="checkbox"/> A9000701020      € 464.10</p> <p><b>for SPEEDY 4teen</b>  <input type="checkbox"/> A9000701020      € 464.10</p> <p><b>for BUDDY</b>  <input type="checkbox"/> A9000701020      € 464.10</p>