



## BUDDY 4all

### Active easy running wheelchair

Maximum user weight limit 100 kg

A9047100000 € 3,690.00

Quotation date: \_\_\_\_\_

Consultation was carried out by  
Staff at the medical supplies shop: \_\_\_\_\_

PRO ACTIV field representative: \_\_\_\_\_

Order date: \_\_\_\_\_

according to GTC PRO ACTIV GmbH

#### Medical suppliers/company address:

Cust. no.: \_\_\_\_\_  
Name: \_\_\_\_\_  
Street: \_\_\_\_\_  
Postcode / city: \_\_\_\_\_  
Phone: \_\_\_\_\_  
Fax: \_\_\_\_\_  
E-Mail: \_\_\_\_\_  
Contact: \_\_\_\_\_

Invoicing address:  same as company address

other (see below)

Cust. no.: \_\_\_\_\_  
Name: \_\_\_\_\_  
Street: \_\_\_\_\_  
Postcode / city: \_\_\_\_\_  
Phone: \_\_\_\_\_  
Fax: \_\_\_\_\_  
E-Mail: \_\_\_\_\_  
Contact: \_\_\_\_\_

Shipping address:  same as company  
address


Other (see below)

Name: \_\_\_\_\_  
Street: \_\_\_\_\_  
Postcode / city: \_\_\_\_\_

#### Customer name/commission:

Surname,  
forename: \_\_\_\_\_  
Street: \_\_\_\_\_  
Postcode / city: \_\_\_\_\_  
Phone: \_\_\_\_\_  
Fax: \_\_\_\_\_  
E-Mail: \_\_\_\_\_  
D.O.B: \_\_\_\_\_  
Paying  
organisation: \_\_\_\_\_

#### Body size:



Body height: \_\_\_\_\_ cm

Weight: \_\_\_\_\_ kg

Lower leg length: \_\_\_\_\_ cm  
(anatomic dimension: hollow of the knee to the foot sole):

Indication: \_\_\_\_\_

This information is mandatory for further processing.

# Basic size

<b>Seat width SB</b>	<input type="checkbox"/> 20	<input type="checkbox"/> 21	<input type="checkbox"/> 22	<input type="checkbox"/> 23	<input type="checkbox"/> 24	<input type="checkbox"/> 25	<input type="checkbox"/> 26	<input type="checkbox"/> 27	<input type="checkbox"/> 28	<input type="checkbox"/> 29	<input type="checkbox"/> 30	<input type="checkbox"/> 31
	<input type="checkbox"/> 32	<input type="checkbox"/> 33	<input type="checkbox"/> 34	<input type="checkbox"/> 35	<input type="checkbox"/> 36	<input type="checkbox"/> 37	<input type="checkbox"/> 38					
	The following dimensions are optionally available subject to surcharge:				<input type="checkbox"/> 39	<input type="checkbox"/> 40	<input type="checkbox"/> 41	<input type="checkbox"/> 42	<input type="checkbox"/> 43	<input type="checkbox"/> 44	<input type="checkbox"/> 45	<input type="checkbox"/> 46

	Size M (Frame length 35 cm (to the bend) and frame height 21 cm) <i>Not in conjunction with G-shape frame stem</i>										
	<input type="checkbox"/> 22	<input type="checkbox"/> 23	<input type="checkbox"/> 24	<input type="checkbox"/> 25	<input type="checkbox"/> 26	<input type="checkbox"/> 27	<input type="checkbox"/> 28	<input type="checkbox"/> 29	<input type="checkbox"/> 30	<input type="checkbox"/> 31	

<b>Seat depth ST</b>	Size L (Frame length 38 cm (to the bend) and frame height 21 cm)										
	<input type="checkbox"/> 26	<input type="checkbox"/> 27	<input type="checkbox"/> 28	<input type="checkbox"/> 29	<input type="checkbox"/> 30	<input type="checkbox"/> 31	<input type="checkbox"/> 32	<input type="checkbox"/> 33	<input type="checkbox"/> 34	<input type="checkbox"/> 35	
	Size XL (Frame length 43 cm (to the bend) and frame height 21 cm)										
	<input type="checkbox"/> 31	<input type="checkbox"/> 32	<input type="checkbox"/> 33	<input type="checkbox"/> 34	<input type="checkbox"/> 35	<input type="checkbox"/> 36	<input type="checkbox"/> 37	<input type="checkbox"/> 38	<input type="checkbox"/> 39	<input type="checkbox"/> 40	

<b>Back height RH</b>	<input type="checkbox"/> 20	<input type="checkbox"/> 21	<input type="checkbox"/> 22	<input type="checkbox"/> 23	<input type="checkbox"/> 24	<input type="checkbox"/> 25	<input type="checkbox"/> 26	<input type="checkbox"/> 27	<input type="checkbox"/> 28	<input type="checkbox"/> 29
	<input type="checkbox"/> 30	<input type="checkbox"/> 31	<input type="checkbox"/> 32	<input type="checkbox"/> 33	<input type="checkbox"/> 34	<input type="checkbox"/> 35	<input type="checkbox"/> 36	<input type="checkbox"/> 37	<input type="checkbox"/> 38	<input type="checkbox"/> 39
	<input type="checkbox"/> 40	<input type="checkbox"/> 41	<input type="checkbox"/> 42	<input type="checkbox"/> 43	<input type="checkbox"/> 44	<input type="checkbox"/> 45				
	Optional in size XL:			<input type="checkbox"/> 46	<input type="checkbox"/> 47	<input type="checkbox"/> 48				

<b>Lower leg length USL</b>	<input type="checkbox"/> 20	<input type="checkbox"/> 21	<input type="checkbox"/> 22	<input type="checkbox"/> 23	<input type="checkbox"/> 24	<input type="checkbox"/> 25	<input type="checkbox"/> 26	<input type="checkbox"/> 27	<input type="checkbox"/> 28	<input type="checkbox"/> 29	<input type="checkbox"/> 30	
	<input type="checkbox"/> 31	<input type="checkbox"/> 32	<input type="checkbox"/> 33	<input type="checkbox"/> 34	<input type="checkbox"/> 35	<input type="checkbox"/> 36	<input type="checkbox"/> 37	<input type="checkbox"/> 38	<input type="checkbox"/> 39	<input type="checkbox"/> 40	<input type="checkbox"/> 41	
	<input type="checkbox"/> 42	<input type="checkbox"/> 43	<input type="checkbox"/> 44									
	Optional in sizes L + XL:			<input type="checkbox"/> 45	<input type="checkbox"/> 46	<input type="checkbox"/> 47						

<b>Seat height back SHh</b>	for 20" drive wheel										
	<input type="checkbox"/> 31	<input type="checkbox"/> 32	<input type="checkbox"/> 33	<input type="checkbox"/> 34	<input type="checkbox"/> 35	<input type="checkbox"/> 36	<input type="checkbox"/> 37				
	for 22" drive wheel										
	<input type="checkbox"/> 33.5	<input type="checkbox"/> 34.5	<input type="checkbox"/> 35.5	<input type="checkbox"/> 36.5	<input checked="" type="checkbox"/> 37.5	<input type="checkbox"/> 38.5	<input type="checkbox"/> 39.5				
	for 24" drive wheel										
	<input type="checkbox"/> 36	<input type="checkbox"/> 37	<input type="checkbox"/> 38	<input checked="" type="checkbox"/> 39	<input type="checkbox"/> 40	<input type="checkbox"/> 41	<input type="checkbox"/> 42				

<b>Seat height front SHv</b>	<input type="checkbox"/> 40	<input type="checkbox"/> 41	<input type="checkbox"/> 42	<input type="checkbox"/> 43	<input type="checkbox"/> 44	<input type="checkbox"/> 45	<input type="checkbox"/> 46	<input type="checkbox"/> 47
------------------------------	-----------------------------	-----------------------------	-----------------------------	-----------------------------	-----------------------------	-----------------------------	-----------------------------	-----------------------------

<b>Tipping point KP</b> <small>(Front of back tube up to axis middle drive wheel)</small>	<input checked="" type="checkbox"/> 5	<input type="checkbox"/> 6	<input type="checkbox"/> 7	<input type="checkbox"/> 8	<input type="checkbox"/> 9	<input type="checkbox"/> 10
--	---------------------------------------	----------------------------	----------------------------	----------------------------	----------------------------	-----------------------------

Please cross on the appropriate dimension (all dimensions are specified in cm).  
Underlined values are standard settings.  
The values SHh and SHv may deviate a little from the specified values depending on the wheel size, wheel camber as well as the front and rear frame angle!

## Frame colour – one-colour design

Information about current colour designs and example images at: [www.proactiv-gmbh.de/farben.php](http://www.proactiv-gmbh.de/farben.php)

### frame colour: Selection from colour palette

Frame colour is specified when ordering

A9000100170-000

#### Powder coatings

##### White RAL 9016

gloss A9000100170-023

matt A9000100170-123 € 40.10

##### Brilliant silver

gloss A9000100170-011

##### Anthracite RAL 7016

gloss A9000100170-001

matt A9000100170-101 € 40.10

##### Black RAL 9005

gloss A9000100170-003

matt A9000100170-019 € 40.10

##### Pure orange RAL 2004

gloss A9000100170-022

matt A9000100170-122 € 40.10

##### Pink magenta RAL 4010

gloss A9000100170-037

matt A9000100170-137 € 40.10

##### Pastel blue RAL 5024

gloss A9000100170-036

matt A9000100170-136 € 40.10

##### Lemon green

gloss A9000100170-005

matt A9000100170-105 € 40.10

##### Yellow RAL 1023

gloss A9000100170-046

matt A9000100170-146 € 40.10

##### Red RAL 3001

gloss A9000100170-047

matt A9000100170-147 € 40.10

##### Lagoon blue

gloss A9000100170-052

matt A9000100170-152 € 40.10

##### Athletic grey

gloss A9000100170-050

matt A9000100170-150 € 40.10

##### Muddy grey

gloss A9000100170-051

matt A9000100170-151 € 40.10

#### “Dormant” powder coatings

##### Deep black, silver effect

gloss A9000100170-007

##### Sky blue

gloss A9000100170-017

##### Marine blue

gloss A9000100170-045

matt A9000100170-145 € 40.10

##### Electric green

gloss A9000100170-015

matt A9000100170-115 € 40.10

##### Violet

gloss A9000100170-044

##### Fire red

gloss A9000100170-018

matt A9000100170-118 € 40.10

## Frame colour – “liquid coating” one-colour design

€ 323.20

##### Storm (ice blue)

gloss A9000100170-500

matt A9000100170-501

##### Coffee

gloss A9000100170-520

matt A9000100170-521

##### Prairie (rose gold)

gloss A9000100170-504

matt A9000100170-505

##### Pearl (mother of pearl-white)

gloss A9000100170-524

matt A9000100170-525

##### Berry

gloss A9000100170-526

matt A9000100170-527

##### Blue ocean

gloss A9000100170-522

matt A9000100170-523

##### Emerald (emerald green)

gloss A9000100170-514

matt A9000100170-515

##### Mint

gloss A9000100170-528

matt A9000100170-529

##### Chilli

gloss A9000100170-516

matt A9000100170-517

## Frame colour – two colour designs

Note: Frames with two colour designs in the gloss level of the design colour.

**Two-colour design comprising one basic colour (on the back of the frame) and one design colour (on the front of the frame)** € 248.30

Please select from the following colours:

### Base colour (on the rear frame)

#### White RAL 9016

A9000100170-200

#### Brilliant silver

A9000100170-201

#### Anthracite RAL 7016

A9000100170-202

#### Black RAL 9005

A9000100170-203

#### Pastel blue RAL 5024

A9000100170-206

#### Lemon green

A9000100170-207

#### Yellow RAL 1023

A9000100170-208

#### Lagoon blue

A9000100170-212

#### Athletic grey

A9000100170-210

#### Muddy grey

A9000100170-211

### Design colour (on the front frame)

#### Black RAL 9005

gloss A9000100170-259

matt A9000100170-260

#### Lemon green

gloss A9000100170-250

matt A9000100170-251

#### Chilli

gloss A9000100170-253

matt A9000100170-254

#### Pure orange RAL 2004

gloss A9000100170-255

matt A9000100170-256

#### Coffee

gloss A9000100170-278

matt A9000100170-279

#### Mint

gloss A9000100170-286

matt A9000100170-287

#### Emerald (emerald green)

gloss A9000100170-274

matt A9000100170-275

#### Storm (ice blue)

gloss A9000100170-264

matt A9000100170-265

#### Prairie (rose gold)

gloss A9000100170-266

matt A9000100170-267

#### Pearl (mother of pearl-white)

gloss A9000100170-282

matt A9000100170-283

#### Berry

gloss A9000100170-284

matt A9000100170-285

#### Blue ocean

gloss A9000100170-280

matt A9000100170-281

## Colour - Attachment parts

Note: If no specifications have been made, the caster fork and bearing block are mounted in black.

### Drive wheel: Aluminium rim, hub and spokes

Black A9000100171-003

### Caster fork and bearing block

Silver A9000100173-011

Black A9000100173-003

## Attachment parts colour design according to customer desires:

### Drive wheel: Aluminium rim and design inlet<sup>1</sup>

Pure orange RAL 2004 A9000100171-022 € 193.10

White RAL 9016 A9000100171-023 € 193.10

Frame colour A9000100171-999 € 230.80

Selection of single colour coating (deviating from the frame colour) A9000100171-998 € 267.10

Colour: \_\_\_\_\_

A9000100170-

### Caster wheel: 5" aluminium rim<sup>2</sup> (3-spoke design)

Pure orange RAL 2004 A9000100174-022 € 44.40

White RAL 9016 A9000100174-023 € 44.40

Frame colour A9000100174-999 € 69.50

Selection of single colour coating (deviating from the frame colour) A9000100174-998 € 80.60

Colour: \_\_\_\_\_

A9000100170-

## Special colour requests (frame or attachment parts)

### Frame colour individual procurement Selection from RAL colour palette

A9000100170-100 € 383.10

### Degree of gloss

gloss

matt

### Colour designation

\_\_\_\_\_

### RAL number

\_\_\_\_\_

**Note:** Attachment parts **not** possible in design colours

<sup>1</sup> Only in conjunction with Liteline hub (for drive wheel aluminium rim with 18 or 36 spokes, object).


<sup>2</sup> Only in conjunction with 5" caster wheel aluminium rim 3-spoke design easy running A9000501200 or good shock absorption A9000501210

## Frames

### Straight front end


(Only in frame sizes L and XL)

- A9047100010-K\*  
Knee angle 105°  
Bearing block permanently welded



### Front end abducted shape 3.5 cm / side

- A9047100030-K\*  
Knee angle 105°  
Bearing block permanently welded  
Abduction starts at:  
Size M: ST 26 cm  
Size L: ST 29 cm  
Size XL: ST 34 cm



### Custom-made seat width dimension 39 - 46 cm

- A9047100060 € 392.20

With inlet in the upper cross tubes,  
lower cross tube double bolted,  
side panel screw connection via the  
reinforcement inlet

## Push handles

### Safety push handles offset to the rear

- Length: 45 cm  
A9041401075-K001  
max. height = SHh + approx. 53 cm

- Length: 55 cm  
A9041401075-K002  
max. height = SHh + approx. 63 cm  
continuous height adjustment  
Safety push handle holder mounted on the  
back cross bar



### Safety push handles aluminium offset to the rear with side handle offset 145 mm

- A9041401080-K\* € 82.80

continuous height adjustment  
Safety push handle holder mounted on the  
back cross bar  
Length: 55 cm  
Up to seat width of 38 cm or body weight of  
60 kg



### Aluminium central push handle Foldable handle area

- A9041401086-K € 94.90


Height-adjustable and foldable in the handle  
area via toothed joint  
Holder welded on the middle on the back cross  
bar  
Up to seat width of 38 cm or body weight of  
50 kg  
Not compatible with headrest fixture  
Note: Not suitable for carrying the wheelchair -  
to do this, grasp the back cross bar.



### Central push handle for brake handle folding gripping area

- A9041401070-K € 131.20


Closed, with continuous height adjustment,  
with foam rubber handle cover.  
Up to seat width of 38 cm or body weight of  
50 kg  
Not compatible with headrest fixture  
Only in conjunction with parking knee lever  
brake for accompanying person or drum brake



### Push handle bar aluminium, closed, foldable via toothed joint

- A9041401069-K € 360.50

Continuous height adjustable push handle  
holder back-positioned  
Up to seat width of 38 cm or body weight of  
60 kg



## Backrest

### Backrest fixed

- A9047400020  
bolted with aluminium side sections  
Back angle 90° to the ground  
Backrest cannot be folded



### Folding backrest

- A9000400023 € 283.20  
With integrated plastic coated aluminium side  
sections  
7-way inclination adjustment



## Back equipment

### Adaptable back with back padding

A9041802010-K  
Cover: Black, breathable, open surface  
Equipment with back cross bar



### Ergonomic aluminium back shell from RH 24 cm

A9000802210-K € 408.30  
With lumbar support  
Colour: Anthracite matt  
(Where not specified under colour selection)  
Back cushions, if desired, please select separately



### Back cushion for ergonomic back shell

A9000802220-K € 201.80  
Cushion thickness: 2 cm  
Cover:  
 Black, breathable, open surface  
 Black, breathable, closed surface  
If not specified, the cover fabric "black, breathable, open surface" is delivered.



### ERGO Back backrest bar for ergonomic back shell from RH 24 cm and SB 30 cm

A9000802019-K € 224.30  
With recessed gripping area  
Colour: Anthracite matt  
(Where not specified under colour selection)  
Only compatible with central push handle art. no. A9041401082  
Back shell and back cushions, if desired, please select separately



### Equipment without back upholstery and back padding or back system


A9000802011 € -53.00  
Equipment with back cross bar



## Seat upholstery


### Body Contour

A9000801025-K



### Aluminium seat plate

A9000801210-K € 270.70  
Assembly between the frame tubes, on bolted welded seat plate holders  
Colour: Anthracite matt  
(Where not specified under colour selection)



### Equipment without seat upholstery

A9000801050



## Colours of back equipment and seat

Information about current colour designs and availability at [www.proactiv-gmbh.com](http://www.proactiv-gmbh.com)  
Unless otherwise specified, the seat and back system is supplied in the standard colour of anthracite matt.

### Anthracite RAL 7016 gloss

A9000100179-995  
 Back shell  
 ERGO Back backrest bar  
 Aluminium seat plate

### in frame colour

A9000100179-999 € 85.30  
 Back shell  
 ERGO Back backrest bar  
 Aluminium seat plate

**For frame coating in two-colour design**

Please specify whether the desired option is to be coated in the basic colour or the design colour


- A9000100179-999 € 85.30
- Back shell
- Base colour
- Design colour
- ERGO Back backrest bar
- Base colour
- Design colour
- Aluminium seat plate
- Base colour
- Design colour

**If another colour is desired, please select this from the wheelchair colour scheme****In colour from one-colour design:**


- A9000100179-998 € 138.10
- Please enter the last 3 digits of the desired colour:
- A9000100170-
- Back shell
- ERGO Back backrest bar
- Aluminium seat plate

**In colour from two-colour design:**

- A9000100179-997 € 223.60
- Please enter the last 3 digits of the desired colour:
- A9000100170-
- Back shell
- ERGO Back backrest bar
- Aluminium seat plate

**Seat cushion****Seat cushion in rectangular contour  
Inlet: Standard**

- A9000805180-K\* € 106.60

**Seat cushion with radius contour  
Inlet: Standard**

- A9000800300-K\* € 148.00

**Seat cushion height:**

- 3 cm       5 cm
- 4 cm       6 cm

**Seat cushion in rectangular contour  
Inlet: Anti-Decubitus**

- up to 30x30 cm  
A9000805204-K\* € 276.20
- >30x30 up to  
50x50 cm  
A9000805200-K\* € 401.90

**Seat cushion with radius contour  
Inlet: Anti-Decubitus**

- up to 30x30 cm  
A9000800311-K\* € 326.50
- >30x30 up to  
50x50 cm  
A9000800310-K\* € 443.70

**Seat cushion height:**

- 4 cm
- 6 cm

**Cover:**

- Black, breathable, open surface
- Black, breathable, closed surface
- If not specified, the cover fabric "black, breathable, open surface" is delivered.

 **Seat cushion is existing**

Brand/type \_\_\_\_\_

Cushion height \_\_\_\_\_ cm

You can also find the seat cushions on our separate "PRO ACTIV Seat and Back Systems" order sheet

**Clothing guard****Aluminium clothing guard**

Permanently mounted on the side guard  
Only in combination with fixed backrest



In black gloss

- A9000403100-K € 224.90

in frame colour (basic colour with two-colour design)

- A9000403101-K € 287.60

in design colour (only possible with frame colour in two-colour design)

- A9000403102-K € 319.00

### Aluminium clothing guard removable via lock function

Adjustable following wheel arch  
Only in combination with folding backrest



in anthracite matt

A9000403215-K € 224.90

in frame colour (basic colour with two-colour design)

A9000403216-K € 287.60

in design colour (only possible with frame colour in two-colour design)

A9000403217-K € 319.00

### Carbon clothing guard removable via lock function

A9000400013-K € 521.20

Adjustable following wheel arch  
Only in combination with folding backrest



**Knee lever brake system, brake lever integrated in the clothing guard see "Brake systems"**

## Armrests

### Armrest, height adjustable, integrated in the clothing guard

Including aluminium clothing guard  
with arm padding 25 cm long  
Removable via lock function  
Only in combination with folding backrest



in anthracite matt

A9000400014-K € 526.70

In frame colour (basic colour with two-colour design)

A9000400015-K € 575.50

in design colour (only with frame colour in two-colour design)

A9000400016-K € 599.20

### Armrest, height adjustable, integrated in the clothing guard

A9000400017-K € 398.50

Including carbon clothing guard  
with arm padding 25 cm long  
Removable via lock function  
Only in combination with folding backrest

## Axle plate

Standard equipment with 12 mm sockets/axle. Further axle plate on request.


### 2-slot axle plate

A9000201000\*  
Quick-release axle Ø 12 mm  
Rear SH:  
24" wheel: approx. 36 - 42 cm  
22" wheel: approx. 33.5 - 39.5 cm  
20" wheel: approx. 31 - 37 cm




### 1-slot axle plate

A9000201001\*  
Only available in black  
(no colour selection)  
Only in combination with 10° camber  
Rear SH can only be adjusted in some cases



### Extended bike plate

A9000201030\* € 291.40  
+ 4° camber integrated  
Adjustable track  
Rear SH:  
24" wheel: approx. 36 - 42 cm  
22" wheel: approx. 33.5 - 39.5 cm  
20" wheel: approx. 31 - 37 cm



## Wheel camber

\*Total width with tight-mounted handrim

### 1° drive wheel camber

A9000202010\*  
Total width = approx. SB + 20 cm\*

### 4° drive wheel camber

A9000202030\*  
Total width = approx. SB + 26 cm\*

### 6° drive wheel camber

A9000202040\*  
Total width = approx. SB + 30 cm\*

### 8° drive wheel camber

A9000202080\* € 71.20  
Total width = approx. SB + 33 cm\*  
Wheel plate mounted on the outside

### 10° drive wheel camber

A9000202050\* € 126.60  
Total width = approx. SB + 39 cm\*  
Wheel plate mounted on the outside



## E-adaptations

### Preparation for "e-motion" mounting

A9000201200\* € 465.00

Equipment: e-motion torque support, 1/2" axle and socket length 50.8 mm (2"). PRO ACTIV and Alber recommend anti-tipping supports are installed in pairs. Please choose "Adaptations" under the selection or order Alber anti-tipping supports directly from Alber.

### Preparation for "e-fix" mounting

A9000201210\* € 465.00

Equipment: e-fix torque support, 1/2" axle and socket length 50.8 mm (2"). PRO ACTIV and Alber recommend anti-tipping supports are installed in pairs. Please choose "Adaptations" under the selection or order Alber anti-tipping supports directly from Alber, as well as BUDDY mount.

Wheelchair delivery to

Alber GmbH  
(Please also enter on p.1 under "Shipping address, other")

Specialist dealer

## Brake system

### Brake system push to lock Standard brake handle

A9000701010-001



### Brake system push to lock Long brake lever

A9000701010-002



### Brake lever at 20°

A9000700002  
Only in conjunction with knee lever brake system, standard or long brake lever



### Brake system push to lock Brake lever with M8 adapter and plastic ball

A9000706001 € 154.70

Only in conjunction with knee lever brake system, standard or long brake lever

Brake lever:

Standard

Long



### Brake system push to lock Brake handle to flip-down

A9000706008 € 146.70

Length: 9 cm



### Knee lever brake with one-hand operation

A9000707020-K € 362.80

via Bowden cable  
Operation:

Right

Left



### Knee lever brake system, brake lever integrated in the clothing guard

Brake lever mounted on the bottom frame tube  
Operation via wire rope  
Clothing guard cannot be removed.  
Not possible with frame size M and tipping point >7 cm



in black gloss for fixed backrest  
in anthracite matt with folding backrest




A9047701015-K € 594.60

in frame colour (basic colour with two-colour design)








A9047701016-K\* € 656.80

in design colour (only possible with frame colour in two-colour design)

A9047701017-K\* € 688.20

<p><b>Parking knee lever brake for accompanying person</b></p> <p><input type="checkbox"/> A9000701020      € 429.10</p> <p>Push handle operation via Bowden cable.</p> 	<p><b>Drum brake</b></p> <p><input type="checkbox"/> A9000702011-K*      € 438.30</p> <p>Operating from the push handle Total width + 3cm (see wheel camber) Please note: - Only in combination with drive wheel aluminium 36-spoke rim wheel - Not in combination with one-hand drive, e-fix or e-motion</p> 	<p><b>Knee lever brake with lock</b></p> <p><input type="checkbox"/> Standard brake handle A9000700214      € 146.70</p> <p><input type="checkbox"/> Long brake lever A9000700215      € 146.70</p> <p>Locking block for brake lever Prevents braking during pushing by accompanying persons Not possible in combination with knee lever brake with one-hand operation.</p> 
---	---	---

## Footrest

<p><b>Continuous footrest</b></p> <p><input type="checkbox"/> A9041601010-K* with tubular resting bracket height-adjustable from USL 26 cm to USL 47 cm possible For equipment of the footrest with "Rubber coated footplate support" please select under "Footrest accessories".</p> 	<p><b>Footrest, positioned on the upper frame</b></p> <p><input type="checkbox"/> A9041606010-K height-adjustable from USL 20 cm to USL 25 cm possible, Rubber coated footplate support required, please select this under "Footrest accessories".</p> 	<p><b>Variable position foot support</b></p> <p><input type="checkbox"/> A9047600001-K      € 259.90</p> <p>adjustable in height, depth and angle incl. footplate support covered with anti-slip film</p> 
<p><b>Footrest folding back to the rear with spring latch mechanism</b></p> <p><input type="checkbox"/> A9047602005-K*      € 321.00</p> <p>with continuous tubular resting bracket from USL 32 cm including caster wheel track extension 20 mm per side For equipment of the footrest with "Rubber coated footplate support" please select under "Footrest accessories".</p> 	<p><b>Variable position foot support and folds up to the rear</b></p> <p><input type="checkbox"/> A9047600002-K      € 345.00</p> <p>adjustable in height, depth and angle incl. footplate support covered with anti-slip film</p> 	<p><b>Footrest folding up to one side</b></p> <p><input type="checkbox"/> A9041601300-K*      € 160.80</p> <p>from USL 25 cm: with continuous tubular resting bracket height adjustable +/- 3 cm cannot be removed from USL 29 cm: can also be pivoted outwards For equipment of the footrest with "Rubber coated footplate support" please select under "Footrest accessories".</p> 
<p><b>Divided footrest</b></p> <p><input type="checkbox"/> A9041603200-K      € 311.40</p> <p>with rubber coated footplate supports height-adjustable can be pivoted outwards and folded up to the side cannot be removed from USL 29 cm and SB 26 cm possible</p> 		

## Footrest accessories

### Rubber coated footplate, 14 cm deep

- A9000610000-K € 116.40  
bolted, angle-adjustable



### Rubber coated footplate, 18 cm deep

- A9000610005-K € 129.50  
bolted, angle-adjustable



### Footplate side stop

- A9000610100 € 55.20  
adjustable  
only in combination with rubber coated footplate!



### Deflector rollers mounted on the footplate

- A9000600006 € 108.60



## Adaptations

### Anti-tipping support

- right side A9041301160-K001 € 184.00  
 left side A9041301160-K002 € 184.00  
with tension spring  
can be folded up  
removable using Quick-pin



### Outdoor wheel attachment with steering fork

- A9041900100 € 782.00  
Equipped with straight-running lock  
can be removed from the wheelchair via quick release lever  
with swivelling pneumatic tyre 12.1/2 x 2.1/4" caster wheel  
can be removed from the support tube via quick-release axle  
for maximum ground pressure at the steering attachment of 30 kg



### Transport restraint connection integrated on the rear frame with belt eye

- A9047900001 € 337.90  
Frame-welded part for mounting a lap belt and the rear retractor. Lap belt and buckle (A9000900001 / A9000900007), please select separately if required.



### Bearing block intermediate plate with belt eye for mounting the front retractor

- Colour: black A9000900008-001 € 168.10  
 Colour: Silver A9000900008-002 € 168.10

Please note:  
The wheel track guide increases by 12 mm per side.



### Headrest with segment padding

A9000900004 € 384.30

3-part  
and support tube ø 25 mm



### Headrest fixture

pipe dimensions support  
tube ø 25 mm € 258.60  
A9000900005-K001

For mounting the headrest segment padding  
Please select separately (A9000900004).  
Reinforced version with double cross-tube.

pipe dimensions support  
tube 15 x 15 mm € 207.60  
A9000400850-001

pipe dimensions support  
tube  
12.7 x 12.7 mm € 207.60  
A9000400850-002

with square holding fixture for headrests.  
Assembly on the back cross bar.  
Reinforced version with double cross-tube.  
Possible from SB 23 cm.  
Not in combination with the central push  
handle.



You will find our PRO ACTIV transport restraint systems  
on our separate "Transport restraint systems" order sheet



### 3" caster wheel

#### 3" Skate wheel (80 mm)

A9000501110 € 65.40

puncture-proof



### 4" caster wheel

#### 4" caster wheel "Standard" solid rubber black

A9000501011-004  
(82A), puncture-proof



#### 4" Skate wheel (100 mm)

A9000501120 € 65.40  
puncture-proof



#### 4" safety light-up wheel (100 mm)

A9000501220 € 87.60  
puncture-proof



#### 4" Frog Legs caster wheel (100 mm)

A9000501230 € 196.50

with replaceable rubber 6 mm, puncture-proof  
Colour: Black, silver rim  
Incl. hair guard spacer (5.3 mm)



## 5" caster wheel

### 5" caster wheel "Standard" PU black

A9000501036-005  
(82A), puncture-proof



### 5" Caster wheel "Standard" PU grey

A9000501031-005  
(82A), puncture-proof



### 5" Caster wheel "2 components"

A9000501030      € 52.00  
puncture-proof



### 5" safety light-up wheel (125 mm)

A9000501221      € 87.60  
puncture-proof



### 5" Aluminium rim 3-spoke design "easy running"

Silver, anodised      € 160.50  
A9000501200

black      € 160.50  
A9000501201

With exchangeable rubber ring 3 mm  
hard, puncture-proof  
Incl. hair guard spacer (3 mm)



### 5" aluminium rim 3-spoke design "Good shock absorption"

Silver, anodised      € 166.30  
A9000501210

black      € 166.30  
A9000501211

with replaceable rubber 6 mm medium hard,  
puncture-proof  
Incl. hair guard spacer (3 mm)



### 5" Frog Legs caster wheel (125 mm)

A9000501232      € 196.50  
with replaceable rubber 6 mm, puncture-proof  
Colour: Black, silver rim  
Incl. hair guard spacer (5.3 mm)



## 6" caster wheel

### 6" Caster wheel "Standard" PU grey

A9000501010-006  
(82A), puncture-proof



### 6" caster wheel "pneumatic tyre"

A9000501050      € 52.90  
SH front min. 42 cm



## Caster fork

Caster fork colours, see "Component colours"

### Caster fork Standard

A9000500100-K001  
with screwed axle

A9000500100-K002      € 154.80  
with quick-release axle



### Heavy duty caster fork

For high loads or  
when using additional electric drives

A9000510111-K001      € 49.90  
with screwed axle

A9000510111-K002      € 202.40  
with quick-release axle



### Suspension fork

incl. 5" aluminium rim "good shock absorption"

A9000500160-K001      € 466.90  
with screwed axle

A9000500160-K002      € 570.40  
with quick-release axle



### Caster fork spare wheel standard

A9000500300-K      € 274.40  
Standard caster forks with quick-release axles



### Caster fork spare wheel Heavy Duty

A9000500301-K      € 312.50  
Heavy duty caster forks with quick-release  
axles



The spare caster fork does not have a caster  
wheel.  
In case a caster wheel is requested, please put  
a cross on a second selection in the category  
"Caster wheels".

## Track distance

### Caster wheel track extension



A9000500400\*      € 54.70  
10 mm per side

A9000500401\*      € 54.70  
15 mm per side

A9000500410\*      € 54.70  
20 mm per side  
For footbar that folds up to the rear with spring  
locking mechanism A9047602005\* including

## Drive wheels

### Drive wheel, 18 spoke aluminium rim

- 20" (ETRTO 451-19)  
A9000203052
- 22" (ETRTO 489-19)  
A9000203062
- 24" (ETRTO 540-19)  
A9000203122

with Liteline hub, spoked radially

Not in combination with:

- Puncture-proof tyres
- One-hand drive
- Drum brake
- Adaptive bike

only up to a body weight of 70 kg:

Sturz \ Gewicht	1°	4°	6°	8°	10°
bis 50 Kg	✓	✓	✓	✓	✓
51 - 60 Kg	✓	✓	✓	✓	
61 - 70 Kg	✓	✓			



### Drive wheel, 36 spoke aluminium rim

- 20" (ETRTO 451-19)  
A9000203050
- 22" (ETRTO 489-19)  
A9000203060
- 24" (ETRTO 540-19)  
A9000203000

Up to 120 kg body weight with Liteline hub (left image)

Over 120 kg body weight with aluminium high flange (right image) radially spoked



### Drive wheel double rim 36 spokes

- 24" (ETRTO 540)  
A9000203100
- with aluminium front hub radially spoked  
Use with high loads

€ 145.90



### Spinergy running wheel "Spox Sport" with 24 PBO spokes Rim, hub and spokes in black

- 24" (ETRTO 540)  
A9000203550

€ 725.00



### Spinergy "Everyday Light Extreme" with 12 PBO spokes Rim and hub in black, spokes in white

- 24" (ETRTO 540)  
A9000203554

€ 885.00



### Spinergy "Everyday Light Extreme" with 12 PBO spokes Rim, hub and spokes in black


- 24" (ETRTO 540)  
A9000203552

€ 885.00



## Handrims

x = dimension handrim attachment



Handrim attachment

General:

x = 11 mm

x = 19 mm

Handrim attachment

Curve handrim:

x = 10 mm

x = 18 mm

Unless otherwise specified, the handrim is mounted closely (11 mm / 10 mm)

### Stainless steel Ø 20 mm

20" handrim for  
20" wheel  
A9000205010-001

22" handrim for  
22" wheel  
A9000205010-002

24" handrim for  
24" wheel  
A9000205010-004



### Aluminium Ø 19 mm silver anodised

24" handrim for  
24" wheel  
A9000205019-004



### Aluminium Ø 16 mm powder coated in the frame colour

(basic colour with two-colour design)

= Gloss surface

= Matt surface

Note: Matt coating is not recommended due to the low grip.

20" handrim for  
20" wheel  
A9000205021-001 **€ 63.10**

22" handrim for  
22" wheel  
A9000205021-002 **€ 63.10**



### Aluminium ø 19 mm powder coated in frame colour

(basic colour with two-colour design)

= Gloss surface

= Matt surface

Note: Matt coating is not recommended due to the low grip.

20" handrim for  
20" wheel  
A9000205020-001 **€ 63.10**

22" handrim for  
22" wheel  
A9000205020-002 **€ 63.10**

24" handrim for  
24" wheel  
9000205020-004 **€ 63.10**



### Aluminium ø 16 mm polished and hard-anodised

20" handrim for  
20" wheel  
A9000205031-001 **€ 152.50**

22" handrim for  
22" wheel  
A9000205031-002 **€ 152.50**



### Aluminium ø 19 mm polished and hard-anodised










20" handrim for  
20" wheel  
A9000205030-001 **€ 152.50**

22" handrim for  
22" wheel  
A9000205030-002 **€ 152.50**

24" handrim for  
24" wheel  
A9000205030-004 **€ 152.50**





<p><b>Titanium ø 19 mm</b></p> <p><input type="checkbox"/> 24" handrim for 24" wheel A9000205040-004      <b>€ 310.30</b></p> 	<p><b>Silicone handrim cover</b></p> <p><input type="checkbox"/> for 20" handrim A9000205100-001      <b>€ 215.20</b></p> <p><input type="checkbox"/> for 22" handrim A9000205100-002      <b>€ 215.20</b></p> <p><input type="checkbox"/> for 24" handrim A9000205100-004      <b>€ 215.20</b></p> <p>Surface without structure compatible with handrim Ø 19 mm or Ø 20 mm</p> 	<p><b>Ultra-Grip "very grippy"</b></p> <p><input type="checkbox"/> 20" handrim for 20" wheel A9000205070-001      <b>€ 372.10</b></p> <p><input type="checkbox"/> 22" handrim for 22" wheel A9000205070-002      <b>€ 372.10</b></p> <p><input type="checkbox"/> 24" handrim for 24" wheel A9000205070-004      <b>€ 372.10</b></p> 
<p><b>Ergo-Grip</b></p> <p><input type="checkbox"/> 20" handrim for 20" wheel A9000205080-001      <b>€ 377.90</b></p> <p><input type="checkbox"/> 22" handrim for 22" wheel A9000205080-002      <b>€ 377.90</b></p> <p><input type="checkbox"/> 24" handrim for 24" wheel A9000205080-004      <b>€ 377.90</b></p> <p>Attachment fitting to the rim</p> 	<p><b>Maxgripp with cell and natural rubber cover very good grip, soft</b></p> <p><input type="checkbox"/> 20" handrim for 20" wheel A9000200788      <b>€ 377.90</b></p> <p><input type="checkbox"/> 22" handrim for 22" wheel A9000200789      <b>€ 377.90</b></p> <p><input type="checkbox"/> 24" handrim for 24" wheel A9000200787      <b>€ 377.90</b></p> 	<p><b>Curve black handrim anodised, ergonomic handrim profile</b></p> <p><input type="checkbox"/> 24" handrim for 24" wheel H 24 mm A9000205130-004      <b>€ 370.60</b></p> <p><input type="checkbox"/> 24" handrim for 24" wheel H 31 mm size L A9000205140-004      <b>€ 370.60</b></p> 
<p><b>Curve silver anodised ergonomic handrim profile</b></p> <p><input type="checkbox"/> 24" handrim for 24" wheel H 31 mm, size L A9000205160-004      <b>€ 375.50</b></p> 	<p><b>Curve powder-coated, ergonomic handrim profile</b></p> <p><input type="checkbox"/> = Frame colour (gloss surface) <input type="checkbox"/> = black (gloss surface)</p> <p><input type="checkbox"/> 24" handrim for 24" wheel H 24 mm A9000205135-004      <b>€ 432.40</b></p> <p><input type="checkbox"/> 24" handrim for 24" wheel H 31 mm size L A9000205145-004      <b>€ 432.40</b></p>  <p>Profilhöhe H: 24mm oder 31mm</p>	<p><b>Curve Grip, with Tetra-Grip coating ergonomic handrim profile</b></p> <p><input type="checkbox"/> 24" handrim for 24" wheel H 24 mm A9000205155-004      <b>€ 481.40</b></p> 

<p><b>Quadro with Tetra-Grip coating</b></p> <p><input type="checkbox"/> 24" handrim for 24" wheel A9000205170-004      <b>€ 492.90</b></p> 	<p><b>Gekko black anodised</b></p> <p><input type="checkbox"/> 24" handrim for 24" wheel A9000205175-004      <b>€ 477.30</b></p> <p>with a robust and replaceable silicone lip, fully integrated into the profile Profile height 31 mm</p> 	<p><b>Gekko silver anodised</b></p> <p><input type="checkbox"/> 24" handrim for 24" wheel A9000205185-004      <b>€ 488.80</b></p> <p>with a robust and replaceable silicone lip, fully integrated into the profile Profile height 31 mm</p> 
<p><b>Gekko S black anodised</b></p> <p><input type="checkbox"/> 22" handrim for 22" wheel A9000205195-002      <b>€ 477.30</b></p> <p><input type="checkbox"/> 24" handrim for 24" wheel A9000205195-004      <b>€ 477.30</b></p> <p>with a robust and replaceable silicone lip, fully integrated into the profile Profile height 24 mm</p> 	<p><b>Equipment without handrims</b></p> <p><input checked="" type="checkbox"/> A9000205990</p>	
<p><b>Spoke guard</b></p>		
<p><b>Finger- and spoke guard "transparent"</b></p> <p><input type="checkbox"/> for 20" wheel A9000206300-001      <b>€ 191.50</b></p> <p><input type="checkbox"/> for 22" wheel A9000206300-002      <b>€ 191.50</b></p> <p><input type="checkbox"/> for 24" wheel A9000206300-004      <b>€ 191.50</b></p> 		

## Tyres

### Schwalbe "RightRun" High-pressure tyre

22x1" (489-25)  
A9000204001-002

24x1" (540-25)  
A9000204001-004

Colour: black  
7.5 bar



### Schwalbe "RightRun" High-pressure tyre

20x1" (451-25)  
A9000204000-001

22x1" (489-25)  
A9000204000-002

24x1" (540-25)  
A9000204000-004

Colour: Grey, black wall  
7.5 bar



### Schwalbe "DownTown" Profiled tyres

24x1" (540-25)  
A9000204010-004

7.5 bar



### Schwalbe "MarathonPlus Evolution" with "smart guard" puncture protection and reflector strips

22x1" (489-25) € 45.50  
A9000204030-002

24x1" (540-25) € 45.50  
A9000204031-004

Colour: black  
10 bar



### Schwalbe "MarathonPlus Evolution" with "smart guard" puncture protection

24x1" (540-25) € 45.50  
A9000204030-014

Colour: black  
10 bar



### Schwalbe "MarathonPlus" with "smart guard" puncture protection

24x1" (540-25) € 45.50  
A9000204030-004

Colour: Grey, black wall  
10 bar



### Schwalbe "One"

24x1" (540-25) € 86.00  
A9000204123-004

Colour: black  
9 bar



### Puncture-proof tyre "solid rubber tyres"

24x1" (540-25) € 113.40  
A9000204201-004

Colour: black  
Not in combination with the 18 spoke drive  
wheel, Spinergy and one-hand drive.



## Accessories

### 24" MTB spare wheel set

with double rim (ETRTO 507), aluminium front hub, quick-release axle and stainless steel handrim  
Mountain bike tyres, selectively:

A9000203410 € 650.40

24x1.75" Schwalbe "Marathon Plus"  
Brake setting compatible with 24x1" wheel



A9000203413 € 650.40

24x2.1" Schwalbe "SMART SAM Sport"  
Brake setting not compatible with 24x1" wheel



NOTE: We charge a surcharge for the use of drum brakes

A9000200004 € 210.60

### Through axle with assistance, Large push button

A9000200096 € 97.80

Interlock dimension: 122 / 110 mm  
Axle diameter: 12 mm  
Not in conjunction with lockable drive wheel axle clip "Tetra Clip",  
Not in combination with one-hand drive



### Hand air compressor with pressure gauge

A9000901100 € 98.30

up to 14 bar



### Length-adjustable safety belt with buckle

A9000902100 € 74.60



### Adjustable safety belt with hook and fastener

A9000902110 € 59.60



### Crutch holder with nylon cup

A9000902310 € 127.00

flexible






### SIGMA lighting set

A9000903270 € 99.80

Consisting of front light and rear lights, chargeable via Micro USB socket, incl. Micro USB charging cable. Lighting set compliant with road traffic regulations



<p><b>Transport bag for drive wheels</b></p> <p><input type="checkbox"/> for 24" wheel          Ø 76 cm x 2 cm          A9000903530-004</p> <p style="text-align: right;"><b>€ 71.90</b></p> 	<p><b>Accessory bag</b></p> <p><input type="checkbox"/> A9000903600          17 cm x 13cm x 5 cm</p> <p style="text-align: right;"><b>€ 70.30</b></p> 	<p><b>Therapy table</b></p> <p><input type="checkbox"/> Width 350 mm          (for SB up to 37 cm) <b>€ 444.00</b>          A9000902800-001</p> <p><input type="checkbox"/> Width 400 mm          (SB 38 - 41 cm) <b>€ 444.00</b>          A9000902800-002</p> <p><input type="checkbox"/> Width 420 mm          (from SB 42 cm) <b>€ 444.00</b>          A9000902800-003</p> <p>removable table plate from Lexan, transparent,          with front edging, with 2 items Positioning stops          Ø 25 mm          (for adaptation to frame tube diameter          Ø 25 mm)</p> 
<p><b>Further accessories, e.g., adaptive bikes, steering attachments, etc. see separate order sheet or <a href="http://www.proactiv-gmbh.com">www.proactiv-gmbh.com</a></b></p>		

**Notes and remarks:**