



Quotation date: _____
Consultation was carried out by
Staff at the medical supplies shop: _____
PRO ACTIV field representative: _____

Order date: _____
according to GTC PRO ACTIV Reha-Technik GmbH

Medical suppliers/company address:

Cust. no.: _____
Name: _____
Street: _____
Postcode,
city: _____
Phone: _____
Fax: _____
E-Mail: _____
Contact: _____

Shipping address: **same as company address**
 Other (see below)

Name: _____
Street: _____
Postcode,
city: _____

Customer name/commission:

Surname, _____
forename: _____
Street: _____
Postcode, city: _____
Phone: _____
Fax: _____
E-Mail: _____
D.O.B: _____
Paying
organisation: _____
Indication: _____

One-hand drive with a double handrim


Maximum user weight limit 120 kg

Invoicing address: **same as company address**
 other (see below)

Cust. no.: _____
Name: _____
Street: _____
Postcode,
city: _____
Phone: _____
Fax: _____
E-Mail: _____
Contact: _____

Version, mounting with new wheelchair orders

Please note: In the case of a wheelchair with one-hand drive, the wheelchair attachment parts are only available in **black**.



Wheel camber: 1°, 4°, 6°, 8°, 10°, 12° for TRAVELER, SPEEDY A1 / A2, BUDDY classic / 4all
 Wheel camber: 2° for SPEEDY 4all / 4you / 4all Ergo / 4you Ergo, SPEEDY 4teen

TRAVELER incl. detachable drive shaft <input type="checkbox"/> A9030208010* € 1,391.50	SPEEDY 4all / 4you / 4all Ergo / 4you Ergo / SPEEDY 4teen <input type="checkbox"/> A9016208010* € 1,225.90	SPEEDY A1 / A2 / BUDDY classic / BUDDY 4all <input type="checkbox"/> A9010208010* € 1,225.90
--	---	---

Please select the desired operating side:

Operation: Right Left

Version, mounting with retrofitting

Retrofit kit for PRO ACTIV wheelchair

incl. stainless steel handrim, high-pressure smooth-running or high-pressure profiled tyres.
 For further choices at extra cost, see wheelchair order sheet.

<input type="checkbox"/> On SPEEDY A1 / A2 <input type="checkbox"/> On BUDDY classic <input type="checkbox"/> E8000208100 € 1,942.80	On TRAVELER incl. detachable drive shaft <input type="checkbox"/> E8000208200 € 2,225.30	On BUDDY 4all <input type="checkbox"/> E8000208100 € 1,942.80
---	---	---

On SPEEDY 4all / 4you / 4all Ergo / 4you Ergo, SPEEDY 4teen
 Axle tube one-handed powder-coated

E8016200100-K € 132.80

Colour: _____

Seat width SB: _____ cm

Product number: _____

Handrim A

Attachment: tight x = 11 mm wide x = 19 mm

Stainless steel handrim

for 20" wheel, Ø 452 mm
 A9000205010-101
 for 22" wheel, Ø 488 mm
 A9000205010-102
 for 24" and 25" wheel Ø 539 mm
 A9000205010-104
 for 26" wheel, Ø 590 mm
 A9000205010-106

Aluminium handrim powder coated in frame colour (basic colour with two-colour design)

= Gloss surface
 = Matt surface
 Note: Matt coating is not recommended due to the low grip.

for 20" wheel, Ø 452 mm
 A9000205020-101 € 31.60
 for 22" wheel, Ø 488 mm
 A9000205020-102 € 31.60
 for 24" and 25" wheel, Ø 538 mm
 A9000205020-104 € 31.60
 for 26" wheel, Ø 590 mm
 A9000205020-106 € 31.60

Aluminium handrim polished and hard-anodised

for 20" wheel, Ø 452 mm
 A9000205030-101 € 76.20
 for 22" wheel, Ø 488 mm
 A9000205030-102 € 76.20
 for 24" and 25" wheel, Ø 538 mm
 A9000205030-104 € 76.20
 for 26" wheel, Ø 590 mm
 A9000205030-106 € 76.20

Ultra-Grip handrim

for 20" wheel, Ø 452 mm
 A9000205070-101 € 186.10
 for 22" wheel, Ø 486 mm
 A9000205070-102 € 186.10
 for 24" and 25" wheel, Ø 548 mm
 A9000205070-104 € 186.10
 for 26" wheel, Ø 587 mm
 A9000205070-106 € 186.10

Handrim B with 3 spoke fastening

Please select handrim B min. 2" smaller than handrim A

Attachment: tight y = 28 mm wide y = 38 mm

Stainless steel handrim

18", Ø 398 mm
 A9000208010-100
 20", Ø 452 mm
 A9000208010-101
 22", Ø 488 mm
 A9000208010-102
 24", Ø 539 mm
 A9000208010-104

Aluminium handrim powder coated in frame colour (basic colour with two-colour design)

= Gloss surface
 = Matt surface
 Note: Matt coating is not recommended due to the low grip.

18", Ø 397 mm
 A9000208020-100 € 31.60
 20", Ø 452 mm
 A9000208020-101 € 31.60
 22", Ø 488 mm
 A9000208020-102 € 31.60
 24", Ø 538 mm
 A9000208020-104 € 31.60

Aluminium handrim polished and hard-anodised

18", Ø 397 mm
 A9000208030-100 € 76.20
 20", Ø 452 mm
 A9000208030-101 € 76.20
 22", Ø 488 mm
 A9000208030-102 € 76.20
 24", Ø 538 mm
 A9000208030-104 € 76.20

Ultra-Grip handrim

18", Ø 397 mm
 A9000208070-100 € 186.10
 20", Ø 452 mm
 A9000208070-101 € 186.10
 22", Ø 486 mm
 A9000208070-102 € 186.10
 24", Ø 548 mm
 A9000208070-104 € 186.10

Handrim C

Attachment: tight x = 11 mm wide x = 19 mm

Stainless steel handrim

- for 20" wheel, Ø 448 mm
 A9000205010-101
for 22" wheel, Ø 488 mm
 A9000205010-102
for 24" and 25" wheel, Ø 538 mm
 A9000205010-104
for 26" wheel, Ø 585 mm
 A9000205010-106

Aluminium handrim powder coated in frame colour (basic colour with two-colour design)

- = Gloss surface
 = Matt surface
Note: Matt coating is not recommended due to the low grip.
- for 20" wheel, Ø 448 mm
 A9000205020-101 € 31.60
for 22" wheel, Ø 488 mm
 A9000205020-102 € 31.60
for 24" and 25" wheel, Ø 538 mm
 A9000205020-104 € 31.60
for 26" wheel, Ø 585 mm
 A9000205020-106 € 31.60

Aluminium handrim polished and hard-anodised

- for 20" wheel, Ø 448 mm
 A9000205030-101 € 76.20
for 22" wheel, Ø 488 mm
 A9000205030-102 € 76.20
for 24" and 25" wheel, Ø 538 mm
 A9000205030-104 € 76.20
for 26" wheel, Ø 585 mm
 A9000205030-106 € 76.20

Ultra-Grip handrim

- for 20" wheel, Ø 448 mm
 A9000205070-101 € 186.10
for 22" wheel, Ø 488 mm
 A9000205070-102 € 186.10
for 24" and 25" wheel, Ø 538 mm
 A9000205070-104 € 186.10
for 26" wheel, Ø 585 mm
 A9000205070-106 € 186.10

without handrim

Other

Silicone handrim cover without structure

- for 20" handrim, Ø 448 mm
 A9000205100-001 € 107.50 / Item
for 22" handrim, Ø 488 mm
 A9000205100-002 € 107.50 / Item
for 24" handrim, Ø 538 mm
 A9000205100-004 € 107.50 / Item
for 26" handrim, Ø 585 mm
 A9000205100-006 € 107.50 / Item

Please select which of the two handrims the silicone handrim cover should be mounted on:

- Handrim A
 Handrim B
 Handrim C

Knee lever brake with one-hand operation via Bowden cable

- for TRAVELER
 A9000707020-K € 362.80
- for SPEEDY 4all / 4you / 4all Ergo /
4you Ergo
 A9000707020-K € 362.80
- for SPEEDY A1 / A2
 A9000707020-K € 362.80
- for SPEEDY 4teen
 A9000707020-K € 362.80
- for BUDDY classic / 4all
 A9000707020-K € 362.80
- Operation right side
 Operation left side

Knee lever brake for accompanying persons

Push handle operation via Bowden cable.

- for TRAVELER
 A9000701020 € 429.10
- for SPEEDY A1 / A2
 A9000701020 € 429.10
- for SPEEDY 4teen
 A9000701020 € 429.10
- for BUDDY classic / 4all
 A9000701020 € 429.10

Location, Date

Signature End customer/dealer

The signature field is intended for the consultant and insured person. This confirms that the intended equipment is intended for the insured person.
- 4 - 2023-02-16 / The prices are plus statutory value-added tax and shipping costs = Included in the basic price = Position selected

Order sheet one-hand drive Rev. 07 UK