



## SPEEDY F3

**Easy running rigid frame wheelchair  
 with individual frame geometry made to measure  
 permanently welded chassis**

Maximum user weight limit 120 kg

A9014100000 € 5,560.00

**Quotation date:** \_\_\_\_\_  
**Consultation was carried out by**  
**Staff at the medical supplies shop:** \_\_\_\_\_  
**Enquiry no.:** \_\_\_\_\_

**Order date:** \_\_\_\_\_  
 according to GTC PRO ACTIV Reha-Technik GmbH  
**Order no.:** \_\_\_\_\_

**Medical suppliers / company address:**

Cust. no.: \_\_\_\_\_  
 Name: \_\_\_\_\_  
 Street: \_\_\_\_\_  
 Postcode / city: \_\_\_\_\_  
 Phone: \_\_\_\_\_  
 Fax: \_\_\_\_\_  
 E-Mail: \_\_\_\_\_  
 Contact: \_\_\_\_\_

**Shipping address:**  same as company address  
 Other (see below)

Name: \_\_\_\_\_  
 Street: \_\_\_\_\_  
 Postcode / city: \_\_\_\_\_

**Customer name / commission:**


ID / Comm. No.: \_\_\_\_\_  
 Surname, forename: \_\_\_\_\_  
 Street: \_\_\_\_\_  
 Postcode / city: \_\_\_\_\_  
 Phone: \_\_\_\_\_  
 Fax: \_\_\_\_\_  
 E-Mail: \_\_\_\_\_  
 D.o.B.: \_\_\_\_\_  
 Paying organisation: \_\_\_\_\_

**Invoicing address:**  same as company address  
 Other (see below)

Cust. no.: \_\_\_\_\_  
 Name: \_\_\_\_\_  
 Street: \_\_\_\_\_  
 Postcode / city: \_\_\_\_\_  
 Phone: \_\_\_\_\_  
 Fax: \_\_\_\_\_  
 E-Mail: \_\_\_\_\_  
 Contact: \_\_\_\_\_

**Note: We kindly would like to receive a picture of the customer in the wheelchair (the face can be blurred).**

**Body size:**





Body height: \_\_\_\_\_ cm  
 Weight: \_\_\_\_\_ kg  
 Lower leg length: \_\_\_\_\_ cm  
 (anatomic dimension: hollow of the knee to the foot sole)  
 Indication: \_\_\_\_\_

This information is necessary for checking. If the information is missing, a check of the relevant items cannot be made!

## Frame geometry


Customised production possible on request and after examination.



A9014100010  
Steering wheel bearing blocks  
angle-adjustable



A9014100020  
Steering wheel mount welded



### Wheelchair dimensions:


<b>A:</b> SB Seat width	_____ cm	<b>J:</b> Knee angle	_____ °
<b>B:</b> ST Seat depth (Seat upholstery)	_____ cm	<b>J1:</b> Angle of the front Framework tube to vertical	_____ °
<b>C:</b> RH Back height	_____ cm	<b>K:</b> Position of the steering wheel bearing block	_____ cm
<b>D:</b> USL Lower leg length	_____ cm	<b>L:</b> Distance between steering wheel bearing blocks	_____ cm
<b>E:</b> SHh Seat height back	_____ cm	<b>L1:</b> Clearance between the steering wheels	_____ cm
<b>F:</b> SHv Seat height front	_____ cm	<b>M:</b> Lower leg width top	_____ cm
<b>G:</b> KP Tipping point (inner edge of rear tube to centre of drive wheel axle)	_____ cm	<b>N:</b> Lower leg width at the bottom	_____ cm
<b>H:</b> Frame length to the front frame bend	_____ cm	<b>Alternative for specifying the knee angle (J):</b>	
<b>I:</b> Start frame retraction (max. H - 6 cm)	_____ cm	<b>X:</b> Inner edge of back tube to front edge of footrest	_____ cm
		<b>X1:</b> Inner edge of back tube to front edge of footrest	_____ cm

The frame is fitted with a base support plate as standard.

### Equipment without base plate

A9015100001

ATTENTION: If the frame is designed without a base support plate, it is NOT possible to fit adaptations (anti-tipping support, transit wheels, wheelbase extension, etc.)!




### Seat kink

A9014100001

On request, the frame can be designed with a customised seat kink. (If a seat kink is selected without specifying dimensions, the following standard dimensions are used: X = 16 cm, Y = rear seat height)

X: \_\_\_\_\_ Y: \_\_\_\_\_



The dimensional tolerances according to DIN EN ISO 13920 apply to the design (dimensions in mm)

Accuracy level	up to 30	over 30 to 120	over 120 to 400	over 400 to 1000
C	+/- 1	+/- 3	+/- 4	+/- 6

## Frame colour - One-colour design

Information on current colour designs and example images can be found at: <https://www.proactiv-gmbh.com/colours/>

### Frame colour: Selection from colour palette

Frame colour will be specified in the order

A9000100170-000

Raw, blasted, not coated

A9000100170-020

Note: In combination with moisture e.g. hand perspiration, can affect dark aluminum discoloration.

### Powder colours

**white RAL 9016**

gloss A9000100170-023

matt A9000100170-123 € 43.40

**brilliant silver**

gloss A9000100170-011

**anthracite RAL 7016**

gloss A9000100170-001

matt A9000100170-101 € 43.40

**black RAL 9005**

gloss A9000100170-003

matt A9000100170-019 € 43.40

**pure orange RAL 2004**

gloss A9000100170-022

matt A9000100170-122 € 43.40

**pink magenta RAL 4010**

gloss A9000100170-037

matt A9000100170-137 € 43.40

**pastel blue RAL 5024**

gloss A9000100170-036

matt A9000100170-136 € 43.40

**lemon green**

gloss A9000100170-005

matt A9000100170-105 € 43.40

**yellow RAL 1023**

gloss A9000100170-046

matt A9000100170-146 € 43.40

**red RAL 3001**

gloss A9000100170-047

matt A9000100170-147 € 43.40

**lagoon blue**

gloss A9000100170-052

matt A9000100170-152 € 43.40

**athletic grey RAL 7046**

gloss A9000100170-050

matt A9000100170-150 € 43.40

**muddy grey RAL 7033**

gloss A9000100170-051

matt A9000100170-151 € 43.40

**peach RAL 3022**

gloss A9000100170-038

matt A9000100170-138 € 43.40

### Powder colours „dormant“

**deep black effect silver**

gloss A9000100170-007

**sky blue**

gloss A9000100170-017

**navy blue**

gloss A9000100170-045

matt A9000100170-145 € 43.40

**electric green**

gloss A9000100170-015

matt A9000100170-115 € 43.40

**violet**

gloss A9000100170-044

**fire red**

gloss A9000100170-018

matt A9000100170-118 € 43.40

### Frame colour - One-colour design special colours „liquid coating“

€ 349.60

**storm (ice blue)**

gloss A9000100170-500

matt A9000100170-501

**prairie (rose gold)**

gloss A9000100170-504

matt A9000100170-505

**berry**

gloss A9000100170-526

matt A9000100170-527

**blue ocean**

gloss A9000100170-522

matt A9000100170-523

**emerald**

gloss A9000100170-514

matt A9000100170-515

**chili**

gloss A9000100170-516

matt A9000100170-517

## Frame colour - Two-colour design

Note: In the two-colour design, both colours are the same gloss level as the design colour.

**Two-colour design with one base colour (on the rear frame) and one design colour (on the front frame)** € 268.50

Please select a colour below:

### Base colour (on the rear frame)

<b>white RAL 9016</b> <input type="checkbox"/> A9000100170-200	<b>brilliant silver</b> <input type="checkbox"/> A9000100170-201	<b>anthracite RAL 7016</b> <input type="checkbox"/> A9000100170-202
<b>black RAL 9005</b> <input type="checkbox"/> A9000100170-203	<b>pastel blue RAL 5024</b> <input type="checkbox"/> A9000100170-206	<b>lemon green</b> <input type="checkbox"/> A9000100170-207
<b>yellow RAL 1023</b> <input type="checkbox"/> A9000100170-208	<b>lagoon blue</b> <input type="checkbox"/> A9000100170-212	<b>athletic grey RAL 7046</b> <input type="checkbox"/> A9000100170-210
<b>muddy grey RAL 7033</b> <input type="checkbox"/> A9000100170-211		

### Design colour (on the front frame)

<b>black RAL 9005</b> <input type="checkbox"/> gloss A9000100170-259 <input type="checkbox"/> matt A9000100170-260	<b>lemon green</b> <input type="checkbox"/> gloss A9000100170-250 <input type="checkbox"/> matt A9000100170-251	<b>chili</b> <input type="checkbox"/> gloss A9000100170-253 <input type="checkbox"/> matt A9000100170-254
<b>pure orange RAL 2004</b> <input type="checkbox"/> gloss A9000100170-255 <input type="checkbox"/> matt A9000100170-256	<b>emerald</b> <input type="checkbox"/> gloss A9000100170-274 <input type="checkbox"/> matt A9000100170-275	<b>storm (ice blue)</b> <input type="checkbox"/> gloss A9000100170-264 <input type="checkbox"/> matt A9000100170-265
<b>prairie (rose gold)</b> <input type="checkbox"/> gloss A9000100170-266 <input type="checkbox"/> matt A9000100170-267	<b>berry</b> <input type="checkbox"/> gloss A9000100170-284 <input type="checkbox"/> matt A9000100170-285	<b>blue ocean</b> <input type="checkbox"/> gloss A9000100170-280 <input type="checkbox"/> matt A9000100170-281

## Colour - Attachment parts

Note: If not specified, the caster fork and bearing block will be mounted in black.

**Drive wheel: aluminum rim, hub and spokes**

Black A9000100171-003

**Caster fork and bearing block**

Silver A9000100173-011  
 Black A9000100173-003

## Colour design of attachment parts according to customer request:

**Drive wheel: aluminum rim and design inlet<sup>1</sup> according to colour selection / hub and spokes black**

frame colour A9000100171-999 € 249.60

Selection from one-colour design (different from the frame colour) A9000100171-998 € 288.90

Colour: \_\_\_\_\_  
 A9000100170-

**Caster wheel: 5" aluminum rim<sup>2</sup> (6-spoke design)**

frame colour A9000100174-999 € 75.20

Selection from one-colour design (different from the frame colour) A9000100174-998 € 87.10

Colour: \_\_\_\_\_  
 A9000100170-

## Special colour request (Frame colour - One-colour design)

**Frame colour for single purchasing Selection from RAL colour palette**

A9000100170-100 € 414.40

**gloss level**

gloss  
 matt

**colour designation**

\_\_\_\_\_

**RAL number**

\_\_\_\_\_

**Note:** Attachments **not** available in design colour options

<sup>1</sup> Only in combination with Liteline hub

<sup>2</sup> Only in combination with caster wheel 5" aluminum rim 6-spoke-design good shock absorption A900050121\*, as well as with suspension fork 5-spoke-design A9000500160\*

## Folding backrest


A9000400021  
7-way angle adjustable,  
with integrated plastic-coated aluminum side  
parts



## Clothing guard for folding backrest

### Selecting the clothing protection size

A9000400001  
If not specified, a 2/3 wheel overlap will be set where possible.  
The number of notches determines the size of the clothing guard.



Size 1 (30 mm)

Size 2 (36 mm)

Size 3 (46 mm)

### Clothing guard removable via lock function

Cover can be adjusted to follow the wheel arch

#### Aluminium version

in anthracite matt

A9000403215-K € 243.20

in frame colour (base colour for two-colour  
design)

A9000403216-K € 311.10

in design colour (only possible with frame  
colour in two-colour design)

A9000403217-K € 345.10



#### Carbon version

A9000400013-K € 547.30



## Screwed backrest

Backrest not foldable.

### Aluminium side guard

A9015402470-002  
Back angle 90° to the floor.  
Note: Please note that the tyre may not be  
completely covered.



### Aluminium side guard, Customized contour

A9015402470-001 € 167.70


Measurement adjustment:

A = \_\_\_\_ cm / B = \_\_\_\_ cm

Adjustment to the upper edge of the tyre

Adjustment to the upper edge of the rim

If not specified, the side guard will be adjusted  
to the upper edge of the tyre.



### Aluminium side guard, with bolted clothing guard

A9015400025 € 361.20

Wheel cover runs parallel to the contour of the  
drive wheel tyres at a distance of approx.  
8 mm - 12 mm



### Side guard with integrated clothing guard

Individually manufactured from one piece, wheel cover runs parallel to the contour of the drive  
wheel tyres at a distance of approx. 8 mm - 12 mm

#### Aluminium version


A9015400026 € 656.80

in frame colour (base colour for two-colour  
design)



#### Carbon version

A9015400028 € 904.20



## Back tubes & push handles

### Equipment without push handles

A9000401000



### Back tube with integrated handles

A9000401051-K € 130.40  
straight.



### Aluminium push handles fixed in back tube

A9000401050\* € 130.40  
from RH 23 cm



### Horizontal push handles screwed into the back tube

A9000401076 € 130.40  
can be unscrewed



### Safety push handles with continuous height adjustment

A9000401060\* € 216.40  
from RH 30 cm with adaptable back  
from RH 35 cm with lordosis support  
from RH 43 cm with ergonomic backshell



### Safety push handles back-positioned length 45 cm

A9000401075\* € 372.90  
continuous height adjustment,  
Safety push handle holder mounted on the  
back cross bar  
not in combination with swing-away armrest



## Additional option

### Back tube with lordosis support

A9000400005-K € 114.40

- 1 cm
- 2 cm
- 3 cm



## Back equipment

If not specified, the „Adaptable back with back padding“ will be delivered.

### Adaptable back with back padding

A9000802010-K001  
Cover: Black, breathable, open surface  
Equipment with back cross bar



### Adaptable back with back padding „light“

A9000802320-K001  
Cover: Black, breathable, open surface  
Cushion thickness: 9 mm  
Equipment with back cross bar

A9000802320-K002  
Cover: Black, breathable, open surface  
Cushion thickness: 9 mm  
Equipment without back cross bar



### Adaptable back with back padding „Comfort“

A9000815101-K € 131.40  
Cover: Black, breathable, open surface  
Front padding: 20 mm  
Rear padding: 5 mm  
Equipment with back cross bar

### Equipment without back upholstery and back padding or back system


A9000802011 € - 53.00  
Equipment with back cross bar

## Seat upholstery

If not specified, the „Body Contour“ seat upholstery will be delivered.


### Body Contour

A9000801025-K  
approved up to a body weight of 120 kg



### Open belt system

A9000801010-K € 79.90




### Equipment without seat upholstery

A9000801050




## Seat cushion



**Seat cushion for individual frame geometry in rectangular contour**  
Inlet: Standard

A9000805190-K € 164.40



**Seat cushion for individual frame geometry with radius contour**  
Inlet: Standard

A9000805320-K € 206.10

**Seat cushion height:**

3 cm       5 cm  
 4 cm       6 cm

**Seat cushion for individual frame geometry in rectangular contour**  
Inlet: Anti-Decubitus

A9000805210-K € 484.40

**Seat cushion for individual frame geometry with radius contour**  
Inlet: Anti-Decubitus

A9000805310-K € 528.70

**Seat cushion height:**

4 cm       6 cm

### Cover:

Black, breathable, open surface

Black, breathable, closed surface

If not specified, the cover material „Black, breathable, open surface“ will be delivered.

**Seat cushion is provided**

## Wheel camber

\*Total width with tight-mounted handrim.

### 1° drive wheel camber

A9016202010  
Total width = approx. SB + 20 cm\*

### 2,5° drive wheel camber

A9016202051  
Total width = approx. SB + 23 cm\*

### 4° drive wheel camber

A9016202030  
Total width = approx. SB + 26 cm\*

### 6° drive wheel camber

A9016202040  
Total width = approx. SB + 30 cm\*

## E-Adaptations

**Preparation for the „e-motion“ / „e-motion DuoDrive“ mounting**

A9012201200 € 849.00

Equipment: Custom frame distances depending on frame geometry, e-motion torque arm, 1/2" axle, locking length 114 mm, 1/2" drive wheel bushing M27 x 1.5; length 52 mm

PRO ACTIV and Alber recommend mounting two anti-tipping supports. Please select under „Adaptations“ or order Alber tipping supports directly from Alber.

## Brake system

If not specified, the „Brake system push to lock with standard brake handle" will be delivered.

**Recommendation:** For KP 7 cm or lower, we recommend the „Brake system push to lock with brake handle 20° outwards“.

### Integral parking brake positioned on upper frame tube

A9015705000-K € 529.80

Please note:

In order to maintain the function of the foldable backrest, the plastic-coated side part version is mounted "-5", guided low by the contour, thus preventing collision with the integral brake when the backrest is folded down.



### Knee lever brake with standard brake lever

A9000701010-001



### Knee lever brake with long brake lever

A9000701010-002



### Knee lever brake with low brake lever

A9000701010-003



### Knee lever brake system Brake lever angled 20°

standard brake handle A9000700002

long brake handle A9000701010-004



### Brake system push to lock Brake handle to flip-down

A9000706008 € 158.60  
length: 9 cm



## Footrest

If not specified, the „continuous footrest" will be delivered.

Alle Fußbügelrohre inklusive verschweißstem Auflagebügel.

Bitte beachten: Bei Auswahl der „Fußbrettauflage gummiert" entfällt der Auflagebügel.



with support bracket



without support bracket

### Continuous footrest

A9016601010-K\*  
height adjustable +/- 3 cm



### Footrest folds up to the rear with spring locking mechanism

A9016602005-K\* € 287.40

from USL approx. 36 cm, depending on the individual frame geometry



### Footrest folds up to the rear

A9016602010-K\* € 307.90  
with 20 mm operating lever for manual release



## Footrest accessories

### Rubber coated footplate

A9000610000-K € 125.90

bolted, angle-adjustable  
Depth: 14 cm



### Footplate side stop

A9000610100 € 59.70

adjustable  
only in combination with rubber coated footplate



## Wheelbase extension

### Wheelbase extension angle-adjustable and removable

A9015100005 € 746.80


with continuous axle tube,  
height adjustment via gear joint  
adjustable in 6° angular increments  
wheel position displacement approx. 14 cm  
adjustable track  
Attachment possibility must be checked by  
PRO ACTIV.  
Only in combination with base plate.  
Up to 120 kg load capacity.



### Wheelbase extension permanently welded

A9015100002 € 631.60

with second axle tube,  
standard distance between axle tube 150 mm,  
other distances according to drawing set.  
Track width independently adjustable.  
Configuration must be checked by PRO ACTIV.  
Only possible without base plate.



Please select the wheel camber below:

1° drive wheel camber

4° drive wheel camber

6° drive wheel camber

Note: The camber of the wheelbase extension cannot be smaller than in the everyday position. The maximum possible camber is 6°.

If not specified, the next higher camber angle than in the everyday position will be mounted.

## Adaptations

### Anti-tipping support

right side  
A9000300100-K003 € 187.90

left side  
A9000300100-K004 € 187.90

swiveling,  
screwed to the welded base plate with a  
tension spring



### Transit wheels

right side  
A9000301100-001 € 164.80

left side  
A9000301100-002 € 164.80

with universal adapter, removable  
not in combination with anti-tipping support,  
right side or left side  
not in combination with tilting clamp, right side  
or left side



### Tilting clamp

right side  
A9000302101-001 € 140.10

left side  
A9000302101-002 € 140.10

with universal adapter and Quick pin  
not in combination with anti-tipping support,  
right side or left side  
not in combination with transit wheels, right  
side or left side



### Combination of adaptations with wheelbase extension

Also required if a wheelbase extension is  
desired in conjunction with anti-tipping support,  
transit wheels or tilting clamp.

A9015300000 € 840.70

Bracket for implementing compatibility:  
Wheelbase extension combined with anti-  
tippers, transit rollers or tilting brackets.  
Drive axle with radial gearing  
incl. anti-tipping support adapter.

### 4 Dance middle anti-tipping support

A9000300010 € 517.40

removable via quick release lever with rotating  
caster fork and skate wheel Ø 80 mm  
adjustable ground clearance and angle  
Colour: black



### 3" caster wheel

#### 3" skate wheel (80 mm)

A9000501110 € 70.80  
puncture-proof



### 4" caster wheel

#### 4" caster wheel „Standard“ solid rubber black

A9000501011-004 (82A), puncture-proof



#### 4" skate wheel (100 mm)

A9000501120 € 70.80  
puncture-proof



#### 4" safety light-up wheel (100 mm)

A9000501220 € 94.80  
puncture-proof



#### 4" CorMORun caster wheel (100 mm)

A9000501042 € 173.00  
with replaceable rubber 6 mm,  
puncture-proof  
Rim colour: Black  
Incl. hair guard spacer (5.3 mm)



#### 4" Frog Legs caster wheel (100 mm)

A9000501230 € 206.30  
with replaceable rubber 6 mm,  
puncture-proof  
Colour: black, silver rim  
Incl. hair guard spacer (5.3 mm)



### 5" caster wheel

#### 5" caster wheel „Standard“ PU black

A9000501036-005 (82A), puncture-proof



#### 5" caster wheel „Standard“ PU grey

A9000501031-005 (82A), puncture-proof



#### 5" safety light-up wheel (125 mm)

A9000501221 € 94.80  
puncture-proof



#### 5" aluminum rim 3-spoke design „good shock absorption“

silver anodised A9000501210 € 179.80

black A9000501211 € 179.80

with replaceable rubber 6 mm  
medium hard, puncture-proof  
Incl. hair guard spacer (3 mm)



#### 5" CorMORun caster wheel (125 mm)

A9000501041 € 181.30  
with replaceable rubber 6 mm,  
puncture-proof  
Rim colour: Black  
Incl. hair guard spacer (5.3 mm)



#### 5" Frog Legs caster wheel (125 mm)

A9000501232 € 206.30  
with replaceable rubber 6 mm,  
puncture-proof  
Colour: black, silver rim  
Incl. hair guard spacer (5.3 mm)



## 6" caster wheel

### 6" caster wheel „Standard“ PU grey

A9000501010-006  
(82A), puncture-proof



### 6" caster wheel „pneumatic tyre“

A9000501050 € 57.20



### 6" caster wheel „Flexel-System“

A9000501070 € 57.20  
puncture-proof



## Caster fork

For caster fork colours, see „Colour of add-on parts“

If not specified, the „Standard caster fork with screw axle“ will be delivered.

### Wish for a caster fork size

A9014500010

If not specified, a size will be selected depending on the frame geometry.

Please note that it is not possible to select a size for the short caster fork.

Size 1

Size 2

Size 3

### Caster fork Standard

A9000500100-K001  
with screwed axle

A9000500100-K002 € 167.40  
with quick-release axle



### Caster fork Heavy Duty for high loads

Extra weight approx. 100 - 150 gr.

A9000510111-K001 € 54.00  
with screwed axle

A9000510111-K002 € 218.90  
with quick-release axle



### Suspension fork

inkl. 5" aluminum rim „good shock absorption“

A9000500160-K002 € 616.90  
with quick-release axle



### Caster fork short

A9000500110-K001  
with screwed axle

Not in combination with 6" steering wheels

A9000500110-K002 € 167.40  
with quick-release axle

Not in combination with 5" or 6" caster wheels



### Caster fork Heavy duty short for high loads

A9000510120-K001 € 54.00  
with screwed axle

Not in combination with 6" steering wheels

A9000510120-K002 € 218.90  
with quick-release axle

Not in combination with 5" or 6" caster wheels



### Caster fork spare wheel „Standard“

A9000500300-K € 296.70  
Caster fork Standard with quick-release axles



### Caster fork spare wheel „Heavy Duty“

A9000500301-K € 337.90  
Caster fork Heavy Duty with quick-release axles



The spare caster fork does not have a caster wheel.  
In case a caster wheel is requested, please put a cross on a second selection in the category „Caster wheels“.

## Track distance

### Caster wheel track extension



A9000500400\*  
10 mm per side

€ 59.10

A9000500401\*  
15 mm per side

€ 59.10

A9000500410\*  
20 mm per side

€ 59.10

## Drive wheels

If not specified, the „drive wheel aluminum rim 36 spokes“ will be delivered.

### Drive wheel 18 spoke aluminum rim

24" (ETRTO 540)  
A9000203122  
with Liteline hub, radially spoked  
Not in combination with:  
- Puncture-proof tires  
- Handbike  
up to 70 kg body weight:

Sturz	1°	4°	6°	8°	10°
Gewicht					
bis 50 Kg	✓	✓	✓	✓	✓
51 - 60 Kg	✓	✓	✓	✓	
61 - 70 Kg	✓	✓			



### Drive wheel 36 spoke aluminum rim

24" (ETRTO 540)  
A9000203000-010  
with Liteline hub (left image),  
radially spoked  
up to 120 kg body weight

24" (ETRTO 540)  
A9000203000-020  
with aluminum front hub (right image),  
radially spoked  
over 120 kg body weight



### Drive wheel double rim 36-spokes

24" (ETRTO 540)  
A9000203100 € 157.80

25" (ETRTO 559)  
A9000203200 € 157.80

26" (ETRTO 590)  
A9000203300 € 157.80  
with aluminum front hub,  
radially spoked  
For use with high loads and body weight  
>120 kg



### Spinergy wheel „Spox Sport“ with 24 PBO spokes; rim, hub and spokes in black

24" (ETRTO 540)  
A9000203550 € 784.10

25" (ETRTO 559)  
A9000201559 € 784.10

26" (ETRTO 590)  
A9000203555 € 784.10



### Spinergy „Everyday Light Extreme“ with 12 PBO spokes; rim and hub in black, spokes in white

24" (ETRTO 540)  
A9000203554 € 957.20

Recommendation: For heavy loads  
(e.g. higher body weight), we recommend the  
Spinergy „Spox Sport“ wheel with 24 PBO  
spokes.



### Spinergy „Everyday Light Extreme“ with 12 PBO spokes; rim, hub and spokes in black


24" (ETRTO 540)  
A9000203552 € 957.20

Recommendation: For heavy loads  
(e.g. higher body weight), we recommend the  
Spinergy „Spox Sport“ wheel with 24 PBO  
spokes.



# Handrims

x = dimension handrim attachment



Handrim attachment general:

- x = 11 mm
- x = 19 mm

Handrim attachment Curve + Gekko handrim:

- x = 10 mm
- x = 18 mm

If not specified, the handrim is mounted wide (19 mm / 18 mm).

## Stainless steel Ø 20 mm

- 24" handrim for 24" wheel and 25" wheel A9000205010-004
- 26" handrim for 26" wheel A9000205010-006



## Aluminum Ø 19 mm silver anodised

- 24" handrim for 24" wheel and 25" wheel A9000205019-004



## Aluminum Ø 19 mm powder coated in frame colour

(basic colour with two-colour design)

- = Gloss surface
- = Matt surface

Note: Matt coating is not recommended due to the low grip.

- 24" handrim for 24" wheel and 25" wheel A9000205020-004 **€ 68.30**
- 26" handrim for 26" wheel A9000205020-006 **€ 68.30**



## Aluminum Ø 19 mm polished and hard-anodised

- 24" handrim for 24" wheel and 25" wheel **€ 164.90** A9000205030-004
- 26" handrim for 26" wheel **€ 164.90** A9000205030-006



## Titanium Ø 19 mm

- 24" handrim for 24" wheel and 25" wheel **€ 335.60** A9000205040-004
- 26" handrim for 26" wheel **€ 335.60** A9000205040-006



## Silicone handrim cover

- for handrim 24" A9000205100-004 **€ 232.80**
  - for handrim 26" A9000205100-006 **€ 232.80**
- Surface without structure, compatible with handrim Ø 19 mm or Ø 20 mm.



## Ultra-Grip „very grippy“

- 24" handrim for 24" wheel and 25" wheel **€ 402.40** A9000205070-004
- 26" handrim for 26" wheel **€ 402.40** A9000205070-006



## Ergo-Grip

- 24" handrim for 24" wheel and 25" wheel **€ 402.40** A9000205080-004
  - 26" handrim for 26" wheel **€ 402.40** A9000205080-006
- Attachment fitting to the rim

Note: Due to the complex geometry, concentricity deviations of +/- 2-3 mm are possible.



## Curve black anodised, ergonomic handrim profile

- 24" handrim for 24" wheel and 25" wheel Profile height 24 mm **€ 402.40** A9000205130-004
- 24" handrim for 24" wheel and 25" wheel Profile height 31 mm Size L **€ 402.40** A9000205140-004
- 26" handrim for 26" wheel Profile height 31 mm Size L **€ 402.40** A9000205140-006



<p><b>Curve silver anodised ergonomic handrim profile</b></p> <p><input type="checkbox"/> 24" handrim for 24" wheel and 25" wheel Profile height 31 mm <b>€ 402.40</b> Size L A9000205160-004</p> 	<p><b>Curve powder-coated, ergonomic handrim profile</b></p> <p><input type="checkbox"/> = frame colour (gloss surface) <input type="checkbox"/> = black (gloss surface)</p> <p><input type="checkbox"/> 24" handrim for 24" wheel and 25" wheel Profile height 24 mm <b>€ 467.60</b> A9000205135-004</p> <p><input type="checkbox"/> 24" handrim for 24" wheel and 25" wheel Profile height 31 mm <b>€ 467.60</b> Size L A9000205145-004</p> <p><input type="checkbox"/> 26" handrim for 26" wheel Profile height 31 mm <b>€ 467.60</b> Size L A9000205145-006</p> 	<p><b>Curve Grip, with Tetra-Grip coating, ergonomic handrim profile</b></p> <p><input type="checkbox"/> 24" handrim for 24" wheel and 25" wheel Profile height 31 mm <b>€ 520.70</b> Size L A9000205150-004</p> 
<p><b>Quadro with Tetra-Grip coating</b></p> <p><input type="checkbox"/> 24" handrim for 24" wheel and 25" wheel <b>€ 520.70</b> A9000205170-004</p> <p><input type="checkbox"/> 25" handrim for 25" wheel <b>€ 520.70</b> A9000205170-005</p> <p><input type="checkbox"/> 26" handrim for 26" wheel <b>€ 520.70</b> A9000205170-006</p> 	<p><b>Gekko black anodised</b></p> <p><input type="checkbox"/> 24" handrim for 24" wheel and 25" wheel <b>€ 520.70</b> A9000205175-004</p> <p><input type="checkbox"/> 25" handrim for 25" wheel <b>€ 520.70</b> A9000205175-005</p> <p><input type="checkbox"/> 26" handrim for 26" wheel <b>€ 520.70</b> A9000205175-006 with a robust and replaceable silicone lip, fully integrated into the profile Profile height 31 mm</p> 	<p><b>Gekko silver anodised</b></p> <p><input type="checkbox"/> 24" handrim for 24" wheel and 25" wheel <b>€ 520.70</b> A9000205185-004 with a robust and replaceable silicone lip, fully integrated into the profile Profile height 31 mm</p> 
<p><b>Gekko S black anodised</b></p> <p><input type="checkbox"/> 24" handrim for 24" wheel and 25" wheel <b>€ 520.70</b> A9000205195-004 with a robust and replaceable silicone lip, fully integrated into the profile Profile height 24 mm</p> 	<p><b>Equipment without handrims</b></p> <p><input checked="" type="checkbox"/> A9000205990</p>	

## Spoke guard

### Finger- and spoke guard „transparent“

for 24" wheel  
A9000206300-004 € 207.10

for 26" wheel  
A9000206300-006 € 207.10

not in combination with Spinerger drive wheels.



## Tyres

If not specified, the „Schwalbe RightRun in black“ tyres will be delivered.

### Schwalbe „RightRun“ High-pressure tyres

24x1" (540-25)  
A9000204001-004

Colour: Black  
7.5 bar



### Schwalbe „RightRun“ High-pressure tyres

24x1" (540-25)  
A9000204000-004

25x1" (559-25)  
A9000204000-005

26x1" (590-25)  
A9000204000-006

Colour: Grey, black wall  
7.5 bar



### Schwalbe „DownTown“ Profiled tyres

24x1" (540-25)  
A9000204010-004

7.5 bar



### Schwalbe „Marathon Plus Evolution“ with „smart guard“ puncture protection

24x1" (540-25)  
A9000204030-014 € 49.20

25x1" (559-25)  
A9000204030-015 € 49.20

26x1" (590-25)  
A9000204030-016 € 49.20

Colour: Black  
10 bar



### Schwalbe „Marathon Plus“ with „smart guard“ puncture protection

24x1" (540-25)  
A9000204030-004 € 49.20

Colour: Grey, black wall  
10 bar



### Schwalbe „One“

24x1" (540-25)  
A9000204123-004 € 93.00

Colour: Black  
9 bar



### Puncture-proof „solid rubber tyres“

24x1" (540-25)  
A9000204201-004 € 122.70

Colour: Black  
not in combination with Spinerger and 18 spoke  
drive wheels.



## Accessories

### 24" MTB spare wheel set

with double rim (ETRTO 507), aluminum front hub, quick-release axle and stainless steel handrim  
Mountain bike tyres, selectively:

A9000203410      € 703.40  
24x1.75" Schwalbe „Marathon Plus“  
Brake setting compatible with 24x1" wheel




A9000203413      € 703.40  
24x2.1" Schwalbe „SMART SAM Sport“  
Brake setting not compatible with 24x1" wheel



### Through axle with assistance, Large push button

A9000200096      € 105.80  
Interlock dimension: 122 / 110 mm  
Axle diameter: 12 mm



### Hand air compressor with pressure gauge

A9000901100      € 106.30  
up to 14 bar



### Crutch holder with flexible nylon cup

A9000902310\*      € 137.40



### Length-adjustable safety belt with buckle

A9000902100      € 80.60



### Adjustable safety belt with hook and fastener

A9000902110      € 64.50



### SIGMA lighting set

A9000903270      € 107.90

Consisting of front light and rear lights,  
chargeable via Micro USB socket,  
incl. Micro USB charging cable  
Lighting set compliant with road traffic  
regulations



### Transport bag for drive wheels

for 24" wheel  
A9000903530-004      € 77.80  
Ø 76 cm x 2 cm

for 26" wheel  
A9000903530-006      € 77.80  
Ø 83 cm x 2 cm



## Separate order sheets

### Category „Rehabilitation technology accessories“

e.g. Seat & backrest systems  
can be found at:

[www.proactiv-gmbh.com/downloads/](http://www.proactiv-gmbh.com/downloads/)

### Additional products

such as hand bikes or traction &  
attachment devices can be found at:

[www.proactiv-gmbh.com/downloads/](http://www.proactiv-gmbh.com/downloads/)

## Services

Consultation, clarification of the installation situation and measurement through PRO ACTIV

A9990000000 € 354.50

Handover and instruction on the wheelchair by PRO ACTIV

A9990000006 € 354.50

Notes and remarks:

### Location, Date

### Signature End customer/dealer

The signature field is intended for the consultant and insured person. This confirms that the intended equipment is intended for the insured person.  
- 17 - 2026-03-03 / The prices are plus statutory value-added tax and shipping costs  = included in the basic price  = Position selected