

Im Hofstätt 11  
D-72359 Dotternhausen  
Tel +49 (0)7427 9480-0  
Fax +49 (0)7427 9480-7025  
e-mail info@proactiv-gmbh.de  
www.proactiv-gmbh.de



## LIFT solid

### with electric seat height adjustment

With limited arm or finger function, the seat lifts upwards up to 30 cm with an electronic lifting function. The user reaches the desired work height and new action areas at the push of a button.

Maximum user weight limit 120 kg

A9061100000\* € 11,760.00

Quotation date: \_\_\_\_\_

Consultation was carried out by  
Staff at the medical supplies shop: \_\_\_\_\_

Enquiry no.: \_\_\_\_\_

Order date: \_\_\_\_\_

according to GTC PRO ACTIV Reha-Technik GmbH

Order no.: \_\_\_\_\_

### Medical suppliers / company address:

Cust. no.: \_\_\_\_\_  
Name: \_\_\_\_\_  
Street: \_\_\_\_\_  
Postcode / city: \_\_\_\_\_  
Phone: \_\_\_\_\_  
Fax: \_\_\_\_\_  
E-Mail: \_\_\_\_\_  
Contact: \_\_\_\_\_

Shipping address:  same as company address

Other (see below)

Name: \_\_\_\_\_  
Street: \_\_\_\_\_  
Postcode / city: \_\_\_\_\_

### Customer name / commission:

ID / Comm. No.: \_\_\_\_\_  
Surname, forename: \_\_\_\_\_  
Street: \_\_\_\_\_  
Postcode / city: \_\_\_\_\_  
Phone: \_\_\_\_\_  
Fax: \_\_\_\_\_  
E-Mail: \_\_\_\_\_  
D.o.B.: \_\_\_\_\_  
Paying organisation: \_\_\_\_\_

Invoicing address:  same as company address


Other (see below)

Cust. no.: \_\_\_\_\_  
Name: \_\_\_\_\_  
Street: \_\_\_\_\_  
Postcode / city: \_\_\_\_\_  
Phone: \_\_\_\_\_  
Fax: \_\_\_\_\_  
E-Mail: \_\_\_\_\_  
Contact: \_\_\_\_\_

# Basic size


## Body size:

Body height: \_\_\_\_\_ cm  
 Weight: \_\_\_\_\_ kg  
 Lower leg length: \_\_\_\_\_ cm  
 (anatomic dimension: hollow of the knee to the foot sole)  
 Indication: \_\_\_\_\_



This information is necessary for checking. If the information is missing, a check of the relevant items cannot be made!

## Wheelchair dimensions:



<b>Seat width SB</b>	<input type="checkbox"/> 35	<input type="checkbox"/> 36	<input type="checkbox"/> 37	<input type="checkbox"/> 38	<input type="checkbox"/> 39	<input type="checkbox"/> 40					
	<input type="checkbox"/> 41	<input type="checkbox"/> 42	<input type="checkbox"/> 43	<input type="checkbox"/> 44	<input type="checkbox"/> 45	<input type="checkbox"/> 46					
	The following dimensions are optionally available subject to surcharge:						<input type="checkbox"/> 30	<input type="checkbox"/> 31	<input type="checkbox"/> 32	<input type="checkbox"/> 33	<input type="checkbox"/> 34
<b>Seat depth ST</b>	Frame size 1		Frame size 2			Frame size 3					
	<input type="checkbox"/> 35	<input type="checkbox"/> 36	<input type="checkbox"/> 39	<input type="checkbox"/> 40	<input type="checkbox"/> 41	<input type="checkbox"/> 44	<input type="checkbox"/> 45	<input type="checkbox"/> 46			
	<input type="checkbox"/> 37	<input type="checkbox"/> 38	<input type="checkbox"/> 42	<input type="checkbox"/> 43		<input type="checkbox"/> 47	<input type="checkbox"/> 48				
	The following dimensions are optionally available subject to surcharge:						<input type="checkbox"/> 25	<input type="checkbox"/> 26	<input type="checkbox"/> 27	<input type="checkbox"/> 28	<input type="checkbox"/> 29
<b>Back height RH</b>	<input type="checkbox"/> 20	<input type="checkbox"/> 21	<input type="checkbox"/> 22	<input type="checkbox"/> 23	<input type="checkbox"/> 24	<input type="checkbox"/> 25	<input type="checkbox"/> 26	<input type="checkbox"/> 27	<input type="checkbox"/> 28	<input type="checkbox"/> 29	
	<input type="checkbox"/> 30	<input type="checkbox"/> 31	<input type="checkbox"/> 32	<input type="checkbox"/> 33	<input type="checkbox"/> 34	<input type="checkbox"/> 35	<input type="checkbox"/> 36	<input type="checkbox"/> 37	<input type="checkbox"/> 38	<input type="checkbox"/> 39	
	<input type="checkbox"/> 40	<input type="checkbox"/> 41	<input type="checkbox"/> 42	<input type="checkbox"/> 43	<input type="checkbox"/> 44	<input type="checkbox"/> 45	<input type="checkbox"/> 46	<input type="checkbox"/> 47	<input type="checkbox"/> 48		
<b>Lower leg length USL</b>	<input type="checkbox"/> 30	<input type="checkbox"/> 31	<input type="checkbox"/> 32	<input type="checkbox"/> 33	<input type="checkbox"/> 34	<input type="checkbox"/> 35	<input type="checkbox"/> 36	<input type="checkbox"/> 37	<input type="checkbox"/> 38	<input type="checkbox"/> 39	<input type="checkbox"/> 40
	<input type="checkbox"/> 41	<input type="checkbox"/> 42	<input type="checkbox"/> 43	<input type="checkbox"/> 44	<input type="checkbox"/> 45	<input type="checkbox"/> 46	<input type="checkbox"/> 47	<input type="checkbox"/> 48	<input type="checkbox"/> 49	<input type="checkbox"/> 50	
	The following dimensions are optionally available subject to surcharge:						<input type="checkbox"/> 20	<input type="checkbox"/> 21	<input type="checkbox"/> 22	<input type="checkbox"/> 23	<input type="checkbox"/> 24
<b>Seat height back SHh</b>	<input type="checkbox"/> 44	<input type="checkbox"/> 47									
<b>Seat height front SHv</b>	with 4" caster wheel:	<input type="checkbox"/> 46	<input type="checkbox"/> 48	<input type="checkbox"/> 50	<input type="checkbox"/> 52	<input type="checkbox"/> 52					
	with 5" caster wheel:	<input type="checkbox"/> 48	<input type="checkbox"/> 50	<input type="checkbox"/> 52	<input type="checkbox"/> 54	<input type="checkbox"/> 54					
	with 6" caster wheel:	<input type="checkbox"/> 50	<input type="checkbox"/> 52	<input type="checkbox"/> 54	<input type="checkbox"/> 56	<input type="checkbox"/> 56					

Please tick the required dimensions (all dimensions are in cm)

We recommend a seat height difference (difference between SHv and SHh) between 5 cm and 12 cm.

The SHh and SHv values might be a bit different from the specified values,

depending on the wheel size and camber.

Individual dimensions available on request.

## Frame colour - One-colour design

Information on current colour designs and example images can be found at: <https://www.proactiv-gmbh.com/colours/>

### Frame colour: Selection from colour palette

Frame colour will be specified in the order

A9000100170-000

Raw, blasted, not coated

A9000100170-020

Note: In combination with moisture e.g. hand perspiration, can affect dark aluminum discoloration.

### Powder colours

**white RAL 9016**

gloss A9000100170-023

matt A9000100170-123 € 43.40

**brilliant silver**

gloss A9000100170-011

**anthracite RAL 7016**

gloss A9000100170-001

matt A9000100170-101 € 43.40

**black RAL 9005**

gloss A9000100170-003

matt A9000100170-019 € 43.40

**pure orange RAL 2004**

gloss A9000100170-022

matt A9000100170-122 € 43.40

**pink magenta RAL 4010**

gloss A9000100170-037

matt A9000100170-137 € 43.40

**pastel blue RAL 5024**

gloss A9000100170-036

matt A9000100170-136 € 43.40

**lemon green**

gloss A9000100170-005

matt A9000100170-105 € 43.40

**yellow RAL 1023**

gloss A9000100170-046

matt A9000100170-146 € 43.40

**red RAL 3001**

gloss A9000100170-047

matt A9000100170-147 € 43.40

**lagoon blue**

gloss A9000100170-052

matt A9000100170-152 € 43.40

**athletic grey RAL 7046**

gloss A9000100170-050

matt A9000100170-150 € 43.40

**muddy grey RAL 7033**

gloss A9000100170-051

matt A9000100170-151 € 43.40

**peach RAL 3022**

gloss A9000100170-038

matt A9000100170-138 € 43.40

### Powder colours „dormant“

**deep black effect silver**

gloss A9000100170-007

**sky blue**

gloss A9000100170-017

**navy blue**

gloss A9000100170-045

matt A9000100170-145 € 43.40

**electric green**

gloss A9000100170-015

matt A9000100170-115 € 43.40

**violet**

gloss A9000100170-044

**fire red**

gloss A9000100170-018

matt A9000100170-118 € 43.40

### Frame colour - One-colour design special colours „liquid coating“

€ 349.60

**storm (ice blue)**

gloss A9000100170-500

matt A9000100170-501

**prairie (rose gold)**

gloss A9000100170-504

matt A9000100170-505

**berry**

gloss A9000100170-526

matt A9000100170-527

**blue ocean**

gloss A9000100170-522

matt A9000100170-523

**emerald**

gloss A9000100170-514

matt A9000100170-515

**chili**

gloss A9000100170-516

matt A9000100170-517

## Frame colour - Two-colour design

Note: In the two-colour design, both colours are the same gloss level as the design colour.

**Two-colour design with one base colour (on the rear frame) and one design colour (on the front frame)** € 268.50

Please select a colour below:

### Base colour (on the rear frame)

<b>white RAL 9016</b> <input type="checkbox"/> A9000100170-200	<b>brilliant silver</b> <input type="checkbox"/> A9000100170-201	<b>anthracite RAL 7016</b> <input type="checkbox"/> A9000100170-202
<b>black RAL 9005</b> <input type="checkbox"/> A9000100170-203	<b>pastel blue RAL 5024</b> <input type="checkbox"/> A9000100170-206	<b>lemon green</b> <input type="checkbox"/> A9000100170-207
<b>yellow RAL 1023</b> <input type="checkbox"/> A9000100170-208	<b>lagoon blue</b> <input type="checkbox"/> A9000100170-212	<b>athletic grey RAL 7046</b> <input type="checkbox"/> A9000100170-210
<b>muddy grey RAL 7033</b> <input type="checkbox"/> A9000100170-211		

### Design colour (on the front frame)

<b>black RAL 9005</b> <input type="checkbox"/> gloss A9000100170-259 <input type="checkbox"/> matt A9000100170-260	<b>lemon green</b> <input type="checkbox"/> gloss A9000100170-250 <input type="checkbox"/> matt A9000100170-251	<b>chili</b> <input type="checkbox"/> gloss A9000100170-253 <input type="checkbox"/> matt A9000100170-254
<b>pure orange RAL 2004</b> <input type="checkbox"/> gloss A9000100170-255 <input type="checkbox"/> matt A9000100170-256	<b>emerald</b> <input type="checkbox"/> gloss A9000100170-274 <input type="checkbox"/> matt A9000100170-275	<b>storm (ice blue)</b> <input type="checkbox"/> gloss A9000100170-264 <input type="checkbox"/> matt A9000100170-265
<b>prairie (rose gold)</b> <input type="checkbox"/> gloss A9000100170-266 <input type="checkbox"/> matt A9000100170-267	<b>berry</b> <input type="checkbox"/> gloss A9000100170-284 <input type="checkbox"/> matt A9000100170-285	<b>blue ocean</b> <input type="checkbox"/> gloss A9000100170-280 <input type="checkbox"/> matt A9000100170-281

## Colour - Attachment parts

Note: The attachments are always mounted in black.

**Drive wheel: aluminum rim, hub and spokes**

Black A9000100171-003

**Caster fork and bearing block**

Black A9000100173-003

### Colour design of attachment parts according to customer request:

**Drive wheel: aluminum rim and design inlet<sup>1</sup> according to colour selection / hub and spokes black**

frame colour A9000100171-999 € 249.60

Selection from one-colour design (different from the frame colour) A9000100171-998 € 288.90

Colour: \_\_\_\_\_

A9000100170-

**Caster wheel: 5" aluminum rim<sup>2</sup> (6-spoke design)**

frame colour A9000100174-999 € 75.20

Selection from one-colour design (different from the frame colour) A9000100174-998 € 87.10

Colour: \_\_\_\_\_

A9000100170-

## Special colour request (Frame colour - One-colour design)

**Frame colour for single purchasing Selection from RAL colour palette**

A9000100170-100 € 414.40

**gloss level**

gloss

matt

**colour designation**

\_\_\_\_\_

**RAL number**

\_\_\_\_\_

**Note:** Attachments **not** available in design colour options

<sup>1</sup> Only in combination with Liteline hub

<sup>2</sup> Only in combination with caster wheel 5" aluminum rim 6-spoke-design good shock absorption A900050121\*, as well as with suspension fork 5-spoke-design A9000500160\*

## Frame


### Straight front end

A9060100010-K\*  
Knee angle 105°



### Front end V-shape

A9060100020-K\*  
Knee angle 105°  
Note: Footrest area with V-shape  
 $x = SB - (0.25 \times USL + 6 \text{ cm})$



## Individual production

On request and after checking

### Seat width and Seat depth

A9060100100-K € 1,249.50

Seat width 30 cm - 34 cm  
and / or

Seat depth 25 cm - 34 cm

If necessary, with forward displacement of the backrest

### Seat width 47 cm - 50 cm

A9060100200-K € 1,230.10

### Lower leg leght 20 cm - 29 cm

A9060100300-K € 237.60

only for continuous footrest

## Back tubes & push handles

### Equipment without push handles

A9000401000



### Back tube with integrated handles

A9000401051-K € 130.40  
straight.



### Aluminium push handles fixed in back tube

A9000401050\* € 130.40

from RH 23 cm



### Horizontal push handles screwed into the back tube

A9000401076 € 130.40

can be unscrewed



### Safety push handles with continuous height adjustment

A9000401060\* € 216.40

mount in the back tube  
from RH 30 cm with adaptable back  
from RH 35 cm with lordosis support  
from RH 43 cm with ergonomic backshell



### Safety push handles back-positioned length 45 cm

A9000401075\* € 372.90

continuous height adjustment,  
Safety push handle holder mounted on the back cross bar  
not in combination with flip-down armrest



## Additional option

### Back tube with lordosis support

A9000400005-K € 114.40

1 cm

2 cm

3 cm



## Backrest

### Folding backrest

A9060400011  
7-way angle adjustable,  
with integrated plastic-coated aluminum side parts



## Back equipment

If not specified, the „Adaptable back with back padding“ will be delivered.

### Adaptable back with back padding

A9000802010-K001  
Cover: Black, breathable, open surface  
Equipment with back cross bar



### Adaptable back with back padding „light“

A9000802320-K001  
Cover: Black, breathable, open surface  
Cushion thickness: 9 mm  
Equipment with back cross bar

A9000802320-K002  
Cover: Black, breathable, open surface  
Cushion thickness: 9 mm  
Equipment without back cross bar



### Adaptable back with back padding „Comfort“

A9000815101-K **€ 131.40**  
Cover: Black, breathable, open surface  
Front padding: 20 mm  
Rear padding: 5 mm  
Equipment with back cross bar

### Back rise + 3 - 5 cm

A9000802100 **€ 159.10**  
for the nominal back height dimension  
Only in combination with height-adjustable  
push handles



### Equipment without back upholstery and back padding or back system


A9000802011 **€ - 53.00**  
Equipment with back cross bar

## Seat upholstery

If not specified, the „Body Contour“ seat upholstery will be delivered.


### Body Contour

A9000801025-K



### Open belt system

A9000801010-K **€ 79.90**





### Equipment without seat upholstery

A9000801050



## Seat cushion



### Seat cushion in rectangular contour Inlet: Standard

A9000805180-K\* **€ 118.30**

### Seat cushion with radius contour Inlet: Standard

A9000800300-K\* **€ 163.10**

### Seat cushion height:

3 cm  5 cm  
 4 cm  6 cm

### Seat cushion in rectangular contour Inlet: Anti-Decubitus

up to 30x30 cm  
A9000805204-K\* **€ 301.70**

>30x30 up to 50x50 cm  
A9000805200-K\* **€ 437.70**

### Seat cushion with radius contour Inlet: Anti-Decubitus

up to 30x30 cm  
A9000800311-K\* **€ 356.10**

>30x30 up to 50x50 cm  
A9000800310-K\* **€ 482.90**

### Seat cushion height:

4 cm  6 cm

### Cover:

- Black, breathable, open surface  
 Black, breathable, closed surface

If not specified, the cover material „Black, breathable, open surface“ will be delivered.

**Seat cushion is provided**

## Lift mechanism

### Electric seat height adjustment by 30 cm

Release on the right secured to the frame  
A9061300120-001

Release on the left secured to the frame  
A9061300120-002

Note: On request and after testing, the release is also possible on both sides (right and left)

### Battery charging socket moved upwards

A9061300020 € 253.90

positioned under the rear seat  
With toggle switch for switching between charging the battery and operating the lift function

## Clothing guard

### Plug in clothing guard

A9061403110 € 250.20

## Armrests

### Arm padding 25 cm length mounted on clothing guard

A9061403150 € 368.20

only in combination with plug in clothing guard

### Flip-down armrest

A9061403140 € 346.40

with arm padding 25 cm long  
not in combination with plug in clothing guard

## Axle plate

Standard equipment with 12 mm sockets/axles

### Drive wheel plate

A9061201100  
Rear seat height approx. 47 cm  
(for 24" wheel)

### Axle plate lowered


A9061201120-K € 260.90

Rear seat height approx. 44 cm  
(for 24" wheel)

### Extended drive wheel plate

A9060201110-K € 292.00

Wheelbase extension by 7.5 cm



## Wheel camber

\*Total width with tight-mounted handrim.

### 1° drive wheel camber

A9060202000  
Total width = approx. SB + 20 cm\*

## E-Adaptations

### Preparation for the „e-motion" mounting

A9060201200-K € 502.90

Equipment: e-motion torque arm, 1/2" axle and socket length 50.8 mm (2")  
Not allowed in combination with „e-motion DuoDrive".

### Preparation for the „e-fix" mounting

E35/E36 (current models) € 502.90  
A9061201230-K

E25/E26 (predecessor models for reuse) € 502.90  
A9061201220-K

Equipment: 1/2" axle, socket length 50.8 mm (2")  
Torque support e-fix and battery holder.

The e-fix attachment must be carried out by Fa. Alber, as an additional safety device is required in accordance with Alber installation guidelines.

## Brake system

If not specified, the „Knee lever brake system mounted at the rear” will be delivered.

### Knee lever brake system mounted at the rear

A9061701010



### Knee lever brake system mounted low down

A9061701015 € 156.10

Brake holding rail mounted on the axle plate, straight brake lever version



### Dual knee lever brake system

A9061701020 € 373.20

One brake mounted at the rear and another brake mounted at the front on the wheel plate with long brake lever in straight design

## Footrest

If not specified, the „continuous footrest” will be delivered.

### Continuous footrest

A9016601010-K\* with tubular resting bracket height adjustable +/- 3 cm If you want to equip the footrest with a „rubber coated footplate” please select this under „Footrest accessories”.



### Footrest folds up to the rear

A9016602010-K\* € 307.90

with tubular resting bracket and 20 mm operating lever from frame size 2 from USL 42 cm V-shape from SB 38 cm If you want to equip the footrest with a „rubber coated footplate” please select this under „Footrest accessories”.



### Divided footrest

A9016603200-K € 336.80

with rubber coated footplate, height adjustable +/- 3 cm, angle adjustment, can be pivoted outwards and folded up to the side Including frame reinforcement V-shape from SB 40 cm



### „Swing Away“ footrest

A9016604100-K € 404.70

with rubber coated footplate, separated in the middle, can be folded to the side, removable, can be folded up, height-adjustable, angle-adjustable Including frame reinforcement from USL 34 cm



### Removable front frame, foldable with divided footrest

A9065630000-K € 845.70

height-adjustable, with rubber plate Only with straight front end From 100 kg user weight only on request and after checking



## Footrest accessories

### Rubber coated footplate

A9000610000-K € 125.90

bolted, angle-adjustable Depth: 14 cm



### Footplate side stop

A9000610100 € 59.70

adjustable only in combination with rubber coated footplate



## Adaptations

### Anti-tipping support, removable

A9061300100  
adjustable with Quickpin



### 3" caster wheel

#### 3" skate wheel (80 mm)

A9000501110      € 70.80  
puncture-proof



### 4" caster wheel

#### 4" caster wheel „Standard“ solid rubber black

A9000501011-004  
(82A), puncture-proof



#### 4" skate wheel (100 mm)

A9000501120      € 70.80  
puncture-proof



#### 4" safety light-up wheel (100 mm)

A9000501220      € 94.80  
puncture-proof



#### 4" CorMORun caster wheel (100 mm)

A9000501042      € 173.00  
with replaceable rubber 6 mm,  
puncture-proof  
Rim colour: Black  
Incl. hair guard spacer (5.3 mm)



#### 4" Frog Legs caster wheel (100 mm)

A9000501230      € 206.30  
with replaceable rubber 6 mm,  
puncture-proof  
Colour: black, silver rim  
Incl. hair guard spacer (5.3 mm)



### 5" caster wheel

#### 5" caster wheel „Standard“ PU black

A9000501036-005  
(82A), puncture-proof



#### 5" caster wheel „Standard“ PU grey

A9000501031-005  
(82A), puncture-proof



#### 5" safety light-up wheel (125 mm)

A9000501221      € 94.80  
puncture-proof



<p><b>5" aluminum rim 3-spoke design „good shock absorption“</b></p> <p><input type="checkbox"/> silver anodised A9000501210 € 179.80</p> <p><input type="checkbox"/> black A9000501211 € 179.80</p> <p>with replaceable rubber 6 mm medium hard, puncture-proof Incl. hair guard spacer (3 mm)</p> 	<p><b>5" CorMORun caster wheel (125 mm)</b></p> <p><input type="checkbox"/> A9000501041 € 181.30</p> <p>with replaceable rubber 6 mm, puncture-proof Rim colour: Black Incl. hair guard spacer (5.3 mm)</p> 	<p><b>5" Frog Legs caster wheel (125 mm)</b></p> <p><input type="checkbox"/> A9000501232 € 206.30</p> <p>with replaceable rubber 6 mm, puncture-proof Colour: black, silver rim Incl. hair guard spacer (5.3 mm)</p> 
--	---	--

6" castor wheels available on request after checking, only with V-shape front end.

### Caster fork

If not specified, the „Standard caster fork with screw axle“ will be delivered.

#### Caster fork Standard

- A9000500100-K001 with screwed axle
- A9000500100-K002 € 167.40 with quick-release axle



#### Caster fork Heavy Duty

- for high loads or for use with additional electric drives
- A9000510111-K001 € 54.00 with screwed axle
- A9000510111-K002 € 218.90 with quick-release axle



#### Caster fork spare wheel „Standard“

- A9000500300-K € 296.70
- Caster fork Standard with quick-release axles



#### Caster fork spare wheel „Heavy Duty“

- A9000500301-K € 337.90
- Caster fork Heavy Duty with quick-release axles



The spare caster fork does not have a caster wheel.  
In case a caster wheel is requested, please put a cross on a second selection in the category „Caster wheels“.

### Track distance

#### Caster wheel track extension

- A9060500400-K € 73.60
- 13 mm per side



## Drive wheels

If not specified, the „drive wheel aluminum rim 36 spokes“ will be delivered.

### Drive wheel 36 spoke aluminum rim

24" (ETRTO 540)  
A9000203000-010  
with Liteline hub (left image),  
radially spoked  
up to 120 kg body weight

24" (ETRTO 540)  
A9000203000-020  
with aluminum front hub (right image),  
radially spoked  
over 120 kg body weight



### Drive wheel double rim 36-spokes

24" (ETRTO 540)  
A9000203100 € 157.80

25" (ETRTO 559)  
A9000203200 € 157.80

26" (ETRTO 590)  
A9000203300 € 157.80  
with aluminum front hub,  
radially spoked  
For use with high loads and body weight  
>120 kg



### Spinergy wheel „Spox Sport“ with 24 PBO spokes; rim, hub and spokes in black

24" (ETRTO 540)  
A9000203550 € 784.10


25" (ETRTO 559)  
A9000201559 € 784.10

26" (ETRTO 590)  
A9000203555 € 784.10



## Handrims

x = dimension handrim attachment



Handrim attachment  
general:

- x = 11 mm
- x = 19 mm

Handrim attachment  
Curve + Gekko handrim:

- x = 10 mm
- x = 18 mm

If not specified, the handrim is mounted wide (19 mm / 18 mm).

### Stainless steel Ø 20 mm

24" handrim for  
24" wheel and 25" wheel  
A9000205010-004

26" handrim for  
26" wheel  
A9000205010-006



### Aluminum Ø 19 mm silver anodised

24" handrim for  
24" wheel and 25" wheel  
A9000205019-004



### Aluminum Ø 19 mm powder coated in frame colour

(basic colour with two-colour design)

- = Gloss surface
- = Matt surface

Note: Matt coating is not recommended due to the low grip.

24" handrim for  
24" wheel and 25" wheel € 68.30  
A9000205020-004

26" handrim for  
26" wheel € 68.30  
A9000205020-006



### Aluminum Ø 19 mm polished and hard-anodised

24" handrim for  
24" wheel and 25" wheel € 164.90  
A9000205030-004

26" handrim for  
26" wheel € 164.90  
A9000205030-006



### Titanium Ø 19 mm

24" handrim for  
24" wheel and 25" wheel € 335.60  
A9000205040-004

26" handrim for  
26" wheel € 335.60  
A9000205040-006



### Silicone handrim cover

for handrim 24"  
A9000205100-004 € 232.80

for handrim 26"  
A9000205100-006 € 232.80

Surface without structure,  
compatible with handrim Ø 19 mm or  
Ø 20 mm.



<p><b>Ultra-Grip „very grippy“</b></p> <p><input type="checkbox"/> 24" handrim for 24" wheel and 25" wheel <b>€ 402.40</b> A9000205070-004</p> <p><input type="checkbox"/> 26" handrim for 26" wheel <b>€ 402.40</b> A9000205070-006</p> 	<p><b>Ergo-Grip</b></p> <p><input type="checkbox"/> 24" handrim for 24" wheel and 25" wheel <b>€ 402.40</b> A9000205080-004</p> <p><input type="checkbox"/> 26" handrim for 26" wheel <b>€ 402.40</b> A9000205080-006</p> <p>Attachment fitting to the rim</p> <p>Note: Due to the complex geometry, concentricity deviations of +/- 2-3 mm are possible.</p> 	<p><b>Curve black anodised, ergonomic handrim profile</b></p> <p><input type="checkbox"/> 24" handrim for 24" wheel and 25" wheel Profile height 24 mm <b>€ 402.40</b> A9000205130-004</p> <p><input type="checkbox"/> 24" handrim for 24" wheel and 25" wheel Profile height 31 mm <b>€ 402.40</b> Size L A9000205140-004</p> <p><input type="checkbox"/> 26" handrim for 26" wheel Profile height 31 mm <b>€ 402.40</b> Size L A9000205140-006</p> 
<p><b>Curve silver anodised ergonomic handrim profile</b></p> <p><input type="checkbox"/> 24" handrim for 24" wheel and 25" wheel Profile height 31 mm <b>€ 402.40</b> Size L A9000205160-004</p> 	<p><b>Curve powder-coated, ergonomic handrim profile</b></p> <p><input type="checkbox"/> = frame colour (gloss surface)</p> <p><input type="checkbox"/> = black (gloss surface)</p> <p><input type="checkbox"/> 24" handrim for 24" wheel and 25" wheel Profile height 24 mm <b>€ 467.60</b> A9000205135-004</p> <p><input type="checkbox"/> 24" handrim for 24" wheel and 25" wheel Profile height 31 mm <b>€ 467.60</b> Size L A9000205145-004</p> <p><input type="checkbox"/> 26" handrim for 26" wheel Profile height 31 mm <b>€ 467.60</b> Size L A9000205145-006</p> 	<p><b>Curve Grip, with Tetra-Grip coating, ergonomic handrim profile</b></p> <p><input type="checkbox"/> 24" handrim for 24" wheel and 25" wheel Profile height 31 mm <b>€ 520.70</b> Size L A9000205150-004</p> 
<p><b>Quadro with Tetra-Grip coating</b></p> <p><input type="checkbox"/> 24" handrim for 24" wheel and 25" wheel <b>€ 520.70</b> A9000205170-004</p> <p><input type="checkbox"/> 25" handrim for 25" wheel <b>€ 520.70</b> A9000205170-005</p> <p><input type="checkbox"/> 26" handrim for 26" wheel <b>€ 520.70</b> A9000205170-006</p> 	<p><b>Gekko black anodised</b></p> <p><input type="checkbox"/> 24" handrim for 24" wheel and 25" wheel <b>€ 520.70</b> A9000205175-004</p> <p><input type="checkbox"/> 25" handrim for 25" wheel <b>€ 520.70</b> A9000205175-005</p> <p><input type="checkbox"/> 26" handrim for 26" wheel <b>€ 520.70</b> A9000205175-006</p> <p>with a robust and replaceable silicone lip, fully integrated into the profile Profile height 31 mm</p> 	<p><b>Gekko silver anodised</b></p> <p><input type="checkbox"/> 24" handrim for 24" wheel and 25" wheel <b>€ 520.70</b> A9000205185-004</p> <p>with a robust and replaceable silicone lip, fully integrated into the profile Profile height 31 mm</p> 

### Gekko S black anodised

- 24" handrim for 24" wheel and 25" wheel **€ 520.70**  
A9000205195-004

with a robust and replaceable silicone lip, fully integrated into the profile  
Profile height 24 mm



### Equipment without handrims

- A9000205990

## Spoke guard

### Finger- and spoke guard „transparent“

- for 24" wheel **€ 207.10**  
A9000206300-004

- for 26" wheel **€ 207.10**  
A9000206300-006

not in combination with Spinergy drive wheels.



## Tyres

If not specified, the „Schwalbe RightRun in black“ tyres will be delivered.

### Schwalbe „RightRun“ High-pressure tyres

- 24x1" (540-25)  
A9000204001-004

Colour: Black  
7.5 bar



### Schwalbe „RightRun“ High-pressure tyres

- 24x1" (540-25)  
A9000204000-004

- 25x1" (559-25)  
A9000204000-005

- 26x1" (590-25)  
A9000204000-006

Colour: Grey, black wall  
7.5 bar



### Schwalbe „DownTown“ Profiled tyres

- 24x1" (540-25)  
A9000204010-004

7.5 bar



### Schwalbe „Marathon Plus Evolution“ with „smart guard“ puncture protection

- 24x1" (540-25) **€ 49.20**  
A9000204030-014

- 25x1" (559-25) **€ 49.20**  
A9000204030-015

- 26x1" (590-25) **€ 49.20**  
A9000204030-016

Colour: Black  
10 bar



### Schwalbe „Marathon Plus“ with „smart guard“ puncture protection

- 24x1" (540-25) **€ 49.20**  
A9000204030-004

Colour: Grey, black wall  
10 bar



### Schwalbe „One“


- 24x1" (540-25) **€ 93.00**  
A9000204123-004

Colour: Black  
9 bar



<p><b>Puncture-proof „solid rubber tyres”</b></p> <p><input type="checkbox"/> 24x1" (540-25) <b>€ 122.70</b> A9000204201-004 Colour: Black Not in combination with Spinergy drive wheels.</p> 	
--	--

## Accessories

<p><b>Through axle with assistance, Large push button</b></p> <p><input type="checkbox"/> A9000200096 <b>€ 105.80</b> Interlock dimension: 122 / 110 mm Axle diameter: 12 mm</p> 	<p><b>Hand air compressor with pressure gauge</b></p> <p><input type="checkbox"/> A9000901100 <b>€ 106.30</b> up to 14 bar</p> 	<p><b>Crutch holder with flexible nylon cup</b></p> <p><input type="checkbox"/> A9000902310* <b>€ 137.40</b></p> 
<p><b>Length-adjustable safety belt with buckle</b></p> <p><input type="checkbox"/> A9000902100 <b>€ 80.60</b></p> 	<p><b>Adjustable safety belt with hook and fastener</b></p> <p><input type="checkbox"/> A9000902110 <b>€ 64.50</b></p> 	<p><b>SIGMA lighting set</b></p> <p><input type="checkbox"/> A9000903270 <b>€ 107.90</b> Consisting of front light and rear lights, chargeable via Micro USB socket, incl. Micro USB charging cable Lighting set compliant with road traffic regulations</p> 
<p><b>Transport bag for drive wheels</b></p> <p><input type="checkbox"/> for 24" wheel A9000903530-004 <b>€ 77.80</b> Ø 76 cm x 2 cm</p> <p><input type="checkbox"/> for 26" wheel A9000903530-006 <b>€ 77.80</b> Ø 83 cm x 2 cm</p> 	<p><b>Stand-up aid handles via locking pin</b></p> <p><input type="checkbox"/> A9060403200 <b>€ 470.30</b> removable, height-adjustable mounted on the seat tube depending on the seat width, a collision with the frame and/or brakes is possible and may restrict height adjustability.</p> 	

## Separate order sheets

**Category „Rehabilitation technology accessories”**  
e.g. Seat & backrest systems  
can be found at:  
[www.proactiv-gmbh.com/downloads/](http://www.proactiv-gmbh.com/downloads/)

## Services

**Consultation, clarification of the installation situation and measurement through PRO ACTIV**

A9990000000 **€ 354.50**

**Handover and instruction on the wheelchair by PRO ACTIV**

A9990000006 **€ 354.50**

**Notes and remarks:**

---

**Location, Date**

**Signature End customer/dealer**

The signature field is intended for the consultant and insured person. This confirms that the intended equipment is intended for the insured person.

- 15 - 2026-03-03 / The prices are plus statutory value-added tax and shipping costs     = included in the basic price     = Position selected

Order sheet LIFT solid - electric Rev. 07 UK