



## **FREEWAY trekking steering attachment**

### **up to 120 kg load capacity**

The FREEWAY trekking is a steering attachment to facilitate handling of uneven surfaces and longer tours, even with longer downhill slopes. The adjustable steering column allows the handle bar position to be adjusted to a different angle based on changing requirements.

### **basic equipment**

- Quick-Snap docking plate for connecting to the frame adapter
- Handle bar with adjustable angle and foam handles and brake handle on the right side
- 220 mm steering column length
- Axle position can be adjusted using clamping fork attachment for easy height adjustment for adaptation to the (existing) Quick-Snap adapter
- Hydraulic disc brake
- Tyre Schwalbe Marathon 20"
- user manual

A9122000030

€ 1,020.00

**Quotation date:** \_\_\_\_\_

**Consultation was carried out by**  
**Staff at the medical supplies shop:** \_\_\_\_\_

**Enquiry no.:** \_\_\_\_\_

**Order date:** \_\_\_\_\_

according to GTC PRO ACTIV Reha-Technik GmbH

**Order no.:** \_\_\_\_\_

### **Medical suppliers / company address:**

**Cust. no.:** \_\_\_\_\_

**Name:** \_\_\_\_\_

**Street:** \_\_\_\_\_

**Postcode / city:** \_\_\_\_\_

**Phone:** \_\_\_\_\_

**Fax:** \_\_\_\_\_

**E-Mail:** \_\_\_\_\_

**Contact:** \_\_\_\_\_

**Different invoicing address (see below)**

**Cust. no.:** \_\_\_\_\_

**Name:** \_\_\_\_\_

**Street:** \_\_\_\_\_

**Postcode / city:** \_\_\_\_\_

**Different shipping address (see below)**

**Cust. no.:** \_\_\_\_\_

**Name:** \_\_\_\_\_

**Street:** \_\_\_\_\_

**Postcode / city:** \_\_\_\_\_

### **Customer name / commission:**

**ID / Comm. No.:** \_\_\_\_\_

**Surname, forename:** \_\_\_\_\_

**Street:** \_\_\_\_\_

**Postcode / city:** \_\_\_\_\_

**Phone:** \_\_\_\_\_

**E-Mail:** \_\_\_\_\_

**D.o.B.:** \_\_\_\_\_

**Paying organisation:** \_\_\_\_\_

**Indication:** \_\_\_\_\_

### **Body size**

Body height: \_\_\_\_\_ cm

Upper body height: \_\_\_\_\_ cm

Weight: \_\_\_\_\_ kg

This information is necessary for checking!

## Fork colour - One-Colour design

Information on current colour designs and example images can be found at: <https://www.proactiv-gmbh.com/colours/>

### Fork colour: Selection from colour palette

Frame colour will be specified in the order

A9000100170-000

#### Powder colours

##### white RAL 9016

gloss A9000100170-023

matt A9000100170-123 € 43.40

##### brilliant silver

gloss A9000100170-011

##### anthracite RAL 7016

gloss A9000100170-001

matt A9000100170-101 € 43.40

##### black RAL 9005

gloss A9000100170-003

matt A9000100170-019 € 43.40

##### pure orange RAL 2004

gloss A9000100170-022

matt A9000100170-122 € 43.40

##### pink magenta RAL 4010

gloss A9000100170-037

matt A9000100170-137 € 43.40

##### pastel blue RAL 5024

gloss A9000100170-036

matt A9000100170-136 € 43.40

##### lemon green

gloss A9000100170-005

matt A9000100170-105 € 43.40

##### yellow RAL 1023

gloss A9000100170-046

matt A9000100170-146 € 43.40

##### red RAL 3001

gloss A9000100170-047

matt A9000100170-147 € 43.40

##### lagoon blue

gloss A9000100170-052

matt A9000100170-152 € 43.40

##### athletic grey RAL 7046

gloss A9000100170-050

matt A9000100170-150 € 43.40

##### muddy grey RAL 7033

gloss A9000100170-051

matt A9000100170-151 € 43.40

##### peach RAL 3022

gloss A9000100170-038

matt A9000100170-138 € 43.40

#### Powder colours „dormant“

##### deep black effect silver

gloss A9000100170-007

##### sky blue

gloss A9000100170-017

##### navy blue

gloss A9000100170-045

matt A9000100170-145 € 43.40

##### electric green

gloss A9000100170-015

matt A9000100170-115 € 43.40

##### violet

gloss A9000100170-044

##### fire red

gloss A9000100170-018

matt A9000100170-118 € 43.40

## Fork colour – One-colour design special colours „liquid coating“

€ 349.60

##### storm (ice blue)

gloss A9000100170-500

matt A9000100170-501

##### prairie (rose gold)

gloss A9000100170-504

matt A9000100170-505

##### berry

gloss A9000100170-526

matt A9000100170-527

##### blue ocean

gloss A9000100170-522

matt A9000100170-523

##### emerald

gloss A9000100170-514

matt A9000100170-515

##### chili

gloss A9000100170-516

matt A9000100170-517

## Accessories

### Parking stand FREEWAY swiveling

A9120900001 € 205.90

for easy adaptation in combination with „Quick-Snap“ adapter  
Colour: black



### Parking disc brake

A9122300000 € 168.00

operating lever mounted on the handlebar support tube.  
Designed as a tension spreader.

### Plug-in mudguard

A9122502100 € 109.40



<p><b>Mudguard aluminium black</b></p> <p><input type="checkbox"/> A9122500005      € 153.50</p> 	<p><b>Aluminum luggage rack „cargo systems“</b></p> <p><input type="checkbox"/> A9120900000      € 210.00</p> <p>Colour: black Maximum load weight 10 kg</p> 	<p><b>SIGMA BC 8.0 ATS bike computer</b></p> <p><input type="checkbox"/> A9110901000      € 70.10</p> <p>wireless, with mount</p> 
<p><b>SIGMA lighting set</b></p> <p><input type="checkbox"/> A9120900019      € 123.90</p> <p>consisting of front light and rear lights, chargeable via Micro USB socket, incl. Micro USB charging cable. Lighting set compliant with road traffic regulations</p> 		
<b>Frame adapter</b>		
<p><b>SPEEDY rigid frame adapter left-hand operation with „Quick-Snap-Adaption“</b></p> <p><input type="checkbox"/> A9110700141      € 537.20</p> <p>with quick-release clamp and centering pin holder Colour: black</p> <p>Please select accordingly:</p> <p><b>Centring pin for Ø 25 mm</b> for A1 and A2</p> <p><input type="checkbox"/> A9110700022</p> <p><b>Centring pin for Ø 40 mm</b> for 4all, 4all Ergo, 4you Ergo, F2, F3 and 4teen</p> <p><input type="checkbox"/> A9110700023</p> 	<p><b>TRAVELER folding frame adapter with „Quick-Snap-Adaption“</b></p> <p><input type="checkbox"/> A9110701300      € 537.20</p> <p>with quick-release clamp and centering pin holder on the folding mechanism Colour: black</p> 	<p><b>BUDDY rigid frame adapter with „Quick-Snap-Adaption“</b></p> <p><input type="checkbox"/> A9110701400      € 537.20</p> <p>with quick-release clamp and centering pin mount on the frame Colour: black</p> 

**Wheelchair front adapter with universal clamp and „Quick-Snap-Adaption“**

A9110701500 € 572.30

Adapter attachment to the wheelchair stem. Operation on the right-hand side.



Mounting situation for PA products:  
Straight front end from SB 31 cm to 54 cm  
V-shaped front end from SB 30 cm to 50 cm

**Please select the frame clamp below:**

clamping range 95 mm  
A9110701525  
Tube diameter 23-27 mm

clamping range 95 mm  
A9110701530  
Tube diameter 28-32 mm

clamping range 95 mm  
A9110701535  
Tube diameter 33-37 mm

**PLEASE NOTE:**

If the exact tube diameter is not specified, an adapter for a tube diameter of 23 - 27 mm will be delivered as standard.

**Please select the clamping sleeve below:**

clamping sleeve [80 mm] included as standard  
for close clearance between footrest and drive wheel.

clamping sleeve [110 mm] € 46.70  
A9199000117  
for medium distance between footrest and drive wheel.

clamping sleeve [170 mm] € 56.60  
A9199000154  
for a wide distance between footrest and drive wheel.

**PLEASE NOTE:**

If not specified, the clamping sleeve with 80 mm length will be delivered.

**Combination adapter PA-Triride, steerable Quick-Snap version**

A9199000083 € 659.70

with Quick-Snap Adaptation for adapting PRO ACTIV handbikes and pulling devices as well as Triride products that are equipped with the PRO ACTIV adapter attachment package (A9140700000).

Frame clamps are NOT included in the delivery. Please select these separately.



**Please select the frame clamp below:**

Ø 25 mm  
E8140BS-C-106 € 205.20

Ø 30 mm  
E8140BS-C-107 € 205.20

Ø 32 mm  
E8140BS-C-108 € 205.20

**Options for users with reduced functions**

**Lever extension excenter lock for Quick Snap connection**

Makes it easier to open and close the Quick Snap adapter. Tetra version

A9110700005 € 73.90

straight design  
only in combination with rigid frame adapter

A9110700006 € 73.90

curved design  
only in combination with front adapter



**Adaptation help for centering pin**

A9110700002 € 55.00  
to simplify the insertion of the central adapter.



**Parking stand with operating handle for high comfort of use**

A9110900004 € 58.70  
with longer operating handle



**Chest strap for upper body fixation**

consisting of elastic webbing material in combination with neoprene with Velcro and operating tab.  
Colour: black

size S (105 cm)  
A9000800016 € 241.80

size M (125 cm)  
A9000800013 € 241.80



size L (145 cm)  
A9000800014 € 241.80

size XL (165 cm)  
A9000800017 € 241.80

## Services

### Consultation, clarification of the installation situation and measurement through PRO ACTIV

A9990000000 € 354.50

### Handover and instruction of the handbike / traction device by PRO ACTIV

A9990000011 € 354.50

### Attachment of the handbike / traction device to the wheelchair by PRO ACTIV

at specialist dealer

on a PRO ACTIV wheelchair

A9990000012 € 230.00

product no.: \_\_\_\_\_

on a third-party product

A9990000013 € 366.90

manufacturer: \_\_\_\_\_

model: \_\_\_\_\_

frame tube Ø \_\_\_\_ mm

### Attachment of the handbike / traction device to the wheelchair through specialist dealer

including mounting of the frame adapter

on a PRO ACTIV wheelchair

A9990000008

product no.: \_\_\_\_\_

on a third-party product

A9990000005

manufacturer: \_\_\_\_\_

model: \_\_\_\_\_

frame tube Ø \_\_\_\_ mm

Notes and remarks:

### Location, Date

The signature field is intended for the consultant and insured person. This confirms that the intended equipment is intended for the insured person

- 5 - 2026-03-03 / The prices are plus statutory value-added tax and shipping costs

= included in the basic price  = Position selected

Order sheet FREEWAY Rev. 08 UK

### Signature End customer/dealer