



Quotation date: _____
Consultation was carried out by
Staff at the medical supplies shop: _____
Enquiry no.: _____

Order date: _____
 according to GTC PRO ACTIV Reha-Technik GmbH
Order no.: _____

TRAVELER 4you Ergo Equipment „Edition PURE“ Extra-lightweight equipment with a total wheelchair weight starting at 10 kg with customized frame geometry

A9030105000 € 5,690.00

Medical suppliers / company address:
 Cust. no.: _____
 Name: _____
 Street: _____
 Postcode / city: _____
 Phone: _____
 Fax: _____
 E-Mail: _____
 Contact: _____

Invoicing address: same as company address
 Other (see below)


Cust. no.: _____
 Name: _____
 Street: _____
 Postcode / city: _____
 Phone: _____
 Fax: _____
 E-Mail: _____
 Contact: _____

Shipping address: same as company address
 Other (see below)

Name: _____
 Street: _____
 Postcode / city: _____

Customer name / commission:
 ID / Comm. No.: _____
 Surname, forename: _____
 Street: _____
 Postcode / city: _____
 Phone: _____
 Fax: _____
 E-Mail: _____
 D.o.B.: _____
 Paying organisation: _____


Body size:



Body height: _____ cm
 Weight: _____ kg
 Lower leg length: _____ cm
 (anatomic dimension: hollow of the knee to the foot sole)
 Indication: _____

This information is necessary for checking. If the information is missing, a check of the relevant items cannot be made!

Basic size



Wheelchair dimensions:

- | | | | |
|--|----------|--|----------|
| A: SB Seat width
(min. 33 cm – max. 52 cm) | _____ cm | H: Start frame retraction
(possible from ST 34 cm)
(I – min. 6 cm) | _____ cm |
| B: ST Seat depth
(min. 35 cm) | _____ cm | I: Frame length to the front frame bend
(min. 41 cm)
(min. ST - 1 cm) | _____ cm |
| C: RH Back height
(from RH 22 cm) | _____ cm | J: KW Knee angle
(from 90°) | _____ ° |
| D: USL Lower leg length | _____ cm | J1: Angle of the front Framework tube to vertical
(recommended > 10°) | _____ ° |
| E: SHh Seat height back
(min. 39 cm – max. 47 cm)
with 24" drive wheels | _____ cm | K: Position of the steering wheel bearing block | _____ cm |
| F: SHv Seat height front
(min. 42 cm – max. 52 cm)
with 5" caster wheels | _____ cm | Front end | |
| G: KP Tipping point
(inner edge of rear tube to centre of drive wheel axle)
(KP min. 3 cm – max. 11 cm) | _____ cm | <input type="checkbox"/> Straight (without retraction)
USL-width = SB – 5 cm | |
| Attachment wheelbase extension possible from KP 7 cm | | <input type="checkbox"/> V-shape (7°) | |
| | | <input type="checkbox"/> with frame retraction
parallel front end, straight shape
USL-width min. SB - 11 cm
Please specify the USL-width: | |
| | | M: Lower leg width | _____ cm |

The dimensional tolerances according to DIN EN ISO 13920 apply to the design (dimensions in mm)

Accuracy level	up to 30	over 30 to 120	over 120 to 400	over 400 to 1000
C	+/- 1	+/- 3	+/- 4	+/- 6

Frame colour - One-colour design

Information on current colour designs and example images can be found at: <https://www.proactiv-gmbh.com/colours/>

Frame colour: Selection from colour palette

Frame colour will be specified in the order

A9000100170-000

Raw, blasted, not coated

A9000100170-020

Note: In combination with moisture e.g. hand perspiration, can affect dark aluminum discoloration.

Powder colours

white RAL 9016

gloss A9000100170-023

matt A9000100170-123

€ 43.40

brilliant silver

gloss A9000100170-011

anthracite RAL 7016

gloss A9000100170-001

matt A9000100170-101

€ 43.40

black RAL 9005

gloss A9000100170-003

matt A9000100170-019

€ 43.40

pure orange RAL 2004

gloss A9000100170-022

matt A9000100170-122

€ 43.40

pink magenta RAL 4010

gloss A9000100170-037

matt A9000100170-137

€ 43.40

pastel blue RAL 5024

gloss A9000100170-036

matt A9000100170-136

€ 43.40

lemon green

gloss A9000100170-005

matt A9000100170-105

€ 43.40

yellow RAL 1023

gloss A9000100170-046

matt A9000100170-146

€ 43.40

red RAL 3001

gloss A9000100170-047

matt A9000100170-147

€ 43.40

lagoon blue

gloss A9000100170-052

matt A9000100170-152

€ 43.40

athletic grey RAL 7046

gloss A9000100170-050

matt A9000100170-150

€ 43.40

muddy grey RAL 7033

gloss A9000100170-051

matt A9000100170-151

€ 43.40

peach RAL 3022

gloss A9000100170-038

matt A9000100170-138

€ 43.40

Powder colours „dormant“

deep black effect silver

gloss A9000100170-007

sky blue

gloss A9000100170-017

navy blue

gloss A9000100170-045

matt A9000100170-145

€ 43.40

electric green

gloss A9000100170-015

matt A9000100170-115

€ 43.40

violet

gloss A9000100170-044

fire red

gloss A9000100170-018

matt A9000100170-118

€ 43.40

Frame colour - One-colour design special colours „liquid coating“

€ 349.60

storm (ice blue)

gloss A9000100170-500

matt A9000100170-501

prairie (rose gold)

gloss A9000100170-504

matt A9000100170-505

berry

gloss A9000100170-526

matt A9000100170-527

blue ocean

gloss A9000100170-522

matt A9000100170-523

emerald

gloss A9000100170-514

matt A9000100170-515

chili

gloss A9000100170-516

matt A9000100170-517

Frame colour - Two-colour design

Note: In the two-colour design, both colours are the same gloss level as the design colour.

Two-colour design with one base colour (on the rear frame) and one design colour (on the front frame) € 268.50

Please select a colour below:

Base colour (on the rear frame)

white RAL 9016 <input type="checkbox"/> A9000100170-200	brilliant silver <input type="checkbox"/> A9000100170-201	anthracite RAL 7016 <input type="checkbox"/> A9000100170-202
black RAL 9005 <input type="checkbox"/> A9000100170-203	pastel blue RAL 5024 <input type="checkbox"/> A9000100170-206	lemon green <input type="checkbox"/> A9000100170-207
yellow RAL 1023 <input type="checkbox"/> A9000100170-208	lagoon blue <input type="checkbox"/> A9000100170-212	athletic grey RAL 7046 <input type="checkbox"/> A9000100170-210
muddy grey RAL 7033 <input type="checkbox"/> A9000100170-211		

Design colour (on the front frame)

black RAL 9005 <input type="checkbox"/> gloss A9000100170-259 <input type="checkbox"/> matt A9000100170-260	lemon green <input type="checkbox"/> gloss A9000100170-250 <input type="checkbox"/> matt A9000100170-251	chili <input type="checkbox"/> gloss A9000100170-253 <input type="checkbox"/> matt A9000100170-254
pure orange RAL 2004 <input type="checkbox"/> gloss A9000100170-255 <input type="checkbox"/> matt A9000100170-256	emerald <input type="checkbox"/> gloss A9000100170-274 <input type="checkbox"/> matt A9000100170-275	storm (ice blue) <input type="checkbox"/> gloss A9000100170-264 <input type="checkbox"/> matt A9000100170-265
prairie (rose gold) <input type="checkbox"/> gloss A9000100170-266 <input type="checkbox"/> matt A9000100170-267	berry <input type="checkbox"/> gloss A9000100170-284 <input type="checkbox"/> matt A9000100170-285	blue ocean <input type="checkbox"/> gloss A9000100170-280 <input type="checkbox"/> matt A9000100170-281

Colour - Attachment parts

Note: If not specified, the caster fork and bearing block will be mounted in black.

Drive wheel: aluminum rim, hub and spokes

Black A9000100171-003

Caster fork and bearing block

Silver A9000100173-011
 Black A9000100173-003

Colour design of attachment parts according to customer request:

Drive wheel: aluminum rim and design inlet¹ according to colour selection / hub and spokes black

frame colour A9000100171-999 € 249.60

Selection from one-colour design (different from the frame colour) A9000100171-998 € 288.90

Colour: _____
A9000100170-

Caster wheel: 5" aluminum rim² (6-spoke design)

frame colour A9000100174-999 € 75.20

Selection from one-colour design (different from the frame colour) A9000100174-998 € 87.10

Colour: _____
A9000100170-

Special colour request (Frame colour - One-colour design)

Frame colour for single purchasing Selection from RAL colour palette

A9000100170-100 € 414.40

gloss level

gloss
 matt

colour designation

RAL number

Note: Attachments **not** available in design colour options

¹ Only in combination with Liteline hub

² Only in combination with caster wheel 5" aluminum rim 6-spoke-design good shock absorption A900050121*, as well as with suspension fork 5-spoke-design A9000500160*

- Equipment package

A9030105000-010

€ 1,503.30

with the „PURE“ equipment package, the product is approved for a maximum load capacity of 100 kg.

The package includes the following backrest equipment, clothing guard, seat upholstery, brake system and footrest for TRAVELER pure:

Side guard Carbon with integrated clothing guard

A9015400028

Clothing guard and side guard in one piece, individually manufactured. Backrest permanently bolted, cannot be folded, backrest angle 90° to the floor, adjustable in 5 positions in 4° steps, wheel cover runs contour-parallel to the drive wheel tyres at a distance of approx. 8 mm - 12 mm.



Adaptable back with back padding „light“

A9000802320-K002

Back padding with extra thin padding for easy, space-saving folding of the wheelchair. Cover: Black, breathable, open surface
Cushion thickness: 9 mm
Equipment without back cross bar



Open belt system

A9000801011-K



Integral parking brake

A9035705000

Brake lever fixed under the upper frame tube



Footrest without protrusion, rounded

A9000610015

Lightweight design with radius end (R 12 mm) in the front area and without protrusion to reduce the overall wheelchair length by 30 mm to improve vehicle loading. Depth: 15 cm



Push handles

Equipment without push handles

A9000401000



Back tube with integrated handles

A9000401051-K € 130.40
extra weight: ca. 80 g
straight.




Aluminium push handles fixed in back tube

A9000401050* € 130.40
extra weight: ca. 230 g
from RH 23 cm




Seat cushion



Seat cushion for individual frame geometry in rectangular contour Inlet: Standard

A9000805190-K € 164.40



Seat cushion for individual frame geometry with radius contour Inlet: Standard

A9000805320-K € 206.10

Seat cushion for individual frame geometry in rectangular contour Inlet: Anti-Decubitus

A9000805210-K € 484.40

Seat cushion for individual frame geometry with radius contour Inlet: Anti-Decubitus

A9000805310-K € 528.70

Seat cushion height:

3 cm 5 cm
 4 cm 6 cm

Seat cushion height:

4 cm 6 cm

Cover:

- Black, breathable, open surface
 Black, breathable, closed surface

If not specified, the cover material „Black, breathable, open surface“ will be delivered.

Seat cushion is provided

Wheel camber

*Total width with tight-mounted handrim.

1° drive wheel camber

A9030202010
Total width = approx. SB + 20 cm*

2,5° drive wheel camber

A9016202051
Total width = approx. SB + 23 cm*

4° drive wheel camber

A9030202030
Total width = approx. SB + 26 cm*

Footrest

If not specified, the „continuous footrest“ will be delivered.

Continuous footrest

A9030601000-K
with tubular resting bracket, foldable, removable
When folding the wheelchair, the footrest folds down completely from USL 37 cm.
For V-shape: If the footrest area x is smaller than 22 cm, the footrest will be parallel.



Footrest folds up to the rear with spring locking mechanism

A9035602005-K* € 353.20
extra weight: ca. 94 g
removable, with tubular resting bracket
for KW 105° from USL 39 cm
for KW 95° from USL 38 cm



Caster wheel & Caster fork

For caster fork colours, see „Colour of add-on parts“

5" aluminum rim 3-spoke design „good shock absorption“

silver anodised
A9000501210 € 179.80
weight: 252 g/Pcs.

black
A9000501211 € 179.80
weight: 252 g/Pcs.

with replaceable rubber 6 mm
medium hard, puncture-proof
Incl. hair guard spacer (3 mm)



CorMORun caster wheel

4" (100 mm)
A9000501042 € 173.00
weight: 200 g/Pcs.

5" (125 mm)
A9000501041 € 181.30
weight: 278 g/Pcs.

with replaceable rubber 6 mm,
puncture-proof
Rim colour: Black
Incl. hair guard spacer (5.3 mm)



Frog Legs caster wheel

4" (100 mm)
A9000501230 € 206.30
weight: 179 g/Pcs.

5" (125 mm)
A9000501232 € 206.30
weight: 278 g/Pcs.

with replaceable rubber 6 mm,
puncture-proof
Colour: black, silver rim
Incl. hair guard spacer (5.3 mm)



Caster fork Standard

A9000500100-K001
with screwed axle € 167.40

A9000500100-K002
with quick-release axle



Drive wheels

Drive wheel 18 spoke aluminum rim

24" (ETRTO 540)
A9000203123
weight: 715 g/Pcs.

with Liteline hub, radially spoked
Inlet in frame color

Not in combination with:

- Transport restraint connection
- Handbike

up to 70 kg body weight:

Sturz	1°	4°	6°	8°	10°
Gewicht					
bis 50 Kg	✓	✓	✓	✓	✓
51 - 60 Kg	✓	✓	✓	✓	
61 - 70 Kg	✓	✓			



Spinergy „Everyday Light Extreme“ with 12 PBO spokes; rim and hub in black, spokes in white

24" (ETRTO 540)
A9000203554 € 957.20
weight: 710 g/Pcs.

Recommendation: For heavy loads
(e.g. higher body weight), we recommend the
Spinergy „Spox Sport“ wheel with 24 PBO
spokes.



Spinergy „Everyday Light Extreme“ with 12 PBO spokes; rim, hub and spokes in black

24" (ETRTO 540)
A9000203552 € 957.20
weight: 710 g/Pcs.

Recommendation: For heavy loads
(e.g. higher body weight), we recommend the
Spinergy „Spox Sport“ wheel with 24 PBO
spokes.



Drive wheel 36 spoke aluminum rim

24" (ETRTO 540)
A9000203000-010
weight: 877 g/Pcs.

with Liteline hub




Spinergy wheel „Spox Sport“ with 24 PBO spokes; rim, hub and spokes in black

24" (ETRTO 540)
A9000203550 € 784.10
weight: 730 g/Pcs.



Handrims

x = dimension handrim attachment



Handrim attachment
general:

- x = 11 mm
 x = 19 mm

Handrim attachment
Curve + Gekko handrim:

- x = 10 mm
 x = 18 mm

If not specified, the handrim is mounted wide (19 mm / 18 mm).

Aluminum Ø 19 mm silver anodised



- 24" handrim for
24" wheel and 25" wheel
A9000205019-004
weight: 295 g/Pcs.



Titanium Ø 19 mm

- 24" handrim for
24" wheel and 25" wheel
A9000205040-004
weight: 418 g/Pcs. € 335.60



Curve black anodised, ergonomic handrim profile

- 24" handrim for
24" wheel and 25" wheel
Profile height 24 mm
Size S
A9000205130-004
weight: 330 g/Pcs. € 402.40



Tyres

Schwalbe „One“



- 24x1" (540-25)
A9000204123-004
weight: 205 g/Pcs.

Colour: Black
9 bar

€ 93.00



Schwalbe „Marathon Plus Evolution“ with „smart guard“ puncture protection

- 24x1" (540-25)
A9000204030-014
weight: 512 g/Pcs. € 49.20

Colour: Black
10 bar



Schwalbe „RightRun“ High-pressure tyres

- 24x1" (540-25)
A9000204001-004
weight: 294 g/Pcs.

Colour: Black
7.5 bar



Services

Consultation, clarification of the installation situation and measurement through PRO ACTIV

- A9990000000 € 354.50

Handover and instruction on the wheelchair by PRO ACTIV

- A9990000006 € 354.50

Notes and remarks:

(Note: If you would like to choose other options, please note the items with item number from the „TRAVELER 4you Ergo“ order sheet below)

Location, Date

Signature End customer/dealer

The signature field is intended for the consultant and insured person. This confirms that the intended equipment is intended for the insured person.

- 8 - 2026-03-03 / The prices are plus statutory value-added tax and shipping costs = included in the basic price = Position selected

Order sheet TRAVELER 4you Ergo Edition PURE Rev. 02 UK