



WHEEL-E

up to 150 kg trailer load and 10 kg payload
equipped with:

- Docking plate for quick connection to the „Quick-Snap adapter“
- Display for driving speed, assistance level and charge level display
- Driving speed up to 6 km/h with back gear, adjustable via drive signal transmitter on the right
- with hydraulic disc brake system and additional mechanical disc brake with parking function
- Parking stand with maneuvering rollers
- Handlebar, attached to the angle-adjustable stem
- Sigma LED battery lighting set with front and rear light
- Bell
- Drive wheel 20" (ETRTO 406), rim black
- Drive unit black

„The powerful, high-performance drive
for demanding conditions“

Climbing power up to 15%

A9130100031-020

€ 6,530.00

Standard motor

55 Nm torque / rated power 1000 W

Included in the base price

Driving comfort control with more powerful HT motor

70 Nm torque / rated power 1000 W

A9130100050 € 511.20

incl. downhill brake in cruise control function via active engine brake. One-hand operation possible. Speed up and brake with the thumb throttle.

Quotation date: _____

**Consultation was carried out by
Staff at the medical supplies shop:** _____

Enquiry no.: _____

Order date: _____

according to GTC PRO ACTIV Reha-Technik GmbH

Order no.: _____

Medical suppliers / company address:

Cust. no.: _____

Name: _____

Street: _____

Postcode / city: _____

Phone: _____

Fax: _____

E-Mail: _____

Contact: _____

Different invoicing address (see below)

Cust. no.: _____

Name: _____

Street: _____

Postcode / city: _____

Different shipping address (see below)

Cust. no.: _____

Name: _____

Street: _____

Postcode / city: _____

Customer name / commission:

ID / Comm. No.: _____

Surname, forename: _____

Street: _____

Postcode / city: _____

Phone: _____

E-Mail: _____

D.o.B.: _____

Paying organisation: _____

Indication: _____

Brake system

Parking disc brake, operating lever mounted on the handlebar support tube

A9130300106 € 164.30

designed as a pull-spreader.
Not in combination with double disc brake system.

Double disc brake system „Magura“ brake disc Ø 203 mm

A9130300105 € 203.00

recommended from 80 kg trailer load and steep inclines, braking and locking function via locking bracket



Battery and display

Battery 48V 10,5Ah 504Wh

A9130900016
Range approx. 25 km - 45 km
Battery mounted on the right, balance weight mounted on the left



Double Battery 48V 10.5Ah 504Wh

A9130900017 € 1,558.60

Range approx. 50 km - 90 km
positioned with plug-in connection via magnetic plug on the steering column
incl. 2nd charger



4.3" LCD colour display with USB connection

A9130200030 € 232.50

mounting on the right side of the handlebar clamp (replaces the 2.6" standard display)

Tyres

20x2.15" Schwalbe „Marathon 365“

A9117105103
(ETRTO 406-55)
profile tyre
puncture protection level 5



20x2.40" Schwalbe „Super Moto-X“

A9130105132 € 28.70
(ETRTO 406-62)
puncture protection level 6



20x4.0" tyres BFT

A9130105130 € 44.00

(ETRTO 406-100)
Profile tyres, Fat Tire




Handlebar

Steering column

A9130600001
length of the steering column 220 mm

A9130600002
length of the steering column 280 mm

A9130600003
length of the steering column 340 mm
angle adjustable via clamped pivot point on the fork bridge



Accessories

SKS mudguard

A9130500003
removable



Mudguard MX Colour: orange

A9130500002
permanently mounted



Mudguard MX Colour: black

A9130500004
permanently mounted



Aluminum luggage rack „cargo systems”

A9110903290 € 188.80

Colour: black
max. load weight on the loading platform 10 kg,
max. 15 kg in combination with luggage bags
at lower loading height.



Bottle holder incl. drinking bottle

A9130900001 € 23.90

mounted on the left side of the battery holder
rail. Not in combination with double battery



LED headlights

A9130900004 € 121.80

LED front light powered by the drive system
battery. Rear light as plug-in light with battery
power supply. Lighting set not required.

Rear-view mirror on the left

A9130900002 € 93.30

Rear-view mirror on the right

A9130900003 € 93.30

Options for users with reduced functions

Twist grip Triride

A9130200019 € 46.70

instead of the thumb throttle handle the twist
grip is installed.
Mounted on the right side of the handlebar.



Gas assistance for WHEEL-E

A9130200017 € 19.20

tool-free clamping on the twist grip of the
WHEEL-E



Gas operation Tetra

A9130400100 € 344.40

right side



Brake operation Tetra

A9130400101 € 375.80

Brake operating on the left with brake operation
via rotary function with hand rest



Chest strap for upper body fixation

size S (105 cm) € 241.80
A9000800016

size M (125 cm) € 241.80
A9000800013

size L (145 cm) € 241.80
A9000800014

size XL (165 cm) € 241.80
A9000800017

consisting of elastic webbing material in
combination with neoprene with Velcro and
operating tab.
Colour: black



Lever extension excenter lock for Quick Snap connection


Makes it easier to open and close the Quick
Snap adapter. Tetra version

A9110700005 € 73.90

straight design
only in combination with rigid frame adapter

A9110700006 € 73.90

curved design
only in combination with front adapter



<p>Adaptation help for centering pin</p> <p><input type="checkbox"/> A9110700002 € 55.00 to simplify the insertion of the central adapter.</p> 	<p>Parking stand with operating lever for high operating comfort</p> <p><input type="checkbox"/> A9110900004 € 58.70 with longer operating lever</p>	
Frame adapter		
<p>SPEEDY rigid frame adapter with „Quick-Snap-Adaption“ with quick-release clamp and centering pin holder Colour: black</p> <p><input type="checkbox"/> A9110700141 € 537.20</p> <p>Heavy Duty version</p> <p><input type="checkbox"/> A9110700143 € 816.70 up to a maximum body weight of 200 kg</p> <p>Please select accordingly:</p> <p>Centring pin for Ø 25 mm for A1 and A2</p> <p><input type="checkbox"/> A9110700022</p> <p>Centring pin for Ø 40 mm for 4all, 4all Ergo, 4you Ergo, F2, F3 and 4teen</p> <p><input type="checkbox"/> A9110700023</p> 	<p>TRAVELER folding frame adapter with „Quick-Snap-Adaption“</p> <p><input type="checkbox"/> A9110701300 € 537.20 with quick-release clamp and centering pin holder on the folding mechanism Colour: black</p> 	<p>BUDDY rigid frame adapter with „Quick-Snap-Adaption“</p> <p><input type="checkbox"/> A9110701400 € 537.20 with quick-release clamp and centering pin mount on the frame Colour: black</p> 
<p>Wheelchair front adapter with universal clamp and „Quick-Snap-Adaption“</p> <p><input type="checkbox"/> A9110701500 € 572.30 Adapter attachment to the wheelchair stem.</p> <p>Mounting situation for PA products: Straight front end from SB 31 cm to 54 cm V-shaped front end from SB 30 cm to 50 cm</p> <p>Please select the frame clamp below:</p> <p><input checked="" type="checkbox"/> A9110701525 € 572.30 Tube diameter 23-27 mm</p> <p>PLEASE NOTE: If the exact tube diameter is not specified, an adapter for a tube diameter of 23 - 27 mm will be delivered as standard.</p> <p><input checked="" type="checkbox"/> clamping sleeve [80 mm] included as standard for close clearance between footrest and drive wheel.</p> <p>PLEASE NOTE: If not specified, the clamping sleeve with 80 mm length will be delivered.</p>	<p>Wheelchair front adapter Heavy Duty > 120 kg - 200 kg with universal clamp and „Quick-Snap-Adaption“</p> <p><input type="checkbox"/> A9199701510 € 871.60 Adapter with reinforced double-tube steering head connection for adaptations with wheelchairs for user weights > 120 kg - 200 kg Adapter attachment to the wheelchair frame stem</p>  <p><input checked="" type="checkbox"/> A9110701530 € 572.30 Tube diameter 28-32 mm</p> <p><input type="checkbox"/> clamping sleeve [110 mm] € 46.70 A9199000117 for medium distance between footrest and drive wheel.</p>	<p><input checked="" type="checkbox"/> A9110701535 € 572.30 Tube diameter 33-37 mm</p> <p><input type="checkbox"/> clamping sleeve [170 mm] € 56.60 A9199000154 for a wide distance between footrest and drive wheel.</p>

Clamp lever with 3-way star handle for actuating the frame clamping

A9199701507 € 106.10

mounted on the right and left for easy terminal operation with limited manual function
Only in combination with front adapter



Wheelchair-dependent accessories


Extended bike plate

for SPEEDY A1 / A2 A9110100100-002 € 370.00

for BUDDY A9110100100-003 € 370.00

for TRAVELER A9110100100-001 € 370.00

4° camber integrated
Track independently adjustable



Wheelbase extension via „grid mechanism with pulling operation“

for SPEEDY 4all / 4teen / 4you / TRAVELER 4you A9110102001 € 564.70

removable, with continuous axle tube, track width independently adjustable

Please select the wheel camber:

1° drive wheel camber

4° drive wheel camber

6° drive wheel camber



Wheelbase extension angle-adjustable and removable

for SPEEDY F2 / F3 A9015100005 € 746.80

with continuous axle tube, height adjustment via gear joint adjustable in 6° angular increments wheel position displacement approx. 14 cm adjustable track

Attachment possibility must be checked by PRO ACTIV.

Only in combination with base plate. Up to 120 kg load capacity.

Please select the wheel camber below:

1° drive wheel camber

4° drive wheel camber

6° drive wheel camber



Technical vehicle inspection and accessories

Activation max. speed 15 km/h

A9130900101 € 145.00

use only on private property, outside the scope of application of road traffic regulations

Front brake

A9125500120 € 152.50

as hill start assist with anti-roll back device, compatible for installation in combination with the wheelchair knee lever brake.

Required for WHEEL-E with 15 km/h individual operating permit.



Activation max. speed 24 km/h

A9130900106 € 145.00

use only on private property, outside the scope of application of road traffic regulations

Required in combination with wheelbase extension:

retaining strips adjustment range 295 mm 4all

A9000700000 € 216.70

retaining strips adjustment range 295 mm Ergo

A9000700173 € 430.30

Please select:

Excentric lever (2 pcs.) for removable knee lever brake

A9000700207 € 152.50

Knurled screws for adjustable knee lever brake with 2 position stops

A9010708100 € 129.50

LEGAL WARNING for the WHEEL-E towing device.

Equipping the WHEEL-E with a console for a maximum speed of 6 km/h.
 The factory settings of the WHEEL-E PRO ACTIV are designed for an electric drive system that supports a maximum speed of 6 km/h.
 Consoles for a top speed of 15 km/h or 24 km/h are optionally available. PRO ACTIV points out that the specific laws of your country must be observed for the use of the 15 km/h and 24 km/h consoles. When ordering, please confirm with your signature that you have taken note of this fact.

Services**Consultation, clarification of the installation situation and measurement through PRO ACTIV**

A9990000000 € 354.50

Handover and instruction of the handbike / traction device by PRO ACTIV

A9990000011 € 354.50

Attachment of the handbike / traction device to the wheelchair by PRO ACTIV at specialist dealer**on a PRO ACTIV wheelchair**

A9990000012 € 230.00

product no.: _____

on a third-party product

A9990000013 € 366.90

manufacturer: _____

model: _____

frame tube Ø ____ mm

Attachment of the handbike / traction device to the wheelchair through specialist dealer

including mounting of the frame adapter

on a PRO ACTIV wheelchair

A9990000008

product no.: _____

on a third-party product

A9990000005

manufacturer: _____

model: _____

frame tube Ø ____ mm

Notes and remarks:**Location, Date**

The signature field is intended for the consultant and insured person. This confirms that the intended equipment is intended for the insured person.
 - 6 - 2026-03-03 / The prices are plus statutory value-added tax and shipping costs

Signature End customer/dealer

= included in the basic price = Position selected