



PRO ACTIV Reha-Technik GmbH
 Im Hofstätt 11
 D-72359 Dotternhausen
 Phone +49 (0)7427 9480-0
 Fax +49 (0)7427 9480-7025
 E-mail: info@proactiv-gmbh.de
 www.proactiv-gmbh.com



☐ Offer date: _____

☐ Order date: _____
 according to GTC PRO ACTIV Reha-Technik GmbH

Consultant: _____

*** Required fields, please fill in completely!**

Order address / customer:

Cust. no.:* _____

Last name:* _____

First name:* _____

Date of birth: _____

Street:* _____

Postcode / city:* _____

Phone:* (home) _____

Phone: (work) _____

E-mail:* _____

Collection: ☐ Yes ☐ No

Biometric data:

Body height: _____

Weight: _____

Inside leg: _____

M-Porter A9590000002 **€ 1,794.12**

- Aluminium hydroformed frame, Frame height: 52 cm
- Control set with steering stop 11/8 inch
- Fixed mudguard at front and rear
- Ursus "Jumbo" two-leg stand
- Bell
- Front and rear lighting
- Stainless steel load carrier
- Compliant with road traffic regulations

Ready-to-ride weight: starting at 21.90 kg

Permitted total weight: max. 200 kg

Defined components:

Rigid aluminium fork

Colour selection: see colour design at the back

Drive wheel 24"

Continental "Air King" 24x2.4"

(60-507), with reflective strips, load capacity LI38 = 132 kg per tyre
 Colour: black

Nuvinci N360 hub gear

robust, infinitely adjustable planetary gear

Magura HS22 hydraulic rim brake

IDM at the front, with self-manufactured brake support at the rear

Ergotec Viper saddle support 31.6 mm

Length: 400 mm; colour: black

Selle Royal ARIEL "Strengtex" saddle

Width: 16 cm; length: 27 cm

Colour: black

Riser handlebar

Width: 680 mm

Herrmans "Diamond lock" grips

2-component grip with double clamping

SHIELD 5650L Abus frame lock

Security level 9

Front end

Front end 6°

☒ A9500400501-110
Length: 110 mm
Clamping: Oversize 31.8 mm
Colour: black



Ergotec Swell RE front end

☐ A9500400550 € 39.58
Length: 100 mm
Clamping: Oversize 31.8 mm
angle-adjustable
Colour: black



Pedals

Offroad pedals

☒ A9500500300
CR-MO axle,
with reflector (K-mark)
Platform with jagged profile



Touring pedals

☒ A9500500310
CR-MO axle,
with reflector (K-mark)



Ergo Grip pedals

☐ A9500500320 € 26.13
CR-MO axle, black, plastic / sand paper on
both sides, integrated reflectors



Colour design – frame and attachment parts

For the following components, a colour can be chosen from the information sheet “Pedelec Colour Design”:

RAL powder paint or “Dormant” powder paint or “Liquid coating” custom paint

Frame colour:

☒ RAL powder paint and “Dormant” powder paint:

Please enter the last 3 digits of the colour:

A9500100170 - in words (colour): _____

☐ “Liquid coating” custom paint:

Please enter the last 2 digits of the custom colour:

A9500100170 - 5 in words (colour): _____

€ 161.68

Fork colour for the “Rigid aluminium fork” version:

☒ Fork with frame colour:

A9500400010

**Individual colour selection for the “Rigid aluminium fork”
different from the frame colour:**

☐ RAL powder paint and “Dormant” powder paint:

A9500400020

Please enter the last 3 digits of the colour:

A9500100170 - in words (colour): _____

€ 57.56

☐ “Liquid coating” custom paint:

A9500400030

Please enter the last 2 digits of the custom colour:

A9500100172 - 5 in words (colour): _____

€ 98.32

Notes and remarks:

Leasing

☐ Yes

☐ Still being discussed

☐ No

Leasing companies: _____ (please fill in the “Pedelec leasing form”!)

I hereby confirm the information provided above and agree to receiving a written offer without obligation.

Location, date End customer signature (offer)

Location, date Consultant signature (offer)

Location, date End customer signature (order)

Location, date Consultant signature (order)

Info sheet: Colour Designs Pedelecs

Colour Designs - Powder Coatings

Information on current colour designs and sample images at: www.proactiv-gmbh.com/colours/

Powder Coatings

white RAL 9016 <input checked="" type="checkbox"/> gloss A9500100170-023	anthracite RAL 7016 <input checked="" type="checkbox"/> gloss A9500100170-001 <input checked="" type="checkbox"/> matt A9500100170-101	black RAL 9005 <input checked="" type="checkbox"/> gloss A9500100170-003 <input checked="" type="checkbox"/> matt A9500100170-019
pure orange RAL 2004 <input checked="" type="checkbox"/> gloss A9500100170-022 <input checked="" type="checkbox"/> matt A9500100170-122	yellow RAL 1023 <input checked="" type="checkbox"/> gloss A9500100170-046 <input checked="" type="checkbox"/> matt A9500100170-146	red RAL 3001 <input checked="" type="checkbox"/> gloss A9500100170-047 <input checked="" type="checkbox"/> matt A9500100170-147
„Dormant“ Powder Coatings		
sky blue <input checked="" type="checkbox"/> gloss A9500100170-017	electric green <input checked="" type="checkbox"/> gloss A9500100170-015	violet <input checked="" type="checkbox"/> gloss A9500100170-044
fire red <input checked="" type="checkbox"/> gloss A9500100170-018		
Powder Coatings “Trend colours”		
athletic grey <input checked="" type="checkbox"/> gloss A9500100170-303 € 30.59 <input checked="" type="checkbox"/> matt A9500100170-304 € 30.59	muddy grey <input checked="" type="checkbox"/> gloss A9500100170-305 € 30.59 <input checked="" type="checkbox"/> matt A9500100170-306 € 30.59	lagoon blue <input checked="" type="checkbox"/> gloss A9500100170-301 € 30.59 <input checked="" type="checkbox"/> matt A9500100170-302 € 30.59

Colour Designs - Custom Colours ‘Liquid Coating’

€ 161.68

storm (ice blue) <input type="checkbox"/> gloss A9500100170-500 <input type="checkbox"/> matt A9500100170-501	coffee <input type="checkbox"/> gloss A9500100170-520 <input type="checkbox"/> matt A9500100170-521	prairie (rose gold) <input type="checkbox"/> gloss A9500100170-504 <input type="checkbox"/> matt A9500100170-505
blue ocean <input type="checkbox"/> gloss A9500100170-522 <input type="checkbox"/> matt A9500100170-523	mint <input type="checkbox"/> gloss A9500100170-528 <input type="checkbox"/> matt A9500100170-529	