

Im Hofstätt 11  
D-72359 Dotternhausen  
Phone +49 (0)7427 9480-0  
Fax +49 (0)7427 9480-7025  
e-mail info@proactiv-gmbh.de  
www.proactiv-gmbh.com



## NJ1 compact bike

with *dérailleurs* or free-wheel hub  
and disk brake

Maximum user weight limit 100 kg

- girder chassis with reinforcing brace
- drive unit, can be decoupled via the eccentric quick-lock system
- A-shaped handles with integrated switching and brake fittings
- infinitely adjustable backrest
- 6° wheel camber
- hydraulic disc brake system "Magura"
- mechanical parking brake with parking lock
- safety flag
- black drive unit

A9100100001 € 5,100.00

Quotation date: \_\_\_\_\_

Consultation was carried out by  
Staff at the medical supplies shop: \_\_\_\_\_

PRO ACTIV field representative: \_\_\_\_\_

Order date: \_\_\_\_\_

According to GTC PRO ACTIV Reha-Technik GmbH

### Medical suppliers/company address:

Cust. no.: \_\_\_\_\_

Name: \_\_\_\_\_

Street: \_\_\_\_\_

Postcode / city: \_\_\_\_\_

Phone: \_\_\_\_\_

Fax: \_\_\_\_\_

E-Mail: \_\_\_\_\_

Contact: \_\_\_\_\_

Invoicing address, other (see below)

Cust. no.: \_\_\_\_\_

Name: \_\_\_\_\_

Street: \_\_\_\_\_

Postcode / city: \_\_\_\_\_

Shipping address, other (see below)

Cust. no.: \_\_\_\_\_

Name: \_\_\_\_\_

Street: \_\_\_\_\_

Postcode / city: \_\_\_\_\_

### Customer name/commission:

Surname, forename: \_\_\_\_\_

Street: \_\_\_\_\_

Postcode / city: \_\_\_\_\_

Phone: \_\_\_\_\_

E-Mail: \_\_\_\_\_

D.O.B: \_\_\_\_\_

Paying organisation: \_\_\_\_\_

Indication: \_\_\_\_\_

### Body size

Body height: \_\_\_\_\_ cm

Weight: \_\_\_\_\_ kg

## Chassis colour

Information on current colour designs and example images under: [www.proactiv-gmbh.de/farben.php](http://www.proactiv-gmbh.de/farben.php)

### White RAL 9016

- gloss A9000100170-023  
 matt A9000100170-123 € 40.10

### Brilliant silver

- gloss A9000100170-011

### Anthracite RAL 7016

- gloss A9000100170-001  
 matt A9000100170-101 € 40.10

### Black RAL 9005

- gloss A9000100170-003  
 matt A9000100170-019 € 40.10

### Deep black, silver effect

- gloss A9000100170-007

### Pure orange RAL 2004

- gloss A9000100170-022

### Fire red

- gloss A9000100170-018  
 matt A9000100170-118 € 40.10

### Pink magenta RAL 4010

- gloss A9000100170-037  
 matt A9000100170-137 € 40.10

### Sky blue

- gloss A9000100170-017

### Pastel blue RAL 5024

- gloss A9000100170-036  
 matt A9000100170-136 € 40.10

### Marine blue

- gloss A9000100170-045  
 matt A9000100170-145 € 40.10

### Electric green

- gloss A9000100170-015  
 matt A9000100170-115 € 40.10

### Lemon green

- gloss A9000100170-005  
 matt A9000100170-105 € 40.10

### Violet

- gloss A9000100170-044

### Powder coating transparent

- gloss A9000100170-029

### Yellow RAL 1023

- gloss A9000100170-046  
 matt A9000100170-146 € 40.10

### Red RAL 3001

- gloss A9000100170-047  
 matt A9000100170-147 € 40.10

### Lagoon blue

- gloss A9000100170-052  
 matt A9000100170-152 € 40.10

### Athletic grey

- gloss A9000100170-050  
 matt A9000100170-150 € 40.10

### Muddy grey

- gloss A9000100170-051  
 matt A9000100170-151 € 40.10

## Dérailleur

### Shimano Deore 1x11 Rapid Fire 11-speed dérailleur

- A9110200032  
 Mechanical dérailleur gear with operation on the handle. 11-speed cassette 11-42 teeth

#### Chain ring selection:

- 32 teeth  
 38 teeth



### Shimano Deore XT 1x11 Rapid Fire 11-speed dérailleur

- A9110200034 € 182.90  
 Mechanical dérailleur gear with operation on the handle. 11-speed cassette 11-42 teeth top chain ring with chain guards

#### Chain ring selection:

- 32 teeth  
 38 teeth

### Shimano Deore XT 1x11 DI2 with button operation 11-speed dérailleur

- A9110200035 € 1,414.50  
 electronic dérailleur with operation on the handle. 11-speed cassette 11-42 teeth

#### Chain ring selection:

- 32 teeth  
 38 teeth



### Shimano Deore XT 1x11 DI2 With chin operation 11-speed dérailleur

- A9110200033 € 1,552.50  
 Electronic dérailleur gear with chin operation (tetra equipment). 11-speed cassette 11-42 teeth

#### Chain ring selection:

- 32 teeth  
 38 teeth



## Hub gear

### Shimano "Alfine" 11-speed hub gear

A9100200026 € 563.50

"Rapid Fire" handle  
Gear ratio: Sprocket 32 teeth,  
Cassette 19 teeth



### Rohloff "Speedhub 500/14" 14-speed hub gear

A9100200027 € 1,426.00

"Grip Shift" twist grip shifter  
Gear ratio: Sprocket 32 teeth,  
Cassette 16 teeth



## Gear change accessories

### Mountain drive 2.5: 1 gear reduction for hills in combination with hub gear

A9100200003 € 759.00

Planetary gear set for doubling of the number  
of gears  
permits switching to the gear reduction for hills  
when stationary.  
With 38 teeth chain wheel as standard.  
Only in combination with grip size 40, 45, 48  
and 53 cm.



## Brake system

### Powerful brake system with brake disc Ø 203 mm

A9110300105 € 79.10

Recommended from a body weight of 80 kg.



### Double disc brake system "Magura" brake disc Ø 203 mm

A9100300100 € 222.10

No need for mechanical disk brake.  
Only in combination with the powerful  
Brake system  
Recommended from a body weight of 100 kg

## Dimensions

### Crank type and grip width

#### Bottom bracket shaft with Standard crank

Grip width 40 cm  
A9110500100-040

Grip width 45 cm  
A9110500100-045

Grip width 50 cm  
A9110500100-050









#### Bottom bracket shaft with Tetra pedal crank (flared)

With flared crank arms for max. space for  
upper body

Grip width 48 cm (176 mm)  
A9110500100-048

Grip width 53 cm (226 mm)  
A9110500100-053

Grip width 58 cm (276 mm)  
A9110500100-058

<p><b>Standard pedal crank</b></p> <ul style="list-style-type: none"> <li><input checked="" type="checkbox"/> Length 155 mm A9110600100-155</li> <li><input checked="" type="checkbox"/> Length 175 mm A9110600100-175</li> <li><input checked="" type="checkbox"/> Length 195 mm A9110600100-195</li> </ul> 	<p><b>Tetra pedal crank</b></p> <ul style="list-style-type: none"> <li><input type="checkbox"/> Length 155 mm A9100600105-155      <b>€ 129.50</b></li> <li><input type="checkbox"/> Length 175 mm A9100600105-175      <b>€ 129.50</b></li> <li><input type="checkbox"/> Individual length 195 mm to 140 mm      <b>€ 246.30</b> A9100600100</li> </ul> 	
<p><b>Fixed bottom bracket support</b></p> <ul style="list-style-type: none"> <li><input checked="" type="checkbox"/> S (X/Y 67/153) A9100501100-003</li> <li><input checked="" type="checkbox"/> M (X/Y 90/250) A9100501100-001</li> </ul> 	<p><b>Adjustable bottom bracket support</b></p> <p>the angle and length can be adjusted via the clamped rotary point on the fork bridge.</p> <p>For the pre-setting, please specify the bottom bracket support size that comes as close to your desired position.</p> <ul style="list-style-type: none"> <li><input type="checkbox"/> S (180 mm)      <b>€ 220.90</b> A9100500001-002</li> <li><input type="checkbox"/> M (240 mm)      <b>€ 220.90</b> A9100500001-003</li> </ul> <p>If no details are provided, the adjustable bottom bracket support shall be supplied in S.</p> 	
<b>Leg rest</b>		
<p><b>Leg rest separated</b></p> <ul style="list-style-type: none"> <li><input checked="" type="checkbox"/> A9100100008 Adaptation on the fork for short legs Only possible up to body size 165 cm</li> </ul> 	<p><b>Leg rest wrapping</b></p> <ul style="list-style-type: none"> <li><input checked="" type="checkbox"/> A9100100009 Adaptation on the fork</li> </ul> 	<p><b>Foot stopper for leg rest</b></p> <ul style="list-style-type: none"> <li><input type="checkbox"/> A9100100012      <b>€ 86.80</b></li> </ul> 
<p><b>Thigh anchorage</b></p> <ul style="list-style-type: none"> <li><input type="checkbox"/> A9100100013      <b>€ 204.70</b></li> </ul> 		

## Back equipment

### Backrest with filling plate and back cushion

A9100800002  
Cover: Black, breathable, open surface

### Body contour backrest upholstery

A9100800010 € 85.10  
with back cushion 3 cm  
Cover: Black, breathable, open surface




## Chassis

### Chassis geometry

Length 116 cm  
Seat width approx.  
39 cm  
A9100100000-002

Length 116 cm  
Seat width approx.  
43 cm  
A9100100000-003 € 194.70

Heavy Duty version  
up to 120 kg body  
weight  
A9100800013 € 405.10  
Recommended for off-road  
use



## Seat upholstery

### Open belt system seat upholstery

Chassis: 116 cm  
SB 39 cm, ST 40 cm

A9100801010-002

Chassis: 116 cm  
SB 43 cm, ST 40 cm

A9100801010-003



### Body Contour upholstery


breathable, cushioning effect

Chassis: 116 cm  
SB 39 cm, ST 40 cm

A9100801025-002

Chassis: 116 cm  
SB 43 cm, ST 40 cm

A9100801025-003



## Seat cushion

### Standard seat cushion

Cover: Black, breathable, open surface  
Seat cushion height: 4 cm

Chassis: 116 cm  
SB 39 cm, ST 40 cm

A9100000190-002 € 106.60

Chassis: 116 cm  
SB 43 cm, ST 40 cm

A9100000190-003 € 106.60



### Seat cushion with thigh guide

Cover: Black, breathable, open surface  
Seat cushion height: 4 cm  
Only in combination with thigh anchorage.

Chassis: 116 cm  
SB 39 cm, ST 40 cm

A9100800103-002 € 194.60

Chassis: 116 cm  
SB 43 cm, ST 40 cm

A9100800103-003 € 194.60



## Wheels

### 26x1" Euroline running wheels with double-wall rim

A9100104260  
(ETRTO 559)  
Running wheels front and back with 36 spokes

### 26x1" Mavic XC717 running wheels with double-wall rim

A9100104262 **€ 283.60**  
(ETRTO 559)  
extra light version: rim weight 395 g,  
Running wheels front with 32 spokes  
Running wheels back with 18 spokes  
(Weight saving on the bike approx. 600 g)

## Tyres

### 26x1" Schwalbe "One"

A9100105160  
(ETRTO 559-25)  
V-Guard WheelStar



### 26x2" Schwalbe "Marathon Plus Tour"

A9100105140 **€ 67.60**  
(ETRTO 559-50)  
Touring tyres



## Bike accessories

### Plug-in protection plate

A9100502100 **€ 55.30**



### Leg guard

A9100100231 **€ 162.70**  
Guard for MTB tyres  
Screwed onto the fork tube  
Colour: Matt black

### Collision guard clamped on axle tube

A9100800100 **€ 207.50**  
with 40 mm mounting clip  
Colour: black



### Neck rest

A9100800150 **€ 248.40**  
height-adjustable  
angle-adjustable



### Hip safety belt

A9100800200 **€ 112.90**



### SIGMA BC 9.16 ATS Bike Computer

A9110901000 **€ 66.80**  
wireless, with bracket



### SIGMA BC 14.16 STS CAD Bike Computer

A9110900018 **€ 163.20**  
with pedal frequency set  
wireless



### Bracket ø 22 mm x 42 mm for fixing of accessories (speedometer, lighting or similar)

A9110903271 **Item price € 22.10**  
Colour: Natural anodised

\_\_\_\_\_ Items



### SIGMA lighting set

A9120900019 **€ 117.90**  
Consisting of front light and rear lights,  
chargeable via Micro USB socket,  
incl. Micro USB charging cable.  
Lighting set compliant with road traffic  
regulations



## Options for users with reduced function

### Palm rest on the handle

A9110601200 € 97.80

Padded, bolted, angle-adjustable.  
Not in combination with other handles.



### Tetra hand support, pivotable, adjustable grip distance

A9110601101 € 161.30



### Tetra hand support, pivotable, adjustable grip width

A9110601102 € 182.40

Very stiff outer support as round tube ø 20 mm with foam rubber cover  
Cannot be pivoted in conjunction with Tetra hand support



### Tetra shell handles

A9110601100 € 164.10



### Notes and remarks:

### Location, Date

The signature field is intended for the consultant and insured person. This confirms that the intended equipment is intended for the insured person.  
- 7 - 2023-02-16 / The prices are plus statutory value-added tax and shipping costs  = Included in the basic price  = Position required

### Signature End customer/dealer