

Mobilität wird einfach!

Im Hofstätt 11
D- 72359 Dotternhausen
Phone +49 (0)7427 9480-0
Fax +49 (0)7427 9480-7025
E-Mail info@proactiv-gmbh.de
www.proactiv-gmbh.com



Quotation date: _____

Order date: _____

according to GTC PRO ACTIV GmbH

Medical suppliers/company address:

Cust. no.: _____
Name: _____
Street: _____
Postcode, city: _____
Phone: _____
Fax: _____
E-Mail: _____
Contact: _____

Invoicing address, other (see below)

Cust. no.: _____
Name: _____
Street: _____
Postcode, city: _____

Shipping address, other (see below)

Cust. no.: _____
Name: _____
Street: _____
Postcode, city: _____

Customer name/commission:

Name: _____
Street: _____
Postcode, city: _____
Phone: _____
Fax: _____
E-Mail: _____
DoB: _____
Paying organisation: _____
Indication: _____

NJ1 Compact bike RS

with dérailleurs or free-wheel hub and disk brake

Maximum user weight limit 100 kg

- girder chassis with reinforcing brace
- drive unit, can be decoupled via the eccentric quick-lock system
- A-shaped handles with integrated switching and brake fittings
- infinitely adjustable backrest
- 6° wheel camber
- hydraulic disc brake system "Magura"
- Shimano V-brake (rim brake) incl. locking mechanism
- safety flag

9100100000 € 3,940.00

Services:

Delivery, instruction and handing over at the customer by PRO ACTIV

9990000003 € 425.00

plus shipping costs

Body size

Body height: _____ cm

Weight: _____ kg

Chassis colour

Information on current colour designs and example images under:
www.proactiv-gmbh.de/farben.php

White RAL 9016

- gloss 9000100170-023
 matt 9000100170-123 € 33.30

Brilliant silver

- gloss 9000100170-011

Anthracite RAL 7016

- gloss 9000100170-001
 matt 9000100170-101 € 33.30

Black RAL 9005

- gloss 9000100170-003
 matt 9000100170-019 € 33.30

Deep black, silver effect

- gloss 9000100170-007

Pure orange RAL 2004

- gloss 9000100170-022

Fire Red

- gloss 9000100170-018
 matt 9000100170-118 € 33.30

Pink magenta RAL 4010

- gloss 9000100170-037
 matt 9000100170-137 € 33.30

Sky blue

- gloss 9000100170-017

Pastel blue RAL 5024

- gloss 9000100170-036
 matt 9000100170-136 € 33.30

Marine blue

- gloss 9000100170-045
 matt 9000100170-145 € 33.30

Apple green

- gloss 9000100170-014
 matt 9000100170-114 € 33.30

Lemon green

- gloss 9000100170-005
 matt 9000100170-105 € 33.30

Violet

- gloss 9000100170-044

Yellow RAL 1023

- gloss 9000100170-046
 matt 9000100170-146 € 33.30

Red RAL 3001

- gloss 9000100170-047
 matt 9000100170-147 € 33.30

Transparent, powder coated

- gloss 9000100170-029

Drive unit colour

Silver, anodised drive unit

- 9110100170-000

Black Shadow, deep black drive unit RAL 9005, smooth gloss

- 9110100170-001 € 103.80

- fork tube R and L,
- fork bridge, top and bottom
- bottom bracket support

Dérailleur

Shimano Deore 3x9 Rapide Fire 27-speed dérailleur

- 9100200018
Gear ratio:
sprocket 22-32-44 teeth,
cassette 11-13-15-17-20-23-26-30-32 teeth



Shimano Deore XT 1x11 Rapid Fire 11- speed chain shift derailleur

- 9110200034 € 151.00
mechanical derailleur operation from the
handle. 11-way cassette 11-42 teeth



Shimano Deore XT 1x11 Di2 Button operation 11- speed chain shift derailleur

- 9110200035 € 1,175.00
electronic derailleur operation from the
handle. 11-way cassette 11-42 teeth



Shimano Deore XT 1x11 Di2 Rapid Fire 11- speed chain shift dérailleur

- 9110200032 € 1,138.00
electronic derailleur operation from the
handle. 11-way cassette 11-42 teeth



Shimano Deore XT 1x11 Di2 chin operation 11- speed chain shift derailleur

- 9110200033 € 1,245.00
electronic derailleur with chin operation
(Execution Tetra).
11-way cassette 11-42 teeth



Hub drive system

Shimano "Alfine" 11-speed hub gear

- 9100200026 € 88.60
Handle "Rapid Fire"
Gear ratio: Sprocket 32 teeth,
Cassette 20 teeth
(with bottom bracket support L – 19 teeth)



Rohloff "Speedhub 500/14" 14-speed hub gear

- 9100200027 € 1,037.90
"Grip Shift" twist grip shifter
Gear ratio: Sprocket 32 teeth,
Cassette 16 teeth



Gear change accessories

Mountain drive 2.5: 1 gear reduction for hills in combination with hub gear

9100200003 € 499.80

Planetary gear set for doubling of the number of gears

All gears (the hub gear and the gear reduction for hills) can be changed when stationary!

In combination with grip width 40 or 45 cm!



Mountain drive 2.5: 1 gear reduction for hills in combination with chain shift

9100200006 € 371.30

Planetary gear set

Replaces the 3-fold sprocket and allows the gear reduction for hills to be switched when stationary!

In combination with grip width 40 or 45 cm!



Gear accessories for chain shift

"Grip Shift" twist grip shifter 3-way, mounting left side

9110200017 € 43.30

SRAM Centera for 3-way dérailleur

In combination with Shimano Deore or XT chain shift!



"Grip Shift" twist grip shifter 9-way, mounting right side

9110200016 € 43.30

SRAM Centera for 9-fold dérailleur

In combination with Shimano Deore or XT chain shift!



Hand lever to directly switch the 3-way dérailleur

9100200020 € 55.80

direct hand operation of the dérailleur on the 3-way sprocket



Brake system

Powerful HD brake system "Magura"

9100300101 € 78.60

Recommended from a body weight of 80 kg



Second independent disk brake system "Magura"

9100300100 € 183.90

Shimano V-brake rim brake omitted

Only in combination with the powerful HD brake system "Magura"!

Chassis


Chassis geometry

Length 104 cm
Seat width approx. 39 cm
9100100000-001

Length 116 cm
Seat width approx. 39 cm
9100100000-002

Length 116 cm € 161.20
Seat width approx. 43 cm
9100100000-003

Heavy Duty version € 355.50
up to 120 kg body weight
9100800013




Dimensions

Pedal bearing shaft

Grip width 40 cm
9100500100-040

Grip width 45 cm
9100500100-045

Grip width 50 cm
9100500100-050




Fixed bottom bracket support

XS (X/Y 113/130)
9100501100-004

S (X/Y 67/153)
9100501100-003

M (X/Y 90/250)
9100501100-001



Adjustable bottom bracket support

9100500001 € 182.90

Adjustable angle via clamped hinge on the fork bridge and in their length.



For the pre-setting, please specify the bottom bracket support size that comes as close to your desired position.

Size: _____

If no details are provided, the adjustable bottom bracket support shall be supplied in S.



Pedal crank

Length 155 mm
9100600100-155

Length 175 mm
9100600100-175

Length 195 mm
9100600100-195



Customer-specific crank length according to dimension

_____ mm (minimum length 125 mm)

9100600101 € 155.80

Leg rest

Leg rest separated

910010008
Adaptation on the fork
for short legs
Only possible up to body size 165 cm!




Leg rest wrapping

910010009
Adaptation on the fork



Leg rest adapted to the chassis

910010010 € 400.40
Legs not influenced by steering movement.
Increases the vehicle width!
Not in combination with the leg shield.



Foot stopper for leg rest

910010012 € 71.90



Thigh anchorage

910010013 € 127.90



Back equipment

Backrest with filler and 1.5 cm foam padding

910080000



Body contour back padding

910080010 € 123.80
with back padding 3 cm
Cover: black/breathable



Seat upholstery


Open belt system seat upholstery

9100801010



Body contour seat upholstery

9100801025



Seat cushion

shape	seat cushion "Standard" 9100000190	seat cushion with upper leg guidance 9100800103
dimensions in cm		
Chassis: 104 SB: 39 ST: 30	-001 <input type="checkbox"/>	-001 <input type="checkbox"/>
Chassis: 116 SB: 39 ST: 40	-002 <input type="checkbox"/>	-002 <input type="checkbox"/>
Chassis: 116 SB: 43 ST: 40	-003 <input type="checkbox"/>	-003 <input type="checkbox"/>
price	€ 88,30	€ 161,10

Cover: black, breathable, open fabric surface
height: 4cm

seat cushion with upper leg guidance only in combination with hip safety belt.

Wheels

Running wheels with double-wall rim

26" (559)
9100104260

Tyres

Schwalbe Durano

26x1.1" (559-28)
9100105000
Racing-Slic



Schwalbe Marathon Plus Tour

26x2,0" (559-50)
9100105140 € 58,80



Bike accessories

Plug-in protection plate

9110502100 € 45.80



Leg guard

9100100231 € 141.50
Guard for MTB tyres
Screwed onto the fork tube
Colour: black matt

Collision guard clamped on the axle tube

9100800100 € 171.80



Neck rest

9100800150 € 206.00
height-adjustable
angle-adjustable



Hip safety belt

9100800200 € 93.50



**SIGMA BC 1200
Bike Computer**

9110901000 € 55.30
wireless



Your additional notes:

**SIGMA BC 1606L DTS
Bike Computer**

9110900018 € 135.20
with pedal frequency set
wireless



**Bracket ø 22 mm x 42 mm for fixing
of accessories
(speedometer, lighting or similar)**

9110903271 Item price €
18.30

Colour: natural anodised



_____ Items

**Lighting set SIGMA Lightster LED
with bracket**

9100900019 € 102.50
with battery charger and set of batteries
(4 pcs. NICD) front and rear lights



Options for users with reduced function

**Tetra hand support, pivotable,
adjustable grip distance**

9110601101 € 140.30



Tetra shell handles

9110601100 € 135.90



Tetra backhand support

9110601102 € 158.60
pivotable, adjustable distance to the bike
handle



I herewith confirm the details listed above

Location, Date

Signature (End customer)

Location, Date

Signature (Dealer)