



Quotation date: _____
Consultation was carried out by
Staff at the medical supplies shop: _____

Enquiry no.: _____

Order date: _____
according to GTC PRO ACTIV Reha-Technik GmbH

Order no.: _____

Medical suppliers / company address:

Cust. no.: _____

Name: _____

Street: _____

Postcode / city: _____

Phone: _____

Fax: _____

E-Mail: _____

Contact: _____

Different invoicing address (see below)

Cust. no.: _____

Name: _____

Street: _____

Postcode / city: _____

Different shipping address (see below)

Cust. no.: _____

Name: _____

Street: _____

Postcode / city: _____

Customer name / commission:

ID / Comm. No.: _____

Surname, forename: _____

Street: _____

Postcode / city: _____

Phone: _____

E-Mail: _____

D.o.B.: _____

Paying organisation: _____

Indication: _____

NJ1 el. compact bike

approved up to 100 kg payload

- girder chassis with integrated rechargeable battery holder
- drive unit, can be decoupled via the eccentric quick-lock system
- A-shaped handles with integrated shift and brake controls
- Backrest continuously angle-adjustable
- 6° wheel camber
- Hydraulic disc brake system
- additional mechanical disc brake with parking lock
- Safety flag with holder
- Leg protection guard
- Drive unit black

Neodrives Z15 drive unit

„The powerful hub drive for maximum performance with TFT display and 11-speed Shimano derailleur gears”
torque 12/40 Nm, power 250 watt

Wheel size: 26"

A9100100011-026 € 8,780.00

>> Expected to be available from April 2026 <<

Body size

Body height: _____ cm

Weight: _____ kg

Chassis Colour

Information on current colour designs and example images can be found at: <https://www.proactiv-gmbh.com/colours/>

white RAL 9016

- gloss A9000100170-023
 matt A9000100170-123 € 43.40

brilliant silver

- gloss A9000100170-011

anthracite RAL 7016

- gloss A9000100170-001
 matt A9000100170-101 € 43.40

black RAL 9005

- gloss A9000100170-003
 matt A9000100170-019 € 43.40

deep black effect silver

- gloss A9000100170-007

pure orange RAL 2004

- gloss A9000100170-022

fire red

- gloss A9000100170-018
 matt A9000100170-118 € 43.40

pink magenta RAL 4010

- gloss A9000100170-037
 matt A9000100170-137 € 43.40

sky blue

- gloss A9000100170-017

pastel blue RAL 5024

- gloss A9000100170-036
 matt A9000100170-136 € 43.40

navy blue

- gloss A9000100170-045
 matt A9000100170-145 € 43.40

electric green

- gloss A9000100170-015
 matt A9000100170-115 € 43.40

lemon green

- gloss A9000100170-005
 matt A9000100170-105 € 43.40

violet

- gloss A9000100170-044

blasted, transparent powder coated

- gloss A9000100170-029

yellow RAL 1023

- gloss A9000100170-046
 matt A9000100170-146 € 43.40

red RAL 3001

- gloss A9000100170-047
 matt A9000100170-147 € 43.40

lagoon blue

- gloss A9000100170-052
 matt A9000100170-152 € 43.40

athletic grey RAL 7046

- gloss A9000100170-050
 matt A9000100170-150 € 43.40

muddy grey RAL 7033

- gloss A9000100170-051
 matt A9000100170-151 € 43.40

peach RAL 3022

- gloss A9000100170-038
 matt A9000100170-138 € 43.40

Battery Neodrives

The specified ranges are dependent on the route profile, user weight and speed.

Battery neodrives DT

36V 16Ah 576Wh
range approx. 75 km

- A9100900109



Double battery neodrives DT

36V 16Ah 576Wh
range approx. 150 km

- A9100900110 € 1,412.40

plug-in device for easy switching between batteries via the magnetic connector

Derailleur gear

Shimano Deore 1x11 Rapid Fire 11-speed derailleur

- A9110200032
mechanical derailleur gear with operation on the handle. 11-speed cassette 11-42 teeth

Selection of the chainring:

- 32 teeth
 38 teeth



Shimano Deore XT 1x11 Rapid Fire 11-speed derailleur

- A9110200034 € 192.20
mechanical derailleur gear with control on the handle. 11-speed cassette 11-42 teeth top chain ring with chain guards

Selection of the chainring:

- 32 teeth
 38 teeth

Shimano Deore XT 1x11 DI2 with push button operation 11-speed derailleur

- A9110200035 € 1,486.10
electronic derailleur gear with operation on the handle. 11-speed cassette 11-42 teeth

Operation: right side left side

Selection of the chainring:

- 32 teeth
 38 teeth



**Shimano Deore XT 1x11 DI2
with chin control
11-speed derailleur gear system**

A9110200033 € 1,631.10

electronic derailleur with chin control
(tetra equipment)
11-speed cassette 11-42 teeth

Selection of the chainring:

32 teeth

38 teeth



Gearshift accessories

**Mountaindrive 2,5 : 1
Hill reduction
in combination with hub gear**

A9100200003 € 724.90

Planetary gearbox for doubling the gears.
Allows you to shift downhill while standing still.
As standard with 38 teeth chainring.
In combination with grip width 40 cm, 45 cm,
48 cm or 53 cm



Brake system

**Powerful braking system with brake
disc Ø 203 mm**

A9110300105 € 83.10

recommended from 80 kg body weight.



**Double disc brake system „Magura”
brake disc Ø 203 mm**

A9100300100 € 233.30

eliminates the mechanical disk brake.
Only in combination with the powerful braking
system.
recommended from 100 kg body weight.

**Mechanical back-pedal disc brake
with crank-free function**

A9110300004 € 738.40

compatible with derailleur and freewheel hub.
In combination with:
- Bottom bracket shaft reach 45 cm to 55 cm
- pedal crank length 155 mm or 175 mm
up to 120 kg body weight
Not in combination with Mountaindrive 2.5 : 1
hill reduction!

Dimensions

Crank type and grip width

Standard crank

- grip width 40 cm
A9110500100-040
- grip width 45 cm
A9110500100-045
- grip width 50 cm
A9110500100-050

Tetra crank (flared)

with extended crank arms for maximum upper body space

- grip width 48 cm (176 mm)
A9110500100-048
- grip width 53 cm (226 mm)
A9110500100-053
- grip width 58 cm (276 mm)
A9110500100-058

Standard pedal crank

- length 155 mm
A9110600100-155
- length 175 mm
A9110600100-175
- length 195 mm
A9110600100-195



Tetra pedal crank


- length 155 mm
A9100600105-155 € 136.10
- length 175 mm
A9100600105-175 € 136.10
- individual length
195 mm to 140 mm € 258.70
A9100600100



Fixed bottom bracket support

(Limited availability – only as long as stocks last. Alternative: Adjustable bottom bracket support)

- S (X/Y 67/153)
A9100501100-003
- M (X/Y 90/250)
A9100501100-001



Adjustable bottom bracket support

the angle and length can be adjusted via the clamped rotary point on the fork bridge and length.

For the pre-setting, please specify the bottom bracket support size that comes as close to your desired position.

- S (180 mm)
A9100500001-002 € 232.10
- M (240 mm)
A9100500001-003 € 232.10

If not specified, the adjustable bottom bracket support will be delivered in S.



Legrest

Divided legrest

- A9100100008
adaptation on the fork
for short legs
Only possible up to body size 165 cm



All-round legrest

- A9100100009
adaptation on the fork



Foot stop for legrest

- A9100100012 € 91.20




Upper thigh fixation <input type="checkbox"/> A9100100013 € 215.10 	
---	--

Back equipment

Backrest with filling plate and back cushion <input checked="" type="checkbox"/> A9100800002 Cover: Black, breathable, open surface	Backrest upholstery Body Contour <input type="checkbox"/> A9100800010 € 92.40 with back cushion 3 cm Cover: Black, breathable, open surface 
--	---

Chassis

Chassis geometry <input checked="" type="checkbox"/> length 116 cm SB approx. 39 cm A9100100000-002 <input type="checkbox"/> length 116 cm SB approx. 43 cm € 204.60 A9100100000-003	<input type="checkbox"/> Heavy Duty version up to 120 kg body weight A9100800013 € 425.60 (Recommended for off-road use) Required in addition to A9100100000*.	
--	---	--

Seat upholstery

Seat upholstery open belt system Chassis: 116 cm SB 39 cm, ST 40 cm <input checked="" type="checkbox"/> A9100801010-002	Chassis: 116 cm SB 43 cm, ST 40 cm <input type="checkbox"/> A9100801010-003	
---	---	---

Seat upholstery Body Contour Chassis: 116 cm SB 39 cm, ST 40 cm <input type="checkbox"/> A9100801025-002	Chassis: 116 cm SB 43 cm, ST 40 cm <input type="checkbox"/> A9100801025-003	
--	---	---

Seat cushion

Seat cushion standard Cover: Black, breathable, open surface Seat cushion height: 4 cm Chassis: 116 cm SB 39 cm, ST 40 cm <input type="checkbox"/> A9100000190-002 € 115.00	Chassis: 116 cm SB 43 cm, ST 40 cm <input type="checkbox"/> A9100000190-003 € 115.00	
--	---	---

Seat cushion with upper thigh guide Cover: Black, breathable, open surface Seat cushion height: 4 cm Only in combination with upper thigh fixation Chassis: 116 cm SB 39 cm, ST 40 cm <input type="checkbox"/> A9100800103-002 € 207.50	Chassis: 116 cm SB 43 cm, ST 40 cm <input type="checkbox"/> A9100800103-003 € 207.50	
---	---	---

Tyres

26x1.5" Schwalbe „Marathon“

A9100105100
(ETRTO 559-40)
touring tyre



26x2" Schwalbe „Marathon Plus Tour“

A9100105140 € 71.10
(ETRTO 559-50)
touring tyre



26x2.1" Schwalbe „Smart Sam“

A9100105150 € 74.60
(ETRTO 559-54)



Accessories Bike

Plug-in mudguard

A9100502100 € 58.10



Rear guard clamped on the axle tube

A9100800100 € 218.10



Neck rest

A9100800150 € 261.00
height-adjustable, angle-adjustable



Pelvic fixation belt

A9100800200 € 121.70



Bracket Ø 22 mm x 42 mm for mounting accessories (tacho, lighting, etc.)

A9110903271 pcs. á € 23.20
Colour: Natural anodised
_____ piece



SIGMA lighting set

A9120900019 € 123.90

consisting of front light and rear lights,
chargeable via Micro USB socket,
incl. Micro USB charging cable.
Lighting set compliant with road traffic
regulations



Options for users with reduced functions

Palm rest on the handle

A9110601200 € 102.80
upholstered, bolted, angle-adjustable.
Not in combination with other handles



Tetra hand support, pivotable, adjustable grip distance

A9110601101 € 169.40



Tetra backrest swivel and adjustable grip width

A9110601102 € 191.60

very stiff external support as round tube
Ø 20 mm with foam rubber coating



Tetra shell handles

A9110601100 € 172.40



Operating satellite mounted on handle

Mounting: right side left side

with existing support line (placed along a
brake line)

A9116200000 € 103.10

without existing support line (an additional
blind cable is required)

A9116200001 € 146.20

Services

Consultation, clarification of the installation situation and measurement through PRO ACTIV

A9990000000 € 354.50

Handover and instruction of the handbike / traction device by PRO ACTIV

A9990000011 € 354.50

Notes and remarks:

Location, Date

The signature field is intended for the consultant and insured person. This confirms that the intended equipment is intended for the insured person.

Signature End customer/dealer

- 7 - 2026-03-03 / The prices are plus statutory value-added tax and shipping costs = included in the basic price = Position selected