

Im Hofstätt 11
D-72359 Dotternhausen
Tel +49 (0)7427 9480-0
Fax +49 (0)7427 9480-7025
e-mail info@proactiv-gmbh.de
www.proactiv-gmbh.de



SPEEDY 4all Ergo **Active easy running rigid frame** **wheelchair**

Maximum user weight limit 120 kg

A9017100001 € 4,120.00

☐ **Quotation date:** _____

Consultation was carried out by
Staff at the medical supplies shop: _____

Enquiry no.: _____

☐ **Order date:** _____

according to GTC PRO ACTIV Reha-Technik GmbH

Order no.: _____

Medical suppliers / company address:

Cust. no.: _____
Name: _____
Street: _____
Postcode / city: _____
Phone: _____
Fax: _____
E-Mail: _____
Contact: _____

Shipping address: ☐ **same as company address**
☐ **Other (see below)**

Name: _____
Street: _____
Postcode / city: _____

Customer name / commission:

ID / Comm. No.: _____
Surname, forename: _____
Street: _____
Postcode / city: _____
Phone: _____
Fax: _____
E-Mail: _____
D.o.B.: _____
Paying organisation: _____

Invoicing address: ☐ **same as company address**
☐ **Other (see below)**

Cust. no.: _____
Name: _____
Street: _____
Postcode / city: _____
Phone: _____
Fax: _____
E-Mail: _____
Contact: _____

Basic size


Body size:

Body height: _____ cm

Weight: _____ kg


Lower leg length: _____ cm
(anatomic dimension: hollow of the knee to the foot sole)

Indication: _____



This information is necessary for checking. If the information is missing, a check of the relevant items cannot be made!

Wheelchair dimensions:



Seat width SB	☐ 30	☐ 31	☐ 32	☐ 33	☐ 34	☐ 35	☐ 36	☐ 37	☐ 38	☐ 39	
	☐ 40	☐ 41	☐ 42	☐ 43	☐ 44	☐ 45	☐ 46				
	The following dimensions are optionally available subject to surcharge:							☐ 47	☐ 48	☐ 49	☐ 50
Seat depth ST	Short frame (43 cm)				Long frame (47 cm)						
	☐ 36	☐ 37	☐ 38	☐ 39	☐ 40	☐ 41	☐ 42	☐ 43	☐ 44		
	☐ 40	☐ 41	☐ 42	☐ 43	☐ 45	☐ 46	☐ 47	☐ 48			
Back height RH	☐ 20	☐ 21	☐ 22	☐ 23	☐ 24	☐ 25	☐ 26	☐ 27	☐ 28	☐ 29	
	☐ 30	☐ 31	☐ 32	☐ 33	☐ 34	☐ 35	☐ 36	☐ 37	☐ 38	☐ 39	
	☐ 40	☐ 41	☐ 42	☐ 43	☐ 44	☐ 45	☐ 46	☐ 47	☐ 48		
Lower leg length USL	☐ 37	☐ 38	☐ 39	☐ 40	☐ 41	☐ 42	☐ 43				
	☐ 44	☐ 45	☐ 46	☐ 47	☐ 48	☐ 49	☐ 50				
Seat height back SHh (with 24" drive wheels)	☐ 38	☐ 39	☐ 40	☐ 41	☐ 42	☐ 43					
	The following dimensions are optionally available subject to surcharge:						☐ 44	☐ 45			
							☐ 46	☐ 47			
Seat height front SHv (with 5" caster wheels)	☐ 44	☐ 45	☐ 46	☐ 47	☐ 48	☐ 49	☐ 50	☐ 51	☐ 52		
Tipping point KP (inner edge of rear tube to centre of drive wheel axle)	☐ 3	☐ 4	☐ 5	☐ 6	☐ 7	☒ 8	☐ 9	☐ 10	☐ 11		

Please tick the required dimensions (all dimensions are in cm)

The underlined values are the standard settings. If not specified, the standard value will be delivered.

The SHh and SHv values might be a bit different from the specified values, depending on the wheel size and camber. Individual dimensions available on request.

We recommend a seat height difference (difference between SHv and SHh) between 5 cm and 12 cm.

Frame colour - One-colour design

Information on current colour designs and example images can be found at: <https://www.proactiv-gmbh.com/colours/>

Frame colour: Selection from colour palette

Frame colour will be specified in the order

☒ A9000100170-000

Raw, blasted, not coated

☒ A9000100170-020

Note: In combination with moisture e.g. hand perspiration, can affect dark aluminum discoloration.

Powder colours

white RAL 9016

☒ gloss A9000100170-023

☐ matt A9000100170-123 € 43.40

brilliant silver

☒ gloss A9000100170-011

anthracite RAL 7016

☒ gloss A9000100170-001

☐ matt A9000100170-101 € 43.40

black RAL 9005

☒ gloss A9000100170-003

☐ matt A9000100170-019 € 43.40

pure orange RAL 2004

☒ gloss A9000100170-022

☐ matt A9000100170-122 € 43.40

pink magenta RAL 4010

☒ gloss A9000100170-037

☐ matt A9000100170-137 € 43.40

pastel blue RAL 5024

☒ gloss A9000100170-036

☐ matt A9000100170-136 € 43.40

lemon green

☒ gloss A9000100170-005

☐ matt A9000100170-105 € 43.40

yellow RAL 1023

☒ gloss A9000100170-046

☐ matt A9000100170-146 € 43.40

red RAL 3001

☒ gloss A9000100170-047

☐ matt A9000100170-147 € 43.40

lagoon blue

☒ gloss A9000100170-052

☐ matt A9000100170-152 € 43.40

athletic grey RAL 7046

☒ gloss A9000100170-050

☐ matt A9000100170-150 € 43.40

muddy grey RAL 7033

☒ gloss A9000100170-051

☐ matt A9000100170-151 € 43.40

peach RAL 3022

☒ gloss A9000100170-038

☐ matt A9000100170-138 € 43.40

Powder colours „dormant“

deep black effect silver

☒ gloss A9000100170-007

sky blue

☒ gloss A9000100170-017

navy blue

☒ gloss A9000100170-045

☐ matt A9000100170-145 € 43.40

electric green

☒ gloss A9000100170-015

☐ matt A9000100170-115 € 43.40

violet

☒ gloss A9000100170-044

fire red

☒ gloss A9000100170-018

☐ matt A9000100170-118 € 43.40

Frame colour - One-colour design special colours „liquid coating“

€ 349.60

storm (ice blue)

☐ gloss A9000100170-500

☐ matt A9000100170-501

prairie (rose gold)

☐ gloss A9000100170-504

☐ matt A9000100170-505

berry

☐ gloss A9000100170-526

☐ matt A9000100170-527

blue ocean

☐ gloss A9000100170-522

☐ matt A9000100170-523

emerald

☐ gloss A9000100170-514

☐ matt A9000100170-515

chili

☐ gloss A9000100170-516

☐ matt A9000100170-517

Frame colour - Two-colour design

Note: In the two-colour design, both colours are the same gloss level as the design colour.

☐ **Two-colour design with one base colour (on the rear frame) and one design colour (on the front frame)**

€ 268.50

Please select a colour below:

Base colour (on the rear frame)

white RAL 9016

☐ A9000100170-200

brilliant silver

☐ A9000100170-201

anthracite RAL 7016

☐ A9000100170-202

black RAL 9005

☐ A9000100170-203

pastel blue RAL 5024

☐ A9000100170-206

lemon green

☐ A9000100170-207

yellow RAL 1023

☐ A9000100170-208

lagoon blue

☐ A9000100170-212

athletic grey RAL 7046

☐ A9000100170-210

muddy grey RAL 7033

☐ A9000100170-211

Design colour (on the front frame)

black RAL 9005

☐ gloss A9000100170-259

☐ matt A9000100170-260

lemon green

☐ gloss A9000100170-250

☐ matt A9000100170-251

chili

☐ gloss A9000100170-253

☐ matt A9000100170-254

pure orange RAL 2004

☐ gloss A9000100170-255

☐ matt A9000100170-256

emerald

☐ gloss A9000100170-274

☐ matt A9000100170-275

storm (ice blue)

☐ gloss A9000100170-264

☐ matt A9000100170-265

prairie (rose gold)

☐ gloss A9000100170-266

☐ matt A9000100170-267

berry

☐ gloss A9000100170-284

☐ matt A9000100170-285

blue ocean

☐ gloss A9000100170-280

☐ matt A9000100170-281

Colour - Attachment parts

Note: If not specified, the caster fork and bearing block will be mounted in black.

Drive wheel: aluminum rim, hub and spokes

☐ Black A9000100171-003

Caster fork and bearing block

☐ Silver A9000100173-011

☐ Black A9000100173-003

Colour design of attachment parts according to customer request:

Drive wheel: aluminum rim and design inlet¹ according to colour selection / hub and spokes black

☐ frame colour A9000100171-999 € 249.60

☐ Selection from one-colour design (different from the frame colour) A9000100171-998 € 288.90

Colour: _____

A9000100170- ☐ ☐ ☐

Caster wheel: 5" aluminum rim² (6-spoke design)

☐ frame colour A9000100174-999 € 75.20

☐ Selection from one-colour design (different from the frame colour) A9000100174-998 € 87.10

Colour: _____

A9000100170- ☐ ☐ ☐

Special colour request (Frame colour - One-colour design)

Frame colour for single purchasing Selection from RAL colour palette

☐ A9000100170-100 € 414.40

gloss level

☐ gloss

☐ matt

colour designation

RAL number

Note: Attachments **not** available in design colour options

¹ Only in combination with Liteline hub

² Only in combination with caster wheel 5" aluminum rim 6-spoke-design good shock absorption A9000500121*, as well as with suspension fork 5-spoke-design A9000500160*


Frame

Note: Footrest area with V-shape $x = SB - (0.25 \times USL + 6 \text{ cm})$

Straight front end

☐ A9017100012-K*
Knee angle 95°
Min. 4 cm seat slope and
max. caster wheel size 5"


☐ A9017100010-K*
Knee angle 105°



Front end V-shape

☐ A9017100020-K* € 257.00
Knee angle 95°
from SB 33 cm

☐ A9017100023-K* € 257.00
Knee angle 105°
from SB 33 cm



Seat width 47 cm - 52 cm

☐ A9016100060-K € 305.90
with frame reinforcement
rated for a load capacity of 120 kg

Back tubes & push handles

Equipment without push handles

☐ A9000401000



Back tube with integrated handles

☐ A9000401051-K € 130.40
straight.



Aluminium push handles fixed in back tube

☐ A9000401050* € 130.40
from RH 23 cm



Horizontal push handles screwed into the back tube

☐ A9000401076 € 130.40
can be unscrewed



Safety push handles with continuous height adjustment

☐ A9000401060* € 216.40
mount in the back tube
from RH 30 cm with adaptable back
from RH 35 cm with lordosis support
from RH 43 cm with ergonomic backshell



Safety push handles back-positioned length 45 cm

☐ A9000401075* € 372.90
continuous height adjustment,
Safety push handle holder mounted on the
back cross bar
not in combination with swing-away armrest



Additional option

Back tube with lordosis support

☐ A9000400005-K € 114.40

- ☐ 1 cm
- ☐ 2 cm
- ☐ 3 cm



Backrest

If not specified, the „folding backrest“ will be delivered.


Folding backrest

☐ A9000400021
7-way angle adjustable,
with integrated plastic-coated aluminum side
parts



Backrest folding function blocked for Ergo models

☐ A9000400036
Locking bolt for angle adjustment replaced by
screw connection.
Back angle 90° to the floor.



Fixed backrest

☐ A9000400020 € 99.10
screwed to aluminium side parts
Back angle 90° to the floor
Backrest cannot fold down

Back equipment

If not specified, the „Adaptable back with back padding“ will be delivered.

Adaptable back with back padding

☒ A9000802010-K001

Cover: Black, breathable, open surface
Equipment with back cross bar



Adaptable back with back padding „light“

☒ A9000802320-K001

Cover: Black, breathable, open surface
Cushion thickness: 9 mm
Equipment with back cross bar

☒ A9000802320-K002

Cover: Black, breathable, open surface
Cushion thickness: 9 mm
Equipment without back cross bar



Adaptable back with back padding „Comfort“

☒ A9000815101-K € 131.40

Cover: Black, breathable, open surface
Front padding: 20 mm
Rear padding: 5 mm
Equipment with back cross bar

Back rise + 3 - 5 cm

☒ A9000802100 € 159.10

for the nominal back height dimension
Only in combination with height-adjustable push handles



Equipment without back upholstery and back padding or back system

☒ A9000802011 € - 53.00

Equipment with back cross bar


Seat upholstery

If not specified, the „Body Contour“ seat upholstery will be delivered.

Body Contour


☒ A9000801025-K

approved up to a body weight of 120 kg



Open belt system

☒ A9000801010-K € 79.90





Equipment without seat upholstery

☒ A9000801050



Seat cushion



Seat cushion in rectangular contour Inlet: Standard

☒ A9000805180-K* € 118.30

Seat cushion with radius contour Inlet: Standard

☒ A9000800300-K* € 163.10

Seat cushion height:

☐ 3 cm ☐ 5 cm
☐ 4 cm ☐ 6 cm

Seat cushion in rectangular contour Inlet: Anti-Decubitus

☒ up to 30x30 cm A9000805204-K* € 301.70

☒ >30x30 up to 50x50 cm A9000805200-K* € 437.70

Seat cushion with radius contour Inlet: Anti-Decubitus

☒ up to 30x30 cm A9000800311-K* € 356.10

☒ >30x30 up to 50x50 cm A9000800310-K* € 482.90

Seat cushion height:

☐ 4 cm ☐ 6 cm

Cover:

☐ Black, breathable, open surface

☐ Black, breathable, closed surface

If not specified, the cover material „Black, breathable, open surface“ will be delivered.

☐ Seat cushion is provided

Clothing guard

Aluminum clothing guard removable via lock function

Adjustable following wheel arch
Only in combination with folding backrest



in anthracite matt

☐ A9000403215-K € 243.20

in frame colour (base colour for two-colour design)

☐ A9000403216-K € 311.10

in design colour (only possible with frame colour in two-colour design)

☐ A9000403217-K € 345.10

Carbon clothing guard removable via lock function

☐ A9000400013-K € 547.30

Adjustable following wheel arch
Only in combination with folding backrest



Armrests

Armrest, height adjustable, integrated in the clothing guard

Including aluminum clothing guard
with arm padding 25 cm long
Removable via lock function
Only in combination with folding backrest



in anthracite matt

☐ A9000400014-K € 569.60

in frame colour (base colour for two-colour design)

☐ A9000400015-K € 621.90

in design colour (only possible with frame colour in two-colour design)

☐ A9000400016-K € 648.10

Armrest, height adjustable, integrated in the clothing guard

☐ A9000400017-K € 862.00

Including carbon clothing guard
with arm padding 25 cm long
Removable via lock function
Only in combination with folding backrest

Armrests folding up, removable

☐ A9000403141-K € 358.50

with arm padding 25 cm long
Positioning: upper arm padding edge approx.
25 cm to 31 cm above the seat tube,
adjustable in 1.5 cm steps
Not possible in combination with back-
positioned safety push handles




Axle plate

Standard equipment with 12 mm sockets / axles


Standard axle plate

☐ A9016201000
Rear seat height approx. 38 cm - 43 cm
(for 24" wheel)



Long axle plate

☐ A9016201020 € 234.70
Rear seat height approx. 43 cm - 47 cm
(for 24" wheel)



Wheel camber

*Total width with tight-mounted handrim.

1° drive wheel camber

☐ A9016202010

Total width = approx. SB + 20 cm*

2,5° drive wheel camber

☐ A9016202051

Total width = approx. SB + 23 cm*

4° drive wheel camber

☐ A9016202030

Total width = approx. SB + 26 cm*

6° drive wheel camber

☐ A9016202040

Total width = approx. SB + 30 cm*

E-Adaptations

Preparation for the „e-motion“ / „e-motion DuoDrive“ mounting

☐ A9016201200-K € 502.90

Equipment: e-motion torque arm, 1/2" axle and socket length 50.8 mm (2")
PRO ACTIV and Alber recommend mounting two anti-tipping supports. Please select under „Adaptations“.

Preparation for the „e-fix“ mounting

☐ A9016201210-K € 502.90

Equipment: 1/2" axle, socket length 50.8 mm (2") and e-fix torque support.
PRO ACTIV and Alber recommend mounting two anti-tipping supports. Please select under „Adaptations“.
Please order Alber anti-tipping supports and corresponding components direct from Alber.

Brake system

If not specified, the „Brake system push to lock with standard brake handle“ will be delivered.

Recommendation: For KP 7 cm or lower, we recommend the „Brake system push to lock with brake handle 20° outwards“.

Integral parking brake

☐ A9035705000 € 248.50

Brake lever fixed under the upper frame tube
Only in combination with 4° drive wheel camber



Knee lever brake with standard brake lever

☐ A9000701010-001



Knee lever brake with long brake lever

☐ A9000701010-002



Knee lever brake with low brake lever

☐ A9000701010-003



Knee lever brake system Brake lever angled 20°

☐ standard brake handle A9000700002

☐ long brake handle A9000701010-004



Brake system push to lock Brake handle to flip-down

☐ A9000706008 € 158.60
length: 9 cm



Drum brake

☐ A9000702010-K* € 474.00

Operating from the push handle
Total wheelchair width (see selection „Wheel camber“) + 3 cm

Depending on the seat height, the long wheel plate is mounted in order to protect the drum brake anchors.

Please note:

- Only in combination with drive wheel aluminum 36-spoke rim wheel
- Not in combination with one-hand drive, e-fix or e-motion



Footrest

If not specified, the „continuous footrest“ will be delivered.

Continuous footrest

☒ A9016601010-K*
with tubular resting bracket
height adjustable +/- 3 cm
for KW 95° from USL 32 cm
for KW 105° from USL 37 cm
If you want to equip the footrest with a „rubber coated footplate“ please select this under „Footrest accessories“.



Footrest folding up to one side

☐ A9016601300-K* € 173.90
with tubular resting bracket
height adjustable +/- 3 cm
can be swung away,
cannot be removed
for KW 95° from USL 35 cm
for KW 105° from USL 40 cm
V-shape from SB 40 cm
If you want to equip the footrest with a „rubber coated footplate“ please select this under „Footrest accessories“.



Footrest folds up to the rear with spring locking mechanism

☐ A9016602005-K* € 287.40
with tubular resting bracket
for KW 95° from USL 36 cm
for KW 105° from USL 37 cm
If you want to equip the footrest with a „rubber coated footplate“ please select this under „Footrest accessories“.



Footrest folds up to the rear

☐ A9016602010-K* € 307.90
with tubular resting bracket and 20 mm operating lever
for KW 95° from USL 38 cm
for KW 105° from USL 42 cm
straight front end from SB 33 cm
V-shape from SB 38 cm
Support width from 20 cm
If you want to equip the footrest with a „rubber coated footplate“ please select this under „Footrest accessories“.



Divided footrest

☐ A9016603200-K € 336.80
with rubber coated footplate,
height adjustable +/- 3 cm, angle adjustment,
can be pivoted outwards and folded up to the side
Including frame reinforcement
for KW 95° from USL 37 cm
for KW 105° from USL 42 cm
V-shape from SB 40 cm



„Swing Away“ footrest

☐ A9016604100-K € 404.70
with rubber coated footplate,
separated in the middle,
can be folded to the side, removable,
can be folded up, height-adjustable,
angle-adjustable
Including frame reinforcement
from USL 34 cm



Footrest accessories

Rubber coated footplate

☐ A9000610000-K € 125.90
bolted, angle-adjustable
Depth: 14 cm



Footplate side stop

☐ A9000610100 € 59.70
adjustable
only in combination with rubber coated footplate



Adaptations

Anti-tipping support

☐ right side A9000300100-K001 € 169.10
☐ left side A9000300100-K002 € 169.10
with tension spring and universal adapter pivotable
For optimum pivoting behaviour the seat slope should be > 3 cm



Preparation for mounting the anti-tipping support with lifting function from Alber (e-fix, e-motion)

☐ A9016300105 € 296.70
Universal adapter is prepared for mounting the anti-tipping support with lifting function from Alber.
Supplied without anti-tipping support and anti-tipping support mount.
On the customer side, please order this from Alber



Transit wheels

☐ right side A9000301100-003 € 164.80
☐ left side A9000301100-004 € 164.80
with universal adapter, removable
not in combination with anti-tipping support,
right side or left side



Tilting clamp

☐ right side
A9000302100-K001 € 140.10

☐ left side
A9000302100-K002 € 140.10

with universal adapter and Quick pin
not in combination with anti-tipping support,
right side or left side



Removable wheelbase extension SPEEDY 4all / 4you

☐ A9016102001 € 581.30

wheel position displacement approx. 14 cm
adjustable track

Please select the wheel camber:

☐ 1° drive wheel camber

☐ 4° drive wheel camber

☐ 6° drive wheel camber

Note: The camber of the wheelbase extension
cannot be less than in the everyday position.
A max. camber of 6° is possible.

If not specified, the next higher wheel camber
will be mounted than in the everyday position.



3" caster wheel

3" skate wheel (80 mm)

☐ A9000501110 € 70.80

puncture-proof



4" caster wheel

4" caster wheel „Standard“ solid rubber black

☒ A9000501011-004
(82A), puncture-proof



4" skate wheel (100 mm)

☐ A9000501120 € 70.80

puncture-proof



4" safety light-up wheel (100 mm)

☐ A9000501220 € 94.80

puncture-proof



4" CorMORun caster wheel (100 mm)

☐ A9000501042 € 173.00

with replaceable rubber 6 mm,
puncture-proof
Rim colour: Black
Incl. hair guard spacer (5.3 mm)



4" Frog Legs caster wheel (100 mm)

☐ A9000501230 € 206.30

with replaceable rubber 6 mm,
puncture-proof
Colour: black, silver rim
Incl. hair guard spacer (5.3 mm)



5" caster wheel

5" caster wheel „Standard“ PU black

☒ A9000501036-005
(82A), puncture-proof



5" caster wheel „Standard“ PU grey

☒ A9000501031-005
(82A), puncture-proof



5" safety light-up wheel (125 mm)

☐ A9000501221 € 94.80
puncture-proof



5" aluminum rim 3-spoke design „good shock absorption“

☐ silver anodised A9000501210 € 179.80

☐ black A9000501211 € 179.80

with replaceable rubber 6 mm
medium hard, puncture-proof
Incl. hair guard spacer (3 mm)



5" CorMORun caster wheel (125 mm)

☐ A9000501041 € 181.30

with replaceable rubber 6 mm,
puncture-proof
Rim colour: Black
Incl. hair guard spacer (5.3 mm)



5" Frog Legs caster wheel (125 mm)

☐ A9000501232 € 206.30

with replaceable rubber 6 mm,
puncture-proof
Colour: black, silver rim
Incl. hair guard spacer (5.3 mm)



6" caster wheel

6" caster wheel „Standard“ PU grey

☒ A9000501010-006
(82A), puncture-proof



6" caster wheel „pneumatic tyre“

☐ A9000501050 € 57.20



6" caster wheel „Flexel-System“

☐ A9000501070 € 57.20
puncture-proof



Caster fork

For caster fork colours, see „Colour of add-on parts“

If not specified, the „Standard caster fork with screw axle“ will be delivered.

Caster fork Standard

☒ A9000500100-K001
with screwed axle

☐ A9000500100-K002 € 167.40
with quick-release axle



Caster fork Heavy Duty for high loads

Extra weight approx. 100 - 150 gr.

☐ A9000510111-K001 € 54.00
with screwed axle

☐ A9000510111-K002 € 218.90
with quick-release axle



Suspension fork

inkl. 5" aluminum rim „good shock absorption“

☐ A9000500160-K002 € 616.90
with quick-release axle



Caster fork spare wheel „Standard“

☐ A9000500300-K € 296.70
Caster fork Standard with quick-release axles



Caster fork spare wheel „Heavy Duty“


☐ A9000500301-K € 337.90
Caster fork Heavy Duty with quick-release axles



The spare caster fork does not have a caster wheel.
In case a caster wheel is requested, please put a cross on a second selection in the category „Caster wheels“.

Track distance

Caster wheel track extension



☐ A9000500400*
10 mm per side

€ 59.10

☐ A9000500401*
15 mm per side

€ 59.10

☐ A9000500410*
20 mm per side

€ 59.10

Drive wheels

If not specified, the „drive wheel aluminum rim 36 spokes“ will be delivered.

Drive wheel 18 spoke aluminum rim

☐ 24" (ETRTO 540)
A9000203122
with Liteline hub, radially spoked
Not in combination with:
- Puncture-proof tires
- One-hand drive
- Drum brake
- Transport restraint connection
- Handbike
up to 70 kg body weight:

Sturz	1°	4°	6°	8°	10°
Gewicht					
bis 50 Kg	✓	✓	✓	✓	✓
51 - 60 Kg	✓	✓	✓	✓	
61 - 70 Kg	✓	✓			



Drive wheel 36 spoke aluminum rim

☐ 24" (ETRTO 540)
A9000203000-010
with Liteline hub (left image),
radially spoked
up to 120 kg body weight
☐ 24" (ETRTO 540)
A9000203000-020
with aluminum front hub (right image),
radially spoked
over 120 kg body weight



Drive wheel double rim 36-spokes

☐ 24" (ETRTO 540)
A9000203100 € 157.80
☐ 25" (ETRTO 559)
A9000203200 € 157.80
☐ 26" (ETRTO 590)
A9000203300 € 157.80
with aluminum front hub,
radially spoked
For use with high loads and body weight
>120 kg



Spinergy wheel „Spox Sport“ with 24 PBO spokes; rim, hub and spokes in black

☐ 24" (ETRTO 540)
A9000203550 € 784.10
☐ 25" (ETRTO 559)
A9000201559 € 784.10
☐ 26" (ETRTO 590)
A9000203555 € 784.10



Spinergy „Everyday Light Extreme“ with 12 PBO spokes; rim and hub in black, spokes in white

☐ 24" (ETRTO 540)
A9000203554 € 957.20
Recommendation: For heavy loads
(e.g. higher body weight), we recommend the
Spinergy „Spox Sport“ wheel with 24 PBO
spokes.




Spinergy „Everyday Light Extreme“ with 12 PBO spokes; rim, hub and spokes in black

☐ 24" (ETRTO 540)
A9000203552 € 957.20
Recommendation: For heavy loads
(e.g. higher body weight), we recommend the
Spinergy „Spox Sport“ wheel with 24 PBO
spokes.



Handrims

x = dimension handrim attachment



Handrim attachment
general:

- ☐ x = 11 mm
☐ x = 19 mm

Handrim attachment
Curve + Gekko handrim:

- ☐ x = 10 mm
☐ x = 18 mm

If not specified, the handrim is mounted wide (19 mm / 18 mm).

Stainless steel Ø 20 mm

- ☐ 24" handrim for
24" wheel and 25" wheel
A9000205010-004
- ☐ 26" handrim for
26" wheel
A9000205010-006



Aluminum Ø 19 mm silver anodised

- ☐ 24" handrim for
24" wheel and 25" wheel
A9000205019-004



Aluminum Ø 19 mm powder coated in frame colour

(basic colour with two-colour design)

- ☐ = Gloss surface
☐ = Matt surface

Note: Matt coating is not recommended due to the low grip.

- ☐ 24" handrim for
24" wheel and 25" wheel € 68.30
A9000205020-004
- ☐ 26" handrim for
26" wheel € 68.30
A9000205020-006



Aluminum Ø 19 mm polished and hard-anodised

- ☐ 24" handrim for
24" wheel and 25" wheel € 164.90
A9000205030-004
- ☐ 26" handrim for
26" wheel € 164.90
A9000205030-006




Titanium Ø 19 mm

- ☐ 24" handrim for
24" wheel and 25" wheel € 335.60
A9000205040-004
- ☐ 26" handrim for
26" wheel € 335.60
A9000205040-006



Silicone handrim cover

- ☐ for handrim 24"
A9000205100-004 € 232.80
- ☐ for handrim 26"
A9000205100-006 € 232.80
- Surface without structure,
compatible with handrim Ø 19 mm or
Ø 20 mm.



Ultra-Grip „very grippy“

- ☐ 24" handrim for
24" wheel and 25" wheel € 402.40
A9000205070-004
- ☐ 26" handrim for
26" wheel € 402.40
A9000205070-006



Ergo-Grip

- ☐ 24" handrim for
24" wheel and 25" wheel € 402.40
A9000205080-004
- ☐ 26" handrim for
26" wheel € 402.40
A9000205080-006
- Attachment fitting to the rim

Note: Due to the complex geometry,
concentricity deviations of +/- 2-3 mm are
possible.



Curve black anodised, ergonomic handrim profile

- ☐ 24" handrim for
24" wheel and 25" wheel
Profile height 24 mm € 402.40
A9000205130-004
- ☐ 24" handrim for
24" wheel and 25" wheel
Profile height 31 mm € 402.40
Size L
A9000205140-004
- ☐ 26" handrim for
26" wheel
Profile height 31 mm € 402.40
Size L
A9000205140-006



Curve silver anodised ergonomic handrim profile ☐ 24" handrim for 24" wheel and 25" wheel Profile height 31 mm € 402.40 Size L A9000205160-004 	Curve powder-coated, ergonomic handrim profile ☐ = frame colour (gloss surface) ☐ = black (gloss surface) ☐ 24" handrim for 24" wheel and 25" wheel Profile height 24 mm € 467.60 A9000205135-004 ☐ 24" handrim for 24" wheel and 25" wheel Profile height 31 mm € 467.60 Size L A9000205145-004 ☐ 26" handrim for 26" wheel Profile height 31 mm € 467.60 Size L A9000205145-006 	Curve Grip, with Tetra-Grip coating, ergonomic handrim profile ☐ 24" handrim for 24" wheel and 25" wheel Profile height 31 mm € 520.70 Size L A9000205150-004 
Quadro with Tetra-Grip coating ☐ 24" handrim for 24" wheel and 25" wheel € 520.70 A9000205170-004 ☐ 25" handrim for 25" wheel € 520.70 A9000205170-005 ☐ 26" handrim for 26" wheel € 520.70 A9000205170-006 	Gekko black anodised ☐ 24" handrim for 24" wheel and 25" wheel € 520.70 A9000205175-004 ☐ 25" handrim for 25" wheel € 520.70 A9000205175-005 ☐ 26" handrim for 26" wheel € 520.70 A9000205175-006 with a robust and replaceable silicone lip, fully integrated into the profile Profile height 31 mm 	Gekko silver anodised ☐ 24" handrim for 24" wheel and 25" wheel € 520.70 A9000205185-004 with a robust and replaceable silicone lip, fully integrated into the profile Profile height 31 mm 
Gekko S black anodised ☐ 24" handrim for 24" wheel and 25" wheel € 520.70 A9000205195-004 with a robust and replaceable silicone lip, fully integrated into the profile Profile height 24 mm 	Equipment without handrims ☒ A9000205990	

Spoke guard

Finger- and spoke guard „transparent“

☐ for 24" wheel
A9000206300-004 € 207.10

☐ for 26" wheel
A9000206300-006 € 207.10

not in combination with Spinergy drive wheels.



Tyres

If not specified, the „Schwalbe RightRun in black“ tyres will be delivered.

Schwalbe „RightRun“ High-pressure tyres

☐ 24x1" (540-25)
A9000204001-004
Colour: Black
7.5 bar



Schwalbe „RightRun“ High-pressure tyres

☐ 24x1" (540-25)
A9000204000-004

☐ 25x1" (559-25)
A9000204000-005

☐ 26x1" (590-25)
A9000204000-006
Colour: Grey, black wall
7.5 bar



Schwalbe „DownTown“ Profiled tyres

☐ 24x1" (540-25)
A9000204010-004
7.5 bar



Schwalbe „Marathon Plus Evolution“ with „smart guard“ puncture protection

☐ 24x1" (540-25)
A9000204030-014 € 49.20

☐ 25x1" (559-25)
A9000204030-015 € 49.20

☐ 26x1" (590-25)
A9000204030-016 € 49.20

Colour: Black
10 bar



Schwalbe „Marathon Plus“ with „smart guard“ puncture protection

☐ 24x1" (540-25)
A9000204030-004 € 49.20

Colour: Grey, black wall
10 bar



Schwalbe „One“

☐ 24x1" (540-25)
A9000204123-004 € 93.00
Colour: Black
9 bar



Puncture-proof „solid rubber tyres“

☐ 24x1" (540-25)
A9000204201-004 € 122.70

Colour: Black
not in combination with Spinergy drive wheels.



Accessories

24" MTB spare wheel set

with double rim (ETRTO 507), aluminum front hub, quick-release axle and stainless steel handrim
Mountain bike tyres, selectively:

☐ A9000203410 € 703.40

24x1.75" Schwalbe "Marathon Plus"
Brake setting compatible with 24x1" wheel



☐ A9000203413 € 703.40

24x2.1" Schwalbe "SMART SAM Sport"
Brake setting not compatible with 24x1" wheel




NOTE: We charge a surcharge for the use of drum brakes

☐ A9000200004 € 232.20

Through axle with assistance, Large push button

☐ A9000200096 € 105.80

Interlock dimension: 122 / 110 mm
Axle diameter: 12 mm
Not in combination with one-hand drive



Hand air compressor with pressure gauge

☐ A9000901100 € 106.30

up to 14 bar



Crutch holder with flexible nylon cup

☐ A9000902310* € 137.40



Length-adjustable safety belt with buckle

☐ A9000902100 € 80.60



Adjustable safety belt with hook and fastener

☐ A9000902110 € 64.50



SIGMA lighting set

☐ A9000903270 € 107.90

Consisting of front light and rear lights, chargeable via Micro USB socket, incl. Micro USB charging cable
Lighting set compliant with road traffic regulations



Transport bag for drive wheels

☐ for 24" wheel
A9000903530-004 € 77.80
Ø 76 cm x 2 cm

☐ for 26" wheel
A9000903530-006 € 77.80
Ø 83 cm x 2 cm



Separate order sheets

Category „Rehabilitation technology accessories“

e.g. Seat & backrest systems, one-hand drive or transport restraint systems can be found at:

www.proactiv-gmbh.com/downloads/

Additional products

such as hand bikes or traction & attachment devices can be found at:

www.proactiv-gmbh.com/downloads/

Additional sheets

such as heavy-duty versions etc. can be found at:

www.proactiv-gmbh.com/downloads/

Services

Consultation, clarification of the installation situation and measurement through PRO ACTIV

☐ A9990000000 € 354.50

Handover and instruction on the wheelchair by PRO ACTIV

☐ A9990000006 € 354.50

Notes and remarks:

Location, Date

Signature End customer/dealer

The signature field is intended for the consultant and insured person. This confirms that the intended equipment is intended for the insured person.
- 17 - 2026-03-03 / The prices are plus statutory value-added tax and shipping costs ☐ = included in the basic price ☒ = Position selected

Order sheet SPEEDY 4all Ergo Rev. 04 UK