



NJ1 e-assistant

up to 120 kg trailer load and 10 kg payload
4 different basic models for different purposes,
equipped with:

- Docking plate for quick connection to the „Quick-Snap Adapter“
- A-shaped handles with integrated switching and brake fittings
- hydraulic disk brake system, right side operation
- additional brake, operation on the handle on the left with locking bar for use as parking brake
- chain protection tube
- starting assistant (throttle) for travelling without moving the crank up to 6 km/h
- Parking stand
- Lighting set
- Mudguard
- Colour display sMMi Z15 Connect
- Drinking bottle with holder
- Drive unit black

Drive Neodrives

Torque 12/40 Nm, power 250 Watt

Wheel size:

- | | |
|---|------------|
| <input type="checkbox"/> 20" (ETRTO 406)
A9116100011-020 | € 6,880.00 |
| <input type="checkbox"/> 24" (ETRTO 507)
A9116100011-024 | € 6,880.00 |

>> Expected to be available from April 2026 <<

Quotation date: _____

Consultation was carried out by
Staff at the medical supplies shop: _____

Enquiry no.: _____

Order date: _____

according to GTC PRO ACTIV Reha-Technik GmbH

Order no.: _____

Medical suppliers / company address:

Cust. no.: _____

Name: _____

Street: _____

Postcode / city: _____

Phone: _____

Fax: _____

E-Mail: _____

Contact: _____

Different invoicing address (see below)

Cust. no.: _____

Name: _____

Street: _____

Postcode / city: _____

Different shipping address (see below)

Cust. no.: _____

Name: _____

Street: _____

Postcode / city: _____

Customer name / commission:

ID / Comm. No.: _____

Surname, forename: _____

Street: _____

Postcode / city: _____

Phone: _____

E-Mail: _____

D.o.B.: _____

Paying organisation: _____

Indication: _____

Batteries Neodrives

The specified ranges depend on the distance profile, user weight and speed.

Battery Neodrives DT 36V 16Ah 576Wh

Range approx. 45 km - 75 km

- A9116900116
mounted on the left of the fork



Double Battery Neodrives DT 36V 16Ah 576Wh

Range approx. 85 km - 150 km

- A9116900118 € 1,290.40

Montage beidseitig an der Gabel,
Umsteckvorrichtung an der Tretlagerstütze
montiert zum einfachen Wechsel zwischen den
Akkus über den Magnetstecker



Derailleur gears

Shimano Deore 1x11 Rapid Fire 11-speed derailleur

- A9110200032
mechanical derailleur gear with operation on
the handle. 11-speed cassette 11-42 teeth

Selection of the chainring:

- 32 teeth recommended for
24" drive wheel
- 38 teeth recommended for
20" drive wheel
- 42 teeth for lower crank
frequency or driving
at higher speeds



Shimano Deore XT 1x11 Rapid Fire 11-speed derailleur

- A9110200034 € 192.20

mechanical derailleur gear with control on the
handle. 11-speed cassette 11-42 teeth
top chain ring with chain guards

Selection of the chainring:

- 32 teeth recommended for
24" drive wheel
- 38 teeth recommended for
20" drive wheel
- 42 teeth for lower crank
frequency or driving
at higher speeds



Shimano Deore XT 1x11 DI2 with push button operation 11-speed derailleur

- A9110200035 € 1,486.10

electronic derailleur gear with operation on the
handle. 11-speed cassette 11-42 teeth

Operation: right side left side

If not specified, the operation will be mounted
on the „right side“.

Selection of the chainring:

- 32 teeth recommended for
24" drive wheel
- 38 teeth recommended for
20" drive wheel
- 42 teeth for lower crank
frequency or driving
at higher speeds



Shimano Deore XT 1x11 DI2 with chin control 11-speed derailleur

- A9110200033 € 1,631.10

electronic derailleur gear with chin control
(tetra equipment).

11-speed cassette 11-42 teeth

Selection of the chainring:

- 32 teeth recommended for
24" drive wheel
- 38 teeth recommended for
20" drive wheel
- 42 teeth for lower crank
frequency or driving
at higher speeds



Gear change accessories

Mountain drive 2.5 : 1 Hill reduction

A9100200003 € 797.40

Planetary gearbox for doubling the gears.
Allows you to shift downhill while standing still.
As standard with 38 teeth chainring.
In combination with grip width 40 cm, 45 cm,
48 cm or 53 cm



Brake system

Mechanical parking disk brake can be operated via the spreader

A9110300002 € 84.60

mounted on the bottom bracket support,
brake handle on the handle, left side is omitted



Driving brake - supplement for the mechanical parking brake

A9117300120 € 116.10

The parking disc brake is supplemented with a
brake handle on the left. Thus, the mechanical
disc brake can be used as a driving brake and
parking brake.


Powerful braking system with brake disc Ø 203 mm

A9110300105 € 83.10

Recommended from 80 kg trailer load or body
weight.

Not required if the hydraulic coaster disc brake
is selected

Not in combination with a Rohloff gear system



Double powerful brake system brake disc Ø 203 mm

for 20"
A9110300100-020 € 172.50

locking function can be activated using the
locking bracket on the handle.

for 24"
A9110300100-024 € 347.10

can be operated as a parking brake together
with V-brake rim brake using the locking lever.
Recommended from 100 kg trailer load

Mechanical backpedal disk brake with crank release function

A9110300004 € 738.40

compatible with derailleur and free-running
shifting hub.

In combination with:

- Grip width 45 cm, 50 cm, 53 cm or 58 cm
 - Pedal crank crank lengths 155 mm or 175 mm
up to 120 kg body weight
- Not in combination with Mountaindrive 2.5 : 1
gear reduction for hills

Tyres

Schwalbe „Marathon Plus“

20x1,75" (406-47)
A9130105128

24x1,75" (507-47) € 28.70
A9110105135

unpuncturable, with anti-puncture belt
Smart Guard



Dimensions

Crank type and grip width

Standard crank

- grip width 40 cm
A9110500100-040
- grip width 45 cm
A9110500100-045
- grip width 50 cm
A9110500100-050

Tetra crank (flared)

with extended crank arms for maximum upper body space

- grip width 48 cm (176 mm)
A9110500100-048
- grip width 53 cm (226 mm)
A9110500100-053
- grip width 58 cm (276 mm)
A9110500100-058

Standard pedal crank

- length 155 mm
A9110600100-155
- length 175 mm
A9110600100-175
- length 195 mm
A9110600100-195




Tetra pedal crank

- length 155 mm
A9100600105-155 € 136.10
- length 175 mm
A9100600105-175 € 136.10
- individual length
195 mm to 140 mm € 258.70
A9100600100



Fixed bottom bracket support

- L (X/Y 158/362 mm)
A9110501100-002
- XL (X/Y 158/410 mm) € 66.70
A9110501100-006
- XXL (X/Y 170/465 mm) € 66.70
A9110501100-007



Adjustable bottom bracket support

the angle and length can be adjusted via the clamped rotary point on the fork bridge and length.

For the pre-setting, please specify the bottom bracket support size that comes as close to your desired position.

- M (300 mm) € 232.10
A9110500001-001
- L (370 mm) € 232.10
A9110500001-002
- XL (420 mm) € 232.10
A9110500001-003

If not specified, the adjustable bottom bracket support will be delivered in XL. Available up to size XL.



Handles

Palm rest on the handle

- A9110601200 € 102.80
upholstered, bolted, angle-adjustable.



Tetra hand support, pivotable, adjustable grip distance

- A9110601101 € 169.40



Tetra backrest swivel and adjustable grip width

- A9110601102 € 191.60
very stiff external support as round tube
Ø 20 mm with foam rubber coating
Not in combination with Tetra hand support
pivotable



Tetra ergonomic handle, freestanding

- A9100600001 € 164.60
crank handles arranged vertically
can be combined with all other handle options



Tetra shell handles

- A9110601100 € 172.40



Accessories Bike

Black aluminum mudguard

- A9116500005 € 112.40
Only in combination with 20" drive wheel

Schutzblech und Kettenschutzscheibe individuelle Farbauswahl

- A9116500010 € 248.60
from colour palette one-colour design
powder colours wheelchair
A9000100170- _ _ _ _
Colour: _____
Only in combination with 20" drive wheel

LED lighting set for Neodrives System

- A9116900002 € 151.70
LED front headlight powered by the drive
system battery
Rear light as plug-in light with battery power
supply
Lighting set from the basic equipment is not
included



Bracket Ø 22 mm x 42 mm for mounting accessories (tacho, lighting, etc.)

- A9110903271 pcs. á € 23.20
Colour: natural anodized
_____ piece



Maneuvering castors for parking stand

- A9110900017 € 122.10
Maneuvering castors mounted on the parking
stand via inlet



Maneuvering handle mounted on the chain guard

- A9110500000 € 64.00
with foam rubber coating

Aluminum luggage rack „cargo systems”

A9110903290 € 188.80

Colour: black

max. load weight on the loading platform 10 kg, max. 15 kg in conjunction with luggage bags at lower loading height.



Balancing weight for battery

A9116900000 € 128.30

with drinking bottle and holder.

For compensating the weight of the rechargeable battery when using the NJ1 with only one rechargeable battery and for improving traction.



Operating satellite on the handle

Mounting: Right side Left side

If not specified, the mounting is on the left side.

with existing support line (laid along a brake line)

A9116200000 € 103.10

without existing support line (an additional blind line is required)

A9116200001 € 146.20

Options for users with reduced functions

Lever extension excenter lock for Quick Snap connection

Makes it easier to open and close the Quick Snap adapter. Tetra version

A9110700005 € 73.90

straight design

only in combination with rigid frame adapter

A9110700006 € 73.90

curved design

only in combination with front adapter



Adaptation help for centering pin

A9110700002 € 55.00

to simplify the insertion of the central adapter.



Parking stand with operating lever for high operating comfort

A9110900004 € 58.70

with longer operating lever

Chest strap for upper body fixation

size S (105 cm) € 241.80
A9000800016

size M (125 cm) € 241.80
A9000800013

size L (145 cm) € 241.80
A9000800014

size XL (165 cm) € 241.80
A9000800017

consisting of elastic webbing material in combination with neoprene with Velcro and operating tab.

Colour: black



Frame adapter

SPEEDY rigid frame adapter with „Quick-Snap-Adaption“

A9110700141 € 537.20

with quick-release clamp and centering pin holder
Colour: black

Please select accordingly:

Centring pin for Ø 25 mm
for A1 and A2

A9110700022

Centring pin for Ø 40 mm
for 4all, 4all Ergo, 4you Ergo, F2, F3 and 4teen

A9110700023



TRAVELER folding frame adapter with „Quick-Snap-Adaption“

A9110701300 € 537.20

with quick-release clamp and centering pin holder on the folding mechanism
Colour: black



BUDDY rigid frame adapter with „Quick-Snap-Adaption“

A9110701400 € 537.20

with quick-release clamp and centering pin mount on the frame
Colour: black



Wheelchair front adapter with universal clamp and „Quick-Snap-Adaption“

A9110701500 € 572.30

Adapter attachment to the wheelchair stem.
Operation on the right-hand side.

Mounting situation for PA products:
Straight front end from SB 31 cm to 54 cm
V-shaped front end from SB 30 cm to 50 cm

Please select the frame clamp below:

clamping range 95 mm
A9110701525

Tube diameter 23-27 mm

clamping range 95 mm
A9110701530

Tube diameter 28-32 mm

clamping range 95 mm
A9110701535

Tube diameter 33-37 mm

PLEASE NOTE:

If the exact tube diameter is not specified, an adapter for a tube diameter of 23 - 27 mm will be delivered as standard.

Please select the clamping sleeve below:

clamping sleeve [80 mm]
included as standard

for close clearance between footrest and drive wheel.

clamping sleeve [110 mm] € 46.70
A9199000117

for medium distance between footrest and drive wheel.

clamping sleeve [170 mm] € 56.60
A9199000154

for a wide distance between footrest and drive wheel.

PLEASE NOTE:

If not specified, the clamping sleeve with 80 mm length will be delivered.

Clamp lever with 3-way star handle for actuating the frame clamping

A9199701507 € 106.10


mounted on the right and left for easy terminal operation with limited manual function
Only in combination with front adapter



Wheelchair-dependent accessories

Extended bike plate

- for SPEEDY A1 / A2
A9110100100-002 € 370.00
- for BUDDY
A9110100100-003 € 370.00
- for TRAVELER
A9110100100-001 € 370.00
- 4° camber integrated
Track independently adjustable



Wheelbase extension via „grid mechanism with pulling operation“

- for SPEEDY 4all /
4teen / 4you /
TRAVELER 4you
A9110102001 € 564.70
- removable, with continuous axle tube, track
width independently adjustable
- Please select the wheel camber:
- 1° drive wheel camber
- 4° drive wheel camber
- 6° drive wheel camber



Wheelbase extension angle-adjustable and removable

- for SPEEDY F2 / F3
A9015100005 € 746.80
- with continuous axle tube,
height adjustment via gear joint
adjustable in 6° angular increments
wheel position displacement approx. 14 cm
adjustable track
- Attachment possibility must be checked by
PRO ACTIV.
- Only in combination with base plate.
Up to 120 kg load capacity.
- Please select the wheel camber below:
- 1° drive wheel camber
- 4° drive wheel camber
- 6° drive wheel camber



Services

Consultation, clarification of the installation situation and measurement through PRO ACTIV

- A9990000000 € 354.50

Handover and instruction of the handbike / traction device by PRO ACTIV

- A9990000011 € 354.50

Attachment of the handbike / traction device to the wheelchair by PRO ACTIV at specialist dealer

on a PRO ACTIV wheelchair

- A9990000012 € 230.00

product no.: _____

on a third-party product

- A9990000013 € 366.90

manufacturer: _____

model: _____

frame tube Ø ____ mm

Attachment of the handbike / traction device to the wheelchair through specialist dealer

including mounting of the frame adapter

on a PRO ACTIV wheelchair

- A9990000008

product no.: _____

on a third-party product

- A9990000005

manufacturer: _____

model: _____

frame tube Ø ____ mm

Notes and remarks:

Location, Date

The signature field is intended for the consultant and insured person. This confirms that the intended equipment is intended for the insured person.

Signature End customer/dealer

- 8 - 2026-03-03 / The prices are plus statutory value-added tax and shipping costs = included in the basic price = Position selected

Order sheet NJ1 e-assistant Rev. 05 UK