



## LITTY 4you Wheelchair for children

Maximum user weight limit 50 kg  
Weight from 6.7 kg

A9042100000 € 4,560.00

**Quotation date:** \_\_\_\_\_  
**Consultation was carried out by**  
**Staff at the medical supplies shop:** \_\_\_\_\_  
**Enquiry no.:** \_\_\_\_\_

**Order date:** \_\_\_\_\_  
according to GTC PRO ACTIV Reha-Technik GmbH  
**Order no.:** \_\_\_\_\_

### Medical suppliers / company address:

Cust. no.: \_\_\_\_\_  
Name: \_\_\_\_\_  
Street: \_\_\_\_\_  
Postcode / city: \_\_\_\_\_  
Phone: \_\_\_\_\_  
Fax: \_\_\_\_\_  
E-Mail: \_\_\_\_\_  
Contact: \_\_\_\_\_

**Invoicing address:**  same as company address  
 Other (see below)

Cust. no.: \_\_\_\_\_  
Name: \_\_\_\_\_  
Street: \_\_\_\_\_  
Postcode / city: \_\_\_\_\_  
Phone: \_\_\_\_\_  
Fax: \_\_\_\_\_  
E-Mail: \_\_\_\_\_  
Contact: \_\_\_\_\_


**Shipping address:**  same as company address  
 Other (see below)

Name: \_\_\_\_\_  
Street: \_\_\_\_\_  
Postcode / city: \_\_\_\_\_

### Customer name / commission:

ID / Comm. No.: \_\_\_\_\_  
Surname, forename: \_\_\_\_\_  
Street: \_\_\_\_\_  
Postcode / city: \_\_\_\_\_  
Phone: \_\_\_\_\_  
Fax: \_\_\_\_\_  
E-Mail: \_\_\_\_\_  
D.o.B.: \_\_\_\_\_  
Paying organisation: \_\_\_\_\_

### Body size:



Body height: \_\_\_\_\_ cm  
Weight: \_\_\_\_\_ kg  
Lower leg length: \_\_\_\_\_ cm  
(anatomic dimension: hollow of the knee to the foot sole)  
Indication: \_\_\_\_\_


This information is necessary for checking. If the information is missing, a check of the relevant items cannot be made!

# Basic size

## Seat plate and back shell screwed to the side part

Please note: Please select the seat plate and back shell separately.

### Wheelchair dimensions:



Seat width **SB**

<input type="checkbox"/> 16	<input type="checkbox"/> 17	<input type="checkbox"/> 18	<input type="checkbox"/> 19	<input type="checkbox"/> 20	<input type="checkbox"/> 21	<input type="checkbox"/> 22	<input type="checkbox"/> 23
<input type="checkbox"/> 24	<input type="checkbox"/> 25	<input type="checkbox"/> 26	<input type="checkbox"/> 27	<input type="checkbox"/> 28	<input type="checkbox"/> 29	<input type="checkbox"/> 30	<input type="checkbox"/> 31
<input type="checkbox"/> 32	<input type="checkbox"/> 33	<input type="checkbox"/> 34	<input type="checkbox"/> 35	<input type="checkbox"/> 36	<input type="checkbox"/> 37	<input type="checkbox"/> 38	

Seat depth **ST**

<input type="checkbox"/> 18	<input type="checkbox"/> 19	<input type="checkbox"/> 20	<input type="checkbox"/> 21	<input type="checkbox"/> 22	<input type="checkbox"/> 23	<input type="checkbox"/> 24	<input type="checkbox"/> 25
<input type="checkbox"/> 26	<input type="checkbox"/> 27	<input type="checkbox"/> 28	<input type="checkbox"/> 29	<input type="checkbox"/> 30	<input type="checkbox"/> 31	<input type="checkbox"/> 32	<input type="checkbox"/> 33
<input type="checkbox"/> 34	<input type="checkbox"/> 35	<input type="checkbox"/> 36	<input type="checkbox"/> 37	<input type="checkbox"/> 38			

Back height **RH**

<input type="checkbox"/> 20	<input type="checkbox"/> 21	<input type="checkbox"/> 22	<input type="checkbox"/> 23	<input type="checkbox"/> 24	<input type="checkbox"/> 25	<input type="checkbox"/> 26	<input type="checkbox"/> 27	<input type="checkbox"/> 28	<input type="checkbox"/> 29
<input type="checkbox"/> 30	<input type="checkbox"/> 31	<input type="checkbox"/> 32	<input type="checkbox"/> 33	<input type="checkbox"/> 34	<input type="checkbox"/> 35	<input type="checkbox"/> 36	<input type="checkbox"/> 37	<input type="checkbox"/> 38	<input type="checkbox"/> 39
<input type="checkbox"/> 40	<input type="checkbox"/> 41	<input type="checkbox"/> 42	<input type="checkbox"/> 43	<input type="checkbox"/> 44	<input type="checkbox"/> 45				

Lower leg length **USL**

<input type="checkbox"/> 20	<input type="checkbox"/> 21	<input type="checkbox"/> 22	<input type="checkbox"/> 23	<input type="checkbox"/> 24	<input type="checkbox"/> 25	<input type="checkbox"/> 26	<input type="checkbox"/> 27	<input type="checkbox"/> 28	<input type="checkbox"/> 29
<input type="checkbox"/> 30	<input type="checkbox"/> 31	<input type="checkbox"/> 32	<input type="checkbox"/> 33	<input type="checkbox"/> 34	<input type="checkbox"/> 35	<input type="checkbox"/> 36	<input type="checkbox"/> 37	<input type="checkbox"/> 38	<input type="checkbox"/> 39
<input type="checkbox"/> 40	<input type="checkbox"/> 41	<input type="checkbox"/> 42	<input type="checkbox"/> 43	<input type="checkbox"/> 44	<input type="checkbox"/> 45				

Seat height back **SHh**

<input type="checkbox"/> 31	<input type="checkbox"/> 32	<input type="checkbox"/> 33	<input type="checkbox"/> 34	<input type="checkbox"/> 35	<input type="checkbox"/> 36
<input type="checkbox"/> 37	<input type="checkbox"/> 38	<input type="checkbox"/> 39	<input type="checkbox"/> 40	<input type="checkbox"/> 41	<input type="checkbox"/> 42

Seat height front **SHv**

<input type="checkbox"/> 34	<input type="checkbox"/> 35	<input type="checkbox"/> 36	<input type="checkbox"/> 37	<input type="checkbox"/> 38	<input type="checkbox"/> 39	<input type="checkbox"/> 40
<input type="checkbox"/> 41	<input type="checkbox"/> 42	<input type="checkbox"/> 43	<input type="checkbox"/> 44	<input type="checkbox"/> 45	<input type="checkbox"/> 46	<input type="checkbox"/> 47

Tipping point **KP**

(inner edge of rear tube to centre of drive wheel axle)

\_\_\_\_\_ cm (If not specified, a standard value of 8 cm will be used)

**H:** Position of the steering wheel bearing block

\_\_\_\_\_ cm (If not specified: 20" and 22" = 35 cm or 24" = 38 cm)

**I:** Distance between steering wheel bearing blocks

\_\_\_\_\_ cm (If not specified: Seat width + 19 cm)

**J:** Distance front edge footrest to back

\_\_\_\_\_ cm (If not specified: Seat depth + 15 cm)

Different dimensions are available for an additional charge

**The dimensional tolerances according to DIN EN ISO 13920 apply to the design (dimensions in mm)**

Accuracy level	up to 30	over 30 to 120	over 120 to 400	over 400 to 1000
<b>C</b>	<b>+/- 1</b>	<b>+/- 3</b>	<b>+/- 4</b>	<b>+/- 6</b>

## Frame colour - One-colour design

### Frame, side parts and seat system

Information on current colour designs and example images can be found at: <https://www.proactiv-gmbh.com/colours/>

#### Frame colour: Selection from colour palette

##### Frame colour will be specified in the order

A9000100170-000

##### Raw, blasted, not coated

A9000100170-020

Note: In combination with moisture e.g. hand perspiration, can affect dark aluminum discoloration.

#### Powder colours

##### white RAL 9016

gloss A9000100170-023

matt A9000100170-123 € 43.40

##### brilliant silver

gloss A9000100170-011

##### anthracite RAL 7016

gloss A9000100170-001

matt A9000100170-101 € 43.40

##### black RAL 9005

gloss A9000100170-003

matt A9000100170-019 € 43.40

##### pure orange RAL 2004

gloss A9000100170-022

matt A9000100170-122 € 43.40

##### pink magenta RAL 4010

gloss A9000100170-037

matt A9000100170-137 € 43.40

##### pastel blue RAL 5024

gloss A9000100170-036

matt A9000100170-136 € 43.40

##### lemon green

gloss A9000100170-005

matt A9000100170-105 € 43.40

##### yellow RAL 1023

gloss A9000100170-046

matt A9000100170-146 € 43.40

##### red RAL 3001

gloss A9000100170-047

matt A9000100170-147 € 43.40

##### lagoon blue

gloss A9000100170-052

matt A9000100170-152 € 43.40

##### athletic grey RAL 7046

gloss A9000100170-050

matt A9000100170-150 € 43.40

##### muddy grey RAL 7033

gloss A9000100170-051

matt A9000100170-151 € 43.40

##### peach RAL 3022

gloss A9000100170-038

matt A9000100170-138 € 43.40

#### Powder colours „dormant“

##### deep black effect silver

gloss A9000100170-007

##### sky blue

gloss A9000100170-017

##### navy blue

gloss A9000100170-045

matt A9000100170-145 € 43.40

##### electric green

gloss A9000100170-015

matt A9000100170-115 € 43.40

##### violet

gloss A9000100170-044

##### fire red

gloss A9000100170-018

matt A9000100170-118 € 43.40

## Frame colour - One-colour design special colours „liquid coating“

### Frame, side parts and seat system

€ 349.60

##### storm (ice blue)

gloss A9000100170-500

matt A9000100170-501

##### prairie (rose gold)

gloss A9000100170-504

matt A9000100170-505

##### berry

gloss A9000100170-526

matt A9000100170-527

##### blue ocean

gloss A9000100170-522

matt A9000100170-523

##### emerald

gloss A9000100170-514

matt A9000100170-515

##### chili

gloss A9000100170-516

matt A9000100170-517

## Frame colour - Two-colour design

Note: In the two-colour design, both colours are the same gloss level as the design colour.

- Two-colour design consisting of one frame colour and a second colour for the side parts resp. seat plate and back shell** € 268.50  
A9042100170-200

Please select colours below:

Two-colour design LITTY consisting of a base colour at the frame with axle tube and a second colour for the side parts and if desired, for the seat system. Select from one-colour design:

A: Frame colour: A9000100170-

B: Side part colour: A9000100170-

Colour seat plate and back shell:  A or  B

## Colour - Attachment parts

Note: If not specified, the caster fork and bearing block will be mounted in black.

**Drive wheel: aluminum rim, hub and spokes**

- Black A9000100171-003

**Caster fork and bearing block**

- Silver A9000100173-011  
 Black A9000100173-003

## Colour design of attachment parts according to customer request:

**Drive wheel: aluminum rim and design inlet according to colour selection / hub and spokes black**

- frame colour A9000100171-999 € 249.60  
 Selection from one-colour design (different from the frame colour) A9000100171-998 € 288.90  
Colour: \_\_\_\_\_  
A9000100170-

**Caster wheel: 5" aluminum rim<sup>1</sup>**  
(6-spoke design)

- frame colour A9000100174-999 € 75.20  
 Selection from one-colour design (different from the frame colour) A9000100174-998 € 87.10  
Colour: \_\_\_\_\_  
A9000100170-

## Special colour request (Frame colour - One-colour design)

**Frame colour for single purchasing Selection from RAL colour palette**

- A9000100170-100 € 414.40

**gloss level**

- gloss  
 matt

**colour designation**

\_\_\_\_\_

**RAL number**

\_\_\_\_\_

**Note:** Attachments **not** available in design colour options

<sup>1</sup> Only in combination with caster wheel 5" aluminum rim 6-spoke-design good shock absorption A9000501210 or A9000501211

## Push handles

### Equipment without push handles

A9000401000



### Central push handle aluminum

A9042401085-K  
height-adjustable  
Holder welded to the center of the back crossbar



### Safety push handles back-positioned

Length: 45 cm  
max. height = SHh + approx. 53 cm  
A9042401075-K001

Length: 55 cm  
max. height = SHh + approx. 63 cm  
A9042401075-K002  
continuous height-adjustable  
push handle holder back-positioned



### Aluminium central push handle for brake handles Foldable handle area

A9042401086-K € 141.90


height-adjustable and foldable in the handle area via toothed joint, holder welded on the middle on the back cross bar.  
Approved up to a body weight of 50 kg  
Only in combination with A9042701020 Parking knee lever brake for accompanying person



### Aluminium central push handle Foldable handle area

A9042401084-K € 102.60

height-adjustable and foldable in the handle area via toothed joint, holder welded on the middle on the back cross bar.  
Approved up to a body weight of 50 kg



### Aluminium push handle closed, foldable via toothed joint

A9042401069-K € 389.90

continuous height-adjustable  
push handles shifted to the back



## Back equipment

If not specified, the „ergonomic back shell" will be delivered.

### Ergonomic back shell aluminum

A9042802210-K € 441.60

in frame colour  
from RH 20 cm and SB 18 cm



### Back cushion for ergonomic back shell

A9042802220-K € 221.30

Cover:

- Black, breathable, open surface
- Black, breathable, closed surface

If not specified, the cover material „Black, breathable, open surface" will be delivered.



### Equipment without back shell and back cushion

A9042802011



## Seat system

### Seat plate aluminum

A9042801210-K  
in frame colour  
bolted between the aluminum side parts  
adjustable in height, angle and depth



## Seat cushion



### Seat cushion in rectangular contour Inlet: Standard

A9000805180-K\* € 118.30

### Seat cushion with radius contour Inlet: Standard

A9000800300-K\* € 163.10

### Seat cushion height:

3 cm       5 cm  
 4 cm       6 cm

### Seat cushion in rectangular contour Inlet: Anti-Decubitus

up to 30x30 cm  
A9000805204-K\* € 301.70

>30x30 up to 50x50 cm  
A9000805200-K\* € 437.70

### Seat cushion with radius contour Inlet: Anti-Decubitus

up to 30x30 cm  
A9000800311-K\* € 356.10

>30x30 up to 50x50 cm  
A9000800310-K\* € 482.90

### Seat cushion height:

4 cm       6 cm

### Cover:

Black, breathable, open surface

Black, breathable, closed surface

If not specified, the cover material „Black, breathable, open surface“ will be delivered.

Seat cushion is provided

## Wheel camber

### 6° drive wheel camber

A9042202010

### 8° drive wheel camber

A9042202020

### 9,5° drive wheel camber

A9042202060

### 11° drive wheel camber

A9042202030

### 13° drive wheel camber

A9042202040

### 15° drive wheel camber

A9042202050

## E-Adaptations

### Preparation for the „e-motion“ mounting

A9042201200-K € 502.90

Equipment: e-motion torque arm,  
1/2" axle and socket length 50.8 mm (2")  
Up to 40 kg: 9,5° wheel camber,  
over 40 kg: 6° and 8° wheel camber  
Only in combination with 22" or 24"  
PRO ACTIV and Alber recommend mounting  
two anti-tipping supports. Please select under  
„Adaptations“



## Brake system

If not specified, the „Knee lever brake“ will be delivered.

### Knee lever brake

A9042701010  
Brake lever integrated in side part



### Knee lever brake with position lock for brake lever

A9042701030 € 304.30

integrated into the side part when braked and  
unbraked.  
Operated by an accompanying person.



### Parking knee lever brake for accompanying person

A9042701020 € 464.10

Push handle operation via Bowden cable.  
Only in combination with knee lever brake and  
push handles



## Footrest

If not specified, the „fixed footrest“ will be delivered.

### Fixed footrest

A9042601010-K  
adjustable in height, depth and angle  
including footplate cover with anti-slip film



### Footrest folds up to the rear

A9042602010-K € 286.80  
including footplate cover with anti-slip film



## Adaptations

### Anti-tipping support - tiltable upwards

right side A9042301160-12\* € 207.10  
 left side A9042301160-22\* € 207.10  
removable via Quick-pin



### Outdoor wheel attachment with steering fork

A9042900200-K € 845.70  
Equipped with straight-running lock  
can be removed from the wheelchair via quick  
release lever with swivelling pneumatic tyre  
12.1/2 x 2.1/4" caster wheel  
can be removed from the support tube via  
quick-release axle  
for maximum ground pressure at the  
steering attachment of 30 kg



## 3" caster wheel

### 3" skate wheel (80 mm)

A9000501110 € 70.80  
puncture-proof



## 4" caster wheel

### 4" caster wheel „Standard“ solid rubber black

A9000501011-004  
(82A), puncture-proof



### 4" skate wheel (100 mm)

A9000501120 € 70.80  
puncture-proof



### 4" safety light-up wheel (100 mm)

A9000501220 € 94.80  
puncture-proof



### 4" CorMORun caster wheel (100 mm)

A9000501042 € 173.00  
with replaceable rubber 6 mm,  
puncture-proof  
Rim colour: Black  
Incl. hair guard spacer (5.3 mm)



### 4" Frog Legs caster wheel (100 mm)

A9000501230 € 206.30  
with replaceable rubber 6 mm,  
puncture-proof  
Colour: black, silver rim  
Incl. hair guard spacer (5.3 mm)



## 5" caster wheel

### 5" caster wheel „Standard“ PU black

A9000501036-005  
(82A), puncture-proof



### 5" caster wheel „Standard“ PU grey

A9000501031-005  
(82A), puncture-proof



### 5" safety light-up wheel (125 mm)

A9000501221 **€ 94.80**  
puncture-proof



### 5" aluminum rim 3-spoke design „good shock absorption“

silver anodised  
A9000501210 **€ 179.80**

black  
A9000501211 **€ 179.80**

with replaceable rubber 6 mm  
medium hard, puncture-proof  
Incl. hair guard spacer (3 mm)



### 5" CorMORun caster wheel (125 mm)

A9000501041 **€ 181.30**  
with replaceable rubber 6 mm,  
puncture-proof  
Rim colour: Black  
Incl. hair guard spacer (5.3 mm)



### 5" Frog Legs caster wheel (125 mm)

A9000501232 **€ 206.30**  
with replaceable rubber 6 mm,  
puncture-proof  
Colour: black, silver rim  
Incl. hair guard spacer (5.3 mm)



## Caster fork

For caster fork colours, see „Colour of add-on parts“  
If not specified, the „Standard caster fork with screw axle“ will be delivered.

### Caster fork Standard

A9000500100-K001  
with screwed axle

A9000500100-K002 **€ 167.40**  
with quick-release axle



### Caster fork short

A9000500110-K001  
with screwed axle

A9000500110-K002 **€ 167.40**  
with quick-release axle  
Not in combination with 5" caster wheels



### Lenkradgabel-Wechselsatz LITTY

A9000500300-K **€ 296.70**  
with quick-release axle

The article is delivered without caster wheels.  
If caster wheels are desired, please tick a  
second option under the section „Caster  
wheels.“



## Drive wheels

### Drive wheel aluminum rim, 18 spokes with Liteline hub

radially spoked  
Design inlet in frame colour

20" (ETRTO 451)  
A9000203053


22" (ETRTO 489)  
A9000203063

24" (ETRTO 540)  
A9000203123



# Handrims

x = dimension handrim attachment



Handrim attachment  
general:

- x = 11 mm
- x = 19 mm

Handrim attachment  
Curve + Gekko handrim:

- x = 10 mm
- x = 18 mm

If not specified, the handrim is mounted tightly (11 mm / 10 mm).

## Stainless steel Ø 20 mm

- 20" handrim for 20" wheel  
A9000205010-001
- 22" handrim for 22" wheel  
A9000205010-002
- 24" handrim for 24" wheel  
A9000205010-004



## Aluminum Ø 19 mm silver anodised

- 22" handrim for 22" wheel  
A9000205019-002
- 24" handrim for 24" wheel  
A9000205019-004



## Aluminum Ø 16 mm powder coated in frame colour

(basic colour with two-colour design)

- = Gloss surface
- = Matt surface

Note: Matt coating is not recommended due to the low grip.

- 20" handrim for 20" wheel  
A9000205021-001 **€ 68.30**
- 22" handrim for 22" wheel  
A9000205021-002 **€ 68.30**



## Aluminum Ø 19 mm powder coated in frame colour

(basic colour with two-colour design)

- = Gloss surface
- = Matt surface

Note: Matt coating is not recommended due to the low grip.

- 20" handrim for 20" wheel **€ 68.30**  
A9000205020-001
- 22" handrim for 22" wheel **€ 68.30**  
A9000205020-002
- 24" handrim for 24" wheel **€ 68.30**  
A9000205020-004



## Aluminum Ø 16 mm polished and hard-anodised

- 20" handrim for 20" wheel **€ 164.90**  
A9000205031-001
- 22" handrim for 22" wheel **€ 164.90**  
A9000205031-002



## Aluminum Ø 19 mm polished and hard-anodised

- 20" handrim for 20" wheel **€ 164.90**  
A9000205030-001
- 22" handrim for 22" wheel **€ 164.90**  
A9000205030-002
- 24" handrim for 24" wheel **€ 164.90**  
A9000205030-004



## Titanium Ø 19 mm

- 24" handrim for 24" wheel **€ 335.60**  
A9000205040-004



## Silicone handrim cover

- for handrim 20" **€ 232.80**  
A9000205100-001
  - for handrim 22" **€ 232.80**  
A9000205100-002
  - for handrim 24" **€ 232.80**  
A9000205100-004
- Surface without structure, compatible with handrim Ø 19 mm or Ø 20 mm.



## Ultra-Grip „very grippy“

- 20" handrim for 20" wheel **€ 402.40**  
A9000205070-001
- 22" handrim for 22" wheel **€ 402.40**  
A9000205070-002
- 24" handrim for 24" wheel **€ 402.40**  
A9000205070-004



<p><b>Ergo-Grip</b></p> <p><input type="checkbox"/> 20" handrim for 20" wheel A9000205080-001 € 402.40</p> <p><input type="checkbox"/> 22" handrim for 22" wheel A9000205080-002 € 402.40</p> <p><input type="checkbox"/> 24" handrim for 24" wheel A9000205080-004 € 402.40</p> <p>Attachment fitting to the rim</p> <p>Note: Due to the complex geometry, concentricity deviations of +/- 2-3 mm are possible.</p> 	<p><b>Maxgripp with cellular and natural rubber coating very grippy, soft</b></p> <p><input type="checkbox"/> 20" handrim for 20" wheel A9000200788 € 402.40</p> <p><input type="checkbox"/> 22" handrim for 22" wheel A9000200789 € 402.40</p> <p><input type="checkbox"/> 24" handrim for 24" wheel A9000200787 € 402.40</p> 	<p><b>Curve black anodised, ergonomic handrim profile</b></p> <p><input type="checkbox"/> 22" handrim for 22" wheel Profile height 24 mm A9000205130-002 € 402.40</p> <p><input type="checkbox"/> 24" handrim for 24" wheel Profile height 24 mm A9000205130-004 € 402.40</p> 	
<p><b>Curve powder-coated, ergonomic handrim profile</b></p> <p><input type="checkbox"/> = frame colour (gloss surface)</p> <p><input type="checkbox"/> = black (gloss surface)</p> <p><input type="checkbox"/> 24" handrim for 24" wheel Profile height 24 mm A9000205135-004 € 467.60</p> 	<p><b>Curve Grip, with Tetra-Grip coating, ergonomic handrim profile</b></p> <p><input type="checkbox"/> 24" handrim for 24" wheel Profile height 24 mm A9000205155-004 € 520.70</p> 	<p><b>Quadro with Tetra-Grip coating</b></p> <p><input type="checkbox"/> 24" handrim for 24" wheel A9000205170-004 € 520.70</p> 	
<p><b>Gekko black anodised</b></p> <p><input type="checkbox"/> 22" handrim for 22" wheel A9000205175-002 € 520.70</p> <p><input type="checkbox"/> 24" handrim for 24" wheel A9000205175-004 € 520.70</p> <p>with a robust and replaceable silicone lip, fully integrated into the profile Profile height 31 mm</p> 	<p><b>Gekko silver anodised</b></p> <p><input type="checkbox"/> 24" handrim for 24" wheel A9000205185-004 € 520.70</p> <p>with a robust and replaceable silicone lip, fully integrated into the profile Profile height 31 mm</p> 	<p><b>Gekko S black anodised</b></p> <p><input type="checkbox"/> 22" handrim for 22" wheel A9000205195-002 € 520.70</p> <p><input type="checkbox"/> 24" handrim for 24" wheel and 25" wheel A9000205195-004 € 520.70</p> <p>with a robust and replaceable silicone lip, fully integrated into the profile Profile height 24 mm</p> 	
<p><b>Equipment without handrims</b></p> <p><input checked="" type="checkbox"/> A9000205990</p>			
<p><b>Spoke guard</b></p>			
<p><b>Finger- and spoke guard „transparent“</b></p> <p><input type="checkbox"/> for 20" wheel A9000206300-001 € 207.10</p> <p><input type="checkbox"/> for 22" wheel A9000206300-002 € 207.10</p>			<p><input type="checkbox"/> for 24" wheel A9000206300-004 € 207.10</p> 

## Tyres

If not specified, the „Schwalbe RightRun in black“ tyres will be delivered.

### Schwalbe „RightRun“ High-pressure tyres

22x1" (489-25)  
A9000204001-002

24x1" (540-25)  
A9000204001-004

Colour: Black  
7.5 bar



### Schwalbe „RightRun“ High-pressure tyres

20x1" (451-25)  
A9000204000-001

22x1" (489-25)  
A9000204000-002

24x1" (540-25)  
A9000204000-004

Colour: Grey, black wall  
7.5 bar



### Schwalbe „Marathon Plus Evolution“ with „smart guard“ puncture protection and reflector strips

22x1" (489-25)  
A9000204030-002

€ 49.20

Colour: Black  
10 bar



### Schwalbe „Marathon Plus Evolution“ with „smart guard“ puncture protection

24x1" (540-25)  
A9000204030-014

€ 49.20

Colour: Black  
10 bar



### Schwalbe „MarathonPlus“ Schwalbe „Marathon Plus“ with „smart guard“ puncture protection

24x1" (540-25)  
A9000204030-004

€ 49.20

Colour: Grey, black wall  
10 bar




## Accessories

### Through axle with assistance, Large push button

A9000200096 € 105.80

Interlock dimension: 122 / 110 mm  
Axle diameter: 12 mm



### Hand air compressor with pressure gauge

A9000901100 € 106.30

up to 14 bar



### Length-adjustable safety belt with buckle

A9000902100 € 80.60



### Adjustable safety belt with hook and fastener

A9000902110 € 64.50



### Transport bag for drive wheels

for 24" wheel  
A9000903530-004 € 77.80  
Ø 76 cm x 2 cm



### Therapy table

width 350 mm  
A9000902800-001 € 480.20

removable table plate from Lexan, transparent,  
with front edging, with 2 positioning stops  
Ø 25 mm (for adaptation to frame tube  
diameter Ø 25 mm)



## Services

Consultation, clarification of the installation situation and measurement through PRO ACTIV

A9990000000 € 354.50

Handover and instruction on the wheelchair by PRO ACTIV

A9990000006 € 354.50

Notes and remarks:

### Location, Date

### Signature End customer/dealer

The signature field is intended for the consultant and insured person. This confirms that the intended equipment is intended for the insured person.  
- 12 - 2026-03-03 / The prices are plus statutory value-added tax and shipping costs  = included in the basic price  = Position selected  
Order sheet LITTY 4you Rev. 06 UK