

Customized production - model „SPEEDY 4all Ergo“

Frame design

Frame design Heavy Duty AF5 up to 150 kg body weight

A9017100011 € 730.60

from 121 kg to 150 kg payload resp.
SB 47 cm to 52 cm

- Frame reinforcement with 2 additional diagonal frame cross tubes and reinforcements
 - Reinforced drive wheel mount
 - Stainless steel backrest locking pin
 - Back tube reinforcement inlay on side part screw thread
 - Push handle permanently mounted, steel Ø 20 x 1.5 mm
 - Push handles in stainless steel design
 - Side plate aluminium 3 mm
- Please note: Aluminium clothing guard removable via lock function (A9000403210-K) additionally required!
- Footrest or foot bar, continuous, reinforced design
 - Knee lever brake bolt, stainless steel

Note:

Please select the „Heavy Duty“ steering wheel fork under Steering wheel fork!

Services

Consulting and customization by PRO ACTIV according to individual customer requirements

A9990000017 € 612.00

Additional requirements