



## FOXX offroad compact bike

approved up to 120 kg payload

### basic equipment

- Front wheels 24" (507-30) tubeless ready
- Rear wheel 27.5" (587-35) tubeless ready
- full suspension with 100 mm travel at the front and 140 mm travel at the rear
- A-shaped handles with integrated shift and brake controls
- Hydraulic disc brake system at front and rear
- Backrest angle and distance adjustable
- Adjustable bottom bracket position
- Mudguard for the rear wheel
- Underride guard
- Battery lighting set

A9105100000

€ 16,975.00

**Quotation date:** \_\_\_\_\_

**Consultation was carried out by**  
**Staff at the medical supplies shop:** \_\_\_\_\_

**Enquiry no.:** \_\_\_\_\_

**Order date:** \_\_\_\_\_

according to GTC PRO ACTIV Reha-Technik GmbH

**Order no.:** \_\_\_\_\_

### Medical suppliers / company address:

**Cust. no.:** \_\_\_\_\_

**Name:** \_\_\_\_\_

**Street:** \_\_\_\_\_

**Postcode / city:** \_\_\_\_\_

**Phone:** \_\_\_\_\_

**Fax:** \_\_\_\_\_

**E-Mail:** \_\_\_\_\_

**Contact:** \_\_\_\_\_

**Different invoicing address (see below)**

**Cust. no.:** \_\_\_\_\_

**Name:** \_\_\_\_\_

**Street:** \_\_\_\_\_

**Postcode / city:** \_\_\_\_\_

**Different shipping address (see below)**

**Cust. no.:** \_\_\_\_\_

**Name:** \_\_\_\_\_

**Street:** \_\_\_\_\_

**Postcode / city:** \_\_\_\_\_

### Customer name / commission:

**ID / Comm. No.:** \_\_\_\_\_

**Surname, forename:** \_\_\_\_\_

**Street:** \_\_\_\_\_

**Postcode / city:** \_\_\_\_\_

**Phone:** \_\_\_\_\_

**E-Mail:** \_\_\_\_\_

**D.o.B.:** \_\_\_\_\_

**Paying organisation:** \_\_\_\_\_

### Body size

Body height: \_\_\_\_\_ cm

Upper body height: \_\_\_\_\_ cm

Weight: \_\_\_\_\_ kg

Indication: \_\_\_\_\_

This information is necessary for checking.

## Drive system Bafang

The specified ranges depend on the route profile, user weight and speed.

### Bafang M560 motor

Maximum power for all off-road adventures

A9105100010  
500W, 130Nm torque,  
assistance up to 25 km/h,  
starting assistance up to 15 km/h



### Battery Bafang 48V 696Wh

Range: approx. 50 km - 75 km

A9105900100  
Battery mounted in the main frame under the  
seat plate



### Double battery Bafang 48V 696Wh

Range: approx. 100 km - 150 km

A9105900105      **€ 880.70**  
Batteries mounted in the main frame under the  
seat plate.  
Interchangeable plug-in connection mounted  
on the seat plate for easy switching between  
batteries.



## Drive system Shimano

The specified ranges depend on the route profile, user weight and speed.

### Shimano EP6 Cargo motor

Drive system for off-road and on-road use

A9105100020      **€ 334.80**  
250W, 85Nm torque,  
assistance up to 25 km/h,  
starting aid up to 6 km/h  
Pedelec system with road approval



### Battery Shimano 36V 630Wh

Range: approx. 45 km - 70 km

A9105900110      **€ 628.70**  
Battery mounted in the main frame under the  
seat plate



### Double battery Shimano 36V 630Wh

Range: approx. 90 km - 140 km

A9105900115      **€ 1,447.50**  
Batteries mounted in the main frame under the  
seat plate.  
Interchangeable plug-in connection mounted  
on the seat plate for easy switching between  
batteries.



## Colour attachments

If not specified, the legrest, rear mudguard and chain guard disc will be delivered in UDS black.

### Leg rest, rear mudguard and chain guard disc UDS black

A9105100170-001

### Leg rest, rear mudguard and chain guard disc Pure orange RAL 2004 gloss

A9105100170-002

### Leg rest, rear mudguard and chain guard disc customized colour selection

A9105100170-003      **€ 225.60**

from colour palette one-colour design  
powder colours wheelchair

A9000100170- \_\_\_\_\_

Colour: \_\_\_\_\_

## Leg rest

### Leg rest size M

for body heights from 175 cm to 200 cm

Included in the base price



### Leg rest size S

up to a body height of 175 cm

A9105100180



## Derailleur gears

### Shimano Cues U8000 11-speed derailleur gear system

A9105200010  
with Linkglide technology  
Cassette 11-50 T



### Sram GX-Eagle AXS 12-speed derailleur gear system

A9105200011      € 790.50  
electric radio control  
Cassette 11-50 T



## Gear hub

### Rohloff 14-speed gear hub with electric operation

A9105200020      € 2,047.10  
powered by the Bafang system



## Chassis

### Standard chassis

A9105400010  
Fox Float DPS Performance front and Manitou  
Mara Inline rear with three modes each (Open,  
Medium, Firm)  
Rebound and compression adjustable

### High-quality off-road chassis

A9105400020      € 865.20  
Fox Float DPS Factory front and Manitou Mara  
Pro rear with three modes each (Open, Medium,  
Firm)  
Rebound and compression adjustable  
High- and low-speed compression adjustable

## Dimensions

### Crank type and grip width

#### Standard crank

- grip width 40 cm  
A9110500100-040
- grip width 45 cm  
A9110500100-045
- grip width 50 cm  
A9110500100-050

#### Tetra crank (flared)

with extended crank arms for maximum  
upper body space

- grip width 48 cm (176 mm)  
A9110500100-048
- grip width 53 cm (226 mm)  
A9110500100-053
- grip width 58 cm (276 mm)  
A9110500100-058

#### Standard pedal crank

- length 155 mm  
A9110600100-155
- length 175 mm  
A9110600100-175
- length 195 mm  
A9110600100-195



#### Tetra pedal crank

- length 155 mm      € 136.10  
A9100600105-155
- length 175 mm      € 136.10  
A9100600105-175
- individual length  
195 mm to 140 mm      € 258.70  
A9100600100



### Adjustable bottom bracket support

the angle and length can be adjusted via the clamped rotary point on the fork bridge and length.

For the pre-setting, please specify the bottom bracket support size that comes as close to your desired position.

S (180 mm)  
A9105500010-01

M (240 mm)  
A9105500010-02

If not specified, the adjustable bottom bracket support will be delivered in S.



## Handle options

### Palm rest on the handle

A9110601200 € 102.80

upholstered, bolted, angle-adjustable.  
Not in combination with other handles.



### Tetra hand support, pivotable, adjustable grip distance

A9110601101 € 169.40



### Tetra backrest swivel and adjustable grip width

A9110601102 € 191.60

very stiff external support as round tube  
Ø 20 mm with foam rubber coating



## Upper body fixation

### Chest strap for upper body fixation

size S (105 cm) € 241.80  
A9000800016

size M (125 cm) € 241.80  
A9000800013

size L (145 cm) € 241.80  
A9000800014

consisting of elastic webbing material in combination with neoprene with Velcro and operating tab.  
Colour: black



### Pelvi.Loc thigh fixation

A9105800000 € 339.90



## Tyres

### Offroad tyres

A9105105110  
Front: Smart Sam 24 x 2.1"  
Rear: Magic Marry 27.5 x 2.6" Soft



### Offroad tyres 2

A9105105111 € 61.00  
Front: Hans Dampf 24 x 2.35"  
Rear: Magic Marry 27.5 x 2.6" Soft



### Trekking tyres

A9105105100 € 48.20  
Front: Land Cruiser Plus 24 x 2.0"  
Rear: Johnny Watts 27.5 x 2.6"



<b>Puncture protection band</b> <input type="checkbox"/> A9000206100-K      € 55.80		
<b>Accessories</b>		
<b>Lighting for front and rear</b> <input type="checkbox"/> A9105900000      € 225.60 incl. reflectors for StVZO 	<b>Fidlock TWIST water bottle</b> <input type="checkbox"/> A9105900003      € 52.60 mounted on the head tube bottle with a capacity of 590 ml 	<b>Spacers for thigh fixation with a height of 20 mm</b> <input type="checkbox"/> A9105900010      € 38.90 each 20 mm high (max. 60 mm high / 3 pairs) ____ pair 
<b>Mudguards for the front wheels</b> <input type="checkbox"/> A9105900001      € 213.20 UDS black coated 	<b>Damper pump</b> <input type="checkbox"/> A9105900015      € 40.90 up to 20 bar 	<b>Bags and tool set</b> <input type="checkbox"/> A9105900004      € 142.00 Set consisting of a bag placed on the frame next to the damper (1.1 liter capacity) which for example offers space for a cell phone, wallet and tubes and another bag placed between the tubes of the head tube mount (0.6 liter capacity) which is ideal for storing tools. 
<b>Services</b>		
<b>Consultation, clarification of the installation situation and measurement through PRO ACTIV</b> <input type="checkbox"/> A9990000000      € 354.50	<b>Handover and instruction of the handbike / traction device by PRO ACTIV</b> <input type="checkbox"/> A9990000011      € 354.50	

Notes and remarks:

Location, Date

Signature End customer/dealer

The signature field is intended for the consultant and insured person. This confirms that the intended equipment is intended for the insured person.

## LEGAL NOTICE - FOXX offroad handbike with SHIMANO or BAFANG drive system

PRO ACTIV offers the FOXX offroad handbike with a choice of two different drive systems. The utilisation options differ depending on the equipment and the applicable legal regulations.

### SHIMANO drive system (pedelec system according to EN 15194)

The FOXX with SHIMANO drive system fulfils the requirements of a pedelec in accordance with EU standard EN 15194. With this equipment, the vehicle is permitted for use on public roads and public paths in Germany without the need for an operating permit or an insurance licence plate. This also applies to use on cycle paths.

### BAFANG drive system (more powerful off-road system)

The FOXX with BAFANG M560 drive system has a motor output of 500 W and a torque of 130 Nm. PRO ACTIV expressly points out that the FOXX with BAFANG drive system may not be operated on public roads or on public paths, squares or cycle paths in Germany without an individual operating licence and without an insurance licence plate. In this case, operation is only permitted on private property and off-road terrain (e.g. on non-public trails or bike parks). Operation in public areas without a valid operating licence constitutes an infringement of the Compulsory Insurance Act (§ 1 PflVG) and can be prosecuted in accordance with § 6 Para. 1 PflVG.

Please also note: If you use the FOXX with BAFANG drive outside Germany, it is your responsibility to inform yourself in advance about the legal requirements and traffic regulations applicable in the respective country.

### Confirmation of knowledge

I confirm with my signature that I have read, understood and taken note of the above information. I undertake to use the FOXX with BAFANG drive system in Germany only with a valid operating licence on public roads and only on private property or off-road terrain without an operating licence. When using the vehicle in other countries, I will inform myself on my own responsibility about the legal requirements applicable there and comply with them.

---

**Customer name in block letters**

---

**Place, date**

**Signature (end customer)**