



SPEEDY F1

Rigid frame sports wheelchair
 with customized frame geometry and custom-
 made equipment,
permanently welded chassis
 Maximum user weight limit 120 kg

A9011100000 € 5,240.00

Quotation date: _____
 Consultation was carried out by
 Staff at the medical supplies shop: _____
Enquiry no.: _____

Order date: _____
 according to GTC PRO ACTIV Reha-Technik GmbH
Order no.: _____

Medical suppliers / company address:

Cust. no.: _____
Name: _____
Street: _____
Postcode / city: _____
Phone: _____
Fax: _____
E-Mail: _____
Contact: _____

Invoicing address: same as company address
 Other (see below)

Cust. no.: _____
Name: _____
Street: _____
Postcode / city: _____
Phone: _____
Fax: _____
E-Mail: _____
Contact: _____

Shipping address: same as company address
 Other (see below)

Name: _____
Street: _____
Postcode / city: _____

Customer name / commission:

ID / Comm. No.: _____
Surname, forename: _____
Street: _____
Postcode / city: _____
Phone: _____
Fax: _____
E-Mail: _____
D.o.B.: _____
Paying organisation: _____

SPEEDY F1 basket



SPEEDY F1 tennis



SPEEDY F1 badminton



SPEEDY F1 4dance



Basic size


Körpermaße:

Body height: _____ cm

Weight: _____ kg


Lower leg length: _____ cm
(anatomic dimension: hollow of the knee to the foot sole)

Indication: _____



This information is necessary for checking. If the information is missing, a check of the relevant items cannot be made!

Wheelchair dimensions:



- | | |
|--|---|
| <p>A: SB Seat width _____ cm</p> <p>B: ST Seat depth _____ cm</p> <p>C: RH Back height _____ cm</p> <p>D: USL Lower leg length _____ cm</p> <p>E: SHh Seat height back _____ cm</p> <p>F: SHv Seat height front _____ cm</p> <p>G: KP Tipping point _____ cm
(inner edge of rear tube to centre of drive wheel axle)</p> <p>H: Frame length to the front frame bend _____ cm</p> | <p>I: Start frame retraction _____ cm
(I - min. 6 cm)</p> <p>J: KW Knee angle _____ °</p> <p>or</p> <p>J1: Angle of the front Framework tube to vertical _____ °</p> <p>K: Position of the steering wheel bearing block _____ cm
(inner edge of backrest to center of bearing block)</p> <p>L: Distance between steering wheel bearing blocks _____ cm</p> <p>M: Lower leg width top _____ cm</p> <p>N: Lower leg width at the bottom _____ cm</p> <p>O: Position of the tip support bearing block _____ cm
(in basketball, according to the rules, the drive wheel forms the rearmost point)</p> |
|--|---|

The dimensional tolerances according to DIN EN ISO 13920 apply to the design (dimensions in mm)

Accuracy level	up to 30	over 30 to 120	over 120 to 400	over 400 to 1000
C	+/- 1	+/- 3	+/- 4	+/- 6

Frame colour - One-colour design

Information on current colour designs and example images can be found at: <https://www.proactiv-gmbh.com/colours/>

Frame colour: Selection from colour palette

Frame colour will be specified in the order

A9000100170-000

Raw, blasted, not coated

A9000100170-020

Note: In combination with moisture e.g. hand perspiration, can affect dark aluminum discoloration.

Powder colours

white RAL 9016

gloss A9000100170-023

matt A9000100170-123 € 43.40

brilliant silver

gloss A9000100170-011

anthracite RAL 7016

gloss A9000100170-001

matt A9000100170-101 € 43.40

black RAL 9005

gloss A9000100170-003

matt A9000100170-019 € 43.40

pure orange RAL 2004

gloss A9000100170-022

matt A9000100170-122 € 43.40

pink magenta RAL 4010

gloss A9000100170-037

matt A9000100170-137 € 43.40

pastel blue RAL 5024

gloss A9000100170-036

matt A9000100170-136 € 43.40

lemon green

gloss A9000100170-005

matt A9000100170-105 € 43.40

yellow RAL 1023

gloss A9000100170-046

matt A9000100170-146 € 43.40

red RAL 3001

gloss A9000100170-047

matt A9000100170-147 € 43.40

lagoon blue

gloss A9000100170-052

matt A9000100170-152 € 43.40

athletic grey RAL 7046

gloss A9000100170-050

matt A9000100170-150 € 43.40

muddy grey RAL 7033

gloss A9000100170-051

matt A9000100170-151 € 43.40

peach RAL 3022

gloss A9000100170-038

matt A9000100170-138 € 43.40

Powder colours „dormant“

deep black effect silver

gloss A9000100170-007

sky blue

gloss A9000100170-017

navy blue

gloss A9000100170-045

matt A9000100170-145 € 43.40

electric green

gloss A9000100170-015

matt A9000100170-115 € 43.40

violet

gloss A9000100170-044

fire red

gloss A9000100170-018

matt A9000100170-118 € 43.40

Frame colour - One-colour design special colours „liquid coating“

€ 349.60

storm (ice blue)

gloss A9000100170-500

matt A9000100170-501

prairie (rose gold)

gloss A9000100170-504

matt A9000100170-505

berry

gloss A9000100170-526

matt A9000100170-527

blue ocean

gloss A9000100170-522

matt A9000100170-523

emerald

gloss A9000100170-514

matt A9000100170-515

chili

gloss A9000100170-516

matt A9000100170-517

Frame colour - Two-colour design

Note: In the two-colour design, both colours are the same gloss level as the design colour.

- Two-colour design with one base colour (on the rear frame) and one design colour (on the front frame)** € 268.50

Please select a colour below:

Base colour (on the rear frame)

white RAL 9016 <input type="checkbox"/> A9000100170-200	brilliant silver <input type="checkbox"/> A9000100170-201	anthracite RAL 7016 <input type="checkbox"/> A9000100170-202
black RAL 9005 <input type="checkbox"/> A9000100170-203	pastel blue RAL 5024 <input type="checkbox"/> A9000100170-206	lemon green <input type="checkbox"/> A9000100170-207
yellow RAL 1023 <input type="checkbox"/> A9000100170-208	lagoon blue <input type="checkbox"/> A9000100170-212	athletic grey RAL 7046 <input type="checkbox"/> A9000100170-210
muddy grey RAL 7033 <input type="checkbox"/> A9000100170-211		

Design colour (on the front frame)

black RAL 9005 <input type="checkbox"/> gloss A9000100170-259 <input type="checkbox"/> matt A9000100170-260	lemon green <input type="checkbox"/> gloss A9000100170-250 <input type="checkbox"/> matt A9000100170-251	chili <input type="checkbox"/> gloss A9000100170-253 <input type="checkbox"/> matt A9000100170-254
pure orange RAL 2004 <input type="checkbox"/> gloss A9000100170-255 <input type="checkbox"/> matt A9000100170-256	emerald <input type="checkbox"/> gloss A9000100170-274 <input type="checkbox"/> matt A9000100170-275	storm (ice blue) <input type="checkbox"/> gloss A9000100170-264 <input type="checkbox"/> matt A9000100170-265
prairie (rose gold) <input type="checkbox"/> gloss A9000100170-266 <input type="checkbox"/> matt A9000100170-267	berry <input type="checkbox"/> gloss A9000100170-284 <input type="checkbox"/> matt A9000100170-285	blue ocean <input type="checkbox"/> gloss A9000100170-280 <input type="checkbox"/> matt A9000100170-281

Colour - Attachment parts

Note: If not specified, the caster fork and bearing block will be mounted in black.

Drive wheel: aluminum rim, hub and spokes

- Black A9000100171-003

Caster fork and bearing block

- Silver A9000100173-011
 Black A9000100173-003

Colour design of attachment parts according to customer request:

Drive wheel: aluminum rim and design inlet¹ according to colour selection / hub and spokes black

- frame colour A9000100171-999 € 249.60
- Selection from one-colour design (different from the frame colour) A9000100171-998 € 288.90
- Colour: _____
A9000100170-

Caster wheel: 5" aluminum rim² (6-spoke design)

- frame colour A9000100174-999 € 75.20
- Selection from one-colour design (different from the frame colour) A9000100174-998 € 87.10
- Colour: _____
A9000100170-

Special colour request (Frame colour - One-colour design)

Frame colour for single purchasing Selection from RAL colour palette

- A9000100170-100 € 414.40

gloss level

- gloss
 matt

colour designation

RAL number

Note: Attachments **not** available in design colour options

¹ Only in combination with Liteline hub

² Only in combination with caster wheel 5" aluminum rim 6-spoke-design good shock absorption A900050121*, as well as with suspension fork 5-spoke-design A9000500160*


Backrest

Aluminium side guard

A9015402470-002

Backrest screwed, not foldable, back angle 90° to the floor.

Note: Please note that the tyre may not be completely covered.




Aluminium side guard, customized contour

A9015402470-001 € 167.70

Measurement adjustment:

- A = ____ cm / B = ____ cm
 Adjustment to the upper edge of the tyre
 Adjustment to the upper edge of the rim

If not specified, the side guard will be adjusted to the upper edge of the tyre.



Side part with welded tube bracket, contour following the wheel arch

A9011400400-002 € 454.20

with filling plate

Backrest screwed, not foldable, Adjustment to the upper edge of the tyre in frame colour

Not in combination with Body Contour



Side part with welded tube bracket, customized contour

A9011400400-001 € 524.60


with filling plate

Backrest screwed, not foldable, in frame colour

Measurement adjustment:

A = ____ cm / B = ____ cm

Not in combination with Body Contour



Folding backrest

A9011400010 € 323.20

7-way angle adjustable, with integrated aluminum side parts



Back equipment

If not specified, the „Adaptable back with back padding“ will be delivered.

Adaptable back with back padding

A9000802010-K001

Cover: Black, breathable, open surface
Equipment with back cross bar



Adaptable back with back padding „light“

A9000802320-K001

Cover: Black, breathable, open surface
Cushion thickness: 9 mm
Equipment with back cross bar

A9000802320-K002

Cover: Black, breathable, open surface
Cushion thickness: 9 mm
Equipment without back cross bar



Adaptable back with back padding „Comfort“

A9000815101-K € 131.40

Cover: Black, breathable, open surface
Front padding: 20 mm
Rear padding: 5 mm
Equipment with back cross bar

Seat upholstery

If not specified, the „Body Contour“ seat upholstery will be delivered.

Body Contour


A9000801025-K

approved up to a body weight of 120 kg
not in combination with side part with welded tube bracket (item no. A9011400400-001 and A9011400400-002)





Open belt system

A9000801010-K € 79.90



Seat cushion



Seat cushion for individual frame geometry in rectangular contour
Inlet: Standard

A9000805190-K € 164.40

Seat cushion for individual frame geometry with radius contour
Inlet: Standard

A9000805320-K € 206.10

Seat cushion height:

3 cm 5 cm
 4 cm 6 cm

Seat cushion for individual frame geometry in rectangular contour
Inlet: Anti-Decubitus

A9000805210-K € 484.40

Seat cushion for individual frame geometry with radius contour
Inlet: Anti-Decubitus

A9000805310-K € 528.70

Seat cushion height:

4 cm 6 cm

Cover:

Black, breathable, open surface

Black, breathable, closed surface

If not specified, the cover material „Black, breathable, open surface“ will be delivered.

Seat cushion is provided

Wheel camber

12° drive wheel camber

A9011202000-012

15° drive wheel camber

A9011202000-015

18° drive wheel camber

A9011202000-018

20° drive wheel camber

A9011202000-020

Footrest

If not specified, the „continuous footrest“ will be delivered.

Continuous footrest

A9016601010-K001
with tubular resting bracket
height adjustable +/- 3 cm



Continuous footrest with rubber coated footplate

A9011610000-K € 120.90
height adjustable +/- 3 cm
Cover: Edged at the front, screwed
Depth: 27 cm



Continuous footrest with footrest cover „box shape“

A9011610020 € 214.50
rubber coated, height adjustable +/- 3 cm,
Cover:
Width: 20 cm,
Length: 31,5 cm + 1,5 cm adjustable



Ram bar

Offence bar with radius arch

A9011100500 € 351.50



Offence bar with trapezoid shape

A9011100501 € 351.50



Defence bar, removable

A9011100502 € 351.50



Rear ram bar

- A9011100600 € 292.00
designed with a curved tube



Adaptations

Central anti-tipping support oversized permanently welded

- Anti-tipping fork with screwed axle A9011300378 € 358.50
- Anti-tipping fork removable via floating axle A9011300374 € 401.90




2 anti-tipping supports with oversized tube

- Anti-tipping fork with screwed axle A9011300170 € 464.50
- Anti-tipping fork removable via floating axle A9011300270 € 586.80



Central anti-tipping support with oversized tube, length adjustable (tennis)

- A9011300370 € 441.60
with permanently mounted tipping wheel rotary fork



central anti-tipping support 4 Dance

- A9000300010 € 517.40
removable via quick release lever with rotating caster fork and skate wheel Ø 80 mm adjustable ground clearance and angle
Colour: black



3" caster wheel

3" skate wheel (80 mm)

- A9000501110 € 70.80
puncture-proof



Drive wheels

If not specified, the „drive wheel aluminum rim 36 spokes“ will be delivered.

Drive wheel 36 spoke aluminum rim

24" (ETRTO 540)
A9000203000-010
with Liteline hub (left image),
radially spoked
up to 120 kg body weight

24" (ETRTO 540)
A9000203000-020
with aluminum front hub (right image),
radially spoked
over 120 kg body weight



Drive wheel double rim 36-spokes


24" (ETRTO 540)
A9000203100 € 157.80

25" (ETRTO 559)
A9000203200 € 157.80

26" (ETRTO 590)
A9000203300 € 157.80

28" (ETRTO 622)
A9000203310 € 157.80

with aluminum front hub, radially spoked
For use with high loads and body weight
>120 kg



Spinergy wheel „Spox Sport“ with 24 PBO spokes; rim, hub and spokes in black

24" (ETRTO 540)
A9000203550 € 784.10


25" (ETRTO 559)
A9000201559 € 784.10

26" (ETRTO 590)
A9000203555 € 784.10



Handrims

x = dimension handrim attachment



If not specified, the handrim is mounted wide (19 mm).

Stainless steel Ø 20 mm

24" handrim for
24" wheel and 25" wheel
A9000205010-004

26" handrim for
26" wheel
A9000205010-006

28" handrim for
28" wheel
A9000205010-008



Aluminum Ø 19 mm silver anodised

24" handrim for
24" wheel and 25" wheel
A9000205019-004



Aluminum Ø 19 mm polished and hard-anodised

24" handrim for
24" wheel and 25" wheel € 164.90
A9000205030-004

26" handrim for
26" wheel € 164.90
A9000205030-006



Titanium Ø 19 mm

24" handrim for
24" wheel and 25" wheel € 335.60
A9000205040-004

26" handrim for
26" wheel € 335.60
A9000205040-006



Silicone handrim cover

for handrim 24" € 232.80
A9000205100-004

for handrim 26" € 232.80
A9000205100-006

Surface without structure,
compatible with handrim Ø 19 mm or
Ø 20 mm.



Spoke guard

Finger- and spoke guard „transparent“

for 24" wheel
A9000206300-004 € 207.10

for 26" wheel
A9000206300-006 € 207.10

Not in combination with Spinergy drive wheels.



Tyres

If not specified, the „Schwalbe RightRun in black“ tyres will be delivered.

Schwalbe „RightRun“ High-pressure tyres

24x1" (540-25)
A9000204001-004

Colour: Black
7.5 bar



Schwalbe „RightRun“ High-pressure tyres

24x1" (540-25)
A9000204000-004

25x1" (559-25)
A9000204000-005

26x1" (590-25)
A9000204000-006

Colour: Grey, black wall
7.5 bar



28x1" high-pressure lightweight

A9000204000-008
(ETRTO 622)
Colour: Grey, black wall
7.5 bar

Accessories

Hand air compressor with pressure gauge

A9000901100 € 106.30
up to 14 bar



Transport bag for drive wheels

for 24" wheel
A9000903530-004 € 77.80
Ø 76 cm x 2 cm

for 26" wheel
A9000903530-006 € 77.80
Ø 83 cm x 2 cm



Separate order sheets

Category „Rehabilitation technology accessories“

e.g. Seat & backrest systems
can be found at:

www.proactiv-gmbh.com/downloads/

Services

Consultation, clarification of the installation situation and measurement through PRO ACTIV

A9990000000 € 354.50

Handover and instruction on the wheelchair by PRO ACTIV

A9990000006 € 354.50

Notes and remarks:

Location, Date

Signature End customer/dealer

The signature field is intended for the consultant and insured person. This confirms that the intended equipment is intended for the insured person.
- 10 - 2026-03-03 / The prices are plus statutory value-added tax and shipping costs = included in the basic price = Position selected
Order sheet SPEEDY F1 Rev. 10 UK