



☐ **Quotation date:** _____

Consultation was carried out by

Staff at the medical supplies shop: _____

PRO ACTIV field representative: _____

☐ **Order date:** _____

According to GTC PRO ACTIV Reha-Technik GmbH

Medical suppliers/company address:

Cust. no.: _____

Name: _____

Street: _____

Postcode / city: _____

Phone: _____

Fax: _____

E-Mail: _____

Contact: _____

☐ **Invoicing address, other**(see below)

Cust. no.: _____

Name: _____

Street: _____

Postcode / city: _____

☐ **Shipping address, other** (see below)

Cust. no.: _____

Name: _____

Street: _____

Postcode / city: _____

Customer name/commission:

Surname, forename: _____

Street: _____

Postcode / city: _____

Phone: _____

E-Mail: _____

D.O.B: _____

Paying organisation: _____

Indication: _____

NJ1 e-assistant

Up to 120 kg trailer load and 10 kg load capacity

4 different basic models for different purposes, equipped with:

- Docking plate for quick coupling to the "Quick-Snap adapter"
- A-shaped handles with integrated switching and brake fittings
- hydraulic disk brake system, right side operation
- additional brake, operation on the handle on the left with locking bar for use as parking brake
- chain protection tube
- starting assistant (throttle) for travelling without moving the crank Up to 6 km/h
- parking stand
- Lighting set
- mudguard
- colour display sMMi Z15 Connect
- drinking bottle with holder
- black drive unit

Neodrives

Torque 12/40 Nm, power 250 Watt

Wheel size:

☐ 20" (ETRTO 406) A9116100011-020 € 6,750.00

☐ 24" (ETRTO 507) A9116100011-024 € 6,750.00

Please note: When purchasing an NJ1 e-assistant in combination with a wheelchair, the discounted VAT rate applies for the complete unit.

Body size

Body height: _____ cm

Upper body height: _____ cm

Weight: _____ kg

Information essential to process the order.

Neodrives, rechargeable batteries

The ranges specified depend on the range profile, user weight and speed.

Neodrives G2 RR rechargeable battery

36V 14.25 Ah 513 Wh

Range approx. 40-70 km

☐ A9116900106 € 80.00

Assembly to the left on the fork



Currently not available

Neodrives G2 RR double rechargeable battery

36V 14.25 Ah 513 Wh

Range approx. 80-140 km

☐ A9116900108 € 1,326.10

Assembly on both sides of the fork, re-plugging unit mounted on the bottom bracket support to allow the simple replacement of the rechargeable batteries via a magnetic connector



Currently not available

Neodrives DT rechargeable battery

36V 16Ah 576Wh

Range approx. 45-75 km

☐ A9116900116

Assembly to the left on the fork



Neodrives DT double rechargeable battery

36V 16Ah 576Wh

Range approx. 85-150 km

☐ A9116900118 € 1,228.20

Assembly on both sides of the fork, re-plugging unit mounted on the bottom bracket support to allow the simple replacement of the rechargeable batteries via a magnetic connector



Dérailleur

Shimano Deore 1x11 Rapid Fire 11-speed dérailleur

☐ A9110200032

Mechanical dérailleur gear with operation on the handle. 11-speed cassette 11-42 teeth

Chain ring selection:

- ☐ 32 teeth recommended for 24" drive wheel
- ☐ 38 teeth recommended for 20" drive wheel
- ☐ 42 teeth for lower crank frequency or driving at higher speeds



Shimano Deore XT 1x11 Rapid Fire 11-speed dérailleur

☐ A9110200034 € 182.90

Mechanical dérailleur gear with operation on the handle. 11-speed cassette 11-42 teeth

Chain ring selection:

- ☐ 32 teeth recommended for 24" drive wheel
- ☐ 38 teeth recommended for 20" drive wheel
- ☐ 42 teeth for lower crank frequency or driving at higher speeds




**Shimano Deore XT 1x11 DI2
with button operation
11-speed dérailleur**

☐ A9110200035 € 1,414.50

Electronic dérailleur with operation on the handle. 11-speed cassette 11-42 teeth

Chain ring selection:

- ☐ 32 teeth recommended for 24" drive wheel
☐ 38 teeth recommended for 20" drive wheel
☐ 42 teeth for lower crank frequency or driving at higher speeds



**Shimano Deore XT 1x11 DI2
With chin operation
11-speed dérailleur**

☐ A9110200033 € 1,552.50

Electronic dérailleur gear with chin operation (tetra equipment). 11-speed cassette 11-42 teeth

Chain ring selection:

- ☐ 32 teeth recommended for 24" drive wheel
☐ 38 teeth recommended for 20" drive wheel
☐ 42 teeth for lower crank frequency or driving at higher speeds



Gear change accessories

**Mountain drive 2.5: 1 gear reduction
for hills**

☐ A9100200003 € 759.00

Planetary gear set for doubling of the number of gears permits switching to the gear reduction for hills when stationary. With 38 teeth chain wheel as standard. Only in combination with grip size 40, 45, 48 and 53 cm.



Brake system

**V-brake parking brake operated via
locking lever**

☐ A9110300002 € 80.50

mounted on the bottom bracket support, brake handle on the handle, left side is omitted



**Powerful brake system
with brake disc Ø 203 mm**

☐ A9110300105 € 79.10

Recommended from a towed weight of 80 kg or body weight. Not required with hydraulic back pedal disk brake selected



**Double powerful brake system disc
brake Ø 203 mm**

☐ A9110300100* € 330.40

For 24": together with V-brake rim brake as parking brake, can be operated via locking lever. Recommended from a towed weight of 100 kg

**Mechanical backpedal disk brake with
crank release function**

☐ A9110300004 € 702.80

Compatible with dérailleur and free-running shifting hub.

In combination with:

- Grip width 45 cm, 50 cm, 53 cm or 58 cm
 - Pedal crank crank lengths 155 or 175 mm up to 120 kg body weight
- Not in combination with Mountaindrive 2.5:1 gear reduction for hills

Tyres

20x1.75" Schwalbe "Marathon"

☐ A9110105121
(ETRTO 406-47)
Touring tyres



20x2" Schwalbe "Knobby"

(Low availability – only as long as stocks last.
Alternative option: Schwalbe "Marathon")

☐ A9110105126
(ETRTO 406-54)
Profiled tyres



Schwalbe "Marathon Plus"

☐ 20x1.75" (406-47)
A9130105128 € 27.40

☐ 24x1.75" (507-47)
A9110105135 € 27.40
unpuncturable, with anti-puncture belt
Smart Guard



Dimensions

Crank type and grip width

Standard crank

- ☐ Grip width 40 cm
A9110500100-040
- ☐ Grip width 45 cm
A9110500100-045
- ☐ Grip width 50 cm
A9110500100-050

Tetra pedal crank (flared)

With flared crank arms for max. space for upper body

- ☐ Grip width 48 cm (176 mm)
A9110500100-048
- ☐ Grip width 53 cm (226 mm)
A9110500100-053
- ☐ Grip width 58 cm (276 mm)
A9110500100-058

Standard pedal crank

- ☐ Length 155 mm
A9110600100-155
- ☐ Length 175 mm
A9110600100-175
- ☐ Length 195 mm
A9110600100-195



Tetra pedal crank

- ☐ Length 155 mm € 129.50
A9100600105-155
- ☐ Length 175 mm € 129.50
A9100600105-175
- ☐ Individual length
140 mm to 195 mm € 246.30
A9100600100



Fixed bottom bracket support

- ☐ L (X/Y 158/362 mm)
A9110501100-002
- ☐ XL (X/Y 158/410 mm) € 63.50
A9110501100-006
- ☐ XXL (X/Y 170/465 mm) € 63.50
A9110501100-007



Adjustable bottom bracket support

the angle and length can be adjusted
via the clamped rotary point on the fork bridge
and length.

For the pre-setting, please specify the
bottom bracket support size that comes
as close to your desired position.

- ☐ M (300 mm) € 220.90
A9110500001-001
- ☐ L (370 mm) € 220.90
A9110500001-002
- ☐ XL (420 mm) € 220.90
A9110500001-003

If no details are provided, the adjustable
bottom bracket support shall be supplied in XL.
Deliverable up to size XL



Handles

Palm rest on the handle

☐ A9110601200 € 97.80

Padded, bolted, angle-adjustable.
Not in combination with other handles



Tetra hand support, pivotable, adjustable grip distance

☐ A9110601101 € 161.30



Tetra hand support, pivotable, adjustable grip width

☐ A9110601102 € 182.40

Very stiff outer support as round tube ø 20 mm with foam rubber cover
Cannot be pivoted in conjunction with Tetra hand support



Tetra Ergo handle, free-standing

☐ A9100600001 € 156.70

Crank handles arranged vertically
Can be combined with Tetra hand rest, pivotable



Tetra shell handles

☐ A9110601100 € 164.10



Bike accessories

Mudguard and chain guard, black

☐ A9116500005 € 107.00

Only in combination with 20" drive wheel

Mudguard and chain guard unique colour selection

☐ A9116500010 € 236.70

Made from a solid colour palette
Powder paint wheelchair

A9000100170-__ __ __

Colour: _____

Only in combination with 20" drive wheel

LED lighting set for Neodrives System

☐ A9116900002 € 144.40

Supply of the LED front light via the rechargeable battery of the drive system
Back light as plug-on light supplied by a battery
No need for lighting set from basic equipment

Neodrives rapid charger

☐ A9116900006 € 80.20

Replaces the standard charger



Bracket ø 22 mm x 42 mm for fixing of accessories (speedometer, lighting or similar)

☐ A9110903271 Item price € 22.10

Colour: Natural anodised

_____ Items



Manoeuvring rollers for parking stands

☐ A9110900017 € 116.20

Manoeuvring rollers mounted to the parking stand via inlet



Manoeuvring bar mounted on the chain guard

☐ A9110500000 € 160.80

with foam rubber cover

Aluminium luggage rack "cargo systems"

☐ A9110903290 € 179.70

Colour: black
Max. load weight 10 kg on the loading platform,
Max. 15 kg in combination with luggage bags with low loading height.



Compensating weight for rechargeable battery

☐ A9116900000 € 122.20

with drinking bottle and holder. For compensating the weight of the rechargeable battery with only one rechargeable battery, and for improving traction.



Operating satellite mounted on the handle

Assembly: Right ☐ Left ☐

with existing support line (laid along a brake line)

☐ A9116200000 € 96.30

without an existing support line (an additional dummy line is required)

☐ A9116200001 € 136.50

Options for users with reduced function**Quick release clip lever extension for fixed frame adapter**

☐ A9110700004 € 98.90

to make the opening and closing of the quick release clips easier.

Only in combination with rigid frame adapter

**Securing the lever extension with locking bolt**

☐ A9199000084 € 82.90

Can only be used for the lever extension on the quick release clip.

Eccentric connection lever extension for Quick Snap connection

simplifies opening or closing the Quick Snap Adapter. Tetra version

☐ A9110700005 € 70.30


Straight version

Only in combination with rigid frame adapter

☐ A9110700006 € 70.30

Curved version

Only in conjunction with front adapter

**Centring column adaptation assistance**

☐ A9110700002 € 52.40

in order to make inserting the central adapter easier

**Parking stand with operating lever for high operating comfort**

☐ A9110900004 € 55.90

with longer actuating lever and reduced spring force

Chest strap for upper body fixation

☐ Size S (105 cm) A9000800016 € 220.80

☐ Size M (125 cm) A9000800013 € 220.80

☐ Size L (145 cm) A9000800014 € 220.80

☐ Size XL (165 cm) A9000800017 € 220.80

Consisting of elastic webbing material in combination with neoprene with Velcro and operating tab.
Colour: black

**Frame adapters****SPEEDY fixed frame adapter with "Quick Snap Adaptation"**

☐ A9110700141 € 492.50

with quick release clips and centring column mount

Colour: black

**TRAVELER folding frame adapter with "Quick Snap Adaptation"**

☐ A9110701300 € 526.70

with quick release clips and centring column mount on the folding mechanism

Colour: black

**BUDDY fixed frame adapter with "Quick Snap Adaptation"**

☐ A9110701400 € 492.50

with quick release clips and centring column mount on the frame

Colour: black



Wheelchair front adapter with universal clamping and "Quick Snap Adaptation"

☐ A9110701500 € 561.10

Adapter mount on the wheelchair front end
Operation on the right side

Installation situation with PA products:

Frame stem G-shape from SB 31 cm to 54 cm

Frame stem V-shape from SB 30 cm to 50 cm

PLEASE NOTE:

If the tube diameter is not specified, as
standard we shall deliver an adapter for a tube
diameter of 23 - 27 mm.



Please make the following selection:

☐ A9110701525

Tube diameter 23-27 mm

☐ A9110701530

Tube diameter 28-32 mm

☐ A9110701535

Tube diameter 33-37 mm

Clamp lever with 3-way star handle for actuating the frame clamping

☐ A9199701507 € 101.00

mounted on the right and left for easy terminal
operation with limited manual function
Only in conjunction with front adapter



Wheelchair-dependent accessories

Extended bike plate

☐ For SPEEDY A1/A2
A9110100100-002 € 352.20

☐ For BUDDY
A9110100100-003 € 352.20

☐ For TRAVELER
A9110100100-001 € 352.20

4° camber integrated
Track can be adjusted independently



Wheelbase extension via "grid mechanism with pulling operation"

☐ For SPEEDY 4all /
4teen / 4you /
TRAVELER 4you € 537.50
A9110102001

removable, with continuous axle tube
4°/ 6° camber integrated
Track can be adjusted independently



Wheelbase extension, angle-adjustable and removable

☐ For SPEEDY F2/F4
A9015100005 € 677.30

with continuous axle tube,
height adjustment via gear joint
adjustable in 6° angular increments
wheel position displacement approx. 14 cm
4°/6° camber integrated, adjustable track
Attachment possibility must be checked by
PRO ACTIV.

Only in combination with base plate.
Up to 120 kg load capacity.



Services

Handover and instruction for the
handbike/traction device by the dealer
incl. mounting the frame adapter

on a PRO ACTIV wheelchair

☐ A9990000008

Prod. no.: _____

on a third-party product

☐ A9990000005

Manufacturer: _____

Model: _____

Frame tube ø ____ mm

Notes and remarks:

Location, Date

Signature End customer/dealer

The signature field is intended for the consultant and insured person. This confirms that the intended equipment is intended for the insured person.

- 7 - 2023-02-16 / The prices are plus statutory value-added tax and shipping costs ☐ = Included in the basic price ☒ = Position required
Order sheet NJ1 e-assistant Rev. 03 UK