

Im Hofstätt 11
 D-72359 Dotternhausen
 Tel +49 (0)7427 9480-0
 Fax +49 (0)7427 9480-7025
 e-mail info@proactiv-gmbh.de
 www.proactiv-gmbh.de



SPEEDY 4teen

Active easy running rigid frame wheelchair

Maximum user weight limit 100 kg

A9050100001 € 4,080.00

Quotation date: _____
Consultation was carried out by
Staff at the medical supplies shop: _____
Enquiry no.: _____

Order date: _____
 according to GTC PRO ACTIV Reha-Technik GmbH
Order no.: _____

Medical suppliers / company address:

Cust. no.: _____
 Name: _____
 Street: _____
 Postcode / city: _____
 Phone: _____
 Fax: _____
 E-Mail: _____
 Contact: _____

Shipping address: same as company address
 Other (see below)

Name: _____
 Street: _____
 Postcode / city: _____


Customer name / commission:

ID / Comm. No.: _____
 Surname, forename: _____
 Street: _____
 Postcode / city: _____
 Phone: _____
 Fax: _____
 E-Mail: _____
 D.o.B.: _____
 Paying organisation: _____

Invoicing address: same as company address
 Other (see below)

Cust. no.: _____
 Name: _____
 Street: _____
 Postcode / city: _____
 Phone: _____
 Fax: _____
 E-Mail: _____
 Contact: _____


Body size:



Body height: _____ cm
 Weight: _____ kg
 Lower leg length: _____ cm
 (anatomic dimension: hollow of the knee to the foot sole)
 Indication: _____

This information is necessary for checking. If the information is missing, a check of the relevant items cannot be made!

Basic size



Size L (Frame length R 39 cm to arch)

Seat width SB	<input type="checkbox"/> 26	<input type="checkbox"/> 27	<input type="checkbox"/> 28	<input type="checkbox"/> 29	<input type="checkbox"/> 30	<input type="checkbox"/> 31	<input type="checkbox"/> 32				
	<input type="checkbox"/> 33	<input type="checkbox"/> 34	<input type="checkbox"/> 35	<input type="checkbox"/> 36	<input type="checkbox"/> 37	<input type="checkbox"/> 38					
Seat depth ST	<input type="checkbox"/> 32	<input type="checkbox"/> 33	<input type="checkbox"/> 34	<input type="checkbox"/> 35	<input type="checkbox"/> 36	<input type="checkbox"/> 37	<input type="checkbox"/> 38	<input type="checkbox"/> 39	<input type="checkbox"/> 40		
	<input type="checkbox"/> 20	<input type="checkbox"/> 21	<input type="checkbox"/> 22	<input type="checkbox"/> 23	<input type="checkbox"/> 24	<input type="checkbox"/> 25	<input type="checkbox"/> 26	<input type="checkbox"/> 27	<input type="checkbox"/> 28	<input type="checkbox"/> 29	
Back height RH	<input type="checkbox"/> 30	<input type="checkbox"/> 31	<input type="checkbox"/> 32	<input type="checkbox"/> 33	<input type="checkbox"/> 34	<input type="checkbox"/> 35	<input type="checkbox"/> 36	<input type="checkbox"/> 37	<input type="checkbox"/> 38	<input type="checkbox"/> 39	
	<input type="checkbox"/> 40	<input type="checkbox"/> 41	<input type="checkbox"/> 42	<input type="checkbox"/> 43	<input type="checkbox"/> 44	<input type="checkbox"/> 45					
	<input type="checkbox"/> 20	<input type="checkbox"/> 21	<input type="checkbox"/> 22	<input type="checkbox"/> 23	<input type="checkbox"/> 24	<input type="checkbox"/> 25	<input type="checkbox"/> 26	<input type="checkbox"/> 27	<input type="checkbox"/> 28	<input type="checkbox"/> 29	
Lower leg length USL	<input type="checkbox"/> 30	<input type="checkbox"/> 31	<input type="checkbox"/> 32	<input type="checkbox"/> 33	<input type="checkbox"/> 34	<input type="checkbox"/> 35	<input type="checkbox"/> 36	<input type="checkbox"/> 37	<input type="checkbox"/> 38	<input type="checkbox"/> 39	
	<input type="checkbox"/> 40	<input type="checkbox"/> 41	<input type="checkbox"/> 42	<input type="checkbox"/> 43	<input type="checkbox"/> 44	<input type="checkbox"/> 45	<input type="checkbox"/> 46	<input type="checkbox"/> 47	<input type="checkbox"/> 48	<input type="checkbox"/> 49	
	<input type="checkbox"/> 20	<input type="checkbox"/> 21	<input type="checkbox"/> 22	<input type="checkbox"/> 23	<input type="checkbox"/> 24	<input type="checkbox"/> 25	<input type="checkbox"/> 26	<input type="checkbox"/> 27	<input type="checkbox"/> 28	<input type="checkbox"/> 29	
Seat height back SHh	<input type="checkbox"/> 36	<input type="checkbox"/> 37	<input type="checkbox"/> 38	<input type="checkbox"/> 39	<input type="checkbox"/> 40	<input type="checkbox"/> 41	<input type="checkbox"/> 42	<input type="checkbox"/> 43	<input type="checkbox"/> 44	with 24" wheels	
Seat height front SHv	<input type="checkbox"/> 40	<input type="checkbox"/> 41	<input type="checkbox"/> 42	<input type="checkbox"/> 43	<input type="checkbox"/> 44	<input type="checkbox"/> 45	<input type="checkbox"/> 46	<input type="checkbox"/> 47	with 5" caster wheels		
Tipping point KP (inner edge of rear tube to centre of drive wheel axle)	<input type="checkbox"/> 1	<input type="checkbox"/> 2	<input type="checkbox"/> 3	<input type="checkbox"/> 4	<input type="checkbox"/> 5	<input type="checkbox"/> 6	<input type="checkbox"/> 7	<input checked="" type="checkbox"/> 8	<input type="checkbox"/> 9	<input type="checkbox"/> 10	

Size XL (Frame length R 43 cm to arch)

Seat width SB	<input type="checkbox"/> 26	<input type="checkbox"/> 27	<input type="checkbox"/> 28	<input type="checkbox"/> 29	<input type="checkbox"/> 30	<input type="checkbox"/> 31	<input type="checkbox"/> 32				
	<input type="checkbox"/> 33	<input type="checkbox"/> 34	<input type="checkbox"/> 35	<input type="checkbox"/> 36	<input type="checkbox"/> 37	<input type="checkbox"/> 38					
Seat depth ST	<input type="checkbox"/> 36	<input type="checkbox"/> 37	<input type="checkbox"/> 38	<input type="checkbox"/> 39	<input type="checkbox"/> 40	<input type="checkbox"/> 41	<input type="checkbox"/> 42	<input type="checkbox"/> 43			
	<input type="checkbox"/> 20	<input type="checkbox"/> 21	<input type="checkbox"/> 22	<input type="checkbox"/> 23	<input type="checkbox"/> 24	<input type="checkbox"/> 25	<input type="checkbox"/> 26	<input type="checkbox"/> 27	<input type="checkbox"/> 28	<input type="checkbox"/> 29	
Back height RH	<input type="checkbox"/> 30	<input type="checkbox"/> 31	<input type="checkbox"/> 32	<input type="checkbox"/> 33	<input type="checkbox"/> 34	<input type="checkbox"/> 35	<input type="checkbox"/> 36	<input type="checkbox"/> 37	<input type="checkbox"/> 38	<input type="checkbox"/> 39	
	<input type="checkbox"/> 40	<input type="checkbox"/> 41	<input type="checkbox"/> 42	<input type="checkbox"/> 43	<input type="checkbox"/> 44	<input type="checkbox"/> 45	<input type="checkbox"/> 46	<input type="checkbox"/> 47	<input type="checkbox"/> 48		
	<input type="checkbox"/> 20	<input type="checkbox"/> 21	<input type="checkbox"/> 22	<input type="checkbox"/> 23	<input type="checkbox"/> 24	<input type="checkbox"/> 25	<input type="checkbox"/> 26	<input type="checkbox"/> 27	<input type="checkbox"/> 28	<input type="checkbox"/> 29	<input type="checkbox"/> 30
Lower leg length USL	<input type="checkbox"/> 31	<input type="checkbox"/> 32	<input type="checkbox"/> 33	<input type="checkbox"/> 34	<input type="checkbox"/> 35	<input type="checkbox"/> 36	<input type="checkbox"/> 37	<input type="checkbox"/> 38	<input type="checkbox"/> 39	<input type="checkbox"/> 40	<input type="checkbox"/> 41
	<input type="checkbox"/> 42	<input type="checkbox"/> 43	<input type="checkbox"/> 44	<input type="checkbox"/> 45	<input type="checkbox"/> 46	<input type="checkbox"/> 47	<input type="checkbox"/> 48	<input type="checkbox"/> 49	<input type="checkbox"/> 50	<input type="checkbox"/> 51	<input type="checkbox"/> 52
	<input type="checkbox"/> 20	<input type="checkbox"/> 21	<input type="checkbox"/> 22	<input type="checkbox"/> 23	<input type="checkbox"/> 24	<input type="checkbox"/> 25	<input type="checkbox"/> 26	<input type="checkbox"/> 27	<input type="checkbox"/> 28	<input type="checkbox"/> 29	<input type="checkbox"/> 30
Seat height back SHh	<input type="checkbox"/> 36	<input type="checkbox"/> 37	<input type="checkbox"/> 38	<input type="checkbox"/> 39	<input type="checkbox"/> 40	<input type="checkbox"/> 41	<input type="checkbox"/> 42	<input type="checkbox"/> 43	<input type="checkbox"/> 44	with 24" wheels	
Seat height front SHv	<input type="checkbox"/> 40	<input type="checkbox"/> 41	<input type="checkbox"/> 42	<input type="checkbox"/> 43	<input type="checkbox"/> 44	<input type="checkbox"/> 45	<input type="checkbox"/> 46	<input type="checkbox"/> 47	with 5" caster wheels		
Tipping point KP (inner edge of rear tube to centre of drive wheel axle)	<input type="checkbox"/> 1	<input type="checkbox"/> 2	<input type="checkbox"/> 3	<input type="checkbox"/> 4	<input type="checkbox"/> 5	<input type="checkbox"/> 6	<input type="checkbox"/> 7	<input checked="" type="checkbox"/> 8	<input type="checkbox"/> 9	<input type="checkbox"/> 10	

Size XXL (Frame length R 47 cm to arch)

Seat width SB	<input type="checkbox"/> 26	<input type="checkbox"/> 27	<input type="checkbox"/> 28	<input type="checkbox"/> 29	<input type="checkbox"/> 30	<input type="checkbox"/> 31	<input type="checkbox"/> 32				
	<input type="checkbox"/> 33	<input type="checkbox"/> 34	<input type="checkbox"/> 35	<input type="checkbox"/> 36	<input type="checkbox"/> 37	<input type="checkbox"/> 38					
Seat depth ST	<input type="checkbox"/> 40	<input type="checkbox"/> 41	<input type="checkbox"/> 42	<input type="checkbox"/> 43	<input type="checkbox"/> 44	<input type="checkbox"/> 45	<input type="checkbox"/> 46	<input type="checkbox"/> 47	<input type="checkbox"/> 48		
	<input type="checkbox"/> 20	<input type="checkbox"/> 21	<input type="checkbox"/> 22	<input type="checkbox"/> 23	<input type="checkbox"/> 24	<input type="checkbox"/> 25	<input type="checkbox"/> 26	<input type="checkbox"/> 27	<input type="checkbox"/> 28	<input type="checkbox"/> 29	
Back height RH	<input type="checkbox"/> 30	<input type="checkbox"/> 31	<input type="checkbox"/> 32	<input type="checkbox"/> 33	<input type="checkbox"/> 34	<input type="checkbox"/> 35	<input type="checkbox"/> 36	<input type="checkbox"/> 37	<input type="checkbox"/> 38	<input type="checkbox"/> 39	
	<input type="checkbox"/> 40	<input type="checkbox"/> 41	<input type="checkbox"/> 42	<input type="checkbox"/> 43	<input type="checkbox"/> 44	<input type="checkbox"/> 45	<input type="checkbox"/> 46	<input type="checkbox"/> 47	<input type="checkbox"/> 48		
	<input type="checkbox"/> 20	<input type="checkbox"/> 21	<input type="checkbox"/> 22	<input type="checkbox"/> 23	<input type="checkbox"/> 24	<input type="checkbox"/> 25	<input type="checkbox"/> 26	<input type="checkbox"/> 27	<input type="checkbox"/> 28	<input type="checkbox"/> 29	<input type="checkbox"/> 30
Lower leg length USL	<input type="checkbox"/> 31	<input type="checkbox"/> 32	<input type="checkbox"/> 33	<input type="checkbox"/> 34	<input type="checkbox"/> 35	<input type="checkbox"/> 36	<input type="checkbox"/> 37	<input type="checkbox"/> 38	<input type="checkbox"/> 39	<input type="checkbox"/> 40	<input type="checkbox"/> 41
	<input type="checkbox"/> 42	<input type="checkbox"/> 43	<input type="checkbox"/> 44	<input type="checkbox"/> 45	<input type="checkbox"/> 46	<input type="checkbox"/> 47	<input type="checkbox"/> 48	<input type="checkbox"/> 49	<input type="checkbox"/> 50	<input type="checkbox"/> 51	<input type="checkbox"/> 52
	<input type="checkbox"/> 20	<input type="checkbox"/> 21	<input type="checkbox"/> 22	<input type="checkbox"/> 23	<input type="checkbox"/> 24	<input type="checkbox"/> 25	<input type="checkbox"/> 26	<input type="checkbox"/> 27	<input type="checkbox"/> 28	<input type="checkbox"/> 29	<input type="checkbox"/> 30
Seat height back SHh	<input type="checkbox"/> 39	<input type="checkbox"/> 40	<input type="checkbox"/> 41	<input type="checkbox"/> 42	<input type="checkbox"/> 43	<input type="checkbox"/> 44	<input type="checkbox"/> 45	<input type="checkbox"/> 46	<input type="checkbox"/> 47	with 24" wheels	
Seat height front SHv	<input type="checkbox"/> 43	<input type="checkbox"/> 44	<input type="checkbox"/> 45	<input type="checkbox"/> 46	<input type="checkbox"/> 47	<input type="checkbox"/> 48	<input type="checkbox"/> 49	<input type="checkbox"/> 50	with 5" caster wheels		
Tipping point KP (inner edge of rear tube to centre of drive wheel axle)	<input type="checkbox"/> 1	<input type="checkbox"/> 2	<input type="checkbox"/> 3	<input type="checkbox"/> 4	<input type="checkbox"/> 5	<input type="checkbox"/> 6	<input type="checkbox"/> 7	<input checked="" type="checkbox"/> 8	<input type="checkbox"/> 9	<input type="checkbox"/> 10	

Frame colour - One-colour design

Information on current colour designs and example images can be found at: <https://www.proactiv-gmbh.com/colours/>

Frame colour: Selection from colour palette

Frame colour will be specified in the order

A9000100170-000

Raw, blasted, not coated

A9000100170-020

Note: In combination with moisture e.g. hand perspiration, can affect dark aluminum discoloration.

Powder colours

white RAL 9016

gloss A9000100170-023

matt A9000100170-123

€ 43.40

brilliant silver

gloss A9000100170-011

anthracite RAL 7016

gloss A9000100170-001

matt A9000100170-101

€ 43.40

black RAL 9005

gloss A9000100170-003

matt A9000100170-019

€ 43.40

pure orange RAL 2004

gloss A9000100170-022

matt A9000100170-122

€ 43.40

pink magenta RAL 4010

gloss A9000100170-037

matt A9000100170-137

€ 43.40

pastel blue RAL 5024

gloss A9000100170-036

matt A9000100170-136

€ 43.40

lemon green

gloss A9000100170-005

matt A9000100170-105

€ 43.40

yellow RAL 1023

gloss A9000100170-046

matt A9000100170-146

€ 43.40

red RAL 3001

gloss A9000100170-047

matt A9000100170-147

€ 43.40

lagoon blue

gloss A9000100170-052

matt A9000100170-152

€ 43.40

athletic grey RAL 7046

gloss A9000100170-050

matt A9000100170-150

€ 43.40

muddy grey RAL 7033

gloss A9000100170-051

matt A9000100170-151

€ 43.40

peach RAL 3022

gloss A9000100170-038

matt A9000100170-138

€ 43.40

Powder colours „dormant“

deep black effect silver

gloss A9000100170-007

sky blue

gloss A9000100170-017

navy blue

gloss A9000100170-045

matt A9000100170-145

€ 43.40

electric green

gloss A9000100170-015

matt A9000100170-115

€ 43.40

violet

gloss A9000100170-044

fire red

gloss A9000100170-018

matt A9000100170-118

€ 43.40

Frame colour - One-colour design special colours „liquid coating“

€ 349.60

storm (ice blue)

gloss A9000100170-500

matt A9000100170-501

prairie (rose gold)

gloss A9000100170-504

matt A9000100170-505

berry

gloss A9000100170-526

matt A9000100170-527

blue ocean

gloss A9000100170-522

matt A9000100170-523

emerald

gloss A9000100170-514

matt A9000100170-515

chili

gloss A9000100170-516

matt A9000100170-517

Frame colour - Two-colour design

Note: In the two-colour design, both colours are the same gloss level as the design colour.

Two-colour design with one base colour (on the rear frame) and one design colour (on the front frame) € 268.50

Please select a colour below:

Base colour (on the rear frame)

white RAL 9016 <input type="checkbox"/> A9000100170-200	brilliant silver <input type="checkbox"/> A9000100170-201	anthracite RAL 7016 <input type="checkbox"/> A9000100170-202
black RAL 9005 <input type="checkbox"/> A9000100170-203	pastel blue RAL 5024 <input type="checkbox"/> A9000100170-206	lemon green <input type="checkbox"/> A9000100170-207
yellow RAL 1023 <input type="checkbox"/> A9000100170-208	lagoon blue <input type="checkbox"/> A9000100170-212	athletic grey RAL 7046 <input type="checkbox"/> A9000100170-210
muddy grey RAL 7033 <input type="checkbox"/> A9000100170-211		

Design colour (on the front frame)

black RAL 9005 <input type="checkbox"/> gloss A9000100170-259 <input type="checkbox"/> matt A9000100170-260	lemon green <input type="checkbox"/> gloss A9000100170-250 <input type="checkbox"/> matt A9000100170-251	chili <input type="checkbox"/> gloss A9000100170-253 <input type="checkbox"/> matt A9000100170-254
pure orange RAL 2004 <input type="checkbox"/> gloss A9000100170-255 <input type="checkbox"/> matt A9000100170-256	emerald <input type="checkbox"/> gloss A9000100170-274 <input type="checkbox"/> matt A9000100170-275	storm (ice blue) <input type="checkbox"/> gloss A9000100170-264 <input type="checkbox"/> matt A9000100170-265
prairie (rose gold) <input type="checkbox"/> gloss A9000100170-266 <input type="checkbox"/> matt A9000100170-267	berry <input type="checkbox"/> gloss A9000100170-284 <input type="checkbox"/> matt A9000100170-285	blue ocean <input type="checkbox"/> gloss A9000100170-280 <input type="checkbox"/> matt A9000100170-281

Colour - Attachment parts

Note: If not specified, the caster fork and bearing block will be mounted in black.

Drive wheel: aluminum rim, hub and spokes

Black A9000100171-003

Caster fork and bearing block

Silver A9000100173-011
 Black A9000100173-003

Colour design of attachment parts according to customer request:

Drive wheel: aluminum rim and design inlet¹ according to colour selection / hub and spokes black

frame colour A9000100171-999 € 249.60

Selection from one-colour design (different from the frame colour) A9000100171-998 € 288.90

Colour: _____
A9000100170-

Caster wheel: 5" aluminum rim² (6-spoke design)

frame colour A9000100174-999 € 75.20

Selection from one-colour design (different from the frame colour) A9000100174-998 € 87.10

Colour: _____
A9000100170-

Special colour request (Frame colour - One-colour design)

Frame colour for single purchasing Selection from RAL colour palette

A9000100170-100 € 414.40

gloss level

gloss
 matt

colour designation

RAL number

Note: Attachments **not** available in design colour options


¹ Only in combination with Liteline hub

² Only in combination with caster wheel 5" aluminum rim 6-spoke-design good shock absorption A900050121*, as well as with suspension fork 5-spoke-design A9000500160*

Frame


Straight front end

A9050100010-K*
Knee angle 100°



Front end abducted shape

A9050100030-K* € 383.10
Knee angle 100°
2.5 cm per side
Abduction starts at:
ST 31 cm, Size L
ST 35 cm, Size XL
ST 39 cm, Size XXL
Length of the abduction up to opened
dimensions approx. 55 mm



Caster wheel bearing block

Base plate for bearing block welded

A9050500000
X (see drawing) for frame size
Size L: 63 mm
Size XL: 74 mm
Size XXL: 80 mm



Push handles

Aluminium push handles fixed in back tube

A9050401050
from RH 23 cm



Safety push handles with continuous height adjustment

A9050401060* € 86.30
mount in the back tube
from RH 30 cm with adaptable back
from RH 35 cm with lordosis support
from RH 43 cm with ergonomic backshell



Safety push handles back-positioned

A9050401075* € 86.30


continuous height-adjustable
push handle holder back-positioned

Length 45 cm

Length 55 cm

If not specified, the „length 45 cm” will be
delivered.



<p>Aluminum push handles, back-positioned with side handle offset 145 mm</p> <p><input type="checkbox"/> A9050401080-K* € 163.90</p> <p>continuous height-adjustable Safety push handle holder mounted on the back cross bar Length: 55 cm Up to SB 38 cm or body weight 60 kg</p> 	<p>Aluminium central push handle Foldable handle area</p> <p><input type="checkbox"/> A9050401086-K € 188.90</p> <p>height-adjustable and foldable in the handle area via toothed joint, holder welded on the middle on the back cross bar. Up to SB 38 cm or body weight 50 kg Not in combination with headrest mount Please note: Not suitable for carrying the wheelchair — use the back crossbar for this purpose.</p> 	<p>Central push handle for brake handles Foldable grip area</p> <p><input type="checkbox"/> A9050401070-K € 232.80</p> <p>closed, with continuous height adjustment, with foam rubber handle cover Up to SB 38 cm or body weight 50 kg Not in combination with headrest mount Only in combination with parking lever brake for accompanying person or drum brake.</p> 
<p>Aluminium push handle closed, foldable via toothed joint</p> <p><input type="checkbox"/> A9050401069-K € 476.20</p> <p>continuous height-adjustable push handles shifted to the back Up to SB 38 cm or body weight 60 kg</p> 	<p>Equipment without push handles</p> <p><input checked="" type="checkbox"/> A9000401000</p> 	
<p>Additional option</p>		
<p>Back tube with lordosis support</p> <p><input type="checkbox"/> A9000400005-K € 114.40</p> <p><input type="checkbox"/> 1 cm <input type="checkbox"/> 2 cm <input type="checkbox"/> 3 cm</p>		
<p>Backrest</p>		
<p>Fixed backrest</p> <p><input checked="" type="checkbox"/> A9050400020</p> <p>screwed to aluminium side parts Back angle 90° to the floor Backrest cannot fold down</p>	<p>Folding backrest</p> <p><input type="checkbox"/> A9000400024* € 306.30</p> <p>7-way angle adjustable, with integrated plastic-coated aluminum side parts</p> 	<p>Backrest folding function blocked</p> <p><input type="checkbox"/> A9000400038* € 312.30</p> <p>Locking bolt for angle adjustment replaced by screw connection Back angle 90° to the floor with integrated plastic-coated aluminum side parts</p> 

Back equipment

If not specified, the „Adaptable back with back padding“ will be delivered.

Adaptable back with back padding

A9000802010-K001

Cover: Black, breathable, open surface
Equipment with back cross bar



Adaptable back with back padding „light“

A9000802320-K001

Cover: Black, breathable, open surface
Cushion thickness: 9 mm
Equipment with back cross bar

A9000802320-K002

Cover: Black, breathable, open surface
Cushion thickness: 9 mm
Equipment without back cross bar



Adaptable back with back padding „Comfort“

A9000815101-K **€ 131.40**

Cover: Black, breathable, open surface
Front padding: 20 mm
Rear padding: 5 mm
Equipment with back cross bar

Back crossbar flared

by 20 mm
A9000400835-020 **€ 109.30**

by 30 mm
A9000400835-030 **€ 109.30**

by 50 mm
A9000400835-050 **€ 140.80**

by 80 mm
A9000400835-080 **€ 155.30**

Only in combination with back-positioned push handles



Equipment without back upholstery and back padding or back system

A9000802011 **€ - 53.00**


Equipment with back cross bar

Seat upholstery

If not specified, the „Body Contour“ seat upholstery will be delivered.

Body Contour

A9000801025-K



Seat plate aluminum

A9000801212-K **€ 292.70**



screwed onto the frame tubes
Colour: anthracite matt
(If no colour is specified under „Select colour“)



Equipment without seat upholstery

A9000801050

Seat cushion



Seat cushion in rectangular contour Inlet: Standard

A9000805180-K* € 118.30

Seat cushion with radius contour Inlet: Standard

A9000800300-K* € 163.10

Seat cushion height:

3 cm 5 cm
 4 cm 6 cm

Seat cushion in rectangular contour Inlet: Anti-Decubitus

up to 30x30 cm
A9000805204-K* € 301.70

>30x30 up to 50x50 cm
A9000805200-K* € 437.70

Seat cushion with radius contour Inlet: Anti-Decubitus

up to 30x30 cm
A9000800311-K* € 356.10

>30x30 up to 50x50 cm
A9000800310-K* € 482.90

Seat cushion height:

4 cm 6 cm

Cover:

- Black, breathable, open surface
 Black, breathable, closed surface

If not specified, the cover material „Black, breathable, open surface“ will be delivered.

Seat cushion is provided

Clothing guard

Aluminum clothing guard fixed to the side panel

Only in combination with fixed backrest



in black gloss

A9000403100-K € 243.20

in frame colour (base colour for two-colour design)

A9000403101-K € 311.10

in design colour (only possible with frame colour in two-colour design)

A9000403102-K € 345.10

Aluminum clothing guard removable via lock function

Adjustable following wheel arch

Only in combination with folding backrest



in anthracite matt

A9000403210-K € 243.20

in frame colour (base colour for two-colour design)

A9000403211-K € 311.10

in design colour (only possible with frame colour in two-colour design)

A9000403212-K € 345.10

Carbon clothing guard removable via lock function

A9000400012-K € 547.30

Adjustable following wheel arch

Only in combination with folding backrest



Armrests

Armrest, height adjustable, integrated in the clothing guard

Including aluminum clothing guard
with arm padding 25 cm long
Removable via lock function
Only in combination with folding backrest



in anthracite matt

A9000400008-K € 569.60

in frame colour (base colour for two-colour design)

A9000400009-K € 621.90

in design colour (only possible with frame colour in two-colour design)

A9000400010-K € 648.10

Armrest, height adjustable, integrated in the clothing guard

A9000400018-K € 862.00

Including carbon clothing guard
with arm padding 25 cm long
Removable via lock function
Only in combination with folding backrest

Arm padding 25 cm length mounted on clothing guard

A9000403151 € 165.80

Upper arm padding edge depending on the clothing guard setting mounted approx. 20 cm to 24 cm above the seat tube.
Only in combination with clothing guard



Armrests folding up, removable

A9000403141-K € 358.50

with arm padding 25 cm long
Positioning: upper arm padding edge approx. 25 cm to 31 cm above the seat tube, adjustable in 1.5 cm steps
Not possible in combination with back-positioned safety push handles




Axle plate

Standard equipment with 12 mm sockets / axles

Standard axle plate


A9016201000
Rear seat height:
Size L (22") approx. 33,5 cm - 37,5 cm
Size L/XL (24") approx. 36 cm - 40 cm
Size XXL (24") approx. 39 cm - 43 cm



Long axle plate

A9016201020 € 234.70

Rear seat height:
Size L (22") approx. 35,5 cm - 41,5 cm
Size L/XL (24") approx. 38 cm - 44 cm
Size XXL (24") approx. 41 cm - 47 cm



Wheel camber

*Total width with tight-mounted handrim.

1° drive wheel camber

A9016202010
Total width = approx. SB + 20 cm*

2,5° drive wheel camber

A9016202051
Total width = approx. SB + 23 cm*

4° drive wheel camber

A9016202030
Total width = approx. SB + 26 cm*

6° drive wheel camber

A9016202040
Total width = approx. SB + 30 cm*

8° drive wheel camber

A9016202050 € 77.00
Total width = approx. SB + 33 cm*
Only possible up to a body weight of 50 kg

E-Adaptations

Preparation for the „e-motion“ / „e-motion DuoDrive“ mounting

A9016201200-K € 502.90

Equipment: e-motion torque arm, 1/2" axle and socket length 50.8 mm (2")
PRO ACTIV and Alber recommend mounting two anti-tipping supports. Please select under „Adaptations“.

Preparation for the „e-fix“ mounting

A9016201210-K € 502.90

Equipment: 1/2" axle, socket length 50.8 mm (2") and e-fix torque support.
PRO ACTIV and Alber recommend mounting two anti-tipping supports. Please select under „Adaptations“.
Please order Alber anti-tipping supports and corresponding components direct from Alber.

Brake system

If not specified, the „Brake system push to lock with standard brake handle" will be delivered.

Recommendation: For KP 7 cm or lower, we recommend the „Brake system push to lock with brake handle 20° outwards“.

Knee lever brake with standard brake lever

A9000701010-001



Knee lever brake with long brake lever

A9000701010-002



Knee lever brake with low brake lever

A9000701010-003



Knee lever brake system Brake lever angled 20°

standard brake handle A9000700002

long brake handle A9000701010-004



Brake lever with M8 adapter and plastic ball

A9000706001 € 170.50

only in combination with knee lever brake system, standard or long brake handle



Brake system push to lock Brake handle to flip-down

A9000706008 € 158.60

Length: 9 cm



Knee lever brake with one-hand operation

right side A9000707020-001 € 392.30

left side A9000707020-002 € 392.30

via Bowden cable



Parking knee lever brake for accompanying person

A9000701020 € 464.10

Push handle operation via Bowden cable.



Knee lever brake with lock

standard brake handle A9000700214 € 158.60

long brake handle A9000700215 € 158.60

Locking block for brake lever
Prevents braking during pushing by accompanying persons
Not in combination with knee lever brake with one-hand operation, brake lever angled 20° as well as parking knee lever brake for accompanying person



Drum brake

A9000702010-K* € 474.00

Operating from the push handle

Total wheelchair width (see selection „Wheel camber“) + 3 cm

Depending on the seat height, the long wheel plate is mounted in order to protect the drum brake anchors.

Please note:

- Only in combination with drive wheel aluminum 36-spoke rim wheel
- Not in combination with one-hand drive, e-fix or e-motion



Footrest

If not specified, the „continuous footrest“ will be delivered.

Continuous footrest

A9016601010-K*
with tubular resting bracket
height adjustable +/- 3 cm
Size L / XL from USL 26 cm
Size XXL from USL 29 cm
If you want to equip the footrest with a „rubber coated footplate“ please select this under „Footrest accessories“.



Footrest folding up to one side

A9016601300-K* € 173.90
with tubular resting bracket
height adjustable +/- 3 cm
cannot be removed
for Size L / XL from USL 25 cm
for Size XXL from USL 28 cm
can also be swivelled outwards:
for Size L / XL from USL 29 cm
for Size XXL from USL 32 cm
If you want to equip the footrest with a „rubber coated footplate“ please select this under „Footrest accessories“.



Footrest folds up to the rear with spring locking mechanism

A9016602005-K* € 287.40
with tubular resting bracket
for Size L / XL from USL 31 cm
for Size XXL from USL 35 cm
If you want to equip the footrest with a „rubber coated footplate“ please select this under „Footrest accessories“.



Footrest folds up to the rear

A9016602010-K* € 307.90
with tubular resting bracket and 20 mm operating lever
Please note:
Size L / XL from USL 39 cm
Size XXL from USL 42 cm
straight front end from SB 33 cm
If you want to equip the footrest with a „rubber coated footplate“ please select this under „Footrest accessories“.

Footrest folds up to the rear for USL < 39 cm and from SB 33 cm

A9000600007-K € 352.50
Footrest joint clamped to the frame,
with tubular resting bracket and 20 mm operating lever
If you want to equip the footrest with a „rubber coated footplate“ please select this under „Footrest accessories“.
Not in combination with therapy table



Divided footrest

A9016603200-K € 336.80
with rubber coated footplate,
height adjustable +/- 3 cm, angle adjustment,
can be pivoted outwards and folded up to the side
Including frame reinforcement
Please note:
Size L / XL from USL 31 cm
Size XXL from USL 34 cm
cannot be removed



„Swing Away“ footrest

A9016604100-K € 404.70
with rubber coated footplate,
separated in the middle,
can be folded to the side, removable,
can be folded up, height-adjustable,
angle-adjustable
Including frame reinforcement
Please note:
Possible with straight shape Size L / XL / XXL
and abducted shape Size XXL



High-distanced footrest

A9041606010-K
height adjustable
Possible from USL 20 cm to 25 cm.
Rubber coated footplate necessary, please select this under „Footrest accessories“.



Footrest accessories

Rubber coated footplate with 14 cm depth

- A9000610000-K € 125.90
bolted, angle-adjustable



Rubber coated footplate with 18 cm depth

- A9000610005-K € 140.10
bolted, angle-adjustable



Footplate side stop

- A9000610100 € 59.70
adjustable
only in combination with rubber coated footplate



Deflector rollers mounted on the footplate

- A9000600006 € 117.40
only in combination with rubber coated footplate



Adaptations

Anti-tipping support

- right side A9000300100-K001 € 169.10
 left side A9000300100-K002 € 169.10
with tension spring and universal adapter pivotable
For optimum pivoting behaviour the seat slope should be > 3 cm



Preparation for mounting the anti-tipping support with lifting function from Alber (e-fix, e-motion)

- A9016300105 € 296.70
Universal adapter is prepared for mounting the anti-tipping support with lifting function from Alber.
Supplied without anti-tipping support and anti-tipping support mount.
On the customer side, please order this from Alber



Transit wheels

- right side A9000301100-003 € 164.80
 left side A9000301100-004 € 164.80
with universal adapter, removable
not in combination with anti-tipping support, right side or left side



Tilting clamp

- right side A9000302100-K001 € 140.10
 left side A9000302100-K002 € 140.10
with universal adapter and Quick pin
not in combination with anti-tipping support, right side or left side



Removable wheelbase extension SPEEDY 4all / 4you

- A9016102001 € 581.30
wheel position displacement approx. 14 cm adjustable track
Please select the wheel camber:
 1° drive wheel camber
 4° drive wheel camber
 6° drive wheel camber

Note: The camber of the wheelbase extension cannot be less than in the everyday position. A max. camber of 6° is possible.

If not specified, the next higher wheel camber will be mounted than in the everyday position.



Outdoor wheel attachment with steering fork

- A9000900100 € 845.70
Equipped with straight-running lock can be removed from the wheelchair via quick release lever with swivelling pneumatic tyre 12.1/2 x 2.1/4" caster wheel can be removed from the support tube via quick-release axle for maximum ground pressure at the steering attachment of 30 kg



3" caster wheel

3" skate wheel (80 mm)

A9000501110 € 70.80
puncture-proof



4" caster wheel

4" caster wheel „Standard“ solid rubber black

A9000501011-004 (82A), puncture-proof



4" skate wheel (100 mm)

A9000501120 € 70.80
puncture-proof



4" safety light-up wheel (100 mm)

A9000501220 € 94.80
puncture-proof



4" CorMORun caster wheel (100 mm)

A9000501042 € 173.00
with replaceable rubber 6 mm,
puncture-proof
Rim colour: Black
Incl. hair guard spacer (5.3 mm)



4" Frog Legs caster wheel (100 mm)

A9000501230 € 206.30
with replaceable rubber 6 mm,
puncture-proof
Colour: black, silver rim
Incl. hair guard spacer (5.3 mm)



5" caster wheel

5" caster wheel „Standard“ PU black

A9000501036-005 (82A), puncture-proof



5" caster wheel „Standard“ PU grey

A9000501031-005 (82A), puncture-proof



5" safety light-up wheel (125 mm)

A9000501221 € 94.80
puncture-proof



5" aluminum rim 3-spoke design „good shock absorption“

silver anodised A9000501210 € 179.80

black A9000501211 € 179.80

with replaceable rubber 6 mm
medium hard, puncture-proof
Incl. hair guard spacer (3 mm)



5" CorMORun caster wheel (125 mm)

A9000501041 € 181.30
with replaceable rubber 6 mm,
puncture-proof
Rim colour: Black
Incl. hair guard spacer (5.3 mm)



5" Frog Legs caster wheel (125 mm)

A9000501232 € 206.30
with replaceable rubber 6 mm,
puncture-proof
Colour: black, silver rim
Incl. hair guard spacer (5.3 mm)



6" caster wheel

6" caster wheel „Standard“ PU grey

A9000501010-006
(82A), puncture-proof



6" caster wheel „pneumatic tyre“

A9000501050 € 57.20
Front seat height min. 42 cm



6" caster wheel „Flexel-System“

A9000501070 € 57.20
puncture-proof



Caster fork

For caster fork colours, see „Colour of add-on parts“

If not specified, the „Standard caster fork with screw axle“ will be delivered.

Caster fork Standard

A9000500100-K001
with screwed axle

A9000500100-K002 € 167.40
with quick-release axle



Caster fork Heavy Duty

for high loads or
with the use of additional electric drives

A9000510111-K001 € 54.00
with screwed axle

A9000510111-K002 € 218.90
with quick-release axle



Suspension fork

inkl. 5" aluminum rim „good shock absorption“

A9000500160-K002 € 616.90
with quick-release axle



Caster fork spare wheel „Standard“

A9000500300-K € 296.70
Caster fork Standard with quick-release axles



Caster fork spare wheel „Heavy Duty“

A9000500301-K € 337.90
Caster fork Heavy Duty with quick-release axles



The spare caster fork does not have a caster wheel.
In case a caster wheel is requested, please put a cross on a second selection in the category „Caster wheels“.

Track distance

Caster wheel track extension



A9000500400* € 59.10
10 mm per side

A9000500401* € 59.10
15 mm per side

A9000500410* € 59.10
20 mm per side

Drive wheels

If not specified, the „drive wheel 36 spoke aluminum rim“ will be delivered.

Drive wheel 18 spoke aluminum rim

22" (ETRTO 489-19)
A9000203062

24" (ETRTO 540)
A9000203122

with Liteline hub, radially spoked

Not in combination with:

- Puncture-proof tires
- One-hand drive
- Drum brake
- Transport restraint connection
- Handbike

up to 70 kg body weight:

Sturz	1°	4°	6°	8°	10°
Gewicht					
bis 50 Kg	✓	✓	✓	✓	✓
51 - 60 Kg	✓	✓	✓	✓	
61 - 70 Kg	✓	✓			



Drive wheel 36 spoke aluminum rim

22" (ETRTO 489-19)
A9000203060-010

24" (ETRTO 540)
A9000203000-010

with Liteline hub (left image),
radially spoked

22" (ETRTO 451-19)
A9000203060-020

24" (ETRTO 540)
A9000203000-020

with aluminum front hub (right image),
radially spoked
in combination with one-hand drive



Drive wheel double rim 36-spokes

24" (ETRTO 540)
A9000203100

€ 157.80

25" (ETRTO 559)
A9000203200

€ 157.80

26" (ETRTO 590)
A9000203300

€ 157.80

with aluminum front hub,
radially spoked
for use with high loads



Spinergy wheel „Spox Sport“ with 24 PBO spokes; rim, hub and spokes in black

24" (ETRTO 540)
A9000203550

€ 784.10

25" (ETRTO 559)
A9000201559

€ 784.10

26" (ETRTO 590)
A9000203555

€ 784.10



Spinergy „Everyday Light Extreme“ with 12 PBO spokes; rim and hub in black, spokes in white

24" (ETRTO 540)
A9000203554

€ 957.20

Recommendation: For heavy loads
(e.g. higher body weight), we recommend the
Spinergy „Spox Sport“ wheel with 24 PBO
spokes.



Spinergy „Everyday Light Extreme“ with 12 PBO spokes; rim, hub and spokes in black

24" (ETRTO 540)
A9000203552


€ 957.20

Recommendation: For heavy loads
(e.g. higher body weight), we recommend the
Spinergy „Spox Sport“ wheel with 24 PBO
spokes.



Handrims

x = dimension handrim attachment



Handrim attachment
general:

- x = 11 mm
- x = 19 mm

Handrim attachment
Curve + Gekko handrim:

- x = 10 mm
- x = 18 mm

If not specified, the handrim is mounted wide (19 mm / 18 mm).

Stainless steel Ø 20 mm

- 22" handrim for 22" wheel
A9000205010-002
- 24" handrim for 24" wheel and 25" wheel
A9000205010-004
- 26" handrim for 26" wheel
A9000205010-006



Aluminum Ø 19 mm silver anodised

- 22" handrim for 22" wheel
A9000205019-002
- 24" handrim for 24" wheel and 25" wheel
A9000205019-004



Aluminum Ø 19 mm powder coated in frame colour

(basic colour with two-colour design)

- = Gloss surface
- = Matt surface

Note: Matt coating is not recommended due to the low grip.

- 22" handrim for 22" wheel **€ 68.30**
A9000205020-002
- 24" handrim for 24" wheel and 25" wheel **€ 68.30**
A9000205020-004
- 26" handrim for 26" wheel **€ 68.30**
A9000205020-006



Aluminum Ø 19 mm polished and hard-anodised

- 22" handrim for 22" wheel **€ 164.90**
A9000205030-002
- 24" handrim for 24" wheel and 25" wheel **€ 164.90**
A9000205030-004
- 26" handrim for 26" wheel **€ 164.90**
A9000205030-006



Titanium Ø 19 mm

- 24" handrim for 24" wheel and 25" wheel **€ 335.60**
A9000205040-004
- 26" handrim for 26" wheel **€ 335.60**
A9000205040-006



Silicone handrim cover

- for handrim 22" **€ 232.80**
A9000205100-002
- for handrim 24" **€ 232.80**
A9000205100-004
- for handrim 26" **€ 232.80**
A9000205100-006

Surface without structure, compatible with handrim Ø 19 mm or Ø 20 mm



<p>Ultra-Grip „very grippy“</p> <p><input type="checkbox"/> 22" handrim for 22" wheel € 402.40 A9000205070-002</p> <p><input type="checkbox"/> 24" handrim for 24" wheel and 25" wheel € 402.40 A9000205070-004</p> <p><input type="checkbox"/> 26" handrim for 26" wheel € 402.40 A9000205070-006</p> 	<p>Ergo-Grip</p> <p><input type="checkbox"/> 22" handrim for 22" wheel € 402.40 A9000205080-002</p> <p><input type="checkbox"/> 24" handrim for 24" wheel and 25" wheel € 402.40 A9000205080-004</p> <p><input type="checkbox"/> 26" handrim for 26" wheel € 402.40 A9000205080-006</p> <p>Attachment fitting to the rim</p> <p>Note: Due to the complex geometry, concentricity deviations of +/- 2-3 mm are possible.</p> 	<p>Curve black anodised, ergonomic handrim profile</p> <p><input type="checkbox"/> 24" handrim for 24" wheel and 25" wheel Profile height 24 mm € 402.40 A9000205130-004</p> <p><input type="checkbox"/> 24" handrim for 24" wheel and 25" wheel Profile height 31 mm € 402.40 Size L A9000205140-004</p> <p><input type="checkbox"/> 26" handrim for 26" wheel Profile height 31 mm € 402.40 Size L A9000205140-006</p> 
<p>Curve silver anodised ergonomic handrim profile</p> <p><input type="checkbox"/> 24" handrim for 24" wheel and 25" wheel Profile height 31 mm € 402.40 Size L A9000205160-004</p> 	<p>Curve powder-coated, ergonomic handrim profile</p> <p><input type="checkbox"/> = frame colour (gloss surface) <input type="checkbox"/> = black (gloss surface)</p> <p><input type="checkbox"/> 24" handrim for 24" wheel and 25" wheel Profile height 24 mm € 467.60 A9000205135-004</p> <p><input type="checkbox"/> 24" handrim for 24" wheel and 25" wheel Profile height 31 mm € 467.60 Size L A9000205145-004</p> <p><input type="checkbox"/> 26" handrim for 26" wheel Profile height 31 mm € 467.60 Size L A9000205145-006</p> 	<p>Curve Grip, with Tetra-Grip coating, ergonomic handrim profile</p> <p><input type="checkbox"/> 24" handrim for 24" wheel and 25" wheel Profile height 31 mm € 520.70 Size L A9000205150-004</p> 
<p>Quadro with Tetra-Grip coating</p> <p><input type="checkbox"/> 24" handrim for 24" wheel and 25" wheel € 520.70 A9000205170-004</p> <p><input type="checkbox"/> 25" handrim for 25" wheel € 520.70 A9000205170-005</p> <p><input type="checkbox"/> 26" handrim for 26" wheel € 520.70 A9000205170-006</p> 	<p>Gekko black anodised</p> <p><input type="checkbox"/> 24" handrim for 24" wheel and 25" wheel € 520.70 A9000205175-004</p> <p><input type="checkbox"/> 25" handrim for 25" wheel € 520.70 A9000205175-005</p> <p><input type="checkbox"/> 26" handrim for 26" wheel € 520.70 A9000205175-006</p> <p>with a robust and replaceable silicone lip, fully integrated into the profile Profile height 31 mm</p> 	<p>Gekko silver anodised</p> <p><input type="checkbox"/> 24" handrim for 24" wheel and 25" wheel € 520.70 A9000205185-004</p> <p>with a robust and replaceable silicone lip, fully integrated into the profile Profile height 31 mm</p> 

<p>Gekko S black anodised</p> <p><input type="checkbox"/> 22" handrim for 22" wheel € 520.70 A9000205195-002</p> <p><input type="checkbox"/> 24" handrim for 24" wheel and 25" wheel € 520.70 A9000205195-004</p> <p>with a robust and replaceable silicone lip, fully integrated into the profile Profile height 24 mm</p> 	<p>Equipment without handrims</p> <p><input checked="" type="checkbox"/> A9000205990</p>	
---	---	--

Spoke guard

<p>Finger- and spoke guard „transparent“ Not in combination with Spinergy drive wheels</p>			
<p><input type="checkbox"/> for 22" wheel € 207.10 A9000206300-002</p>	<p><input type="checkbox"/> for 24" wheel € 207.10 A9000206300-004</p>	<p><input type="checkbox"/> for 26" wheel € 207.10 A9000206300-006</p>	

Tyres

If not specified, the „Schwalbe RightRun in black“ tyres will be delivered.

<p>Schwalbe „RightRun“ High-pressure tyres</p> <p><input checked="" type="checkbox"/> 22x1" (489-25) A9000204001-002</p> <p><input checked="" type="checkbox"/> 24x1" (540-25) A9000204001-004</p> <p>Colour: Black 7.5 bar</p> 	<p>Schwalbe „RightRun“ High-pressure tyres</p> <p><input type="checkbox"/> 22x1" (489-25) A9000204000-002</p> <p><input type="checkbox"/> 24x1" (540-25) A9000204000-004</p> <p><input type="checkbox"/> 25x1" (559-23) A9000204000-005</p> <p><input type="checkbox"/> 26x1" (590-25) A9000204000-006</p> <p>Colour: Grey, black wall 7.5 bar</p> 	<p>Schwalbe „DownTown“ Profiled tyres</p> <p><input type="checkbox"/> 24x1" (540-25) A9000204010-004</p> <p>7.5 bar</p> 
--	--	---

<p>Schwalbe „Marathon Plus Evolution“ with „smart guard“ puncture protection</p> <p><input type="checkbox"/> 24x1" (540-25) € 49.20 A9000204030-014</p> <p><input type="checkbox"/> 25x1" (559-25) € 49.20 A9000204030-015</p> <p><input type="checkbox"/> 26x1" (590-25) € 49.20 A9000204030-016</p> <p>Colour: Black 10 bar</p> 	<p>Schwalbe „Marathon Plus“ with „smart guard“ puncture protection</p> <p><input type="checkbox"/> 24x1" (540-25) € 49.20 A9000204030-004</p> <p>Colour: Grey, black wall 10 bar</p> 	<p>Schwalbe „One“</p> <p><input type="checkbox"/> 24x1" (540-25) € 93.00 A9000204123-004</p> <p>Colour: Black 9 bar</p> 
---	---	--

Puncture-proof „solid rubber tyres” <input type="checkbox"/> 22x1" (489-25) € 122.70 A9000204200-002 Colour: Grey <input type="checkbox"/> 24x1" (540-25) € 122.70 A9000204201-004 Colour: Black Not in combination with 18 spoke or Spinergy drive wheels.	
	

Accessories

24" MTB spare wheel set

with double rim (ETRTO 507), aluminum front hub, quick-release axle and stainless steel handrim
Mountain bike tyres, selectively:

A9000203410 € 703.40
24x1.75" Schwalbe "Marathon Plus"
Brake setting compatible with 24x1" wheel



A9000203413 € 703.40
24x2.1" Schwalbe "SMART SAM Sport"
Brake setting not compatible with 24x1" wheel




NOTE: We charge a surcharge for the use of drum brakes

A9000200004 € 232.20

Through axle with assistance, Large push button

A9000200096 € 105.80
Interlock dimension: 122 / 110 mm
Axle diameter: 12 mm
Not in combination with one-hand drive



Hand air compressor with pressure gauge

A9000901100 € 106.30
up to 14 bar



Crutch holder with flexible nylon cup

A9000902310* € 137.40



Length-adjustable safety belt with buckle

A9000902100 € 80.60



Adjustable safety belt with hook and fastener

A9000902110 € 64.50



SIGMA lighting set

A9000903270 € 107.90

Consisting of front light and rear lights, chargeable via Micro USB socket, incl. Micro USB charging cable
Lighting set compliant with road traffic regulations



Transport bag for drive wheels

for 24" wheel
A9000903530-004 € 77.80
Ø 76 cm x 2 cm
 for 26" wheel
A9000903530-006 € 77.80
Ø 83 cm x 2 cm



Separate order sheets

Category „Rehabilitation technology accessories”

e.g. Seat & backrest systems, one-hand drive or transport restraint systems can be found at:
www.proactiv-gmbh.com/downloads/

Additional products

such as hand bikes or traction & attachment devices can be found at:
www.proactiv-gmbh.com/downloads/

Additional sheets

such as individual knee angle etc. can be found at:
www.proactiv-gmbh.com/downloads/

Services

Consultation, clarification of the installation situation and measurement through PRO ACTIV

A9990000000 € 354.50

Handover and instruction on the wheelchair by PRO ACTIV

A9990000006 € 354.50

Notes and remarks:

Location, Date

The signature field is intended for the consultant and insured person. This confirms that the intended equipment is intended for the insured person.
- 20 - 2026-03-03 / The prices are plus statutory value-added tax and shipping costs = included in the basic price = Position selected

Order sheet SPEEDY 4teen Rev. 06 UK

Signature End customer/dealer