

☐ **Quotation date:** _____

Consultation was carried out by
Staff at the medical supplies shop: _____

Enquiry no.: _____

☐ **Order date:** _____

according to GTC PRO ACTIV Reha-Technik GmbH

Order no.: _____

Medical suppliers / company address:

Cust. no.: _____

Name: _____

Street: _____

Postcode / city: _____

Phone: _____

Fax: _____

E-Mail: _____

Contact: _____

☐ **Different invoicing address (see below)**

Cust. no.: _____

Name: _____

Street: _____

Postcode / city: _____

☐ **Different shipping address (see below)**

Cust. no.: _____

Name: _____

Street: _____

Postcode / city: _____

Customer name / commission:

ID/Com. No.: _____

Surname, forename: _____

Street: _____

Postcode / city: _____

Phone: _____

E-Mail: _____

D.o.B.: _____

Paying organisation: _____

Indication: _____



HUSK-E adaptive bike mechanical

manual without power-assisted drive

up to 140 kg load capacity*

*load capacity = user weight + wheelchair + items transported on or in the wheelchair

basic equipment

- Docking plate for quick connection to the „Quick-Snap adapter“
- Adjustable and foldable bottom bracket support to minimise pack size
- Disc brake Ø 180 mm with brake lever or coaster disc brake possible
- Mechanical parking disc brake
- Parking stand with extended operating lever
- Lighting set
- Mudguard
- 20" drive wheel

A9117100050-020

€ 4,230.00

Note: The HUSK-E mechanical can be upgraded to an electrically assisted hand bike with a retrofit kit if required. The retrofit kits can be found at shop.proactiv-gmbh.com.

Body size

Body height: _____ cm

Upper body height: _____ cm

User weight: _____ kg

This information is necessary for checking!

Heavy Duty version

Extended trailer load capacity*
from 140 kg to 170 kg

☐ A9117100100 € 507.00

Only possible in combination with an HD adapter.

Transmission

Chainring on bottom bracket with 28 teeth

☐ A9117200050
for a higher top speed

Chainring on bottom bracket with 22 teeth

☐ A9117200051
for a very light response

Derailleur gears

Shimano Deore 1x11 Rapid Fire 11-speed derailleur gears

☐ A9117200029
mechanical derailleur gear system with control on the handle.
11-speed cassette 11-42 teeth
and a gear range of 381 %



Shimano Deore XT 1x11 Rapid Fire 11-speed derailleur gears

☐ A9117200019 € 223.30
mechanical derailleur gear system with control on the handle.
11-speed cassette 11-42 teeth
and a gear range of 418 %



Shimano Deore XT 1x11 Di2 with push-button operation 11-speed derailleur gears

☐ A9117200028 € 1,486.10
electronic derailleur gear with operation on the handle via button.
11-speed cassette 11-42 teeth
and a gear range of 418 %

Operation: right side ☐ left side ☐
If not specified, the operation will be mounted on the „right side“.



Shimano Deore XT 1x11 Di2 with chin control 11-speed derailleur gears

☐ A9117200027 € 1,631.10
electronic derailleur gear with chin control (tetra-equipment).
11-speed cassette 11-42 teeth
and a gear range of 418 %



Brake system

Double hydraulic disc brake system „Magura“

☐ A9117300100 € 168.10
Two Magura 2-piston callipers
Parking brake via locking bar on the handle
Omitting the mechanical parking disk brake.

Mechanical back-pedal disc brake with crank release function

☐ A9117300110 € 738.40
compatible with derailleur gears and freewheel hub.
Incl. second mechanical disc brake designed as a parking brake.
The hydraulic disc brake is omitted.
In combination with:
- grip width 45 cm, 50 cm, 53 cm or 58 cm
- crank length 155 mm or 175 mm
up to 120 kg body weight



Driving brake - supplement for the mechanical parking brake

☐ A9117300120 € 116.10
The parking disc brake is supplemented with a brake handle on the left-hand handle.
This means that the mechanical disc brake can be used as a driving brake and parking brake.

Colour mudguard

If not specified, the mudguard and chain guard will be delivered in anthracite RAL 7016 matt.

Mudguard, chain guard and motor cover anthracite RAL 7016 matt

☒ A9117500004

Mudguard, chain guard and motor cover pure orange RAL 2004 gloss

☒ A9117500006



Mudguard, chain guard and motor cover lemon green gloss

☒ A9117500007



Mudguard, chain guard and motor cover customized colour selection

☐ A9117500010 **€ 153.50**

from colour palette one-colour design
powder colours wheelchair

A9000100170-____

Colour: _____

Tyres

20x2.15" Schwalbe „Marathon 365“

☒ A9117105103
(ETRTO 406-55)
profile tyre
puncture protection level 5



20x2.15" Schwalbe „Big Ben Plus“

☐ A9117105102 **€ 28.70**
(ETRTO 406-55)
profile tyre
puncture protection level 6



Dimensions

Crank type and grip width

Standard crank

☐ grip width 40 cm
A9117500100-040

☐ grip width 45 cm
A9117500100-045

☐ grip width 50 cm
A9117500100-050

Tetra crank (flared)

With extended crank arms for maximum upper body space

☐ grip width 48 cm
A9117500100-048

☐ grip width 53 cm
A9117500100-053

☐ grip width 58 cm
A9117500100-058

Bottom bracket support

maximum length from the pivot point in the fork to the bottom bracket shaft

☐ length 520 mm (short)
using shortened holding tubes **€ 155.30**
A9117500001

Recommendation: up to approx. 165 cm body height

☒ length 570 mm (standard)
A9117500008

Recommendation: from approx. 165 cm to 185 cm body height

☐ length 620 mm (long)
using a longer bottom bracket distance **€ 204.50**
A9117500003

Recommendation: from approx. 185 cm body height

Standard pedal crank ☐ length 155 mm A9110600100-155 ☐ length 175 mm A9110600100-175 ☐ length 195 mm A9110600100-195 	Tetra pedal crank ☐ length 155 mm A9100600105-155 € 136.10 ☐ length 175 mm A9100600105-175 € 136.10 ☐ Individual length 140 mm to 195 mm A9100600100 € 258.70 	
--	---	--

Handles

Palm rest on the handle ☐ A9110601200 € 102.80 upholstered, bolted, angle-adjustable. 	Tetra hand support, pivotable, adjustable grip distance ☐ A9110601101 € 169.40 	Tetra backrest swivel and adjustable grip width ☐ A9110601102 € 191.60 very stiff external support as round tube Ø 20 mm with foam rubber coating 
Tetra ergonomic handle, freestanding ☐ A9100600001 € 164.60 crank handles arranged vertically extends the reach by 8.6 cm can be combined with all other handle options 	Tetra shell handles ☐ A9110601100 € 172.40 	

Accessories Bike

Maneuvering castors for parking stand ☐ A9110900017 € 122.10 Maneuvering castors mounted on the parking stand via inlet 	Panniers incl. holders ☐ A9117900220 € 315.40 Brackets mounted on the fork arms max. load weight 5 kg each 	SIGMA BC 14.0 WL CAD bike computer ☐ A9110900018 € 171.50 with pedal frequency set, wireless 
SIGMA BC 8.0 ATS bike computer ☐ A9110901000 € 70.10 wireless, with mount 		

Options for users with reduced functions

Lever extension excenter lock for Quick Snap connection

☐ A9110700006 € 73.90

simplifies the opening and closing of the Quick Snap adapter
Tetra version, curved



Adaptation help for centering pin

☐ A9110700002 € 55.00
to simplify the insertion of the central adapter.



Chest strap for upper body fixation

☐ size S (105 cm) A9000800016 € 241.80

☐ size M (125 cm) A9000800013 € 241.80

☐ size L (145 cm) A9000800014 € 241.80

☐ size XL (165 cm) A9000800017 € 241.80

consisting of elastic webbing material in combination with neoprene with Velcro and operating tab.
Colour: black



Frame adapter

SPEEDY rigid frame adapter right-hand operation with „Quick-Snap-Adaption“

☐ A9110700142 € 537.20

with quick-release clamp and centering pin holder
Colour: black

Heavy Duty version

☐ A9110700144 € 816.70

for a load capacity* of 140 kg – 220 kg

Please select accordingly:

Centring pin for Ø 25 mm
for A1 and A2

☐ A9110700022

Centring pin for Ø 40 mm
for 4all, 4all Ergo, 4you Ergo, F2, F3 and 4teen

☐ A9110700023



TRAVELER folding frame adapter right-hand operation with „Quick-Snap-Adaption“

☐ A9110701302 € 537.20

with quick-release clamp and centering pin holder on the folding mechanism
Colour: black

Heavy Duty version

☐ A9110700146 € 816.70

for a load capacity* of 140 kg – 200 kg



BUDDY rigid frame adapter right-hand operation with „Quick-Snap-Adaption“

☐ A9110701401 € 537.20

with quick-release clamp and centering pin mount on the frame
Colour: black



Wheelchair front adapter right-hand operation with universal clamp and „Quick-Snap-Adaption“ ☐ A9110701502 € 572.30 Adapter attachment to the wheelchair stem. Operation on the right-hand side.	Wheelchair front adapter HD right-hand operation with universal clamp and „Quick-Snap-Adaption“ ☐ A9199701511 € 871,60 Adapter with reinforced double-tube steering head connection for wheelchair adaptations with a load capacity* of 140 kg – 220 kg Adapter mounting on the wheelchair frame extension 	
Mounting situation for PA products: Straight front end from SB 31 cm to 54 cm V-shaped front end from SB 30 cm to 50 cm		
Please select the frame clamp below:		
☐ clamping range 95 mm A9110701525 Tube diameter 23-27 mm	☐ clamping range 95 mm A9110701530 Tube diameter 28-32 mm	☐ clamping range 95 mm A9110701535 Tube diameter 33-37 mm
☐ clamping range 70 mm € 76.50 A9110900072 Tube diameter 23-27 mm	☐ clamping range 70 mm € 76.50 A9110900078 Tube diameter 28-32 mm	☐ clamping range 70 mm € 76.50 A9110900079 Tube diameter 33-37 mm
☐ clamping range 75 mm € 268.10 A9199000080 Tube diameter Küschall 2.0		
PLEASE NOTE: If the exact tube diameter is not specified, an adapter for a tube diameter of 23 - 27 mm and a clamping range of 95 mm will be delivered as standard.		
Please select the clamping sleeve below:		
☐ clamping sleeve [80 mm] included as standard for close clearance between footrest and drive wheel.	☐ clamping sleeve [110 mm] € 46.70 A9199000117 for medium distance between footrest and drive wheel.	☐ clamping sleeve [170 mm] € 56.60 A9199000154 for a wide distance between footrest and drive wheel.
PLEASE NOTE: If not specified, the clamping sleeve with 80 mm length will be delivered.		
Clamp lever with 3-way star handle for actuating the frame clamping ☐ A9199701507 € 106.10 mounted on the right and left for easy terminal operation with limited manual function Only in conjunction with front adapter 		

Wheelchair-dependent accessories

Extended bike plate

- ☐ for SPEEDY A1 / A2
A9110100100-002 € 370.00
- ☐ for BUDDY
A9110100100-003 € 370.00
- ☐ for TRAVELER
A9110100100-001 € 370.00
- 4° camber integrated
Track independently adjustable



Wheelbase extension via „grid mechanism with pulling operation“

- ☐ for SPEEDY 4all /
4teen / 4you /
TRAVELER 4you
A9110102001 € 564.70

removable, with continuous axle tube, track
width independently adjustable
Up to 150 kg user weight.

Please select the wheel camber:

- ☐ 1° drive wheel camber
- ☐ 4° drive wheel camber
- ☐ 6° drive wheel camber

Note: The camber of the wheelbase extension
must not be less than that in the normal driving
position. A maximum camber of 6° is permitted.



Wheelbase extension angle-adjustable and removable

- ☐ for SPEEDY F2 / F3
A9015100005 € 746.80

with continuous axle tube,
height adjustment via gear joint
adjustable in 6° angular increments
wheel position displacement approx. 14 cm
adjustable track

Attachment possibility must be checked by
PRO ACTIV.

Only in combination with base plate.
Up to 120 kg user weight.

Please select the wheel camber below:

- ☐ 1° drive wheel camber
- ☐ 4° drive wheel camber
- ☐ 6° drive wheel camber

Note: The camber of the wheelbase extension
must not be less than that in the normal driving
position. A maximum camber of 6° is permitted.



Services

Consultation, clarification of the installation situation and measurement through PRO ACTIV

- ☐ A9990000000 € 354.50

Handover and instruction of the handbike / traction device by PRO ACTIV

- ☐ A9990000011 € 354.50

Attachment of the handbike / traction device to the wheelchair by PRO ACTIV at specialist dealer

on a PRO ACTIV wheelchair

- ☐ A9990000012 € 230.00

product no.: _____

on a third-party product

- ☐ A9990000013 € 366.90

manufacturer: _____

model: _____

frame tube Ø ____ mm

Attachment of the handbike / traction device to the wheelchair through specialist dealer

including mounting of the frame adapter

on a PRO ACTIV wheelchair

- ☒ A9990000008

product no.: _____

on a third-party product

- ☒ A9990000005

manufacturer: _____

model: _____

frame tube Ø ____ mm

Notes and remarks: