

Im Hofstätt 11  
D-72359 Dotternhausen  
Tel +49 (0)7427 9480-0  
Fax +49 (0)7427 9480-7025  
e-mail info@proactiv-gmbh.de  
www.proactiv-gmbh.de



## **SPEEDY A2 for stair climber devices** **Active easy running rigid frame** **wheelchair**

Maximum user weight limit 120 kg

A9020100001 € 4,010.00

**Quotation date:** \_\_\_\_\_

**Consultation was carried out by**  
**Staff at the medical supplies shop:** \_\_\_\_\_

**Enquiry no.:** \_\_\_\_\_

**Order date:** \_\_\_\_\_

according to GTC PRO ACTIV Reha-Technik GmbH

**Order no.:** \_\_\_\_\_

### **Medical suppliers / company address:**

**Cust. no.:** \_\_\_\_\_  
**Name:** \_\_\_\_\_  
**Street:** \_\_\_\_\_  
**Postcode / city:** \_\_\_\_\_  
**Phone:** \_\_\_\_\_  
**Fax:** \_\_\_\_\_  
**E-Mail:** \_\_\_\_\_  
**Contact:** \_\_\_\_\_

**Shipping address:**  **same as company address**  
 **Other (see below)**

**Name:** \_\_\_\_\_  
**Street:** \_\_\_\_\_  
**Postcode / city:** \_\_\_\_\_

### **Customer name / commission:**


**ID / Comm. No.:** \_\_\_\_\_  
**Surname, forename:** \_\_\_\_\_  
**Street:** \_\_\_\_\_  
**Postcode / city:** \_\_\_\_\_  
**Phone:** \_\_\_\_\_  
**Fax:** \_\_\_\_\_  
**E-Mail:** \_\_\_\_\_  
**D.o.B.:** \_\_\_\_\_  
**Paying organisation:** \_\_\_\_\_

**Invoicing address:**  **same as company address**  
 **Other (see below)**

**Cust. no.:** \_\_\_\_\_  
**Name:** \_\_\_\_\_  
**Street:** \_\_\_\_\_  
**Postcode / city:** \_\_\_\_\_  
**Phone:** \_\_\_\_\_  
**Fax:** \_\_\_\_\_  
**E-Mail:** \_\_\_\_\_  
**Contact:** \_\_\_\_\_

# Basic size

## Body size:



Body height: \_\_\_\_\_ cm


Weight: \_\_\_\_\_ kg

Lower leg length: \_\_\_\_\_ cm  
(anatomic dimension: hollow of the knee to the foot sole)

Indication: \_\_\_\_\_

This information is necessary for checking. If the information is missing, a check of the relevant items cannot be made!

## Wheelchair dimensions:



<b>Seat width SB</b>	<input type="checkbox"/> 38 <input type="checkbox"/> 39 <input type="checkbox"/> 40 <input type="checkbox"/> 41 <input type="checkbox"/> 42 <input type="checkbox"/> 43 <input type="checkbox"/> 44 <input type="checkbox"/> 45 <input type="checkbox"/> 46 <input type="checkbox"/> 47 <input type="checkbox"/> 48 <input type="checkbox"/> 49 <input type="checkbox"/> 50									
<b>Seat depth ST</b>	<b>Short frame (43 cm)</b> <input type="checkbox"/> 36 <input type="checkbox"/> 37 <input type="checkbox"/> 38 <input type="checkbox"/> 39 <input type="checkbox"/> 40 <input type="checkbox"/> 41 <input type="checkbox"/> 42 <input type="checkbox"/> 43				<b>Long frame (47 cm)</b> <input type="checkbox"/> 40 <input type="checkbox"/> 41 <input type="checkbox"/> 42 <input type="checkbox"/> 43 <input type="checkbox"/> 44 <input type="checkbox"/> 45 <input type="checkbox"/> 46 <input type="checkbox"/> 47 <input type="checkbox"/> 48					
<b>Back height RH</b>	<input type="checkbox"/> 20 <input type="checkbox"/> 21 <input type="checkbox"/> 22 <input type="checkbox"/> 23 <input type="checkbox"/> 24 <input type="checkbox"/> 25 <input type="checkbox"/> 26 <input type="checkbox"/> 27 <input type="checkbox"/> 28 <input type="checkbox"/> 29 <input type="checkbox"/> 30 <input type="checkbox"/> 31 <input type="checkbox"/> 32 <input type="checkbox"/> 33 <input type="checkbox"/> 34 <input type="checkbox"/> 35 <input type="checkbox"/> 36 <input type="checkbox"/> 37 <input type="checkbox"/> 38 <input type="checkbox"/> 39 <input type="checkbox"/> 40 <input type="checkbox"/> 41 <input type="checkbox"/> 42 <input type="checkbox"/> 43 <input type="checkbox"/> 44 <input type="checkbox"/> 45 <input type="checkbox"/> 46 <input type="checkbox"/> 47 <input type="checkbox"/> 48									
<b>Lower leg length USL</b>	<b>Short frame</b> <input type="checkbox"/> 33 <input type="checkbox"/> 34 <input type="checkbox"/> 35 <input type="checkbox"/> 36 <input type="checkbox"/> 37 <input type="checkbox"/> 38 <input type="checkbox"/> 39 <input type="checkbox"/> 40 <input type="checkbox"/> 41 <input type="checkbox"/> 42 <input type="checkbox"/> 43 <input type="checkbox"/> 44 <input type="checkbox"/> 45 <input type="checkbox"/> 46 <input type="checkbox"/> 47 <input type="checkbox"/> 48 <input type="checkbox"/> 49 <input type="checkbox"/> 50					<b>Long frame</b> <input type="checkbox"/> 35 <input type="checkbox"/> 36 <input type="checkbox"/> 37 <input type="checkbox"/> 38 <input type="checkbox"/> 39 <input type="checkbox"/> 40 <input type="checkbox"/> 41 <input type="checkbox"/> 42 <input type="checkbox"/> 43 <input type="checkbox"/> 44 <input type="checkbox"/> 45 <input type="checkbox"/> 46 <input type="checkbox"/> 47 <input type="checkbox"/> 48 <input type="checkbox"/> 49 <input type="checkbox"/> 50 <input type="checkbox"/> 51 <input type="checkbox"/> 52				
<b>Seat height back SHh</b> (with 24" drive wheels)	<b>Short frame</b> <input type="checkbox"/> 37 <input type="checkbox"/> 38 <input type="checkbox"/> 39 <input type="checkbox"/> 40 <input type="checkbox"/> 41 <input type="checkbox"/> 42 <input type="checkbox"/> 43				<b>Long frame</b> <input type="checkbox"/> 39 <input type="checkbox"/> 40 <input type="checkbox"/> 41 <input type="checkbox"/> 42 <input type="checkbox"/> 43 <input type="checkbox"/> 44 <input type="checkbox"/> 45					
<b>Seat height front SHv</b> (with 5" caster wheels)	<b>Short frame</b> <input type="checkbox"/> 44 <input type="checkbox"/> 45 <input type="checkbox"/> 46 <input type="checkbox"/> 47 <input type="checkbox"/> 48 <input type="checkbox"/> 49 <input type="checkbox"/> 50 <input type="checkbox"/> 51				<b>Long frame</b> <input type="checkbox"/> 46 <input type="checkbox"/> 47 <input type="checkbox"/> 48 <input type="checkbox"/> 49 <input type="checkbox"/> 50 <input type="checkbox"/> 51 <input type="checkbox"/> 52 <input type="checkbox"/> 53 <input type="checkbox"/> 54					
<b>Tipping point KP</b> (inner edge of rear tube to centre of drive wheel axle)	<input type="checkbox"/> 6 <input type="checkbox"/> 7 <input checked="" type="checkbox"/> 8 <input type="checkbox"/> 9 <input type="checkbox"/> 10 <input type="checkbox"/> 11 <input type="checkbox"/> 12 <input type="checkbox"/> 13									

Please tick the required dimensions (all dimensions are in cm)

The underlined values are the standard settings. If not specified, the standard value will be delivered.

The SHh and SHv values might be a bit different from the specified values, depending on the wheel size and camber. Individual dimensions available on request.

We recommend a seat height difference (difference between SHv and SHh) between 5 cm and 12 cm.

**scalamobil attachment is possible from SB 38 cm.**

## Frame colour - One-colour design

Information on current colour designs and example images can be found at: <https://www.proactiv-gmbh.com/colours/>

### Frame colour: Selection from colour palette

Frame colour will be specified in the order

A9000100170-000

Raw, blasted, not coated

A9000100170-020

Note: In combination with moisture e.g. hand perspiration, can affect dark aluminum discoloration.

### Powder colours

**white RAL 9016**

gloss A9000100170-023

matt A9000100170-123 € 43.40

**brilliant silver**

gloss A9000100170-011

**anthracite RAL 7016**

gloss A9000100170-001

matt A9000100170-101 € 43.40

**black RAL 9005**

gloss A9000100170-003

matt A9000100170-019 € 43.40

**pure orange RAL 2004**

gloss A9000100170-022

matt A9000100170-122 € 43.40

**pink magenta RAL 4010**

gloss A9000100170-037

matt A9000100170-137 € 43.40

**pastel blue RAL 5024**

gloss A9000100170-036

matt A9000100170-136 € 43.40

**lemon green**

gloss A9000100170-005

matt A9000100170-105 € 43.40

**yellow RAL 1023**

gloss A9000100170-046

matt A9000100170-146 € 43.40

**red RAL 3001**

gloss A9000100170-047

matt A9000100170-147 € 43.40

**lagoon blue**

gloss A9000100170-052

matt A9000100170-152 € 43.40

**athletic grey RAL 7046**

gloss A9000100170-050

matt A9000100170-150 € 43.40

**muddy grey RAL 7033**

gloss A9000100170-051

matt A9000100170-151 € 43.40

**peach RAL 3022**

gloss A9000100170-038

matt A9000100170-138 € 43.40

### Powder colours „dormant“

**deep black effect silver**

gloss A9000100170-007

**sky blue**

gloss A9000100170-017

**navy blue**

gloss A9000100170-045

matt A9000100170-145 € 43.40

**electric green**

gloss A9000100170-015

matt A9000100170-115 € 43.40

**violet**

gloss A9000100170-044

**fire red**

gloss A9000100170-018

matt A9000100170-118 € 43.40

### Frame colour - One-colour design special colours „liquid coating“

€ 349.60

**storm (ice blue)**

gloss A9000100170-500

matt A9000100170-501

**prairie (rose gold)**

gloss A9000100170-504

matt A9000100170-505

**berry**

gloss A9000100170-526

matt A9000100170-527

**blue ocean**

gloss A9000100170-522

matt A9000100170-523

**emerald**

gloss A9000100170-514

matt A9000100170-515

**chili**

gloss A9000100170-516

matt A9000100170-517

## Frame colour - Two-colour design

Note: In the two-colour design, both colours are the same gloss level as the design colour.

**Two-colour design with one base colour (on the rear frame) and one design colour (on the front frame)** € 268.50

Please select a colour below:

### Base colour (on the rear frame)

<b>white RAL 9016</b> <input type="checkbox"/> A9000100170-200	<b>brilliant silver</b> <input type="checkbox"/> A9000100170-201	<b>anthracite RAL 7016</b> <input type="checkbox"/> A9000100170-202
<b>black RAL 9005</b> <input type="checkbox"/> A9000100170-203	<b>pastel blue RAL 5024</b> <input type="checkbox"/> A9000100170-206	<b>lemon green</b> <input type="checkbox"/> A9000100170-207
<b>yellow RAL 1023</b> <input type="checkbox"/> A9000100170-208	<b>lagoon blue</b> <input type="checkbox"/> A9000100170-212	<b>athletic grey RAL 7046</b> <input type="checkbox"/> A9000100170-210
<b>muddy grey RAL 7033</b> <input type="checkbox"/> A9000100170-211		

### Design colour (on the front frame)

<b>black RAL 9005</b> <input type="checkbox"/> gloss A9000100170-259 <input type="checkbox"/> matt A9000100170-260	<b>lemon green</b> <input type="checkbox"/> gloss A9000100170-250 <input type="checkbox"/> matt A9000100170-251	<b>chili</b> <input type="checkbox"/> gloss A9000100170-253 <input type="checkbox"/> matt A9000100170-254
<b>pure orange RAL 2004</b> <input type="checkbox"/> gloss A9000100170-255 <input type="checkbox"/> matt A9000100170-256	<b>emerald</b> <input type="checkbox"/> gloss A9000100170-274 <input type="checkbox"/> matt A9000100170-275	<b>storm (ice blue)</b> <input type="checkbox"/> gloss A9000100170-264 <input type="checkbox"/> matt A9000100170-265
<b>prairie (rose gold)</b> <input type="checkbox"/> gloss A9000100170-266 <input type="checkbox"/> matt A9000100170-267	<b>berry</b> <input type="checkbox"/> gloss A9000100170-284 <input type="checkbox"/> matt A9000100170-285	<b>blue ocean</b> <input type="checkbox"/> gloss A9000100170-280 <input type="checkbox"/> matt A9000100170-281

## Colour - Attachment parts

Note: If not specified, the caster fork and bearing block will be mounted in black.

**Drive wheel: aluminum rim, hub and spokes**

Black A9000100171-003

**Caster fork and bearing block**

Silver A9000100173-011  
 Black A9000100173-003

## Colour design of attachment parts according to customer request:

**Drive wheel: aluminum rim and design inlet<sup>1</sup> according to colour selection / hub and spokes black**

frame colour A9000100171-999 € 249.60

Selection from one-colour design (different from the frame colour) A9000100171-998 € 288.90

Colour: \_\_\_\_\_  
A9000100170-

**Caster wheel: 5" aluminum rim<sup>2</sup> (6-spoke design)**

frame colour A9000100174-999 € 75.20

Selection from one-colour design (different from the frame colour) A9000100174-998 € 87.10

Colour: \_\_\_\_\_  
A9000100170-

## Special colour request (Frame colour - One-colour design)

**Frame colour for single purchasing Selection from RAL colour palette**

A9000100170-100 € 414.40

**gloss level**

gloss  
 matt

**colour designation**

\_\_\_\_\_

**RAL number**

\_\_\_\_\_

**Note:** Attachments **not** available in design colour options

<sup>1</sup> Only in combination with Liteline hub

<sup>2</sup> Only in combination with caster wheel 5" aluminum rim 6-spoke-design good shock absorption A900050121\*, as well as with suspension fork 5-spoke-design A9000500160\*

## Frame

### Straight front end


A9020100010-K\*



### Front end V-shape

A9020100020-K\* € 257.00

from SB 33 cm  
Note: Footrest area with V-shape  
 $x = SB - (0.25 \times USL + 6 \text{ cm})$



## Knee angle

### Knee angle 105°

A9000100050-001

### Knee angle 95°

A9000100050-002 € 339.00

## Rear frame angle

### Rear frame angle 105°

A9010100100-001  
KP from 6 cm to 13 cm

## Back tubes & push handles

### Equipment without push handles

A9000401000



### Back tube with integrated handles

A9000401051-K € 130.40  
straight.



### Aluminium push handles fixed in back tube

A9000401050\* € 130.40  
from RH 23 cm



### Horizontal push handles screwed into the back tube

A9000401076 € 130.40  
can be unscrewed



### Safety push handles with continuous height adjustment

A9000401060\* € 216.40  
mount in the back tube  
from RH 30 cm with adaptable back  
from RH 35 cm with lordosis support  
from RH 43 cm with ergonomic backshell



### Safety push handles back-positioned length 45 cm

A9000401075\* € 372.90  
continuous height adjustment,  
Safety push handle holder mounted on the  
back cross bar  
not in combination with swing-away armrest



## Additional option

### Back tube with lordosis support

A9000400005-K € 114.40

- 1 cm
- 2 cm
- 3 cm



## Backrest

If not specified, the „folding backrest“ will be delivered.


### Folding backrest

A9000400011  
7-way angle adjustable,  
with integrated plastic-coated aluminum side  
parts



### Backrest folding function blocked

A9000400035  
Locking bolt for angle adjustment replaced by  
screw connection.  
Back angle 90° to the floor.



### Fixed backrest

A9000400020 **€ 99.10**  
screwed to aluminium side parts  
Back angle 90° to the floor  
Backrest cannot fold down

## Back equipment

If not specified, the „Adaptable back with back padding“ will be delivered.

### Adaptable back with back padding

A9000802010-K001  
Cover: Black, breathable, open surface  
Equipment with back cross bar



### Adaptable back with back padding „light“

A9000802320-K001  
Cover: Black, breathable, open surface  
Cushion thickness: 9 mm  
Equipment with back cross bar

A9000802320-K002  
Cover: Black, breathable, open surface  
Cushion thickness: 9 mm  
Equipment without back cross bar



### Adaptable back with back padding „Comfort“

A9000815101-K **€ 131.40**  
Cover: Black, breathable, open surface  
Front padding: 20 mm  
Rear padding: 5 mm  
Equipment with back cross bar

### Back rise + 3 - 5 cm

A9000802100 **€ 159.10**  
for the nominal back height dimension  
only in combination with height-adjustable push  
handles



### Equipment without back upholstery and back padding or back system


A9000802011 **€ - 53.00**  
Equipment with back cross bar

## Seat upholstery

If not specified, the „Body Contour“ seat upholstery will be delivered.


### Body Contour

A9000801025-K  
approved up to a body weight of 120 kg



### Open belt system

A9000801010-K **€ 79.90**





### Equipment without seat upholstery

A9000801050



## Seat cushion



### Seat cushion in rectangular contour Inlet: Standard

A9000805180-K\* € 118.30

### Seat cushion with radius contour Inlet: Standard

A9000800300-K\* € 163.10

### Seat cushion height:

3 cm       5 cm  
 4 cm       6 cm

### Seat cushion in rectangular contour Inlet: Anti-Decubitus

up to 30x30 cm  
A9000805204-K\* € 301.70

>30x30 up to 50x50 cm  
A9000805200-K\* € 437.70

### Seat cushion with radius contour Inlet: Anti-Decubitus

up to 30x30 cm  
A9000800311-K\* € 356.10

>30x30 up to 50x50 cm  
A9000800310-K\* € 482.90

### Seat cushion height:

4 cm       6 cm

### Cover:

- Black, breathable, open surface  
 Black, breathable, closed surface

If not specified, the cover material „Black, breathable, open surface“ will be delivered.

Seat cushion is provided

## Clothing guard

### Aluminum clothing guard removable via lock function

Adjustable following wheel arch  
Only in combination with folding backrest



in anthracite matt

A9000403210-K € 243.20

in frame colour (base colour for two-colour design)

A9000403211-K € 311.10

in design colour (only possible with frame colour in two-colour design)

A9000403212-K € 345.10

### Carbon clothing guard removable via lock function

A9000400012-K € 547.30

Adjustable following wheel arch  
Only in combination with folding backrest



## Armrests

### Armrest, height adjustable, integrated in the clothing guard

Including aluminum clothing guard  
with arm padding 25 cm long  
Removable via lock function  
Only in combination with folding backrest



in anthracite matt

A9000400008-K € 569.60

in frame colour (base colour for two-colour design)

A9000400009-K € 621.90

in design colour (only possible with frame colour in two-colour design)

A9000400010-K € 648.10

<p><b>Armrest, height adjustable, integrated in the clothing guard</b></p> <p><input type="checkbox"/> A9000400018-K      € 862.00</p> <p>Including carbon clothing guard with arm padding 25 cm long Removable via lock function Only in combination with folding backrest</p>	<p><b>Armrests folding up, removable</b></p> <p><input type="checkbox"/> A9000403141-K      € 358.50</p> <p>with arm padding 25 cm long Positioning: upper arm padding edge approx. 25 cm to 31 cm above the seat tube, adjustable in 1.5 cm steps not possible in combination with back-positioned safety push handles</p> 	<p><b>DESK side part with armrest 38 cm</b></p> <p><input type="checkbox"/> A9010403160      € 479.50</p> <p>removable, height-adjustable (Overall width + 4 cm) not possible in combination with aluminum / carbon clothing guard</p>
---	---	--


## Axle plate

Standard equipment with 12 mm sockets / axles

### 2-slot axle plate

A9000201000\*

Seat height back (with 24" wheel):  
short frame = approx. 37 cm - 43 cm  
long frame = approx. 39 cm - 45 cm



### 3-slot axle plate

A9000201020\*      € 174.70

Seat height back (with 24" wheel):  
short frame = approx. 40 cm - 46 cm  
long frame = approx. 42 cm - 48 cm



### Tipping point fixation wheel plate

A9000200006      € 163.80

Consisting of  
slotted nut tilt point fixing  
holding plate tilt point fixing  
including fastening material

## Mounting Scala-Mobil

Possible from SB 38 cm.  
For SB 38 cm - 40 cm the attachment must be made by Alber GmbH.  
The Scala-Mobil and the corresponding brackets can be ordered from Alber GmbH.

Wheelchair delivery to:

Alber GmbH       Fachhändler

## Wheel camber

\*Total width with tight-mounted handrim.

### 1° drive wheel camber

A9000202010\*

Total width = approx. SB + 20 cm\*

### 4° drive wheel camber

A9000202030\*

Total width = approx. SB + 26 cm\*

### 6° drive wheel camber

A9000202040\*

Total width = approx. SB + 30 cm\*

## E-Adaptations

### Preparation for the „e-motion“ / „e-motion DuoDrive“ mounting

A9000201200      € 502.90

Equipment: e-motion torque arm, 1/2" axle and socket length 50.8 mm (2")  
PRO ACTIV and Alber recommend mounting two anti-tipping supports. Please select under „Adaptations“.

### Preparation for the „e-fix“ mounting

A9000201210      € 502.90

Equipment: 1/2" axle, socket length 50.8 mm (2") and e-fix torque support.  
PRO ACTIV and Alber recommend mounting two anti-tipping supports. Please select under „Adaptations“.  
Please order Alber anti-tipping supports and corresponding components direct from Alber.

## Brake system

If not specified, the „Brake system push to lock with standard brake handle“ will be delivered.

**Recommendation:** For KP 7 cm or lower, we recommend the „Brake system push to lock with brake handle 20° outwards“.

### Knee lever brake with standard brake lever

A9000701010-001



### Knee lever brake with long brake lever

A9000701010-002



### Knee lever brake with low brake lever

A9000701010-003



<p><b>Knee lever brake system</b> <b>Brake lever angled 20°</b></p> <p><input checked="" type="checkbox"/> standard brake handle A9000700002</p> <p><input type="checkbox"/> long brake handle A9000701010-004</p> 	<p><b>Brake system push to lock</b> <b>Brake handle to flip-down</b></p> <p><input type="checkbox"/> A9000706008      <b>€ 158.60</b></p> <p>length: 9 cm</p> 	<p><b>Drum brake</b></p> <p><input type="checkbox"/> A9000702011-K*      <b>€ 474.00</b></p> <p>Operating from the push handle Total wheelchair width (see selection „Wheel camber“) + 3 cm Depending on the seat height, the long wheel plate is mounted in order to protect the drum brake anchors. Please note: - Only in combination with drive wheel aluminum 36-spoke rim wheel - Not in combination with one-hand drive, e-fix or e-motion</p> 
---	---	---

## Footrest

If not specified, the „continuous footrest“ will be delivered.

### Continuous footrest

A9016601010-K\*  
with tubular resting bracket,  
height adjustable +/- 3 cm  
If you want to equip the footrest with a „rubber coated footplate“ please select this under „Footrest accessories“.  
If x < 22 cm, the footrest is made parallel.



### Footrest folding up to one side

A9016601300-K\*      **€ 173.90**

with tubular resting bracket  
height adjustable +/- 3 cm  
can be swung away,  
cannot be removed  
V-shape from SB 40 cm  
If you want to equip the footrest with a „rubber coated footplate“ please select this under „Footrest accessories“.



### Footrest folds up to the rear with spring locking mechanism

A9016602005-K\*      **€ 287.40**

with tubular resting bracket  
for KW 105° short frame from USL 35 cm,  
long frame from USL 37 cm  
for KW 95° short frame from USL 34 cm,  
long frame from USL 36 cm  
If you want to equip the footrest with a „rubber coated footplate“ please select this under „Footrest accessories“.



### Footrest folds up to the rear

A9016602010-K\*      **€ 307.90**


with tubular resting bracket and 20 mm operating lever  
for KW 95° from USL 38 cm  
for KW 105° from USL 42 cm  
straight front end from SB 33 cm  
V-shape from SB 38 cm  
Support width from 20 cm  
If you want to equip the footrest with a „rubber coated footplate“ please select this under „Footrest accessories“.



### Divided footrest

A9016603200-K      **€ 336.80**

with rubber coated footplate,  
height adjustable +/- 3 cm, angle adjustment,  
can be pivoted outwards and folded up to the side, cannot be removed  
V-shape from SB 40 cm



### „Swing Away“ footrest

A9016604100-K      **€ 404.70**

with rubber coated footplate,  
separated in the middle,  
can be folded to the side, removable,  
can be folded up, height-adjustable,  
angle-adjustable  
Including frame reinforcement  
from USL 34 cm



## Footrest accessories

### Rubber coated footplate

A9000610000-K      **€ 125.90**

bolted, angle-adjustable  
Depth: 14 cm



### Footplate side stop

A9000610100      **€ 59.70**

adjustable  
only in combination with rubber coated footplate



## Adaptations

### Anti-tipping support

right side  
A9000300110-K\* € 169.10

left side  
A9000300111-K\* € 169.10

with tension spring and universal adapter  
pivotable

For optimum pivoting behaviour the  
seat slope should be > 3 cm



### 3" caster wheel

#### 3" skate wheel (80 mm)

A9000501110 € 70.80  
puncture-proof



### 4" caster wheel

#### 4" caster wheel „Standard“ solid rubber black

A9000501011-004  
(82A), puncture-proof



#### 4" skate wheel (100 mm)

A9000501120 € 70.80  
puncture-proof



#### 4" safety light-up wheel (100 mm)

A9000501220 € 94.80  
puncture-proof



#### 4" CorMORun caster wheel (100 mm)

A9000501042 € 173.00  
with replaceable rubber 6 mm,  
puncture-proof  
Rim colour: Black  
Incl. hair guard spacer (5.3 mm)



#### 4" Frog Legs caster wheel (100 mm)

A9000501230 € 206.30  
with replaceable rubber 6 mm,  
puncture-proof  
Colour: black, silver rim  
Incl. hair guard spacer (5.3 mm)



### 5" caster wheel

#### 5" caster wheel „Standard“ PU black

A9000501036-005  
(82A), puncture-proof



#### 5" caster wheel „Standard“ PU grey

A9000501031-005  
(82A), puncture-proof



#### 5" safety light-up wheel (125 mm)

A9000501221 € 94.80  
puncture-proof




<p><b>5" aluminum rim 3-spoke design „good shock absorption“</b></p> <p><input type="checkbox"/> silver anodised A9000501210 € 179.80</p> <p><input type="checkbox"/> black A9000501211 € 179.80</p> <p>with replaceable rubber 6 mm medium hard, puncture-proof Incl. hair guard spacer (3 mm)</p> 	<p><b>5" CorMORun caster wheel (125 mm)</b></p> <p><input type="checkbox"/> A9000501041 € 181.30</p> <p>with replaceable rubber 6 mm, puncture-proof Rim colour: Black Incl. hair guard spacer (5.3 mm)</p> 	<p><b>5" Frog Legs caster wheel (125 mm)</b></p> <p><input type="checkbox"/> A9000501232 € 206.30</p> <p>with replaceable rubber 6 mm, puncture-proof Colour: black, silver rim Incl. hair guard spacer (5.3 mm)</p> 
--	---	--

<p><b>6" caster wheel</b></p>		
<p><b>6" caster wheel „Standard“ PU grey</b></p> <p><input checked="" type="checkbox"/> A9000501010-006 (82A), puncture-proof</p> 	<p><b>6" caster wheel „pneumatic tyre“</b></p> <p><input type="checkbox"/> A9000501050 € 57.20</p> 	<p><b>6" caster wheel „Flexel-System“</b></p> <p><input type="checkbox"/> A9000501070 € 57.20</p> <p>puncture-proof</p> 

<p><b>Caster fork</b></p> <p>For caster fork colours, see „Colour of add-on parts“ If not specified, the „Standard caster fork with screw axle“ will be delivered.</p>		
<p><b>Caster fork Standard</b></p> <p><input checked="" type="checkbox"/> A9000500100-K001 with screwed axle</p> <p><input type="checkbox"/> A9000500100-K002 € 167.40 with quick-release axle</p> 	<p><b>Caster fork Heavy Duty for high loads</b> Extra weight approx. 100 - 150 gr.</p> <p><input type="checkbox"/> A9000510111-K001 € 54.00 with screwed axle</p> <p><input type="checkbox"/> A9000510111-K002 € 218.90 with quick-release axle</p> 	<p><b>Suspension fork</b> inkl. 5" aluminum rim „good shock absorption“</p> <p><input type="checkbox"/> A9000500160-K002 € 616.90 with quick-release axle</p> 
<p><b>Caster fork spare wheel „Standard“</b></p> <p><input type="checkbox"/> A9000500300-K € 296.70</p> <p>Caster fork Standard with quick-release axles</p> 	<p><b>Caster fork spare wheel „Heavy Duty“</b></p> <p><input type="checkbox"/> A9000500301-K € 337.90</p> <p>Caster fork Heavy Duty with quick-release axles</p> 	<p>The spare caster fork does not have a caster wheel. In case a caster wheel is requested, please put a cross on a second selection in the category „Caster wheels“.</p>

## Track distance

### Caster wheel track extension



A9000500400\*  
10 mm per side

€ 59.10

A9000500401\*  
15 mm per side

€ 59.10

A9000500410\*  
20 mm per side

€ 59.10

## Drive wheels

If not specified, the „drive wheel 36 spoke aluminum rim“ will be delivered.

### Drive wheel 18 spoke aluminum rim

24" (ETRTO 540)  
A9000203122  
with Liteline hub, radially spoked  
Not in combination with:  
- Puncture-proof tires  
- One-hand drive  
- Drum brake  
- Transport restraint connection  
- Handbike  
up to 70 kg body weight:

Sturz	1°	4°	6°	8°	10°
Gewicht					
bis 50 Kg	✓	✓	✓	✓	✓
51 - 60 Kg	✓	✓	✓	✓	
61 - 70 Kg	✓	✓			



### Drive wheel 36 spoke aluminum rim

24" (ETRTO 540)  
A9000203000-010  
with Liteline hub (left image),  
radially spoked  
up to 120 kg body weight  
 24" (ETRTO 540)  
A9000203000-020  
with aluminum front hub (right image),  
radially spoked  
over 120 kg body weight



### Drive wheel double rim 36-spokes

24" (ETRTO 540)  
A9000203100 € 157.80  
 25" (ETRTO 559)  
A9000203200 € 157.80  
 26" (ETRTO 590)  
A9000203300 € 157.80  
with aluminum front hub,  
radially spoked  
For use with high loads and body weight  
>120 kg



### Spinergy wheel „Spox Sport“ with 24 PBO spokes; rim, hub and spokes in black

24" (ETRTO 540)  
A9000203550 € 784.10  
 25" (ETRTO 559)  
A9000201559 € 784.10  
 26" (ETRTO 590)  
A9000203555 € 784.10



### Spinergy „Everyday Light Extreme“ with 12 PBO spokes; rim and hub in black, spokes in white

24" (ETRTO 540)  
A9000203554 € 957.20  
Recommendation: For heavy loads  
(e.g. higher body weight), we recommend the  
Spinergy „Spox Sport“ wheel with 24 PBO  
spokes.




### Spinergy „Everyday Light Extreme“ with 12 PBO spokes; rim, hub and spokes in black

24" (ETRTO 540)  
A9000203552 € 957.20  
Recommendation: For heavy loads  
(e.g. higher body weight), we recommend the  
Spinergy „Spox Sport“ wheel with 24 PBO  
spokes.



# Handrims

x = dimension handrim attachment



Handrim attachment general:

- x = 11 mm
- x = 19 mm

Handrim attachment Curve + Gekko handrim:

- x = 10 mm
- x = 18 mm

If not specified, the handrim is mounted wide (19 mm / 18 mm).

## Stainless steel Ø 20 mm

- 24" handrim for 24" wheel and 25" wheel A9000205010-004
- 26" handrim for 26" wheel A9000205010-006



## Aluminum Ø 19 mm silver anodised

- 24" handrim for 24" wheel and 25" wheel A9000205019-004



## Aluminum Ø 19 mm powder coated in frame colour

(basic colour with two-colour design)

- = Gloss surface
- = Matt surface

Note: Matt coating is not recommended due to the low grip.

- 24" handrim for 24" wheel and 25" wheel A9000205020-004 **€ 68.30**
- 26" handrim for 26" wheel A9000205020-006 **€ 68.30**



## Aluminum Ø 19 mm polished and hard-anodised

- 24" handrim for 24" wheel and 25" wheel **€ 164.90** A9000205030-004
- 26" handrim for 26" wheel **€ 164.90** A9000205030-006



## Titanium Ø 19 mm

- 24" handrim for 24" wheel and 25" wheel **€ 335.60** A9000205040-004
- 26" handrim for 26" wheel **€ 335.60** A9000205040-006



## Silicone handrim cover

- for handrim 24" A9000205100-004 **€ 232.80**
  - for handrim 26" A9000205100-006 **€ 232.80**
- Surface without structure, compatible with handrim Ø 19 mm or Ø 20 mm.



## Ultra-Grip „very grippy“

- 24" handrim for 24" wheel and 25" wheel **€ 402.40** A9000205070-004
- 26" handrim for 26" wheel **€ 402.40** A9000205070-006



## Ergo-Grip

- 24" handrim for 24" wheel and 25" wheel **€ 402.40** A9000205080-004
  - 26" handrim for 26" wheel **€ 402.40** A9000205080-006
- Attachment fitting to the rim

Note: Due to the complex geometry, concentricity deviations of +/- 2-3 mm are possible.



## Curve black anodised, ergonomic handrim profile

- 24" handrim for 24" wheel and 25" wheel Profile height 24 mm **€ 402.40** A9000205130-004
- 24" handrim for 24" wheel and 25" wheel Profile height 31 mm Size L **€ 402.40** A9000205140-004
- 26" handrim for 26" wheel Profile height 31 mm Size L **€ 402.40** A9000205140-006



<p><b>Curve silver anodised ergonomic handrim profile</b></p> <p><input type="checkbox"/> 24" handrim for 24" wheel and 25" wheel Profile height 31 mm    <b>€ 402.40</b> Size L A9000205160-004</p> 	<p><b>Curve powder-coated, ergonomic handrim profile</b></p> <p><input type="checkbox"/> = frame colour (gloss surface) <input type="checkbox"/> = black (gloss surface)</p> <p><input type="checkbox"/> 24" handrim for 24" wheel and 25" wheel Profile height 24 mm    <b>€ 467.60</b> A9000205135-004</p> <p><input type="checkbox"/> 24" handrim for 24" wheel and 25" wheel Profile height 31 mm    <b>€ 467.60</b> Size L A9000205145-004</p> <p><input type="checkbox"/> 26" handrim for 26" wheel Profile height 31 mm    <b>€ 467.60</b> Size L A9000205145-006</p> 	<p><b>Curve Grip, with Tetra-Grip coating, ergonomic handrim profile</b></p> <p><input type="checkbox"/> 24" handrim for 24" wheel and 25" wheel Profile height 31 mm    <b>€ 520.70</b> Size L A9000205150-004</p> 
<p><b>Quadro with Tetra-Grip coating</b></p> <p><input type="checkbox"/> 24" handrim for 24" wheel and 25" wheel    <b>€ 520.70</b> A9000205170-004</p> <p><input type="checkbox"/> 25" handrim for 25" wheel    <b>€ 520.70</b> A9000205170-005</p> <p><input type="checkbox"/> 26" handrim for 26" wheel    <b>€ 520.70</b> A9000205170-006</p> 	<p><b>Gekko black anodised</b></p> <p><input type="checkbox"/> 24" handrim for 24" wheel and 25" wheel    <b>€ 520.70</b> A9000205175-004</p> <p><input type="checkbox"/> 25" handrim for 25" wheel    <b>€ 520.70</b> A9000205175-005</p> <p><input type="checkbox"/> 26" handrim for 26" wheel    <b>€ 520.70</b> A9000205175-006 with a robust and replaceable silicone lip, fully integrated into the profile Profile height 31 mm</p> 	<p><b>Gekko silver anodised</b></p> <p><input type="checkbox"/> 24" handrim for 24" wheel and 25" wheel    <b>€ 520.70</b> A9000205185-004 with a robust and replaceable silicone lip, fully integrated into the profile Profile height 31 mm</p> 
<p><b>Gekko S black anodised</b></p> <p><input type="checkbox"/> 24" handrim for 24" wheel and 25" wheel    <b>€ 520.70</b> A9000205195-004 with a robust and replaceable silicone lip, fully integrated into the profile Profile height 24 mm</p> 	<p><b>Equipment without handrims</b></p> <p><input checked="" type="checkbox"/> A9000205990</p>	

## Spoke guard

### Finger- and spoke guard „transparent“

for 24" wheel  
A9000206300-004      € 207.10

for 26" wheel  
A9000206300-006      € 207.10

Not in combination with Spinergy drive wheels.



## Tyres

If not specified, the „Schwalbe RightRun in black“ tyres will be delivered.

### Schwalbe „RightRun“ High-pressure tyres

24x1" (540-25)  
A9000204001-004

Colour: Black  
7.5 bar



### Schwalbe „RightRun“ High-pressure tyres

24x1" (540-25)  
A9000204000-004

25x1" (559-25)  
A9000204000-005

26x1" (590-25)  
A9000204000-006

Colour: Grey, black wall  
7.5 bar



### Schwalbe „DownTown“ Profiled tyres

24x1" (540-25)  
A9000204010-004

7.5 bar



### Schwalbe „Marathon Plus Evolution“ with „smart guard“ puncture protection

24x1" (540-25)  
A9000204030-014      € 49.20

25x1" (559-25)  
A9000204030-015      € 49.20

26x1" (590-25)  
A9000204030-016      € 49.20

Colour: Black  
10 bar



### Schwalbe „Marathon Plus“ with „smart guard“ puncture protection

24x1" (540-25)  
A9000204030-004      € 49.20

Colour: Grey, black wall  
10 bar



### Schwalbe „One“

24x1" (540-25)  
A9000204123-004      € 93.00

Colour: Black  
9 bar



### Puncture-proof „solid rubber tyres“

24x1" (540-25)  
A9000204201-004      € 122.70

Colour: Black

Not in combination with Spinergy drive wheels.



## Accessories

### 24" MTB spare wheel set

with double rim (ETRTO 507), aluminum front hub, quick-release axle and stainless steel handrim  
Mountain bike tyres, selectively:

A9000203410 € 703.40  
24x1.75" Schwalbe „Marathon Plus"  
Brake setting compatible with 24x1" wheel



A9000203413 € 703.40  
24x2.1" Schwalbe „SMART SAM Sport"  
Brake setting not compatible with 24x1" wheel




NOTE: We charge a surcharge for the use of drum brakes

A9000200004 € 232.20

### Through axle with assistance, Large push button

A9000200096 € 105.80  
Interlock dimension: 122 / 110 mm  
Axle diameter: 12 mm  
Not in combination with one-hand drive



### Hand air compressor with pressure gauge

A9000901100 € 106.30  
up to 14 bar



### Crutch holder with flexible nylon cup

A9000902310\* € 137.40



### Length-adjustable safety belt with buckle

A9000902100 € 80.60



### Adjustable safety belt with hook and fastener

A9000902110 € 64.50



### SIGMA lighting set

A9000903270 € 107.90

Consisting of front light and rear lights, chargeable via Micro USB socket, incl. Micro USB charging cable  
Lighting set compliant with road traffic regulations



### Transport bag for drive wheels

for 24" wheel  
A9000903530-004 € 77.80  
Ø 76 cm x 2 cm

for 26" wheel  
A9000903530-006 € 77.80  
Ø 83 cm x 2 cm



## Separate order sheets

### Category „Rehabilitation technology accessories”

e.g. Seat & backrest systems, one-hand drive or transport restraint systems can be found at:  
[www.proactiv-gmbh.com/downloads/](http://www.proactiv-gmbh.com/downloads/)

### Additional products

such as hand bikes or traction & attachment devices can be found at:  
[www.proactiv-gmbh.com/downloads/](http://www.proactiv-gmbh.com/downloads/)

### Additional sheets

such as heavy-duty versions etc. can be found at:  
[www.proactiv-gmbh.com/downloads/](http://www.proactiv-gmbh.com/downloads/)

## Services

Consultation, clarification of the installation situation and measurement through PRO ACTIV

A9990000000 € 354.50

Handover and instruction on the wheelchair by PRO ACTIV

A9990000006 € 354.50

Notes and remarks:

**Location, Date**

**Signature End customer/dealer**

The signature field is intended for the consultant and insured person. This confirms that the intended equipment is intended for the insured person.  
- 17 - 2026-03-03 / The prices are plus statutory value-added tax and shipping costs  = included in the basic price  = Position selected  
Order sheet SPEEDY A2 Rev. 07 UK