

Im Hofstätt 11
D-72359 Dotternhausen
Tel +49 (0)7427 9480-0
Fax +49 (0)7427 9480-7025
e-mail info@proactiv-gmbh.de
www.proactiv-gmbh.de



SPEEDY A1

Active easy running rigid frame wheelchair

Maximum user weight limit 120 kg

A9010100001 € 4,010.00

☐ Quotation date: _____

Consultation was carried out by
Staff at the medical supplies shop: _____

Enquiry no.: _____

☐ Order date: _____

according to GTC PRO ACTIV Reha-Technik GmbH

Order no.: _____

Medical suppliers / company address:

Cust. no.: _____
Name: _____
Street: _____
Postcode / city: _____
Phone: _____
Fax: _____
E-Mail: _____
Contact: _____

Shipping address: ☐ same as company address
☐ Other (see below)

Name: _____
Street: _____
Postcode / city: _____

Customer name / commission:

ID / Comm. No.: _____
Surname, forename: _____
Street: _____
Postcode / city: _____
Phone: _____
Fax: _____
E-Mail: _____
D.o.B.: _____
Paying organisation: _____

Invoicing address: ☐ same as company address
☐ Other (see below)

Cust. no.: _____
Name: _____
Street: _____
Postcode / city: _____
Phone: _____
Fax: _____
E-Mail: _____
Contact: _____

Basic size


Body size:

Body height: _____ cm

Weight: _____ kg


Lower leg length: _____ cm
(anatomic dimension: hollow of the knee to the foot sole)

Indication: _____



This information is necessary for checking. If the information is missing, a check of the relevant items cannot be made!

Wheelchair dimensions:



Seat width **SB**

☐ 30	☐ 31	☐ 32	☐ 33	☐ 34	☐ 35	☐ 36	☐ 37	☐ 38	☐ 39
☐ 40	☐ 41	☐ 42	☐ 43	☐ 44	☐ 45	☐ 46	☐ 47	☐ 48	☐ 49
☐ 50									

Seat depth **ST**

Short frame (43 cm)					Long frame (47 cm)				
☐ 36	☐ 37	☐ 38	☐ 39		☐ 40	☐ 41	☐ 42	☐ 43	☐ 44
☐ 40	☐ 41	☐ 42	☐ 43		☐ 45	☐ 46	☐ 47	☐ 48	

Back height **RH**

☐ 20	☐ 21	☐ 22	☐ 23	☐ 24	☐ 25	☐ 26	☐ 27	☐ 28	☐ 29
☐ 30	☐ 31	☐ 32	☐ 33	☐ 34	☐ 35	☐ 36	☐ 37	☐ 38	☐ 39
☐ 40	☐ 41	☐ 42	☐ 43	☐ 44	☐ 45	☐ 46	☐ 47	☐ 48	

Lower leg length **USL**

Short frame					Long frame				
☐ 33	☐ 34	☐ 35	☐ 36	☐ 37	☐ 35	☐ 36	☐ 37	☐ 38	☐ 39
☐ 38	☐ 39	☐ 40	☐ 41	☐ 42	☐ 40	☐ 41	☐ 42	☐ 43	☐ 44
☐ 43	☐ 44	☐ 45	☐ 46	☐ 47	☐ 45	☐ 46	☐ 47	☐ 48	☐ 49
☐ 48	☐ 49	☐ 50			☐ 50	☐ 51	☐ 52		

Seat height back **SHh**
(with 24" drive wheels)

Short frame				Long frame			
☐ 37	☐ 38	☐ 39	☐ 40	☐ 39	☐ 40	☐ 41	☐ 42
☐ 41	☐ 42	☐ 43		☐ 43	☐ 44	☐ 45	

Seat height front **SHv**
(with 5" caster wheels)

Short frame					Long frame				
☐ 44	☐ 45	☐ 46	☐ 47	☐ 48	☐ 46	☐ 47	☐ 48	☐ 49	☐ 50
☐ 49	☐ 50	☐ 51	☐ 52		☐ 51	☐ 52	☐ 53	☐ 54	

Tipping point **KP**
(inner edge of rear tube to centre of drive wheel axle)

☐ 4	☐ 5	☐ 6	☐ 7	☒ 8	☐ 9	☐ 10	☐ 11	☐ 12	☐ 13
----------------------------	----------------------------	----------------------------	----------------------------	---------------------------------------	----------------------------	-----------------------------	-----------------------------	-----------------------------	-----------------------------

Please tick the required dimensions (all dimensions are in cm)

The underlined values are the standard settings. If not specified, the standard value will be delivered.

The SHh and SHv values might be a bit different from the specified values, depending on the wheel size and camber. Individual dimensions available on request.

We recommend a seat height difference (difference between SHv and SHh) between 5 cm and 12 cm.

Frame colour - One-colour design		
Information on current colour designs and example images can be found at: https://www.proactiv-gmbh.com/colours/		
Frame colour: Selection from colour palette		
Frame colour will be specified in the order ☒ A9000100170-000	Raw, blasted, not coated ☒ A9000100170-020 Note: In combination with moisture e.g. hand perspiration, can affect dark aluminum discoloration.	
Powder colours		
white RAL 9016 ☒ gloss A9000100170-023 ☐ matt A9000100170-123 € 43.40	brillant silver ☒ gloss A9000100170-011	anthracite RAL 7016 ☒ gloss A9000100170-001 ☐ matt A9000100170-101 € 43.40
black RAL 9005 ☒ gloss A9000100170-003 ☐ matt A9000100170-019 € 43.40	pure orange RAL 2004 ☒ gloss A9000100170-022 ☐ matt A9000100170-122 € 43.40	pink magenta RAL 4010 ☒ gloss A9000100170-037 ☐ matt A9000100170-137 € 43.40
pastel blue RAL 5024 ☒ gloss A9000100170-036 ☐ matt A9000100170-136 € 43.40	lemon green ☒ gloss A9000100170-005 ☐ matt A9000100170-105 € 43.40	yellow RAL 1023 ☒ gloss A9000100170-046 ☐ matt A9000100170-146 € 43.40
red RAL 3001 ☒ gloss A9000100170-047 ☐ matt A9000100170-147 € 43.40	lagoon blue ☒ gloss A9000100170-052 ☐ matt A9000100170-152 € 43.40	athletic grey RAL 7046 ☒ gloss A9000100170-050 ☐ matt A9000100170-150 € 43.40
muddy grey RAL 7033 ☒ gloss A9000100170-051 ☐ matt A9000100170-151 € 43.40	peach RAL 3022 ☒ gloss A9000100170-038 ☐ matt A9000100170-138 € 43.40	
Powder colours „dormant“		
deep black effect silver ☒ gloss A9000100170-007	sky blue ☒ gloss A9000100170-017	navy blue ☒ gloss A9000100170-045 ☐ matt A9000100170-145 € 43.40
electric green ☒ gloss A9000100170-015 ☐ matt A9000100170-115 € 43.40	violet ☒ gloss A9000100170-044	fire red ☒ gloss A9000100170-018 ☐ matt A9000100170-118 € 43.40
Frame colour - One-colour design special colours „liquid coating“		
€ 349.60		
storm (ice blue) ☐ gloss A9000100170-500 ☐ matt A9000100170-501	prairie (rose gold) ☐ gloss A9000100170-504 ☐ matt A9000100170-505	berry ☐ gloss A9000100170-526 ☐ matt A9000100170-527
blue ocean ☐ gloss A9000100170-522 ☐ matt A9000100170-523	emerald ☐ gloss A9000100170-514 ☐ matt A9000100170-515	chili ☐ gloss A9000100170-516 ☐ matt A9000100170-517

Frame colour - Two-colour design

Note: In the two-colour design, both colours are the same gloss level as the design colour.

☐ **Two-colour design with one base colour (on the rear frame) and one design colour (on the front frame)**

€ 268.50

Please select a colour below:

Base colour (on the rear frame)

white RAL 9016

☐ A9000100170-200

brilliant silver

☐ A9000100170-201

anthracite RAL 7016

☐ A9000100170-202

black RAL 9005

☐ A9000100170-203

pastel blue RAL 5024

☐ A9000100170-206

lemon green

☐ A9000100170-207

yellow RAL 1023

☐ A9000100170-208

lagoon blue

☐ A9000100170-212

athletic grey RAL 7046

☐ A9000100170-210

muddy grey RAL 7033

☐ A9000100170-211

Design colour (on the front frame)

black RAL 9005

☐ gloss A9000100170-259

☐ matt A9000100170-260

lemon green

☐ gloss A9000100170-250

☐ matt A9000100170-251

chili

☐ gloss A9000100170-253

☐ matt A9000100170-254

pure orange RAL 2004

☐ gloss A9000100170-255

☐ matt A9000100170-256

emerald

☐ gloss A9000100170-274

☐ matt A9000100170-275

storm (ice blue)

☐ gloss A9000100170-264

☐ matt A9000100170-265

prairie (rose gold)

☐ gloss A9000100170-266

☐ matt A9000100170-267

berry

☐ gloss A9000100170-284

☐ matt A9000100170-285

blue ocean

☐ gloss A9000100170-280

☐ matt A9000100170-281

Colour - Attachment parts

Note: If not specified, the caster fork and bearing block will be mounted in black.

Drive wheel: aluminum rim, hub and spokes

☐ Black A9000100171-003

Caster fork and bearing block

☐ Silver A9000100173-011

☐ Black A9000100173-003

Colour design of attachment parts according to customer request:

Drive wheel: aluminum rim and design inlet¹ according to colour selection / hub and spokes black

☐ frame colour A9000100171-999 € 249.60

☐ Selection from one-colour design (different from the frame colour) A9000100171-998 € 288.90

Colour: _____

A9000100170- ☐ ☐ ☐

Caster wheel: 5" aluminum rim² (6-spoke design)

☐ frame colour A9000100174-999 € 75.20

☐ Selection from one-colour design (different from the frame colour) A9000100174-998 € 87.10

Colour: _____

A9000100170- ☐ ☐ ☐

Special colour request (Frame colour - One-colour design)

Frame colour for single purchasing Selection from RAL colour palette

☐ A9000100170-100 € 414.40

gloss level

☐ gloss

☐ matt

colour designation

RAL number

Note: Attachments **not** available in design colour options


¹ Only in combination with Liteline hub

² Only in combination with caster wheel 5" aluminum rim 6-spoke-design good shock absorption A9000500121*, as well as with suspension fork 5-spoke-design A9000500160*

Frame

Straight front end


☒ A9010100010-K*



Front end V-shape

☐ A9010100020-K* € 257.00

from SB 33 cm
Note: Footrest area with V-shape
 $x = SB - (0.25 \times USL + 6 \text{ cm})$



Knee angle

Knee angle 105°

☒ A9000100050-001

Knee angle 95°

☐ A9000100050-002 € 339.00

Rear frame angle

Rear frame angle 105°

☒ A9010100100-001
KP from 6 cm to 13 cm

Rear frame angle 95°

☐ A9010100100-002 € 129.20
KP from 4 cm to 10 cm

Back tubes & push handles

Equipment without push handles

☒ A9000401000



Back tube with integrated handles

☐ A9000401051-K € 130.40
straight.



Aluminium push handles fixed in back tube

☐ A9000401050* € 130.40
from RH 23 cm



Horizontal push handles screwed into the back tube

☐ A9000401076 € 130.40
can be unscrewed



Safety push handles with continuous height adjustment

☐ A9000401060* € 216.40
mount in the back tube
from RH 30 cm with adaptable back
from RH 35 cm with lordosis support
from RH 43 cm with ergonomic backshell



Safety push handles back-positioned length 45 cm

☐ A9000401075* € 372.90
continuous height adjustment,
Safety push handle holder mounted on the back cross bar
Not in combination with swing-away armrest



Additional option

Back tube with lordosis support

☐ A9000400005-K € 114.40

- ☐ 1 cm
☐ 2 cm
☐ 3 cm



Backrest

If not specified, the „folding backrest“ will be delivered.


Folding backrest

☒ A9000400011
7-way angle adjustable,
with integrated plastic-coated aluminum side
parts



Backrest folding function blocked

☒ A9000400035
Locking bolt for angle adjustment replaced by
screw connection.
Back angle 90° to the floor.



Fixed backrest

☐ A9000400020 **€ 99.10**
screwed to aluminium side parts
Back angle 90° to the floor
Backrest cannot fold down

Back equipment

If not specified, the „Adaptable back with back padding“ will be delivered.

Adaptable back with back padding

☒ A9000802010-K001
Cover: Black, breathable, open surface
Equipment with back cross bar



Adaptable back with back padding „light“

☒ A9000802320-K001
Cover: Black, breathable, open surface
Cushion thickness: 9 mm
Equipment with back cross bar

☒ A9000802320-K002
Cover: Black, breathable, open surface
Cushion thickness: 9 mm
Equipment without back cross bar



Adaptable back with back padding „Comfort“

☐ A9000815101-K **€ 131.40**
Cover: Black, breathable, open surface
Front padding: 20 mm
Rear padding: 5 mm
Equipment with back cross bar

Back rise + 3 - 5 cm

☐ A9000802100 **€ 159.10**
for the nominal back height dimension
Only in combination with height-adjustable
push handles



Equipment without back upholstery and back padding or back system


☐ A9000802011 **€ - 53.00**
Equipment with back cross bar

Seat upholstery

If not specified, the „Body Contour“ seat upholstery will be delivered.


Body Contour

☒ A9000801025-K
approved up to a body weight of 120 kg



Open belt system

☐ A9000801010-K **€ 79.90**





Equipment without seat upholstery

☒ A9000801050



Seat cushion



Seat cushion in rectangular contour Inlet: Standard

☐ A9000805180-K* € 118.30

Seat cushion with radius contour Inlet: Standard

☐ A9000800300-K* € 163.10

Seat cushion height:

☐ 3 cm ☐ 5 cm
☐ 4 cm ☐ 6 cm

Seat cushion in rectangular contour Inlet: Anti-Decubitus

☐ up to 30x30 cm A9000805204-K* € 301.70

☐ >30x30 up to 50x50 cm A9000805200-K* € 437.70

Seat cushion with radius contour Inlet: Anti-Decubitus

☐ up to 30x30 cm A9000800311-K* € 356.10

☐ >30x30 up to 50x50 cm A9000800310-K* € 482.90

Seat cushion height:

☐ 4 cm ☐ 6 cm

Cover:

☐ Black, breathable, open surface

☐ Black, breathable, closed surface

If not specified, the cover material „Black, breathable, open surface“ will be delivered.

☐ Seat cushion is provided

Clothing guard

Aluminum clothing guard removable via lock function

Adjustable following wheel arch

Only in combination with folding backrest



in anthracite matt

☐ A9000403210-K € 243.20

in frame colour (base colour for two-colour design)

☐ A9000403211-K € 311.10

in design colour (only possible with frame colour in two-colour design)

☐ A9000403212-K € 345.10

Carbon clothing guard removable via lock function

☐ A9000400012-K € 547.30

Adjustable following wheel arch

Only in combination with folding backrest



Armrests

Armrest, height adjustable, integrated in the clothing guard

Including aluminum clothing guard

with arm padding 25 cm long

Removable via lock function

Only in combination with folding backrest



in anthracite matt

☐ A9000400008-K € 569.60

in frame colour (base colour for two-colour design)

☐ A9000400009-K € 621.90

in design colour (only possible with frame colour in two-colour design)

☐ A9000400010-K € 648.10


Armrest, height adjustable, integrated in the clothing guard ☐ A9000400018-K € 862.00 Including carbon clothing guard with arm padding 25 cm long Removable via lock function Only in combination with folding backrest	Armrests folding up, removable ☐ A9000403141-K € 358.50 with arm padding 25 cm long Positioning: upper arm padding edge approx. 25 cm to 31 cm above the seat tube, adjustable in 1.5 cm steps Not possible in combination with back-positioned safety push handles 	DESK side part with armrest 38 cm ☐ A9010403160 € 479.50 removable, height-adjustable (Overall width + 4 cm) Not possible in combination with aluminum / carbon clothing guard
--	---	--

Axle plate

Standard equipment with 12 mm sockets / axles


2-slot axle plate

☐ A9000201000*
Seat height back (with 24" wheel):
short frame = approx. 37 cm - 43 cm
long frame = approx. 39 cm - 45 cm



3-slot axle plate

☐ A9000201020* € 174.70
Seat height back (with 24" wheel):
short frame = approx. 40 cm - 46 cm
long frame = approx. 42 cm - 48 cm




Extended wheel plate

☐ A9000201010* € 174.70
Wheelbase extension by 7.5 cm
short frame = approx. 37 cm - 43 cm
long frame = approx. 39 cm - 45 cm



Extended bike plate

☐ A9000201030* € 315.20
+ 4° camber integrated, adjustable track
short frame = approx. 37 cm - 43 cm
long frame = approx. 39 cm - 45 cm



Tipping point fixation wheel plate

☐ A9000200006 € 163.80
Consisting of
slotted nut tilt point fixing
holding plate tilt point fixing
including fastening material

Wheel camber

*Total width with tight-mounted handrim.

1° drive wheel camber

☐ A9000202010*
Total width = approx. SB + 20 cm*

4° drive wheel camber

☐ A9000202030*
Total width = approx. SB + 26 cm*

6° drive wheel camber

☐ A9000202040*
Total width = approx. SB + 30 cm*

E-Adaptations

Preparation for the „e-motion“ / „e-motion DuoDrive“ mounting

☐ A9000201200 € 502.90
Equipment: e-motion torque arm, 1/2" axle and socket length 50.8 mm (2")
PRO ACTIV and Alber recommend mounting two anti-tipping supports. Please select under „Adaptations“.

Preparation for the „e-fix" mounting

☐ A9000201210 € 502.90
Equipment: 1/2" axle, socket length 50.8 mm (2") and e-fix torque support.
PRO ACTIV and Alber recommend mounting two anti-tipping supports. Please select under „Adaptations“.
Please order Alber anti-tipping supports and corresponding components direct from Alber.

Brake system

If not specified, the „Brake system push to lock with standard brake handle" will be delivered.

Recommendation: For KP 7 cm or lower, we recommend the „Brake system push to lock with brake handle 20° outwards“.

Knee lever brake with standard brake lever

☐ A9000701010-001



Knee lever brake with long brake lever

☐ A9000701010-002



Knee lever brake with low brake lever

☐ A9000701010-003



Knee lever brake system Brake lever angled 20° ☐ standard brake handle A9000700002 ☐ long brake handle A9000701010-004	Brake system push to lock Brake handle to flip-down ☐ A9000706008 € 158.60 length: 9 cm	Drum brake ☐ A9000702010-K* € 474.00 Operating from the push handle Total wheelchair width (see selection „Wheel camber“) + 3 cm Depending on the seat height, the long wheel plate is mounted in order to protect the drum brake anchors. Please note: - Only in combination with drive wheel aluminum 36-spoke rim wheel - Not in combination with one-hand drive, e-fix or e-motion 
---	--	--

Footrest

If not specified, the „continuous footrest“ will be delivered.

Continuous footrest

☐ A9016601010-K*
with tubular resting bracket,
height adjustable +/- 3 cm
If you want to equip the footrest with a „rubber
coated footplate“ please select this under
„Footrest accessories“.
If x < 22 cm, the footrest is made parallel.



Footrest folding up to one side

☐ A9016601300-K* **€ 173.90**
with tubular resting bracket
height adjustable +/- 3 cm
can be swung away,
cannot be removed
V-shape from SB 40 cm
If you want to equip the footrest with a „rubber
coated footplate“ please select this under
„Footrest accessories“.



Footrest folds up to the rear with spring locking mechanism

☐ A9016602005-K* **€ 287.40**
with tubular resting bracket
for KW 105° short frame from USL 35 cm,
long frame from USL 37 cm
for KW 95° short frame from USL 34 cm,
long frame from USL 36 cm
If you want to equip the footrest with a „rubber
coated footplate“ please select this under
„Footrest accessories“.



Footrest folds up to the rear

☐ A9016602010-K* **€ 307.90**
with tubular resting bracket and 20 mm
operating lever
for KW 95° from USL 38 cm
for KW 105° from USL 42 cm
straight front end from SB 33 cm
V-shape from SB 38 cm
If you want to equip the footrest with a „rubber
coated footplate“ please select this under
„Footrest accessories“.



Divided footrest

☐ A9016603200-K **€ 336.80**
with rubber coated footplate,
height adjustable +/- 3 cm, angle adjustment,
can be pivoted outwards and folded up to the
side, cannot be removed
V-shape from SB 40 cm



„Swing Away“ footrest

☐ A9016604100-K **€ 404.70**
with rubber coated footplate,
separated in the middle,
can be folded to the side, removable,
can be folded up, height-adjustable,
angle-adjustable
Including frame reinforcement
from USL 34 cm



Footrest accessories

Rubber coated footplate

☐ A9000610000-K **€ 125.90**
bolted, angle-adjustable
Depth: 14 cm



Footplate side stop

☐ A9000610100 **€ 59.70**
adjustable
only in combination with rubber coated
footplate



Adaptations

Anti-tipping support

☐ right side A9000300110-K* € 169.10

☐ left side A9000300111-K* € 169.10

with tension spring and universal adapter
pivotable

For optimum pivoting behaviour the
seat slope should be > 3 cm



Transit wheels

☐ A9000301110-K* € 329.60

with universal adapter, removable
not in combination with anti-tipping support or
tilting clamp



Tilting clamp

☐ right side A9000302110* € 140.10

☐ left side A9000302111* € 140.10

with universal adapter and Quick pin
not in combination with anti-tipping support,
right side or left side



3" caster wheel

3" skate wheel (80 mm)

☐ A9000501110 € 70.80

puncture-proof



4" caster wheel

4" caster wheel „Standard“ solid rubber black

☐ A9000501011-004
(82A), puncture-proof



4" skate wheel (100 mm)

☐ A9000501120 € 70.80

puncture-proof



4" safety light-up wheel (100 mm)

☐ A9000501220 € 94.80

puncture-proof



4" CorMORun caster wheel (100 mm)

☐ A9000501042 € 173.00

with replaceable rubber 6 mm,
puncture-proof
Rim colour: Black
Incl. hair guard spacer (5.3 mm)



4" Frog Legs caster wheel (100 mm)

☐ A9000501230 € 206.30

with replaceable rubber 6 mm,
puncture-proof
Colour: black, silver rim
Incl. hair guard spacer (5.3 mm)



5" caster wheel

5" caster wheel „Standard“ PU black

☐ A9000501036-005
(82A), puncture-proof



5" caster wheel „Standard“ PU grey

☐ A9000501031-005
(82A), puncture-proof



5" safety light-up wheel (125 mm)

☐ A9000501221 € 94.80
puncture-proof



5" aluminum rim 3-spoke design „good shock absorption“ ☐ silver anodised A9000501210 € 179.80 ☐ black A9000501211 € 179.80 with replaceable rubber 6 mm medium hard, puncture-proof Incl. hair guard spacer (3 mm) 	5" CorMORun caster wheel (125 mm) ☐ A9000501041 € 181.30 with replaceable rubber 6 mm, puncture-proof Rim colour: Black Incl. hair guard spacer (5.3 mm) 	5" Frog Legs caster wheel (125 mm) ☐ A9000501232 € 206.30 with replaceable rubber 6 mm, puncture-proof Colour: black, silver rim Incl. hair guard spacer (5.3 mm) 
6" caster wheel		
6" caster wheel „Standard“ PU grey ☒ A9000501010-006 (82A), puncture-proof 	6" caster wheel „pneumatic tyre“ ☐ A9000501050 € 57.20 	6" caster wheel „Flexel-System“ ☐ A9000501070 € 57.20 puncture-proof 
Caster fork		
For caster fork colours, see „Colour of add-on parts“ If not specified, the „Standard caster fork with screw axle“ will be delivered.		
Caster fork Standard ☒ A9000500100-K001 with screwed axle ☐ A9000500100-K002 € 167.40 with quick-release axle 	Caster fork Heavy Duty for high loads Extra weight approx. 100 - 150 gr. ☐ A9000510111-K001 € 54.00 with screwed axle ☐ A9000510111-K002 € 218.90 with quick-release axle 	Suspension fork inkl. 5" aluminum rim „good shock absorption“ ☐ A9000500160-K002 € 616.90 with quick-release axle 
Caster fork spare wheel „Standard“ ☐ A9000500300-K € 296.70 Caster fork Standard with quick-release axles 	Caster fork spare wheel „Heavy Duty“ ☐ A9000500301-K € 337.90 Caster fork Heavy Duty with quick-release axles 	The spare caster fork does not have a caster wheel. In case a caster wheel is requested, please put a cross on a second selection in the category „Caster wheels“.

Track distance

Caster wheel track extension



☐ A9000500400*
10 mm per side

€ 59.10

☐ A9000500401*
15 mm per side

€ 59.10

☐ A9000500410*
20 mm per side

€ 59.10

Drive wheels

If not specified, the „drive wheel 36 spoke aluminum rim“ will be delivered.

Drive wheel 18 spoke aluminum rim

☐ 24" (ETRTO 540)
A9000203122
with Liteline hub, radially spoked
Not in combination with:
- Puncture-proof tires
- One-hand drive
- Drum brake
- Transport restraint connection
- Handbike
up to 70 kg body weight:

Sturz	1°	4°	6°	8°	10°
Gewicht					
bis 50 Kg	✓	✓	✓	✓	✓
51 - 60 Kg	✓	✓	✓	✓	
61 - 70 Kg	✓	✓			



Drive wheel 36 spoke aluminum rim

☐ 24" (ETRTO 540)
A9000203000-010
with Liteline hub (left image),
radially spoked
up to 120 kg body weight
☐ 24" (ETRTO 540)
A9000203000-020
with aluminum front hub (right image),
radially spoked
over 120 kg body weight



Drive wheel double rim 36-spokes

☐ 24" (ETRTO 540)
A9000203100 € 157.80
☐ 25" (ETRTO 559)
A9000203200 € 157.80
☐ 26" (ETRTO 590)
A9000203300 € 157.80
with aluminum front hub,
radially spoked
For use with high loads and body weight
>120 kg



Spinergy wheel „Spox Sport“ with 24 PBO spokes; rim, hub and spokes in black

☐ 24" (ETRTO 540)
A9000203550 € 784.10
☐ 25" (ETRTO 559)
A9000201559 € 784.10
☐ 26" (ETRTO 590)
A9000203555 € 784.10



Spinergy „Everyday Light Extreme“ with 12 PBO spokes; rim and hub in black, spokes in white

☐ 24" (ETRTO 540)
A9000203554 € 957.20
Recommendation: For heavy loads
(e.g. higher body weight), we recommend the
Spinergy „Spox Sport“ wheel with 24 PBO
spokes.




Spinergy „Everyday Light Extreme“ with 12 PBO spokes; rim, hub and spokes in black

☐ 24" (ETRTO 540)
A9000203552 € 957.20
Recommendation: For heavy loads
(e.g. higher body weight), we recommend the
Spinergy „Spox Sport“ wheel with 24 PBO
spokes.



Handrims

x = dimension handrim attachment



Handrim attachment
general:

- ☐ x = 11 mm
☐ x = 19 mm

Handrim attachment
Curve + Gekko handrim:

- ☐ x = 10 mm
☐ x = 18 mm

If not specified, the handrim is mounted wide (19 mm / 18 mm).

Stainless steel Ø 20 mm

- ☐ 24" handrim for
24" wheel and 25" wheel
A9000205010-004
- ☐ 26" handrim for
26" wheel
A9000205010-006



Aluminum Ø 19 mm silver anodised

- ☐ 24" handrim for
24" wheel and 25" wheel
A9000205019-004



Aluminum Ø 19 mm powder coated in frame colour

(basic colour with two-colour design)

- ☐ = Gloss surface
☐ = Matt surface

Note: Matt coating is not recommended due to the low grip.

- ☐ 24" handrim for
24" wheel and 25" wheel € 68.30
A9000205020-004
- ☐ 26" handrim for
26" wheel € 68.30
A9000205020-006



Aluminum Ø 19 mm polished and hard-anodised

- ☐ 24" handrim for
24" wheel and 25" wheel € 164.90
A9000205030-004
- ☐ 26" handrim for
26" wheel € 164.90
A9000205030-006




Titanium Ø 19 mm

- ☐ 24" handrim for
24" wheel and 25" wheel € 335.60
A9000205040-004
- ☐ 26" handrim for
26" wheel € 335.60
A9000205040-006



Silicone handrim cover

- ☐ for handrim 24"
A9000205100-004 € 232.80
- ☐ for handrim 26"
A9000205100-006 € 232.80
- Surface without structure,
compatible with handrim Ø 19 mm or
Ø 20 mm.



Ultra-Grip „very grippy“

- ☐ 24" handrim for
24" wheel and 25" wheel € 402.40
A9000205070-004
- ☐ 26" handrim for
26" wheel € 402.40
A9000205070-006



Ergo-Grip

- ☐ 24" handrim for
24" wheel and 25" wheel € 402.40
A9000205080-004
- ☐ 26" handrim for
26" wheel € 402.40
A9000205080-006
- Attachment fitting to the rim

Note: Due to the complex geometry,
concentricity deviations of +/- 2-3 mm are
possible.



Curve black anodised, ergonomic handrim profile

- ☐ 24" handrim for
24" wheel and 25" wheel
Profile height 24 mm € 402.40
A9000205130-004
- ☐ 24" handrim for
24" wheel and 25" wheel
Profile height 31 mm € 402.40
Size L
A9000205140-004
- ☐ 26" handrim for
26" wheel
Profile height 31 mm € 402.40
Size L
A9000205140-006



Curve silver anodised ergonomic handrim profile ☐ 24" handrim for 24" wheel and 25" wheel Profile height 31 mm € 402.40 Size L A9000205160-004 	Curve powder-coated, ergonomic handrim profile ☐ = frame colour (gloss surface) ☐ = black (gloss surface) ☐ 24" handrim for 24" wheel and 25" wheel Profile height 24 mm € 467.60 A9000205135-004 ☐ 24" handrim for 24" wheel and 25" wheel Profile height 31 mm € 467.60 Size L A9000205145-004 ☐ 26" handrim for 26" wheel Profile height 31 mm € 467.60 Size L A9000205145-006 	Curve Grip, with Tetra-Grip coating, ergonomic handrim profile ☐ 24" handrim for 24" wheel and 25" wheel Profile height 31 mm € 520.70 Size L A9000205150-004 
Quadro with Tetra-Grip coating ☐ 24" handrim for 24" wheel and 25" wheel € 520.70 A9000205170-004 ☐ 25" handrim for 25" wheel € 520.70 A9000205170-005 ☐ 26" handrim for 26" wheel € 520.70 A9000205170-006 	Gekko black anodised ☐ 24" handrim for 24" wheel and 25" wheel € 520.70 A9000205175-004 ☐ 25" handrim for 25" wheel € 520.70 A9000205175-005 ☐ 26" handrim for 26" wheel € 520.70 A9000205175-006 with a robust and replaceable silicone lip, fully integrated into the profile Profile height 31 mm 	Gekko silver anodised ☐ 24" handrim for 24" wheel and 25" wheel € 520.70 A9000205185-004 with a robust and replaceable silicone lip, fully integrated into the profile Profile height 31 mm 
Gekko S black anodised ☐ 24" handrim for 24" wheel and 25" wheel € 520.70 A9000205195-004 with a robust and replaceable silicone lip, fully integrated into the profile Profile height 24 mm 	Equipment without handrims ☒ A9000205990	

Spoke guard

Finger- and spoke guard „transparent“

☐ for 24" wheel
A9000206300-004 € 207.10

☐ for 26" wheel
A9000206300-006 € 207.10

not in combination with Spinergy drive wheels.



Tyres

If not specified, the „Schwalbe RightRun in black“ tyres will be delivered.

Schwalbe „RightRun“ High-pressure tyres

☐ 24x1" (540-25)
A9000204001-004

Colour: Black
7.5 bar



Schwalbe „RightRun“ High-pressure tyres

☐ 24x1" (540-25)
A9000204000-004

☐ 25x1" (559-25)
A9000204000-005

☐ 26x1" (590-25)
A9000204000-006

Colour: Grey, black wall
7.5 bar



Schwalbe „DownTown“ Profiled tyres

☐ 24x1" (540-25)
A9000204010-004

7.5 bar



Schwalbe „Marathon Plus Evolution“ with „smart guard“ puncture protection

☐ 24x1" (540-25)
A9000204030-014 € 49.20

☐ 25x1" (559-25)
A9000204030-015 € 49.20

☐ 26x1" (590-25)
A9000204030-016 € 49.20

Colour: Black
10 bar



Schwalbe „Marathon Plus“ with „smart guard“ puncture protection

☐ 24x1" (540-25)
A9000204030-004 € 49.20


Colour: Grey, black wall
10 bar



Schwalbe „One“

☐ 24x1" (540-25)
A9000204123-004 € 93.00

Colour: Black
9 bar



Puncture-proof „solid rubber tyres“

☐ 24x1" (540-25)
A9000204201-004 € 122.70

Colour: Black
not in combination with Spinergy drive wheels.



Accessories

24" MTB spare wheel set

with double rim (ETRTO 507), aluminum front hub, quick-release axle and stainless steel handrim
Mountain bike tyres, selectively:

☐ A9000203410 € 703.40

24x1.75" Schwalbe „Marathon Plus"
Brake setting compatible with 24x1" wheel



☐ A9000203413 € 703.40

24x2.1" Schwalbe „SMART SAM Sport"
Brake setting not compatible with 24x1" wheel




NOTE: We charge a surcharge for the use of drum brakes

☐ A9000200004 € 232.20

Through axle with assistance, Large push button

☐ A9000200096 € 105.80

Interlock dimension: 122 / 110 mm
Axle diameter: 12 mm
Not in combination with one-hand drive



Hand air compressor with pressure gauge

☐ A9000901100 € 106.30

up to 14 bar



Crutch holder with flexible nylon cup

☐ A9000902310* € 137.40



Length-adjustable safety belt with buckle

☐ A9000902100 € 80.60



Adjustable safety belt with hook and fastener

☐ A9000902110 € 64.50



SIGMA lighting set

☐ A9000903270 € 107.90

Consisting of front light and rear lights,
chargeable via Micro USB socket,
incl. Micro USB charging cable
Lighting set compliant with road traffic
regulations



Transport bag for drive wheels

☐ for 24" wheel
A9000903530-004 € 77.80
Ø 76 cm x 2 cm

☐ for 26" wheel
A9000903530-006 € 77.80
Ø 83 cm x 2 cm



Separate order sheets

Category „Rehabilitation technology accessories"

e.g. Seat & backrest systems, one-hand
drive or transport restraint systems
can be found at:

Additional products

such as hand bikes or traction &
attachment devices can be found at:
www.proactiv-gmbh.com/downloads/

Additional sheets

such as heavy-duty versions etc.
can be found at:
www.proactiv-gmbh.com/downloads/

Services

Consultation, clarification of the
installation situation and measurement
through PRO ACTIV

☐ A9990000000 € 354.50

Handover and instruction on the
wheelchair by PRO ACTIV

☐ A9990000006 € 354.50

Notes and remarks:

Location, Date

The signature field is intended for the consultant and insured person. This confirms that the intended equipment is intended for the insured person.
- 17 - 2026-03-03 / The prices are plus statutory value-added tax and shipping costs ☐ = included in the basic price ☒ = Position selected

Signature End customer/dealer