

Im Hofstätt 11
 D-72359 Dotternhausen
 Tel +49 (0)7427 9480-0
 Fax +49 (0)7427 9480-7025
 e-mail info@proactiv-gmbh.de
 www.proactiv-gmbh.de



LITTY 4all Wheelchair for children

Maximum user weight limit 50 kg
 Weight from 6.7 kg

A9043100000-* € 4,060.00

Quotation date: _____
Consultation was carried out by
Staff at the medical supplies shop: _____
Enquiry no.: _____

Order date: _____
 according to GTC PRO ACTIV Reha-Technik GmbH
Order no.: _____

Medical suppliers / company address:

Cust. no.: _____
 Name: _____
 Street: _____
 Postcode / city: _____
 Phone: _____
 Fax: _____
 E-Mail: _____
 Contact: _____

Invoicing address: same as company address
 Other (see below)

Cust. no.: _____
 Name: _____
 Street: _____
 Postcode / city: _____
 Phone: _____
 Fax: _____
 E-Mail: _____
 Contact: _____


Shipping address: same as company address
 Other (see below)

Name: _____
 Street: _____
 Postcode / city: _____

Customer name / commission:

ID / Comm. No.: _____
 Surname, forename: _____
 Street: _____
 Postcode / city: _____
 Phone: _____
 Fax: _____
 E-Mail: _____
 D.o.B.: _____
 Paying organisation: _____

Body size:



Body height: _____ cm
 Weight: _____ kg
 Lower leg length: _____ cm
 (anatomic dimension: hollow of the knee to the foot sole)
 Indication: _____


This information is necessary for checking. If the information is missing, a check of the relevant items cannot be made!

Basic size

Seat plate and back shell screwed to the side part

Please note: Please select the seat plate and back shell separately.

Wheelchair dimensions:



| | | | | | | | | | |
|----------------------|-----------------------------|-----------------------------|-----------------------------|-----------------------------|-----------------------------|-----------------------------|-----------------------------|-----------------------------|-----------------------------|
| Seat width SB | <input type="checkbox"/> 18 | <input type="checkbox"/> 20 | <input type="checkbox"/> 22 | <input type="checkbox"/> 24 | <input type="checkbox"/> 26 | <input type="checkbox"/> 28 | <input type="checkbox"/> 30 | <input type="checkbox"/> 32 | <input type="checkbox"/> 34 |
|----------------------|-----------------------------|-----------------------------|-----------------------------|-----------------------------|-----------------------------|-----------------------------|-----------------------------|-----------------------------|-----------------------------|

| | | | | |
|----------------------|-----------------------------|-----------------------------|-----------------------------|-----------------------------|
| Seat depth ST | <input type="checkbox"/> 20 | <input type="checkbox"/> 25 | <input type="checkbox"/> 30 | <input type="checkbox"/> 35 |
|----------------------|-----------------------------|-----------------------------|-----------------------------|-----------------------------|

Note: height adjustment via perforated grid +/- 2 cm

| | | | | |
|-----------------------|-----------------------------|-----------------------------|-----------------------------|-----------------------------|
| Back height RH | <input type="checkbox"/> 20 | <input type="checkbox"/> 25 | <input type="checkbox"/> 30 | <input type="checkbox"/> 35 |
|-----------------------|-----------------------------|-----------------------------|-----------------------------|-----------------------------|

Note: height adjustment via perforated grid +/- 2 cm

| | | | | | | | | | | |
|-----------------------------|-----------------------------|-----------------------------|-----------------------------|-----------------------------|-----------------------------|-----------------------------|-----------------------------|-----------------------------|-----------------------------|-----------------------------|
| Lower leg length USL | <input type="checkbox"/> 20 | <input type="checkbox"/> 21 | <input type="checkbox"/> 22 | <input type="checkbox"/> 23 | <input type="checkbox"/> 24 | <input type="checkbox"/> 25 | <input type="checkbox"/> 26 | <input type="checkbox"/> 27 | <input type="checkbox"/> 28 | <input type="checkbox"/> 29 |
| | <input type="checkbox"/> 30 | <input type="checkbox"/> 31 | <input type="checkbox"/> 32 | <input type="checkbox"/> 33 | <input type="checkbox"/> 34 | <input type="checkbox"/> 35 | <input type="checkbox"/> 36 | <input type="checkbox"/> 37 | <input type="checkbox"/> 38 | <input type="checkbox"/> 39 |
| | <input type="checkbox"/> 40 | <input type="checkbox"/> 41 | <input type="checkbox"/> 42 | <input type="checkbox"/> 43 | <input type="checkbox"/> 44 | | | | | |

| | | | |
|---|------------------------------|------------------------------|------------------------------|
| Selecting the drive wheel size (required for further dimensions) | <input type="checkbox"/> 20" | <input type="checkbox"/> 22" | <input type="checkbox"/> 24" |
|---|------------------------------|------------------------------|------------------------------|

| | | | | | | |
|-----------------------------|------|-----------------------------|-----------------------------|-----------------------------|-----------------------------|-----------------------------|
| Seat height back SHh | 20": | <input type="checkbox"/> 31 | <input type="checkbox"/> 32 | <input type="checkbox"/> 33 | <input type="checkbox"/> 34 | <input type="checkbox"/> 35 |
| | 22": | <input type="checkbox"/> 34 | <input type="checkbox"/> 35 | <input type="checkbox"/> 36 | <input type="checkbox"/> 37 | <input type="checkbox"/> 38 |
| | 24": | <input type="checkbox"/> 37 | <input type="checkbox"/> 38 | <input type="checkbox"/> 39 | <input type="checkbox"/> 40 | <input type="checkbox"/> 41 |

| | | | | | | |
|------------------------------|------|-----------------------------|-----------------------------|-----------------------------|-----------------------------|-----------------------------|
| Seat height front SHv | 20": | <input type="checkbox"/> 34 | <input type="checkbox"/> 35 | <input type="checkbox"/> 36 | <input type="checkbox"/> 37 | <input type="checkbox"/> 38 |
| | 22": | <input type="checkbox"/> 37 | <input type="checkbox"/> 38 | <input type="checkbox"/> 39 | <input type="checkbox"/> 40 | <input type="checkbox"/> 41 |
| | 24": | <input type="checkbox"/> 40 | <input type="checkbox"/> 41 | <input type="checkbox"/> 42 | <input type="checkbox"/> 43 | <input type="checkbox"/> 44 |

| | | | | |
|--|------|----------------------------|----------------------------|----------------------------|
| Tipping point KP (inner edge of rear tube to centre of drive wheel axle) | 20": | <input type="checkbox"/> 5 | <input type="checkbox"/> 6 | <input type="checkbox"/> 7 |
| | 22": | <input type="checkbox"/> 6 | <input type="checkbox"/> 7 | <input type="checkbox"/> 8 |
| | 24": | <input type="checkbox"/> 7 | <input type="checkbox"/> 8 | <input type="checkbox"/> 9 |

| | | | |
|--|------|------|------|
| H: Position of the steering wheel bearing block (only for information) | 20": | 22": | 24": |
| | 37 | 40 | 43 |

I: Distance front edge footrest to back _____ cm (If not specified: Seat depth + 15 cm)

Different dimensions are available for an additional charge.

The dimensional tolerances according to DIN EN ISO 13920 apply to the design (dimensions in mm)

| Accuracy level | up to 30 | over 30 to 120 | over 120 to 400 | over 400 to 1000 |
|----------------|--------------|----------------|-----------------|------------------|
| C | +/- 1 | +/- 3 | +/- 4 | +/- 6 |

Frame colour - One-colour design

Frame, side parts and seat system

Information on current colour designs and example images can be found at: <https://www.proactiv-gmbh.com/colours/>

Frame colour: Selection from colour palette

Frame colour will be specified in the order

A9000100170-000

Raw, blasted, not coated

A9000100170-020

Note: In combination with moisture e.g. hand perspiration, can affect dark aluminum discoloration.

Powder colours

white RAL 9016

gloss A9000100170-023

matt A9000100170-123 € 43.40

brilliant silver

gloss A9000100170-011

anthracite RAL 7016

gloss A9000100170-001

matt A9000100170-101 € 43.40

black RAL 9005

gloss A9000100170-003

matt A9000100170-019 € 43.40

pure orange RAL 2004

gloss A9000100170-022

matt A9000100170-122 € 43.40

pink magenta RAL 4010

gloss A9000100170-037

matt A9000100170-137 € 43.40

pastel blue RAL 5024

gloss A9000100170-036

matt A9000100170-136 € 43.40

lemon green

gloss A9000100170-005

matt A9000100170-105 € 43.40

yellow RAL 1023

gloss A9000100170-046

matt A9000100170-146 € 43.40

red RAL 3001

gloss A9000100170-047

matt A9000100170-147 € 43.40

lagoon blue

gloss A9000100170-052

matt A9000100170-152 € 43.40

athletic grey RAL 7046

gloss A9000100170-050

matt A9000100170-150 € 43.40

muddy grey RAL 7033

gloss A9000100170-051

matt A9000100170-151 € 43.40

peach RAL 3022

gloss A9000100170-038

matt A9000100170-138 € 43.40

Powder colours „dormant“

deep black effect silver

gloss A9000100170-007

sky blue

gloss A9000100170-017

navy blue

gloss A9000100170-045

matt A9000100170-145 € 43.40

electric green

gloss A9000100170-015

matt A9000100170-115 € 43.40

violet

gloss A9000100170-044

fire red

gloss A9000100170-018

matt A9000100170-118 € 43.40

Frame colour - One-colour design special colours „liquid coating“

Frame, side parts and seat system

€ 349.60

storm (ice blue)

gloss A9000100170-500

matt A9000100170-501

prairie (rose gold)

gloss A9000100170-504

matt A9000100170-505

berry

gloss A9000100170-526

matt A9000100170-527

blue ocean

gloss A9000100170-522

matt A9000100170-523

emerald

gloss A9000100170-514

matt A9000100170-515

chili

gloss A9000100170-516

matt A9000100170-517

Frame colour - Two-colour design

Note: In the two-colour design, both colours are the same gloss level as the design colour.

- Two-colour design consisting of one frame colour and a second colour for the side parts resp. seat plate and back shell** € 268.50
A9042100170-200

Please select colours below:

Two-colour design LITTY consisting of a base colour at the frame with axle tube and a second colour for the side parts and if desired, for the seat system. Select from one-colour design:

A: Frame colour: A9000100170-

B: Side part colour: A9000100170-

Colour seat plate and back shell: A or B

Colour - Attachment parts

Note: If not specified, the caster fork and bearing block will be mounted in black.

Drive wheel: aluminum rim, hub and spokes

- Black A9000100171-003

Caster fork and bearing block

- Silver A9000100173-011
 Black A9000100173-003

Colour design of attachment parts according to customer request:

Drive wheel: aluminum rim and design inlet according to colour selection / hub and spokes black

- frame colour A9000100171-999 € 249.60
 Selection from one-colour design (different from the frame colour) A9000100171-998 € 288.90
Colour: _____
A9000100170-

Caster wheel: 5" aluminum rim¹
(6-spoke design)

- frame colour A9000100174-999 € 75.20
 Selection from one-colour design (different from the frame colour) A9000100174-998 € 87.10
Colour: _____
A9000100170-

Special colour request (Frame colour - One-colour design)

Frame colour for single purchasing Selection from RAL colour palette

- A9000100170-100 € 414.40

gloss level

- gloss
 matt

colour designation

RAL number

Note: Attachments **not** available in design colour options

¹ Only in combination with caster wheel 5" aluminum rim 6-spoke-design good shock absorption A9000501210 or A9000501211

Push handles

Equipment without push handles

A9000401000



Central push handle aluminum

A9042401085-K
height-adjustable
Holder welded to the center of the back crossbar



Safety push handles back-positioned

Length: 45 cm
max. height = SHh + approx. 53 cm
A9042401075-K001

Length: 55 cm
max. height = SHh + approx. 63 cm
A9042401075-K002
continuous height-adjustable
push handle holder back-positioned



Aluminium central push handle for brake handles Foldable handle area

A9042401086-K € 141.90

height-adjustable and foldable in the handle area via toothed joint, holder welded on the middle on the back cross bar.
Approved up to a body weight of 50 kg
Only in combination with A9042701020 Parking knee lever brake for accompanying person



Aluminium central push handle Foldable handle area

A9042401084-K € 102.60

height-adjustable and foldable in the handle area via toothed joint, holder welded on the middle on the back cross bar.
Approved up to a body weight of 50 kg



Aluminium push handle closed, foldable via toothed joint

A9042401069-K € 389.90

ontinuous height-adjustable
push handles shifted to the back



Back equipment

If not specified, the „ergonomic back shell" will be delivered.

Ergonomic back shell aluminum

A9042802210-K € 441.60

in frame colour
from RH 20 cm and SB 18 cm



Back cushion for ergonomic back shell

A9042802220-K € 221.30

Cover:

- Black, breathable, open surface
- Black, breathable, closed surface

If not specified, the cover material „Black, breathable, open surface" will be delivered.



Equipment without back shell and back cushion

A9042802011



Seat system

Seat plate aluminum

A9042801210-K
in frame colour
bolted between the aluminum side parts
adjustable in height, angle and depth



Seat cushion



Seat cushion in rectangular contour Inlet: Standard

A9000805180-K* € 118.30

Seat cushion with radius contour Inlet: Standard

A9000800300-K* € 163.10

Seat cushion height:

3 cm 5 cm
 4 cm 6 cm

Seat cushion in rectangular contour Inlet: Anti-Decubitus

up to 30x30 cm
A9000805204-K* € 301.70

>30x30 up to 50x50 cm
A9000805200-K* € 437.70

Seat cushion with radius contour Inlet: Anti-Decubitus

up to 30x30 cm
A9000800311-K* € 356.10

>30x30 up to 50x50 cm
A9000800310-K* € 482.90

Seat cushion height:

4 cm 6 cm

Cover:

- Black, breathable, open surface
 Black, breathable, closed surface

If not specified, the cover material „Black, breathable, open surface“ will be delivered.

Seat cushion is provided

Wheel camber

11° drive wheel camber

A9042202030

Brake system

If not specified, the „Knee lever brake“ will be delivered.

Knee lever brake

A9042701010
Brake lever integrated in side part



Knee lever brake with position lock for brake lever

A9042701030 € 304.30
integrated into the side part when braked and unbraked.
Operated by an accompanying person.



Parking knee lever brake for accompanying person

A9042701020 € 464.10
Push handle operation via Bowden cable.
Only in combination with knee lever brake and
push handles



Footrest

If not specified, the „fixed footrest“ will be delivered.

Fixed footrest

A9042601010-K
adjustable in height, depth and angle
including footplate cover with anti-slip film



Footrest folds up to the rear

A9042602010-K € 286.80
including footplate cover with anti-slip film



Adaptations

Anti-tipping support - tiltable upwards

- right side
A9042301160-12* € 207.10
- left side
A9042301160-22* € 207.10
removable via Quick-pin



Outdoor wheel attachment with steering fork

- A9042900200-K € 845.70

Equipped with straight-running lock can be removed from the wheelchair via quick release lever with swivelling pneumatic tyre 12.1/2 x 2.1/4" caster wheel can be removed from the support tube via quick-release axle for maximum ground pressure at the steering attachment of 30 kg



3" caster wheel

3" skate wheel (80 mm)

- A9000501110 € 70.80
puncture-proof



4" caster wheel

4" caster wheel „Standard“ solid rubber black

- A9000501011-004 (82A), puncture-proof



4" skate wheel (100 mm)

- A9000501120 € 70.80
puncture-proof



4" safety light-up wheel (100 mm)

- A9000501220 € 94.80
puncture-proof



4" CorMORun caster wheel (100 mm)

- A9000501042 € 173.00
with replaceable rubber 6 mm, puncture-proof
Rim colour: Black
Incl. hair guard spacer (5.3 mm)



4" Frog Legs caster wheel (100 mm)

- A9000501230 € 206.30
with replaceable rubber 6 mm, puncture-proof
Colour: black, silver rim
Incl. hair guard spacer (5.3 mm)



5" caster wheel

5" caster wheel „Standard“ PU black

- A9000501036-005 (82A), puncture-proof



5" caster wheel „Standard“ PU grey

- A9000501031-005 (82A), puncture-proof



5" safety light-up wheel (125 mm)

- A9000501221 € 94.80
puncture-proof



5" aluminum rim 3-spoke design „good shock absorption“

silver anodised
A9000501210 € 179.80

black
A9000501211 € 179.80

with replaceable rubber 6 mm
medium hard, puncture-proof
Incl. hair guard spacer (3 mm)



5" CorMORun caster wheel (125 mm)

A9000501041 € 181.30

with replaceable rubber 6 mm,
puncture-proof
Rim colour: Black
Incl. hair guard spacer (5.3 mm)



5" Frog Legs caster wheel (125 mm)

A9000501232 € 206.30

with replaceable rubber 6 mm,
puncture-proof
Colour: black, silver rim
Incl. hair guard spacer (5.3 mm)



Caster fork

For caster fork colours, see „Colour of add-on parts“

If not specified, the „Standard caster fork with screw axle“ will be delivered.

Caster fork Standard

A9000500100-K001
with screwed axle

A9000500100-K002 € 167.40
with quick-release axle



Lenkradgabel-Wechselsatz LITTY

A9000500300-K € 296.70

with quick-release axle

The article is delivered without caster wheels.
If caster wheels are desired, please tick a
second option under the section „Caster
wheels.“



Drive wheels

Drive wheel aluminum rim, 18 spokes with Liteline hub

radially spoked

Design inlet in frame colour

Please select the drive wheel size according to page 2.

20" (ETRTO 451)
A9000203053


22" (ETRTO 489)
A9000203063

24" (ETRTO 540)
A9000203123



Handrims

x = dimension handrim attachment



Handrim attachment general:

- x = 11 mm
- x = 19 mm

Handrim attachment Curve + Gekko handrim:

- x = 10 mm
- x = 18 mm

If not specified, the handrim is mounted tightly (11 mm / 10 mm).

Stainless steel Ø 20 mm

- 20" handrim for 20" wheel
A9000205010-001
- 22" handrim for 22" wheel
A9000205010-002
- 24" handrim for 24" wheel
A9000205010-004



Aluminum Ø 19 mm silver anodised

- 22" handrim for 22" wheel
A9000205019-002
- 24" handrim for 24" wheel
A9000205019-004



Aluminum Ø 16 mm powder coated in frame colour

(basic colour with two-colour design)

- = Gloss surface
- = Matt surface

Note: Matt coating is not recommended due to the low grip.

- 20" handrim for 20" wheel
A9000205021-001 **€ 68.30**
- 22" handrim for 22" wheel
A9000205021-002 **€ 68.30**



Aluminum Ø 19 mm powder coated in frame colour

(basic colour with two-colour design)

- = Gloss surface
- = Matt surface

Note: Matt coating is not recommended due to the low grip.

- 20" handrim for 20" wheel
A9000205020-001 **€ 68.30**
- 22" handrim for 22" wheel
A9000205020-002 **€ 68.30**
- 24" handrim for 24" wheel
A9000205020-004 **€ 68.30**



Aluminum Ø 16 mm polished and hard-anodised

- 20" handrim for 20" wheel
A9000205031-001 **€ 164.90**
- 22" handrim for 22" wheel
A9000205031-002 **€ 164.90**



Aluminum Ø 19 mm polished and hard-anodised

- 20" handrim for 20" wheel
A9000205030-001 **€ 164.90**
- 22" handrim for 22" wheel
A9000205030-002 **€ 164.90**
- 24" handrim for 24" wheel
A9000205030-004 **€ 164.90**



Titanium Ø 19 mm

- 24" handrim for 24" wheel
A9000205040-004 **€ 335.60**



Silicone handrim cover

- for handrim 20" A9000205100-001 **€ 232.80**
 - for handrim 22" A9000205100-002 **€ 232.80**
 - for handrim 24" A9000205100-004 **€ 232.80**
- Surface without structure, compatible with handrim Ø 19 mm or Ø 20 mm.



Ultra-Grip „very grippy“

- 20" handrim for 20" wheel
A9000205070-001 **€ 402.40**
- 22" handrim for 22" wheel
A9000205070-002 **€ 402.40**
- 24" handrim for 24" wheel
A9000205070-004 **€ 402.40**



| | | | |
|---|--|--|--|
| <p>Ergo-Grip</p> <p><input type="checkbox"/> 20" handrim for 20" wheel A9000205080-001 € 402.40</p> <p><input type="checkbox"/> 22" handrim for 22" wheel A9000205080-002 € 402.40</p> <p><input type="checkbox"/> 24" handrim for 24" wheel A9000205080-004 € 402.40</p> <p>Attachment fitting to the rim</p> <p>Note: Due to the complex geometry, concentricity deviations of +/- 2-3 mm are possible.</p>  | <p>Maxgrepp with cellular and natural rubber coating very grippy, soft</p> <p><input type="checkbox"/> 20" handrim for 20" wheel A9000200788 € 402.40</p> <p><input type="checkbox"/> 22" handrim for 22" wheel A9000200789 € 402.40</p> <p><input type="checkbox"/> 24" handrim for 24" wheel A9000200787 € 402.40</p>  | <p>Curve black anodised, ergonomic handrim profile</p> <p><input type="checkbox"/> 22" handrim for 22" wheel Profile height 24 mm A9000205130-002 € 402.40</p> <p><input type="checkbox"/> 24" handrim for 24" wheel Profile height 24 mm A9000205130-004 € 402.40</p>  | |
| <p>Curve powder-coated, ergonomic handrim profile</p> <p><input type="checkbox"/> = frame colour (gloss surface)</p> <p><input type="checkbox"/> = black (gloss surface)</p> <p><input type="checkbox"/> 24" handrim for 24" wheel Profile height 24 mm A9000205135-004 € 467.60</p>  | <p>Curve Grip, with Tetra-Grip coating, ergonomic handrim profile</p> <p><input type="checkbox"/> 24" handrim for 24" wheel Profile height 24 mm A9000205155-004 € 520.70</p>  | <p>Quadro with Tetra-Grip coating</p> <p><input type="checkbox"/> 24" handrim for 24" wheel A9000205170-004 € 520.70</p>  | |
| <p>Gekko black anodised</p> <p><input type="checkbox"/> 22" handrim for 22" wheel A9000205175-002 € 520.70</p> <p><input type="checkbox"/> 24" handrim for 24" wheel A9000205175-004 € 520.70</p> <p>with a robust and replaceable silicone lip, fully integrated into the profile Profile height 31 mm</p>  | <p>Gekko silver anodised</p> <p><input type="checkbox"/> 24" handrim for 24" wheel A9000205185-004 € 520.70</p> <p>with a robust and replaceable silicone lip, fully integrated into the profile Profile height 31 mm</p>  | <p>Gekko S black anodised</p> <p><input type="checkbox"/> 22" handrim for 22" wheel A9000205195-002 € 520.70</p> <p><input type="checkbox"/> 24" handrim for 24" wheel and 25" wheel A9000205195-004 € 520.70</p> <p>with a robust and replaceable silicone lip, fully integrated into the profile Profile height 24 mm</p>  | |
| <p>Equipment without handrims</p> <p><input checked="" type="checkbox"/> A9000205990</p> | | | |
| <p>Spoke guard</p> | | | |
| <p>Finger- and spoke guard „transparent“</p> <p><input type="checkbox"/> for 20" wheel A9000206300-001 € 207.10</p> <p><input type="checkbox"/> for 22" wheel A9000206300-002 € 207.10</p> | | | <p><input type="checkbox"/> for 24" wheel A9000206300-004 € 207.10</p>  |

Tyres

If not specified, the „Schwalbe RightRun in black” tyres will be delivered.

Schwalbe „RightRun“ High-pressure tyres

22x1" (489-25)
A9000204001-002

24x1" (540-25)
A9000204001-004

Colour: Black
7.5 bar



Schwalbe „RightRun“ High-pressure tyres

20x1" (451-25)
A9000204000-001

22x1" (489-25)
A9000204000-002

24x1" (540-25)
A9000204000-004
Colour: Grey, black wall
7.5 bar



Schwalbe „Marathon Plus Evolution“ with „smart guard“ puncture protection and reflector strips

22x1" (489-25)
A9000204030-002

€ 49.20

Colour: Black
10 bar



Schwalbe „Marathon Plus Evolution“ with „smart guard“ puncture protection

24x1" (540-25)
A9000204030-014

€ 49.20

Colour: Black
10 bar



Schwalbe „MarathonPlus“ Schwalbe „Marathon Plus“ with „smart guard“ puncture protection

24x1" (540-25)
A9000204030-004

€ 49.20

Colour: Grey, black wall
10 bar




Accessories

Through axle with assistance, Large push button

A9000200096 € 105.80

Interlock dimension: 122 / 110 mm
Axle diameter: 12 mm



Hand air compressor with pressure gauge

A9000901100 € 106.30

up to 14 bar



Length-adjustable safety belt with buckle

A9000902100 € 80.60



Adjustable safety belt with hook and fastener

A9000902110 € 64.50



Transport bag for drive wheels

for 24" wheel
A9000903530-004 € 77.80
Ø 76 cm x 2 cm



Therapy table

width 350 mm
A9000902800-001 € 480.20
removable table plate from Lexan, transparent,
with front edging, with 2 positioning stops
Ø 25 mm (for adaptation to frame tube
diameter Ø 25 mm)



Services

Consultation, clarification of the installation situation and measurement through PRO ACTIV

A9990000000 € 354.50

Handover and instruction on the wheelchair by PRO ACTIV

A9990000006 € 354.50

Notes and remarks:

Location, Date

Signature End customer/dealer

The signature field is intended for the consultant and insured person. This confirms that the intended equipment is intended for the insured person.
- 12 - 2026-03-03 / The prices are plus statutory value-added tax and shipping costs = included in the basic price = Position selected

Order sheet LITTY 4all Rev. 04 UK