



**Quotation date:** \_\_\_\_\_  
**Consultation was carried out by**  
**Staff at the medical supplies shop:** \_\_\_\_\_

**Enquiry no.:** \_\_\_\_\_

**Order date:** \_\_\_\_\_  
according to GTC PRO ACTIV Reha-Technik GmbH

**Order no.:** \_\_\_\_\_

**Medical suppliers / company address:**

**Cust. no.:** \_\_\_\_\_  
**Name:** \_\_\_\_\_  
**Street:** \_\_\_\_\_  
**Postcode / city:** \_\_\_\_\_  
**Phone:** \_\_\_\_\_  
**Fax:** \_\_\_\_\_  
**E-Mail:** \_\_\_\_\_  
**Contact:** \_\_\_\_\_

**Different invoicing address (see below)**

**Cust. no.:** \_\_\_\_\_  
**Name:** \_\_\_\_\_  
**Street:** \_\_\_\_\_  
**Postcode / city:** \_\_\_\_\_

**Different shipping address (see below)**

**Cust. no.:** \_\_\_\_\_  
**Name:** \_\_\_\_\_  
**Street:** \_\_\_\_\_  
**Postcode / city:** \_\_\_\_\_

**Customer name / commission:**

**ID / Comm. No.:** \_\_\_\_\_  
**Surname, forename:** \_\_\_\_\_  
**Street:** \_\_\_\_\_  
**Postcode / city:** \_\_\_\_\_  
**Phone:** \_\_\_\_\_  
**E-Mail:** \_\_\_\_\_  
**D.o.B.:** \_\_\_\_\_  
**Paying organisation:** \_\_\_\_\_  
**Indication:** \_\_\_\_\_

## NJ1 el. compact bike

approved up to 100 kg payload

- girder chassis with integrated rechargeable battery holder
- drive unit, can be decoupled via the eccentric quick-lock system
- A-shaped handles with integrated shift and brake controls
- Backrest continuously angle-adjustable
- 6° wheel camber
- Hydraulic disc brake system
- additional mechanical disc brake with parking lock
- Safety flag with holder
- Leg protection guard
- Drive unit black

### Neodrives Z15 drive unit

„The powerful hub drive for maximum performance  
with TFT display and 11-speed Shimano derailleur gears”  
torque 12/40 Nm, power 250 watt

Wheel size: 26"

A9100100011-026 € 8,780.00

>> Expected to be available from August 2026 <<

### Body size

Body height: \_\_\_\_\_ cm  
Weight: \_\_\_\_\_ kg

## Chassis Colour

Information on current colour designs and example images can be found at: <https://www.proactiv-gmbh.com/colours/>

### white RAL 9016

- gloss A9000100170-023  
 matt A9000100170-123 € 43.40

### brilliant silver

- gloss A9000100170-011

### anthracite RAL 7016

- gloss A9000100170-001  
 matt A9000100170-101 € 43.40

### black RAL 9005

- gloss A9000100170-003  
 matt A9000100170-019 € 43.40

### deep black effect silver

- gloss A9000100170-007

### pure orange RAL 2004

- gloss A9000100170-022

### fire red

- gloss A9000100170-018  
 matt A9000100170-118 € 43.40

### pink magenta RAL 4010

- gloss A9000100170-037  
 matt A9000100170-137 € 43.40

### sky blue

- gloss A9000100170-017

### pastel blue RAL 5024

- gloss A9000100170-036  
 matt A9000100170-136 € 43.40

### navy blue

- gloss A9000100170-045  
 matt A9000100170-145 € 43.40

### electric green

- gloss A9000100170-015  
 matt A9000100170-115 € 43.40

### lemon green

- gloss A9000100170-005  
 matt A9000100170-105 € 43.40

### violet

- gloss A9000100170-044

### blasted, transparent powder coated

- gloss A9000100170-029

### yellow RAL 1023

- gloss A9000100170-046  
 matt A9000100170-146 € 43.40

### red RAL 3001

- gloss A9000100170-047  
 matt A9000100170-147 € 43.40

### lagoon blue

- gloss A9000100170-052  
 matt A9000100170-152 € 43.40

### athletic grey RAL 7046

- gloss A9000100170-050  
 matt A9000100170-150 € 43.40

### muddy grey RAL 7033

- gloss A9000100170-051  
 matt A9000100170-151 € 43.40

### peach RAL 3022

- gloss A9000100170-038  
 matt A9000100170-138 € 43.40

## Battery Neodrives

The specified ranges are dependent on the route profile, user weight and speed.

### Battery neodrives DT

**36V 16Ah 576Wh**  
range approx. 75 km

- A9100900109



### Double battery neodrives DT

**36V 16Ah 576Wh**  
range approx. 150 km

- A9100900110 € 1,412.40

plug-in device for easy switching between batteries via the magnetic connector

## Derailleur gear

### Shimano Deore 1x11 Rapid Fire 11-speed derailleur

- A9110200032  
mechanical derailleur gear with operation on the handle. 11-speed cassette 11-42 teeth

#### Selection of the chainring:

- 32 teeth  
 38 teeth



### Shimano Deore XT 1x11 Rapid Fire 11-speed derailleur

- A9110200034 € 192.20  
mechanical derailleur gear with control on the handle. 11-speed cassette 11-42 teeth top chain ring with chain guards

#### Selection of the chainring:

- 32 teeth  
 38 teeth

### Shimano Deore XT 1x11 DI2 with push button operation 11-speed derailleur

- A9110200035 € 1,486.10  
electronic derailleur gear with operation on the handle. 11-speed cassette 11-42 teeth

Operation: right side  left side

#### Selection of the chainring:

- 32 teeth  
 38 teeth



**Shimano Deore XT 1x11 DI2  
with chin control  
11-speed derailleur gear system**

A9110200033      € 1,631.10

electronic derailleur with chin control  
(tetra equipment)  
11-speed cassette 11-42 teeth

Selection of the chainring:

32 teeth

38 teeth



## Gearshift accessories

**Mountaindrive 2,5 : 1  
Hill reduction  
in combination with hub gear**

A9100200003      € 724.90

Planetary gearbox for doubling the gears.  
Allows you to shift downhill while standing still.  
As standard with 38 teeth chainring.  
In combination with grip width 40 cm, 45 cm,  
48 cm or 53 cm



## Brake system

**Powerful braking system with brake  
disc Ø 203 mm**

A9110300105      € 83.10

recommended from 80 kg body weight.



**Double disc brake system „Magura”  
brake disc Ø 203 mm**

A9100300100      € 233.30

eliminates the mechanical disk brake.  
Only in combination with the powerful braking  
system.  
recommended from 100 kg body weight.

**Mechanical back-pedal disc brake  
with crank-free function**

A9110300004      € 738.40

compatible with derailleur and freewheel hub.  
In combination with:  
- Bottom bracket shaft reach 45 cm to 55 cm  
- pedal crank length 155 mm or 175 mm  
up to 120 kg body weight  
Not in combination with Mountaindrive 2.5 : 1  
hill reduction!

## Dimensions

### Crank type and grip width

#### Standard crank

- grip width 40 cm  
A9110500100-040
- grip width 45 cm  
A9110500100-045
- grip width 50 cm  
A9110500100-050

#### Tetra crank (flared)

with extended crank arms for maximum upper body space

- grip width 48 cm (176 mm)  
A9110500100-048
- grip width 53 cm (226 mm)  
A9110500100-053
- grip width 58 cm (276 mm)  
A9110500100-058

#### Standard pedal crank

- length 155 mm  
A9110600100-155
- length 175 mm  
A9110600100-175
- length 195 mm  
A9110600100-195



#### Tetra pedal crank


- length 155 mm  
A9100600105-155      € 136.10
- length 175 mm  
A9100600105-175      € 136.10
- individual length  
195 mm to 140 mm      € 258.70  
A9100600100



#### Fixed bottom bracket support

(Limited availability – only as long as stocks last. Alternative: Adjustable bottom bracket support)

- S (X/Y 67/153)  
A9100501100-003
- M (X/Y 90/250)  
A9100501100-001



#### Adjustable bottom bracket support

the angle and length can be adjusted via the clamped rotary point on the fork bridge and length.

For the pre-setting, please specify the bottom bracket support size that comes as close to your desired position.

- S (180 mm)  
A9100500001-002      € 232.10
- M (240 mm)  
A9100500001-003      € 232.10

If not specified, the adjustable bottom bracket support will be delivered in S.



## Legrest

#### Divided legrest

- A9100100008  
adaptation on the fork  
for short legs  
Only possible up to body size 165 cm



#### All-round legrest


- A9100100009  
adaptation on the fork



#### Foot stop for legrest

- A9100100012      € 91.20




<b>Upper thigh fixation</b> <input type="checkbox"/> A9100100013      € 215.10 	
---	--


## Back equipment

<b>Backrest with filling plate and back cushion</b> <input checked="" type="checkbox"/> A9100800002 Cover: Black, breathable, open surface	<b>Backrest upholstery Body Contour</b> <input type="checkbox"/> A9100800010      € 92.40 with back cushion 3 cm Cover: Black, breathable, open surface 
--	---

## Chassis


<b>Chassis geometry</b> <input checked="" type="checkbox"/> length 116 cm SB approx. 39 cm A9100100000-002 <input type="checkbox"/> length 116 cm SB approx. 43 cm      € 204.60 A9100100000-003	<input type="checkbox"/> Heavy Duty version up to 120 kg body weight A9100800013      € 425.60 (Recommended for off-road use) Required in addition to A9100100000*.	
--	---	--

## Seat upholstery

<b>Seat upholstery open belt system</b> Chassis: 116 cm SB 39 cm, ST 40 cm <input checked="" type="checkbox"/> A9100801010-002	Chassis: 116 cm SB 43 cm, ST 40 cm <input type="checkbox"/> A9100801010-003	
---	---	---

<b>Seat upholstery Body Contour</b> Chassis: 116 cm SB 39 cm, ST 40 cm <input type="checkbox"/> A9100801025-002	Chassis: 116 cm SB 43 cm, ST 40 cm <input type="checkbox"/> A9100801025-003	
--	---	---

## Seat cushion

<b>Seat cushion standard</b> Cover: Black, breathable, open surface Seat cushion height: 4 cm Chassis: 116 cm SB 39 cm, ST 40 cm <input type="checkbox"/> A9100000190-002      € 115.00	Chassis: 116 cm SB 43 cm, ST 40 cm <input type="checkbox"/> A9100000190-003      € 115.00	
--	---	---

<b>Seat cushion with upper thigh guide</b> Cover: Black, breathable, open surface Seat cushion height: 4 cm Only in combination with upper thigh fixation Chassis: 116 cm SB 39 cm, ST 40 cm <input type="checkbox"/> A9100800103-002      € 207.50	Chassis: 116 cm SB 43 cm, ST 40 cm <input type="checkbox"/> A9100800103-003      € 207.50	
---	---	---

## Tyres

### 26x1.5" Schwalbe „Marathon“

A9100105100  
(ETRTO 559-40)  
touring tyre



### 26x2" Schwalbe „Marathon Plus Tour“

A9100105140 € 71.10  
(ETRTO 559-50)  
touring tyre



### 26x2.1" Schwalbe „Smart Sam“

A9100105150 € 74.60  
(ETRTO 559-54)



## Accessories Bike

### Plug-in mudguard

A9100502100 € 58.10



### Rear guard clamped on the axle tube

A9100800100 € 218.10



### Neck rest

A9100800150 € 261.00  
height-adjustable, angle-adjustable



### Pelvic fixation belt

A9100800200 € 121.70



### Bracket Ø 22 mm x 42 mm for mounting accessories (tacho, lighting, etc.)

A9110903271 pcs. á € 23.20  
Colour: Natural anodised  
\_\_\_\_\_ piece



### SIGMA lighting set

A9120900019 € 123.90

consisting of front light and rear lights,  
chargeable via Micro USB socket,  
incl. Micro USB charging cable.  
Lighting set compliant with road traffic  
regulations



## Options for users with reduced functions

### Palm rest on the handle

A9110601200 € 102.80  
upholstered, bolted, angle-adjustable.  
Not in combination with other handles



### Tetra hand support, pivotable, adjustable grip distance

A9110601101 € 169.40



### Tetra backrest swivel and adjustable grip width

A9110601102 € 191.60

very stiff external support as round tube  
Ø 20 mm with foam rubber coating



### Tetra shell handles

A9110601100 € 172.40



### Operating satellite mounted on handle

Mounting: right side  left side

with existing support line (placed along a  
brake line)

A9116200000 € 103.10

without existing support line (an additional  
blind cable is required)

A9116200001 € 146.20

## Services

Consultation, clarification of the installation situation and measurement through PRO ACTIV

A9990000000 € 354.50

Handover and instruction of the handbike / traction device by PRO ACTIV

A9990000011 € 354.50

Notes and remarks:

### Location, Date

The signature field is intended for the consultant and insured person. This confirms that the intended equipment is intended for the insured person.

### Signature End customer/dealer

- 7 - 2026-03-03 / The prices are plus statutory value-added tax and shipping costs  = included in the basic price  = Position selected