

Customized production - model „SPEEDY 4all“

Frame design

Frame design Heavy Duty AF4-4 121 kg - 190 kg body weight

A9016100091 € 766.30

- Upper frame tube with inner tube Ø 20 x 2 mm
 - Inner frame tube Ø 20 x 2 mm on the lower frame tube in the area of the steering wheel bearing block
 - Frame reinforcement with 2 additional diagonal frame cross tubes and reinforcements
 - Reinforced drive wheel mount
 - Stainless steel backrest locking bolts
 - Back tube reinforcement inlay on the side part screw thread
 - Push handles in stainless steel design
 - Side part aluminium 3 mm
- Please note: Aluminium clothing guard removable via lock function (A9000403210-K) additionally required!
- Footrest or foot bar, continuous, reinforced design
 - Knee lever brake bolt, stainless steel

Note:

Please select the „Heavy Duty“ steering wheel fork under Steering wheel fork!

Seat width

Seat width 28 cm - 29 cm

A9016100040 € 213.30

SB 28 cm SB 29 cm

Services

Consulting and customization by PRO ACTIV according to individual customer requirements

A9990000017 € 612.00

Additional requirements