

Im Hofstätt 11
D-72359 Dotternhausen
Phone +49 (0)7427 9480-0
Fax +49 (0)7427 9480-7025
e-mail info@proactiv-gmbh.de
www.proactiv-gmbh.com



FREAK Compact bike

Maximum user weight limit 100 kg

- girder chassis with reinforcement brace
- Chassis length 102 cm, SB 35 cm
- Drive unit, can be decoupled via the eccentric quick-lock system
- Black drive unit with 20" drive wheel
- A-shaped handles with integrated switching fittings
- Shimano V-brake rim brake
- infinitely adjustable backrest
- wheel camber 6°, rear wheel size 20"
- safety flag

A9113000013

€ 4,550.00

Quotation date: _____

Consultation was carried out by
Staff at the medical supplies shop: _____

PRO ACTIV field representative: _____

Order date: _____

According to GTC PRO ACTIV Reha-Technik GmbH

Medical suppliers/company address:

Cust. no.: _____

Name: _____

Street: _____

Postcode / city: _____

Phone: _____

Fax: _____

E-Mail: _____

Contact: _____

Invoicing address, other(see below)

Cust. no.: _____

Name: _____

Street: _____

Postcode / city: _____

Shipping address, other (see below)

Cust. no.: _____

Name: _____

Street: _____

Postcode / city: _____

Customer name/commission:

Surname, forename: _____

Street: _____

Postcode / city: _____

Phone: _____

E-Mail: _____

D.O.B: _____

Paying organisation: _____

Indication: _____

Body size

Body height: _____ cm

Weight: _____ kg

Frame colour – one-colour design

Information about current colour designs and example images at: www.proactiv-gmbh.de/farben.php

White RAL 9016

- gloss A9000100170-023
 matt A9000100170-123 € 40.10

Brilliant silver

- gloss A9000100170-011

Anthracite RAL 7016

- gloss A9000100170-001
 matt A9000100170-101 € 40.10

Black RAL 9005

- gloss A9000100170-003
 matt A9000100170-019 € 40.10

Deep black, silver effect

- gloss A9000100170-007

Pure orange RAL 2004

- gloss A9000100170-022

Fire red

- gloss A9000100170-018
 matt A9000100170-118 € 40.10

Pink magenta RAL 4010

- gloss A9000100170-037
 matt A9000100170-137 € 40.10

Sky blue

- gloss A9000100170-017

Pastel blue RAL 5024

- gloss A9000100170-036
 matt A9000100170-136 € 40.10

Marine blue

- gloss A9000100170-045
 matt A9000100170-145 € 40.10

Electric green

- gloss A9000100170-015
 matt A9000100170-115 € 40.10

Lemon green

- gloss A9000100170-005
 matt A9000100170-105 € 40.10

Violet

- gloss A9000100170-044

Yellow RAL 1023

- gloss A9000100170-046
 matt A9000100170-146 € 40.10

Red RAL 3001

- gloss A9000100170-047
 matt A9000100170-147 € 40.10

Lagoon blue

- gloss A9000100170-052
 matt A9000100170-152 € 40.10

Athletic grey

- gloss A9000100170-050
 matt A9000100170-150 € 40.10

Muddy grey

- gloss A9000100170-051
 matt A9000100170-151 € 40.10

Neodrives

E-assistant version with Neodrives

- A9113200000 € 3,093.50

- Neodrives Z10 motor
 - torque 12/40 Nm, output 250 W
 - colour display sMMi Z15 Connect
 - pushing aid for travelling without pedal assistance up to 6 km/h
 - Neodrives RR G2 rechargeable battery with 14.25 Ah / 513 Wh
 - Range approx. 40-70 km
- Not in combination with hub gears!



Hub gear

Not in combination with e-assistant version

Shimano "Nexus" 8-speed hub gear

- A9113000015
 "Grip Shift" twist grip shifter
 Gear ratio: Sprocket 32 teeth, Cassette 20 teeth
 with rear pedal brake
 Shimano V-brake brake lever, left side

Shimano "Alfine" 11-speed hub gear

- A9113200026 € 793.50
 "Rapid Fire" handle
 Gear ratio: Sprocket 32 teeth, Cassette 20 teeth with hydraulic disk brake (Magura MT2) brake lever, right side and Shimano V-brake brake handle, left side



Rohloff "Speedhub 500/14" 14-speed hub gear

- A9113000016 € 1,679.00
 "Grip Shift" twist grip shifter
 Gear ratio: Sprocket 32 teeth, Cassette 16 teeth with hydraulic disk brake (Magura MT2) brake lever, right side and Shimano V-brake brake handle, left side

Dérailleur

Shimano Deore 1x11 Rapid Fire 11-speed dérailleur

A9113200034 € 492.50

mechanical dérailleur gear with operation on the handle
11-speed cassette 11-42 teeth
with hydraulic disk brake (Magura MT4) brake lever, right side
and Shimano V-brake brake handle, left side
top chain ring with chain guards

Chain ring selection:

- 32 teeth
 38 teeth

Shimano Deore XT 1x11 Di2 with button operation 11-speed dérailleur

A9113200035 € 1,713.50

electronic dérailleur with operation on the handle
11-speed cassette 11-42 teeth
with hydraulic disk brake (Magura MT4) brake lever, right side
and Shimano V-brake brake handle, left side
top chain ring with chain guards

Chain ring selection:

- 32 teeth
 38 teeth



Shimano Deore XT 1x11 Di2 With chin operation 11-speed dérailleur

A9113200033 € 1,863.00

Electronic dérailleur gear with chin operation (tetra equipment).
11-speed cassette 11-42 teeth
with hydraulic disk brake (Magura MT4) brake lever, right side
and Shimano V-brake brake handle, left side
top chain ring with chain guards

Chain ring selection:

- 32 teeth
 38 teeth



Gear change accessories

Mountain drive 2.5: 1 Gear reduction for hills

A9100200003 € 759.00

Planetary gear set for doubling of the number of gears
permits switching to the gear reduction for hills when stationary. With 38 teeth chain wheel as standard.
Only in combination with grip size 40 or 45 cm



Brake system

Powerful brake system with brake disc Ø 203 mm

A9110300105 € 79.10

Recommended from a body weight of 80 kg



Dimensions

Crank type and grip width










Bottom bracket shaft with Standard crank

- Grip width 40 cm
A9110500100-040
 Grip width 45 cm
A9110500100-045

Bottom bracket shaft with Tetra pedal crank (flared)

With flared crank arms for max. space for upper body

- Grip width 48 cm (176 mm)
A9110500100-048
 Grip width 53 cm (226 mm)
A9110500100-053

<p>Standard pedal crank</p> <ul style="list-style-type: none"> <input type="checkbox"/> Length 155 mm A9110600100-155 <input type="checkbox"/> Length 175 mm A9110600100-175 <input type="checkbox"/> Length 195 mm A9110600100-195 	<p>Tetra pedal crank</p> <ul style="list-style-type: none"> <input type="checkbox"/> Length 155 mm A9100600105-155 € 129.50 <input type="checkbox"/> Length 175 mm A9100600105-175 € 129.50 <input type="checkbox"/> Individual length 140 mm to 195 mm € 246.30 A9100600100 	
<p>Fixed bottom bracket support</p> <ul style="list-style-type: none"> <input type="checkbox"/> XS (X/Y 113/130 mm) A9100501100-004 <input type="checkbox"/> S (X/Y 67/153 mm) A9100501100-003 	<p>Adjustable bottom bracket support</p> <p>the angle and length can be adjusted via the clamped rotary point on the fork bridge and length</p> <ul style="list-style-type: none"> <input type="checkbox"/> S (180 mm) A9100500001-002 € 220.90 <input type="checkbox"/> M (240 mm) A9100500001-003 € 220.90 <p>If no details are provided, the adjustable bottom bracket support shall be supplied in M.</p> 	
Leg rest		
<p>Leg rest separated</p> <ul style="list-style-type: none"> <input type="checkbox"/> A9100100008 Adaptation on the fork for short legs Only possible up to body size 165 cm 	<p>Foot stopper for leg rest</p> <ul style="list-style-type: none"> <input type="checkbox"/> A9100100012 € 86.80 	<p>Thigh anchorage</p> <ul style="list-style-type: none"> <input type="checkbox"/> A9100100013 € 204.70 
Back equipment		
<p>Backrest with filling plate and back cushion</p> <ul style="list-style-type: none"> <input type="checkbox"/> A9100800002 Cover: Black, breathable, open surface 	<p>Body contour backrest upholstery</p> <ul style="list-style-type: none"> <input type="checkbox"/> A9100800010 € 85.10 with back cushion 3 cm Cover: Black, breathable, open surface 	

Seat upholstery

Open belt system seat upholstery

A9100801010



Body Contour upholstery

A9100801025
breathable, cushioning effect



Seat cushion

“Standard” seat cushion

A9113000190 **€ 106.60**

For 102 cm chassis
Seat width: 35 cm
Height: 4 cm
Cover: Black, breathable, open surface



Tyres

20x1.1" Schwalbe “Durano”

A9113105001
(ETRTO: 406-28)
Race bike tyres



20x1.75" Schwalbe “Marathon Plus”

A9113105120 **€ 67.60**
(ETRTO: 406-47)
unpuncturable
with Smart Guard anti-puncture belt



Accessories

Plug-in protection plate

A9100502100 **€ 55.30**



Leg guard

A9100100231 **€ 162.70**
Guard for MTB tyres
Screwed onto the fork tube
Colour: Matt black

Collision guard clamped on axle tube

A9100800100 **€ 207.50**
with 40 mm mounting clip



Neck rest

A9100800150 **€ 248.40**
height-adjustable
angle-adjustable



Hip safety belt

A9100800200 **€ 112.90**



SIGMA BC 9.16 ATS Bike Computer

A9110901000 **€ 66.80**
wireless, with bracket
Not in combination with e-assistant version



**SIGMA BC 14.16. STS CAD
Bike Computer**

A9110900018 **€ 163.20**

with pedal frequency set
wireless
Not in combination with e-assistant version



**Bracket ø 22 mm x 42 mm for fixing of
accessories
(speedometer, lighting or similar)**

A9110903271 **Item price €
22.10**

Colour: Natural anodised

_____ Items



SIGMA lighting set

A9120900019 **€ 117.90**

Consisting of front light and rear lights,
chargeable via Micro USB socket,
incl. Micro USB charging cable.
Lighting set compliant with road traffic
regulations



Options for users with reduced function

**Tetra hand support, pivotable,
adjustable grip distance**

A9110601101 **€ 161.30**



Tetra shell handles

A9110601100 **€ 164.10**



**Tetra hand support, pivotable,
adjustable grip width**

A9110601102 **€ 182.40**

Very stiff outer support as round tube ø 20 mm
with foam rubber cover



Notes and remarks:

Location, Date

Signature End customer/dealer

The signature field is intended for the consultant and insured person. This confirms that the intended equipment is intended for the insured person.

- 6 - 2023-02-16 / The prices are plus statutory value-added tax and shipping costs = Included in the basic price = Position required

Order sheet FREAK compact bike Rev. 05 UK