



NJ1 adaptive bike

up to 120 kg trailer load and 10 kg payload

- Docking plate for quick coupling to the „Quick-Snap“ adapter
- A-shaped handles with integrated shift and brake fittings
- Additional V-brake (rim brake) with parking lock
- Parking stand
- Chain guard tube
- Black chain guard
- Mudguard
- Bottle cage incl. water bottle
- Sigma lighting set with battery charging unit, front and rear lights

Quotation date: _____

Consultation was carried out by
Staff at the medical supplies shop: _____

Enquiry no.: _____

Order date: _____

according to GTC PRO ACTIV Reha-Technik GmbH

Order no.: _____

Medical suppliers / company address:

Cust. no.: _____

Name: _____

Street: _____

Postcode / city: _____

Phone: _____

Fax: _____

E-Mail: _____

Contact: _____

Different invoicing address (see below)

Cust. no.: _____

Name: _____

Street: _____

Postcode / city: _____

Different shipping address (see below)

Cust. no.: _____

Name: _____

Street: _____

Postcode / city: _____

Customer name / commission:

Surname, forename: _____

Street: _____

Postcode / city: _____

Phone: _____

E-Mail: _____

D.o.B.: _____

Paying organisation: _____

Indication: _____

Body size

Body height: _____ cm

Upper body height: _____ cm

Weight: _____ kg

This information is necessary for checking.

Drive unit

Drive unit black

- 20" (ETRTO 406) A9110100001-020 € 3,890.00
- 24" (ETRTO 507) A9110100001-024 € 3,890.00

Hub gear with integrated back-pedal brake

Shimano „Nexus“ 8-speed hub gear

- A9110200023
307% gear range
„Grip-Shift“ twist grip shifter
Gear ratio: 32-tooth chainring,
20-tooth cassette
(19 teeth for size L bottom bracket)



Derailleur gears with disc brake

Shimano Deore 1x11 Rapid Fire 11-speed derailleur

- A9110200032
mechanical derailleur gear with operation on
the handle. 11-speed cassette 11-42 teeth

Selection of the chainring:

- 32 teeth recommended for
24" drive wheel
- 38 teeth recommended for
20" drive wheel
- 42 teeth for lower crank
frequency or driving
at higher speeds



Shimano Deore XT 1x11 Rapid Fire 11-speed derailleur

- A9110200034 € 192.20
mechanical derailleur gear with control on the
handle. 11-speed cassette 11-42 teeth
top chain ring with chain guards

Selection of the chainring:

- 32 teeth recommended for
24" drive wheel
- 38 teeth recommended for
20" drive wheel
- 42 teeth for lower crank
frequency or driving
at higher speeds

Shimano Deore XT 1x11 D12 with push button operation 11-speed derailleur

- A9110200035 € 1,486.10
electronic derailleur gear with operation on the
handle. 11-speed cassette 11-42 teeth

Operation: right side left side
If not specified, the operation will be mounted
on the „right side“.

Selection of the chainring:

- 32 teeth recommended for
24" drive wheel
- 38 teeth recommended for
20" drive wheel
- 42 teeth for lower crank
frequency or driving
at higher speeds



Hub gear with disc brake

Shimano „Alfine“ 11-speed hub gear

- A9100200026 € 592.00
„Rapid Fire“ handle
Gear ratio: Sprocket 32 teeth,
Cassette 19 teeth



Rohloff „Speedhub 500/14“ 14-speed hub gear

- A9100200027 € 1,498.10
„Grip Shift“ twist grip shifter
Gear ratio: Sprocket 32 teeth,
Cassette 16 teeth



Gear change accessories

Mountain drive 2.5 : 1 Hill reduction

A9100200003 € 797.40

Planetary gearbox for doubling the gears.
Allows you to shift downhill while standing still.
As standard with 38 teeth chainring.
In combination with grip width 40 cm, 45 cm,
48 cm or 53 cm



Brake system

Mechanical parking disk brake can be operated via the spreader

A9110300002 € 84.60

mounted on the bottom bracket support,
brake handle on the handle, left side is omitted



Driving brake - supplement for the mechanical parking brake


A9117300120 € 116.10

The parking disc brake is supplemented with a
brake handle on the left. Thus, the mechanical
disc brake can be used as a driving brake and
parking brake.

Powerful braking system with brake disc Ø 203 mm

A9110300105 € 83.10

Recommended from 80 kg trailer load or body
weight.
Not required if the hydraulic coaster disc brake
is selected
Not in combination with a Rohloff gear system



Double powerful brake system brake disc Ø 203 mm

A9110300100-024 € 347.10

including V-Brake rim brake can be operated
as a parking brake via towing spreader
Recommended from 100 kg trailer load

Mechanical backpedal disk brake with crank release function

A9110300004 € 738.40

compatible with derailleur and free-running
shifting hub.
In combination with:
- Grip width 45 cm, 50 cm, 53 cm or 58 cm
- Pedal crank crank lengths 155 mm or 175 mm
up to 120 kg body weight
Not in combination with Mountaindrive 2.5 : 1
gear reduction for hills

Tyres

Schwalbe „Marathon Plus“

20x1,75" (406-47)
A9130105128

24x1,75" (507-47) € 28.70
A9110105135

unpuncturable, with anti-puncture belt
Smart Guard



Dimensions

Crank type and grip width

Bottom bracket shaft with standard crank

- grip width 40 cm
A9110500100-040
- grip width 45 cm
A9110500100-045
- grip width 50 cm
A9110500100-050

Bottom bracket shaft with tetra crank (flared)

with extended crank arms for maximum upper body space

- grip width 48 cm (176 mm)
A9110500100-048
- grip width 53 cm (226 mm)
A9110500100-053
- grip width 58 cm (276 mm)
A9110500100-058

Standard pedal crank

- length 155 mm
A9110600100-155
- length 175 mm
A9110600100-175
- length 195 mm
A9110600100-195



Tetra pedal crank

- length 155 mm
A9100600105-155 **€ 136.10**
- length 175 mm
A9100600105-175 **€ 136.10**
- individual length
195 mm to 140 mm **€ 258.70**
A9100600100



Fixed bottom bracket support

- L (X/Y 158/362 mm)
A9110501100-002
- XL (X/Y 158/410 mm) **€ 66.70**
A9110501100-006
- XXL (X/Y 170/465 mm) **€ 66.70**
A9110501100-007



Adjustable bottom bracket support

the angle and length can be adjusted via the clamped rotary point on the fork bridge and length.

For the pre-setting, please specify the bottom bracket support size that comes as close to your desired position.

- M (300 mm) **€ 232.10**
A9110500001-001
- L (370 mm) **€ 232.10**
A9110500001-002
- XL (420 mm) **€ 232.10**
A9110500001-003

If not specified, the adjustable bottom bracket support will be delivered in XL. Available up to size XL.



Handles

Palm rest on the handle

- A9110601200 € 102.80
upholstered, bolted, angle-adjustable.



Tetra hand support, pivotable, adjustable grip distance

- A9110601101 € 169.40



Tetra backrest swivel and adjustable grip width

- A9110601102 € 191.60

very stiff external support as round tube
Ø 20 mm with foam rubber coating
Not in combination with Tetra hand support
pivotable



Tetra ergonomic handle, freestanding

- A9100600001 € 164.60
crank handles arranged vertically
can be combined with all other handle options



Tetra shell handles

- A9110601100 € 172.40



Accessories Bike

Black aluminum mudguard

- A9116500005 € 112.40
Only in combination with 20" drive wheel

Mudguard and chain guard customized colour selection

- A9116500010 € 248.60
from colour palette one-colour design
powder colours wheelchair

A9000100170- _ _ _

Colour: _____

Only in combination with 20" drive wheel

Bracket Ø 22 mm x 42 mm for mounting accessories (tacho, lighting, etc.)

- A9110903271 pcs. á € 23.20

Colour: natural anodized

_____ piece



Maneuvering castors for parking stand

- A9110900017 € 122.10

Maneuvering castors mounted on the parking
stand via inlet



Maneuvering handle mounted on the chain guard

- A9110500000 € 64.00

with foam rubber coating



Aluminum luggage rack „cargo systems”

- A9110903290 € 188.80

Colour: black

max. load weight on the loading platform 10 kg,
max. 15 kg in conjunction with luggage bags at
lower loading height.



SIGMA BC 8.0 ATS bike computer

- A9110901000 € 70.10
wireless, with mount



SIGMA BC 14.0 WL CAD bike computer

- A9110900018 € 171.50
with pedal frequency set, wireless



Options for users with reduced functions

Lever extension excenter lock for Quick Snap connection


Makes it easier to open and close the Quick Snap adapter. Tetra version

A9110700005 € 73.90

straight design
only in combination with rigid frame adapter

A9110700006 € 73.90

curved design
only in combination with front adapter



Adaptation help for centering pin

A9110700002 € 55.00
to simplify the insertion of the central adapter.



Parking stand with operating lever for high operating comfort

A9110900004 € 58.70
with longer operating lever

Chest strap for upper body fixation

size S (105 cm) A9000800016 € 238.80

size M (125 cm) A9000800013 € 238.80

size L (145 cm) A9000800014 € 238.80

size XL (165 cm) A9000800017 € 238.80

consisting of elastic webbing material in combination with neoprene with Velcro and operating tab.
Colour: black



Frame adapter

SPEEDY rigid frame adapter with „Quick-Snap-Adaption“

A9110700141 € 537.20

with quick-release clamp and centering pin holder
Colour: black

Please select accordingly:

Centring pin for Ø 25 mm
for A1 and A2

A9110700022

Centring pin for Ø 40 mm
for 4all, 4all Ergo, 4you Ergo, F2, F3 and 4teen

A9110700023



TRAVELER folding frame adapter with „Quick-Snap-Adaption“

A9110701300 € 537.20

with quick-release clamp and centering pin holder on the folding mechanism
Colour: black



BUDDY rigid frame adapter with „Quick-Snap-Adaption“

A9110701400 € 537.20

with quick-release clamp and centering pin mount on the frame
Colour: black



Wheelchair front adapter with universal clamp and „Quick-Snap-Adaption“

A9110701500 € 572.30

Adapter attachment to the wheelchair stem. Operation on the right-hand side.

Mounting situation for PA products:
Straight front end from SB 31 cm to 54 cm
V-shaped front end from SB 30 cm to 50 cm

Please select the frame clamp below:

clamping range 95 mm
A9110701525
Tube diameter 23-27 mm

clamping range 95 mm
A9110701530
Tube diameter 28-32 mm

clamping range 95 mm
A9110701535
Tube diameter 33-37 mm

PLEASE NOTE:

If the exact tube diameter is not specified, an adapter for a tube diameter of 23 - 27 mm will be delivered as standard.

Please select the clamping sleeve below:

clamping sleeve [80 mm] included as standard
for close clearance between footrest and drive wheel.

clamping sleeve [110 mm] € 46.70
A9199000117
for medium distance between footrest and drive wheel.

clamping sleeve [170 mm] € 56.60
A9199000154
for a wide distance between footrest and drive wheel.

PLEASE NOTE:

If not specified, the clamping sleeve with 80 mm length will be delivered.

Clamp lever with 3-way star handle for actuating the frame clamping

A9199701507 € 106.10

mounted on the right and left for easy terminal operation with limited manual function
Only in combination with front adapter



Wheelchair-dependent accessories

Extended bike plate

for SPEEDY A1 / A2 A9110100100-002 € 370.00

for BUDDY A9110100100-003 € 370.00

for TRAVELER A9110100100-001 € 370.00

4° camber integrated
Track independently adjustable



Wheelbase extension via „grid mechanism with pulling operation“

for SPEEDY 4all / 4teen / 4you / TRAVELER 4you A9110102001 € 564.70

removable, with continuous axle tube, track width independently adjustable

Please select the wheel camber:

- 1° drive wheel camber
- 4° drive wheel camber
- 6° drive wheel camber



Wheelbase extension angle-adjustable and removable

for SPEEDY F2 / F3 A9015100005 € 746.80

with continuous axle tube, height adjustment via gear joint adjustable in 6° angular increments
wheel position displacement approx. 14 cm adjustable track
Attachment possibility must be checked by PRO ACTIV.

Only in combination with base plate. Up to 120 kg load capacity.

Please select the wheel camber below:

- 1° drive wheel camber
- 4° drive wheel camber
- 6° drive wheel camber



Services

Consultation, clarification of the installation situation and measurement through PRO ACTIV

A9990000000 € 354.50

Handover and instruction of the handbike / traction device by PRO ACTIV

A9990000011 € 354.50

Attachment of the handbike / traction device to the wheelchair by PRO ACTIV at specialist dealer

on a PRO ACTIV wheelchair

A9990000012 € 230.00

product no.: _____

on a third-party product

A9990000013 € 366.90

manufacturer: _____

model: _____

frame tube Ø ____ mm

Attachment of the handbike / traction device to the wheelchair through specialist dealer

including mounting of the frame adapter

on a PRO ACTIV wheelchair

A9990000008

product no.: _____

on a third-party product

A9990000005

manufacturer: _____

model: _____

frame tube Ø ____ mm

Notes and remarks:

Location, Date

Signature End customer/dealer

The signature field is intended for the consultant and insured person. This confirms that the intended equipment is intended for the insured person.

- 8 - 2026-03-03 / The prices are plus statutory value-added tax and shipping costs = included in the basic price = Position selected