



Quotation date: _____
Consultation was carried out by
Staff at the medical supplies shop: _____
PRO ACTIV field representative: _____

Order date: _____
According to GTC PRO ACTIV Reha-Technik GmbH

Medical suppliers/company address:

Cust. no.: _____
Name: _____
Street: _____
Postcode / city: _____
Phone: _____
Fax: _____
E-Mail: _____
Contact: _____

Invoicing address, other(see below)

Cust. no.: _____
Name: _____
Street: _____
Postcode / city: _____

Shipping address, other (see below)

Cust. no.: _____
Name: _____
Street: _____
Postcode / city: _____

Customer name/commission:

Surname,
forename: _____
Street: _____
Postcode / city: _____
Phone: _____
E-Mail: _____
D.O.B: _____
Paying
organisation: _____
Indication: _____



SPIKE adaptive bike

approved up to 75 kg trailer load and 5 kg load capacity

- docking plate for quick coupling to the "Quick-Snap" adapter
- drive unit with steering fork
- A-shaped handles with integrated switching and brake fittings
- chain guard in black on both sides
- without gear shift
- bottom bracket support can be adjusted in height and angle
- parking stand
- Sigma lighting set with rechargeable battery charger, front and rear lights

A911400000 € 2,550.00

Basic equipment:

- Free-wheel hub with single speed kit, 15 teeth
- Hydraulic disc brakes with Ø 180 mm disc
- Mech. Shimano disk brake, left brake handle

Please note: When purchasing a SPIKE adaptive bike in combination with a wheelchair, the discounted VAT rate applies for the complete unit.

Body size

Body height: _____ cm
Upper body height: _____ cm
Weight: _____ kg

Information essential to process the order.

Fork colour

Information about current colour designs and example images at: www.proactiv-gmbh.de/farben.php

White RAL 9016

- gloss A9000100170-023
 matt A9000100170-123 € 40.10

Brilliant silver

- gloss A9000100170-011

Anthracite RAL 7016

- gloss A9000100170-001
 matt A9000100170-101 € 40.10

Black RAL 9005

- gloss A9000100170-003
 matt A9000100170-019 € 40.10

Deep black, silver effect

- gloss A9000100170-007

Pure orange RAL 2004

- gloss A9000100170-022

Fire Red

- gloss A9000100170-018
 matt A9000100170-118 € 40.10

Pink magenta RAL 4010

- gloss A9000100170-037
 matt A9000100170-137 € 40.10

Sky blue

- gloss A9000100170-017

Pastel blue RAL 5024

- gloss A9000100170-036
 matt A9000100170-136 € 40.10

Marine blue

- gloss A9000100170-045
 matt A9000100170-145 € 40.10

Electric green

- gloss A9000100170-015
 matt A9000100170-115 € 40.10

Lemon green

- gloss A9000100170-005
 matt A9000100170-105 € 40.10

Violet

- gloss A9000100170-044

Yellow RAL 1023

- gloss A9000100170-046
 matt A9000100170-146 € 40.10

Red RAL 3001

- gloss A9000100170-047
 matt A9000100170-147 € 40.10

Lagoon blue

- gloss A9000100170-052
 matt A9000100170-152 € 40.10

Athletic grey

- gloss A9000100170-050
 matt A9000100170-150 € 40.10

Muddy grey

- gloss A9000100170-051
 matt A9000100170-151 € 40.10

Equipment

Shimano Deore 1x11 Rapid Fire 11-speed dérailleur

- A9114200032 € 384.20
 mechanical dérailleur gear with operation on the handle
 11-speed cassette 11-42 teeth

Shimano "Nexus" 7-speed hub gear "Grip Shift" twist grip shifter

- A9114000007 € 516.70
 with V-Brake rim brake, brake handle, left side



Shimano "Nexus" 8-speed hub gear "Grip Shift" twist grip shifter

- A9114000008 € 684.10
 with V-Brake rim brake, brake handle, left side



Mountain drive 2.5: 1 gear reduction for hills In combination with hub gear


- A9100200003 € 759.00
 Planetary gear set for doubling of the number of gears
 All gears (the hub gear and the gear reduction for hills) can be changed when stationary!
 In combination with grip width 40 or 45 cm (for 7-speed or 8-speed hub gears)



Dimensions

Bottom bracket shaft with Standard crank

- Grip width 36 cm
A9114500100-036
- Grip width 40 cm
A9110500100-040
- Grip width 45 cm
A9110500100-045



Standard pedal crank

- Length 135 mm
A9110600100-135
- Length 155 mm
A9110600100-155
- Length 175 mm
A9110600100-175

Tyres

20x2" Schwalbe "Big Apple"

- A9114105125
(ETRTO 406-50)
Balloon tyres



20x1.75" Schwalbe "Marathon"

- A9110105121
(ETRTO 406-47)
Touring tyres



Puncture protection insert for pneumatic tyres

- A9000206100-K **€ 51.60**

Wheelchair-dependent accessories

Extended bike plate

- For SPEEDY A1/A2 **€ 352.20**
A9110100100-002
- For BUDDY **€ 352.20**
A9110100100-003
- For TRAVELER **€ 352.20**
A9110100100-001

4° camber integrated
Track can be adjusted independently



Wheelbase extension via "grid mechanism with pulling operation"

- For SPEEDY 4all /
4teen / 4you /
TRAVELER 4you **€ 537.50**
A9110102001

removable, with continuous axle tube
4°/ 6° camber integrated
Track can be adjusted independently



Frame adapters

SPEEDY fixed frame adapter Left operation with "Quick Snap Adaptation"

- A9110700141 **€ 492.50**
with quick release clips and centring column
mount
Colour: black



TRAVELER folding frame adapter with "Quick Snap Adaptation"

- A9110701300 **€ 526.70**
with quick release clips and centring column
mount on the folding mechanism
Colour: black



BUDDY fixed frame adapter with "Quick Snap Adaptation"

- A9110701400 **€ 492.50**
with quick release clips and centring column
mount on the frame
Colour: black



Wheelchair front adapter with universal clamping and "Quick Snap Adaptation"

A9110701500 € 561.10

Adapter mount on the wheelchair front end
Operation on the right side

Installation situation with PA products:
Frame stem G-shape from SB 31 cm to 54 cm
Frame stem V-shape from SB 30 cm to 50 cm

PLEASE NOTE:

If the tube diameter is not specified, as standard we shall deliver an adapter for a tube diameter of 23 - 27 mm.



Please make the following selection:

A9110701525
Frame tube diameter 23-27 mm

A9110701530
Frame tube diameter 28-32 mm

A9110701535
Frame tube diameter 33-37 mm

Clamp lever with 3-way star handle for actuating the frame clamping

A9199701507 € 101.00

mounted on the right and left for easy terminal operation with limited manual function
Only in conjunction with front adapter



Accessories

Plug-in protection plate

made from plastic
incl. required mounting components

A9122502100 € 104.10

in combination with luggage rack

A9122502101 € 55.30



Aluminium mudguard

incl. required mounting components

A9122500005 € 146.10

in combination with luggage rack

A9122500006 € 107.00

SIGMA BC 9.16 ATS Bike Computer

A9110901000 € 66.80

wireless, with bracket



Chain guard in frame colour

A9114500005 € 81.70

Aluminium luggage rack "cargo systems"

A9120900000 € 199.90

Colour: black
Max. 10 kg load weight
max. 15 kg in combination with luggage bags with low loading height



Manoeuvring rollers for parking stands

A9110900017 € 116.20

Manoeuvring rollers mounted to the parking stand via inlet



Options for users with reduced function

Tetra hand support pivotable and adjustable grip width

A9110601101 € 161.30



Tetra shell handles

A9110601100 € 164.10




Lever extension quick release clip for rigid frame adapters

A9110700004 € 98.90

to make the opening and closing of the quick release clips easier



<p>Securing the lever extension with locking bolt</p> <p><input type="checkbox"/> A9199000084 € 82.90</p> <p>Can only be used for the lever extension on the quick release clip.</p>	<p>Lever extension eccentric lock for Quick Snap connection simplifies opening or closing the Quick Snap Adapter. Tetra version</p> <p><input type="checkbox"/> A9110700005 € 70.30 Straight version Only in combination with rigid frame adapter</p> <p><input type="checkbox"/> A9110700006 € 70.30 Curved version Only in conjunction with front adapter</p> 	<p>Centring column adaptation assistance</p> <p><input type="checkbox"/> A9110700002 € 52.40 in order to make inserting the central adapter easier</p> 
---	--	---

<p>Chest strap for upper body fixation</p> <p><input type="checkbox"/> Size S (105 cm) € 220.80 A9000800016</p> <p><input type="checkbox"/> Size M (125 cm) € 220.80 A9000800013</p> <p><input type="checkbox"/> Size L (145 cm) € 220.80 A9000800014</p> <p><input type="checkbox"/> Size XL (165 cm) € 220.80 A9000800017</p> <p>Consisting of elastic webbing material in combination with neoprene with Velcro and operating tab. Colour: black</p> 	
--	--

Services

<p>Consultation, clarification of the installation situation and measurement by PRO ACTIV</p> <p><input type="checkbox"/> A9990000000 € 327.80</p>	<p>Handover and instruction for the handbike/traction device by PRO ACTIV incl. mounting the frame adapter</p> <p>on a PRO ACTIV wheelchair</p> <p><input type="checkbox"/> A9990000007 € 540.50 Prod. no.: _____</p> <p>on a third-party product</p> <p><input type="checkbox"/> A9990000009 € 667.00 Manufacturer: _____ Model: _____ Frame tube ø ____ mm</p>	<p>Handover and instruction for the handbike/traction device by the dealer incl. mounting the frame adapter</p> <p>on a PRO ACTIV wheelchair</p> <p><input checked="" type="checkbox"/> A9990000008 Prod. no.: _____</p> <p>on a third-party product</p> <p><input checked="" type="checkbox"/> A9990000005 Manufacturer: _____ Model: _____ Frame tube ø ____ mm</p>
---	---	--

Notes and remarks: