

Customized production - model „TRAVELER“

Frame design

Frame design Heavy Duty T1-4 up to 150 kg body weight

A9030100300-K € 634.30

- Reinforced drive wheel plate, 2 slots, 12 mm thick
- Back tube reinforcement insert on side part screw thread
- Stainless steel backrest locking bolts
- Stainless steel push handles
- Side part thickness 3 mm aluminium
- Footrest or foot bar, continuous, reinforced design
- Knee lever brake bolt, stainless steel

Note:

Please select the „Heavy Duty“ steering wheel fork under Steering wheel fork!

Heavy Duty reinforcement measures on frame and add-on components for payloads of 150 kg - 180 kg

A9030100301-K € 170.50

Inner frame tube in seat tube
Frame gusset plates in the area of the folding arm attachment

Is required additionally to A9030100300-K.

Frame

Frame extension abducted

A9030100071 € 1,869.00

Frame extension according to customer specifications
Frame extension abducted, special frame length abduction: 2,5 cm / side
Frame length to frame arch: ____ cm

The folding function of the backrest might be limited.

Customized frame length


A9030100070 € 503.10

Frame length to frame arch ____ cm
(for ST ____ cm)

Customized frame height

A9030100072 € 3,255.00

Frame geometry in high version:
Frame height + 6 cm (not variable)



Individual seat width

Seat width 54 cm

A9030100061-054 € 2,282.20

Seat width 56 cm

A9030100061-056 € 2,282.20

Knee angle

Individual knee angle

A9000100050-003 € 534.50

90° 100° 110° 115°

Services

Consulting and customization by PRO ACTIV according to individual customer requirements

A9990000017 € 612.00

Additional requirements