



NJ1 el. compactbike

Maximum user weight limit 100 kg

- girder chassis with integrated rechargeable battery holder
- drive unit, can be decoupled via the eccentric quick-lock system
- A-shaped handles with integrated shift and brake controls
- Backrest continuously angle-adjustable
- 6° wheel camber
- Hydraulic disc brake system
- additional mechanical disc brake with parking lock
- Safety flag with holder
- Leg protection guard
- Drive unit black

Neodrives Z15 drive unit

„The powerful hub drive for maximum performance with TFT display and 11-speed Shimano derailleur gears”

torque 12/40 Nm, power 250 watt

Wheel size: 26”

A9100100011-026 € 8,780.00

Quotation date: _____

Consultation was carried out by
Staff at the medical supplies shop: _____

Order date: _____

according to GTC PRO ACTIV Reha-Technik GmbH

Medical suppliers / company address:

Cust. no.: _____

Name: _____

Street: _____

Postcode / city: _____

Phone: _____

Fax: _____

E-Mail: _____

Contact: _____

Different invoice address (see below)

Cust. no.: _____

Name: _____

Street: _____

Postcode / city: _____

Different shipping address (see below)

Cust. no.: _____

Name: _____

Street: _____

Postcode / city: _____

Customer name / commission:

Surname, forename: _____

Street: _____

Postcode / city: _____

Phone: _____

E-Mail: _____

D.O.B: _____

Paying organisation: _____

Indication: _____

Body size

Body height: _____ cm

Weight: _____ kg

Chassis colour

Information on current colour designs and example images at: www.proactiv-gmbh.de/farben.php

White RAL 9016

- gloss A9000100170-023
 matt A9000100170-123

€ 42.10

Brilliant silver

- gloss A9000100170-011

Anthracite RAL 7016

- gloss A9000100170-001
 matt A9000100170-101

€ 42.10

Black RAL 9005

- gloss A9000100170-003
 matt A9000100170-019

€ 42.10

Deep black effect silver

- gloss A9000100170-007

Pure orange RAL 2004

- gloss A9000100170-022

Fire red

- gloss A9000100170-018
 matt A9000100170-118

€ 42.10

Pink magenta RAL 4010

- gloss A9000100170-037
 matt A9000100170-137

€ 42.10

Sky blue

- gloss A9000100170-017

Pastel blue RAL 5024

- gloss A9000100170-036
 matt A9000100170-136

€ 42.10

Navy blue

- gloss A9000100170-045
 matt A9000100170-145

€ 42.10

Electric green

- gloss A9000100170-015
 matt A9000100170-115

€ 42.10

Lemon green

- gloss A9000100170-005
 matt A9000100170-105

€ 42.10

Violet

- gloss A9000100170-044

Blasted transparent coated

- gloss A9000100170-029

Yellow RAL 1023

- gloss A9000100170-046
 matt A9000100170-146

€ 42.10

Red RAL 3001

- gloss A9000100170-047
 matt A9000100170-147

€ 42.10

Lagoon blue

- gloss A9000100170-052
 matt A9000100170-152

€ 42.10

Athletic grey

- gloss A9000100170-050
 matt A9000100170-150

€ 42.10

Muddy grey

- gloss A9000100170-051
 matt A9000100170-151

€ 42.10

Battery Neodrives

The specified ranges are dependent on the route profile, user weight and speed.

Battery neodrives DT 36V 16Ah 576Wh

range approx. 75 km

- A9100900109



Double battery neodrives DT 36V 16Ah 576Wh

range approx. 150 km

- A9100900110 € 1,371.30

Plug-in device for easy switching between batteries via the magnetic connector

Derailleur gear

Shimano Deore 1x11 Rapid Fire 11-speed derailleur

- A9110200032

mechanical derailleur gear with operation on the handle. 11-speed cassette 11-42 teeth
Top chainring with chain guards

Selection of the chainring:

- 32 teeth
 38 teeth



Shimano Deore XT 1x11 Rapid Fire 11-speed derailleur

- A9110200034 € 186.60

mechanical derailleur gear with operation on the handle. 11-speed cassette 11-42 teeth

Selection of the chainring:

- 32 teeth
 38 teeth

Shimano Deore XT 1x11 Di2 with push button operation 11-speed derailleur

- A9110200035 € 1,442.80

electronic derailleur gear with operation on the handle. 11-speed cassette 11-42 teeth

Operation: right side left side

Selection of the chainring:

- 32 teeth
 38 teeth



**Shimano Deore XT 1x11 DI2
with chin control
11-speed derailleur**

A9110200033 € 1,583.60

electronic derailleur gear with chin control
(tetra equipment).
11-speed cassette 11-42 teeth

Selection of the chainring:

32 teeth

38 teeth



Gearshift accessories

**Mountaindrive 2.5 : 1
Hill reduction**

A9116200003 € 703.80

for extension to double the number of gears
Planetary gearbox
Allows you to shift downhill while standing still.
In conjunction with grip width 40, 45, 48 or
53 cm



Brake system

**Powerful braking system
with brake disc Ø 203 mm**

A9110300105 € 80.70

Recommended from 80 kg bodyweight
Not required if mechanical back-pedal disc
brake is selected



**Double disc brake system „Magura“
brake disc Ø 203 mm**

A9100300100 € 226.50

eliminates the mechanical disk brake.
Only in conjunction with the powerful braking
system.
Recommended from 100 kg bodyweight

**Mechanical back-pedal disc brake with
crank-free function**

A9110300004 € 716.90

compatible with derailleur and freewheel hub.
In conjunction with:
- Bottom bracket shaft reach 45 to 55 cm
- pedal crank length 155 or 175 mm
up to 120 kg body weight
Not in conjunction with Mountaindrive 2.5:1 hill
reduction!

Dimensions

Crank type and grip width

Standard crank

- Grip width 40 cm
A9110500100-040
- Grip width 45 cm
A9110500100-045
- Grip width 50 cm
A9110500100-050

Tetra crank (flared)

with extended crank arms for maximum upper body space

- Grip width 48 cm (176 mm)
A9110500100-048
- Grip width 53 cm (226 mm)
A9110500100-053
- Grip width 58 cm (276 mm)
A9110500100-058

Standard pedal crank

- Length 155 mm
A9110600100-155
- Length 175 mm
A9110600100-175
- Length 195 mm
A9110600100-195



Tetra pedal crank

- Length 155 mm
A9100600105-155 € 132.10
- Length 175 mm
A9100600105-175 € 132.10
- individual length
195 mm to 140 mm € 251.20
A9100600100




Fixed bottom bracket support

(Limited availability - until stocks last.

Alternative: adjustable bottom bracket support)

- S (X/Y 67/153)
A9100501100-003
- M (X/Y 90/250)
A9100501100-001



Adjustable bottom bracket support

adjustable in angle via clamped pivot point at the fork bridge and in length

Please specify the bottom bracket support size for the presetting that comes closest to your desired position.

- S (180 mm)
A9100500001-002 € 225.30
- M (240 mm)
A9100500001-003 € 225.30

Please note: If not specified, the adjustable bottom bracket support will be delivered in S.



Legrest

Divided legrest

- A9100100008
adaptation on the fork
for short legs
Only possible up to 165 cm body height



All-round legrest

- A9100100009
adaptation on the fork



Foot stop for legrest

- A9100100012 € 88.50



Upper thigh fixation

A9100100013 € 208.80



Back equipment

Backrest with filler panel and back cushion

A9100800002
Cover: Black, breathable, open surface

Back upholstery Body Contour


A9100800010 € 86.80
with back cushion 3 cm
Cover: Black, breathable, open surface



Chassis

Chassis geometry

- Length 116 cm
SB approx. 39 cm
A9100100000-005
- Length 116 cm
SB approx. 43 cm € 198.60
A9100100000-006
- Heavy duty version up
to 120 kg body weight
A9100800013 € 413.20
recommended for off-road
use



Seat upholstery

Seat upholstery open belt system

Chassis: 116 cm
SB 39 cm, ST 40 cm
 A9100801010-002

Chassis: 116 cm
SB 43 cm, ST 40 cm
 A9100801010-003



Seat upholstery Body Contour

Chassis: 116 cm
SB 39 cm, ST 40 cm
 A9100801025-002

Chassis: 116 cm
SB 43 cm, ST 40 cm
 A9100801025-003




Seat cushion

Seat cushion standard

Cover: Black, breathable, open surface
Seat cushion height: 4 cm

Chassis: 116 cm
SB 39 cm, ST 40 cm
 A9100000190-002 € 108.70

Chassis: 116 cm
SB 43 cm, ST 40 cm
 A9100000190-003 € 108.70



Seat cushion with upper thigh guide

Cover: Black, breathable, open surface
 Seat cushion height: 4 cm

Only in conjunction with upper thigh fixation

Chassis: 116 cm
 SB 39 cm, ST 40 cm

A9100800103-002 € 198.50

Chassis: 116 cm
 SB 43 cm, ST 40 cm

A9100800103-003 € 198.50

**Tyres****26x1.5" Schwalbe „Marathon“**

A9100105100
 (ETRTO 559-40)
 Touring tyre

**26x2" Schwalbe „Marathon Plus Tour“**

A9100105140 € 69.00
 (ETRTO 559-50)
 Touring tyre

**26x2.1" Schwalbe „Smart Sam“**

A9100105150 € 72.40
 (ETRTO 559-54)

**Bike accessories****Plug-in mudguard**

A9100502100 € 56.40

**Collision guard clamped on axle tube**

A9100800100 € 211.70

**Neck rest**

A9100800150 € 253.40
 height-adjustable, angle-adjustable

**Pelvic fixation belt**

A9100800200 € 115.20

**Bracket ø 22 mm x 42 mm for mounting accessories (tacho, lighting, etc.)**

A9110903271 pcs. á € 22.50
 Colour: natural anodized
 ___ piece

**SIGMA lighting set**

A9120900019 € 120.30
 consisting of front light and rear lights,
 chargeable via Micro USB socket,
 incl. Micro USB charging cable.
 Lighting set compliant with road traffic
 regulations

**Options for users with reduced functions****Palm rest on the handle**

A9110601200 € 99.80
 upholstered, bolted, angle-adjustable.
 Not in conjunction with other handles.


**Tetra hand support, pivotable, adjustable grip distance**

A9110601101 € 164.50

**Tetra backrest swivel and adjustable grip width**

A9110601102 € 186.00
 very stiff external support as round tube
 ø 20 mm with foam rubber coating



<p>Tetra shell handle</p> <p><input type="checkbox"/> A9110601100 € 167.40</p> 	<p>Operating satellite mounted on handle</p> <p>Mounting: Right side <input type="checkbox"/> Left side <input type="checkbox"/></p> <p>with existing support line (placed along a brake line)</p> <p><input type="checkbox"/> A9116200000 € 98.20</p> <p>without existing support line (an additional blind cable is required)</p> <p><input type="checkbox"/> A9116200001 € 139.20</p>	
---	---	--

Services

Delivery and instruction of the handbike / traction device by the specialist dealer incl. assembly of the frame adapter

on a PRO ACTIV wheelchair

A9990000008

prod. no: _____

on a third-party product

A9990000005

Manufacturer: _____

Model: _____

Frame tube ø ____ mm

Notes and remarks:

Location, Date

The signature field is intended for the consultant and insured person. This confirms that the intended equipment is intended for the insured person.

Signature End customer/dealer

- 7 - 2024-06-01 / The prices are plus statutory value-added tax and shipping costs

= Included in the basic price = Position selected

Order sheet NJ1 el. compactbike Rev. 03 UK